



Microinjection System



The basic Zebrafish Injection system includes: **PZMIII** microscope with lighted base, **PV820** PicoPump, **M3301** micromanipulator, **5430-ALL** PicoNozzle Kit, **M10** magnetic stand, **5052** steel base plate.

Zebrafish (*Danio Rerio*) are rapidly gaining in popularity as bio-medical research subjects because of the ability to generate high resolution, *in vivo* images of the embryos. Zebrafish are easy to maintain and produce a large number of offspring. Additionally, the embryos have a nearly transparent skin, making their development easily visible. These fish are used for a variety of disciplines, including neuroscience, genetics and aging studies.

Serving scientist for over 45 years, WPI offers a variety of instruments for microinjection including pumps, pipettors, microscopes and more. One of our most popular pumps for microinjection is the PV820 Pneumatic PicoPump.

The **PV820** and **PV830**, Pneumatic PicoPumps, were designed to simplify intracellular injection. You get repeatable microinjection in volumes ranging from picoliters to nanoliters. **PV820** offers eject and hold pressure. The hold pressure prevents backfilling of the pipette by capillary action. In addition, the **PV830** also has vacuum pressure which allows you to securely hold a cell with one pipette while you inject it with another. The volume injected is controlled by the inside diameter of the glass tip, the pressure and the time.

Recently, WPI introduced its customizable Microinjection System with everything you need to get started. The basic system is shown here. On the next page you will find many options and accessories you may use to customize your system.



The system depicted on the cover and at left includes components often favored by researchers, indicated with * in the list below: **PV820** Pneumatic PicoPump, **PUL-1000** Micropipette Puller, **M10** magnetic stand on a **5052** steel base plate (not shown), **M3301L** micromanipulator, **5430-ALL** PicoNozzle Kit with a μ Tip, **FD35-100** Fluorodish, **PZMTIII** microscope with lighted base and articulating mirror (optional), **E210** micropipette storage jar, **504156** Regine tweezers #5, **14003** Vannas spring scissors, glass capillaries, **77020** glass tweezers.

Whatever your needs, WPI offers a range of equipment to fill your requirements

Options for Customizing Your System

INJECTOR

- * **PV820** Pneumatic PicoPump with Hold Pressure
- **PV830** Pneumatic PicoPump with Hold Pressure and Vacuum
- **Nanoliter2010**
- **UMPIII** UltraMicroPump



Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's **PicoPumps** use carefully regulated air pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters.



Microprocessor-controlled **Nanoliter 2010** uses direct piston displacement.



The versatile **UMPIII** injector uses microsyringes to deliver picoliter volumes.

MICROSCOPE

- * **PZMIII-MI** Microscope with lighted base and articulating mirror
- **PZMIII** Precision Stereo Zoom Microscope

PZMIII-MI LED Lighted Microscope Base with articulating mirror and variable light intensity. Dual reflection lens/mirror system provides transmitted brightfield/psuedo-darkfield illumination.



PZMIII Precision Stereo Zoom Microscope, available in binocular and trinocular (shown).

PULLERS

- * **PUL1000** Microprocessor-controlled 4-Step micropipette puller
- **PMP-102** Programmable Multipipette Puller

PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, microperfusion or microinjection. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.



For more information, see www.wpiinc.com/microinjection

MANIPULATOR

- * **M3301** Manual Micromanipulator
- **KITE** Manual Micromanipulator
- **DC3001** Motorized Micromanipulator
- **SN-PCZ-50** Miniature Piezo Micromanipulator with controller
- **M4C** Microscope Stage Adapter



Weighing just 550 grams, the **M3301** is a well-built micromanipulator that outsells all others worldwide for high precision experiments where magnification is in the range of up to 250 \times .



M4C Microscope Stage Adapter

LIGHTS

- * **Z-LITE-186** Fiber optic illuminator with (500186) Bifurcated Light Guides
- **ROF-UV** LED ringlight for UV
- **LED-Lite** Modular LED Light Source with Exchangeable LEDs
- **504134** LED Ring light



The **504134** Ring Light Adapter has 72 LEDs. Max opening: 61mm.



Z-LITE-186


ACCESSORIES

- **Nanofil** Microliter syringes
- **MicroFil** for backfilling glass needles
- Glass capillaries
- Pipettors
- **MicroTip** pre-pulled pipettes
- **E2XX** Micropipette Storage Jar
- **801566/801963** Vacuum Pump for use with the PV830
- **Fluorodish** Optical glass bottom dishes
- **M10** or **M-3** Manipulator base
- **Z-MOLDS** Microinjection and Transplantation Molds
- Many surgical instruments

Use the **E2XX** jars to store up to 30 micropipettes, filled or unfilled, up to three inches in length.



MicroFil is a flexible and reusable glass needle as small as 36 gauge for backfilling micropipettes.



Z-MOLDS Microinjection and Transplantation Molds (4 per kit) are designed for zebrafish research. The molds are turned up-side down and placed in liquid agarose gel and are easily removed once it has solidified. Pipette the embryos into the grooves. The embryos self-align.

ELECTROPORATOR

- **CUY21-EDIT-11** Electroporator



CUY21-EDIT-II Electroporator is a combination square wave and exponential wave form generator with single and dual pulse modes. It has single and dual pulse modes for square and decaying pulses. You can set pulses emitting constant current (mA) as well as voltage (V) and generate multiple waveforms. It offers built-in resistance measurement and a measured output recorder.

Microinjection System

ZEBRAFISH REFERENCES

Referenced in 1567 peer reviewed publications

Video: Konantz, J., Antos, C. L. (2014). *Reverse Genetic Morpholino Approach Using Cardiac Ventricular Injection to Transfect Multiple Difficult-to-target Tissues in the Zebrafish Larva* : J. Vis. Exp. (88), e51595, doi:10.3791/51595. <http://www.jove.com/video/51595/reverse-genetic-morpholino-approach-using-cardiac-ventricular>.

Video: Kizil, C., Iltzsch, A., Kaslin, J., Brand, M. (2013). *Micromanipulation of Gene Expression in the Adult Zebrafish Brain Using Cerebroventricular Microinjection of Morpholino Oligonucleotides* : J. Vis. Exp. (75), e50415, doi:10.3791/50415. <http://www.jove.com/video/50415/micromanipulation-gene-expression-adult-zebrafish-brain-using>

Video: Holmes, K. E., Wyatt, M. J., Shen, Y., Thompson, D. A., Barald, K. F. (2011). *Direct Delivery of MIF Morpholinos Into the Zebrafish Otolocyst by Injection and Electroporation Affects Inner Ear Development* : J. Vis. Exp. (47), e2466, doi:10.3791/2466. <http://www.jove.com/video/2466/direct-delivery-mif-morpholinos-into-zebrafish-otocyst-injection>

Video: Cianciolo Cosentino, C., Roman, B. L., Drummond, I. A., Hukriede, N. A. (2010). *Intravenous Microinjections of Zebrafish Larvae to Study Acute Kidney Injury* : J. Vis. Exp. (42), e2079, doi:10.3791/2079. <http://www.jove.com/video/2079/intravenous-microinjections-zebrafish-larvae-to-study-acute-kidney>

Video: Russek-Blum, N., Nabel-Rosen, H., Levkowitz, G. (2010). *Two-Photon-Based Photoactivation in Live Zebrafish Embryos*. J. Vis. Exp. (46), e1902, doi:10.3791/1902.

Video: Kague, E., Weber, C., Fisher, S. (2010). *Mosaic Zebrafish Transgenesis for Evaluating Enhancer Sequences* : J. Vis. Exp. (41), e1722, doi:10.3791/1722. <http://www.jove.com/video/1902/two-photon-based-photoactivation-in-live-zebrafish-embryos>

Yabing Li, Zhihua Han, Xinmei Zheng, Zhiyuan Ma, Hongling Liu, John P. Giesy, Yuwei Xie, Hongxia Yu. (2014). *Comparison of waterborne and in ovo nanoinjection exposures to assess effects of PFOS on zebrafish embryos* : Environ Sci Pollut Res. DOI 10.1007/s11356-014-3527-y. <http://www.usask.ca/toxicology/jgiesy/pdf/publications/JA-799-temp.pdf>

Katarina Vajn, Denis Suler, Jeffery A. Plunkett, Martin Oudega. (2014). *Temporal Profile of Endogenous Anatomical Repair and Functional*

Recovery following Spinal Cord Injury in Adult Zebrafish : Published: August 26, 2014 DOI: 10.1371/journal.pone.0105857. <http://journals.plos.org/plosone/article?id=10.1371/journal.pone.0105857>

Henson HE, Parupalli C, Ju B, Taylor MR. (2014). *Functional and genetic analysis of cho-roid plexus development in zebrafish* : Frontiers in Neuroscience, v.8, 364. doi:10.3389/fnins.2014.00364. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC4226144>.

Wagley S, Hemsley C, Thomas R, et al. (2014). *The twin arginine translocation system is essential for aerobic growth and full virulence of Burkholderia thailandensis* : J Bacteriol. 2014;196(2):407-416. doi:10.1128/JB.01046-13. <http://jlb.asm.org/content/196/2/407.full>

Perathoner S, Daane JM, Henrion U, et al. (2014). *Bioelectric signaling regulates size in zebrafish fins* : PLoS Genet. 2014;10(1):e1004080. doi:10.1371/journal.pgen.1004080. <http://journals.plos.org/plosgenetics/article?id=10.1371/journal.pgen.1004080>

Jo DH, Son D, Na Y, et al. (2013). *Orthotopic transplantation of retinoblastoma cells into vitreous cavity of zebrafish for screening of anticancer drugs*. Mol Cancer. 2013;12:71. doi:10.1186/1476-4598-12-71. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3707771/>

Jo DH, Son D, Na Y, et al. (2013). *Orthotopic transplantation of retinoblastoma cells into vitreous cavity of zebrafish for screening of anticancer drugs*. Mol Cancer. 2013;12(1):71. doi:10.1186/1476-4598-12-71. <http://www.molecular-cancer.com/content/12/1/71>

Gramma A, Engert F. (2012). *Direction selectivity in the larval zebrafish tectum is mediated by asymmetric inhibition* : Front Neural Circuits. 2012;6:59. doi:10.3389/fncir.2012.00059. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3432856/>

Huber-Reggi SP, Chen C-C, Grimm L, Straumann D, Neuhauss SCF, Huang MY-Y. (2012). *Severity of infantile nystagmus syndrome-like ocular motor phenotype is linked to the extent of the underlying optic nerve projection defect in zebrafish bella-donna mutant* : J Neurosci. 2012;32(50):18079-18086. doi:10.1523/JNEUROSCI.4378-12.2012. <http://www.jneurosci.org/content/32/50/18079.full>

Nakaya N, Sultana A, Lee H-S, Tomarev SI. (Oct. 26, 2012) *Olfactomedin 1 interacts with the Nogo A receptor complex to regulate axon growth*. J Biol Chem. 2012;287(44):37171-37184. doi:10.1074/jbc.M112.389916. <http://www.jbc.org/content/287/44/37171.full>

Seo J, Yun C-O, Kwon O-J, et al. (2012). *A proteoliposome containing apolipoprotein A-I mutant (V156K) enhances rapid tumor regression activity of human origin oncolytic adenovirus in tumor-bearing zebrafish and mice* : Mol Cells. 2012;34(2):143-148. doi:10.1007/s10059-012-2291-4. <http://dx.doi.org/10.1007/s10059-012-2291-4>

Cho K-H. (2011) *Enhanced delivery of rapamycin by V156K-apoA-I high-density lipoprotein inhibits cellular proatherogenic effects and senescence and promotes tissue regeneration* : J Gerontol A Biol Sci Med Sci. 2011;66(12):1274-1285. doi:10.1093/gerona/glr169. <http://biomedgerontology.oxfordjournals.org/content/66A/12/1274.full>

Kizil C, Brand M. (2011). *Cerebroventricular microinjection (CVMi) into adult zebrafish brain is an efficient misexpression method for forebrain ventricular cells* : PLoS One. 2011;6(11):e27395. doi:10.1371/journal.pone.0027395. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3208640/>

Russek-Blum N, Nabel-Rosen H, Levkowitz G. (2010). *Two-photon-based photoactivation in live zebrafish embryos* : J Vis Exp. 2010;(46). doi:10.3791/1902. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3159643/>

Wang Y, Nakayama M, Pitulescu ME, et al. (2010). *Ephrin-B2 controls VEGF-induced angiogenesis and lymphangiogenesis* : Nature. Vol 465, 27 May 2010. doi:10.1038/nature09002. <http://www.med.upenn.edu/timm/documents/nature09002.pdf>

Zhu Y, Song D, Tran N-T, Nguyen N. (2007). *The effects of the members of growth hormone family knockdown in zebrafish development* : General and Comparative Endocrinology 150 (2007) 395-404. <http://core.ecu.edu/biol/zhuy/publications/2007Zhu.pdf>

Kopeika J, Zhang T, Rawson D. (2006). *Zebrafish embryos (Danio rerio) using microinjection* : Cryo Letters. 27(5):319-328. <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=1855166&ol=pmcentrez&rendertype=abstract>.



WORLD PRECISION INSTRUMENTS, INC.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA

Tel: 941-371-1003 • **Fax:** 941-377-5428 • **E-mail:** wpi@wpiinc.com • **Internet:** www.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com

Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188870 • E-mail: wpide@wpi-europe.com

China & Hong Kong: Rm 29a, No8 Donfang Rd., Pudong District, Shanghai 200120 PRC • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

Brazil: Av. Conselheiro Nébias, 756 sala 2611, Santos-CEP: 11045-002, São Paulo Brazil • Tel: (013) 406-29703 • E-mail: info@brazil.wpiinc.com



Microinjection System



The basic Zebrafish Injection system includes:
PZMIII microscope with lighted base, **PV820** PicoPump, **M3301** micromanipulator, **5430-ALL** PicoNozzle Kit, **M10** magnetic stand, **5052** steel base plate.

Zebrafish (*Danio Rerio*) are rapidly gaining in popularity as bio-medical research subjects because of the ability to generate high resolution, *in vivo* images of the embryos. Zebrafish are easy to maintain and produce a large number of offspring. Additionally, the embryos have a nearly transparent skin, making their development easily visible. These fish are used for a variety of disciplines, including neuroscience, genetics and aging studies.

Serving scientist for over 45 years, WPI offers a variety of instruments for microinjection including pumps, pipettors, microscopes and more. One of our most popular pumps for microinjection is the PV820 Pneumatic PicoPump.

The **PV820** and **PV830**, Pneumatic PicoPumps, were designed to simplify intracellular injection. You get repeatable microinjection in volumes ranging from picoliters to nanoliters. **PV820** offers eject and hold pressure. The hold pressure prevents backfilling of the pipette by capillary action. In addition, the **PV830** also has vacuum pressure which allows you to securely hold a cell with one pipette while you inject it with another. The volume injected is controlled by the inside diameter of the glass tip, the pressure and the time.

Recently, WPI introduced its customizable Microinjection System with everything you need to get started. The basic system is shown here. On the next page you will find many options and accessories you may use to customize your system.



The system depicted on the cover and at left includes components often favored by researchers, indicated with * in the list below: **PV820** Pneumatic PicoPump, **PUL-1000** Micropipette Puller, **M10** magnetic stand on a **5052** steel base plate (not shown), **M3301L** micromanipulator, **5430-ALL** PicoNozzle Kit with a μ Tip, **FD35-100** Fluorodish, **PZMTIII** microscope with lighted base and articulating mirror (optional), **E210** micropipette storage jar, **504156** Regine tweezers #5, **14003** Vannas spring scissors, glass capillaries, **77020** glass tweezers.

Whatever your needs, WPI offers a range of equipment to fill your requirements

Options for Customizing Your System

INJECTOR

- * **PV820** Pneumatic PicoPump with Hold Pressure
- **PV830** Pneumatic PicoPump with Hold Pressure and Vacuum
- **Nanoliter2010**
- **UMPIII** UltraMicroPump



Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's **PicoPumps** use carefully regulated air pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters.



Microprocessor-controlled **Nanoliter 2010** uses direct piston displacement.



The versatile **UMPIII** injector uses microsyringes to deliver picoliter volumes.

MICROSCOPE

- * **PZMIII-MI** Microscope with lighted base and articulating mirror
- **PZMIII** Precision Stereo Zoom Microscope

PZMIII-MI LED Lighted Microscope Base with articulating mirror and variable light intensity. Dual reflection lens/mirror system provides transmitted brightfield/psuedo-darkfield illumination.



PZMIII Precision Stereo Zoom Microscope, available in binocular and trinocular (shown).

PULLERS

- * **PUL1000** Microprocessor-controlled 4-Step micropipette puller
- **PMP-102** Programmable Multipipette Puller

PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, microperfusion or microinjection. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.




For more information, see www.wpiinc.com/microinjection

MANIPULATOR

- * **M3301** Manual Micromanipulator
- **KITE** Manual Micromanipulator
- **DC3001** Motorized Micromanipulator
- **SN-PCZ-50** Miniature Piezo Micromanipulator with controller
- **M4C** Microscope Stage Adapter



Weighing just 550 grams, the **M3301** is a well-built micromanipulator that outsells all others worldwide for high precision experiments where magnification is in the range of up to 250 \times .



M4C Microscope Stage Adapter

Z-LITE-186

LIGHTS

- * **Z-LITE-186** Fiber optic illuminator with (500186) Bifurcated Light Guides
- **ROF-UV** LED ringlight for UV
- **LED-Lite** Modular LED Light Source with Exchangeable LEDs
- **504134** LED Ring light



The **504134** Ring Light Adapter has 72 LEDs. Max opening: 61mm.


ACCESSORIES

- **Nanofil** Microliter syringes
- **MicroFil** for backfilling glass needles
- Glass capillaries
- Pipettors
- **MicroTip** pre-pulled pipettes
- **E2XX** Micropipette Storage Jar
- **801566/801963** Vacuum Pump for use with the PV830
- **Fluorodish** Optical glass bottom dishes
- **M10** or **M-3** Manipulator base
- **Z-MOLDS** Microinjection and Transplantation Molds
- Many surgical instruments

Use the **E2XX** jars to store up to 30 micropipettes, filled or unfilled, up to three inches in length.



MicroFil is a flexible and reusable glass needle as small as 36 gauge for backfilling micropipettes.



Z-MOLDS Microinjection and Transplantation Molds (4 per kit) are designed for zebrafish research. The molds are turned up-side down and placed in liquid agarose gel and are easily removed once it has solidified. Pipette the embryos into the grooves. The embryos self-align.

ELECTROPORATOR

- **CUY21-EDIT-11** Electroporator



CUY21-EDIT-II Electroporator is a combination square wave and exponential wave form generator with single and dual pulse modes. It has single and dual pulse modes for square and decaying pulses. You can set pulses emitting constant current (mA) as well as voltage (V) and generate multiple waveforms. It offers built-in resistance measurement and a measured output recorder.

Microinjection System

ZEBRAFISH REFERENCES

Referenced in 1567 peer reviewed publications

Video: Konantz, J., Antos, C. L. (2014). *Reverse Genetic Morpholino Approach Using Cardiac Ventricular Injection to Transfect Multiple Difficult-to-target Tissues in the Zebrafish Larva* : J. Vis. Exp. (88), e51595. doi:10.3791/51595. <http://www.jove.com/video/51595/reverse-genetic-morpholino-approach-using-cardiac-ventricular>.

Video: Kizil, C., Iltzsche, A., Kaslin, J., Brand, M. (2013). *Micromanipulation of Gene Expression in the Adult Zebrafish Brain Using Cerebroventricular Microinjection of Morpholino Oligonucleotides* : J. Vis. Exp. (75), e50415, doi:10.3791/50415. <http://www.jove.com/video/50415/micromanipulation-gene-expression-adult-zebrafish-brain-using>

Video: Holmes, K. E., Wyatt, M. J., Shen, Y., Thompson, D. A., Barald, K. F. (2011). *Direct Delivery of MIF Morpholinos Into the Zebrafish Otolarynx by Injection and Electroporation Affects Inner Ear Development* : J. Vis. Exp. (47), e2466, doi:10.3791/2466. <http://www.jove.com/video/2466/direct-delivery-mif-morpholinos-into-zebrafish-otocyst-injection>

Video: Cianciolo Cosentino, C., Roman, B. L., Drummond, I. A., Hukriede, N. A. (2010). *Intravenous Microinjections of Zebrafish Larvae to Study Acute Kidney Injury* : J. Vis. Exp. (42), e2079, doi:10.3791/2079. <http://www.jove.com/video/2079/intravenous-microinjections-zebrafish-larvae-to-study-acute-kidney>

Video: Russek-Blum, N., Nabel-Rosen, H., Levkowitz, G. (2010). *Two-Photon-Based Photoactivation in Live Zebrafish Embryos*. J. Vis. Exp. (46), e1902, doi:10.3791/1902.

Video: Kague, E., Weber, C., Fisher, S. (2010). *Mosaic Zebrafish Transgenesis for Evaluating Enhancer Sequences* : J. Vis. Exp. (41), e1722, doi:10.3791/1722. <http://www.jove.com/video/1902/two-photon-based-photoactivation-in-live-zebrafish-embryos>

Yabing Li, Zhihua Han, Xinmei Zheng, Zhiyuan Ma, Hongling Liu, John P. Giesy, Yuwei Xie, Hongxia Yu. (2014). *Comparison of waterborne and in ovo nanoinjection exposures to assess effects of PFOS on zebrafish embryos* : Environ Sci Pollut Res. DOI 10.1007/s11356-014-3527-y. <http://www.usask.ca/toxicology/jgiesy/pdf/publications/JA-799-temp.pdf>

Katarina Vajn, Denis Suler, Jeffery A. Plunkett, Martin Oudega. (2014). *Temporal Profile of Endogenous Anatomical Repair and Functional Recovery following Spinal Cord Injury in Adult*

Zebrafish : Published: August 26, 2014 DOI: 10.1371/journal.pone.0105857. <http://journals.plos.org/plosone/article?id=10.1371/journal.pone.0105857>

Henson HE, Parupalli C, Ju B, Taylor MR. (2014). *Functional and genetic analysis of chorioid plexus development in zebrafish* : Frontiers in Neuroscience, v.8, 364. doi:10.3389/fnins.2014.00364. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC4226144>.

Wagley S, Hemsley C, Thomas R, et al. (2014). *The twin arginine translocation system is essential for aerobic growth and full virulence of Burkholderia thailandensis* : J Bacteriol. 2014;196(2):407-416. doi:10.1128/JB.01046-13. <http://jlb.asm.org/content/196/2/407.full>

Perathoner S, Daane JM, Henrion U, et al. (2014). *Bioelectric signaling regulates size in zebrafish fins* : PLoS Genet. 2014;10(1):e1004080. doi:10.1371/journal.pgen.1004080. <http://journals.plos.org/plosgenetics/article?id=10.1371/journal.pgen.1004080>

Jo DH, Son D, Na Y, et al. (2013). *Orthotopic transplantation of retinoblastoma cells into vitreous cavity of zebrafish for screening of anticancer drugs*. Mol Cancer. 2013;12:71. doi:10.1186/1476-4598-12-71. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3707771/>

Jo DH, Son D, Na Y, et al. (2013). *Orthotopic transplantation of retinoblastoma cells into vitreous cavity of zebrafish for screening of anticancer drugs*. Mol Cancer. 2013;12(1):71. doi:10.1186/1476-4598-12-71. <http://www.molecular-cancer.com/content/12/1/71>

Gramma A, Engert F. (2012). *Direction selectivity in the larval zebrafish tectum is mediated by asymmetric inhibition* : Front Neural Circuits. 2012;6:59. doi:10.3389/fncir.2012.00059. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3432856/>

Huber-Reggi SP, Chen C-C, Grimm L, Straumann D, Neuhaus SCF, Huang MY-Y. (2012). *Severity of infantile nystagmus syndrome-like ocular motor phenotype is linked to the extent of the underlying optic nerve projection defect in zebrafish belladonna mutant* : J Neurosci. 2012;32(50):18079-18086. doi:10.1523/JNEUROSCI.4378-12.2012. <http://www.jneurosci.org/content/32/50/18079.full>

Nakaya N, Sultana A, Lee H-S, Tomarev SI. (Oct. 26, 2012) *Olfactomedin 1 interacts with*

the Nogo A receptor complex to regulate axon growth. J Biol Chem. 2012;287(44):37171-37184. doi:10.1074/jbc.M112.389916. <http://www.jbc.org/content/287/44/37171.full>

Seo J, Yun C-O, Kwon O-J, et al. (2012). *A proteoliposome containing apolipoprotein A-I mutant (V156K) enhances rapid tumor regression activity of human origin oncolytic adenovirus in tumor-bearing zebrafish and mice* : Mol Cells. 2012;34(2):143-148. doi:10.1007/s10059-012-2291-4. <http://dx.doi.org/10.1007/s10059-012-2291-4>

Cho K-H. (2011) *Enhanced delivery of rapamycin by V156K-apoA-I high-density lipoprotein inhibits cellular proatherogenic effects and senescence and promotes tissue regeneration* : J Gerontol A Biol Sci Med Sci. 2011;66(12):1274-1285. doi:10.1093/geronol/glr169. <http://biomedgerontology.oxfordjournals.org/content/66A/12/1274.full>

Kizil C, Brand M. (2011). *Cerebroventricular microinjection (CVMI) into adult zebrafish brain is an efficient misexpression method for forebrain ventricular cells* : PLoS One. 2011;6(11):e27395. doi:10.1371/journal.pone.0027395. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3208640/>

Russek-Blum N, Nabel-Rosen H, Levkowitz G. (2010). *Two-photon-based photoactivation in live zebrafish embryos* : J Vis Exp. 2010;(46). doi:10.3791/1902. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3159643/>

Wang Y, Nakayama M, Pitulescu ME, et al. (2010). *Ephrin-B2 controls VEGF-induced angiogenesis and lymphangiogenesis* : Nature. Vol 465, 27 May 2010. doi:10.1038/nature09002. <http://www.med.upenn.edu/timm/documents/nature09002.pdf>

Zhu Y, Song D, Tran N-T, Nguyen N. (2007). *The effects of the members of growth hormone family knockdown in zebrafish development* : General and Comparative Endocrinology 150 (2007) 395-404. <http://core.ecu.edu/biol/zhu/publications/2007Zhu.pdf>

Kopeika J, Zhang T, Rawson D. (2006). *Zebrafish embryos (Danio rerio) using microinjection* : Cryo Letters. 27(5):319-328. <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=1855166&tool=pmcentrez&rendertype=abstract>.



WORLD PRECISION INSTRUMENTS, INC.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA
Tel: 941-371-1003 • **Fax:** 941-377-5428 • **E-mail:** wpi@wpiinc.com • **Internet:** www.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com
Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188870 • E-mail: wpide@wpi-europe.com

China & Hong Kong: Rm 29a, No8 Dongfang Rd., Pudong District, Shanghai 200120 PRC • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com

Brazil: Av. Conselheiro Nébias, 756 sala 2611, Santos-CEP: 11045-002, São Paulo Brazil • Tel: (013) 406-29703 • E-mail: info@brazil.wpiinc.com