





Let us Be Your Partner in Discovery

We are at the dawn of an exciting new day of science—improving the lives of all humanity with breakthroughs like pheno-genomic testing and therapy. Daily, researchers like you are discovering new models, pathways and techniques which impact our health and our environment. These discoveries result in life-changing—and in many cases life-prolonging—therapies and revolutionary protocols. We are also working with researchers in monitoring environmental influences in our oceans and streams. New surgical instruments are being created to assist in specific surgeries and diagnoses. This year, for example, we added a new line of tools specifically for work under a microscope. They are non-reflective and are also modified with special tips to enable surgeries in arteries and veins.

At WPI, our solutions enable pioneers on the cusp of discovery to impact the quality of life for all of us for generations to come. WPI provides solutions for a wide range of disciplines, including electrophysiology, molecular biology, cellular biology, surgical procedures and physiology, to advance the life science knowledge base. As researchers continue the quest for answers to life's mysteries, we provide the instruments to make laboratory research effective, reliable and reproducible. We aim to be your partner in the scientific adventure and instill confidence in your journey of discovery. Your research methods and ideas, combined with our scientists and engineers, help us develop novel solutions for your work.

Our customers' journey of discovery, and the promise it brings of life-changing impact, inspires our team. Let us be your research partner with quality, cost-effective instruments for wherever science takes you.



Table of Contents

New Product Summ	nary			2
Cell & Tissue				3
Epithelial Physiology Muscle Physiology Biofluorometer	Micromanipulators Mounting Devices Magnetic Holding Devices	Live Cell Imaging Vibration Isolation		
Microinjection				39
Microinjection Systems Microinjection Accessories	Pinpoint Injector Picoliter Injectors			
Fluid Handling				49
Microperfusion Peristaltic Pumps UltraMicroPumps	Syringes UMP3 Application Nanoliter Injectors	Syringe Pumps Microfluidics Pipetters	Pipetter Tips Vacuum Pumps Luer Fittings	
Electrophysiology				73
Amplifiers Stimulators Stimulus Isolators	Electrophysiology Accessories Microelectrodes Wire, Cables & Connectors	Pipette Holders Capillary Glass Specialty Wire & Half Cells	Micropipettes Pullers, Microforges Micropressure Systems	
Animal Physiology				117
Stereotaxic Anesthesia Blood Pressure Monitoring	Thermometers Temperature Control Neuroscience Cannulas	Analgesia Adhesives		
Biosensing				143
Biosensors Free Radical Analyzers Oxygen Detection Systems	Data Acquisition pH Detection Biosensing Accessories			
Optical Measureme	ent Systems			161
Qualified Optic Cables Optical Detection Flow Cells	Sample Holders Dipping Probes Cuvettes	Spectrometers Light Sources Biofluorometer	Microscopy & Cameras Illumination	
Microdissection Ins	truments			192
Forceps Scissors Needle Holders	Black Instruments Punches Clamps, Scalpels	Rodent Accessories Microdissection Kits Surgical Equipment	Care and Cleaning	
Lab Supplies				212
Adhesives Glass Capillaries Micropipettes	Luer Fittings BNC Cables & Connectors Wire	Bases & Mounts Calipers & Micrometers Pipetters & Tips	Illumination Culture Dishes	
Index				226

Products at WPI




The new MICRO-ePORE™ cell penetrator is a simple and versatile pin-point controlled electroporation system that can be used for efficient microinjection of reagents into oocytes and pre-implantation stage mammalian embryos. Patent pending Flutter Electrode Technology

Patent pending Flutter Electrode Technology assists in small, clean, precise membrane penetration without tearing or damaging the membrane. See page 42.

QUALIFIED OPTICAL FIBERS

The use of qualified optical fibers ensures better fibers, better assemblies, better probes. WPI excels at making optical fibers with high batch-to-batch performance consistency. WPI UV fiber has industry leading degradation resistance to deep UV light. WPI fiber composition and manufacturing is unique so that we can offer you qualified and guaranteed



FLUORESCENCE PROBE

WPI offers fluorescence probes to match your application, both single and double emission. These probes may be used for the detection of the transient response of free ion concentrations, like calcium, potassium, sodium or magnesium. Fluorescence probes can also be designed to detect pH or membrane potential. Auto-fluorescence like the detection of ATPase activity via NADH or FAD is another application. For details, talk with a WPI application specialist. See page 178.



OPTICAL OXYGEN SENSOR



Now you can detect oxygen in small sample in living tissues using the phase shift between a reference signal and a measured signal. **BioOxy** is a new and innovative technology for measuring oxygen in gaseous and aqueous phase. **BioOxy** is an optical oxygen sensor with important advantages over using common Clark type electrodes. These sensors have a rapid response time with the $\mathsf{t}_{90} \leq 3$ seconds (in gas phase). No oxygen is consumed making measurements. These sensors are perfect for bioprocess control and making oxygen measurement in small samples. See page 152.

TEER MEASUREMENT CHAMBERS



ENDOHM-6G

ENDOHM-24GSNAP

Our new EndOhm chambers are now made of glass, and the new insert holder with 120° tri-supports can accommodate three leg inserts! The crystal clear glass chamber allows visualization of apical electrode positioning. See page 5.

STERILIZATION TRAYS AND BASKETS

With top quality structural integrity, these trays are ideal for the handling and storage of all standard microdissection and surgical instruments. The unique "grid" design of the base makes it easy to install the silicone finger mats required to protect the delicate instruments. They have burr-free edges, and we offer a large selection to fit your needs. See www.wpiinc. com/baskets. See page 211.



OPHTHALMIC KITS AND INSTRUMENTS

We offer kits for specific research and veterinary procedures. And we have over 100 ophthalmic instruments each sold separately. See www.wpiinc. com/eyes. See page 198.

CERAMIC COATED INSTRUMENTS

Our black instruments are coated with titanium nitride (TiN), an extremely hard ceramic material. The TiN coating hardens and protects the cutting edge. Ceramic coated anti-reflective instruments are perfect for microscopy and



microsurgical applications. Coating surgical instruments with a black ceramic adds a thin layer to the metal instrument, making the instrument harder and giving you greater precision. This anti-glare surface minimizes reflections off the surface of your instruments. The incredibly smooth coating improves the instruments' resistance to corrosion and minimizes friction. The ceramic coating is virtually impenetrable, because the raw material is bonded to the instrument both physically and chemically. These instruments are much more resilient to the pressure of daily use and chemical processing. Coated instruments last considerably longer. See page 203.

STERI-LITE MICROBEAD STERILIZER

Ideal for the sterilization of small research instruments, the chamber is filled with glass beads and heated up to 300°C. Eliminate bacteria, spores and other microorganisms. Chamber temperature is displayed on the large LED control panel and may be adjusted with the control knob. See page 210.

ST5191

Cell & Tissue



Comprehensive Products for Cell & Tissue Applications

WPI's cell and tissue products cover a broad range of applications. Our products are tested and optimized to provide you with instruments that are reliable in supporting your research requirements. We offer a range of products for TEER measurement, live cell imaging, microscopy studies and muscle physiology applications. WPI's products have been cited in 100s of reference papers and in many cases are the only products available of their kind. WPI's **EVOM2™** is the original instrument designed specifically to perform non-destructive TEER measurement on epithelial monolayer cell cultures.

Epithelial Volt/Ohm (TEER) Meter

Non-destructively test for epithelial monolayer confluence in 2D cell cultures

Features

- Measures trans-epithelial electrical resistance or trans-epithelilal voltage
- Compatible with 12 and 24 well culture plate systems out of the box Includes industry standard STX2 hand held "chopstick" electrodes
- · Analog output for recording resistance or voltage measurements
- Auto ranging from 0-10 KΩ
- · Battery powered

Benefits

- You can verify performance and calibrate the meter for TEER function using provided test resistor
- Battery powered meter is portable

A variety of accessory electrodes are available for measuring TEER in 6and 96-well fixed (HTS) and removable well culture systems (See STX100 series (page 6) and Endohm (page 5) electrodes

Applications

• TEER and trans-epithelial voltage measurements in 2D cell cultures

The **EVOM™** was the first instrument designed specifically to perform routine Trans Epithelial Electrical Resistance (TEER) measurement in tissue culture research. **EVOM2™** is the next generation, redesigned for ease of use. The EVOM2™ not only qualitatively measures cell monolayer health, but also quantitatively measures cellular confluence. The unique electronic circuit of the EVOM2™ and the included STX2 electrode detect the confluence of the cellular monolayer. When combined with WPI's Endohm chamber, the **EVOM2™** can also be used to perform more accurate quantitative measurements or lower resistance measurements like transendothelial electrical resistance measurements.

Isolated battery power for 10 hours of use

The isolated power source of the **EVOM2™** was specifically designed to avoid adverse effects on tissue and the formation of electrode metal $% \left(1\right) =\left(1\right) \left(1\right$ deposits, even when it is plugged into a standard wall outlet. Now, the **EVOM2™** is always on when you need it. In addition, its rechargeable battery allows up to 10 hours of mobile use.

Accurate reading every time

The four-and-a-half digit readout provides a range of 1-10,000 Ω . The included test electrode lets you calibrate the resistance measurements for an accurate reading every time, and the voltage meter never needs calibration. An analog BNC output is standard with the EVOM2™ providing an output port for recording data or remote display of the EVOM2™ output.

Electrode pair to measure voltage, pass current

EVOM2™ comes complete with the popular STX2 "chopstick" electrodes, 4 mm wide and 1 mm thick. Each stick of the electrode pair contains a silver/silver-chloride pellet for measuring voltage and a silver electrode for passing current. The small size of each electrode is designed to facilitate placement of the electrodes into a variety of standard cell culture wells.

References

Di, S., Gujie, M., & Thomas, W. (2016). Magnetic ferri-liposomes for triggered drug release across the blood- brain barrier. Frontiers in Bioengineering and

Biotechnology, 4. http://doi.org/10.3389/conf.FBIOE.2016.01.00061

Paradis, A., Leblanc, D., & Dumais, N. (2016). Optimization of an in vitro human blood-brain barrier model: Application to blood monocyte transmigration assays. MethodsX, 3, 25-34. http://doi.org/10.1016/j.mex.2015.11.009





Meenach, S. A., Tsoras, A. N., McGarry, R. C., Mansour, H. M., & Anderson, K. W. (2016). Development of three-dimensional lung multicellular spheroids in air- and liquid-interface culture for the evaluation of anticancer therapeutics. International Journal of Oncology, 48(4), 1701-9. http://doi.org/10.3892/ iio.2016.3376

Turdalieva, A., Solandt, J., Shambetova, N., Xu, H., Blom, H., Brismar, H., ... Fu, Y. (2016). Bioelectric and Morphological Response of Liquid-Covered Human Airway Epithelial Calu-3 Cell Monolayer to Periodic Deposition of Colloidal 3-Mercaptopropionic-Acid Coated CdSe-CdS/ZnS Core-Multishell Quantum Dots. PloS One, 11(2), e0149915. http://doi.org/10.1371/journal.pone.0149915

Hirano, M., & Hirano, K. (2016). Myosin di-phosphorylation and peripheral actin bundle formation as initial events during endothelial barrier disruption. Scientific Reports, 6(1), 20989. http://doi.org/10.1038/srep20989

Marvin, S. A., et. al. (2016). Type I Interferon Response Limits Astrovirus Replication and Protects against Increased Barrier Permeability In Vitro and In Vivo. Journal of Virology, 90(4), 1988-1996. http://doi.org/10.1128/JVI.02367-15

Jacobs, M. E., et. al. (2016). SGK1 regulation by miR-466g in cortical collecting duct cells. American Journal of Physiology-Renal Physiology, 310(11), F1251–F1257. http://doi.org/10.1152/ajprenal.00024.2016

Yu, C., et. al. (2015). The Effects of Glucagon-like Peptide-2 on the Tight Junction and Barrier Function in IPEC-J2 Cells through Phosphatidylinositol 3-kinase–Protein Kinase B-Mammalian Target of Rapamycin Signaling Pathway. Asian-Australasian Journal of Animal Sciences, 29(5), 731–738. http://doi.org/10.5713/ajas.15.0415

SPECIFICATIONS

EVOM2™

±200 mV

300523

±200 mV

RESOLUTION	0.1 mV	0.1 mV
RESISTANCE RANGE	0 to 9999 Ω	0 to 100 KΩ
RESISTANCE RESOLUTION	1 Ω	10 Ω
AC SQUARE WAVE CURRENT	±10 μA nominal at	1 μA nominal at 12.5
AC SQUARE WAVE CORREIN	12.5 Hz	Hz
		geable 6V NiMH
POWER	2700 mAH batt	ery with external
	12 VDC supply	/ for recharging
NOMINAL BATTERY RUN TIME	10 h	nours
BNC OUTPUT	1-10 V (1 mV/Ω)	1-10 V (1 mV/10 Ω)
DIMENSIONS	19 x 11 x 6 cm (7.	.25" x 4.25" x 2.30")
WEIGHT	1.4 kg	g (3 lb.)
ELECTRODE CONNECTION	RJ-11connector	(telephone style)
TEST RESISTOR	Externa	Ι, 1000 Ω
ENVIDONIMENTAL DANICE	10-38°C	(50-100°F)

ORDERING INFORMATION

ENVIRONMENTAL RANGE

MEMBRANE VOLTAGE RANGE

EVOM2	Epithelial Tissue Voltohmmeter 9999Ω range (includes STX2)
300523	EVOM2™ 100 KΩ Range (includes STX2)
91799	Epithelial Tissue Voltohmmeter 9999Ω range (includes STX3)

0-90% non-condensing relative humidity

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

STX2	Replacement "Chopstick" Electrode Set
STX3	Adjustable Tip Spread "Chopstick" Electrode Set
3993	Electrode Adapter (for electrodes with 2 mm pins)
91736	Replacement Battery, Rechargeable NiMH
91750	EVOM2™ Test Resistor

Reproducible Resistance of Endothelial Tissue

TEER measurement of epithelial and endothelial cell cultures

Features

- Compatible with original EVOM™ and EVOM2™ meters
- · Adjustable apical electrode height
- Crystal clear glass chamber allows visualization of apical electrode positioning
- New insert holder with 120° tri-supports for three leg inserts
- Three sizes cover a range of well cup sizes from a variety of manufacturers

Benefits

- Stability and reproducibility superior to the STX2 electrodes to 1% tolerance
- Can be used with 6, 12 or 24 well plates with removable inserts
- Symmetrical electrode pattern disperses test current uniformly
- Simple test procedure to verify electrode performance

Applications

 TEER measurement for removable culture cup systems using EVOM2™ meters for endothelial and epithelial cell cultures

Using WPI's **EVOM2™** resistance meter, Endohm chambers provide reproducible resistance measurements of endothelial and epithelial monolayers in culture cups. Transfer cups from their culture wells to the Endohm chamber for measurement rather than using hand-held electrodes. The chamber and the cap each contain a pair of concentric electrodes: a voltage-sensing silver/silver chloride pellet in the center plus an annular current electrode. The height of the top electrode can be adjusted to fit cell culture cups of different manufacturers.

Make more precise measurements with Endohms

Endohm's symmetrically opposing circular disc electrodes, situated above and beneath the membrane, allow a more uniform current density to flow across the membrane than with STX2 electrodes. The background resistance of a blank insert is reduced from 150 Ω (when using WPl's hand-held STX2 electrodes) to less than 5 Ω . With Endohm's fixed electrode geometry, variation of readings on a given sample is reduced from 10-30 Ω with STX2 electrodes (depending on your experience) to 1-2 Ω . Compared with other resistance measurement methods, Endohm with EVOM2 $^{\rm IM}$ offers a much more convenient and economic solution to "leaky tissue" measurement. Because of the uniform density of the AC square wave current from EVOM2 $^{\rm IM}$, errors caused by electrode polarization or membrane capacitance are largely eliminated. Endohm together with EVOM2 $^{\rm IM}$ offers the most accurate and economical endothelial ohm meter now available. To date, cups from Costar, Millipore, and Falcon have been tested. Endohm chambers may be sterilized with EtO, alcohol or a bactericide. Do not autoclave.

References

Srinivasan, B., Kolli, A. R., Esch, M. B., Abaci, H. E., Shuler, M. L., & Hickman, J. J. (2015). TEER measurement techniques for in vitro barrier model systems. *Journal of Laboratory Automation*, 20(2), 107–26. http://doi.org/10.1177/2211068214561025

ENDOHM-12 COMPATIBILITY CHART				
Corning	Millipore	Membrane Diameter	Growth Surface Area	Membrane Pore Size
3401		12 mm	1.12 cm ²	0.4 μm
3402	PITP01250	12 mm	1.12 cm ²	3.0 µm
3403	PITT01250	12 mm	1.12 cm ²	3.0 µm
3493		12 mm	1.12 cm ²	0.4 µm
3494		12 mm	1.12 cm ²	3 µm
3460	PIHT15R48* PET Insert	12 mm	1.12 cm ²	0.4 µm
	PIRP15R48* PET Insert	12 mm	1.12 cm ²	1 µm
3462	PISP15R48* PET Insert	12 mm	1.12 cm ²	3 µm
	PIMP15R48* PET Insert	12 mm	1.12 cm ²	5 μm
	PIEP30R48* PIEP15R48* PIEP15R48* PET Insert	12 mm	1.12 cm ²	8 µm



	ENDOHM-24	COMPATIBILITY C	HART
Corning	Millipore	Material	Pore Size (µm)
3407		Polycarbonate	0.4
3801		Polycarbonate	0.4
3802		Polycarbonate	3.0
3412	PIHT30R48*	Polycarbonate	0.4
3414		Polycarbonate	3.0
	PITT03050	Polycarbonate	3.0
3428		Polycarbonate	8.0
3450		Polyester	0.4
3452		Polyester	3.0
3491		Collagen	0.4
3492		Collagen	3.0
	PICMORG50	Organotypic Insert	0.4
	PIHA03050	HA Insert	0.45
	PIHP03050	PCF Insert	0.4
	PICM03050	HA Mixed Cellulose Esters	0.4
	PIHT30R48*	PET Insert	0.4
	PIRP30R48*	PET Insert	1.0
	PISP30R48*	PET Insert	3.0
	PIMP30R48*	PET Insert	5.0
	PIEP30R48*	PET Insert	8

ENDOHM-6 COMPATIBILITY CHART				
Corning	Millipore	Membrane Diameter	Growth Surface Area	Membrane Pore Size
3470		6.5 mm	0.33 cm ²	0.4 µm
3472	PITP01250	6.5 mm	0.33 cm ²	3 µm
3413	PCF Insert	6.5 mm	0.33 cm ²	0.4 µm
3415	PITP 01250 PCF Insert	6.5 mm	0.33 cm ²	3 µm
3421		6.5 mm	0.33 cm ²	5 µm
3422	PIEP 01250 PCF Insert	6.5 mm	0.33 cm ²	8 µm
3495	PISP12R48* PIHT12R48* PET Insert	6.5 mm	0.33 cm ²	0.4 µm
	PIHA012 50 (HA Insert)	6.5 mm	0.33 cm ²	0.45 µm
	PICM012 50 (CM Insert)	6.5 mm	0.33 cm ²	0.4 µm
3496	PISP12R48* PET Insert	6.5 mm	0.33 cm ²	3 µm
	PIRP12R48* PET Insert	6.5 mm	0.33 cm ²	1 µm
	PIMP12R48* PET Insert	6.5 mm	0.33 cm ²	5 μm
	PIEP12R48* PET Insert	6.5 mm	0.33 cm ²	8 µm
	PIXP01250 PCF Insert	6.5 mm	0.33 cm ²	12 µm
	PITT01250 PIHP 01250		_	1.0 µm 3.0 µm

^{*} New Chambers Support Tri-Leg Inserts

	ORDERING INFORMATION
ENDOHM-6	
ENDOHM-1	2G Endohm for 12 mm Culture Cup (12 wells per plate)
ENDOHM-2	4GSNAP Endohm for 24 mm & Costar Snapwell™ Cup (6 wells per plate)
	Requires EVOM2™, EVOM™, EVOMX™ or Millicell ERS-2
53330-01	Replacement Endohm Cable

HTS Electrodes for Use with EVOM™ & EVOM2™

For High Throughput Screening (HTS) cell culture filter plates

Features

- Designed for use with 24-well HTS plate (Corning Costar and BD Falcon) and with 96-well plates (Millipore)
- Improved accuracy down to 5 Ω
- Sterilize with EtO, alcohol or bactericide

Benefits

- Smaller tip size than the STX2 electrode constructed for durability fits neatly into the keyhole shaped filter well
- Electrode design reduces chance of contamination

Applications

 STX-100 Electrodes are designed for TEER measurement in HTS culture plates using the EVOM2™

Semi-permeable HTS (high throughput screening) culture plates have become a standard tool for pharmaceutical and institutional research in epithelial transport. HTS culture plates have well cups which are bonded together into a single assembly and are not removable. This makes the plates ideal for automated applications, but imposes a significant inconvenience when TEER must be measured in the absence of an automated system. WPl's **STX100** series electrodes provide a cost-effective alternative to automation, allowing HTS well plates to be measured manually using a hand-held electrode.

The spatial orientation of an electrode during TEER measurement can have a significant effect on the resistance reading. When compared to the STX-2 "chopsticks" electrode, the STX100 series electrodes offer a technical advantage. The design of the STX100 series electrodes guarantees spatial repeatability. Each STX100 electrode is designed to self-align based on the form factor of the apical and basal access ports of the HTS plate. Each manufacturer of HTS culture systems has a unique form factor, so the STX100 electrodes are manufacturer specific.

WPI developed **STX100** electrodes for 24 and 96 well plates specific to Corning, BD Falcon and Millipore. Refer to the charts for compatibility information or contact your WPI sales representative for assistance. If an automated system is preferable for your application, WPI's automated TEER measurement system (REMS) on page 7 is an affordable solution.

	STX100 C COMPA	TIBILITY CHAR	T
Corning	Description	Pore Size	Membrane
3379	HTS Transwell-24	0.4 µm	PET
3378	HTS Transwell-24	0.4 µm	PET
3396	HTS Transwell-24	0.4 µm	
3397	HTS Transwell-24	0.4 µm	
3398	HTS Transwell-24	3.0 µm	
3399	HTS Transwell-24	3.0 µm	

	STX100F COMPATIBILITY CHART			
Falcon	Description	Pore Size	Membrane	
351180	BD Falcon (24 well) HTS Multiwell Insert System	1.0 µm	PET	
351181	BD Falcon (24 well) HTS Multiwell Insert System	1.0 µm	PET	
351182	BD Falcon (24 well) HTS Multiwell Insert System	3.0 µm	PET	
351183	BD Falcon HTS (24 well) Multiwell Insert System	3.0 µm	PET	
351184	BD Falcon HTS (24 well) Multiwell Insert System	8.0 µm	PET	
351185	BD Falcon (24 well) HTS Multiwell Insert System	8.0 µm	PET	
354803	BD BioCoat (24 well) HTS Fibrillar Collagen Multiwell Insert System	1.0 µm	PET	
354804	BD BioCoat (24 well) HTS Fibrillar Collagen Multiwell Insert System	1.0 µm	PC	



STX100 M COMPATIBILITY CHART				
Millipore Description Pore Size Membrane				
PSRP004R1	96-Well Plate	1.0 µm	PET	
PSHT004R5	96-Well Plate	0.4 µm	PCF	
PSRP004R5	96-Well Plate	1.0 µm	PET	
PSHT004S5	96-Well Plate	0.4 µm	PCF	
PSHT004R1	96-Well Plate	1.0 µm	PCF	

SIX 100 C90 COMPATIBLETT CHART			
Corning	Description	Pore Size	Membrane
3380	HTS Transwell-96 System	1.0 µm	PC
3392	HTS Transwell-96 System	1.0 µm	PET
3381	HTS Transwell-96 System	0.4 µm	PC
3391	HTS Transwell-96 System	0.4 µm	PC
3385	HTS Transwell-96 Well Plate	3.0 µm	PC
3386	HTS Transwell-96 Well Plate	3.0 µm	PC
3387	HTS Transwell-96 Well Plate	5.0 µm	PC
3388	HTS Transwell-96 Well Plate	5.0 µm	PC
3374	HTS Transwell-96 Well Plate	8.0 µm	PET
3384	HTS Transwell-96 Well Plate	8.0 µm	PET



	ORDERING INFORMATION
STX100C	STX100 for Corning Costar HTS Transwell-24
STX100F	STX100 for BD Falcon 24 well HTS Multiwell Insert System
	STX100 for Millipore MultiscreenTM HTS 96-Well Plate
STX100C96	STX100 for Corning HTS 96-Well Plate

Automated TEER Measurement System

PC-controlled high throughput TEER measurement for epithelial monolayer

Features

- Automates TEER measurement and data logging for use with HTS well plates
- · PC controlled positioning
- Data acquisition in LabView
- Manufacturer specific electrodes available for 24- and 96-well HTS plates
- Plate configuration files and sample sequences are user-definable
- Two user-defined rinse locations
- · Manual mode

Benefits

- Speed—capable of acquiring TEER data on a 96-well plate in less than five minutes
- Automation reduces the possibility for human error

Applications

 Automated measurement and data logging of TEER for 24 and 96 well HTS culture plates

The **REMS** AutoSampler automates measurements of TEER epithelial or endothelial monolayers cultured on HTS well plates. It is a PC controlled tissue resistance measurement system that offers reproducibility, accuracy, flexibility and ease-of-operation. Automated measurement of tissue resistance in cell culture microplates provides the advantages of speed, precision, decreased opportunity for contamination and the rapid availability of measured resistance data.

The main components of the **REMS** AutoSampler include: the robotic sampler that moves the electrode over each well of the microplate, the electrode which is located on the robotic arm, a base plate for the 24- and 96-well tray, a Windows-based data acquisition card, the **REMS** electrode interface unit and the **REMS** software to operate the system on a Windows-based computer.

Automate TEER measurements

The **REMS** AutoSampler automates TEER measurements that would otherwise be performed manually with WPl's **EVOM2**TM Epithelial Voltohmmeter. Automated tissue resistance measurements up to 20 k Ω can be performed on 24- or 96-well HTS microplates. See *www.wpiinc.com/REMS* for manufacturer plate compatibility.

Precisely and reproducibly positions electrode

The **REMS** AutoSampler will automatically measure and record tissue resistance from a user-specified matrix of culture wells on the microplate. According to the specified sequence, the robotic arm moves over the identified wells taking TEER measurements. By means of an x-y-z locating system, the electrode is positioned precisely into the well. The ability of the **REMS** AutoSampler to reproducibly position the electrode contributes to consistent TEER measurements. TEER data is incrementally stored as the electrode moves from one well to the next.

Compact electrode pair

The use of AC current to measure resistance provides several advantages over DC current, including:

- Absence of offset voltages on measurements
- There is a zero net current being passed through the membrane and, therefore, it is not adversely affected by a current charge
- No electrochemical deposition of electrode metal.

Rinse and calibration check stations

The **REMS** AutoSampler also features two rinse stations. If occasional rinsing of the **REMS** electrode is required, it may be sent to a rinse station by pressing the rinse station button on the menu bar.



References

Gallego, M., Grootaert, C., Mora, L., Aristoy, M. C., Van Camp, J., & Toldrá, F. (2016). Transepithelial transport of dry-cured ham peptides with ACE inhibitory activity through a Caco-2 cell monolayer. *Journal of Functional Foods*, 21, 388–395. http://doi.org/10.1016/j.jff.2015.11.046

Kamiloglu, S., Grootaert, C., Capanoglu, E., Ozkan, C., Smagghe, G., Raes, K., & Van Camp, J. (2016). Anti-inflammatory potential of black carrot (Daucus carota L.) polyphenols in a co-culture model of intestinal Caco-2 and endothelial EA.hy926 cells. *Molecular Nutrition & Food Research*. http://doi.org/10.1002/mnfr.201600455

Cain, M. D., Salimi, H., Gong, Y., Yang, L., Hamilton, S. L., Heffernan, J. R., ... Klein, R. S. (2017). Virus entry and replication in the brain precedes blood-brain barrier disruption during intranasal alphavirus infection. *Journal of Neuroimmunology*, 308, 118–130. http://doi.org/10.1016/j.jneuroim.2017.04.008

REMS AUTOSAMPLER SPECIFICATIONS

MEMBRANE RESISTANCE RANGE AC SQUARE WAVE CURRENT ELECTRODE POSITIONING ELECTRODE ARM SPEED

TYPICAL MEASUREMENT TIME 24-WELL SCAN PATTERN LINE VOLTAGE (User Specified) DIMENSIONS

WEIGHT

Auto-ranging or 0–2000 Ω and 0–20 k Ω \pm 20 μ A @ 12.5 Hz Resolution in X, Y and Z: \pm 1 mm X- and Y-axis: 250 mm/sec Z-axis: 247.3 mm/sec

1 min., 10 sec. Preset or user-defined 100/120 V or 220/240 V 53.5 × 43.7 × 37.1 cm (21½ × 17½ × 14½ in.) 24 kg (52 lb.)

ORDERING INFORMATION

SYS-REMS Automated Tissue Resistance Measuring System

Includes robot sampler, data acquisition board; computer, display, keyboard, mouse; software for Windows 7 or 10; and electrode (SPECIFY WHEN ORDERING).

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

REMS-24 Replacement Electrode for 24-well HTS Plate
REMS-96 Replacement Electrode for Millipore™ 96-well Plate
REMS-96C Replacement Electrode for Corning 96-well HTS plate
REMS-24M Replacement Electrode for Millipore 24-well HTS plate
See www.wpiinc.com/REMS detailed information.

Ussing System for Epithelial Research

Non-destructive TEER measurement for epithelial tissue

Features

- · Direct connect low-resistance electrodes
- · Simple operation, easy to control temperature and clean after use
- Luer type leak-free attachment of tubing and electrodes
- Recessed electrode ports avoid bubble formation
- Secure membrane holding by sharp stainless steel pins or O-ring
- Specialized chamber adapts cell culture insert (Costar Snapwell) for monolayer cell culture
- Chambers with rectangular openings for tubular tissues from small

Benefits

- Leak free design of Ussing chambers
- · Can be used with monolayer cell culture inserts
- Optional drains for guick evacuation of radioactive or toxic substances
- · Circulation reservoirs available in two sizes
- Control temperature with a circulating water bath (option) available

Applications

- Ion transport studies
- Nutrient transport studies

WPI's Ussing System offers researchers a quick, effective means of making low-resistance electrical connections to the Ussing chamber without need of long agar bridges or Calomel half-cells. Ag/AgCl half-cells screw into short tubes which plug firmly into place in the chamber's Luer ports. These direct-connect electrodes eliminate the inconvenience and expense of Calomel half-cells in open liquids. The system includes one Ussing Chamber (eight sizes available), support stand, electrode kit, glass circulation reservoir (two sizes available), and a tubing start-up kit (25 feet of 0.375-in. tubing, 10 feet of 0.156-in. tubing, plus four male Luer fittings, two compressor clamps, one Y-connector, and one clip). Sixteen possible system configurations are available. Components are also available separately. (Preamplifier in photo not included.)

Leak free design of Ussing chambers


WPI's classical Ussing Chambers are well established perfusion chambers that are easy to operate, easy to control temperature, and easy to clean after use. Ussing Chambers are machined from solid acrylic with eight entry ports for fluid lines, electrodes, or agar bridges. For easy, leak-free attachment of tubing and electrodes, all eight ports are Luer type. The four ports for voltage and current electrodes are recessed to prevent formation of air bubbles in the chamber. The fluid compartments in each side of the chamber are separated by the epithelial membrane being studied. Sharp stainless steel pins on one side of the chamber hold the membrane in position and mate with holes in the opposite chamber interface. (In the CHM4, tissue is held by an O-ring instead of pins.)

Can be used with monolayer cell culture inserts

The CHM5 chamber adapts the Costar 12 mm Snapwell, a cell culture insert for monolayer cell culture, into WPI's "classical" epithelial voltage clamp system. Classical Ussing Chambers have not been widely used for monolayer cell culture inserts, because most inserts have a very deep profile, limiting good fluid perfusion at the surface of the membrane and limiting voltage electrodes from measuring the potential close to the surface of the membrane. CHM5 solves these problems: Perfusion fluid is introduced into the chamber at an angle so that it flows directly to the surface of the membrane. The voltage electrode is also inserted into the chamber at an angle to reduce the distance between the surface of the membrane and the electrode.

Rectangular openings for tubular tissue

Two small chambers with rectangular openings are designed for tubular tissue from small animals such as the mouse intestinal tract membrane



Complete Ussing System includes stand, glass reservoir, electrodes, Ussing chamber and tubing.

(CHM6) and rat intestinal tract membrane (CHM7). The rectangular opening more closely matches the shape of the tissue than would a circular opening, significantly increasing the membrane area available for testing. The larger membrane area increases the transport rate of low permeability chemicals. It also reduces the electrical resistance of the system for easier current clamping.

Optional drains

Drains may be added to Ussing chambers to allow quick and complete evacuation of radioactive or toxic substances. To have drains added at the time of order, add a "D" to the part number (such as "USS1LD"). The cost of the drain will be added to the cost of the chamber or system ordered.

Cartridge electrodes

The Electrode Kit contains four voltage/current electrodes, plus four Luer-tipped cartridges. Electrodes are threaded and screw securely into the end of each cartridge. The Luer tip then plugs securely into the Luer openings of the chamber. The cable from each electrode terminates with a 2 mm pin which may be plugged into mos voltage/current clamps.

The miniature electrode-gel cartridge EKV and EKC Cartridge Electrodes is a small plastic tube with a male

Luer tip identical to those at the tip of hypodermic syringes. The tube may be filled with different gel materials. Agar is commonly used, but other gel materials may also be satisfactory.

CHM8 Chamber Clear acrylic chambers let you see Voltage Fluid out your experiment in progress . . Electrode Suture quality pins Fluid in minimize tissue damage Current Electrode Alignment **Optional drain** Dots for hazardous Guide Luer ports allow easy fitting material also Pins of fluid lines and electrodes available

Assembled chambers are 101.6 mm (4 in.) long.

Circulation reservoirs available in two sizes


Hand-blown borosilicate glass reservoirs with jacketed chambers for temperature control are available in two sizes — **5210** holds 20 mL per side, and **5362** (at left) holds 10 mL per side (useful when expensive chemicals are involved). Reservoir condenser caps prevent air bubbles and turbulence in fluid reservoirs.



The Julabo circulating bath is ideal for controlling temperatures of external systems. With a powerful 15 L/min flow rate, the pump provides optimum heat exchange. The tap water cooling feature is standard with a range of 20-100°C.

g bath is ideal for ures of external erful 15 L/min provides optimum tap water cooling with a range of

Reservoir openings and pins



References

M.Khera, G.T. Somogyi, S. Kiss, T.B. Boone, C.P. Smith "Botulinum toxin A inhibits ATP release from bladder urothelium after chronic spinal cord injury" *Neurochemistry International* 45. 2004: 987-993

McLamb, B. L., Gibson, A. J., Overman, E. L., Stahl, C., & Moeser, A. J. (2013). Early Weaning Stress in Pigs Impairs Innate Mucosal Immune Responses to Enterotoxigenic E. coli Challenge and Exacerbates Intestinal Injury and Clinical Disease. *PLoS ONE*, 8(4), e59838. http://doi.org/10.1371/journal.pone.0059838

Khera, M., Somogyi, G. T., Kiss, S., Boone, T. B., & Smith, C. P. (2004). Botulinum toxin A inhibits ATP release from bladder urothelium after chronic spinal cord injury. *Neurochemistry International*, 45(7), 987–93. http://doi.org/10.1016/j.neuint.2004.06.001

ORDERING INFORMATION

USSING SYSTEMS, LARGE RESERVOIR					
USS1L	1L Medium Chamber, Stand, Reservoir, Electrodes, Tubing				
USS2L	Small Chamber, Stand, Reservoir, Electrodes, Tubing				
USS3L	Large Chamber, Stand, Reservoir, Electrodes, Tubing				
USS4L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing				
USS5L	Snap Chamber, Stand, Reservoir, Electrodes, Tubing				
USS6L	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing				
USS7L	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing				
USS8L	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing				

USSING SYSTEMS, SMALL RESERVOIR

USS1S	Medium Chamber, Stand, Reservoir, Electrodes, Tubing
USS2S	Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS3S	Large Chamber, Stand, Reservoir, Electrodes, Tubing
USS4S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing
USS5S	Snap Chamber, Stand, Reservoir, Electrodes, Tubing
USS6S	Small Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS7S	Large Rectangular Chamber, Stand, Reservoir, Electrodes, Tubing
USS8S	Extra Small Chamber, Stand, Reservoir, Electrodes, Tubing

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

System cor	mponents also available separately:
xxxxD	Drain option (add "D" to part number of chamber or system)
СНМ1	Medium Chamber
CHM2	Small Chamber
СНМЗ	Large Chamber
СНМ4	Extra Small Chamber with O-Ring Seal
СНМ5	Snap Chamber (fits Costar Snapwell cups)
СНМ6	Small Rectangular Chamber
СНМ7	Large Rectangular Chamber
CHM8	Extra Small Chamber with Mounting Pins
EK1	Ussing Electrode Kit (2 voltage, 2 current)
EKC	Extra Ussing Current Electrode (red) (each)
EKV	Extra Ussing Voltage Electrode (blue) (each)
5210	Large Glass Circulation Reservoir, (20 mL per side)
5233	Replacement Condenser for 5210
5362	Small Glass Circulation Reservoir, (10 mL per side)
5361	Replacement Condenser for 5362
3955	EKV Cartridges, 35 mm (pkg. of 12)
3960	EKC Cartridges, 58 mm (pkg. of 12)
3669	Tubing Kit (flexible hose and Luer fittings)
3579-20	Replacement Luer fittings for tubing connections (pkg. of 20)
5153	Support Stand
3485	Post Mounting Kit for Preamp
505060	Julabo Circulating Bath, 5L volume, 115V
505061	Julabo Circulating Bath, 13L volume, 115V
505062	Julabo Circulating Bath, 5L volume, 230V
505063	Julabo Circulating Bath, 13L volume, 230V
·	<u> </u>

SPECIFICATIONS							
СНМ1	CHM2	СНМЗ	CHM4	CHM5	СНМ6	СНМ7	СНМ8
(Medium)	(Small)	(Large)	(Extra Small)	(Snap)	(Rect., Small)	(Rect., Large)	(Extra Small)
12 mm	9 mm	13.5 mm	4 mm	12 mm	5 x 14.5 mm	7 x 30 mm	4 mm
1.0 mL	0.75 mL	1.2 mL	0.5 mL	1.7 mL	0.8 mL	5.5 mL	0.5 mL
17 mm	12 mm	18.5 mm	6 mm*	N/A	7 x 16.5 mm	9 x 32 mm	5.5 mm
113 mm ²	63.5 mm ²	143 mm ²	12.6 mm ²	113 mm ²	67.1 mm ²	199.5 mm ²	12.6 mm ²
	(Medium) 12 mm 1.0 mL 17 mm	(Medium) (Small) 12 mm 9 mm 1.0 mL 0.75 mL 17 mm 12 mm	CHM1 CHM2 CHM3 (Medium) (Small) (Large) 12 mm 9 mm 13.5 mm 1.0 mL 0.75 mL 1.2 mL 17 mm 12 mm 18.5 mm	CHM1 CHM2 CHM3 CHM4 (Medium) (Small) (Large) (Extra Small) 12 mm 9 mm 13.5 mm 4 mm 1.0 mL 0.75 mL 1.2 mL 0.5 mL 17 mm 12 mm 18.5 mm 6 mm*	CHM1 CHM2 CHM3 CHM4 CHM5 (Medium) (Small) (Large) (Extra Small) (Snap) 12 mm 9 mm 13.5 mm 4 mm 12 mm 1.0 mL 0.75 mL 1.2 mL 0.5 mL 1.7 mL 17 mm 12 mm 18.5 mm 6 mm* N/A 113 mm² 63.5 mm² 143 mm² 12.6 mm² 113 mm²	CHM1 CHM2 CHM3 CHM4 CHM5 CHM6 (Medium) (Small) (Large) (Extra Small) (Snap) (Rect., Small) 12 mm 9 mm 13.5 mm 4 mm 12 mm 5 x 14.5 mm 1.0 mL 0.75 mL 1.2 mL 0.5 mL 1.7 mL 0.8 mL 17 mm 12 mm 18.5 mm 6 mm* N/A 7 x 16.5 mm 113 mm² 63.5 mm² 143 mm² 12.6 mm² 113 mm² 67.1 mm²	CHM1 CHM2 CHM3 CHM4 CHM5 CHM6 CHM7 (Medium) (Small) (Large) (Extra Small) (Snap) (Rect., Small) (Rect., Large) 12 mm 9 mm 13.5 mm 4 mm 12 mm 5 x 14.5 mm 7 x 30 mm 1.0 mL 0.75 mL 1.2 mL 0.5 mL 1.7 mL 0.8 mL 5.5 mL 17 mm 12 mm 18.5 mm 6 mm* N/A 7 x 16.5 mm 9 x 32 mm 113 mm² 63.5 mm² 143 mm² 12.6 mm² 113 mm² 67.1 mm² 199.5 mm²

^{*}O-ring diam.

Muscle Testing Platform

Customizable system for muscle physiology research

Features

- · Heated cuvette and micrometer
- Study intact muscle fiber bundles*
- Data recording and analysis included
- · Modular design
- · Corrosion-free

Benefits

• System is completely customizable with a host of accessories

Applications

- Electrical stimulation, including twitch and tetany, to analyze force response in amplitude; and kinetics, like contraction and relaxation times and velocities
- Intracellular calcium concentration/ distribution and muscle force in intact muscle fiber bundles.* (Requires a Biofluorometer.)

These applications may require specific components or other electronics that are compatible with this system.

The SI-MT-L is the standard muscle research platform, which can be configured to study intact muscle fiber bundles, muscle strips and small whole muscles. Modular design allows the system to be configured for turnkey solutions for specific applications. The system is built on a solid platform making precise mechanical and optical measurements easy. Like the Cell Tester and the SI-HTB systems, SI-MT-L and SI-MTM-L use SI-KG optical force transducers and are constructed with corrosion-free materials (stainless

steel, anodized aluminum or plastic).

The system is supplied with a **LabTrax-MDAC** data acquisition system to:


- Record signals from the force transducer, the motor position monitor or any further relevant physiological signal using up to eight analog inputs.
- Control the stimulator, or an external stimulus isolator, through an analog output
- Control the position of the linear motor through a second analog output
- · Control other devices through two digital outputs
- · Program numerous customized protocols

Customize the system to suit your needs

The flexibility of the Muscle Testing Platform allows for complete customization. For example, you can add:

- Linear motor for muscle length perturbation studies
- Isolated constant current stimulator for twitch and tetanus stimulation
- Biofluorometer for intracellular calcium studies
- Additional accessories like force transducers, data acquisition system





MDAC Protocol segment, showing the combination of basic signal patterns to generate an experimental protocol.

ORDERING INFORMATION

SYSTEMS AVAILABLE FOR INTACT MUSCLE SAMPLES INCLUDE:

SI-MTM-L System with a Univette heated cuvette & linear motor
SI-MT-L System with a Univette heated cuvette & digital micrometer

Systems include: Base plate; force transducer and stand, digital micrometer or linear motor (optional), micrometer or motor stand, cuvette and table, oxygenation system, evacuation system, signal conditioning system with transducer amplifier and temperature control module, anti-oscillation module and linear motor amplifier (if applicable), LabTrax-MDAC data acquisiton system data acquisition system with MDAC software for recording force and position signals and for controlling stimulation and the movement of the linear motor.

*Contact WPI for help evaluating your specific experimental requirements.

Horizontal Tissue Bath

Research system for higher throughput of complex pharmacological/

physiological assays

Features

- Two channel system for increased productivity, easily expanded to add channels
- Fully independent heating and fluid control for each channel
- · Low profile/small footprint
- Variable volume, chemically inert Teflon bath with shape configurations from variable to fixed 500 µL–10 mL*
- Modular, space-saving, blade-style electronics
- Large variety of force transducers covering mN–N forces
- Can be combined with automated fluid control systems
- Add an electrometer like the WPI Duo773 or Electro705

Benefits

- Tissue bath for volumes as small as 500 µL
- Low profile, small footprint to reduce bench space requirements
- Versatile system is easy to upgrade with added channels and options

Applications

 Tissue fluorescence (calcium, ATPase via NADH, FAD). (Requires the multi-purpose Biofluorometer* (SI-BF-100). See page 20.)

These applications may require specific components or other electronics that are compatible with this system.

The SI-H Horizontal Tissue/Organ Bath system (SI-HTB) combines the ease of use and productivity of a traditional vertical organ bath with the more advantageous features of single tissue physiology platforms.

Bath options 500 µL up

The **SI-HTB** system breaks through the (large) volume limitations of the traditional organ bath, allowing volumes as low as 500 μ L in an inert, Teflon-based bath. The bath design allows multiple shape options for thick, long, flat and thin tissue. When pharmaceuticals are available in precious, small amounts, you will appreciate this standard feature. A wide range of transducers and tissue mounting supports allow a great variety of tissue shape, volume and size.


Low profile

The low profile and small footprint of the bath system, combined with the modular, space-saving, chassis-mounted design of the electronics, reduces the bench space requirement up to 4-fold when compared with standard 4-channel organ baths.

Versatile system, easy to upgrade

The **SI-HTB** combines advanced physiological techniques with the throughput needed in pharmacological assays in one flexible platform. Upgrades to four or more channels are easy and economical.

The motor option (SI-MOT) turns your system into a tissue work-out station with isotonic and eccentric force measurement capabilities. Nearly all established myo-mechanical tests from stretch-release to work-loops and muscle fatigue are now possible in a single organ bath system. Some of these procedures require automated length changing methods, which can be conveniently programmed using WPI's LabTrax-MDAC data acquisition system.



Options

The solid horizontal tissue bath design is ideal for combination with electrophysiology on the same platform. Intracellular measurements can share the stable solid base of the bath system.

WPI's new fiber-optic based, multi-channel Biofluorometer allows for tissue fluorescence measurements (calcium, NADH, FAD) on the **SI-HTB** platform.

Now, you can design a system to meet your needs and budget. And, it is fully upgradeable in the future.

WPI's 16-bit, full speed, Labview-based Muscle Data Acquisition system **LABTRAX-MDAC** is perfect for this platform.

*Contact WPI for help evaluating your specific experimental requirements.

ORDERING INFORMATION

SI-HTB2 Horizontal Tissue Bath, 2-Channel System

2-Channel SI-HTB platform for isometric force (1), SI-KG Force Trandsucers (2), SI-BAM21-LCB Optical Force Transducer Amplifiers (2), SI-TCM2B 2-Channel Temperature Controller (1), SI-BMFA power frame enclosure (1), LabTrax 8/16 with MDAC Data Acquisition software

SI-HTB4 Horizontal Tissue Bath, 4-Channel System

2-Channel SI-HTB platform for isometric force (2), SI-KG Force Trandsucers (4), SI-BAM21-LCB Optical Force Transducer Amplifiers, (4) SI-TCM2B 2-Channel Temperature Controller (2), SI-BMFA power frame enclosure (1), LabTrax 8/16 with MDAC Data Acquisition software (1)

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

LABTRAX-MDAC LabTrax 8/16 8-Channel Data Acquisition with MDAC Software				
SI-AOSUB	Anti-Oscillation Module			
SI-PF100	Programmable Filter			
MINISTAR	Miniature Peristaltic Pump, 1-channel			
PERIPRO-2HS	Peri-Star™ Pro, 2-channel, High Rate, Large Tubing (110-220V)			
PERIPRO-4HS	Peri-Star™ Pro, 4-channel, High Rate, Large Tubing (110-220V)			
PERIPRO-4LS	Peri-Star™ Pro, 4-channel, Low Rate, Small Tubing (110-220V)			
PERIPRO-8LS	Peri-Star™ Pro, 8-channel, Low Rate, Small Tubing (110-220V)			
801566	Mini Vacuum Pump (110V)			
801963	Mini Vacuum Pump (220V)			
SI-COLUB	Constant Load Module			
97204	Pinger for the Anti-Oscillation Module			

Cell Tester System

Perform muscle physiology tests on single skeletal muscle fibers

Features

- Fits ANY inverted microscope
- Unique rotational stage-improves the experimental throughput
- Two integral piezo manipulators included

Benefits

- Simultaneous use of the multi-purpose biofluorometer* (SI-BF-100) for tissue fluorescence (calcium, ATPase via NADH, FAD)
- For use with small multi-cellular preparations or skinned muscle fibers*
- Stimulate and perform perturbations of muscle
- Options available for customization of the

Applications

- · Stretching and relaxation of cells with nanometer resolution
- Classic cross bridge cycling studies in cardiac
- Skeletal muscle fiber axial stretch and isometric
- Intracellular calcium concentration/distribution in muscle fibers as muscle force is measured
- ATPase activity as muscle force is measured

The SI-CTS200 system, which is the result of blending the latest technologies in electronics, mechanics and optics, permits researchers to investigate living systems at a new level of observation. The Cell Tester provides researchers with the comprehensive ability to investigate and characterize the physiological, biomechanical and biophysical properties of single isolated living cells.

Study multi-cellular or skinned fibers

The **SI-CTS200** systems can be used on small multi-cellular preparations or skinned muscle fibers or strips. It is designed to sit on the stage of any standard, research-level, inverted microscope while maintaining the optical path of the microscope for simultaneous fluorescence or confocal

Stimulate and perform perturbations of cells

The Cell Tester is an integrated system of components needed to maintain and handle muscle fibers, stimulate and perform pertubations of the cells, and detect, amplify, and record signals, like contractile force.

Customize your system with options

- Optical force transducer for measuring with nanonewton sensitivity and integrated cell attachment system.
- Equipped with a nanomotor for stretching and relaxing cells with nanometer resolution.
- A rotating cuvette system for easy alignment of cells increases productivity. It is designed to orient cells in the XY plane so that no physical manipulation of the position of the tissue itself is required.
- Interchangeable bath inserts provide a range of options for the handling of live fibers.
- LabTrax-MDAC data acquisition system records and control your experiments, using up to eight analog inputs.
- Biofluorometer (See page 20.)



References

Prosser BL, Ward CW, Lederer WJ. X-ROS signaling: rapid mechanochemo transduction in heart. Science. 2011 Sep 9;333(6048):1440-5. PMID: 21903813

FORCE TRANSDUCER SPECIFICATIONS

RANGE 0-5 mN (0-0.5 g) FORCE RESOLUTION 20 nN at 10x gain COMPLIANCE 10 µm/mN NOISE $0.3 \mu N$

RESONANCE FREQUENCY 250 Hz (This is eliminated from the

measurement by the AOSUB)

TIME RESOLUTION

Resolutions were determined while using the SI-AOSUB anti-oscillation filter.

NANOMOTOR SPECIFICATIONS

TOTAL TRAVEL ±90 µm RESOLUTION 20 nm SMALLEST STEP

INPUT ±10V (calibrated at 10 μm/V)

CELL TESTER PLATFORM SPECIFICATIONS

DISTANCE (bottom of the base plate to the bottom of glass) 0.508mm

ORDERING INFORMATION

SI-CTS200 Complete Cell Tester System

System Includes: Rotating Cuvette System; Micromanipulator System; Signal Conditioning Amplifier with four modules: Optical Transducer Amplifier; Temperature Controller; Anti-Oscillation Unit; Position Controller; data acquisition system with MDAC Software for recording, controlling stimulation and nanomotor position; Force Transducer of choice; Nanomotor; Glass Fiber Cell Mounts (1 set)

SI-CTS200B Cell Tester, Non-Rotating, no Micromanipulator System SI-CTS200A Cell Tester, Manual Platform, no Micromanipulator System

For additional information, including application notes, go to wpiinc.com/celltester

*Contact WPI for help evaluating your specific experimental requirements.

Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!



SI-BMFA power frame enclosure plus optional modules (pages 14-17).

Features

- Ergonomic design
- 8-Channel
- Small footprint
- Backplane design includes provision for configurable communication between modules

Benefits

- Up to eight modules connected through the backplane of the chassis allows researchers to assemble the set of electronics they need for their own custom applications
- Standard configuration options are available for Muscle Tester platforms

Applications

 Process transduction of physical signals, displacement and optical force transducer outputs

The transduction of physical signals in the last decade has increasingly moved in the direction of the computer with an electronic amplifier. Further signal conditioning and analysis of sampled raw data is then handled efficiently in software, like WPI's **MDAC** package. The software preserves the raw data and is highly reliable. Operations such as integration, differentiation, filtration and even waveform generation are now efficiently handled in software. On the other hand, however, the transduction of physical signals such as bio-potentials, force, temperature, pressure or ionic concentrations must be measured with an electronic amplifier.

SIH/WPI's physiology amplifier system focuses on this idea and provides a flexible electronic platform intended to process the transduction of physical signals, displacement transducer outputs from force transducer signals. This platform simply focuses on the reliable transduction of the electronic signal and provides a convenient passage for the translation of real world signals to a computer for analysis.

Eight modules connected through backplane

The system consists of a chassis with eight slots on the backplane and includes an ultra quiet, shielded power supply. All of the module outputs

are routed to rear panel connectors. If you prefer, outputs may be routed internally to the inputs of other modules. The system has a small footprint and may be stacked to provide as many channels as you need.

The **SI-BMFA** Power Frame is the foundation of the SI modular physiology suite. It incorporates a robust power supply that can accommodate up to eight physiology modules, which can be mixed or matched in any combination. Modules are quick and easy to install, thanks to an innovative and mechanically solid card rail system.

The system is flexible and configurable. A variety of modules are available for the Signal Conditioning Amplifier System, and you can mix and match the modules to suit your requirements.

Optical Transducer Amplifier—SI-BAM21-LCB (page 15)

Programmable Filter Module—SI-PF100B (page 17)

Linear Motor Control Module—SI-MOTDB (page 14)

Temperature Control Module—SI-TCM2B (page 16)

Anti-Oscillation Module—SI-AOSUB (page 17)

Constant Load Module—SI-COLUB (page 16)

Standard options available for muscle testers

When the system is ordered with SI-MT (Muscle Tester) system, the Signal Conditioning Amplifier System (chassis) is configured with an SI-BAM21-LCB. Optional modules include an SI-TCM2B Temperature Control Module, an SI-MOTDB Linear Motor Controller, an SI-PF100B Programmable Filter Mondule, the SI-AOSUB Anti-Oscillation Module and the SI-COLUB Constant Load Unit. The Temperature Control Module and Linear Motor Controller require two slots each on the chassis backplane.

ORDERING INFORMATION				
SI-BMFA	Power Frame Enclosure			
SI-MOT-M1	SI-MOT-MT Linear Motor for SI-HTBM system			
SI-MOTDB Linear Motor Controller				
SI-BAM21-LCB Optical Transducer Amplifier				
SI-TCM2B 2-Channel Temperature Control Module				
SI-COLUB Constant Load Module				
	SI-AOSUB Anti-Oscillation Module			
SI-PF100B	Programmable Filter Module			

Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!

Linear Motor Control Module

Features

- Powers the motor and provides an output indicating the actual motor position
- Connects to Analog to Digital Converter output of the computer or data acquisition system (like LabTrax-MDAC) to allow the control of the programmed waveform and timing of the motor control
- Input range of ±10 VDC
- Linear motor position is determined by a DC value applied from the Position In port
- 2-Slot control module for a linear motor

Benefits

- Over current protection that automatically shuts down when the supply voltage dips below the reference value
- External ±10 VDC position command input for control by a constant load module or data acquisition system



SI-MOTDB

The SI-H Linear Motor Controller is designed for use with the SI-H line of muscle physiology research platforms. For systems that require a linear motor, this unit provides the precision control of the motor. A linear motor is required to perform mechanical muscle testing such as slack-tests, isotonic releases, constant velocity releases, stretch releases, after-loaded contractions and eccentric contractions (intact muscle). The position of the linear motor is determined by a combination of the data from the controller indicating the current position and the DC value applied to the front panel at the Position In port. The applied Position In signal can be provided by a data acquisition system (LabTrax-MDAC to use standard or customized protocols). The data acquisition analog output signal is set to define the waveform and timing pattern of force to be applied to the sample.



The SI-MOT motor can be used with any SI-MTM-L Muscle Tester.

SI-MOTDB SPECIFICATIONS

POWER REQUIREMENTS 12 V DC provided by the chassis

INPLIT ±10 V DC TRAVEL 1 mm/2 VDC current

MAXIMUM TRAVEL 7 mm (±3.5 mm from center of travel)

ORDERING INFORMATION

Perfect for: Ratiometric Calcium & ATPase

SI-MOT-MT Linear Motor for SI-MTM-L system

SI-MOTDB Linear Motor Controller

Applications

- · Slack test
- Isotonic release
- Constant velocity release
- Stretch release
- Eccentric/concentric contractions (intact muscle)

Biofluorometer

Reliable, simplified and affordable LED based fluorometer

Recent advancements in optics and LED technology simplify ratiometric calcium imaging, making this equipment more affordable. A breakthrough in WPI patented technology allows the SI-BF-100 to use

wavelengths below 380 nm and produce more light in those spectra. This technology significantly cuts the cost of photometric calcium imaging without sacrificing resolution or quality.

Streamlined system to reduce errors

Up till now, calcium imaging systems have been required to compensate for errors and noise introduced by the complexity of their design. The systems require mechanical filters and use expensive xenon or mercury light sources. The beauty of the SI-BF-100 is its simplicity. The elegance of its design reduces the noise introduced into the system and the errors inherent in traditional systems.



This single wavelength spectrophotometer can be customized for your specific application. See page 20 for details.

Optical Force Transducer Amplifier Module

Features

- · Designed for use with SI-H Muscle Tester Platforms
- Rapid auto zeroing function with fine offset adjustment
- Offset indicator LED's
- Multiple gain ranges with adjustable fine tuning for precise calibration
- 1X, 2X, 5X, 10X gains. Optional factory setting allows for 10X, 20X, 50X and 100X gains

Benefits

- Single amplifier that spans the entire SI-KG optical force transducer line
- Manual calibration to the ±10 V measurement range, for covering maximum forces ranging from 50 mN to 2N. When using the SI-KG7A, can calibrate as small as 5 mN.

Applications

• Muscle physiology studies, including isotonic studies

The SI-BAM21-LCB amplifier for SI-KG Optical Force Transducers is used in conjunction with the SI-H muscle physiology systems. The SI-BAM21-LCB powers the force transducer and converts the output of the transducer to an amplified analog voltage that is proportional to the force applied to the force transducer. The output signal can be multiplied by a factor of 1, 2, 5 or 10 to provide better resolution for a minimal change in applied force.

NOTE: An optional factory setting increases the multiplier by a factor of 10, allowing the signal to be multiplied by 10, 20, 50 and 100.

The SI-BAM21-LCB amplifier works with SI-KG optical force transducers to:

- Generate an analog output (-10 VDC to +10 VDC) that is proportional to the force applied to the tissue sample.
- Supply a DC voltage that powers the SI-KG force transducer to which it is connected.

Also available in a single stand-alone enclosure, either version provides an incredibly quiet, linear and stable transducer signal to your data recording system.



SI-BAM21-LCB

How the SI-BAM21-LCB amplifier works

In a typical setup, a muscle is held by a force transducer. The force transducer is connected to the SI-BAM21-LCB. As the muscle contracts or releases, the transducer converts the force into an electrical current signal which is proportional to the force applied to the transducer. The SI-BAM21-LCB converts the current signal into a voltage signal that can be displayed on the screen of the recording device.

Before initiating an experiment, the SI-BAM21-LCB must first be zeroed. This sets the baseline for measurements to follow.

The output signal is buffered and multiplied by 1, 2, 5 or 10, depending on the Gain switch setting on the front panel of the amplifier module. The ×10 setting is useful when output signals are extremely small. Finally, the force proportional signal is sent through the output amplifier circuit.

The analog output has a range of -10 V to +10 V that drives the LABTRAX-MDAC data acquisition system, multimeter or oscilloscope.

NOTE: When the Signal Conditioning Amplifier System is configured at the factory for an SI-HTB or SI-MT Muscle Tester system, the signal is routed internally from the SI-BAM21-LCB module to the SI-AOSUB module.

SI-BAM21-LCB SPECIFICATIONS

INPUT CONFIGURATION Current to voltage converter

1×. 2×. 5×. 10× **GAIN**

Optional factory setting: 10×, 20×, 50×, 100× INPUT OFFSET ADJUSTMENT ± 2.0 VDC

OUTPUT IMPEDANCE 470 O

POWFR 12 VDC provided by chassis

OUTPUT RANGE

ORDERING INFORMATION

SI-BAM21-LCB Optical Transducer Amplifier

Stand alone Optical Force Transducer Amplifier

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

LABTRAX-MDAC	LabTrax 8/16 with MDAC software
2851	BNC Cable
SI-KG2	0-2 N Force Transducer
SI-KG2B	0-200 mN Force Transducer
SI-KG4	0-50 mN Force Transducer
SI-KG7a	0-5 mN Force Transducer
SI-KG7b	0-10 mN Force Transducer
SI-KGxx	contact WPI for specialty transducers with different range
	6 0 :: 15 T 1 33

See Optical Force Transducers on page 22.



The stand-alone SI-BAM21-LC Force Transducer Amplifier has all the capabilities of the Signal Conditioning Amplifier module (SI-BAM21-LCB).

Signal Conditioning Amplifier System

Choose the amplifier modules you need to measure nearly anything!

Temperature Control Module

Features

- · Uses PID control to maintain a constant temperature with ±0.1°C tolerance
- Easy to control with simple interface
- Also available as a standalone device
- 2-Slot control module for maintaining temperature

Benefits

- · Controls two cuvettes simultaneously
- · User defined high and low alarm warnings

Applications

· Muscle physiology applications using SI-H muscle tester platforms



SI-TCM2B

The SI-H Temperature Control Unit is designed for use with the SI-H line of muscle physiology research platforms. It maintains the temperature of an SI-H cuvette up to 45°C. It is accurate to 0.1°C. The circuit is appropriate to RTD (resistive temperature device) applications. It linearly converts a temperature reading to a voltage that is displayed as a temperature on the SI-TCM2 and can be recorded. This unit is available in a stand-alone model and as a module for the Signal Conditioning Amplifier System backplane.

Call for details and pricing information.

SI-TCM2 SPECIFICATIONS

INPUT CONFIGURATION Current to voltage converter

POWER REQUIREMENTS 12 V DC at 2.5A 50/60 Hz wall adaptor,

2.5 mm ID/5.5 mm OD with positive center DC barrel (included-WPI

#801513)

OPERATING TEMPERATURE RANGE Room temperature

MAXIMUM TEMPERATURE 45°C DISPLAY PRECISION 0.1°C CONTROLLER RESOLUTION 0.1°C

CUVETTE TEMPERATURE SENSOR 1000 Ω RTD (1000 Ω at 0°C)

ORDERING INFORMATION

2-Channel Temperature Control Module SI-TCM2B 2-Channel Temperature Control, Stand-Alone

Constant Load Module

Features

- Offers three modes including Constant Load, External Loop and Bypass
- Can be configured using LabTrax-MDAC data acquisition system

Benefits

· Precise time resolution for feedback control in isotonic muscle testing


Applications

· Muscle physiology applications using SI-H muscle tester platforms

The SI-COLUB Constant Load Module for performing constant load experiments, has augmented flexibility. In its primary mode (Constant Load) the unit takes an external trigger command from the force transducer to perform a constant load cycle. In addition, the module allows for a different external trigger or you can completely bypass the module without having to switch cabling. The Constant Load Module lets you maintain a constant force, muscle length or sarcomere length rather than keep the total length of the preparation constant during an isotonic test. This is accomplished using a feedback loop.



The SI-COLUB monitors a designated parameter to determine how much force is necessary. It also monitors a feedback signal. The motor position command signal driving the motor is constantly adjusted to drive the feedback signal to the commanded setpoint.



SI-COLUB SPECIFICATIONS

COMMAND REQUEST ±10 V **FEEDBACK** ±10 V MOTOR OUTPUT ±10 V

POWER REQUIREMENTS 12 V DC provided by the chassis

ORDERING INFORMATION

SI-COLUB Constant Load Module

Anti-Oscillation Module

Features

- Neutralizes nearly 100% of the unavoidable transducer oscillation in optical force transducer
- High time-resolution of the corrected force signal
- Additional powerful signal smoothing with high time fidelity for all higher harmonics

Benefits

- Goes beyond traditional low-pass filtering of an oscillated signal
- · No phase-shift in the force signal
- High time fidelity in following the applied step signal

Applications


• All types of optical-based force measurements

Within a force transducer, the applied force displaces an elastic component of the transducer. This displacement is transformed into an electrical signal. If the displacement of the elastic component is small and if the displacement is linearly transformed into an electrical signal, the external force and electrical signal are linearly related to each other.

However, shortly after the force changes rapidly the electrical signal does not follow the applied force

SI-AOSUB linearly, as the elastic element and the force transducer pin have a mass. So, the system acts not only as an elastic device but also as a harmonic oscillator, notably when the force changes rapidly, resulting in oscillation of the system around the new displacement level (ringing phenomenon).

To address this problem, WPI's Anti-Oscillation Unit (SI-AOSUB) uses a genuine electronic approach to remove the unavoidable transducer oscillation, based on electronic differentiation of the force signal and subsequent low-pass filtering of the resulting high-frequency signal.



The upper trace is a force transient obtained directly from the bridge amplifier output, and the lower trace shows the signal after it passes through the anti-oscillation module.



To adjust the anti-oscillation filter properly, the transducer is excited at its resonance frequency using a magnetic driver or pulser (WPI #97204). The Pinger (WPI #97204) is included with the anti-oscillation module.

SI-AOSUB SPECIFICATIONS

ORDERING INFORMATION

12 VDC provided by chassis **POWFR**

INIPI IT + 10 VDC

PULSER OUTPUT 0 - 10 VDC adjustable 0.1 Hz - 4.0 KHz

0.1 Hz - 2.0 KHz DAMPING FREQUENCY RANGE **OUTPUT RANGE**

± 10 VDC

SI-AOSUB Anti-Oscillation Module

Pinger for Anti-Oscillation Module 97204

2851 **BNC Cable** Programmable Filter Module

Features

AOSU

- Low Pass Filter
- · Bessel and Butterworth filters
- Programmable cutoff frequency
- Signal may be routed through the backplane or through the front panel BNCs

Benefits

General purpose filter module for nearly all type of noisy signal

Applications

• Muscle physiology applications using SI-H muscle tester platforms

When you use a motor, an SI-PF100B Programmable Filter is necessary to minimize the natural vibration. It is designed so you can eliminate the resonance frequency without affecting the signal of interest. It is a low pass filter set to pass signals of interest below the specified frequency. It can be calibrated from 5 to 1,000 Hz.

You may select either a Bessel or a Butterworth filter. Then, you must carefully select the cutoff frequency based on the typical resonance frequency of your force transducer and your own experimental setup.

When the Signal Conditioning Amplifier System electronics are configured at the factory with an SI-PF100B Programmable Filter, the signal is routed internally from the amplifier module (SI-BAM21-LCB) to the SI-PF100B. If you prefer, the signal may be routed from the amplifier through the ports on the front panel of the Programmable Filter using a standard BNC cable.



SI-PF100 SPECIFICATIONS

POWER 12 VDC provided by chassis

INPUT ± 10 VDC **CUTOFF FREQUENCY RANGE** 5-1,000 Hz

FILTER TYPES BESSEL, BUTTERWORTH

ORDERING INFORMATION

SI-PF100B Programmable Filter Module

SI-PF100 Programmable Filter in Stand-Alone Enclosure



The Programmable Filter Module as stand-alone unit (**SI-PF100**) can be used for filtering two independent input signals, with the installation of a second SI-PF100B module (not included).

ORDER TOLL-FREE: (866) 606-1974 (*U.S. only*) **⋅ Tel:** (941) 371-1003 **⋅ Fax:** (941) 377-5428 **⋅ E-mail:** sales@wpiinc.com

Optical Sarcomere Spacing System

Measure sarcomere spacing with nanometer accuracy

Features

- Monitors high-speed dynamic sarcomere length changes
- Sarcomere spacing measurement with nanometer accuracy
- NIH open-source µManager* software interface between USB 3.0 camera and microscope
- Synchronous use of Sarcomere Spacing plug-in and µManager software
- User defined Region of Interest (ROI) for image capturing from whole cell down to a minimum of 6 consecutive sarcomeres
- Compatible with Windows 10

Benefits

- Camera image capturing using the standard µManager software familiar to researchers
- · Use Image toolboxes for complementary image processing
- Live image processing and display of length changes in sarcomere spacing of contracting muscle cells and fibers
- Synchronous storage of live images and processed data in userfriendly MS Excel format, for further offline analysis
- User defined ROI setting allows focusing on specific cell regions and changing image capturing speed
- High-speed USB 3.0 camera fits to all standard inverted microscopes with C-Mount adapter


Applications

- · Cardiac muscle physiology
- · Skeletal muscle physiology
- · Quality control of muscle slices in meat industry

WPl's Optical Sarcomere Spacing System quantifies live images to detect length changes in sarcomere spacing in nanometer accuracy during twitch or tetanus contractions of muscle cells/fibers. The NIH μ Manager software interface and WPl's image processing plug-in for sarcomere spacing detection monitor and measure sarcomere length changes of the contracting muscle cells in real time with an optimal frame rate of 500 FPS.

The Optical Sarcomere Spacing System fits perfectly with WPI's Cell Tester system (SI-CTS200) to capture twitch contraction of cardio myocytes, synchronous with the force signal. Extend the quantification by adding measurements of Ca²⁺ or ATPase concentration of contracting muscle cells/fibers using the WPI's Biofluorometer (SI-BF-100) simultaneously.





µManager software interface and Optisarc image processing plug-in for sarcomere spacing detection run simultaneously during live image capturing of contracting cardio myocyte on selected ROI (Image with courtesy of Dr. Michael Kohlhass, Universitätsklinikum des Saarlandes, Homburg).

The Optical Sarcomere Spacing System can also be extended for capturing sarcomere length changes in slow-twitch and fast-twitch skeletal muscle cells/fibers.

Sarcomere spacing detection of any muscle slice is also possible using genuine image processing methods for high reproducible measurements of selected sarcomere regions in static images.

OPTISARC CAMERA SPECIFICATIONS

IMAGE SENSOR1/2.5" 5 megapixel monochrome CMOSACTIVE PICTURE ELEMENTS2,592 (H) x 1,944 (V).MAXIMUM FRAME RATE2,106 FPS @ 32 x 32.OPTIMAL FRAME RATE500 FPS @ 296 X 148CONNECTIONUSB 3.0 Micro-B

ORDERING INFORMATION

SI-OSARC	USB 3.0 CMOS Camera with OptiSarc software plug-in for live image capture to detect sarcomere spacing
INV-101	Inverted Microscope with C-mount Adapter
SI-BF-100**	Biofluorometer
SI-CTS200**	Cell Tester

µManager* was developed at Vale laboratory at UCSF, funded by an NIH grant R01-EB007187 from the National Institute of Biomedical Imaging and Bioengineering (NBIB). https://micro-manager.org/ and https://imagej.nih.gov/ij/

*Arthur Edelstein, Nenad Amodaj, Karl Hoover, Ron Vale and Nico Stuurman (2010): Computer Control of Microscopes using µManager. Current Protocols in Molecular Biology, Chapter 14, Unit 14.20, 22 pages.

** For simultaneously use with SI-BF-100 and SI-CTS200, ask WPI about selecting an appropriate microscope.

Ca²⁺ Detection in Muscle Tissue using Fluorescence Spectroscopy

The use of fluorescent probes in cell physiology has emerged as indispensable tool in the analysis of cell functioning over recent years. The physics underlying fluorescence is illustrated by the electronic-state diagram (so-called Jablonski diagram, see Fig. 1), showing the three-stage process to create the fluorescent signal (Excitation - Excited/State Lifetime - Fluorescence Emission) in a fluorophore/indicator and simplified described below.



Fig. 1– Jablonski diagram illustrating the processes of fluorescence by absorption of higher photon energy by a fluorophore and subsequent emission of lower photon energy, resulting in fluorescence during the

Fluorescence is obtained when an excitation photon (hv_{EX}) from an external source, such as a high-power LED, is absorbed by a fluorophore that elevates its energy (S1'). During the fluorescence-lifetime, the elevated energy (S1') decays to a lower energy state S1. Then, fluorescence results in the emission of a photon with lower energy (hv_{EM}) and at a lower wavelength. Fundamental in spectroscopy is the difference in energy or wavelength represented by (hv_{EX} - hv_{EM}), which is called the Stokes shift. The Stokes shift allows efficient discrimination

of the excitation, making fluorescence a very sensitive technique and able to be detected against a low background, isolated from excitation photons.

fluorescence-lifetime.

Fig. 2– Typical excitation and emission spectra, showing the lower wavelength of the excitation source and higher wavelength of the fluorescence.

wavelength

Four essential elements of fluorescence signaling can be then identified to build up a detection system:

- Excitation light source adapted to the absorption bandwidth of the fluorophore (e.g. high-power LED of specific wavelength)
- A fluorophore/indicator (e.g. Fura-8 for free Ca²⁺ detection in muscle tissue)
- Emission wavelength filters to limit the bandwidth of the emission photons or overlapping bands
- A detector system that registers the fluorescence light and produces a recordable output as an electrical signal (e.g. Photomultiplier tubes)

Regardless of the application, compatibility of these four elements is essential for optimizing fluorescence detection.

Example of free Ca2+ detection in muscle tissue

Typically, a fluorescent dye is introduced into tissue or single cells to obtain a fluorescent response of the labeled molecule. A typical example is the detection of the transient increase in the cytoplasmic/myoplasmic

free calcium concentration (Δ [Ca²⁺]) as the intermediate signaling event of the excitation-contraction coupling. The quantification of Δ [Ca²⁺] is done using a monochromatic light to excite the dye labeled Ca²⁺ molecule in a tissue/cell sample either in a tissue bath or microscopic experimental set-up. The emitted fluorescence signal from the indicator dye can be then used to monitor the amplitude and time-course of

the $\Delta[Ca^{2+}]$ detected by sensitive detectors, such as highly sensitive

The ratiometric indicator dye Fura-8 was selected for the detection of free calcium concentration ($\Delta [\text{Ca}^{2+}]$) in heart muscle tissue slices. Fura-8 was excited at 365nm and 410 nm wavelengths and the emission recorded at 535 nm wavelength in dual excitation/ single emission mode. The advantages of choosing this ratiometric measurement technique with dual excitation/single emission using Fura-8 were minimization of movement artifact, cancellation of possible effects of uneven loading, inhomogeneous distribution of fluorescence indicator in the cells or indicator bleaching in the detection of free calcium concentration ($\Delta [\text{Ca}^{2+}]$) in the muscle tissue.

This allowed quantification and comparison between:

photomultiplier tubes (PMT module) or cameras.

- High spatial versus high time resolution techniques on the human left ventricular slices
- The possibility to measure free calcium concentration (Δ[Ca²⁺])
 transients in a horizontal tissue bath on human left ventricular slices
 or murine slices.




Fig. 3– Qualitative representation some results of free Ca2+ detection in human heart slices using the SI-BF-100 system.

Average fluorescence intensities of Fura-8 loaded human left ventricular slices detected at 525 nm, when excited at 340 nm and 410 nm, respectively, and ratios calculated (lower trace) from the imaging data of a Rolera EM-C2 camera (left). Right, the response of the SI-BF-100 detected at 525 nm, when excited at 365 nm and 410 nm wavelength using two aperture settings and calculated ratios (lower traces). Furthermore, fluorescent data collected with the small aperture setting and low-pass filtered at 50 Hz is shown (from: Belz et al., Proc. SPIE 9702, 2016).

References

Belz M., et. al. Fiber optic Biofluorometer for physiological research on muscle slices. Proc. SPIE 9702, Optical Fibers and Sensors for Medical Diagnostics and Treatment Applications XVI, 2016.

Spectrophotometry. Wikipedia, the free encyclopedia, 2017 (Cross-references).

19

Biofluorometer

Reliable, 2-channel, LED-based fluorometer

Features

- Two photomultiplier inputs
- · Light excitation with high power LEDs
- Modes: single excitation & single emission, single excitation & dual emission, dual excitation & single emission
- 2 channel mode (two single excitation & single emission)
- Optical connections: Liquid Light Guides (LLGs) and SMA terminated fibers
- Sampling rates up to 1 kHz (1000 ratios/second)
- Automatic LED light drift correction for long term measurements
- · Automatic room light correction
- Optional fiber optic probes for horizontal tissue bath applications
- · Optional imaging probes for Langendorff systems
- Optional attachments for direct connection to fluorescence microscopes via epifluorescence port (excitation) and C-Mount (Emission) via liquid light guides

Benefits

- Versatile instrument for horizontal tissue bath, Langendorff and microscope applications
- Warmup time < 1 minute
- Low bleaching mode (5% LED On time)
- Can be combined with imaging based Sarcomere Detection System (OptiSarc)
- SMA and Liquid Light Guide connections
- · Single and dual emission or excitation detection methods
- 2-channel instrument for single excitation/single emission dyes
- Customized analysis techniques in WPI's MDAC data acquisition software

Applications

- Ratiometric calcium measurement and ATPase
- Fluorometric applications in neuroscience and cell biology
- The SI-BF100 is an LED-based fluorometer for life science applications.
 It is ideally suited for ratiometric calcium detection (FURA-8) and ATPase detection (via NADH fluorescence). With up to three LED modules (wavelengths), the SI-BF100 covers many fluorometric applications in neuroscience and cell biology.

The **SI-BF-100** enables the detection and analysis of fluorescence signals in four different modes:

- Single excitation/single emission—In this classical mode, a fluorophore
 is excited at one wavelength and the fluorescence signal is detected
 at a single higher wavelength using one photomultiplier. The
 concentration of the analyte is directly proportional to the intensity of
 the detected signal.
- Dual excitation/single emission–A fluorophore is excited at two
 wavelengths and the fluorescence signal is detected at one
 wavelength using one photomultiplier. The concentration of the
 analyte is proportional to the ratio of the two detected fluorescence
 signals. This ratiometric concept minimizes the effect of indicator dye
 bleaching and motion artifact in experiments. A typical example is
 the detection of free calcium in muscle tissue using the indicator dye
 Fura-8™.
- Single excitation/dual emission–A fluorophore is excited at one wavelength and the fluorescence signal is detected at two wavelength using two photomultipliers.
- Dual excitation/dual emission–Two separate fluorophores are excited at different wavelengths and the fluorescence signal of each fluorophore is detected at two separate wavelengths using two photomultipliers.




Experimental Setups

Possible experimental setups for the **SI-BF-100** include:

A: This **SI-BF-100LLG** setup depicts the attachment to a fluorescent microscope.

B: This **SI-BF-100LLG** setup shows a single fiber based detection system with a dichroic mirror and an imaging system at the distal end of a liquid light guide. Dichroic mirrors are used in A and B to separate excitation and emission light.

C: The SI-BF-100SMA setup has a round shaped fiber bundle that is used to deliver excitation light to and pick up fluorescent light from the sample (adapted from Belz et al., 2016).




Single and dual emission probes

All probes use fibers with a 300 μm core diameter. Excitation fibers have 1000 μm SMA connectors for Excitation and Double Emission Probes and 1500 μm SMA connectors for Single Emission Probes.




This SMA Single Emission Probe has a for muscle strips, where the middle fibers are used for excitation and the rectangular fiber pattern, optimized surrounding fibers are used for detection.



Biofluorescence of calcium use in heart muscle slices with Fura-8

Calcium sensing

A specific target molecule in muscle contraction is $Ca^{2^{+}}$ as key intermediate signaling event between excitation and contraction of muscle fibers, and thus essential for the analysis of the force development in muscles. On the cellular level, force production is therefore directly related to the transient increase in the myoplasmic free calcium concentration. More specifically, the assessment of both parameters simultaneously is therefore critical in the evaluation and interpretation of the force development characteristics. Fluorescence techniques used in conjunction with muscle research systems (like WPl's SI-MT-L or SI-HTB2) to record muscle force is an **innovative technique** in cardiac muscle and skeletal muscle physiology. WPl's Biofluorometer (SI-BF-100) was specifically developed to monitor rapid changes in Ca^{2+} transients, i.e. $\Delta[Ca^{2+}]$.



Organ bath setup: Average fluorescence intensities of Fura-8 loaded murine myocardium slices excited at 365/410 nm wavelength and detected at 525 nm. Left, using two apertures for fluorescence detection via the microscope. Right, same experimental condition using optical fiber fluorescence detection with the SI-BF-100SMA in the organ bath, sampled at 1 kHz. The organ bath data are low-pass filtered at 200 Hz or 50 Hz cut-off frequency. Adapted and reprinted with permission from Belz et al., Proc of SPIE Vol 9702, 97020Q-1-97020Q-11, 2016.

Advantages using Fura-8 dye:

- More sensitive to calcium than Fura-2
- Higher signal-to-noise ratio than Fura-2
- Simplified dye loading by incubation for 1 hour at room temperature
- Emission peak response is shifted to longer wavelength (peak at 525 nm)
- Red-shift dual excitation wavelength (354 nm and 415 nm)

Application in heart muscle slices

The Biofluorometer opens a wide field in functional fluorescence research, by studying the fundamental and/or applied aspects of the underlying energetics and signaling aspects of muscle contraction. This is notably useful in:

- Pre-clinical & Toxicological Studies:
 - screening of potential drugs
 - evaluating the side effects of drugs
 - evaluating models of cardiac disease
- Sports & Rehabilitation:
 - disuse vs. overuse
 - · muscle damage
 - function for heart transplantation



The SI-BF-100 Organ/tissue bath setup with an SMA single emission probe for direct sensing of calcium fluorescence signal. Photo courtesy of Professor Andreas Dendorfer, Walter Brendel Zentrum, München, Germany.

System configuration

The **SI-BF-100** setup for detecting Ca²⁺ transients via Fura-8 depends on the excitation and emission wavelengths of the indicator dye Fura-8. The **SI-BF-100** setup depends on your experimental paradigm and the sample size of interest:

- BF-100SMA-C-Fura 8: direct sensing of calcium via direct measurement in tissue sample
- BF-1002XLG-C-Fura 8: direct sensing of calcium via microscope setup in muscle cell or in small tissue sample
- BF-100CAM-C-Fura 8: calcium imaging via microscope setup with camera in muscle cell or in small tissue sample

BF-100SMA-C-Fura 8	Includes
SI-BF-100 Main Unit with SMA connection	(1) SI-BF-100SMA
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Single Emission Probe	(1) 94650

BF-100CAM-C-Fura 8	Includes
SI-BF-100 Main Unit with LLG connection	(1) SI-BF-100LLG
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Liquid Light Guide Ø 3mm	(3) 802407
CM-CAM-2XLG Packaging	(1) 99259

BF-1002XLG-C-Fura 8	Includes
SI-BF-100 Main Unit with LLG connection	(1) SI-BF-100LLG
HP LED 365 nm	(1) 99209-1
HP LED 420 nm	(1) 99209-4
Optical Filter 535 nm, 43 nm BP	(1) 802238
Liquide Light Guide Ø 3mm	(2) 802407
CM-2XLG Packaging	(1) 99261

SI-BF-100 SPECIFICATIONS

MEASUREMENT PRINCIPLE	Fiber optic fluorometer with 2 inputs and 1 output
DETECTOR INPUTS	2 PMTs
EXCITATION	High Power LED Modules: 365 nm, 420 nm, 470 nm, 530 nm (select any 3 modules, when ordering)
ANALOG OUTPUT RANGE	0-10V
OPTICAL CONNECTIONS	Choice of Liquid Light Guide (LLG) or SMA connections
POWER	12 V/2 A (includes external 100 – 240 V / 50 – 60 Hz power supply)
DIMENSIONS (h x w x d)	3.5 x 17 x 13 in. (88 x 431 x 330 mm)

		ATION

SI-BF-100LLG	Biofluorometer with LLG Optical Connections
SI-BF-100SMA	Biofluorometer with SMA Fiber Optic Connections
99261	C-Mount Microscope Attachment for 2x PMTs
99259	C-Mount microscope attachment for 1x camera & 2x PMTs
Includes	1x camera C-Mount adapter with adjustable aperture

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

SI-BF-SMA-UPGRADE	Biofluorometer Upgrade Kit for SMA optical probes
SI-BF-LLG-UPGRADE Biofluorometer Upgrade Kit for LLG connection	
802407 Liquid Light Guide (LLG), 3mm diameter, 6 ft. long	
M3301 Manual Manipulator for securing the probe	
M10 Magnetic Base	
94650 Single Emission, Small Tissue Probe	
94689 Dual Emission, Small Tissue Probe	
CL 1 1	

Check our website for new LED modules, emission filters and dichroic mirrors for specific applications.

KG Optical Force Transducers

- Simple calibration
- Different models to accommodate a wide range of forces and sensitivities
- Nearly insensitive to changes in temperature/ambient light
- Extremely high level of linearity
- Virtually indestructible with normal use
- KG transducers are required for use with SI-BAM21 amplifiers



	UNLOADED TRANSDUCER WITHOUT TISSUE MOUNTING SUPPORT				IT
	Force Range Noise Compliance Frequency				
SI-KG7	0-5 mN	0-0.5 g	0.2 μΝ	10 nm/mN	250 Hz
SI-KG7a	0-5 mN	0-0.5 g	0.4 μΝ	5 nm/mN	500 Hz
SI-KG7b	0-10 mN	0-1 g	1 μΝ	1.5 nm/mN	550 Hz
SI-KG2	0-2 N	0-200 g	250 μΝ	150 nm/mN	1.3 kHz
SI-KG2B	0-0.2 N	0-20 g	80 μΝ		590 Hz
SI-KG4	0-50 mN	0-5 g	15 μΝ	0.5 nm/mN	1.2 kHz
SI-KGxx	Contact	WPI for spe	ecialty tran	sducers with dif	ferent ranges.

TRANSDUCER SPECIFICATIONS

Tissue Mounting Hooks

Mounting hooks can be used in a variety of combinations, depending on the type of tissue to be examined.

Mounting hooks are sold in kits. Currently, there are 11 kit configurations, each available in four different sizes. The mounting hook size that is required depends on the force transducer used.

Vascular hooks are available for mounting blood vessels (rings). They are normally used with a pair of blunt hooks (SI-TM8).




For larger muscles, screw clamps (SI-TM11) and spring clips (SI-TM9, SI-TM10) are available.

The micrometer and motor receive a large (SI-KG4 size) tissue mount. If a smaller tissue mount is used, the 97909 adapter is required. This adapter is included with every SI-MT or SI-HTB system.

Ordering


When ordering tissue mounts, specify the tissue mount configuration and force transducer to be used.



Force Transducer mount pictured on the left; Motor/micrometer mount pictured on the right.

Tweezer Basket Pointed Hook Blunt Hook Spring Clip Screw Clip Tendon Hook





*The **97909** tissue mount adapter tube (OD:0.096", ID:0.035") allows you to use **SI-KG4** size mounting hooks with **SI-KG2** size force transducers.

ORDERING INFORMATION						
	Force Transducer Mount Micrometer/Motor Mount Force Transducers					
SI-TM1 Papillary Muscle	Basket	Pointed Hook	Available for all force transducers			
SI-TM2 General Purpose	Pointed Hook	Pointed Hook	Available for all force transducers			
SI-TM3 Small Skeletal Muscle	Tweezer	Tendon Hook	Available for all force transducers			
SI-TM4 Trabeculae	Basket	Basket	Available for all force transducers			
SI-TM5 Papillary Muscle	Basket	Tweezer	Available for all force transducers			
SI-TM6 Muscle Rings	Blunt Hook	Blunt Hook	Available for all force transducers			
SI-TM7 General Purpose	Tweezer	Tweezer	Available for all force transducers			
SI-TM8 Muscle Rings	Blunt Hook/Vascular Hook	Blunt Hook/Vascular Hook	Available for all force transducers			
SI-TM9 Strong Skeletal Muscle	Tweezer	Spring Clip	SI-KG2, SI-KG2A, SI-KG2B Only			
SI-TM10 Strong Skeletal Muscle	Pointed Hook	Spring Clip	SI-KG2, SI-KG2A, SI-KG2B Only			
SI-TM11 Very Strong Skeletal Muscle	Pointed Hook	Screw Clamp	SI-KG2, SI-KG2A, SI-KG2B Only			

4-Channel Data Acquisition System with Software

Low noise, high resolution system with 8 analog input and 3 analog output channels

Features

- Powerful low-noise (<1 mV RMS) and high-resolution (16 bits) data acquisition system for sampling up-to 8 analog input channels and 3 analog output channels simultaneously, using standard BNC connections
- MDAC software provides easy to use interface controlling, with extensible standard and customized Data Processing and Analysis Tools



Benefits

- Online Channel Math operations, general purpose Fast Fourier analysis (FFT) and digital filtering of Analog In channels
- Numerous basic signal forms can be combined to design experimental protocols, for most physiological applications
 - · Factory designed standard or customized protocols
 - · Semi-automated data analysis toolbox
- Protocol repeat function to avoid time consuming protocol programming of extended experiments

Applications

- Muscle physiology (Can be used with SI-MTM Muscle Testing Platform, SI-CTS200 Cell Tester System, SI-HTB2 Horizontal Tissue Bath and SI-BF-100 Biofluorometer)
- Stand alone general data recorder for Spectroscopy, Neuroscience and Electrophysiology (Can be used with TBR4100 Free Radical Analyzer, Extracellular Bioamplifiers like SYS-DAM50, SYS-DAM80, SYS-900A, ISO-80, EVOM2™ Volt Ohm Meter, ATC2000 Animal Temperature Controller, BP-1 Blood Pressure Monitor or the BAT-12 Microprobe Thermometer)
- Instrument control for software triggered devices like A365/A385/ A395 Constant Current Stimulators, MPS-2 Perfusion System, SYS-PV820/SYS-PV830 Pneumatic PicoPumps, Duo 773 Intracellular amplifiers, and the SYS-TBM4M Transbridge Transducer Amplifier (e.g. for FORT force transducers)

Knowledge of the physiological characteristics of muscle tissue can be useful to quantify beneficial or adverse effects of drug supply on muscle function in pre-clinical and toxicological studies, evaluating muscle dystrophies, training effects in sports and rehabilitation (disuse vs. overuse) and advanced physiology and biomedical research.

This is usually achieved by quantifying the contractile and/or the elastic properties of muscle tissue. This needs the programming of different and specific experimental protocols (isometric, concentric and eccentric, isokinetic or isotonic), so that the physiological structure of interest can be quantified. **LabTrax-MDAC** data acquisition software was designed for use with WPl's Muscle Physiology line to test physiological characteristics of muscle tissues in various conditions, using factory designed standard or customized protocols. The semi-automated Data Analysis Toolbox of standard protocols gives quick access to user-friendly, readable and interpretable results of the experiments.

Variety of muscle physiology applications

The physiological response of muscle tissue to training, disuse, nutrition, drug supply and others factors may be studied by adding accessories to the system, like:

 Study of the muscle's force production capacities in combination with the Ca²⁺ release from the sarcoplasmic reticulum (SR) and ATPase consumption. The perfect instrument for this is WPI's Biofluorometer (SI-BF-100) in combination with any system of WPI's Muscle Physiology line, controlled via LABTRAX-MDAC. Study of the muscle's force production capacities from direct muscle or peripheral nerve stimulation. For this experiment, use WPI's programmable isolated current stimulators (A365, A380 or A395), controlled via LABTRAX-MDAC.

LABTRAX-MDAC provides easily used continuous stimulation protocols, so that especially cardiac cells/tissue remain intact during experimental resting periods.

LABTRAX-MDAC is also well suited for other software triggered instruments or as a stand-alone general data recorder for selected WPI Instruments



The back panel of the Lab-Trax-8/16 has four analog outputs, digital inputs or outputs, a USB port, power socket and power switch.

LAB-TRAX-8/16 SPECIFICATIONS

ANALOG INPUTS 8 BNC connections INPUT RANGE + 10VSYSTEM NOISE < 1 mV RMS 1.500V **ISOLATION OPERATING CURRENT** 800 mA maximum ANALOG OUTPUTS 4 BNC connections OUTPUT RANGE + 10V IMPLEMENTED FILTER 5th order low-pass Bessel filter with 3dB cut-off frequency **OUTPUT IMPEDANCE** 100Ω **OUTPUT CURRENT** 15 mA DIGITAL I/O 16/16 TTL (BNC or DB-9 Connector) LOGIC HIGH VOLTAGE 3.3 V minimum LOGIC LOW VOLTAGE 1.0 V maximum ANALOG & DIGITAL INPUTS Operating voltage protected to ±30V PC INTERFACE USB 2.0 **RESOLUTION** 16 bits **POWER SOURCE** 12V DC

ORDERING INFORMATION		
LABTRAX-MDAC	Lab-Trax-8/16 with MDAC software	
LAB-TRAX-4	4-Channel General Data Acquisition System	
2851	BNC to BNC Cable	

Micromanipulators


Micromanipulators are used when precision work is conducted under a microscope. A micropipette, electrode or probe can be mounted on a micromanipulator and move as little as a micron at a time. This tool can be used for in vitro fertilization, patch clamp experimentation, extracellular recording, microinjection and any application requiring fine mechanical placement (resolution).

In addition to micromanipulators, WPI offers tilt bases, piezo translators and a variety of stands.

Micromanipulators can be broken out into three broad categories

- Manual/Motorized
- Motorized





MICROMANIPULATOR COMPARISON							
Micromanipulator	Manual or Motorized	Resolution	Travel	Stands	Tilt Base	Piezo Translator	Notes
SU-TRIO235	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm Diagonal: 50 mm	Fixed platform stage with imperial/ standard holes, chamber insert and gantry supports	N/A	N/A	Carries up to a kilogram
SU-TRIO245	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm	Fixed platform stage with imperial/ standard holes, chamber insert and gantry supports	N/A	N/A	Carries up to a kilogram
SU-QUAD	Motorized	< 100 nm	X: 25 mm Y: 25 mm Z: 25 mm Diagonal: 30 mm	Fixed platform stage with imperial/ standard holes, chamber insert and gantry supports	N/A	N/A	Quiet mode eliminates electrical noise
M3301	Manual	0.01 mm (X fine) 0.1 mm (X, Y, Z)	X (fine): 10 mm X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3	MPM20	
KITE	Manual	0.1 mm	X (fine): 10 mm X: 35 mm Y,Z: 20 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
M325	Manual	10 μm	X: 25 mm Y,Z: 10 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		
ммј	Manual	0.1 mm	X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Joystick control
MD4	Manual	10 μm (X fine) 100 μm (X, Y, Z)	X (fine): 10 mm X: 37 mm Y: 20 mm Z: 25 mm	M9, M10, M10L, 501622, 501623	TBS, M-3		Holds two electrodes
MM1	Manual (mini)	1.0 µm	3 mm				225 g load
ММЗ	Manual (mini)	1.5 µm	13 mm				340 g load
MM1-3	Manual (mini)	1.0 µm	3 mm	MB-2			225 g load
ММ3-3	Manual (mini)	1.5 µm	13 mm	MB-2 with MM1A (adapter)			340 g load

Manual Micromanipulators

Popular manual micromanipulator

Economy manual micromanipulator



Features

- The most widely used micromanipulator
- · Lightweight 550 g
- Sure, repeatable movement without drift

Benefits

- Control knobs clustered in 8 cm area in a single vertical plane for quick
- Right and left hand orientation options available

Applications

- Microinjection
- · Electrophysiology recording

Weighing just 550 g and employing a slim space-saving design, this wellbuilt micromanipulator outsells all others worldwide for high precision experiments where magnification is in the range of up to 250×. Its design allows units to stand tightly grouped, since all control knobs project to the rear. And because control knobs are clustered within an 8 cm area in a single vertical plane, resolution is quick. The hand works blindly while the eye monitors the microscopic image. Vernier scales allow readings to 0.1 mm; X-axis fine control allows readings to 10 µm.

The instrument employs rack-and-pinion drive, V-shaped guide ways and cross roller bearings, so all movement is sure and repeatable, without drift, side play, backlash or sticking. Contact parts are milled of hardened steel for high performance and long life.

M3301 SPECIFICATIONS		
	TRAVEL RANGE	RESOLUTION
X-axis Fine	10 mm	0.01 mm
X-axis	37 mm	0.1 mm
Y-axis	20 mm	0.1 mm
Z-axis	25 mm	0.1 mm

	ORDERING INFORMATION
M3301R	Manual Manipulator (right-handed)
M3301L	Manual Manipulator (left-handed)
M3301-M3	G-R Manual Manipulator (right-handed) & Tilting Base*
M3301-M3	I-L Manual Manipulator (left-handed) & Tilting Base*
502105	Axis Adiustment Tool

*Requires 5464 weight



- Vernier scales allow readings to 0.1 mm
- X-axis fine control allows readings to 10 µm
- Choice of optional M3 Tilting base which can be mounted to a table with M6 screws
- Left or right-handed versions of the **KITE** micromanipulator are supplied with a standard 12 mm clamp and electrode holder M3301EH

Benefits

- Control knobs clustered in the same plane for quick resolution
- · Right and left hand orientation options available

Applications

- Microinjection
- · Electrophysiology recording

The KITE is an economical manual micromanipulator that is perfect for student use. It is available in right or left hand versions and is also sold bundled with an optional M3 Tilting Base. Our 5 lb. weight (WPI #5464) is frequently sold with this unit to provide stability.

KITE SPECIFICATIONS					
	TRAVEL RANGE RESOLUTION				
X-axis Fine	10 mm	0.01 mm			
X-axis	35 mm	0.1 mm			
Y-axis	20 mm	0.1 mm			
Z-axis	20 mm	0.1 mm			

	ORDERING INFORMATION
KITE-R Kite Manual Manipulator (right-handed) KITE-L Kite Manual Manipulator (left-handed)	
KITE-M3-L Kite (left-handed) + Tilting Base Combo*	
	*Requires 5464 weight

OPTIONAL ACCESSORIES

M3301EH	Replacement Electrode Holder (14 cm long)
15873	Optional Angled Electrode Holder (13 cm long)
M-3	80° Tilting Base M6 x 1 mm screw
5464	5-lb Weight for Tilting Base (Shipping weight: 7 lb (3 kg))
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder
M4C	Microscope Stage Adapter
	Also see magnetic stands (page 20)

Also see magnetic stands (page 28)

Miniature Micropositioners

One and three axis positioners for small jobs





Features

- Precise and smooth motion
- Less than 1 µm maximum wobble

Benefits

- Compact
- Provides precise and smooth motion with no backlash

Applications

· Oocyte injection with Nanoliter Injector

Single stage measures only $5\times11\times26$ mm with 3 mm travel. Provides precise and smooth motion with no backlash, positive spring loaded carriage, straight within 1 μ m and less than 1 μ m maximum wobble. Features fine 80 TPI screw adjustment. 10 mm square mounting surface has a 3.9 mm tapped center hole for transmission and/or mounting. Available in single X (**MM1**) and X-Y-Z (**MM1-3**) axis configurations.

MM3 and MM3-3

Features

- Opens 0-4mm
- · Precise and smooth motion
- Less than 1 μm maximum wobble

Benefits

- Compact
- Provides precise and smooth motion with no backlash

Applications

• Oocyte injection with Nanoliter Injector



Single stage measures only $7\times17\times44$ mm with 13 mm travel. Offers precise and smooth motion with no backlash, positive spring-loaded carriage, straight within 1.5 µm, and less than 1.5 µm maximum wobble. Features fine 80 TPI screw adjustment. 13 mm square mounting surface has a 7mm tapped center hole for transmission and/or mounting. Available in single X (MM3) and X-Y-Z (MM3-3) axis configurations.

MM1

MINI-MICROPOSITIONER SPECIFICATIONS				
	MM1	MM1-3	MM3	MM3-3
AXIS	Χ	X-Y-Z	Х	X-Y-Z
STRAIGHT LINE ACCURACY	Within 1 µm over 3 mm travel	Within 1 µm over 3 mm travel	Within 1.5 µm over 13 mm travel	Within 1.5 µm over 13 mm travel
CLEAR APERTURE	3.9 mm tapped hole, 8-32 thread	3.9 mm tapped hole, 8-32 thread	7 mm tapped hole, 5/16-16 thread	7 mm tapped hole, 5/16-16 thread
LOAD CAPACITY	255 g Normal	255 g Normal	340 g Normal	340 g Normal
FINISH	Black Anodized	Black Anodized	Black Anodized	Black Anodized
WEIGHT	3 g/axis	12 g/axis	14 g/axis	48 g/axis
TYPE	Fine Screw	Fine Screw	Fine Screw	Fine Screw
TRAVEL	3 mm	3 mm	13 mm	13 mm

	ORDERING INFORMATION
MM1	Mini Micropositioner, one axis, 3 mm travel
MM1-3	Mini Micropositioner, three axes, 3 mm travel
MM1-A	Mounting Adapter for MM1 and MM1-3
MM1-C	Clamp for MM1 and MM1-3
ММ3	Micropositioner, one axis, 13 mm travel
MM3-3	Micropositioner, three axes, 13 mm travel
ММ3-А	Mounting Adapter for MM3 and MM3-3
ММ3-С	Clamp for MM3 and MM3-3
MM3-ALL	Complete 3-Axis Micropositioner & Magnetic Stand
MM1-ALL	Complete 3-Axis Mini Micropositioner & Magnetic Stand



Highly stable, motorized, 3-axis manipulator with 25 mm of movement on each axis





Features

- Three independent axes 25 mm orthogonal travel in X, Y and Z
- Sub-micron (less than 100 nm) resolution
- Carries up to a kilogram
- DIP switches on the Control box select direction of movement produced by turn of the Control box knob
- USB interface for computer control
- The TRIO™ comes standard with a universal mounting system suitable for the most popular headstages or pipette holders.

Benefits

- User selectable angle from 0 90° via Control box input
- Fast movement with a top speed of 3mm/sec (while homing)
- Mechanically robust construction for high stability
- Compact, fanless, user-friendly Control box controller preserves bench and rack space
- Push button control of multiple functions WORK, HOME, LOCK, PULSE and RELATIVE
- Mechanical designs and software algorithms allow the motors to power down completely during recording

Applications

- Suited for in vivo and in vitro electrophysiological recording
- · Applications with noise-sensitive recordings

The new **SU-TRIO245** is a highly stable 3-axis manipulator with 25 mm of travel on each axis. The TRIO's synthetic 4th axis can be set in software as any angle between 0 and 90° for diagonal movement. Based on a lead-screw design with a smaller overall size and footprint than most manipulators, the SU-TRIO245 is ideal for applications that require 2 pipettes in one setup or for setups where space is limited.

The compact design of the integrated control box requires minimal bench space. It provides quiet, fan-free operation, and is easy to use. No rack mounted controller is required. Position coordinates, in relative or absolute values, are displayed directly on the control box. The TRIO manipulators use a logarithmic acceleration algorithm that eliminates the need for speed selection. As the knobs on the control box are turned faster, acceleration ramps up. This allows for smooth and intuitive motion control of electrode position without the need to stop and change speeds or lift your hand from the knobs. A Y-axis lockout function (accessible by DIP switch) is also available, allowing X/Z-only axial movement during HOME and WORK repositioning.

An alternate configuration, the SU-TRIO235, removes the Z axis and replaces it with an adjustable diagonal axis. This configuration pairs a traditional X and Y axis with a diagonal axis that can be adjusted in angle with a set screw. The diagonal features 50 mm of travel with X and Y axes that have 25 mm of travel. On this model, we have added software to create a synthetic "Z-axis" by combining the diagonal and X axes in reverse of how we create a synthetic diagonal axis on our other manipulators. The SU-TRIO235 and SU-TRIO245 both feature the same precision bearing and lead-screw design.

TRIO SPECIFICATIONS

TRAVEL SU-TRIO235: 25 mm on X and Y, 50 mm on diagonal

SU-TRIO245: 25 mm on X, Y and Z axes

CONTROL BOX 5.5 x 7.5 x 4" (14 x 19 x 10.2 cm)

WEIGHT 2.2 lbs. (1.0 kg)
POWER 115/230 V, 50/60 Hz

ORDERING INFORMATION

SU-TRIO235	3-Axis Motorized Micromanipulator with X, Y, D Control
SU-TRIO245-L	3-Axis Motorized Micromanipulator with X, Y, Z Control (left-hand)
SU-TRIO245-R	3-Axis Motorized Micromanipulator with X, Y, Z Control (right-hand)
505069	Z-Axis Vertical Extension for Quad/Trio
505071	Rotating Base for MP225 Mechanical
505073	4-in. Dovetail Extension
505074	Mounting Adapter Plate
505085	Rod Holder for Motorized Manipulator

Basic systems include the manipulator, controller, rod holder, 4-inch dovetail extension, mounting adapter plate, Z-axis vertical extension, cables, and power supply

QUAD® 4-Axis Motorized Manipulator

4-axis motorized manipulator with 25 mm of movement on each axis



Features

- Four independent axes 30 mm travel in diagonal for coaxial pipette movement, 25 mm travel in X, Y and Z
- Sub-micron 100 nm resolution
- · Quiet mode eliminates electrical noise
- True diagonal assures coaxial movement
- USB interface for computer control

Benefits

- Compact, fanless, user-friendly control box preserves bench and rack space
- Push button control of multiple functions WORK, HOME, LOCK, PULSE and RELATIVE
- Robotic HOME and WORK position moves for easy automated pipette exchange
- · Display indicates coordinates in relative or absolute
- Five conveniently located control buttons
- When you are ready to begin collecting data, the motor drive electronics can be suppressed

Applications

- Suited for in vivo and in vitro electrophysiological recording
- Cell penetration

TRAVEL

CONTROL BOX WEIGHT POWER

• Applications with noise-sensitive recordings

The QUAD® motorized micromanipulator is easy to use and features four independent axes. Each axis has a 25 mm range of motion, a digital display of position and a control box. The compact, intuitive controller takes up minimal bench space, is fan free and is easy to use. Three axes provide the X, Y, Z-orthogonal motion typical of most motorized micromanipulators. In the QUAD® a true fourth axis moves the electrode coaxially at exactly the desired angle of approach. The fourth axis has 30 mm of travel that significantly extends the range of travel for the system.


The QUAD® is very easy to operate. It includes a control box with digital display of position. The control box inputs for each axis allow facile manual control of electrode position. Five conveniently located buttons control all of the basic functions that you will need in normal operation. Press and hold the WORK button to quickly store a work position; press WORK again and the manipulator will return to the same location. HOME sends the manipulator to an initial location or any user defined position, which is useful for changing electrodes rapidly. When you are ready to begin collecting data the motor drive electronics can be suppressed by pressing LOCK. Display coordinates can be zeroed at any location by pressing RELATIVE; go back to absolute coordinates by pressing RELATIVE again. Finally, PULSE activates a pulse movement mode that produces small, rapid bursts of motion that can be advantageous for sharp electrode cell penetration.

QUAD SPECIFICATIONS
30 mm on the diagonal 4th axis
25 mm on the X, Y and Z axes
5.5 x 7.5 x 4 in. (14 x 19 x 10.2 cm)
2.2 lb. (1.0 kg)
115/220 // 50/60 //-

	ORDERING INFORMATION
SU-QUAD-L	4-Axis Motorized Micromanipulator (left-hand)
SU-QUAD-R	4-Axis Motorized Micromanipulator (right-hand)
505069	Z-Axis Vertical Extension for Quad/Trio
505070	Mounting Adapter Plate for Quad
505071	Rotating Base for MP225 Mechanical
505073	4-inch Dovetail Extension
505085	Rod Holder for Motorized Manipulator

Specialty Micromanipulators

Joystick Controlled Manual Micromanipulator



Features

- X axis can be tilted 90°
- Easy steering motion that translates normal hand movement into smooth sub-millimeter shifts

Benefits

- Fine adjustment for the X and Y axes can be controlled by the joystick
- Joystick manipulation

Applications

Oocyte injection with Nanoliter Injector

	ORDERING INFORMATION
MMJR	Joystick Micromanipulator (right-handed)
MMJL	Joystick Micromanipulator (left-handed)
	Options and replacement parts listed below

Micrometer Slide Micromanipulator



Features

- · Built of precision micrometer-actuated linear slides
- Can be configured for right or left hand use
- Resolution is 10 µm

Benefits

- Each slide is comprised of a large micrometer head and a springreturn linear slide
- Micrometer head is graduated in 10 µm steps which enables repeatable positioning to an accuracy of 2 µm

Applications

• Electrophysiology recording and injection

ORDERING INFORMATION

M325 3-Axis Fine Controlled Manual Micromanipulator

Options and replacement parts listed below

Dual Tool-Holder Micromanipulator



Features

- Scales allow coarse adjustment readings with an accuracy of 100 μm
- X-axis fine control is achieved with a micrometer screw

Benefits

• Dual electrode holders

Applications

• Electrophysiology differential amplification

	ORDERING INFORMATION
MD4R	Double-Holder Micromanipulator (right)
MD4L	Double-Holder Micromanipulator (left)
MD4-M3-R	Double-Holder Micromanipulator (right) + Tilting Base
MD4-M3-L	Double-Holder Micromanipulator (left) + Tilting Base

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

OFFICIAL	ACCESSORIES/REFEACEMENT FARTS
M3301EH	Replacement Electrode Holder (14 cm 5 Ø 7.2 mm)
15873	Angled Electrode Holder (13 cm long)
M4C	Microscope Stage Adapter
500475	Ball Joint, 7 cm long, for Ø 8 mm Holder
500476	Ball Joint, 4 cm long, for Ø 4 mm Holder



M9

Mechanical clamp tightens three rotatable joints simultaneously with one locking knob. Arm adjusts without distortion. Base exerts 100kg magnetic force for great stability. Fine adjustment for precise operations included.

Magnetic Base:

50 (w) x 60 (l) x 55 (h) mm (2.2 x 2.4 x 2.2 in.)

Vertical Holding Power:

100 kgf (220 lb force)

Arms:

L1: 119 mm (4.7 in.) L2: 106 mm (4.2 in.) L3: 25 mm (0.98 in.) Ø 12 mm (0.472 in.)

Clamp Hole:

none

Weight:

1.8 kg (4 lb.)

ORDERING

И9

Magnetic Stand

M10

Similar to M1 with a 12mm diameter sub pole (fits 12mm clamp supplied with M3301, DC3001, MD4 and MMJ manipulators).

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 14 mm (0.55 in.) length: 178 mm (7 in.)

Sub Pole:

diameter: 12 mm (0.47 in.) length: 165 mm (6.5 in.)

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb.)

ORDERING

M10 Magnetic Stand

M10L

Same as M10, but equipped with a taller (14-inch) vertical main pole.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 14 mm (0.55 in.) length: 356 mm (14 in.)

Sub Pole:

diameter: 12 mm (0.47 in.) length: 165 mm (6.5 in.)

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb.)

ORDERING

M10L Magnetic Stand

M11

Bends freely for maximum flexibility. The connecting arm twists and bends like a snake. Lock the arm in position with a flick of the controlling lever.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb force)

Main Pole:

diameter: 16 mm (0.63 in.) length: 315 mm (12.4 in.)

Sub Pole:

none

Clamp Hole:

Adjustable from 6 mm to 8 mm

Weight:

1.4 kg (3 lb.)

ORDERING

M11 Magnetic Stand

Powerful Ball Joint Rare Earth Magnet

Construct holding devices for small parts/equipment

Small but very powerful. Holds 2 kilograms (~5 pounds)! $^{\rm M}$ Steel ball rotates freely 360° on a 180° axis

M3 mounting screw on ball for attachment to equipment Magnet base threaded (M3) for mounting onto a base or equipment

This novel magnetic ball joint has phenomenal holding power for up to 2kg of attached weight while permitting the ball a full 360° rotation on a 180° axis. You can

freely orient your equipment to an infinite number of positions within this rotation. This is made possible by the combination of a steel ball (10 $\,$

mm diameter) and a powerful rare earth magnet contained in the magnet cylinder (Ø10 x 20 mm). Convenient M3 attachment sites are provided on both the ball (male) and the magnet base (female). For use with micromanipulators for the positioning and holding of optical instruments including various lighting sources and lasers, pipettes and any small parts that would benefit from the flexibility offered by this new magnetic ball joint.




500871

ORDERING INFORMATION

500871 Magnetic Ball Joint

WORLD PRECISION INSTRUMENTS

The base of each stand exerts a powerful magnetic force that holds it solidly on ferrous metal surfaces — even vertically or upside down



MB2

Mechanical clamping type tightens three joints simultaneously. Arm is freely adjustable without distortion. Equipped with fine adjuster and medium size magnet for stabilizing the base. Suitable for performing precision operation. Fits clamp for MM1-ALL and MM3-ALL.

Magnetic Base:

30 (w) x 35 (l) x 35 (h) mm (1.2 x 1.4 x 1.4 in.)

Vertical Holding Power:

17 kgf (37 lb. force)

Arm:

L1: 46 mm (1.8 in.)

L2: 46 mm (1.8 in.)

L3: 39 mm (1.5 in.)

Clamp Hole:

Adjustable from 5 to 6.5 mm

Weight:

0.38 kg (0.83 lb)

ORDERING

MB2 Compact Magnetic Stand

M8

A ball joint at the base of the main post allows 360° rotation, offering considerable versatility. The second arm adopts angles up to 75°.

Magnetic Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb. force)

Main Pole:

diameter: 12 mm (0.47 in.) length: 194 mm (7.6 in.)

Sub Pole:

diameter: 10 mm (0.39 in.) length: 165 mm (6.5 in.)

Clamp Hole:

Adjustable from 4.5 mm to 6.5 mm

Weight:

1.8 kg (4 lb.)

ORDERING

M8 Magnetic Stand

A precision base providing stable support for such devices as electrodes and manipulators. Adjustable second arm adopts a variety of angles

Base:

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb. force)

Main Pole:

diameter: 12 mm (0.47 in.) length: 176 mm (6.9 in.)

Sub Pole:

diameter: 10 mm (0.39 in.) length: 165 mm (6.5 in.)

Clamp Hole:

diameter: 4.5 mm and 6.5 mm

Weight:

1.8 kg (4 lb.)

ORDERING

Magnetic Stand

M₁L

10 mm Ø

Same base and support arm as M1, but equipped with a longer (14-inch) vertical post.

Base:

10 mm ∅

50 (w) x 58 (l) x 55 (h) mm (2.0 x 2.3 x 2.2 in.)

Vertical Holding Power:

80 kgf (176 lb. force)

Main Pole:

diameter: 12 mm (0.47 in.) length: 356 mm (14 in.)

Sub Pole:

diameter: 10 mm (0.39 in.) length: 165 mm (6.5 in.)

Clamp Hole:

diameter: 4.5 mm and 6.5 mm

Weight:

1.8 kg (4 lb.)

ORDERING

M1L Magnetic Stand

Versatile Manipulator Mounting Clamp

Mount your Micromanipulator



Three of the stands above — M1, M1L and M8 — have 10 mm diameter mounting rods. The standard mount on several WPI manipulators (DC3001, KITE, M3301, MMJ, and MD4) accommodates a 12 mm rod. In order to use one of these three stands, you will need to replace the manipulator's standard 12mm mounting clamp with the optional M5 clamp.

ORDERING INFORMATION

M5 Ø10 mm Clamp

Magnetic Holding Devices

A solid platform for mounting your manipulators



501651

Round Base

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

Easy ON/OFF operation using lever Thin and powerful magnetic force Generous array of tap holes

Holding Power:

20 kgf (44 lb. force)

Dimension:

75 (OD) x 20 (h) mm 2.9 (OD) x 0.8 (h) in.

Mounting Hole:

4-M4 x 0.7, depth 6 mm* M8 x 1, depth 6 mm Span 35 mm

Weight:

0.7 kg (1.5 lb.)

ORDERING

501651 Magnetic Base, 75 mm diameter 503568 Magnetic Base, 50 mm diameter

* Posts with M4-threads not available from WPI.

An ideal accessory for optical tables and vibration-free platform. Reduces experimental set-up time by allowing free positioning and instant clamp down of optical components. Switchable ON/OFF magnetic circuit permits fine adjustment and precise positioning.

Easy ON/OFF operation using lever Thin and powerful magnetic force

Generous array of tap holes

Holding Power:

20 kgf (44 lb. force)

Dimension:

65 (w) x 65 (l) x 20 (h) mm 2.6 (w) x 2.6 (l) x 0.8 (h) in.

Mounting Hole:

8-M4 x 0.7, depth 6 mm* M8 x 1, depth 6 mm Span 25 mm

Weight:

0.6 kg (1.3 lb.)

ORDERING 501653 Magnetic Base,

	65x65 mm	
503569	Magnetic Base,	
	45x45 mm	
503570	Magnetic Base,	
	90x90 mm	
503571	Magnetic Base,	
	120x120 mm	

* Posts with M4-threads not available from WPI.



MOBITY

MOBITY™ is a new magnetic clamping system. With its ease of use, only one hand is needed to operate the attractive power. The MOBITY™ has a strong 88lbf pull, yet weighs only 1.5 lbs. MOBITY™ meets various applications with 4 tapped holes on the top surface. Requires (1) 9V alkaline battery (included).

Holding Power:

40 kgf (88 lb. force)

55 (w) x 73 (l) x 50 (h) mm 2.2 (w) x 2.9 (l) x 2.0 (h) in.

Mounting Hole:

3-M4, depth 20 mm* M8, depth 15 mm

Weight:

0.7 kg (1.5 lb.)

ORDERING

501652 MOBITY Magnetic Clamping System * Posts with M4-threads not available from WPI. .

M7

A small holder ideal for use where space is limited. Main post unscrews from base which may then be used alone as a switchable magnetic holder.

Magnetic Base:

30 (w) x 35 (l) x 35 (h) mm 1.2 (w) x 1.4 (l) x 1.4 (h) in.

Vertical Holding Power:

20 kgf (44 lb. force)

Main Pole:

Diameter: 7 mm (0.28 in.) Length: 52 mm (2 in.)

Clamp Hole:

Diameter: 6 mm

Weight:


0.36 kg (0.8 lb.)

ORDERING

Compact Magnetic Stand

Base Plates

A magnetic stand requires a steel mounting surface. WPI's steel base plates have plenty of mass to give stability to your experimental setup. Beveled edges make them easy to handle. Rubber feet hold them off the benchtop, making them easier to grasp when moving. The special black coating provides a durable, protective, non-reflective finish.






ORDERING INFORMATION

Steel Base Plate, 4.5 kg (10 lb) 5052 5479 Steel Base Plate, 14.5 kg (32 lb)

Mounting a Manipulator

Tilting Base with Screw Adjustment

At last! A tilt base you can operate with one hand!



- Holes also allow permanent mounting to your bench top
- · Adjust tilt angle by turning knob
- · Manipulator mounting bracket included

Benefits

- 13x15 cm (5x6 in.) footprint saves space in your work area
- Two sets of mounting holes are pre-drilled for WPI manipulators (M3301R shown) but steel platform may be drilled for mounting other devices

Applications


- Microinjection
- Electrophysiology amplification
- Stimulation

ORDERING INFORMATION

TBS Tilt Base with Screw Adjustment

Shipping Weight 7 lb.

Ball-joint holder attachmentHere's a new angle for mounting



	ORDERING INFORMATION
500475	Ball Joint, 7 cm long, for O.D. 5-9 mm Electrode Holder (shown)
500476	Ball Joint, 4 cm long, for O.D. 2.8-4.5 mm Electrode Holder

Manipulator Clamps Interchangeable clamps allow manipulators to be mounted on a variety of supports. M5 Ø 10 mm clamp M2 Ø 12 mm clamp ORDERING INFORMATION M2 Ø 12 mm Clamp M5 Ø 10 mm Clamp

Ø 12.7 mm (½-in.) Clamp

М6



Environmental Control for Live Cell Microscopy

Incubation and environmental control for long term imaging of cellular growth

Features

- Unique diffusion grid, combined with air input and return vents, provide an air flow pattern for consistent, even heating, with no hot or cold spots in the chamber
- External heater that can be placed far enough from the system to eliminate electrical and vibrational interference from the heater
- High degree of temperature precision and stability
- Minimal focal drift after equilibrium is achieved—accuracy ±0.1°C at the sample itself, and 0.2°C across the microscope stage (allowing for uniform heating of multiwell dishes)
- Airflow pattern and temperature uniformity eliminate dramatic changes in environmental temperature when the incubator door opens
- Ergonomic design for ease of use— The focus and x/y stage controls are outside of the incubator itself.
 Large doors allow easy access to the specimen and small ones for cords, tubing, etc.
- Precision, shielded temperature probe
- Simple, one person setup of the system


Benefits

- Choice of controller to manage air flow, heat, carbon dioxide and oxygen
- Control and monitor an external heater or lens warmer with ECU or control and monitor temperature and monitor humidity with ECU-HOC
- Diffuser grid and proper air venting insures consistent air flow inside the acrylic microscope chamber

Applications

- Live cell imaging under a microscope
- Time lapse video research


This Live Cell Microscope Incubator was extensively tested in laboratories. When compared with other systems, it offers dramatic advantages. For example, other incubators for live cell microscopy rely on passive, random



Incubator on a Nikon Eclipse IE-2000 Microscope. Incubators for Zeiss and Olympus microscopes are also available, as well as versions that accommodate confocal modules. All incubators are compatible with all commercially available cameras, light sources, filter wheels, motorized stages and motorized nose pieces.

diffusion of heated air from a single source to maintain the desired temperature setpoint. With no hot air return vent, the heated air escapes from the system through cracks at the microscope/incubator junction in an uncontrolled, random fashion. These systems offer no temperature uniformity, suffer from focus drift and often experience electrical and vibrational interference from the heater. You will also notice dramatic temperature drifts when the imaging environment is disturbed.

This diagram shows how the Stagetop environment connects with the ECU-HOC to control heat, air flow, carbon dioxide and oxygen.



WORLD PRECISION INSTRUMENTS

Choice of controller to manage heat, CO₂ and O₂

This unique, acrylic Live Cell Imaging chamber, combined with an Environmental Control Unit (**ECU**) and an **AirTherm** controller, ensure precision control of your incubator environment.

The ECU comes in different varieties so you have all the control you require.

- With the ECU-H5, you can control air flow and heating.
- In addition to air flow and temperature control, the ECU-HC lets you control the carbon dioxide level. It has an internal sensor.
- The ECU-HOC adds control of the oxygen level, which is accomplished by displacing the oxygen with nitrogen.

Control external heater and monitor humidity

The first four **ECU** units are capable of controlling a simple, external heater, like the AirTherm Satelite (**AirTherm-SAT**) or a microscope lens warmer. The **AirTherm-SMT** can monitor and control temperature and monitor humidity level inside the microscope chamber.

Diffuser grid and venting for consistent air flow

Air flow affects the temperature uniformity of incubators. The red arrows on Fig. 1 and Fig. 3 indicate air flow. The Live Cell Microscope Incubator uses a diffuser grid and proper venting to insure consistent air flow. Traditional incubators with poor air flow suffer with hot and cold spots in the incubator, as seen in thermal images (Fig. 2 and Fig. 4). Warmer temperatures are indicated by red and cooler temperatures by blue.



Fig. 1—Single air input and no venting causes random air flow in a traditional incubator.



Fig. 2—Hot and cold spots result from inconsistent flow.



Fig. 3—A diffusion grid with air input and exhaust vents yields consistent air flow.



Fig. 4—Consistent air flow means uniform heating.



AIRTHERM-SAT





This humidifier may be used inside a microscope chamber for controlling the humidity inside the stagetop environment chamber.

ORDERING AN ENVIRONMENTAL CHAMBER

Acrylic enclosures are essentially custom-built. When ordering a system, you will need to provide the following information:

Microscope	Stage	Stage-Up
Perfect Focus	Camera	Left Port
Right Port	Analyzer	Fluor Attachment
Tirf	White Light Tirf	Binocular D Head
Tilting Head	Filter Wheels	Excitation
Emission	Dual Lamphouse	Transmitted Light Shutter

Cells 35 mm/60 mm Wells

Coverslips

ORDERING INFORMATION FULL CHAMBER SYSTEM WITH ECU CONTROLLER & NEW

AIRTHERM-SAT/SMT

All Systems Include: Proprietary Humidification Module, Stage

Adapter, Stage Dish with Optical Grade Glass
FC-HC Microscope Environmental Chamber, Heat & CO₂ Controller

Requires 100% CO_2 & ambient air supplies.

FC-HOC Microscope Environmental Chamber, Heat Controller, CO_2 and O_2 Controller

FC-H5 Microscope Environmental Chamber, Heat Controller, Digital Control for Bottle Gas

FC-SMT Microscope Environmental Chamber, Heat Controller (AIRTHERM-SMT instead of ECU)

Requires pre-mixed 5% CO₂ gas supply
Optional Humidifier for AirTherm-SMT

AIR-THERM-SMT



CE

Stagetop Environmental Control

Control temperature and CO₂ in a microscope stagetop environment

Features

- Four programmable digital control loops:
 - Independent incubator base temperature PID control with ±0.1°C
 - Independent incubator lid temperature PID control with ±0.1°C precision
 - CO₂ digital PID control with ±0.1% precision
 - Airflow digital PID control from 0-900 SCCM
- · USB-based remote control and data logging
- · Electronic flow meter
- Programmable alarm for out of tolerance condition on all four
- Compact and lightweight

Benefits

- · Compact housing that fits most inverted microscope stages and holds standard culture well plates
- Control temperature and CO₂, O₂ using the Environmental Control Unit
- Control and monitor system parameters using the ECU

Applications

- For short term or long term studies of living cell cultures under a microscope (Live Cell Imaging)
- Time lapse video research

Time lapse video research requires a microscope stagetop incubator.

Compact housing fits most inverted microscope stages

Perfect for live cell imaging, STEV2 (the stagetop environmental control platform) is a compact environmental case that houses your culture wells and fits on a microscope stage inside the live cell microscope incubator.

Control temperature and CO₂ and O₂

This system offers precision control of both temperature and carbon dioxide, as well as remote control and data logging via a USB connection. The system is flexible and easy to configure for a variety of experimental

Control and monitor system parameters

The system includes the Environmental Control Unit electronics which use four programmable loops to control the temperatures of the case and the lid, CO₂ within the environmental case and airflow within the incubator.





The Environmental Control Unit home screen shows the real time readings for the parameters you want to see.

ECU SPECIFICATIONS

Non-dispersive infrared (NDIR), dual beam, 20s response time 110/240V, 50/60Hz

Zirconium Dioxide, diffusion, 4 sec

10 - 50°C (50 - 122°F) Operating Temperature (ambient) Operating Humidity (ambient) 15 - 70% RH, non-condensing Warm up Time 20 minutes

Computer Interface USB via external USB/RS232 converter

Sensor

Power

CO2 Sensor (ST-HC only) Sensor Range

Control Range Control Precision Control Accuracy

O₂ Sensor (ST-HOC only)

Sensor Type

response time Sensor Range 0-25% 0-25%* Control Range Control Precision 0.1% O₂

Control Accuracy ±0.5% (2% of the full scale)

ORDERING INFORMATION

STAGETOP ENVIRONMENTAL CONTROL SYSTEMS

ST-H5	Chamber, Controller, Heat, Digital Flow Control for Bottle Gas
ST-HC	Chamber, Controller with CO ₂ & Heat
	(using CO ₂ internal sensor)
ST-HOC	Chamber, Controller with CO ₂ , O ₂ & Heat
	(using CO ₂ internal sensor)

ACCESSORIES

IRTHERM-SAT-1W	AirTherm Satellite Heater (110V)
IRTHERM-SAT-2W	AirTherm Satellite Heater (230V)

0 - 20% CO₂

<2.5% reading/year

0 - 20%

0.1% CO₂ 0.1 - 3% of reading

Heat Controller/Humidity Monitor

Smart, electrically quiet air heater for live cell imaging & custom incubators

Features

- Precision heat controller for use in live cell imaging and custom incubators
- Control heat and monitor humidity (optional) with a single controller
- · Electrically and acoustically quiet
- · Quick, precise response to thermal change

Benefits

- PID control algorithms allow for tight control of temperature in the environmental chamber
- System includes hoses and probes so that it is ready to connect with your microscope chamber, incubator or controlled environment

Applications

- Live cell imaging systems fitted with full microscope environmental chamber enclosures
- · Control environment in custom incubators

The AirTherm™ SMT is a new generation of heat control system from WPI designed to be used in Live Cell Imaging applications with microscopes fitted with a full microscope environmental chamber enclosure. The standard AirTherm™ SMT controls temperature and, as an option, monitors humidity.

PID algorithm for maintaining tight control

The **AirTherm™ SMT** uses a PID control algorithm to tightly control the temperature and monitor humidity of the controlled environment


With AirTherm™ SMT, the temperature of the sample and microscope optics can be controlled to within 0.2°C. During operation, air is drawn out of the chamber through a flexible hose, heated by the AirTherm™ SMT heater and re-circulated to the chamber by the return hose.

Ready to connect to your microscope chamber

The system is typically used in a closed loop configuration.

The AirTherm™ SMT system includes:

- Two coil-reinforced heater hose pieces and hose clamps.
- Temperature sensor for remote placement in environmental chamber
- Humidity probe for monitoring chamber humidity available as an option.



A typical AirTherm SMT installation places the heated air inflow at the bottom and the cold air return at the top of the microscope chamber.



AIR-THERM-SMT

References

X. Yin, D.A. Knecht, M.A. Lynes "Meallothionein mediates leukocyte chemotaxis" *BMC Immunology* 6. 2005: 6--21

AIRTHERM SPECIFICATIONS

AIR FLOW RATE 20–50 CFM (0.55–1.4 m³/minute)

CONTROL TEMPERATURE RANGE Ambient to 45°C

TEMPERATURE RESOLUTION 0.1°C
TEMPERATURE ACCURACY 0.2°C

ANALOG OUTPUT 0.5°C resolution;

FOR CHART RE-CORDER 0–10 V represents 0–100°C

HEATING VOLUME Less than 50 CF (1400 L), re-circulating

TEMPERATURE SENSOR TYPE Platinum RTD 1000 Ω

HUMIDIFIER TYPE Ultrasonic
HUMIDIFIER TANK CAPACITY 0.5 gallons
HUMIDIFIER DAILY OUTPUT 2 gallons

FUSE For 120 VAC, 8A 250 V 5x20 mm metric

For 230 VAC, 4A 250 V 5x20 mm metric

POWER 450 W, 95–135 V or 220–240 V, 50/60 Hz DIMENSIONS 6½ x 8 x 7½ in. (15.5 x 21 x 19 cm)

ORDERING INFORMATION

AIRTHERM-SMT-2W AirTherm™ SMT Heater, 230V **AIRTHERM-SAT-2W** AirTherm Satellite Heater, 230V

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

15590	Replacement Hoses, 2.5" diam., 4.5 ft
98727	Replacement Temperature Probe
98728	Humidity Probe
300626	AirTherm SMT Humidifier Kit with Sensor & Hose (110 V)
300411	Replacement Humidifier for AirTherm SMT (110 V)
5389	Inlet and Outlet Hoses
96824	3-Way Hose Assembly (3 hoses and T-connector)

Bench Top Vibration Isolation Platforms

Simple setup and adjustment

Features

- Only 4.6" thin
- Portable
- · Available in a variety of payload ranges

Benefits

- Better performance than an air table
- · Simple setup and adjustment
- Requires no air or electricity

Applications

 Microscopy and other research applications that require the isolation of the preparation from vibrations of any kind

These bench top platforms offers 10-100 times better performance than a full size air table in a package only 4.6 inches tall, and without air or electricity! These vibration isolation platforms are extremely easy to use and offer extreme performance — 1.5 Hz horizontal natural frequency and 0.5 Hz vertical natural frequency. There are only two adjustments. This is the thinnest, most portable and most user-friendly isolator ever offered that is capable of delivering this level of performance.

Weight: Approximately 40 lb. (16 kg)

Dimensions: 18" W x 20" D x 4.6" H (457 x 508 x 117 mm)

MK-BM-8100, 50-105 lb. payload weight range (23-48 kg)

Performance

- · Horizontal frequencies are weight dependent
- Horizontal frequency of 1.5 Hz is achieved at or near the upper limits of the payload range
- At the lower limits of the payload range the horizontal frequency is approximately 2.5 Hz
- Vertical frequency is tunable to 0.5 Hz throughout the payload range

	ORDERING INFORMATION
MODEL	PAYLOAD RANGE
MK-BM-825	Vibration Platform, 10 - 30 lb. (4.5 - 14 kg)
MK-BM-850	Vibration Platform, 25 - 55 lb. (11 - 25 kg)
MK-BM-8100	Vibration Platform, 50 - 105 lb. (23 - 48 kg)
MK-BM-8125	Vibration Platform, 90 - 130 lb. (40 - 59 kg)
MK-BM-8150	Vibration Platform, 125 - 155 lb. (57 - 70 kg)
MK-BM-8175*	Vibration Platform, 150 - 180 lb. (68 - 81.5 kg)
MK-BM-8200*	Vibration Platform, 175 - 205 lb. (79.5 - 93 kg)
MK-BM-8225*	Vibration Platform, 200 - 230 lb. (90.5 - 104 kg)
MK-BM-8250*	Vibration Platform, 225 - 255 lb. (102 - 115.5 kg)
	1. 4

* Weight: Approximately 47 lb. (21 kg)/(same dimensions)

Vibration-Free Tables

Vibration-Free Workstation

Features

- Vertical and horizontal vibration isolation
- High performance active-air suspension
- Automatic leveling
- · VibraDamped steel
- Class 100 clean room compatible
- Leveling feet

Benefits

- Eliminate inconsistent and unreliable performance
- · Reduce excessive wear, maintenance and fatigue failures
- Protect sensitive instruments and equipment from faulty operation or failure using Vibration-Free Platforms and Workstations

Applications

 Patch clamping, cell injection, analytical balances and optical microscopes

All buildings vibrate. Activities of people, machinery, heating and ventilation systems, and nearby truck or rail traffic cause all types of vibrations. These vibrations, though acceptable to occupants, cannot be tolerated by equipment used in patch clamping, cell injection, analytical balances and optical microscopes. Additional tabletop sizes and finishes are available, as well as optional accessories such as side rails and casters. See WPI website for available configurations.

Universal Manipulator Stand

Mount manipulators at angles and heights with infinite flexibility

UMS SPECIFICATIONS

DIMENSIONS

Base Plate 10.0 x 12.5 x 1.5 cm (LxWxH) Stand 4.0 x 4.0 x 30 cm (LxWxH) (501622)

4.0 x 4.0 x 30 cm (LxWxH) (301622) 4.0 x 4.0 x 45 cm (LxWxH) (501623)

Mounting Holes

English: 1/4 20 x 1" (2 bolts supplied) Metric: M6 x 25 mm grid (2 bolts supplied)

SHIPPING WEIGHT

501622 9 lb. (4 kg) **501623** 11 lb. (5 kg)

Vibration-free Platform (VFP) not included

501622

	ORDERING INFORMATION
501622	Universal Micromanipulator Stand 30 cm (includes one clamp)
501623	Universal Micromanipulator Stand 45 cm (includes one clamp)
501624	Additional Rotation Clamp


The M3301, KITE, DC3001 and SM325 (but not M325) can be used with these mounts.

Custom order.

part numbers.

Call for pricing and

Microinjection



We offer a complete microinjection system

WPI has been in the microinjection business since the very beginning. We were first on the market with pressure injection equipment, and we have been engineering and delivering innovative instruments ever since.

For essential fluidic control in neurophysiology, we offer solutions across pumping platforms including syringe, peristaltic, microfluidic and injection pumps. For picoliter to microliter ranges, our tools are designed for applications like patch clamp, oocyte work, infusion, perfusion, cellular injection, continuous flow and non-pulsatile flow.



WPI Microinjection System

The system depicted at left includes components often favored by researchers, indicated with * in the list below: PV820
Pneumatic PicoPump, PUL-1000 Micropipette Puller, M4C
stand, M3301R micromanipulator, 5430-XX PicoNozzle Kit

with a µTip, PZMTIII microscope, lighted base, articulating mirror and optional PRO-300 HDS camera and view screen, E2XX micropipette storage jar, Z-MOLDS Microinjection and Transplantation Molds, 14003-G Vannas spring scissors, glass capillaries, 77020 glass tweezers and Fluorodish optical grade glass bottom culture dishes. Whatever your needs, WPI offers a range of equipment to fill your requirements.

Options for Customizing Your System

INJECTORS FOR GLASS PIPETTES

- **PV820** Pneumatic PicoPump with Hold Pressure
- **PV830** Pneumatic PicoPump with Hold Pressure and Vacuum
- NL2010MC2T Nanoliter Injector with SMARTouch Controller



Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's **PicoPumps** use accurate timing and carefully regulated air pressures for securing cells and injecting. Injected volumes range from picoliters to nanoliters.



Microprocessor-controlled Nanoliter Injector **NL2010MC2T** with SMARTouch controller uses positive displacement injection, eliminating the need for pipette calibration. System uses glass micropipettes.

PINPOINT INJECTOR

****** MICRO-ePORE™ Cell Penetrator



A simple and versatile pinpoint controlled electroporation system for clean efficient injection of CRISPER reagents into 2-stage embryos. The cell membrane is easily penetrated, without cell disruption or tear drop formation. The technique has been used in animal embryos, including, mice, rats and bovine, as well as injecting morpholino oligomers into Zebrafish tails.



INJECTORS FOR GLASS PIPETTES OR METAL NEEDLES

The versatile UltraMicroPump *UMP3* injector is a Syringe-based micropump using metal needles as small as 36 gauge and microsyringes to deliver picoliter volumes.

A glass pipette injection kit is also available

Use with the NANOFIL syringe and 1mm OD glass for your microinjection application."





MICROSCOPES

- **PZMTIII** Precision Stereo Zoom Trinocular Microscope
- PZMIII Precision Stereo Zoom Binocular Microscope
- **PZMIII-MI** Binocular Microscope Head with LED illuminated base (articulating mirror)

PZMTIII-MI Trinocular Microscope Head with LED illuminated base (articulating mirror) and variable light intensity. Dual reflection lens/ mirror system provides transmitted brightfield/ pseudo-darkfield illumintation.



The Trinocular Microscope may be used with the PRO-300 HDS.

Photo Courtesy of Benjamin Dubansky, PhD. University of North Texas

NOTE:#Base can be used with other manufacturer's microscope optics

PULLERS

- **PUL-1000** Microprocessor-Controlled 4-Step Micropipette Puller
- SU-P97 Flaming/Brown Pipette Puller
- SU-P1000 Next Generation Pipette Puller

WPI is an authorized distributor for Sutter Instruments

PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, microperfusion or microinjection. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. Perfect for long taper pipettes.



MANIPULATORS

- 🗱 M3301 Manual Micromanipulator
- KITE Manual Micromanipulator
- M4C Microscope Stage Adapter
- **SU-TRIO** 3-Axis Micromanipulator
- SU-QUAD 4-Axis Micromaniulator



Weighing just 550 grams, the M3301 is a well-built micromanipulator that outsells all others worldwide for high precision experiments where magnification is in the range of up to 250×.



GLASS MICROPIPETTES

- Pre-Pulled Long Taper Glass Pipettes (TIP15FLT, TIP12FLT)
- Glass Capillaries
- µTip Pre-Pulled Pipettes with Calibrated Tip ID



Calibrated Tip ID. Get pre-pulled, pipettes with or without a filament. Available in a variety of sizes. Luer connect option is available.




Eliminate the cost and trouble of pulling your own pipettes. Get new low cost, long taper pipettes with or without a filament. Taper length: 12-14 mm.

ACCESSORIES

- PRO-300 HDS Camera and Monitoring Screen
- 5430-XX PicoNozzle Kit
- NanoFil Microliter Syringes 10 or 100 µL
- NFINHLD-G10 1.0 mm Glass Pipette Holder for NanoFil Syringe
- MicroFil for backfilling glass needles
- **Pipetters**
- E2XX Micropipette Storage Jar
- **801566/801963** Vacuum Pump for use with the **PV830**
- FluoroDish Optical Glass Bottom Dishes
- M10 or M-3 Manipulator Base
- **Z-MOLDS** Microinjection and Transplantation Molds
- Many Surgical Instruments

Pinpoint Injector

The next generation injection method for increased knock-in efficiency



Features

- · Touch-screen display—resistive touch panel for use with gloves
- Injection control with foot switch or manually with touch screen
- Intuitive user interface
- · User adjustable frequency and voltage with touch screen or knobs
- · Small footprint
- Four user-programmable protocols
- Adjustable audio continuity tone indicates active probe
- Injection counter indicates total number of injections
- Electrodes include reference electrode option
- Electrode holder compatible with most common manipulators

Benefits

- High yield (80-85%) of viable, 2-cell stage mouse embryos after injection
- Integrates with WPI's PV820/PV830 or the Eppendorf Femtojet® (Femtojet® 4i is a registered Trademark of Eppendorf AG)
- · Hands-free operation with foot switch control
- Works with most inverted microscopes

MICRO-ePore successes in university labs



Photo courtesy of Dr Janet Rossant, The Hospital for Sick Children in Toronto, Canada. An example of the MICRO-ePORE™ electrode setup on an inverted microscope. MICRO-ePORE™ electrode design ensures compatibility with the leading microscope manufacturers.

Applications

- Microinjection into oocytes and pre-implantation stage mammalian embryos, including microinjection of CRISPR-Cas9 reagents into the cytoplasm of two-cell stage embryos
- · Pronuclear rodent zygote microinjection
- · Gene silencing in zebrafish



MICRO-ePORE SPECIFICATIONS

VOLTAGE PARAMETERS 0–3.0 V, at 1 mV increments FREQUENCY PARAMETERS 50–3000 Hz, at 1 Hz increments

PIPETTE RESISTANCE MAXIMUM 40 MΩ

DIMENSIONS 19.7 × 12.7 × 7.6 cm (7.75 × 5 × 3 in.)

WEIGHT 0.9 kg (2 lb)
CERTIFICATIONS CE, ROHS

ORDERING INFORMATION

MICRO-ePORE MICRO-ePore System™

Includes Electrode Holder • Microinjection Headstage, including adapter for Femtojet™ • Holding Pipette • Platinum Reference Electrode • Vacuum Tubing Set • Foot Switch (Auxiliary Knob Controller sold separately)

98991	Auxiliary Knob Controller
SYS-PV820	Pneumatic PicoPump Microinjector
SYS-PV830	Pneumatic PicoPump Microinjector with Vacuum Pressure
1B100-3*	Single-Barrel Standard Borosilicate Glass Tubing, 3 in.
	(76mm) long, 1.0 mm OD, 0.58 mm ID, package of 500
TIP1TW1*	Pre-Pulled Glass Pipettes, 1 µm tip ID, 1.0 mm glass OD,
	Thin wall, fire polished

*Many varieties and choices are available for glass capillaries and pre-pulled µTip micropipettes.

Reliable Pneumatic PicoPump

Pneumatic microinjector with vacuum pressure



PV830

Microinjection of sub-nanoliter volumes through a very small opening in a glass micropipette is often accomplished best with a controlled burst of gas pressure. Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's PicoPumps use precisely regulated pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters.

The **PV830** model provides separate regulation and pressure gages for the eject and hold pressures. In addition to the features of the **PV820**, the **PV830** model also includes a vacuum pressure regulator and vacuum gauge.

See page 44 for more information.

Microelectrode Puller

A compact, versatile and reliable workhorse




PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, patch clamp studies, microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.

This puller is a reasonably priced, compact, versatile and reliable workhorse. The microprocessor, combined with the LCD display, makes the ${\bf PUL-1000}$ easy to use.


See page 108 for more information.

Glass Capillaries


Quality glass, superior prices for microinjection/microelectrodes




Standard, no Filament



Standard Filament



Thin Wall, no Filament



Thin Wall Filament

Our quality borosilicate glass capillaries are available with a large variety of options, including fire polished, filaments, thin wall, specialty glass and multi-barrel.

Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

For details, see page 47.

Reliable Pneumatic PicoPump

Pneumatic microinjector with vacuum pressure



PV830

Features

- Inject into a single cell with picoliter volumes
- Regulated hold, ejection and vacuum pressure
- Carefully regulated air pressures for securing cells and injecting them with fluid
- · Optional independent vacuum regulation
- Pressure Input: 0-120 psi
- · Pressure Output: 0.3-90 psi
- · Choice of two models

Benefits

- Hold pressure prevents backfilling of the pipette by capillary action
- PV830 includes vacuum pressure for filling pipettes from the tip or securing floating cells during injection

Applications

• Intracellular injection in the picoliter to nanoliter range

Although syringe and piston techniques are popular for microinjection in the nanoliter range, microinjection of sub-nanoliter volumes through a very small opening in a glass micropipette is often better accomplished with a controlled burst of gas pressure. Designed to simplify intracellular injection and a variety of other microinjection tasks, WPI's PicoPumps use precisely regulated pressures for securing cells and injecting them with fluid. Injected volumes range from picoliters to nanoliters. Separate ports supply positive and negative pressure—positive pressure for highpressure ejection, and suction for supporting the cell or for filling the pipette from the tip. In the time between injections, a secondary default "hold" pressure is applied to the pipette to prevent fluid uptake through capillary action or diffusion. The pressure pulses are typically controlled by a precision timer. Timing, ejection pressure, holding pressure and suction are adjusted independently by control knobs and monitored with indicator gauges on the front panel. Injection pressure is controlled by a 20-turn regulator on the front panel. The built-in timing function controls the amount of time that the injection pressure is applied with millisecond resolution. Time intervals can range from 10 s down to 10 ms or less. The pressure burst can also be controlled manually or triggered by an external source.

The **PV830** model provides separate regulation and pressure gauges for the eject and hold pressures. In addition to the features of the **PV820**, the **PV830** model also includes a vacuum pressure regulator and vacuum gauge.

Hold Pressure Prevents Backfilling

Like the **PV830**, the **PV820** was designed to simplify intracellular injection. The hold pressure prevents backfilling of the pipette by capillary action. Simple to use. Reliable every time. Each PicoPump is





- (1) MPH6S microelectrode holder
- (1) Handle for the MPH6S (4" hollow tube with male Luer fitting at both ends—diameter 6.25 mm × 100 mm long)
- (1) 5' tubing (0.060" ID, 0.120" OD, male locking Luer fitting on one end and a female locking Luer fitting at the other end, rated for 200 PSI and 86 durometer shore A)

supplied with a **5430-ALL** kit that includes two PicoNozzles and tubing to connect the holders to the pressure and vacuum ports.

Vacuum Pressure for Additional Applications

Eject pressure supplies a high-pressure pulse for injecting fluid. Hold pressure, which is not sufficient to cause fluid ejection, is used to prevent back filling of the pipette by capillary action or diffusion when the solenoid is inactive. Vacuum can be applied through the ejection port, but vacuum in the **PV820** must be externally regulated (see Mini Vacuum Pump, page 71). To fill pipettes from the tip or through a secondary port to secure a floating cell during microinjection, vacuum may be switched from regulated vacuum to atmosphere by using the switch on the front panel.

References

Di Cara F, Sheshachalam A, Braverman NE, Rachubinski RA, Simmonds AJ (2017). Peroxisome-Mediated Metabolism Is Required for Immune Response to Microbial Infection. *Immunity*. 2017 Jul 18;47(1):93-106.e7. https://www.ncbi.nlm.nih.gov/pubmed/28723556

Carrington, B., Varshney, G. K., Burgess, S. M., & Sood, R. (2015). CRISPR-STAT: an easy and reliable PCR-based method to evaluate target-specific sgRNA activity. *Nucleic Acids Research*, 43(22), e157. http://doi.org/10.1093/nar/gkv802

Pneumatic picoliter injector

WPI's **PV820** Pneumatic PicoPump has been providing scientists with precision and repeatable microinjection in volumes ranging from picoliters to nanoliters for decades. The trusted Pneumatic PicoPump is the preferred Zebrafish microinjection pump. The **PV820** offers the same functionality as the **PV830**, minus the independent vacuum regulator.



	PICOPUMP SPECIFICATIONS	
	PV820	PV830
PRESSURE		
PRESSURE INPUT	0 to 120 psi	0 to 120 psi
PRESSURE OUTPUT	0.3 to 90 psi *	0.3 to 90 psi
PRESSURE BURST TIMER (10-turn dial)	10 ms to 10 s in Timed Mode	10 ms to 10 s in Timed Mode
REGULATOR ACCURACY	0.1% (20-turn dial) *	0.1% (20-turn dial) *
REGULATOR REPEATABILITY	0.05 psi *	0.05 psi *
GAUGE ACCURACY	3% at full scale *	3% at full scale *
INPUT CONNECTOR	Quick Connect (¼ in. OD Tubing)	Quick Connect (¼ in. OD Tubing)
OUTPUT CONNECTOR	Barbed (1/16-in. ID Tubing)	Barbed (⅙-in. ID Tubing)
VACUUM		
VACUUM INPUT	0 to 30.0 in. Hg	0 to 30.0 in. Hg
VACUUM OUTPUT	Unregulated	0.2 to 29.9 in. Hg
LOWEST REGULATED VACUUM	Unregulated	3 in. water
REGULATOR ACCURACY	Unregulated	0.1% (20-turn dial)
REGULATOR REPEATABILITY	Unregulated	0.03 in. Hg
GAUGE ACCURACY	None	3% at full scale
INPUT CONNECTOR	Quick Connect (¼ in. OD Tubing)	Quick Connect (¼ in. OD Tubing)
OUTPUT CONNECTOR	Barbed (1/16 in. ID Tubing)	Barbed (⅓₅ in. ID Tubing)
CONTROL	Manual	Manual
VENT	Atmosphere	Atmosphere
CONNECTIONS INCLUDED		
INPUT KIT	10-ft nylon tubing (0.25-in. OD, 1000 բ	osi), one ½-inch female NPT adapter
OUTPUT KIT	Two PicoNo	ozzle Kits
PHYSICAL SPECIFICATIONS		
POWER	95-135 V or 220-240 V, 50/60 Hz	95-135 V or 220-240 V, 50/60 Hz
DIMENSIONS	17 x 3.5 x 9.5 in. (43 x 9 x 24 cm)	17 x 5.25 x 9.5 in. (43 x 13 x 24 cm)
SHIPPING WEIGHT	11 lb (5 kg)	14 lb (6.3 kg)

^{*}Both Hold and Eject Pressure

Ramezani, T., Laux, D. W., Bravo, I. R., Tada, M., & Feng, Y. (2015). Live Imaging of Innate Immune and Preneoplastic Cell Interactions Using an Inducible Gal4/UAS Expression System in Larval Zebrafish Skin. *Journal of Visualized Experiments*, (96), e52107–e52107. http://doi.org/10.3791/52107

Chokshi, P. (2015). Analysis of scn5Laa and scn5Lab Gene Function in Danio rerio (Zebrafish) Heart Development through TALENs/CRISPR-CAS9-mediated Gene Knockout. Retrieved March 14, 2016, from http://thescholarship.ecu.edu/bitstream/handle/10342/4923/Chokshi_ecu_0600O_11444.pdf?sequence=1

Ewart, M.-A., Kennedy, S., MacMillan, D., Raja, A. L., Watt, I. M., & Currie, S. (2014). Altered vascular smooth muscle function in the ApoE knockout mouse during the progression of atherosclerosis.

http://doi.org/10.1016/j.atherosclerosis.2014.02.014

Henson, H. E., Parupalli, C., Ju, B., & Taylor, M. R. (2014). Functional and genetic analysis of choroid plexus development in zebrafish. Frontiers in *Neuroscience*, 8, 364. http://doi.org/10.3389/fnins.2014.00364

Huang, L., Xiao, A., Wecker, A., McBride, D. A., Choi, S. Y., Zhou, W., & Lipschutz, J. H. (2014). A Possible Zebrafish Model of Polycystic Kidney Disease: Knockdown of wnt5a Causes Cysts in Zebrafish Kidneys. *Journal of Visualized Experiments*, (94), e52156–e52156. http://doi.org/10.3791/52156

Konantz, J., & Antos, C. L. (2014). Reverse Genetic Morpholino Approach Using Cardiac Ventricular Injection to Transfect Multiple Difficult-to-target Tissues in the Zebrafish Larva. *Journal of Visualized Experiments*, (88), e51595–e51595. http://doi.org/10.3791/51595

	ORDERING INFORMATION
SYS-PV820	PicoPump w/ hold pressure
SYS-PV830	PicoPump w/ hold pressure and vacuum
Specify I	ine voltage All PicoPumps require external vacuum source.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

3260	Foot Switch
2932	Rack Mount Kit, 3.5-in. high (PV820)
2933	Rack Mount Kit, 5,25-in. high (PV830)
5430-10	PicoNozzle Kit (MPH6S for 1.0 mm pipette & 5' tubing)
5430-12	PicoNozzle Kit (MPH6S for 1.2 mm pipette & 5' tubing)
5430-15	PicoNozzle Kit (MPH6S for 1.5 mm pipette & 5' tubing)
5430-20	PicoNozzle Kit (MPH6S for 2.0 mm pipette & 5' tubing)
5430-ALL	PicoNozzle Kit (for 1.0, 1.2, 1.5 and 1.65 mm pipettes &
	5-ft. tubing)
75122-110	Gaskets for PicoNozzle , 1.0 mm, green, package of 10
75122-210	Gaskets for PicoNozzle — 1.2 mm, black, package of 10
75122-310	Gaskets for PicoNozzle — 1.5 mm, blue, package of 10
75122-410	Gaskets for PicoNozzle — 1.65 mm, red, package of 10
MPH6S	Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)

Micropipette Holder (specify 1.0, 1.2, 1.5 or 2.0 mm)


Replacement Input Kit

MPH6R

3316

Manual Microsyringe Pumps

Precise manual injection of fluid using glass pipettes



Features

- Accepts syringe sizes: 10 µL-1 mL gas tight Luer tip
- Micrometer head has resolution of 10 µm per division
- · Solid, stainless steel frame

Benefits

- · Choice of models, with or without digital display
- · Cost-effective solution

Applications

• Use for perfusion or withdrawal of liquids

MMP and **DMP** are convenient tools for precise manual injection of fluid using glass pipettes or similar injection devices. The design allows visual feedback of flow at the pipette tip. They can also be used as a manual micro syringe pump for perfusion or withdrawal of liquids. The resolution of the injection volume can be continuously varied from 10 nL to the microliter range, depending on the syringe used.

With or without digital display

DMP comes with a digital micrometer that will allow the reading of piston advancement easily with a 0.001 mm resolution. MMP has the traditional mechanical micrometer head with a resolution of 10 μm per division and advances 500 μm per revolution.

Constructed of stainless steel

The entire frame body of the injector is constructed of stainless steel for excellent stability and durability. The piston of the micrometer can be adjusted to the syringe's plunger position. Small diameter PTFE tubing is used for the accuracy and solution compatibility. The unique design of the pipette holder can securely hold any pipette with an outer diameter of between 1.0 mm and 1.5 mm. All necessary accessories for removing air and filling the syringe and tubing with liquid are included. The system comes complete with a 100 μL gas tight syringe. Other syringe sizes may be purchased separately.

DMP & MMP SPECIFICATIONS

TRAVEL DISTANCE 25 mm

ADVANCES RESOLUTION 0.001 mm for DMP and 0.01 mm for MMP SYRINGE SIZE 10 μ L to 1 mL gas tight Luer tip syringe TUBING 1.5 m of PTFE tubing with 0.5 mm ID

PIPETTE HOLDER 0.24" x 5.2"

PIPETTE HOLDER FITS 1.0 to 1.5 mm OD pipette

ORDERING INFORMATION

MMP Manual Microsyringe Pump

DMP Manual Microsyringe Pump with Digital Display

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

MMP-KIT Injection Assembly Parts Kit

Not including valve—see 14057-10, page 220)

Brazil: 011 55 13 40629703 • info@brazil.wpiinc.com • www.wpiinc.com

Microinjection and Transplantation Molds

A simpler method for handling zebrafish embryos

Features

- · For high throughput microinjection research
- Make impressions in agarose gel to facilitate embryo alignment
- Four molds per kit
- Reusable

Benefits

• Organizes and immobilizes embryos for microinjection

Applications

· Zebrafish research

Mold your agarose, pipette in your embryos, and watch them auto aligning in the grooves. **Z-MOLDS** Microinjection and Transplantation Molds (4 per kit) are designed for zebrafish research.

Proteomics & Large Screening–Inject many embryos, up to 1000. The grooves made by the mold in the agarose gel enable the embryos to self align.





Xenograft and Larval Injection–The sloped ridges make perfect angles in the agarose gel, which then makes it easier to do microinjections in the larvae.





Transplantation–Increase the speed of doing microinjections. Simply turn the petri dish as you are injecting.





Standard Microinjection-Use for blastomere transplantation.





ORDERING INFORMATION

Z-MOLDS Microinjection & Transplantation Molds

References

Kitambi SS, Toledo EM, Usoskin D, Wee S, Harisankar A, Svensson R, Sigmundsson, K, Kalderén C, Niklasson M, Kundu S, Aranda S, Westermark B, Uhrbom L, Andäng M, Damberg P, Nelander S, Arenas E, Artursson P, Walfridsson J, Forsberg Nilsson K, Hammarström LG, Ernfors P. Vulnerability of glioblastoma cells to catastrophic vacuolization and death induced by a small molecule. *Cell*. 2014 Apr 10;157(2):313-28.

Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

Features

- · Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

Benefits


- Superior pricing
- Most glass orders ship within 48 hours

Applications

- Microinjection
- Electrophysiology
- Patch clamp
- · Fluid Handling

Fire Polishing

Fire-polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



Making Uniform, Reproducible Microelectrodes


Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

Filament Glass Capillaries


Single barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.


Thin Wall Glass Capillaries

Thin wall single barrel capillaries are offered both with or without inner filaments.









OPPEDING INFORMATION

		OF	RDERIN	G INFOR	RMATION		
	Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	ltem
	3 in. (76 mm)	1.0	0.58	V		500	1B100F-3
	3 in. (76 mm)	1.0	0.58			500	1B100-3
	3 in. (76 mm)	1.2	0.68	V		350	1B120F-3
SS	3 in. (76 mm)	1.2	0.68			350	1B120-3
<u>ila</u>	3 in. (76 mm)	1.5	0.84	V		225	1B150F-3
e	3 in. (76 mm)	1.5	0.84		~	300	1B150-3
cat	4 in. (100 mm)	1.0	0.58	V	V	500	1B100F-4
=	4 in. (100 mm)	1.0	0.58		~	500	1B100-4
ĕ	4 in. (100 mm)	1.2	0.68	V	V	400	1B120F-4
8	4 in. (100 mm)	1.2	0.68			350	1B120-4
5	4 in. (100 mm)	1.5	0.84	~	V	300	1B150F-4
qa	4 in. (100 mm)	1.5	0.84		V	300	1B150-4
an	4 in. (100 mm)	2.0	1.12	V		125	1B200F-4
Single Barrel Standard Borosilicate Glass	4 in. (100 mm)	2.0	1.12		V	200	1B200-4
ē	6 in. (152 mm)	1.0	0.58	V		500	1B100F-6
3ar	6 in. (152 mm)	1.0	0.58			500	1B100-6
e	6 in. (152 mm)	1.2	0.68	V		350	1B120F-6
쪌	6 in. (152 mm)	1.2	0.68			350	1B120-6
S	6 in. (152 mm)	1.5	0.84	~		225	1B150F-6
	6 in. (152 mm)	1.5	0.84			225	1B150-6
	6 in. (152 mm)	2.0	1.12	~		125	1B200F-6
	6 in. (152 mm)	2.0	1.12			125	1B200-6
	3 in. (76 mm)	1.0	0.75	V		500	TW100F-3
	3 in. (76 mm)	1.0	0.75			500	TW100-3
_	3 in. (76 mm)	1.2	0.90	V	V	400	TW120F-3
arc	3 in. (76 mm)	1.2	0.90			350	TW120-3
D	3 in. (76 mm)	1.5	1.12	~		225	TW150F-3
Sta.	3 in. (76 mm)	1.5	1.12		~	300	TW150-3
()	4 in. (100 mm)	1.0	0.75	~		500	TW100F-4
Ĕ	4 in. (100 mm)	1.0	0.75		~	500	TW100-4
-Ba	4 in. (100 mm)	1.2	0.90	V		350	TW120F-4
<u>6</u>	4 in. (100 mm)	1.2	0.90			350	TW120-4
ing	4 in. (100 mm)	1.5	1.12	V		225	TW150F-4
=	4 in. (100 mm)	1.5	1.12		V	300	TW150-4
Thin-Wall Single-Barrel Standard	6 in. (152 mm)	1.0	0.75	V		500	TW100F-6
ī	6 in. (152 mm)	1.0	0.75		V	500	TW100-6
Ē	6 in. (152 mm)	1.2	0.90	V	V	400	TW120F-6
	6 in. (152 mm)	1.2	0.90			350	TW120-6
	6 in. (152 mm)	1.5	1.12	V		225	TW150F-6
	6 in. (152 mm)	1.5	1.12		V	300	TW150-6

Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330.

Microelectrode Puller

Cost-effective and compact micropipette puller

Features

- · Microprocessor controlled
- · Program sequences up to four steps
- Store up to 95 programs in memory for easy recall
- · Multiple factory programs installed

Benefits

- Tempered glass cover to reduce the effects of humidity on puller reproducibility
- Switchable power supply ensures that line voltage fluctuations don't affect reproducibility

Applications

Pull your own microelectrodes and micropipettes



horizontal puller for making glass micropipettes or microelectrodes used in microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heat, pulling force, pull length and cooling time.

This puller is a reasonably priced, compact, versatile and a reliable workhorse. The microprocessor, combined with the LCD display, makes the **PUL-1000** easy to use.

Tempered Glass Cover

The cover of the pulling chamber is made with tempered glass to minimize the temperature effect on the reproducibility of pulled pipettes.

Switchable Power Supply

PUL-1000 has a high quality switching power supply for use anywhere in the world without worry about the line voltage differences. Pulling reproducibility is unaffected by line voltage fluctuation. Heating voltage can be controlled to within 0.1% accuracy even when line voltage fluctuates from 90 to 240 VAC

References

Plautz, C. Z., Williams, H. C., & Grainger, R. M. (2016). Functional Cloning Using a Xenopus Oocyte Expression System. *Journal of Visualized Experiments* (107), e53518–e53518. http://doi.org/10.3791/53518



Komarova, Y., Peloquin, J., & Borisy, G. (2011). Components of a microinjection system. *Cold Spring Harbor Protocols*, 2011(8), 935–9. http://doi.org/10.1101/pdb.ip27

PUL-1000 SPECIFICATIONS

HEATER ELEMENT Platinum/Iridium **PULLING FORCE** Solenoid, adjustable TAPER LENGTH 1-10 mm CAPILLARY OD RANGE 1.0-2.0 mm MAXIMUM CAPILLARY LENGTH 170 mm MINIMUM CAPILLARY LENGTH 55 mm MEMORY SETS 95 **AUTO SHUT-OFF TIME** POWER

 POWER
 90-240VAC, 50/60 Hz

 DIMENSIONS
 34 x 24 x 12 cm

 SHIPPING WEIGHT
 16 lb

ORDERING INFORMATION

PUL-1000 Micopipette Puller

System includes PUL-1000 Puller; TW100-4 Thin Wall Glass Capillaries (package of 500); MF34G-5 package of MicroFil

13834 Filament, Platinum/Iridium, 2.5mm wide, 2.5mm Square Box 504951 Filament, Platinum/Iridium, 4.5mm wide, Trough

Also see these other pullers available from WPI





SU-P97 See page 111.

SU-P1000 See page 110.

Fluid Handling



A large variety of pumps, from picoliter to microliter


We offer a large variety of pumps for everything from microfluidics to general fluid handling applications. Our peristaltic pumps are easy to setup and clean, offer continuous flow with virtually "infinite" volume (depending only on the capacity of your source), require no contact with metal or the pump and are good for large volume pumping. Our syringe pumps provide accurate volume control and are an excellent choice for lower volume applications. Our pneumatic PicoPumps are designed for delivery of very small volumes.

Which Pump is Right for Me?

MiniStar™ Peripro-2HS Peripro-4HS Peripro-4LS Peripro-8LS MPS2 AL-1000	Up to 250 μL/min. LAB 0 0.73 μL/hr to 1699 mL/hr	1 2 4 4 8 AND WHC 8	Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges
Peripro-2HS Peripro-4HS Peripro-4LS Peripro-8LS MPS2 AL-1000	0.8 - 300 mL/min 0.8 - 300 mL/min 0.01-80 mL/min 0.01-80 mL/min SINGLE-CHANNEL Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	1 2 4 4 8 AND WHC 8	Compact design, remote control Calibrated output, replaceable tubing cartridges OLE-CELL SOLUTION EXCHANGE
Peripro-2HS Peripro-4HS Peripro-4LS Peripro-8LS MPS2 AL-1000	0.8 - 300 mL/min 0.8 - 300 mL/min 0.01-80 mL/min 0.01-80 mL/min SINGLE-CHANNEL Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	2 4 4 8 AND WHC 8	Calibrated output, replaceable tubing cartridges DLE-CELL SOLUTION EXCHANGE
Peripro-4HS Peripro-4LS Peripro-8LS MPS2 AL-1000	0.8 - 300 mL/min 0.01-80 mL/min 0.01-80 mL/min SINGLE-CHANNEL Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	4 4 8 AND WHC 8	Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges DLE-CELL SOLUTION EXCHANGE
Peripro-4LS Peripro-8LS MPS2 AL-1000	0.01-80 mL/min 0.01-80 mL/min SINGLE-CHANNEL Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	4 8 AND WHC 8	Calibrated output, replaceable tubing cartridges Calibrated output, replaceable tubing cartridges PLE-CELL SOLUTION EXCHANGE
Peripro-8LS MPS2 AL-1000	0.01-80 mL/min SINGLE-CHANNEL Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	8 AND WHC 8	Calibrated output, replaceable tubing cartridges DLE-CELL SOLUTION EXCHANGE
MPS2 AL-1000	SINGLE-CHANNEL A Up to 250 μL/min. LAB 0.73 μL/hr to 1699 mL/hr	AND WHC	LE-CELL SOLUTION EXCHANGE
AL-1000	Up to 250 μL/min. LAB 0 0.73 μL/hr to 1699 mL/hr	8	
AL-1000	LAB 0.73 µL/hr to 1699 mL/hr		Programmable control low dead volume
	0.73 µL/hr to 1699 mL/hr		SYRINGE PUMPS
	•	1	
	U /2 U /br to 1600 ml /br	2	Push/pull Push/pull (2 networked pumps)
AL-2000	0.73 µL/hr to 1699 mL/hr	1	
AL-1000 HP	1.459 µL/hr - 6120 mL/hr		Infuse/withdraw, High Pressure
AL-4000	1.459 µL/hr to 6120 µL/hr	2	Infuse/withdraw
AL-6000	0.73 µL/hr to 1699 mL/hr	6	Infuse/withdraw
AL-8000	0.454 µL/hr to 163 mL/hr	8	Infuse/withdraw
SPLG100	1.26 pL/min to 88.32 mL/min	1	Infuse only
SPLG101	1.26 pL/min to 25.99 mL/min	2	Infuse only
SPLG110	1.26 pL/min to 88.28 mL/min	1	Infuse/withdraw
SPLG200	0.5 pL/min to 220.97 mL/min	2	Infuse only
SPLG210	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw
SPLG212	0.5 pL/min to 220.97 mL/min	2	Infuse/withdraw programmable
SPLG270	0.5 pL/min to 220.97 mL/min	2+2	Push-pull
SPLG272	0.5 pL/min to 220.97 mL/min	2+2	Push-pull programmable
SP100i	0.0001-519 mL/hr	1	Basic single channel
SP101i	0.001 μL/hr – 35 mL/min	2	Micro dialysis application
SP120p	0.1 μL/hr – 127 mL/hr	1+1	Push pull, single cycle
SP200i	0.001 μL/hr - 145 mL/min	2	RS232 TTL/foot switch
SP210c	0.001 μL/hr - 86 mL/min	2+2	RS232 push pull, continuous
SP210iw	0.001 μL/hr - 145 mL/min	2	RS232 infuse/withdraw
SP220i	0.001 μL/hr - 21 mL/min	10	RS232 infuse only
SP230iw	0.001 μL/hr - 21 mL/min	10	RS232 infuse/withdraw
SP250i	0.001 μL/hr - 21 mL/min	4	RS232 infuse only
SP260p	0.001 μL/hr - 86 mL/min	2+2	RS232 push pull, single cycle
	MICRO SYRIN	GE PUMP	/ STEREOTAXIC INJECTION
UMP3T-1	0.03 nL/min - 10 μL/sec	1	Ultra micro infuse/withdraw RS232
MMP	Manual 100 μL-1 mL syringe	1	Manual
DMP	Manual 100 μL-1 mL syringe	1	Digital readout micrometer
		MICRO	INJECTION
PV820	Injected volumes from pL to nL	1	Injection pressure and holding pressure
PV830	Injected volumes from pL to nL	1	Injection pressure and holding pressure and vacuum
NANOLITER2010	Bolus, 2.3-69 mL/Injection	1	Oocyte injector, infuse only
Microinjection Sy	ystems (Zebrafish, C. Elegans, Droso _l	ohila, Xenopus	s oocytes)
		MICRO	FLUIDICS
ExiGo	50 nL/min - 10 mL/min	1	Infuse only, feedback via integrated flow sensor, includes iPad mini which can control up to four pumps
Mirus	100nL/min – 10 mL/min ±1%	8	Microchip perfusion. Infuse only, reversible flow, ~600 μL dead volume, PC control
Kima	15 – 35 mL/hr ±4%	1	Microchip perfusion. Infuse only, recirculating pump controlled by iPod Touch, wi-fi communication, <300 µL dead volume
Unigo	1 μL/min–1 mL/min; unidirectional (push)	1	A precision, microfluidic, single-channel, pressure pump. Can expand up to 4 channels.
		SUPPLIE	S & TOOLS

Multichannel Perfusion System

For single ion channel and whole-cell solution exchange




MPS-2 is a programmable 8-channel perfusion system designed for single channel and whole-cell patch preparations. Offering the best combination of performance and value, the **MPS-2** incorporates the same high quality solenoid valves found on similar but much more expensive systems. Unlike other perfusion systems on the market, which often compromise performance to fit every possible application, the **MPS-2** is the only perfusion system designed and optimized specifically for single-channel and whole-cell patch perfusion applications.

Options for Manual and Automatic Control

The system can be controlled manually via membrane switches on the front panel or through a PC. Two different manual control modes are offered. One controls each channel independently and the other mode allows you to assign a master channel that will keep the system flow when all other channels are switched off. User-friendly graphic timing software is included, and the programmed perfusion sequence can be started by computer, a patch clamp amplifier or other external trigger, or manually.

Organized, Color-coded, Ribbon Tubing

The perfusion fluid flows through specially designed color-coded polyurethane ribbon style tubing. The color-coding allows you to easily trace each channel for diagnostic checks or set up, and the ribbon style of tubing keeps the system neat and organized. Unlike PVC based tubing, polyurethane tubing contains no plasticizer, which can cause contamination. The tubing ribbon is designed as an economical disposable item, which is often critical when cleanliness is needed.



Micromanifold closeup: Fluid-filled passages are shown in magenta.

Low Flow Resistance and Dead Volume

The most unique feature of the MPS-2 is its perfusion micromanifold. Using the latest microfluidic techniques, the injection molded micromanifold provides the least flow resistance and dead volume of any product on the market. The flow channel inner diameter is approximately 1.0 mm, except for the last 5 mm before the junction point. This design allows a fast flow rate without using a pressured system. Small channels and a unique design at the merging point further reduce the chance of cross contamination. Dead volume is less than 100 nL at junction.

MPS-2 SPECIFICATIONS

CHANNELS 8
VALVE RESPONSE TIME 2 ms

VALVE CONTROL USB, TTL, external start via software

SYRINGE RESERVOIR VOLUME 10 mL MANIFOLD 8 to 1

DEAD VOLUME

OUTPUT TUBING ID 250 µm and 100 µm. 100 µm ID tip, 8 µL/min.

MAXIMUM FLOW RATES (with 50 mm output tubing, gravity fed)

Tubing cap: 7.85 nL/mm
250 µm ID tip, 250-500 µL/min.
Tubing cap: 49.1 nL/mm

< 100 nL excluding the single outlet tubing

ORDERING INFORMATION

MPS-2 Multichannel Perfusion System & Control Software

System includes: stand, PS-2 controller, valve console, syringe holder, power cord, USB cable, DB9-to-BNC 8-cable assembly, (2) 1A fuses, (10) 10 mL syringes, (10) 3-way stopcocks, (10) Luer fitting with barb for 1/16-in. ID tubing, color coded polyurethane tubing ribbon (5 ft), tubing for making micromanifold cleaning adaptor (3 in.), micromanifold holding rod, (2) micromanifolds with 100 µm id tip, (2) micromanifolds with 250 µm id tip and installation software

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

	7.0000000000000000000000000000000000000
502109-15	Color-coded PU Tubing, 1/16" ID x 8 Channels, 15 ft
502110	Micromanifold, 100 μm ID tip, 2 pcs/pk
502125	Micromanifold, 250 μm ID tip, 2 pcs/pk
Cnac	if I line voltage and Micromanifold tin OD when ordering

Specify line voltage and Micromanifold tip OD when ordering.

Affordable High Performance Peristaltic Pump

Digital Peri-Star™ Pro peristaltic pump is affordably priced!

Features

- · Available in 2-, 4- and 8-channel versions
- Display either rotation speed (RPM) or flow rate (mL/min)
- Wide flow range: 0.01 280 mL/min
- Accuracy of flow rate: 0.5% using self calibration function
- Accuracy of speed: 0.1 rpm
- Large backlit digital LCD display
- Programmable for all tubing sizes between 0.8 mm and 6.4 mm ID
- · Easy and fast tubing replacement using snap-on cartridges
- Membrane keypad allows easy programming while protecting controls from fluid entry
- · Actively driven rollers by planetary gears for long lasting tubing life

Benefits

- · Backlit display with water resistant cover
- · Easy setup and calibration with simple instructions on screen
- Both high and low flow pumps employ planetary gears for minimal pulsations and greater accuracy

Applications

- · Chemical transfer
- · Pharmaceutical processing

Peri-Star™ Pro peristaltic pumps provide accurate and precise pumping with convenience and versatility. Peri-Star Pro can be run in either flow rate mode (mL/min) or rotation speed mode (rpm). For good laboratory practice, pumps must be calibrated after changing the tubing and solution. You can easily calibrate Peri-Star Pro to deliver flow as accurate as 0.5% in a wide flow range from 0.01 mL/min to 280 mL/min. Under rotation speed mode, the digitally controlled stepping motor provides accurate and reproducible operation with 0.1% rpm both forward and in reverse.

Backlit display with water resistant cover

Large backlit digital LCD display provides readouts of rotation direction, flow rate or rotation speed, tubing ID, drive status and remote control mode simultaneously. Water resistant membrane keypad allows easy programming while protecting LCD display and controls from fluid entry.

Easy set up and calibrate the pump

Built-in Human Machine Interface (HMI) with screen instructions in plain English steps you through initial setup, calibration and operating procedures. The user-friendly interface reduces the need to frequently check the printed manual for instruction and reference.

Unique planetary gears for accuracy

Peri-Star Pro is available in two versions: a 4-roller version for high flow and an 8-roller version for lower volumes which provides high pressure with minimal pulsations.

A unique planetary gear design with eight actively driven rollers (four rollers for higher flow rate model), together with independent tubing compression fine adjustment, enables high flow accuracy and prolongs tubing life. Snap-on cartridges allow tubing to be changed quickly without cross contamination of solutions.



PERI-STAR PRO SPECIFICATIONS

	Peri-Star Pro 2H / 4H (High Rate)	Peri-Star Pro 4L / 8L (Low Rate)
NUMBER OF ROLLERS	4	8
NUMBER OF CHANNELS	2–4	4–8
ROTOR SPEED RANGE	1–100 rpm	1–100 rpm
FLUID FLOW RANGE	0.8–280 mL/min #17 Tubing: 3.5–280 mL/min	0.01-80 mL/min #14 Tubing: 0.2-18 mL/min
TUBING RANGE	3.1-6.4 mm ID	0.5-2.4 mm ID
SELF-CALIBRATION	Yes	Yes
WORKING ENVIRONMENT	0-45°C, Humidity < 80%	0-45°C, Humidity < 80%
POWER	110 V or 220 V AC, 50 - 60 Hz	110 V or 220 V AC, 50 - 60 Hz
DIMENSIONS	190 x 162 x 275 mm	190 x 162 x 275 mm
SHIPPING WEIGHT	11 lb. (5 kg)	11 lb. (5 kg)

ORDERING INFORMATION

PERI-STAR PRO PUMPS

PERIPRO-2HS	Peri-Star™ Pro, 2-channel, High Rate, Large Tubing (110-220V)
PERIPRO-4HS	Peri-Star™ Pro, 4-channel, High Rate, Large Tubing (110-220V)
PERIPRO-4LS	Peri-Star™ Pro, 4-channel, Low Rate, Small Tubing (110-220V)
PERIPRO-8LS	Peri-Star™ Pro, 8-channel, Low Rate, Small Tubing (110-220V)


OPTIONAL ACCESSORIES/REPLACEMENT PARTS

503049	Replacement Tubing Cartridge, Large
503050	Replacement Tubing Cartridge, Small
503022	Replacement Silicone Tubing, 1 m, 1.6 mm I.D., #14, with stops
503023	Replacement Silicone Tubing, 1 m, 6.4 mm I.D., #17
503120	TTL Control Module

Pump head for PERIPRO-4HS (4 channels, 4 rollers)







Mini Peristaltic Pump

Lightweight, portable Mini★Star™

Features

- Flow Range: 0.06-14.0 mL/min
- Speed: 1-50.0 rpm, forward/reverse
- Stand and clamp shown in the image are included
- · Includes tubing

Benefits

- · Remote control (wired)
- Small footprint
- Quiet operation
- Compact and easily transported

Applications · Small volume

chemical transfer

TUBING (Two-stop Silicone)

- Pharmaceutical processing
- · Sample perfusion

This compact, lightweight peristaltic pump fits just about anywhere. It can be mounted directly on the bench, in a regular rack or to a post. The speed can be adjusted from 1 to 50 rpm. With recommended silicone tubing, the volume can be set from 0.06 to 14.0 mL/min. The MiniStar™ also features a hand held remote control that allows you to start and stop the pump, purge or adjust its speed and direction.

MINISTAR



CHANNEL **SPEED** 1-50.0 rpm, forward/reverse FLOW RANGE 0.06~14.0 mL/min RESOLUTION 1 rpm (0.1 rpm computer control) SPEED CONTROL Remote control DISPLAY Indicators for status and speed **POWER** 12 V DC (110/220 VAC adapter incl.) WORKING CONDITION Temperature 0-40°C, humidity < 80%

0.8~1.0 mm Wall Thickness Outer Diameter ≤ 4.8 mm

135×72×72 mm (L×W×H) DIMENSION OF DRIVER DIMENSION OF REMOTE CONTROL 105×50×16 mm (L×W×H) WEIGHT OF DRIVER 0.5 kg

	ORDERING INFORMATION
MINISTAR	Miniature Peristaltic Pump, 1-channel
504011	MiniStar™ and Stand (as pictured above)
503120	TTL Control Module
503121	Silicone Tubing with stops, 2.4mm ID x 0.8mm wall x 1 m (5-pk)
503122	Silicone Tubing with stops, 1 mm ID x 1 mm wall x 1 m (5-pk)

Microinjection Syringe Pump

Delivering picoliter volumes precisely

Features

- Graphic display with SMARTouch touch screen controller for "intelligent", easy to use interface controlling up to four syringe pumps
- · Splash proof touch screen
- · User configurable mounting bar
- · Dual mode motor drive
- Compatible with all UMP, UMP2 and UMP3 pumps
- · Optional foot switch available
- 5-Digit display

Benefits

- · Accepts a wide variety of microinjection syringes
- · Manual or automated injections
- · Quiet operation for electrophysiology recordings
- Mounts directly on micromanipulator or stereotaxic frame
- · Nominal injections down to 1.0 nL
- · Rapid setup with intuitive touchscreen controller

Applications

- Microinjection
- Neuroscience
- Microfluidics
- · Micro delivery of biochemical agents or dyes

The UltraMicroPump 3 (**UMP3**) is a versatile pump which uses micro syringes to deliver picoliter to milliliter volumes. The pump is optimum for applications that require injections of precise and small amounts of liquid. With its touchscreen controller, UMP3 can displace as little as 0.53 μ L/step (using 10 μ L syringe with 60 mm scale length).

The new SMARTouch™ controller for the UltraMicroPump features Patent Pending technology which includes:

- Total system calibration Calibrate the syringe and the controller together as a system. This feature eliminates the variability of the syringes and delivers the calibrated volume.
- Smart smoothness The controller can be set to automatically adjust microstepping according to the injection rate to deliver the smoothest flow.
- User defined travel limits Set the limits for a specific syringe in the software. This prevents the pump from over-driving the plunger into the syringe, potentially causing syringe breakage.

The MICRO2T SMARTouch™ controller is feature rich. All operations are controlled through interactive touch screen. It has a graphical indication of the flow and the volume remaining in the syringe. It offers automatic end stop detection that is dependent on the syringe volume. You can control two pumps independently from one controller with its dual display. It also has automatic pump detection and a Pause/Resume feature that allows dosing during infusion/withdrawal The volume accumulated is displayed on screen, as well as the percentage of volume left in the syringe. The SMARTouch controller is fully compatible with all earlier versions of the UltraMicroPump.

Low Fluid Dead Volume

Syringes may be filled externally and then inserted into the pump or filled while mounted in the pump. Fluids injected or withdrawn are held entirely within the micro syringe to maintain a low fluid dead volume.

Flexibility in Mounting

For positioning, the UMP3 may be attached to any of several WPI micropositioners such as the **M3301** (manual), **SU-QUAD** or **SU-TRIO** (motorized) or any manual stereotaxic manipulator.

Rapid Setup with Intuitive Touchscreen Controller

An Integral component in the UMP3 system is the SmarTouch touchscreen controller, which provides an "intelligent" and easy-to-use interface to up to four (or two) syringe pumps. Operating parameters are



set with the touchscreen panel. You can save your parameters for instant recall. An optional foot switch offers "hands free" start/stop operation.

Computer Control—A USB port on the rear of the controller can be used to connect it to a computer for scripted protocols.

NOTE: UMP3 accepts glass syringes with barrel diameters from $5.5\ to$ 9 mm.

References

Augestad, I. L., Nyman, A. K. G., Costa, A. I., Barnett, S. C., Sandvig, A., Håberg, A. K., & Sandvig, I. (2017). Effects of Neural Stem Cell and Olfactory Ensheathing Cell Co-transplants on Tissue Remodelling After Transient Focal Cerebral Ischemia in the Adult Rat. *Neurochemical Research*, 1–11. http://doi.org/10.1007/s11064-016-2098-3

Eleftheriadou, I., Dieringer, M., & Poh, X. (2017). Selective transduction of astrocytic and neuronal CNS subpopulations by lentiviral vectors pseudotyped with Chikungunya virus envelope. *Biomaterials*. Retrieved from http://www.sciencedirect.com/science/article/pii/S0142961217300376

Job, M. O., & Kuhar, M. J. (2017). Cart Peptide in the Nucleus Accumbens Regulates Psychostimulants: Correlations Between Psychostimulant and Cart Peptide Effects. *Neuroscience*. http://doi.org/10.1016/j. neuroscience.2017.02.012

Mosberger, A. C., et al.. (2017). Axotomized Corticospinal Neurons Increase Supra-Lesional Innervation and Remain Crucial for Skilled Reaching after Bilateral Pyramidotomy. *Cerebral Cortex* (New York, 1991). http://doi.org/10.1093/cercor/bhw405



The SMARTouch™ controller can control up to two UMP3 pumps.



UMP3 shown mounted on a standard micromanipulator (not included) and using the new SMARTouch controller.

ULTRAMICROPUMP SPECIFICATIONS

(based on 10 µL syringe)

NORMAL MODE

TRAVEL
MINIMUM DISPENSING VOLUME
LINEAR MOTION PER STEP
WEIGHT
MOUNTING ROD DIAMETERS
MAINS POWER SUPPLY
DIMENSIONS

MAXIMUM SYRINGE SIZE

62 mm 0.58 nL / step (10 μL syringe) 3.175 μm/half step 325 g (11.5 oz) 7.9 mm (0.31 in.) 90-264VAC @ 47-63Hz Ø 32 mm x 190 mm (Ø 1.3 in. x 7.5 in.) 1 mL

MICROSTEPPING MODE

Precision is increased eight-fold

	ORDERING INFORMATION
UMP3T-1	UltraMicroPump III (one) and Micro2T Controller
UMP3T-2	UltraMicroPump III (two) and Micro2T Controller
UMP3-3	UltraMicroPump III (three) and Micro4 Controller
UMP3-4	UltraMicroPump III (four) and Micro4 Controller
UMP3	UltraMicroPump III (without controller)
MICRO2T	SMARTouch Controller, Two-Channel
OPTIONAL	ACCESSORIES/REPLACEMENT PARTS
13142	Foot Switch for MICRO2T*
502201	V-clamp for Stereotaxic Frame
503301	Extension Cable, miniDIN (male-female) 10 ft
503207	Small Base Stand and Clamps
*121/12 ic t	ha foot padal for use with the MICPORT ONLY It is not compat

^{*13142} is the foot pedal for use with the MICRO2T ONLY. It is not compatible with the old MICRO4 controller, which uses the 15867 foot pedal. The two foot pedals are not cross-compatible.

Microvolume Syringes



Syringes with Luer Fitting (no needle)								
Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2		
ILS005LT	5 μL	ILS 5 µL Gas-tight Luer tip	6.5 mm	54.1 mm	Υ	Υ		
ILS010LT	10 µL	ILS 10 μL Gas-tight Luer tip	6.5 mm	54.1 mm	Υ	Υ		
ILS025LT	25 µL	ILS 25 μL Gas-tight Luer tip	8.0 mm	60 mm	Υ	Υ		
SGE050TLL	50 µL	SGE 50 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	Υ		
SGE100TLL	100 µL	SGE 100 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	Υ		
SGE250TLL	250 μL	SGE 250 µL Gas-tight Teflon Luer Lock	8.0 mm	60 mm	Υ	N		

ORDERING INFORMATION

Syringes with Replaceable Beveled Needles						
Order No.	Volume	Description	O.D.	SCALE LENGTH	UMP3	UMP2
SGE0005RN*	0.5 µL	SGE 0.5 μL 23 ga (0.63 mm), 70 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE001RN*	1.0 µL	SGE 1.0 µL 26 ga (0.47 mm), 70 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE005RN	5 μL	SGE 5 µL 23 ga (0.63 mm), 70 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE010RNS	10 μL	SGE 10 μL 26 ga (0.47 mm), 50 mm long needle	8.0 mm	54.1 mm	Υ	Υ
SGE025RN	25 µL	SGE 25 μL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ
SGE050RN	50 μL	SGE 50 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ
SGE100RN	100 μL	SGE 100 µL 25 ga (0.50 mm), 50 mm long needle	8.0 mm	60 mm	Υ	Υ

^{*} The capacity of this syringe is so small that the entire sample is contained within the needle. The plunger extends to the tip of the needle, displacing the full sample during injection — which gives the syringe zero dead volume.

SGE and ILS are respective trademarks of Scientific Glass Engineering and Innovative Labor Systeme.

Replacement Needles

RN0005	For Syringe SGE0005RN, 23 ga (0.63 mm) 70 mm long
RN001	For Syringe SGE001RN, 26 ga (0.47 mm) 70 mm long
RN005	For Syringe SGE005RN, 23 ga (0.63 mm) 50 mm long
RN010	For Syringe SGE010RN(S), 26 ga (0.47 mm) 50 mm long, 5-pack
RN025	For Syringes SGE025RN, SGE050RN, SGE0100RN, 26 ga (0.47 mm) 50 mm long, 5-pack

Sub-Microliter Injection System

Includes the smallest dead volume injection when the 10 μ L syringe is used with WPI needles 34-36g

Features

- · The smallest dead volume injection syringe
- Various needle sizes available from 26 ga. and 33-36 ga.
- · Blunt or sharp needles
- · Compatible with WPI's UMP3 microinjection system

Benefits

- Low dead volume (0.5 μL or less)
- · Switching the syringe tip during an experiment is easy
- · Variety of tips.

Applications

- Animal research
- · Capillary electrophoresis
- · Versatile research applications RPE and IO Kits

NanoFil™ is a specially designed 10 µL syringe developed in response to customer requests for improved microinjection in mice and other small animals. It makes quantitative nanoliter injection much easier and more accurate than any other method currently in use.

Low Dead Volume

NanoFil's low dead volume eliminates the need for oil backfilling, a messy process which risks contamination of the injected sample. Injection is now simpler, and less messy, and there is no possibility of oil contamination in critical applications such as ophthalmology research (see the Retinal Pigment Epithelial (RPE) and Intra Ocular (IO) injection kits listed below).

Easily Switch Syringe Tip

When the inner tip diameter of a conventional syringe is reduced to less than 100 μm , it is very difficult to front fill the solution at a reasonable speed. NanoFil solves this problem by using a tip coupling mechanism that makes it possible to change the syringe tip during the experiment. Simply load the sample using a larger tip, such as the 26 gauge needle provided with the syringe, and then replace it with a micro tip for sample injection. On a conventional 10 μL syringe, a solid ring or bushing is permanently bonded to the tubing. Replacing the tip in the middle of the experiment is not practical. With NanoFil, tips can be exchanged by a simple twist of the brass lock, gently pulling out the tip, and replacing with the desired new tip.

Holds Metal Tips and Quartz Tubing

To secure the tip, NanoFil uses an olive-shaped silicon gasket that is similar to, but much sturdier than, some of the microelectrode holders used for electrophysiology recording. The silicone gasket makes it possible to hold not only metal needles but also **Silflex** tubing. Many types of tubing can be easily connected to the syringe as long as the outer diameter (OD) is close to, but not more than, the barrel inner diameter (ID) of 460 µm. Flexible quartz capillaries used in Gas Chromotography (GC) and Capillary Electrophoresis (CE) can also be easily coupled to the syringe.

Variety of Tips

Specially designed needles as small as 36 gauge (110 μ m OD) are offered in both blunt and beveled styles. Our studies have shown that these needles will cause less trauma to the tissue. NanoFil has a unique coupling mechanism that allows many different forms of small tubing and tips to be coupled with the syringe barrel.

Selecting the correct tip for your application

The replaceable needles used with NanoFil are available with either blunt or beveled tips. The blunt tip is used for injection into soft tissue and when a uniform solution distribution is needed. The beveled style is



ferrule: Ø 8 mm

used for applications that involve the penetration of a tough tissue. One of the main factors that can affect the resolution and accuracy of a microinjection in the low nanoliter range is diffusion from the needle opening. When the tip ID is equal to or larger than 100 μm , the error caused by tip diffusion can be in the nanoliter range level [(100 micron)^3 = 1 nanoliter]. With a 36 gauge needle installed on the NanoFil, the error caused by diffusion will be reduced, making accurate injection of a nanoliter possible.

All of WPI's beveled needles have a unique 25° tri-surface bevel that is optimized for microinjection. A 10 degree single-surface beveled tip penetrates better than one with a 25° angle, however the distance between the upper opening to the tip (the dimension "F" in the table on page 59) is longer. As a result, it requires deeper penetration of the tip. This can be an issue for microinjection into very small areas where the dimensions of the anatomy can't accomodate the required depth of insertion.. WPI's unique 25° beveled tip reduces the required insertion depth by incorporating two extra beveled surfaces. The edge of a single surface beveled tip is actually a blade instead of a point. It dulls very quickly. In contrast, the tri-surfaced tip has a sharp point. It penetrates more easily and is more durable. Our tests show that our 33 gauge, 25 degree beveled tip penetrates easier and lasts longer than other manufacturers' 33 gauge, 10 degree single beveled tips. With a 35 gauge tri-surface beveled tip, the resistance to the penetration becomes even less. Each of our needles undergo a penetration test before leaving the factory to guarantee the best results for our customers.

Available Tips

33 gauge: This tip is similar to Hamilton's 7762 and 7803 series removable needles in both tip length and outer diameter. However, our beveled tip version is shorter, more durable, and penetrates better due to the special tri-surface grinding technique. In the past, 33 gauge tips were the smallest size sold by other manufacturers and were frequently cited in literature. However, our new 35 gauge tip is much better for injections involving small animals, especially mice. Compared with Hamilton's 33 gauge, 10 degree beveled tip, our 35 gauge 25 degree beveled tip can reduce the depth of penetration by almost 80%. The distance between the tip and the upper rim of the opening (dimension F on the drawing) is 348 μ m for the 33 gauge tip. The distance for our 35 gauge tip is only 230 μ m. In addition, the smaller tip size significantly reduces the required penetration force. In nearly all applications, a 33 gauge tip can be replaced with our 35 gauge tip and produce better results.

34 gauge

This is a transitional size between the 33 gauge and 35 gauge. If the 35 gauge is too weak and the 33 gauge is too large, this makes a good alternative.

35 gauge

This was the most popular and preferred tip of most scientists during our field trial. The combination of its strength, length, durability and clogging resistance creates a balance with very little compromising of the individual properties. It is much smaller than the 33 gauge tip offered by other manufacturers. It is only slightly larger than the 36 gauge tip but is much stronger and less likely to be clogged. Samples can be directly loaded with this tip. Its 5 mm length is sufficient enough for almost all injection applications in mice.

36 gauge

This is the smallest tip that is commercially available. The tip is so small that it can be inserted into the opening of the 33 gauge needle tip. Because this is pushing the limits of what current technology can



produce, there are some limitations to consider before using it. Its thin diameter makes it necessary to limit its length to 2.5 to 3 mm in order to maintain a usable strength. Since the tip ID is in the 25 to 50 μm range, it is very easily clogged. Therefore, only well filtered solutions can be used. Depending on the viscosity of the sample, you might also need to pre-load the syringe with a regular tip before switching to this tip for injection. We recommend using the 35 gauge tip instead of the 36 gauge unless it is absolutely necessary.

Flexible Quartz Tubing

The flexible quartz tubing tip is made of 160 µm OD polyimide coated quartz tubing with a special adapter sleeve mounted at the end. It is designed for filling glass capillary electrodes or pipettes, just like WPI's traditional MF34G MicroFil. However, unlike the traditional MicroFil,

which has about 50 µL of dead volume in its Luer hub, the dead volume of this tip is less than 0.589 µL. It is useful for loading electrodes with solutions that have a limited volume or are too expensive to waste.

References

Arriaga, G., Macopson, J. J., & Jarvis, E. D. (2015). Transsynaptic Tracing from Peripheral Targets with Pseudorabies Virus Followed by Cholera Toxin and Biotinylated Dextran Amines Double Labeling. *Journal of Visualized Experiments*, (103), e50672–e50672. http://doi.org/10.3791/50672

Delotterie, D. F., Mathis, C., Cassel, J.-C., Rosenbrock, H., Dorner-Ciossek, C., & Marti, A. (2015). Touchscreen tasks in mice to demonstrate differences between hippocampal and striatal functions. *Neurobiology of Learning and Memory*, 120, 16–27. http://doi.org/10.1016/j.nlm.2015.02.007

Han, Z., Banworth, M. J., Makkia, R., Conley, S. M., Al-Ubaidi, M. R., Cooper, M. J., & Naash, M. I. (2015). Genomic DNA nanoparticles rescue rhodopsin-associated retinitis pigmentosa phenotype. FASEB Journal: Official Publication of the Federation of American Societies for Experimental Biology, 29(6), 2535–44. http://doi.org/10.1096/fj.15-270363

Silflex

NF26BV-2

26

460 µm

Using NanoFil™ in different configurations

Direct injection by hand: This is the simplest and most economical way to inject. Any of our tips can be inserted directly into the NanoFil™ syringe. Even the SilFlex tubing can be inserted to switch from hand injection to the other methods listed below. This method is limited by the accuracy of plunger movement that is achievable with a human hand.

Installed on WPI's UMP3 microsyringe pump: This will allow the user to achieve nanoliter resolution and reproducibility. For neural system injection, mount the UMP3 on a stereotaxic frame.

SilFlex tubing and holder: The needle is mounted on a small plastic holder that is connected to the NanoFil by a 35 cm length of flexible tubing. The NanoFil syringe is mounted on the UMP3 pump. This configuration allows the user to hold the animal in one hand and insert the needle with the other. When the needle reaches the desired location, activate the pump using the foot switch and the preprogrammed injection volume will be delivered. This configuration gives a nanoliter level of accuracy and reproducibility. It is best suited for applications such as the RPE and IO injection.

SPECIFICATIONS									
Tip Order Number	Gauge	Tip O.D. "A"	Tip I.D. "B nominal	Tip Length "C"	Total Length "D"	Shank O.D. "E"	Bevel Length "F"	Total Dead Volume	Tip Material
NF33BV-2	33	210 µm	115 µm	10 mm	40 mm	460 µm	≈348 µm	0.416 µL	Stainless Steel
NF34BV-2	34	185 µm	85 µm	5 mm	35 mm	460 µm	≈290 µm	0.199 µL	Stainless Steel
NF35BV-2	35	135 µm	55 µm	5 mm	35 mm	460 µm	≈204 µm	0.435 μL	Stainless Steel
NF36BV-2	36	120 µm	35 µm	3 mm	33 mm	460 µm	≈156 µm	0.340 μL	Stainless Steel
NFQ34-5	34	160 µm	100 µm	55 mm	75 mm	460 µm	n/a	0.589 μL	Quartz
NF33BL-2	33	210 µm	115 µm	10 mm	34 mm	460 µm	≈0	0.416 μL	Stainless Steel
NF34BL-2	34	185 µm	85 µm	5 mm	29 mm	460 µm	≈0	0.199 µL	Stainless Steel
NF35BL-2	35	135 µm	55 µm	5 mm	29 mm	460 µm	≈0	0.435 μL	Stainless Steel
NF36BL-2	36	120 µm	35 µm	3 mm	27 mm	460 µm	≈0	0.340 µL	Stainless Steel

ORDERING INFORMATION

460 µm

NANOFIL NanoFil™ Syringe, 10 microliter

35 cm

Includes two 28-gauge MicroFil needles MF28G; one 1cc plastic syringe; and one 26-gauge beveled needle NF26BV.

DETAIL A

2.749 µL

0.380 μL

NANOFIL-100 NanoFil™ Syringe, 100 microliter 26-gauge beveled needle NF26BV included.

REPLACEMENT NEEDLE

3 mm

NF26BV-2 26G Beveled Needle, 460 µm nominal diam. (package of 2)

NanoFil™ NEEDLEs

100 µm

110 µm

NF33BL-2	33 G blunt NanoFil™ Needle (pkg. of 2)
NF34BL-2	34 G blunt NanoFil™ Needle (pkg. of 2)
NF35BL-2	35 G blunt NanoFil™ Needle (pkg. of 2)
NF36BL-2	36 G blunt NanoFil™ Needle (pkg. of 2)
NF33BV-2	33 G beveled NanoFil™ Needle (pkg. of 2)
NF34BV-2	34 G beveled NanoFil™ Needle (pkg. of 2)
NF35BV-2	35 G beveled NanoFil™ Needle (pkg. of 2)
NF36BV-2	36 G beveled NanoFil™ Needle (pkg. of 2)
NF33-36BL	Assortment of 4 blunt NanoFil™ Needles
NF33-36BV	Assortment of 4 beveled NanoFil™ Needles

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

OI HORAL	ACCESSORIES/REI EACEMENT I ARTS
NFINHLD	NanoFil™ Injection Holder
SILFLEX-2	SilFlex Tubing 35 cm long (pkg. of 2) (dead volume = 2.74 μL)
NFGSK-5	Spare Silicone Gasket for NanoFil™ & Holder (pkg. of 5)
NFQ34-5	34 Gauge Flexible Quartz Tubing for filling (pkg. of 5)

NanoFil™ Application Kits

Designed for eye research

Features

- · Includes SilFlex tubing, gasket, holder and tip assortment
- NanoFil™ syringe sold separately
- Less than 3 µL dead volume
- · Achieve accurate, repetitive, oil-free injections

Benefits

- · No oil back filling necessary
- · Comes with four needle sizes included in the kit

Applications

- Intraocular injection
- · Retinal pigment epithelium injection
- Mouse brain injection

These kits are specially designed for eye research for injecting retinal pigment epithelium (RPE) and intraocular (IO) in addition to brain injection in mice. They are used exclusively with a **NanoFil™** syringe and **UMP3** to achieve accurate and repeatable oil free injections down to submicroliter ranges. Each kit includes two pieces of Silflex tubing, a holder assembly, spare gaskets, and an assortment of four needles — blunt for the RPE kit and beveled tips for the IO kit. Each kit comes with one each of 33, 34, 35 and 36 gauge needles so that first time users can find the best size for their application.

The Silflex tubing is a very important component of the kit. This 35 cm long, flexible tubing has a precise outer diameter for airtight fitting with the syringe. It also has a small inner diameter to minimize dead volume. The SilFlex is coupled to the injection tip with a seal system similar to that of the NanoFil. The dead volume of the entire kit (including the tubing) is less than 3 microliters. All of the components in the kit are constructed of inert, solvent resistant materials for easy cleaning after viral injection.



RPE-KIT



ORDERING INFORMATION

RPE-KIT Retinal Pigment Epithelium (RPE) Injection Kit SilFlex tubing, gasket, holder, and blunt tip mix

IO-KIT Intraocular (IO) Injection Kit SilFlex tubing, holder, gasket, 4 beveled tips (33g, 34g, 35 g, 36g)

503207 Stand & Clamps

NFINHLD NanoFil Injection Holder

NFINHLD-G10 1.0 mm Glass Pipette Holder for NanoFil Syringe

References

Park, S. W., Kim, J. H., Park, W. J., & Kim, J. H. (2015). Limbal Approach-Subretinal Injection of Viral Vectors for Gene Therapy in Mice Retinal Pigment Epithelium. *Journal of Visualized Experiments*, (102), e53030–e53030. http://doi.org/10.3791/53030



Nanoliter Injector

For oocyte injection and applications in the 2 to 70nL range

Features

- · Microprocessor-controlled injector
- · Direct piston displacement in capillary glass
- Filling and injection speeds: 23nL/sec and 46 nL/sec
- Optional standard or Micro2T controller
- · Optional foot switch

Benefits

- · No syringes required
- · Graphical representation of volume status
- Automatically calculated injections based on syringe volume and scale length
- · Intuitive touchscreen interface

Applications

· Oocyte injection

WPI's microprocessor-controlled **Nanoliter 2010** uses direct piston displacement. By either pushing the injection button on the control box or pressing on the optional foot switch, a discrete volume will be injected. Choice of capillary filling and injection speeds are 23nL/sec or 46nL/sec (emptying speed is 92nL/sec). Maximum fluid ejection is 5 μ L. Each unit comes with sufficient glass to pull at least 300 micropipettes (see **PUL-1000**, page 48). Glass is 1.14mm O.D. (nominal) and 0.5 mm I.D.

By setting the DIP switch, the injection volume can be changed from 2.3 to 69.0nL in 16 steps. Up to 100 injections may be performed from a single filling of the micropipette. Since the volume of a normal Xenopus oocyte is about 500nL, the instrument has the capability to inject from less than 1% to over 10% of the total volume of the oocyte in one preset step increment.

Included: 1 vial 3.5 in. capillaries (300), replacement "O" rings, Allen wrench, MicroFil **MF34G** backfilling needle and two sample **µTip™** prepulled micropipettes.

New SMARTouch controller

Already own a unit with the standard controller (shown below) or a **Micro4** digital controller? Now you can upgrade to touchscreen control with the new **SMARTouch**.

References

Chouchane, M., et. al. (2017). Lineage Reprogramming of Astroglial Cells from Different Origins into Distinct Neuronal Subtypes. *Stem Cell Reports*, 9(1), 162–176. http://doi.org/10.1016/j.stemcr.2017.05.009

Moreno-Juan, V., et. al. (2017). Prenatal thalamic waves regulate cortical area size prior to sensory processing. *Nature Communications*, 8, 14172. http://doi.org/10.1038/ncomms14172

Tappert, L., et. al. (2017). Sublethal doses of imidacloprid disrupt sexual communication and host finding in a parasitoid wasp. Scientific Reports, 7(1), 42756. http://doi.org/10.1038/srep42756

Lynagh, T., et. al. (2017). Unique Contributions of an Arginine Side Chain to

Ligand Recognition in a Glutamate-gated Chloride Channel. *The Journal of Biological Chemistry*, 292(9), 3940–3946. http://doi.org/10.1074/jbc.M116.772939

Li, S., et. al. (2017). Small RNA-Seq analysis reveals microRNA-regulation of the Imd pathway during Escherichia coli infection in Drosophila. *Developmental & Comparative Immunology*, 70, 80–87. http://doi.org/10.1016/j.dci.2017.01.008

The standard configuration of the Nanoliter 2010 includes the small controller which is simple to setup and operate.





NANOLITER 2010 SPECIFICATIONS

REMOTE CONTROL	Yes
GLASS OD	1.14 mm
GLASS ID	0.5 mm
STEP	12.7 µm/step
INJECTION SPEED Slow Fast	23 nL/sec 46 nL/sec
FILL SPEED Slow Fast	23 nL/sec 46 nL/sec
EMPTY SPEED	92 nL/sec
SINGLE STEP VOLUME RANGE	2.3 - 69.0 nL
SMALLEST VOLUME	2.3 nL
LARGEST SINGLE STEP VOLUME	69 nL
TO CHANGE VOLUME	Set switch
INJECTIONS PER FILLING, MAX.	100 injections
SHIPPING WEIGHT	3 lb. (1.1 kg)

ORDERING INFORMATION

NANOLITER2010 Nanoliter 2010 and Standard Controller

System includes Nanoliter 2010 injector, small controller, 1 vial 3.5 in.

capillaries (300), replacement "O" rings, Allen wrench, MicroFil™ MF34G

backfilling needle and two sample μTip™ pre-pulled micropipettes.

NL2010MC2T Nanoliter Injector & SMARTouch® Controller

NL2010MC2T Nanoliter Injector & SMARTouch® Controller (small controller not included)

System includes Nanoliter 2010 injector, MICRO2T controller, 1 vial 3.5 in. capillaries (300), replacement "O" rings, Allen wrench, MicroFil™ MF34G backfilling needle and two sample µTip™ pre-pulled micropipettes.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

13142	Foot Switch for Nanoliter 2010
504127	Replacement Nanoliter 2010 Injector Head
504949	Replacement 3.5-in. Glass Capillaries, 1.14 mm OD (300)
504950	Replacement 7-in. Glass Capillaries , 1.14 mm OD (300)
TIP10XV119	Micropipettes (10) Prepulled 10 μm
40239	Adapter for SMARTouch
300521	Replacement O-rings (five)
500778	Replacement Nanoliter Injector Mounting Adapter
300033	Adapter for Micro4

Touchscreen Programmable Pump Series

Legato syringe pumps for reliable delivery with the ease of a touch screen

Features

- · High resolution color touch screen
- · Real time clock
- · Touch pad "lock" feature
- · Full metal chassis
- Built in syringe table with up to 75 lb. linear force
- Built in RS-485 interface to link multiple pumps
- USB port & RS232 Interface, TTL interface
- · Continuous mode of operation
- Protection with a spill dam
- · Analog control option
- · CE, UL, CSA, CB Scheme, EU RoHS compliance

Benefits

- Automatic dispensing of small volumes
- · Very precise flow rate control
- Hands free operation with foot pedal
- Better flow performance with accuracy ±0.35%

Applications

- · Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids

The large touch screen color display lets you see all of the pump's operating parameters to ensure proper operation during the experiments. Syringe size and flow rate are easily displayed, as well as the volume delivered and elapsed time. Set up is easy using the icondriven software. An LED on the front panel makes it easy to see if the pump is running. Advanced micro stepping techniques are employed to further reduce the step angle to eliminate flow pulsation. Accuracy is ±0.5%. A wide dynamic flow range from picoliters per minute to millimeters per minute can be programmed into the pump. These versatile pumps can be connected through an RS485 interface. Add the new Adagio software to maximize the use of the pump's functions and features. Adagio allows you to configure the pump through the software, as well as operate one or multiple pumps. LabVIEW drivers are available on the National Instruments website.

SPLG100 SERIES FLOW RATES			
Syringe	Diameter	Minimum	Maximum
0.5 μL	0.103 mm	1.260 pL/min	1.325 μL/min
1 μL	0.146 mm	2.520 pL/min	2.651 μL/min
2 μL	0.206 mm	5.100 pL/min	5.299 µL/min
5 μL	0.343 mm	14.100 pL/min	14.690 μL/min
10 µL	0.485 mm	28.260 pL/min	29.380 μL/min
25 µL	0.729 mm	63.900 pL/min	66.370 μL/min
50 µL	1.03 mm	127.600 pL/min	132.500 μL/min
100 μL	1.457 mm	255.20 pL/min	265.100 μL/min
250 μL	2.304 mm	638.300 nL/ min	662.900 μL/min
500 μL	3.256 mm	1.275 nL/min	1.324 mL/min
1000 μL	4.608 mm	2.553 nL/min	2.652 mL/min
1 mL	4.699 mm	2.655 nL/min	2.757 mL/min
3 mL	8.585 mm	8.863 nL/min	9.204 mL/min
5 mL	11.989 mm	17.290 nL/min	17.950 mL/min
10 mL	14.427 mm	25.030 nL/min	25.990 mL/min
20 mL	19.05 mm	43.640 nL/min	45.320 mL/min
30 mL	21.59 mm	56.050 nL/min	58.210 mL/min
60 mL	26.594 mm	85.050 nL/ min	88.320 mL/min

SPLG100 Infuse-Only Syringe Pump

The SPLG100 was the first single-syringe infusion-only pump with a touchscreen interface. The SPLG100 has a wide flow rate range from 1.26 pL/min to 88.32 mL/ CE min, depending on syringe size. It SPLG100 accommodates a single syringe

from 0.5 µL to 60 mL. Any type of syringe, including glass, plastic or stainless steel are held securely in place.

SPLG100 SPECIFICATIONS

SYRINGE SIZE
POWER
MOTOR DRIVE CONTROL
LINEAR FORCE (MAXIMUM)
NUMBER OF MICROSTEPS PER
REVOLUTION OF LEAD SCREW
STEP RATE (MIN.)
STEP RATE (MAX.)
DRIVE MOTOR
PUSHER TRAVEL RATE (MIN.)
PLOW RATE (MIN.)
FLOW RATE (MIN.)
DIMENSIONS
WEIGHT

CONNECTORS


0.5 μL to 60 mL 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse Microprocessor with 1/16 microstepping 13.6 kg (30 lb.) @ 100% force selection

15,360

27.5 sec/μstep 26 μsec/μstep 0.9 degree Stepper Motor 0.15 μm/min 159 mm/min 1.26 pL/ min (0.5 μL syringe) 88.32 mL/min (60 mL syringe) 22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in) 2.66 kg (5.9 lb.) RS485 - IEEE-1394 6 pos, USB Type B

SPLG101 Dual Infuse-Only Syringe Pump

The **SPLG101** is ideal for applications where dual syringes are required with small volumes up to 10 mL. It accommodates two syringes from 0.5 µL to 10 mL. The **SPLG101** has a wide flow rate range from 1.26pL/min to 25.99 mL/min, depending on syringe size.



SPLG101 SPECIFICATIONS

SYRINGE SIZE
POWER
MOTOR DRIVE CONTROL
LINEAR FORCE (MAX.)
NUMBER OF MICROSTEPS PER
REVOLUTION OF LEAD SCREW
STEP RATE (MIN.)
STEP RATE (MAX.)
DRIVE MOTOR
PUSHER TRAVEL RATE (MIN.)
PUSHER TRAVEL RATE (MAX.)
FLOW RATE (MIN.)
FLOW RATE (MIN.)
STEP RATE (MIN.)
PUSHER TRAVEL RATE (MIN.)
FLOW RATE (MIN.)
STEP RATE (MIN.)
FLOW RATE (MIN.)
STEP RATE (MIN.)
FLOW RATE (MIN.)
STEP RATE (MIN.)
FLOW RATE (MIN.)
DIMENSIONS
WEIGHT

 $0.5~\mu L$ to 10 mL 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse Microprocessor with 1/16 microstepping 13.6 kg (30 lbs.) @ 100% force selection 15,360


27.5 sec/µstep
26 µsec/µstep
0.9 degree Stepper Motor
0.15 µm/min
159 mm/min
1.26 pL/ min (0.5 µL syringe)
25.99 mL/min (10 mL syringe)
22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)
2.66 kg (5.9 lbs.)
RS-232 - 9 Pin D-Sub connector,
RS485 - IEEE-1394 6 pos,

CONNECTORS

SPLG110 Infuse/Withdraw Syringe Pump

The **SPLG110** offers infuse/withdraw flow control and programmability for up to two multi-step programs of 50 steps each. The **SPLG110** has a wide flow

of 50 steps each. The SPLG110 has a wide flow rate range from 1.26pL/ min to 88.28mL/min, depending on syringe size. The SPLG110 accommodates a single



syringe from 0.5 μ L to 60 mL. Any type of syringe can be used in the unit including glass, plastic or stainless steel. The pump is ideal for more complex multi-step dosing and has multi-mode operation including infusion only, withdrawal only, infusion and withdrawal and withdrawal/infusion modes.

SPLG110 SPECIFICATIONS

SYRINGE SIZE
POWER
MOTOR DRIVE CONTROL
LINEAR FORCE (MAX.)
NUMBER OF MICROSTEPS PER
REVOLUTION OF LEAD SCREW

STEP RATE (MIN.)
STEP RATE (MAX.)
DRIVE MOTOR
PUSHER TRAVEL RATE (MIN.)
PUSHER TRAVEL RATE (MAX.)
FLOW RATE (MIN.)
FLOW RATE (MAX.)
DIMENSIONS
WEIGHT
CONNECTORS

0.5 µL to 60 mL 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse Microprocessor with 1/16 microstepping 13.6 kg (30 lb.) @ 100% Force Selection 15.360

27.5 sec/μstep
26 μsec/μstep
0.9 degree Stepper Motor
0.15 μm/min
159 mm/min
1.26 pL/ min (0.5 μl syringe)
88.28 mL/min (60 ml syringe)
22.6 x 19.05 x 15 cm (9 x 7.5 x 5 in)
2.66 kg (5.9 lb.)
RS-232 - 9 Pin D-Sub connector,
RS485 - IEEE-1394 6 pos,

SPLG210 Infuse/Withdraw Syringe Pump

USB - Type B

The SPLG210 Infuse/Withdraw syringe pump is easy to use with the high resolution touch screen. The basic model works with one syringe or two (from 0.5 µL to 140 mL) and can be reconfigured in the field to be used with multiple syringes. A protective cover over the display prevents leakage into the display. To optimize your bench space, the SPLG210 can be placed on its side to

reduce the footprint to only 3.5 x 9.75 in. The display also rotates with the change to allow you to operate the pump vertically. The programmable model offers maximum flexibility for configuring and running different programs. Up to 40 programs of 20 steps each can be configured and stored for quick recall with the touch of a button.

SPLG210 SPECIFICATIONS

SYRINGE SIZE **POWER** MOTOR DRIVE CONTROL LINEAR FORCE (MAX.) NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW STEP RATE (MIN.) STEP RATE (MAX.) DRIVE MOTOR PUSHER TRAVEL RATE (MIN.) PUSHER TRAVEL RATE (MAX. FLOW RATE (MIN.) FLOW RATE (MAX.) **DIMENSIONS** WEIGHT CONNECTORS

0.5 µL to 140 mL 100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse Microprocessor with 1/16 microstepping 34 kg (75 lb.) @ 100% force selection 6400

27.5 sec/µstep
26 µsec/µstep
1.8 degree stepper motor
0.36 µm/min
190.80 mm/min
5 pL/min (0.5 µL syringe)
215.803 mL/min (140 mL syringe)
8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in))
4.9 kg (10.75 lb.)
RS-232 - 9 Pin D-Sub connector,
RS-485 - IEEE-1394 6 pos,
USB - Type B,

I/O & TTL - 15 Pin D-Sub connector

SPLG270 Push/Pull Continuous

The SPLG270 is a Push-Pull syringe pump. It accommodates two syringes from 0.5 µL to 140 mL for infusion and two syringes for withdrawal. This model supports infusion and withdrawal and withdrawal



simultaneously at user-defined flow rates and with selectable target volumes to control the total volume pumped. It also supports infuse only, withdraw only, infuse/withdraw, withdraw/infuse and continuous mode. The touch screen interface lets you quickly create configurations and recall them for easy use. The 4.3" TFT color display with touch pad interface presents all the pump operating parameters on one easy-to-view run screen. Protective cover over the display prevents leakage into the display.

SPLG270 SPECIFICATIONS

SYRINGE SIZE

POWER

100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse

MOTOR DRIVE CONTROL

LINEAR FORCE (MAX.

NUMBER OF MICROSTEPS PER
REVOLUTION OF LEAD SCREW

STEP RATE (MIN.)

27.5 sec/µstep

STEP RATE (MAX.)

28 µsec/µstep

100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse

Microprocessor with 1/16 microstepping

34 kg (75 lb.) @ 100% force selection

4400

27.5 sec/µstep

26 µsec/µstep

DRIVE MOTOR

1.8° Stepper Motor
PUSHER TRAVEL RATE (MIN.)

PUSHER TRAVEL RATE (MAX.)

190.80 mm/min

FLOW PATE (MAN.)

 FLOW RATE (MIN.)
 5 pL/min (0.5 μL syringe)

 FLOW RATE (MAX.)
 215.803 mL/min (140 mL syringe)

 DIMENSIONS
 8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in)

WEIGHT 4.9 kg (10.75 lb.)
CONNECTORS RS-232 - 9 Pin D-

S RS-232 - 9 Pin D-Sub connector, RS-485 - IEEE-1394 6 pos, USB - Type B, I/O & TTL - 15 Pin D-Sub connector

Adagio Software

Low cost, simple installation
 Flow evolution graph
 Import and export programs
 Quick and easy manual pump control
 Monitor one or more pumps
 Program data logging

The manual pump control tool allows easy direct control of the pump. Pump commands can be entered directly into the log. Multiple programs can be opened at the same time. The program's progression is tracked, and can be stored in a file for later access.

Intuitive Run Screen — Combining multiple parameters simultaneously with internationally recognizable icons allow the Legato™ Series to provide a new level of intuitive syringe pump operation.

	ORDERING INFORMATION
SPLG100	Legato 100 Syringe Pump, Infuse-Only
SPLG101	Legato 101 Syringe Pump, Dual Infuse-Only
SPLG110	Legato 110 Syringe Pump, Infuse/Withdraw
SPLG200	SPL Syringe Pump, Infuse Only
SPLG210	SPL Syringe Pump, Infuse/Withdraw
SPLG212	SPL Syringe Pump, Infuse/Withdraw Programmable
SPLG270	SPL Syringe Pump, Push-Pull
SPLG272	SPL Syringe Pump, Push-Pull Programmable
504576	Small Syringe Multi Rack (for six 30-60 mL syringes or ten
	0.5 μL-20 mL syringes)
504577	Large Syringe Multi Rack (for up to four 60-140 mL plastic
	syringes)
504578	Software Adagio/USB Key

Touchscreen Programmable Pump Series

Legato syringe pumps for reliable delivery with the ease of a touch screen

Features

- · High resolution color touch screen
- Real time clock
- · Touch pad "lock" feature
- · Full metal chassis
- Built in syringe table with up to 75 lb. linear force
- Built in RS-485 interface to link multiple pumps
- USB port & RS232 Interface, TTL interface
- Continuous mode of operation
- · Protection with a spill dam
- · Analog control option
- CE, UL, CSA, CB Scheme, EU RoHS compliance

Benefits

- Automatic dispensing of small volumes
- · Very precise flow rate control
- · Hands free operation with foot pedal
- Better flow performance with accuracy ±0.35%



Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids

When mounted vertically, the display screen of the SPLG series pumps automatically reorients for ease of use.

	SPLG200	SERIES FLOW R	ATES
Syringe	Diameter	Minimum	Maximum
0.5 μL	0.103 mm	3.12 pL/min	1.589 μL/min
1 μL	0.146 mm	6.18 pL/min	3.180 μL/min
2 μL	0.206 mm	12.301 pL/min	6.358 μL/min
5 μL	0.343 mm	33.96 pL/min	17.630 μL/min
10 μL	0.485 mm	67.72 pL/min	35.249 µL/min
25 μL	0.729 mm	153.42 pL/min	79.640 µL/min
50 µL	1.03 mm	306.24 pL/min	158.984 µL/min
100 μL	1.457 mm	612.72 pL/min	318.126 µL/min
250 μL	2.304 mm	1.533 nL/ min	795.51 μL/min
500 μL	3.256 mm	3.06 nL/min	1.588 mL/min
1000 μL	4.608 mm	6.129 nL/min	3.181 mL/min
1 mL	4.699 mm	6.373 nL/min	3.308 mL/min
3 mL	8.585 mm	21.272 nL/min	11.044 mL/min
5 mL	11.989 mm	41.485 nL/min	21.539 mL/min
10 mL	14.427 mm	60.073 nL/min	31.19 mL/min
20 mL	19.05 mm	104.74 nL/min	54.383 mL/min
30 mL	21.59 mm	134.533 nL/min	69.852 mL/min
50 mL	26.594 mm	204.122 nL/ min	105.985 mL/min
100 mL	35.7 mm	367.839 nL/min	190.992 mL/min
140 mL	38.4 mm	415.623 nL/min	215.803 mL/min



SPLG200 Infuse-Only Syringe Pump

The **SPLG200** Infuse only syringe pump is easy to use with the high resolution touch screen. The basic model works with one syringe or two (from 0.5 μL to 140 mL) and can be reconfigured in the field to be used with multiple syringes. To optimize your bench space, the **SPLG200** can be placed on its side to reduce the footprint to only 3.5 x 9.75 inches. The display automatically reorients itself with the change to allow the user to operate the pump vertically. The pump features user definable flow rates with selectable target volumes and time values to control the infusion rate and the total volume. Up to 40 programs of 20 steps each can be configured and stored in the unit for quick recall with the touch of a button.

SPLG200 SPECIFICTIONS

SYRINGE SIZE	0.5 μL to 140 mL
POWER	100-240 VAC: 50/60 Hz, 50W. 0.5 A fuse
MOTOR DRIVE CONTROL	Microprocessor with 1/16 microstepping
LINEAR FORCE (MAX.)	34 kg (75 lbs.) @ 100% force selection
NUMBER OF MICROSTEPS PER REVOLUTION OF LEAD SCREW	6400
STEP RATE (MIN.)	27.5 sec/µstep
STEP RATE (MAX.)	26 µsec/µstep
DRIVE MOTOR	1.8° stepper motor
PUSHER TRAVEL RATE (MIN.)	0.36 μm/min
PUSHER TRAVEL RATE (MAX.)	190.80 mm/min
FLOW RATE (MIN.)	5 pL/min (0.5 μL syringe)
FLOW RATE (MAX.)	215.803 mL/min (140 mL syringe)
DIMENSIONS	8.89 x 25.4 x 27.94 cm (3.5 x 10 x 11 in.)
WEIGHT	4.9 kg (10.75 lb.)
CONNECTORS	RS-232 - 9 Pin D-Sub connector,
	S485 - IEEE-1394 6 pos,
	USB - Type B,
	I/O & TTL - 15 Pin D-Sub connector

	ORDERING INFORMATION
SPLG200	SPL Syringe Pump, Infuse Only
504576	Small Syringe Multi Rack (for six 30-60 mL syringes
	or ten 0.5 μL-20 mL syringes)
504577	Large Syringe Multi Rack
	(for up to four 60-140 mL plastic syringes)
504578	Software Adagio/USB Key

Affordable Syringe Pumps

Syringe pumps for high metering precision at low, pulse-free rates

Features

- · Sturdy and reliable
- Extremely simple to set up and use
- Surprisingly affordable.

Benefits

- Automatic dispensing of small volumes
- Very precise flow rate control
- Hands-free operation with foot pedal
- Better flow performance with accuracy ±0.35%

Applications

- Drug administration
- Chemical applications with slow incorporation of fixed volumes of


On the SP series syringe pumps, the liquid crystal displays (LCDs) prompt you through setup:

- Select syringe from table stored in the pump's memory and displayed on the LCD.
- Enter the volume to be dispensed.
- 3. Enter the flow rate and press Start.

It's fast and simple. Your settings are permanently stored in memory there's no need to re-enter them each day. SP pumps feature preset rate and volume control. Just set the volume you want dispensed. Volume is tracked continuously on the LCD display. Then, when the preset volume has been dispensed, the pump shuts off automatically. The easy-to-read digital display provides real-time readings using both parameters and values for clearer, mistake-free readings. The SP200 Series pumps offer TTL and RS-232C interfaces and automatic shutoff under stall conditions.

SP210c Continuous Cycle Syringe Pump

The **SP210c** holds up to four syringes and can cycle continuously back and forth in a push-pull action. As two syringes are infusing, two other syringes are withdrawing at the same rate. At the end of the set volume the direction is automatically



· Holds four syringes, 10 mL to 60 mL each



SP220i Multi-Syringe Infusion Pump

Ideal for applications requiring multiple syringes, the

been modified to hold up to syringes. · All features of SP200i Accommodates 10

to 140 mL.

SP220i is an adaptation of the SP200i and has



SP100i Single-Syringe Infusion Pump

This inexpensive single-syringe infusion pump combines precision and simplicity with outstanding ease of use and durability.

- Holds any size syringe, 10 μL to 50 mL
- Automatic volume control and shutoff
- Simple menu-driven setup: dispense volume, dispense flow rate, syringe diameter
- · Last settings stored in permanent memory
- Continuous dispense volume display

CE SP100i

SP120p Two-Syringe Push-Pull Pump

A second syringe mount has been added to the basic SP100i, with both syringes activated by a single pusher block for simultaneous infusion and withdrawal.

- · All the features of SP100i
- · Holds two syringes, from $10 \mu L$ to 10 mL.



	ORDERING INFORMATION
SP100i	Syringe Pump, Infusion (Single) 95-135V
SP100iZ	Syringe Pump, Infusion (Single) 220-240V
SP101i	Syringe Pump, Microdialysis (Double, Slow Speed) 95-135V
SP101iZ	Syringe Pump, Microdialysis (Double, Slow Speed) 220-240V
SP120p	Syringe Pump, Infusion-Withdrawal (Double) 95-135V
SP120pZ	Syringe Pump, Infusion-Withdrawal (Double) 220-240V
SP200i	Syringe Pump, Infusion (Double) 95-135V
SP200iZ	Syringe Pump, Infusion (Double) 220-240V
SP210c	Syringe Pump, Infusion-Withdrawal (Continuous) 95-135V
SP210cZ	Syringe Pump, Infusion-Withdrawal (Continuous) 220-240V
SP210iw	Syringe Pump, Infusion & Withdrawal (Double) 95-135V
SP210iwZ	Syringe Pump, Infusion & Withdrawal (Double) 220-240V
SP220i	Syringe Pump, Infusion (Multiple) 95-135V
SP220iZ	Syringe Pump, Infusion (Multiple) 220-240V
SP230iw	Syringe Pump, Infusion & Withdrawal (Multiple) 95-135V
SP230iwZ	Syringe Pump, Infusion & Withdrawal (Multiple) 220-240V
SP250i	Syringe Pump, Infusion (Multiple, Mixed Volumes) 95-135V
SP250iZ	Syringe Pump, Infusion (Multiple, Mixed Volumes) 220-240V
SP260p	Syringe Pump, I/W (Double) Single Cycle Action, 95-135V
SP260pZ	Syringe Pump, I/W (Double) Single Cycle Action, 220-240V
	All 240-volt pumps are CE-approved.
####-A	Audible Alarm (add "A" to pump part number)
####-P	Programmable Ramp Option (SP200 Series)

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

15623	Serial Cable, SP Pump-to-IBM 9-pin "D" Connector
13685	SP Pump-to-Pump "Daisy-Chain" Linking Cable, 7 ft.
13962	Foot Switch for SP200 Series Pumps

Dual Rate Syringe Pump

Two separate syringe pumps in one instrument

Features

- TTL, USB & RS-232 Communication
- · Compatible with a wide variety of syringes
- · Large touch screen display allows quick and easy setup

Benefits

- Two separate syringe pumps in one instrument
- · High or low-pressure operation
- · Continuous fluid delivery

Applications

- · Continuous injection for long-term toxicology testing
- The injection of dyes, perfumes and flavoring in industrial applications
- Applications with liquids or viscous materials in micromanufacturing
- · Continuous injections of reactants into reactor vessels
- · Simultaneous samplings from two sites

Various applications can benefit from the. All three conditions come with the Gemini 88 Plus. No additional upgrades are required.

The intuitive Gemini 88 Plus graphical user interface controlled with a large 7" LCD color touch screen display allows quick and easy setup. The display run screen presents the user with all key dispensing parameters in real time. Syringe tables containing all major syringe manufacturers allow simple selection of any compatible syringe size. Audible Alarms, Adjustable Force and Screen Lock are all features that are available with a touch of the screen.

The SPLG-G88PLUS Dual Rate Syringe Pump opens whole new pumping possibilities. Three operating conditions named Independent, Reciprocating and Twin are available as standard features to satisfy various applications in life science.

Connectivity & Communication

The Gemini 88 Plus comes standard with USB and RS-232 for PC communication and RS-485 for pump-to-pump communication. An entire suite of ASCII commands is available to control the pump remotely with a PC. The pump contains a foot switch input and digital input/output for each independent pumping channel.

Independent Condition

Independent Condition allows the Gemini 88 Plus to operate as two separate syringe pumps named P1 and P2. Each syringe will operate independently with different syringe types, sizes, force, or target (volume or time, mode dependent).

Reciprocating Condition

In Reciprocating Condition, both syringe channels move in opposite directions at the same rate using the same syringe size and type. When combined with a valve box, the Reciprocating Condition can provide the continuous fluidic delivery or a peristaltic pump with the accurate, pulseless, low flow rates provided by a syringe pump.

Twin Condition

Twin Condition allows both syringes to operate in the same mode using the exact same syringe type, syringe size, force, target (volume or time) and flow rate settings. The pump also allows the user to combine both flows for higher speed and volume infusion applications

ORDERING INFORMATION

Brazil: (013) 406-29703 • info@brazil.wpiinc.com • www.wpiinc.com

SPLG-G88PLUS Dual Rate Syringe Pump



SPLG-G88PLUS SPECIFICTIONS

TYPE Microprocessor dual rate infuse/withdraw/con-

tinuous syringe pump

 $\begin{array}{lll} \text{ACCURACY} & \pm 0.25\% \\ \text{REPRODUCIBILITY} & \pm 0.05\% \\ \end{array}$

CONNECTORS

AVERAGE LINEAR FORCE

SYRINGE TYPE Glass, plastic and stainless steel

SYRINGE SIZE MINIMUM 0.5 µL (0.103 mm minimum inner diameter)
SYRINGE SIZE MAXIMUM 60 mL (32.573 mm maximum inner diameter)

FLOW RATE MINIMUM 1.02 pL/min (0.5 µL syringe, 0.103 mm inner

diameter)

FLOW RATE MAXIMUM 106 mL/min (60 mL syringe, 32.573 mm diam-

eter)

FLOW RATE DISPLAY 7 in. color display with touchscreen

USB - Type B

RS-232 - 9-pin D-Sub connector

RS-485 - IEEE-1394, 6 pos for pump-pump com-

munication

TTL Input/Output - Two 15-pin D-sub connec-

tors, one for reach pump mechanism

Foot switch - Two phonojack inputs, one for each

pump mechanism

70 lb. (31.75 kg) at 100% force setting up to a flow rate of 90 mL/min using up to a 60 mL syringe with a 32.573 mm inner diameter. 50 lb. (22.6 kg) at 100% Force Setting for flow rates

(22.6 kg) at 100% Force Setting for flow rates 90 mL/min to 106 mL/min using the same size

syringe.

POWER SUPPLY Input 100-240 VAC, 50-60 Hz, Output 30 V 1.66A

50 Wall

WEIGHT 21 lb. (9.09 kg)

DIMENSIONS (L x D x H) 11 x 15 x 8 in. (28 x 39 x 21 cm)

CLASSIFICATION Class I
POLLUTION Degree 1
INSTALLATION Category II

REGULATORY
CE, ETL (UL & CSA), CB Scheme, EU RoHS, WEEE

Programmable Syringe Pump

The best value in syringe pumps on the market!

Features

- · Cost-effective
- Program sequences without a computer
- Rich command set
- RS-232 and TTL control

Benefits

- Automatic dispensing of small volumes
- · Very precise flow rate control
- Hands-free operation with foot pedal

Applications

- · Drug administration
- Chemical applications with slow incorporation of fixed volumes of fluids



The Aladdin pump series arguably offers the best VALUE for any syringe pump on the market.

These pumps are available in single, dual, 4, 6, and 8 syringe capacities. Although it is one of the most cost-effective pumps available, the Aladdin pump series boasts a versatility that is unmatched in its price range with features typically found only on pumps costing more than twice as

The Aladdin series provides the capability to both inject and withdraw. These pumps can be used for manually triggered injections, but also have the capability to be programmed in multi-step pumping sequences without the requirement for a computer. Even so, this series also has both RS-232 and TTL inputs for external control and status feedback to a computer, if desired. The pumps can be daisy chained for multi pump applications involving push/pull protocols or multiple independent channels. The performance characteristics are admirable, considering

For the budget minded lab looking for a versatile and reliable pump, the Aladdin series is an excellent option.

The Aladdin pump series will accept syringes from Becton Dickinson, Monoject, Terumo and Air-Tite.

	ORDERING INFORMATION
AL-300	Non-Programmable Syringe Pump, Infuse Only
AL-1000	Programmable Syringe Pump
AL-1000HP	Programmable Syringe Pump, High Pressure (100 lb.)
AL-1000HP2	Programmable Syringe Pump, High Pressure (200 lb.)
AL-2000	Two AL1000 Syringe Pumps
Includes CBL-Dual-3 Interconnecting Cable for push/pull or continuous	
	pumping. Valves not included.

AL-4000	Programmable Multiple (up to 4) Syringe Pump
AL-6000	Programmable Multiple (up to 6) Syringe Pump
AL-6000H	Programmable Multiple (up to 6) Syringe Pump
AL-8000	Programmable Multiple (up to 8) Syringe Pump

Specify line voltage. When ordering 220V models, specify UK, Euro or Australian line cord.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

GN-PC7	PC to Pump Cable, 7 ft.
GN-PC25	PC to Pump Cable, 25 ft.
GN-NET7	Pump-to-Pump Network Cable, 7 ft.
GN-NET25	Pump-to-Pump Network Cable, 25 ft.
GN-TTL	Pump-to-Pump Reciprocating Cable
ADPT2	Foot Switch

		ALA	DDIN SPECIFICA	ATIONS		
	AL-300	AL-1000	AL-1000 HP	AL-4000	AL-6000	AL-8000
SYRINGE SIZES	Plastic syringes up to 60 mL	selected glass micro syringes	Plastic syringes up to 60 mL and selected glass micro syringes from 0.5 to 500 µL.	Plastic syringes up to 60 mL	1-60 mL, or 0.5-5 μL micro syringes	Up to 100 mL
NUMBER OF SYRINGES	1	1	1	2 (may be different sizes)	6	8
MOTOR TYPE	Step Motor, 1/8 to 1/2 step modes	Step Motor, 1/8 to 1/2 step modes		Step Motor, 1/2 to 1/2 step modes	Step Motor, 1/8 to 1/2 step modes	Step Motor
STEPS PER REVOLUTIONS	400	400	200		400	200
STEPPING (max./min.)	0.21 μm to 0.850 μm	0.21 μm to 0.850 μm		0.425 μm to 1.7 μm, depending on motor speed	0.21 μm to 0.850 μm	0.132μm to 0.265 μm, depending on motor speed
MOTOR TO DRIVE SCREW RATIO	15/28	15/28	15/28	15/28	15/28	5/1
SPEED (max./min.)	3.774 cm/min. to 0.004 cm/hr	5.1 cm/min / 0.0042 cm/hr	18.370 cm/min / 0.008 cm/hr	18.370 cm/min / 0.008 cm/hr	5.1 cm/min / 0.0042 cm/hr	3.4917cm/min. / 0.0026cm/hr.
PUMPING RATES	1257 mL/hr with 60 mL syringe, to 0.73 μL/hr with 1 mL syringe	to 0.73 ul /hr with 1 ml syringe,	6120 mL/hr with 60 mL syringe, to 1.459 µL/hr with 1 mL syringe	6120 mL/hr with 60 mL syringe, to 1.459 μ L/hr with 1 mL syringe	1699 mL/hr with 60 mL syringe, to 0.73 μL/hr with 1 mL syringe	0.454 µL/hr. with a B-D 1cc syringe to 1163mL/hr. with a B-D 60 cc syringe
MAXIMUM FORCE	35 lb. at minimum speed, 18 lb. at maximum speed	35 lb. at min. speed, 18 lb. at maximum speed	100 lb. at minimum speed, 18 lb. at maximum speed	100 lb. at minimum speed, 15 lb. at maximum speed	35 lb. at min. speed, 18 lb. at max. speed	160 lbs. at min. speed, 30 lbs. at max. speed
PROGRAM PHASES	n/a	41	41	41	41	41
RS-232 PUMP NETWORK	n/a	100 pumps maximum	100 pumps maximum	100 pumps maximum	100 pumps maximum	100 pumps maximum
POWER SUPPLY	Wall adapter 12V DC @ 850 mA	Wall adapter 12V DC @ 850 mA	Wall adapter 12V DC @ 1000 mA	Wall adapter 12V DC @ 1000 mA	Wall adapter 12V DC @ 850 mA	Unregulated linear external wall adapter, country and power source specific
DIMENSIONS	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	22.9 x 14.6 x 11.4 cm (8.75 x 5.75 x 4.5 in.)	26x38.1x12.7cm (10.25 x15 x5 in.)
WEIGHT	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	1.6 kg (3.6 lb.)	10.125 lbs. (4.595kg)

Microfluidic Pumps

Microfluidic Solenoid Pump



Features

- Controlled by the iKima™ application for use with the iPod Touch (included) and iPhone
- Connects easily to Vena8 Endothelial+ biochip and other manufacturer's flow chambers
- Delivers pulses of fresh media to cells seeded inside the microchannels
- Flow rate: 15 35 mL/hr ± 4%
- Dead volume: < 300 µL

Benefits

- Fits inside standard CO_2 incubators, maintaining temperature and humidity
- · Recirculating long-term perfusion pump
- WiFi communication and control via iPod Touch, even with Kima pump incubated

Applications

- · Cell culture under shear stress/flow
- · Biofilm studies
- Cell culture in biochips with adherent cells (HUVECs)
- · Stem cells
- HepG2 cells (human liver cancer cell line)

Kima pump is a microfluidic pump designed to aid cell culture (e.g., epithelial and endothelial cells) under physiological conditions (shear flow) in various biochips and flow chambers, including Vena8 Endothelial+ biochips where it is possible to culture eight cell monolayers simultaneously over 24-48 hours. The solenoid pump produces a pulsatile flow to mimic the cardiovascular system, delivering a preset volume of liquid as programmed.

	ORDERING INFORMATION
CX-KIMAKIT-IPOD	Kima Kit for Cell Culture/Biofilm Studies Includes Kima Pump, iPod Touch with Controller, 100 mL bottle & cap, tubing kit, Vena Endothelial+ biochips (10-pack).
CX-KIMAKIT-PC	Kima Kit for Cell Culture/Biofilm Studies Includes Kima Pump, PC Controller, 100 mL bottle & cap, tubing kit, Vena Endothelial+ biochips (10-pack).
CX-KIMA-I	Kima Microfluidic Pump with iPod Touch & Controller
CX-KIMA-P	Kima Microfluidic Pump with PC Control Software and Controller
CX-KIMA	Kima Microfluidic Pump

Precision Syringe Pump



Features

- Virtually pulse-free low volume delivery
- Fast response time
- Flow sensor feedback to pump
- Touch control

Benefits

- · Pulse-free flow
- · Precise multichannel mixing
- · Excellent long-term flow stability
- · Rapid flow change (ms range)
- WiFi communication and control via iPad mini or USB communication with LabVIEW based PC control software
- Use standard tubing for connection to any microfluidic biochip

Applications

- · Microfluidics and nanofluidics
- · Droplet generation and manipulation
- · Laminar and multilaminar flow studies
- · Cell-based shear flow studies
- Cell and particle manipulation studies

ExiGo is a precision syringe pump based on the 5-phase stepper motor drive that has more microsteps per revolution of the lead screw vs. standard syringe pumps on the market. ExiGo has 250,000 microsteps/ revolution and even at low rotational speed/low flow rates, it has a very low pulsation and high accuracy. When coupled with the flow sensor and active PID feedback; this results in very fast response times for changing flow rates. A standard syringe pump typically has a smaller number of microsteps and so usually the only way a standard "microfluidic" pump can achieve pulse-free flow control is to use small syringes; e.g. 0.5 μL; 1 μL; 5 μL; etc., to achieve non-pulsatile stable flow rates in the nanoliter/minute range. By comparison, the ExiGo pump with the flow sensor can use a standard 250 µL glass syringe to produce stable nonpulsatile flow rates of 10nL/min - 1 mL/min; or a 5 mL plastic syringe to produce stable non-pulsatile flow rates of 100nL/min - 20 mL/min. ExiGo can be used in conjunction with expandable (flexible) element and fluidic resistance in order to dampen any pulsation occurring during the stepper motor operation. As it employs active feedback, the response time of the pumps still remains fast.

	ORDERING INFORMATION
CX-EXIGO-IM-7	Exigo with iPad, Manifold, 7µL/min FS
CX-EXIGO-PM-7	Exigo with LabView, Manifold, 7µL/min FS
CX-EXIGO-I-7	Exigo with iPad, 7µL/min FS
CX-EXIGO-P-7	Exigo with LabView, 7µL/min FS
CX-EXIGO-IM	Exigo with iPad, Manifold
CX-EXIGO-PM	Exigo with LabView, Manifold
CX-EXIGO-I	Exigo with iPad, SmartFlo app
CX-EXIGO-P	Exigo with LabView
CX-EXIGO-M	Exigo with Manifold, Tubing Kit
CX-EXIGO	Exigo with Tubing Kit

Microfluidic Pumps

Microfluidic Pressure Pump



Features

- Precise flow control with active feedback via plug-and-play flow sensor (required add-on)
- Flow rate: 1 μL/min 1 mL/min; unidirectional (push)
- iPad mini or PC (LabVIEW, Matlab, Python, etc.) control which can control/program up to four pump modules independently
- External compressor (required)

Benefits

- 2 modes of pumping: manual flow rate set or preprogrammed flow rate operation
- Preprogrammed mode includes: constant, ramp, step, sine functions
- Side port connections to dock up to four pumps (combination of UniGo and ExiGo) into one setup
- Use standard tubing for connection to any microfluidic biochip.

Applications

Microfluidics where accurate and stable flow rate delivery is required.
 The pressure pump component is based on controlled air injection.

The **UniGo™** Microfluidic Pump is a precision, microfluidic, single-channel pressure pump for a variety of microfluidic applications, where accurate and stable flow rate delivery is required. The pressure pump component is based on controlled air injection. The UniGo pump requires a plug-and-play flow sensor for active feedback and increased flow control. **SmartFlo** application executed on the iPad Mini or LabView-based interface communicates with up to four UniGo™ Microfluidic Pumps racked together, allowing simultaneous control and independent programming of each pump's flow profile. Uniquely, the UniGo™ may be docked together with the **ExiGo™** microfluidic syringe pump, combining the best features of both instruments in one microfluidic set-up. *Note: it is necessary to purchase a flow sensor with the UniGo™ pump.*

ORDERING INFORMATION

CX-UNIGO-ECI-80 UniGo with iPad, External Compressor, 80 μL/min Flow Sensor CX-UNIGO-ECP-80 UniGo with LabView, External Compressor, 80 μL/min Flow Sensor CX-UNIGO-EC-80 UniGoExternal Compressor, 80 μL/min Flow Sensor CX-UNIGO-I-80 UniGo with iPad, No Compressor, 80 μL/min Flow Sensor CX-UNIGO-P-80 UniGo with LabView, No Compressor, 80 μL/min Flow Sensor CX-UNIGO-80 UniGo, No Compressor, 80 μL/min Flow Sensor CX-UNIGO-80 UniGo, No Compressor, 80 μL/min Flow Sensor

Five flow cells are avilable for your Unigo system. Your choice of one flow sensor is included with your Unigo system. Contact WPI at wpi@wpiinc. com for more information.

Microfluidic syringe pump for shear flow studies



Features

- · Unlimited volume
- · Low per-channel cost
- Flow rate: 100nL/min 10 mL/min ± 1%
- Dead volume: ~600 μL

Benefits

- · Higher throughput enabling eight assays in parallel
- · Flow damper to decrease syringe pump pulses
- PC controlled via VenaFluxAssay software

Applications

- · Microfluidic applications
- Single cell analysis
- Microfluidic syringe pump for cell analysis under shear flow in biochips
- Suitable for cell samples and whole blood samples

Mirus is a precision syringe pump, which uses the combination of an expandable (flexible) element and fluidic resistance in order to dampen the pulsation of the syringe pump stepping. The Mirus has a very stable flow profile. Additionally, Mirus is equipped with 3-way valve allowing automatic recharging and washout of syringe. Mirus is also provided with an 8-way flow splitter, allowing multiple executions (eight parallel experiments) simultaneously.

0	RDERING INFORMATION
CX-MIRUS-PRO	Microfluidic Syringe Pump Includes tubing kit, VenaFlux Assay Software PC control and MultiFlow8
CX-MIRUS	Microfluidic Syringe Pump Includes tubing kit, VenaFlux Assay Software PC control
CX-MIRUS-MULTI8	MultiFlow8 Attachment for Mirus Evo Nanopump

Microfluidic Pumps and Accessories

4-channel Microfluidic Pump



Features

- Precise flow control with active feedback via plug-and-play flow sensor (required add-on)
- Flow rate: 1 µL/min 1 mL/min; unidirectional (push)
- iPad mini or PC (LabVIEW, Matlab, Python, etc.) control which can control/program up to four pump modules independently
- · External compressor (required)
- Working pressure- Max 10 bars/145 psi

Benefits

- Two modes of pumping: manual flow rate set or preprogrammed flow rate operation
- Preprogrammed mode includes: constant, ramp, step, sine functions
- Side port connections to dock up to four pumps (combination of UniGo and ExiGo) into one setup
- · Use standard tubing for connection to any microfluidic biochip

Applications

Microfluidics where accurate and stable flow rate delivery is required.
 The pressure pump component is based on controlled air injection.

4U is a 4-channel compact, precise, economical microfluidic pressure pump for a variety of microfluidic applications, where accurate and stable flow rate delivery is required. **SmartFlo** application executed on LabView-based interface communicates with the four separate channels of the **4U** pump, allowing simultaneous control and independent programming of each pump's flow profile. The **4U** pump provides you with the advantage of using four independent flow rates at the same time.

ORD	ORDERING INFORMATION			
CX-4UECIPAD-4X80FS	4-channel Microfluidic Pump, External			
	Compressor, iPad, 80 µL/min flow sensors (4)			
CX-4UECPC-4X80FS	4-channel Microfluidic Pump, External			
	Compressor, PC, 80 µL/min flow sensors (4)			
CX-4UNCIPAD-4X80FS	4-channel Microfluidic Pump, iPad, 80 μL/min			
	Flow Sensors (4)			
CX-4UNCPC-4X80FS	4-channel Microfluidic Pump, PC, 80 μL/min			
	Flow Sensors (4)			
CX-EC	External Compressor			



Vena8 Fluoro+







VenaT4

ExiGo™ Pump Manifold

Features

- Allows you to direct the fluid to three (3) independent ports
- Plug and play connection to the ExiGo™ pump
- Can be programmed to automatically switch between fluidic channels using SmartFlo PC software

Benefits

· Refill syringes automatically

Applications

- · Microfluidics and Nanofluidics
- Droplet generation and manipulation
- · Laminar and multilaminar flow studies
- · Cell-based shear flow studies
- · Cell and particle manipulation studies

The Manifold is a specialized microfluidic channel selector which allows the **ExiGo** pump to direct fluid to one of three microfluidic channels at a time. Accurate flow switching and low dead volume provide exceptional performance.

CX-MF

-	חם	-51		чест	RMATI	\sim \sim 1
	1740		VI - II	VI E. 115		

CX-MF ExiGo Pump Manifold

Microfluidic Biochips for a Variety of Applications

We offer a range of biochips for many applications. All biochips mimic human capillaries by working in tandem with our different instrumentation platforms.

Our new DropChips are either hydrophilic or hydrophobic based. DropChips comes with 5 different channel sizes from 30 μ m to 80 μ m. Please refer to the chart on page 69 for information on specific biochips.

ORDERING INFORMATION

Optional A	Accessories/Replacement Parts
CX-018	Vena8 with Glass Coverslip
CX-023	Biochip-CONN, Biochip SNGL IN
CX-024	Biochip-CONN, Biochip SNGL OUT
CY-027	Tygon Tuhing for Biochin Conn

CX-027 Tygon Tubing for Biochip Conn
 CX-032 Kima-Connect Tubing Set
 CX-045 Kima Pump Cell Culture Bottle
 CX-067 Flow Sensor Exigo 1 mL/min
 CX-073 Dropchip, Hydrophobic, 40 μm x 40 μm
 CX-084 MF8-Connect Biochip1 Inlet CBL

CX-084 MF8-Connect Biochip'i Inlet CBL
CX-087 MULTIFLOW8 Nozzles, Standard



Vena8 Glass Vena Delta Y1 Cover slip



Vena Delta Y2



Microfluidic Biochips

DROPCHIE	(HYDROPHIL	IC/HYDROPH	OBIC) BIOCHI	PS ORDERING	INFORMATIO	N
Chip Size (µm)	30 x 30	40 × 40	50 x 50	60 x 60	70 x 70	80 x 80
Junction Width (cm)	0.003	0.004	0.005	0.006	0.007	0.008
Height (cm)	0.003	0.004	0.005	0.006	0.007	0.008
Length (cm)	2.8	2.8	2.8	2.8	2.8	2.8
Substrate Thickness	0.5	0.5	0.5	0.5	0.5	0.5
Hydrophilic (Pack of 2) Hydrophobic (Pack of 2)	CX-026 CX-034	CX-028 CX-073	CX-029 CX-035	CX-030 CX-036	CX-031 CX-037	CX-033 CX-038
Hydrophilic (Pack of 5) Hydrophobic (Pack of 5)	CX-039 CX-046	CX-040 CX-047	CX-041 CX-048	CX-042 CX-049	CX-043 CX-050	CX-044 CX-051
Hydrophilic (Pack of 10) Hydrophobic (Pack of 10)	CX-052 CX-058	CX-053 CX-059	CX-054 CX-060	CX-055 CX-061	CX-056 CX-062	CX-057 CX-063

STANDARD BIOCHIP	S ORDE	RING INFO	RMATIO	N			
	Vena8 Fluoro+	Vena8 Endothelial+	VenaT4		Glass erslip	Vena	Delta
				Low Flow Rates	High Flow Rates	Y1	Y2
Channel width, b (cm)	0.04	0.08	0.08	0.16	0.08	0.008	0.008
Channel height, h (cm)	0.01	0.012	0.01	0.016	0.008	0.012	0.012
Channel length, L (cm)	2.8	2.8	2.8	2.8	2.8	2.8	2.8
Channel volume (cm³)	0.00112	0.00269	0.00224	0.00717	0.00179	0.00269	0.00269
Channel volume (μL)	1.12	2.69	2.24	7.17	1.79	2.69	2.69
Thickness of bottom substrate (mm)	0.17	0.5		0.17	0.17	0.5	0.5
# of channels / chip	8	8	4	8	8	4	4
# of assays / pack	40/80	40/80	40	80	80	40	40
Pack of 5	CX-002	CX-004					
Pack of 10	CX-001	CX-003	CX-005	CX-009	CX-010	CX-007	CX-008
SPE	CIFICATI	ONS					
Brightfield, phase contrast, immunostaining	1	✓	1	1	✓	1	1
Confocal microscopy	1			1	√		
Cell types: whole blood (human, animal); PRP; platelets; PBMC; monocytes; T-cells (primary and cell lines); eosinophils; neutrophils etc.	✓	1	1	1	1	√	✓
Cell types: adherent cells <i>e.g.</i> endothelial – HUVECs; HMVECs etc.; HepG2; stem cells; muscle cells <i>etc</i> .	✓	1	1	1	✓	1	✓
Protein coatings: collagen, fibronectin, fibrinogen, vWF, VCAM, ICAM, selectins, MadCAM <i>etc</i> .		1	✓	1	✓	✓	✓
АРІ	PLICATIO	ONS					
Platelet adhesion, aggregation and thrombi formation; leuko- cyte rolling, adhesion and migration; thrombosis; immunology (inflammation); infectious diseases (e.g. malaria); sickle cell disease; respiratory (asthma and COPD)	1			1	✓		
Cell adhesion and culture under perfusion / shear flow; leukocyte cell-cell rolling, adhesion and migration; oncology (melanoma, breast cancer etc.); cardiovascular (atherosclerosis, drug eluting stents); immunology (inflammation); respiratory (asthma and COPD)		1					
Biofilm assays, microbe seeding and culturing; biochips with glass coverslips (attached / not attached; treated / non-treated); biochips for the attachment of coupons for biofilm studies				1	✓		
Chemotaxis, transmigration and invasion assays; 2D and 3D cell culture; mimicking tumour microenvironment with gels (ECM gel, hydrogel, matrigel, collagen gel)			✓				

All biochips are:

- Disposable plastic; some with glass cover slips.
- Require no assembly; unlike many standard perfusion chambers / flow chambers.
- Require no Luer lock connections which increase dead volume. Cellix's biochips have a unique plug and play connection with tubing connections which are autoclaveable and reusable.

Ergonomic Pipetter Design

Making repetitive procedures more efficient and comfortable

Features

- Lightweight and comfortable ergonomic design
- Easy calibration using provided tool
- Easy for cleaning and parts replacement
- · CE and ISO13485 Certified

REG100

Benefits

- Less stress on your hand when you are performing repetitive operations
- Save money by ordering sets of any 5, 6 or 7 pipetters
- Made from biologically inactive and chemical inert polymers

Applications

 Routine laboratory use for accurate and affordable pipetting of liquids and solutions



Regal Pipettes use the air displacement method, where displaced air volume by the motion of the piston is equal to the liquid volume drawn into the pipette. These new pipettes are highly accurate. Within ten complete revolutions of the dial, you can set the minimum and maximum volumes. For ease of use, the dial is fixed to the plunger. Since the light plunger action reduces fatigue, results are more precise.

ORDERING INFORMATION (SINGLES)						
Model	Volume Range μL	Increment µL	Nominal Volume µL	Tolerance %	Repeatability %	
REG2	0.2 ~ 2	0.01	0.2 0.5 2	±12.0 ±5.0 ±2.0	≤6.00 ≤2.50 ≤0.70	
REG10	1 ~ 10	0.1	1 5 10	±3.0 ±1.5 ±1.0	≤1.50 ≤0.60 ≤0.40	
REG20	2 ~ 20	0.1	2 10 20	±3.0 ±1.0 ±1.0	≤1.50 ≤0.50 ≤0.30	
REG50	5 ~ 50	0.5	5 20 50	±2.0 ±1.2 ±1.0	≤1.50 ≤0.40 ≤0.20	
REG100	10 ~ 100	1	10 50 100	±2.0 ±0.8 ±0.8	≤0.50 ≤0.30 ≤0.15	
REG200	20 ~ 200	1	20 100 200	±2.0 ±0.8 ±0.8	≤0.50 ≤0.30 ≤0.15	
REG1000	100 ~ 1000	5	100 500 1000	±1.5 ±0.8 ±0.8	≤0.30 ≤0.30 ≤0.15	
REG5K	1000 ~ 5000	50	1000 2000 5000	±1.0 ±0.7 ±0.7	≤0.50 ≤0.25 ≤0.15	
REG10K	1000 ~ 10000	100	1 mL 5 mL 10 mL	±3.0 ±0.7 ±0.7	≤0.30 ≤0.20 ≤0.15	

ORDERING INFORMATION (SETS)
Regal Pipetters (set of any 5) & Stand
Regal Pipetters (set of any 6) & Stand
Regal Pipetters (set of any 7) & Stand
Stand for Regal Pipetters (holds 8)

Universal Pipette Tips Mini Vacuum Pump

Ultra-clear, certified RNase/DNase-free Small, reliable, durable and accurate



Same as leading brands—

at about half the price!

ORDERING INFORMATION **UNIVERSAL FILTER TIPS (STERILE)**

Tip Volume	For Pipetter	Rack	Part No.
0.1 - 10 μL	REG2 REG10 REG20	960 (10 racks of 96)	500199
10 - 200 μL	REG20 REG50 REG100 REG200	960 (10 racks of 96)	500200

ORDERING INFORMATION **UNIVERSAL TIPS**

Tip Volume	For Pipetter	Bulk	Part No.	Rack	Part No.
0.1 - 10 μL	REG2 REG10 REG20	Bag of 1000	500191	960 (10 racks of 96)	500192
5 - 200 µL	REG20 REG50 REG100 REG200	Bag of 1000	500193	960 (10 racks of 96)	500194
100-1000 μL	REG1K	Bag of 1000	500195	1000 (10 racks of 100)	500196
500 - 5000 μL	REG5K	Bag of 250	500197 *	500 (10 racks of 50)	500198 *



Features

- · Durable aluminum exterior
- · Minimal vibration
- · Low noise
- · Extremely long life time

Benefits

- · Quiet operation
- Compact (18 x 7 x 7cm) unit takes little of your bench space
- · Oil free, maintenance free

Applications

- Commonly sold with Muscle Physiology setups
- Excellent accessory for use with WPI's PV830 Pneumatic PicoPump with vacuum
- Ideal for any application requiring a small, reliable pump that provides vacuum pressure up to 250 mbar

This miniature vacuum pump is durable and accurate. The industrialstrength aluminum exterior, neoprene diaphragm and neoprene/silicone valves ensure this pump will stand up to daily use.

MINI VAC SPECIFICATIONS					
POWER SOURCE	230 (50 Hz)	120 (60 Hz)			
FREE FLOW	4.0 L/min.	3.0 L/min.			
AT -100 MBAR	2.0 L/min.	1.5 L/min.			
MOTOR TYPE	Vibrati	ng			
POWER 4.0 W					
MAXIMUM PRESSURE	-				
MAXIMUM VACUUM	-250 m	bar			
PUMP HEAD CONSTRUCTION	Alumin	um			
DIAPHRAGM	CR-neop	rene			
VALVES	CR-neoprene/FPM	(Viton)/Silicone			
DIMENSIONS	DIMENSIONS 185 x 72 x 72 mm				
WEIGHT 850 g					

	INFORMATION ORDERING				
801566	Mini Vacuum Pump, 110 V, US power plug				
801963	Mini Vacuum Pump, 220 V, EURO mains plug				

Luer Valve Assortment Kit

Build your own liquid flow experiment

Features

- Over 300 assorted parts
- Luer fittings for quick and easy connect and disconnect

Benefits

• Sold individually or in kits

Applications

Liquid flow experimental setups

A useful kit (right) for building your own liquid flow experiment. It provides the means to start, stop, add, divide and control a flow of liquid or gas. Included in the kit are *over 300 assorted parts* such as one-way and threeway stopcocks, manifolds, Y-connectors, injection sites,

male and female luer caps, check valves, syringe-activated check valves, slide clamps, roller clamps, and pinch clamps. All (except clamps) have a luer fitting for quick and easy connecting and disconnecting. Includes assorted luer fittings for use with flexible tubing.

The 14011 component parts are also sold individually. See page 220.



ORDERING INFORMATION

14011 Luer Valve Assortment Kit

Luer-to-Tubing Coupler Assortment Kit

Quick connects in nylon and polypropylene

Features

- · Over 250 assorted parts in each kit
- Valves are polycarbonate, and the valve handles are polyethylene. Do not autoclave those parts.

Benefits

- Polypropylene parts (504954) can be autoclaved repeatedly at 121°C/15PSI, 15 min. cycle
- Polypropylene fittings are chemically inert and resistant to most organic and inorganic solvents
- Nylon fittings are strong and can be bonded with adhesive.

Applications

• Liquid flow experimental setups



504954



ORDERING INFORMATION

504954 Luer-to-Tubing Coupler Assortment Kit (Polypropylene)

504955 Luer-to-Tubing Coupler Assortment Kit (Nylon)

Assemble quick-disconnect luer fittings for use with flexible tubing with internal diameters of 1/16", 3/32" and 1/8". A variety of quick-disconnect connectors can be quickly made for connecting small diameter flexible tubing; 3-way connections can be made with the use of the 3-way luer tee; luer plugs, tees, connectors, bulk-head mounts, color coding rings, locking nuts, male and female luers—are all included to enhance the versatility of this kit. The kit has over 250 assorted parts and is offered in two different types of materials. Nylon parts are not autoclavable.

Electrophysiology



Over 50 years of experience with researchers

From our early beginnings in 1967 working with Yale researchers, electrophysiology has been at the heart of our business. WPI amplifiers, stimulators and isolators are designed with quality components so you get a reliable, low-noise signal every time. Our time-tested designs give you affordable solutions for electrophysiology equipment, electrodes, data acquisition and accessories for applications like:

- ·Intracellular/Extracellular Recording
- Voltage Clamp for Ussing
- Stimulation and Isolation

- Optogenetics
- Digital Filtering

Bioamplifiers for Electrophysiology

A family of very low noise battery-operated amplifiers

Features

- Battery powered to eliminate line noise
- · High pass and low pass filtering
- Single ended or differential operation
- DC/AC amplification
- Variable output positioning
- Constructed of high quality components to ensure minimal intrinsic (shot) noise
- Portable
- Rack mountable

Benefits

- · Very low internal noise
- Ultra quiet DC power supply no AC required
- · Instrinsic low susceptibility to ground loops
- Small footprint
- · Cost effective
- Electrostatic Discharge Protection!

Applications

- Amplifying biopotentials from metal electrodes
- Brain slice field stimulation
- EAG (Electroantennogram)
- ERG (Electroretinogram)

WPI's **DAM** series amplifiers are well known as a standard of the industry for amplification of extracellular potentials. These battery powered bio-amplifiers are designed with a compact chassis profile that enables you to locate the unit closer to the preparation and thereby minimize long lead lengths which contribute to noise. Each amplifier is equipped with selectable high and low filters, and a position control to offset galvanic potentials which may develop during recording.

DAM series amplifiers can be used as stand-alone units on any tabletop or use optional clamp-mounting hardware to locate them conveniently within the work area. Alternatively, a pair of amplifiers can be mounted into a standard equipment rack with a rack mount kit (3484). A variety of hook up accessories are available to configure your application.

DAM80 Overview

DAM80, an AC amplifier only, features a very low noise active headstage probe which can be mounted in micromanipulators for up-close cortical recording, for extracellular recording from high impedance glass or metal microelectrodes. The unit also provides a gated current for tissue marking. Microelectrode holder **MEH7W-XX** (sold separately) is recommended for glass microelectrodes. The DAM80 is perfect for gated or manual current generation for histological marking, iontophoresis or cell stimulation. It includes a very low noise remote active headstage that is useful for very high impedance amplification utilizing glass or metal electrodes.



CBI 102

Included with the DAM80 is a **Startup Kit** containing the following accessories needed for basic metal electrode electrophysiology research:

(2) Cable BNC-to-3.5 mm nlug 6 ft (2m)

CDLIGE	(2) cable, bive to 3.5 min plag, one (2m)
5469	(2) Adapter, mini-banana to 0.031 socket
13388	(2) Adapter, mini-banana to 2mm socket
3294	Cable, ground clip to wire, 3 ft
2033	Mini-banana plug, black
2034	Mini-banana plug, red
2035	(2) Mini-banana plug solderable turrent
EP1	Ag/AgCl pellet (70 mm wire) 1 mm diam x 2.5 mm long
M3301EH	(2) Electrode Holder, 14cm
5470	0.031-inch jack on 12-inch wire (package of 4)

FEATURE COMPARISON						
DAM50 DAM80						
INPUT MODE	AC/DC	AC				
INPUT CONFIGURATION	differential/single ended	differential				
GAIN RANGE	100-10K (AC), 10-1K	100-10K (AC)				
	(DC)					
HIGH/LOW FILTERS	yes	yes				
OFFSET POSITION	yes	yes				
CONTROL						
CURRENT GENERATOR	No	Yes				
REMOTE ACTIVE HEADSTAGE	No	Yes				
OUTPUT CONNECTION	BNC	3.5 mm mini phone				
STANDARD INPUT	unterminated wire	mini banana				
CONNECTION*						
POWER SUPPLY	(2) 9V alkaline batteries	(2) 9V alkaline batteries				

^{*}See optional accessories for additional alternatives

DAM50 References

Liu, Y., Wang, Y., Zhu, G., Sun, J., Bi, X., & Baudry, M. (2016). A calpain-2 selective inhibitor enhances learning & Depril memory by prolonging ERK activation. *Neuropharmacology*, 105, 471–477. http://doi.org/10.1016/j. neuropharm.2016.02.022

Ztaou, S., Maurice, N., Camon, J., Guiraudie-Capraz, G., Kerkerian-Le Goff, L., Beurrier, C., ... Amalric, M. (2016). Involvement of Striatal Cholinergic Interneurons and M1 and M4 Muscarinic Receptors in Motor Symptoms of Parkinson's Disease. Journal of Neuroscience, 36(35).

Kentish, S. S. J., Frisby, C. L., Kritas, S., Li, H., Hatzinikolas, G., O'Donnell, T. A., ... Ahern, G. (2015). TRPV1 Channels and Gastric Vagal Afferent Signalling in Lean and High Fat Diet Induced Obese Mice. PloS One, 10(8), e0135892. http://doi.org/10.1371/journal.pone.0135892

DAM80 References

Donnelly, W. T., Bartlett, D., & Leiter, J. C. (2016). Serotonin in the solitary tract nucleus shortens the laryngeal chemoreflex in anaesthetized neonatal rats. Experimental Physiology, 101(7), 946-961. http://doi.org/10.1113/ FP085716

Feng, B., Joyce, S. C., & Gebhart, G. F. (2016). Optogenetic activation of mechanically insensitive afferents in mouse colorectum reveals chemosensitivity. American Journal of Physiology - Gastrointestinal and Liver Physiology, 310(10).

Kaldenbach, F., Bleckmann, H., & Kohl, T. (2016). Responses of infraredsensitive tectal units of the pit viper Crotalus atrox to moving objects. Journal of Comparative Physiology A, 202(6), 389-398. http://doi.org/10.1007/s00359-016-1076-1

Mustafina, A. N., Koroleva, K. S., Giniatullin, R. A., & Sitdikova, G. F. (2016). Acid Sensitive Ion Channels as Target of Hydrogen Sulfide in Rat Trigeminal Neurons. BioNanoScience, 1-3. http://doi.org/10.1007/s12668-016-

Orton, L. D., Papasavvas, C. A., & Rees, A. (2016). Commissural Gain Control Enhances the Midbrain Representation of Sound Location. The Journal of Neuroscience, 36(16), 4470-81. http://doi.org/10.1523/ JNEUROSCI.3012-15.2016

DAM SERIES SPECIFICATIONS

INPUT IMPEDANCE $10^{12} \Omega$ INPUT I FAKAGE CURRENT 50 pA (typical) MAX. DC DIFFERENTIAL SIGNAL ± 2.5 V (DAM 50) GAIN AC: 100x, 1000x, 10000x DC: 10x, 100x, 1000x (DAM50) 100 dB @ 50/60 Hz COMMON MODE REJECTION RATIO INPUT CAPACITANCE 20 pF AC MODE NOISE $0.4~\mu V$ RMS (2 μV p-p) 0.1-100~HzAC MODE NOISE $2.6 \mu V RMS (10 \mu V p-p) 1 Hz-10 kH$ 7.5 μV RMS (30 μV p-p) 3-10 kHz DC MODE NOISE (DAM50) BANDWIDTH FILTER SETTINGS

AC Mode Low frequency, 0.1, 1, 10, 300 Hz AC Mode (DAM80) High frequency, 0.1, 1, 3, 10 kHz DC Mode (DAM50) High frequency, 0.1, 1, 3, 10 kHz **OUTPUT CONNECTORS** BNC on DAM50; 3.5 mm MiniPhone connector on DAM80

OUTPUT VOLTAGE SWING +8 V **OUTPUT IMPEDANCE** 470 Ω BATTERY TEST Audible tone **CALIBRATOR SIGNAL** 10 Hz square wave Approximately 250 mV **POSITION**

DAM80: DC Generator 0 to ±50 μA, variable

EXTERNAL COMMAND Input Voltage ±10 V commands AC OR DC CURRENT WAVEFORM $\pm 50~\mu A$ max. amplitude @ 200 K Ω

BATTERIES DIMENSIONS

CURRENT SOURCE

DAM50 8 x 4 x 1.75 in. (20.3 x 10.2 x 4.4 cm) DAM80 7 x 4 x 1.75 in. (17.8 x 10.2 x 4.4 cm)

2 x 9V alkaline (included)

SHIPPING WEIGHT 3.5 lb (1.6 kg)



OPD	EDIN	CIN	IEODN	/ATION

SYS-DAM50	Bio-amplifier with shielded electrode cable
SYS-DAM80	Bio-amplifier with active probe (DAM80P)

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

DAM80P	Replacement Probe
3072	6 Replacement Modular Cables (DAM50)
3517	2 Optional Shielded Modular Cables (DAM50)
CBL102	3.5 mm Phone plug-to-BNC Cable
2851	BNC-to-BNC Cable
2033	Black Insulated Mini-Banana Plug
2034	Red Insulated Mini-Banana Plug
2035	Uninsulated Mini-Banana Plug
3484	Rack Mount Kit (for 1 or 2 DAM preamps)
3485	Ringstand Mounting Kit
5447	Electrode Adapter (DAM50)
5469	Metal Microelectrode Adapter for DAM80
	(mini-banana plug to 0.031 in. (0.79 mm) socket)
5489	Adapter for Metal Microelectrode (DAM50)
13388	Adapter, mini-banana plug to 2mm socket
5371	Cable, Low Noise (2 mm pin to 2 mm pin)
3578	Adapter Cable for Ag/AgCl pellets (2 mm pin)
300102	Metal Electrode Extension, 4-in. (2 mm to 0.031 in.)
300647	Shielded Electrode Cable Assembly
MEH7W-XX	Microelectrode Holder- 1.0, 1.5 or 2.0 mm OD

See Cables and Connectors, page 96. See Metal Microelectrodes, page 94.



Optional probe 5489 (non-active) for use with DAM50 also includes microelectrode adapter 5469.

Which Amplifier is Right for You?

			AMPL	IFIER COMPA	RISON CH	ART		
Amplifier	AC/DC	Differential	Active Headstage	Stimulation	Isolated	Multi-channel	Battery Powered	Connectors
Intracellul	ar Bioa	mplifiers						
Electro 705	DC	-	*				♦	2 mm pin
Duo773	DC	*	♦	*		2		2 mm pin
Extracellu	lar Bioa	mplifiers						
ISODAM8A	DC	•	opt		•	4 - 8		Mini Banana or 8-pin DIN
ISO80	AC	*	*	*	*		♦	Mini Banana
DAM50	AC/DC	*					♦	RJ-11
DAM80	AC	*	♦	*			*	Mini Banana
Transducer Amplifiers								
BRIDGE8	DC	+				4 - 8		8-pin DIN WPI transducers
ТВМ4М	DC	•				4		8-pin DIN WPI transducers

WPI's Low-Noise Amplifiers Outperform Cheap Imitations

An amplifier is an electronic device that magnifies an input signal. However, the way an amplifier is designed to handle noise and bandwidth limitations greatly affects the quality and sustainability of the final output signal.

Defining Terms

To knowledgeably discuss amplifiers, let's define a few common terms

Gain – The gain is the multiplier defining how much the amplitude of an input signal is increased. A signal with an ×1 gain is not amplified. A ×10 gain produces an output signal ten times greater than the input signal.

Noise – Any unwanted signal fluctuations are called noise. While noise can also result from external sources, for the purpose of this discussion, we are primarily concerned with the noise resulting from the inner workings of the electronic device, our amplifier. This intrinsic noise is called shot (or schott) noise.

Signal to Noise Ratio (SNR) – The ratio of the output signal to the noise of the amplifier is called the signal to noise ratio. The smaller the shot noise signal in an amplifier in comparison with the output signal, the easier the desired signal is to discriminate. When engineering an amplifier, the SNR may be improved by boosting the first stage gain to yield a larger output signal or by using quality components to minimize the shot noise level of the amplifier.

Output Range – The output range determines the maximum output signal that can be generated with the amplifier. It is determined by the maximum voltage of the power supply. If the amplitude of the



output signal is too large for the output range, part of the signal is cut off (clipped).

Rail – The upper or lower limit of the amplifier range is called a rail. Signals that exceed the rail cannot be faithfully reproduced.

DC Offset – DC offsets can appear in biological preparations. This offset is the amount the output signal is displaced away from a zero reference point. It is usually a result of the potential difference at the electrode's tip.

How Amplifiers Work Power Supply Rails Limit the Range

In a perfect world an input signal can be infinitely multiplied by the gain factor to determine the output signal. For example:

Input Signal	Gain	Output Signal
2 mV	×1	2 mV
2 mV	×2	4 mV
2 mV	×10	20 mV
2 mV	×100	200 mV
2 mV	×10,000	20 V


In the real world, the power supply rails limit the possible output range of the amplifier. For example, a bio-amplifier could have a range of ± 5.0 V. In order for the output signal to be faithfully reproduced, the input signal times the gain factor must fall within the voltage window set by the power rails. Otherwise, the output signal goes off scale, and the input signal is not faithfully reproduced. This is called "hitting the rail."

In our example, a 1.0 μ V input signal at an $\times 10^6$ gain would generate a 1.0V output signal. Since the power supply is rated up to +5.0V, this output signal is clearly visible. If the input signal in this example is greater than 5.0 μ V, the output signal would be greater than +5.0V. Since 5.0V is the top of the range that the power supply is capable of producing, the output signal hits the upper rail and gets cut off. This amplifier gives a +5.0V DC output signal for all input signals greater than or equal to 5.0 μ V. So, a smaller gain factor should be used to bring the output signal back into the dynamic output range of the amplifier.

Noise Limits Amplifier Usability

All electronic devices produce their own internal electronic noise, an unavoidable signal that can mask the output signal. For example, if

Signal vs. Noise



The higher the signal to noise ratio, the more discernible the desired signal.

the input signal is 2mV and the noise is 1 mV, the signal to noise ratio is two to one (2:1), and the output signal would be undetectable. It is nearly impossible to discern which part of this output is generated by noise and which part is the desired signal. (See figure.)

Ideally, the signal to noise ratio should be at least 50 to 1 to produce a quality output signal. A good signal to noise ratio can be achieved in one of two ways:

- Boost the output signal by increasing the gain.
- Reduce the noise.

While increasing the gain is the simplest solution, too much gain can impose a limitation on the dynamic range of the amplifier. Reducing noise is a more complicated solution, but it offers a greater range and more stability.

Two-Stage Amplifiers

Bio-amplifiers usually involve multiple stages of amplification.

Stage One – The unadulterated signal coming into the amplifier is unaffected by the intrinsic noise of the amplifier. Then, it runs through the critical first stage of amplification where the signal is boosted by the primary gain factor to produce an output signal with the desired signal

to noise ratio. The intrinsic noise is not amplified in the first stage. Higher gain factors used in the first stage of amplification can seriously limit the dynamic range available at output stage. Large stage one gains also limit the gain factor available in the second stage of amplification.

Stage 2 – The stage one output signal enters the second stage of amplification where both the signal and the noise from the first stage are amplified together by the second stage gain factor so that the signal is large enough to be seen on a chart recorder or data acquisition system. The second stage amplification is the gain the user controls. It does not change the signal to noise ratio.

Instead of using high gains in the first stage of amplification, a well constructed bio-amplifier that uses high quality components, like WPI's DAM series amplifiers, minimizes the noise in the first stage of amplification so that the dynamic range is retained throughout the amplification process. Poorly designed amplifiers simply increase the gain of the first stage amplification until the desired signal to noise ratio is reached.

Boost the power rails?

Theoretically, increasing the voltage rails powering the amplifier will increase the available dynamic output range. It would seem natural to increase the power supply rails coming into the amplifier in order to provide the capability for greater first stage gains. However, most data acquisition systems are limited to a maximum input signal ranging between ±10.0V. Therefore, it is not practical to increase the power rails of bio-amplifier beyond ±10.0V. Since the industry standard limits us to ±10.0V power supply rails, the only way to improve the signal to noise ratio is to minimize the shot noise in the first stage of amplification. This is why high quality amplifier components are imperative.

Why a Signal Flatlines


Regardless of the amplifier used, biological potentials are often accompanied by a DC offset, because the electrodes polarize over time. The DC offset naturally increases over time. Since the poorly constructed amplifier that utilizes greater first stage gain has restricted its dynamic range, it has limited ability to handle this offset. As the offset continues to increase, the output signal may eventually be forced by the offset into the rail causing the flat line (clipping the signal). (See Figure.)


The amplifier that minimizes the noise in the first stage amplification offers a larger dynamic output range and handles a much greater offset value.


WPI's Amplifiers

The purchase of a low-noise amplifier pays dividends in the end. WPl's amplifiers were engineered for the bio-medical researcher. While 20-30 μV of noise is common in bio-amplifiers, WPl's DAM series amplifiers generate 0.4 μV RMS (root mean squared) at 0.1-100 Hz. (That's equal to 2 μV peak-to-peak.)

Result of Amplifier Using Gain to Control Signal to Noise Ratio







As the offset naturally increases over time, a poorly constructed amplifier will not be able to faithfully reproduce the signal. This offset can also result from gain drift as temperature rises.

Isolated Differential Amplifier

Excellent recording performance for extracellular nerve AP

Features

- Battery powered, rechargeable
- · High pass and low pass filtering
- Active remote headstage
- AC only amplification
- Electrode impedance test function
- Electrode current generation with polarity select
- Variable output positioning

Benefits

- Ultra quiet DC power supply
- Intrinsic low susceptibility to ground loops
- High signal to noise ratio due to remote head stage
- Small footprint
- · Stimulation/histological marking current

Applications

- Amplifying bio-potentials using metal microelectrodes
- · Brain slice field potentials
- EAG (Electroantennogram)
- ERG (Electroretinogram)
- Monitor extracellular nerve action potentials
- Use for cell marking, stimulation or electrode cleaning
- In vivo cortical recording

The **ISO-80** provides low noise AC coupled amplification and offers excellent recording performance for monitoring extracellular nerve action potentials *in vitro* and in living animals. The **ISO-80** is provided with a remote headstage (1 m cable) which incorporates an electrode impedance test function and a constant current stimulator. The constant current stimulator can be used for cell marking, stimulation or electrode cleaning. Typical applications include measuring EMG, EEG, extracellular and action potentials *in vitro* or *in vivo*. The ISO-80 system is DC isolated from the subject ground and employs state of the art electro-magnetic shielding for improved noise rejection. The amplifier employs both high pass and low pass filtering with gain from 100 to 10,000. The lowest low-pass setting is 5 Hz and the upper passband is 10 kHz.

Included with the ISO-80 is a Startup Kit containing the following accessories needed for basic metal electrode electrophysiology research:

CBL102 Cable, BNC-to-3.5 mm plug, 6 ft (2m) (two)
5469 Adapter, mini-banana to 0.031 skt. (two)
13388 Adapter, mini-banana to 2mm skt. (two)
3294 Cable, ground clip to wire, 3 ft
Mini-banana plug, black

2034 Mini-banana plug, red
2035 Mini-banana plug solderable turrent (two)

EP1 Ag/AgCl pellet (70 mm wire) 1 mm diam x 2.5 mm long

M3301EH Electrode Holder, 14cm (two)

5470 0.031-inch jack on 12-inch wire (package of 4)



ISO-80 SPECIFICATIONS

INPUT RESISTANCE $>10^{11} \Omega$, Common mode and differential

FILTER SETTINGS

Low frequency 5, 10, 100, 300 Hz High frequency 100 Hz, 1, 3, 10 kHz

MAX. OUTPUT VOLTAGE SWING ± 8 volts ELECTRODE IMPEDANCE RANGE $100 \text{ k}\Omega$

ELECTRODE IMPEDANCE RANGE 100 k Ω - 10 M Ω @ 300 Hz STIMULATION CURRENT 0 to ±20 μ A (constant current) MAXIMUM STIMULATION VOLTAGE ±15 V

MAXIMUM ELECTRODE VOLTAGE ±40 V
DISPLAY 3½-digit LCD
BATTERY TEST Low battery display

POWER (2) 9V Ni-MH batteries & charger, supplied

SHIPPING WEIGHT 4 lb (1.8 kg)

ORDERING INFORMATION

ISO-80 Isolated Bioamplifier w/ active probe (ISO80P)

Specify line voltage

OPTIONAL ACCESSORIES/REPLACEMENT PARTS


ISO-80P Replacement ISO-80 Probe
CBL102 3.5 mm phone plug-to-BNC cable

References


Lynn, P. A., et. al. (2003). Rectal intraganglionic laminar endings are transduction sites of extrinsic mechanoreceptors in the guinea pig rectum. Gastroenterology, 125(3), 786–94. Retrieved from http://www.ncbi.nlm.nih.gov/pubmed/12949724

Vastani, N., et. al. (2013). Sensitivities of rat primary sensory afferent nerves to magnesium. European Journal of Anaesthesiology, 30(1), 21–28. http://doi.org/10.1097/EJA.0b013e32835949ab

Vastani, N., et. al. (2013). Sensitivities of rat primary sensory afferent nerves to magnesium. European Journal of Anaesthesiology, 30(1), 21–28. http://doi.org/10.1097/EJA.0b013e32835949ab







Single-Ended Application

Optional Differential Application

Low Noise Modular Amplifier System

Isolated, low noise bio-amplifier

Features (Iso-DAM8A)

- · High pass and low pass filtering
- Optional active remote head stage
- AC/DC amplification
- · Variable gain adjustment
- Input is optically isolated
- 50/60 Hz notch filter
- Pre-optical isolation DC offset
- Post-optical isolation zeroing
- Independent module power switch

Benefits (Iso-DAM8A)

- Chassis accepts combination of bioamplifiers and transducer amplifiers
- Flexible channel count (1–8) allows expandability
- Notch filter targets AC line noise sources
- · Variable gain output amplitude
- Wide ±10V output range

ISO-DAMBA STATES STA

74020 **((**

Applications (Iso-DAM8A)

- · Amplifying biopotentials unsing metal microelectrodes
- Brain slice field potentials
- EAG (Electroantennogram)
- ERG (Electroretinogram)

Iso-DAM8A Modules Isolated Low Noise Bio-Amplifier

The ISO-DAM8A is a compact modular standard rack-mountable DC amplifier system. Each channel is electrically isolated from the others and from ground. No current can flow from the input terminals and electrodes, thus, the instrument is intrinsically safe and cannot cause any electrical stimulus or shock to the preparation, in addition ground loop noise is minimized. Systems can be purchased with one, two, three or up to eight preamplifier modules or mixed with Bridge8 transducer amplifier modules (see next page). The user can then select an appropriate low pass filter setting, gain and offset on the channel amplifier panel. A notch filter has been added to reduce line frequency interference. An optional headstage preamplifier (10x gain) allows low noise extracellular (DC) recording with Iso-DAM8A and adds greater signal bandwidth than a shielded cable of the same length. The Iso-DAM8A amplifier and headstage configuration is optimally suited for use with our metal microelectrodes and can be easily configured for many applications. Each amplifier channel has a coaxial (BNC) connector located on the rear

ISO-DAM8A SPECIFICATIONS

EACH CHANNEL
INPUT IMPEDANCE
INPUT LEAKAGE CURRENT
INPUT DC OFFSET
GAIN

COMMON MODE REJECTION EQUIVALENT NOISE SIGNAL INPUT

BANDWIDTH FILTER SETTINGS High Filter (Low Pass) (kHz) Low Filter (High Pass) (Hz) Notch Filter (Hz)

OUTPUT VOLTAGE SWING OUTPUT RESISTANCE ENCLOSURE DIMENSIONS SHIPPING WEIGHT > 10¹² Ohms 10 pA (typical) ± 100 mV

×10, ×100, ×1000, ×10,000 > 100 dB @ 50/60 Hz

< 0.36 μV rms (1.8 μV p-p) 0.1-10 Hz, Gain>10

< 1 μV rms (5 μV P-P) 0.1-10 kHz

10 to 21 lb (4.5 to 9.5 kg)

0.1, 0.5, 1, 3, 10 0.1, 1, 10, 300 50, 60 ± 7.5 Volts 220 Ohms 7 × 17 × 9.2 in. (18 × 43 × 23 cm)

Low noise tranducer amplifier

Features (BRIDGE8)

- Wide range of fixed gains with independent variable gain adjustment
- · Low pass filter
- Single ended or differential transducer compatibility
- Dual range output offset correction
- · Independent module power switch
- Provides ± voltage excitation to transducers

Benefits (BRIDGE8)

- Chassis accepts combination of bioamplifiers and transducer amplifiers
- Flexible channel count (1–8) allows expandability
- Output LEDs confirm transducer output balance

Applications (BRIDGE8)

- WPI force transducers
- Wheatstone bridge transducers
- · Muscle force measurement



BRIDGE8

DGE8 C

Bridge8 Modules Low Noise Transducer Amplifier

Bridge8 is a modular, rack-mountable amplifier system. It is specifically designed for use as a signal conditioning amplifier with strain gages and other powered transducers. **Bridge8** includes differential amplifiers featuring high input impedance, high common mode rejection and low current leakage input terminals for low noise operation. It features a half bridge switch and channel offset A wide variety of WPI transducers are available for force, temperature, pressure and light measurements. The **Bridge8** amplifier is a clear choice for convenience and quality.

BRIDGE8 SPECIFICATIONS

INPUT IMPEDANCE
AMPLIFICATION
INPUT LEAKAGE CURRENT
VOLTAGE OFFSET ADJUSTMENT
AMPLIFIER OUTPUT VOLTAGE
EXCITATION VOLTAGE
EQUIVALENT NOISE SIGNAL INPUT

LOW PASS FILTER BAND (KHz)

> 10¹² Ohms 1, 10, 50, 100, 500, 1000 & Adjustable 0.1 pA at 25 °C ±50 mV (low); ±100 mV (high) ±4.4 V (10 mA, max.) 10 V (±5.0 V) 100 mA, max. <0.4 µV RMS (2 µV p-p) 0.1-100 Hz Gain >10 <3 µV RMS (15 µV p-p) 0.1-100 Hz Gain >10 0.03, 0.1, 0.3, 1, 5, 10,

"Wide Band" R-C Butterworth 6 dB /octave)

74020	Iso-DAM8A Single Channel Module
74030	ISDB chassis and power supply
74040	Iso-DAM8A Active Headstage (separate)
BRIDGE8	Bridge8 Transducer Amplifier Module
OPTIONAL A	ACCESSORIES/REPLACEMENT PARTS
74050	ISDB Blank panels
74016	Replacement Cable, Input bare 5-ft wire
2933	ISDB Rack Mount Kit, 5¼-in. High
FORT10g	Force Transducer 10 g
FORT25	Force Transducer 25 g
FORT100	Force Transducer 100 g
FORT250	Force Transducer 250 g
FORT1000	Force Transducer 1000 g
500184	BNC to BNC 10 foot cable
3161	8-pin DIN plug
3718	Package of 4, 8-pin DIN (startup kit)
3491	Extension Cable (DIN male, DIN female), 5 ft (1.5 m)

ORDERING INFORMATION

Dual Microprobe Intracellular Amplifier

2-Channel intracellular amplifier for dual and differential studies



Features

- · Two channels for differential or intracellular ISE
- Built in DC current generator with external control input
- Built in low pass filter
- Bridge balance circuit to null out electrode voltage drop
- Tickle circuit
- Built in test ports for each channel
- Dual capacitance compensations and output offset controls

Benefits

- Dual channel, single ended recording
- Differential recording
- Bridge circuit nulls electrode voltage drop
- Assign low pass filter to either channel
- Very high impedance channel can be used with intracellular ISE

Applications

- Intracellular electrophysiology using sharp micropipettes
- · Brain slice intracellular recording
- In vivo intracellular recording from brain and spinal cord

For intracellular dual or differential studies, the **Duo773** has separate negative capacity controls and built-in active filtering that allows the precise balancing of time constants for artifact-free differential measurement. Comes complete with two probe headstages, $10^{15}\Omega$ and $10^{11}\Omega$ probes to monitor signals from ion-specific micro-electrodes as well as KCl-filled electrodes.

Headstage for precise positioning

Two gold-plated, epoxy sealed miniature active probes can be positioned directly to the measurement site. Microelectrode holders containing an Ag/AgCl electrochemical half-cells plug directly into the probes. Stray capacitance can be reduced by placing the included driven guard shield over the microelectrode holder at the end of the probe.

Capacity compensation

Channel A can compensate up to 10 pF of electrode shunt capacity and Channel B can compensate up to 50 pF.

Tickler circuit for penetration

A Tickler Circuit assists in cell penetration. The frequency and amplitude of the oscillations may be varied for differences in membrane thickness or cell size. The duration of tickle can be controlled either by using the momentary switch, a foot switch, or by applying a signal to the remote tickler input.

Active filters

Low pass settings on a -40 dB/decade active filter vary the cutoff from 1 to 30 kHz. Either probe or bridge outputs may be selected for filtering.

Current injection

Channel B can eject current through the microelectrode by applying a command signal to the stimulus input connector. The resulting output from the probe will be a constant current replica of the input signal. Two ranges of current delivery are provided: 50 nA and 500 nA or by an external source. This source can be useful for delivering hyperpolarizing currents to stabilize the cell membrane potential and as a holding current for microiontophoresis.


Bridge balance

Subtracts the excess electrode voltage associated with delivering current through the recording micropipette. Electrode resistances up to 1000 $M\Omega$ can be balanced in two ranges. The balanced signal is available from x10 or x50 front panel output connectors.

Independent outputs


The **Duo773** has an output for each probe independent of gain filtering or balancing. In addition the **Duo773** has a 10x and a 50x output for easy integration to most data acquisition programs.

Typical setup



See cables and connectors, page 96 SeeDri-Ref, page 158.

Optional holders for intracellular amplifiers



See Microelectrode Holders, page 98.

References

Huo, Q., Chen, M., He, Q., Zhang, J., & Li, B. (2016). Prefrontal Cortical GABAergic Dysfunction Contributes to Aberrant UP-State Duration in APP Knockout Mice. Cerebral Retrieved from http://cercor.oxfordjournals.org/content/early/2016/08/23/cercor.bhw218.short

Spong, K. E., Rodríguez, E. C., & Robertson, R. M. (2016). Spreading depolarization in the brain of Drosophila is induced by inhibition of the Na+/ K+-ATPase and mitigated by a decrease in activity of protein kinase G. *Journal of Neurophysiology*, jn.00353.2016. http://doi.org/10.1152/jn.00353.2016

Zhang, J., Chen, M., Li, B., Lv, B., Jin, K., & Zheng, S. (2016). Altered striatal rhythmic activity in cylindromatosis knock-out mice due to enhanced GABAergic inhibition. Retrieved from http://www.sciencedirect.com/science/article/pii/S002839081630274X

Bredeloux, P., Finday, I., & Pasqualin, C. (2016). 0194: Functional consequences of-adrenergic receptors activation in the rat pulmonary veins and left atria. Archives of Retrieved from http://www.sciencedirect.com/science/article/pii/S1878648016304311

DUO 773 S	PECIFICATIO	VS
HEADSTAGE (PROBE)	712P (red, port "B")	715P (blue, port "A")
ACTIVE PROBE INPUT IMPEDANCE	>10 ¹¹ Ω	1015 Ω
GAIN	x1, x10	x1
OUTPUT RESISTANCE	100 Ω	100 Ω
OUTPUT VOLTAGE RANGE	±10 V	± 10V
MAXIMUM INPUT VOLTAGE	±15 V	±15 V
PROBE LEAKAGE CURRENT	5 X 10 ⁻¹² A	10 ⁻¹⁴ A
DC POSITION ADJUST RANGE	± 300 mV	± 300 mV
ELECTRODE RESISTANCE TEST CURRENT	1 nA	1 pA, 1 nA selectable
INPUT CAPACITY COMPENSATION	+10 to -50 pF	0 to -10 pF
NOISE		
Input shorted	<50 μV p-p 10kHz bandwidth	<50 μV p-p 10kHz bandwidth
20 MΩ carbon resistor	<200 µV p-p 10kHz bandwidth	<200 µV p-p 10kHz bandwidth

RISE TIME

SYS-773

10-90% direct input small signal 1 µs, typical

10-90% through 20 MΩ (-C "on") 25 μs, typical		
(712P only)**		
± 50 nA low range, ± 500 nA high range		
± 500 nA low range, ± 5 μA high range		
20 mV/nA low range, 2 mV/nA high range		
100 mV/nA low range, 10 mV/nA high range		
3V low range, 10V high range		
0-100 ΜΩ, 0-1000 ΜΩ		
x 10, x 50		
40 dB/decade, continuously variable 1-30 kHz		
3.5-digit LED		
200 mV, 2000 mV, 20 V, 200 nA, 2000 nA		
1 digit		
17 x 5.25 x 10 in. (43 x 13 x 25 cm)		
Diameter: 12 mm Length: 34 mm		
05 425 // - 220 240 // 50/60 //		
95-135 V or 220-240 V, 50/60 Hz		
15 lb (7 kg)		

^{*} Although injected currents are "constant," the maximum current in a given situation will always be limited by the system compliance of 10 V.

ORDERING INFORMATION

OPTIONAL ACCESSORIES/REDI ACEMENT DARTS	
electrodes. Spe	ecify line voltage
	H1SF microelectrode holders for 1.0, 1.2, 1.5 or 2.0 mm glass
Includes two p	robes (712P and 715P or two 712P) with driven guard shields
	Bao 775 Electronnete

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

Duo 773 Electromete

712P	Replacement probe (includes calibration)*
715P	Replacement probe (includes calibration)*
*Instrumer	t should be returned to WPI for free calibration with new probe.
2933	Rack Mount Kit, 5¼-in. high
2547	Driven Guard Shield for 712P & 715P Probes
15790	Replacement Probe Handle
TW100F-4	Glass capillary with filament
TW150F-4	Glass capillary with filament
3259	Foot Switch

^{**}The 712P headstage may be used on either A or B channels, however Current Injection specifications do not apply when used on channel A. The 715P headstage may not be used on the B channel.

Dual Channel Differential Electrometer

Electrochemical measurements with ion specific electrodes

Features

· Dual channel with very high input impedance

 Separate outputs for Channel A, B and A-B (Differential)

- Independent DC offset controls
- Test port
- Standby mode

Benefits

- Measure changes in intracellular ion content electrochemically
- · Stable and drift free
- Excellent amplification with low noise
- Driven guard shield for reduced noise and stray capacitance
- Set probe leakage current

Applications

 Measure intracellular ion concentrations for K+, Ca2+, H+ and other

The **FD223a** electrometer was designed specifically for use with intracellular ion selective electrodes fabricated using glass micropipettes and liquid ion exchangers.

The active head stages allow the researcher to locate the probes directly at the measurement site to minimize noise that would normally be picked up by longer cable runs.

Driven guard shields cover the micropipette holders to further reduce the potential for interference from external sources of electromagnetic noise.

The **FD223a** is equipped with a test resistance port which is used to measure and adjust each probe for minimum leakage current. Each channel has a standby mode which clamps the head stage input voltage to zero, preventing extreme saturation or possible damage to the high impedance input amplifier.

References

E. Ermolayeva, H. Hohmeyer, E. Johannes, D. Sanders

"Calciumdependent membrane depolarisation activated by phytochrome in the moss Physcomitrella patens" *Planta* 199. 1996: 352-358

ORDERING INFORMATION

FD223A FD223a Dual Channel Differential Electrometer

2 probes, driven guard shields and micropipette holder MEH1SF included for all glass microelectrodes O.D. 1.0, 1.2, 1.5 or 2.0 mm.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

M3301L	Left-hand Micromanipulator
M3301R	Right-hand Micromanipulator
M-3	80° Tilting base
RC1T	Reference cell (Ag/AgCl)
2547	Driven guard shield for FD223AP Probe
MEH1SF	Microelectrode holder
FD223AP	Replacement probe (includes calibration)

See cables and connectors, page 96. See microelectrode holders, page 101. See capillary glass, page 104.





FD223A SPECIFICATIONS

RISE TIME (10 TO 90%) 5 µs, small signal

NOISE (0.1 Hz TO 10 KHz) SASELINE STABILITY < 100 μ V p-p, input shorted \pm 0.1 mV/day

POSITION CONTROLS RANGE ±600 mV

PHYSICAL DIMENSIONS Case: 8.8 x 21.0 x 17.5 cm (H x W x D)

Probe: 12.7 x 65 mm (D x L), 1.8 m cable

POWER 90-265 VAC, 50/60 Hz, 10 VA PROBE HANDLE 6.5 x 65 mm (D x L)

SHIPPING WEIGHT 2.5 kg

OPERATING CONDITIONS: Equipment is intended to be operated in a controlled laboratory environment. Temperature: 0-40 °C; altitude: sea level to 2000 m;

relative humidity: 0-95%.



2547 Driven Guard Shield

Ultra Quiet Intracellular Amplifier

High quality, intracellular amplifier perfect for students

Features

- · Driven guard shield
- Test port
- Ground port
- Portable
- Remote headstage

Benefits

- Cost effective
- Battery powered
- Capacitance compensation

Applications

Measure intracellular action potentials

Electro 705, a battery operated, low noise, wide band electrometer preamplifier, is designed for intracellular voltage measurement. Two 705's can be linked together to form a

high impedance differential electrometer pair. Each instrument includes a miniature gold plated active probe to which a microelectrode can be attached using the WPI microelectrode holder supplied.

Remote headstage

Easily mounted in any manipulator, this compact probe, containing the first stage of amplification, includes a microelectrode holder, which plugs directly into the probe input.

Battery power

Four 9V alkaline batteries (included) power the **Electro 705** for approximately 500 hours giving a clean, low noise source of power, making the **Electro 705** the quietest amplifier available. Batteries can be easily tested by the press of a button.

Capacitance Compensation

Corrects for loss of rise time caused by the presence of electrode capacity. Up to 50 pF of electrode shunt capacity may be neutralized.

Driven Guard Shield

Stray capacitance can be further reduced by placing the driven guard shield (included) over the microelectrode holder at the input end of the probe.

Test Features

A Tickler Circuit offers a momentary oscillation that helps achieve cell penetration. The **Electo 705**

provides a 1 nA electrode test current. Electrode resistance is monitored at the 1X output as a voltage (1 mV/M). The Probe Test Port allows the convenience of testing the amplifier's intrinsic noise and gain without cumbersome external test hookups. Head stage leakage current can also be adjusted with minimum effort. The Baseline Position Control adds or subtracts up to 300 mV to the headstage output, allowing artifact voltages such as liquid junction potentials to be nulled prior to recording.

Differential Output

Two **Electro 705** units can be connected in tandem to create an optional differential amplifier probe system.

References

Wan, E., Kushner, J. S., Zakharov, S., Nui, X.-W., Chudasama, N., Kelly, C., ... Marx, S. O. (2013). Reduced vascular smooth muscle BK channel current underlies heart failure-induced vasoconstriction in mice. *FASEB Journal : Official Publication of the Federation of American Societies for Experimental Biology*, 27(5), 1859–67. http://doi.org/10.1096/fj.12-223511

Gokina, N. I., Bonev, A. D., Gokin, A. P., & Goloman, G. (2013). Role of impaired endothelial cell Ca2+ signaling in uteroplacental vascular dysfunction during diabetic rat pregnancy. *American Journal of Physiology - Heart and Circulatory Physiology*, 304(7).



Thomas, R. C., & Bers, D. M. (2013). How to make calcium-sensitive minielectrodes. *Cold Spring Harbor Protocols*, 2013(4), 370–3. http://doi.org/10.1101/pdb.prot072850

ELECTRO 705 SPECIFICATIONS

Probe Input Resistance $10^{12} \, \Omega$ Probe Input Leakage Current $\pm 10 \, \text{pA}$, adjustable to 0
Input Capacitance Compensation 0 to 50 pFNoise level 500 uV peak to peak
Rise time 15 uS, 10 - 90%Output impedance 100 ohms, both outputs
Tickle 9V p-p adjustable from 200 Hz to 3000 Hz

Position Range ±300mV

Electrode R Test 1mV/M ohm

Common Mode Rejection > 104, in pair

Common Mode Rejection > 104, in paired operation

Power Four 9V alkaline cells, supplied

Dimensions 21.6 x 8.9 x 14.3cm (8.5 x 3.5 x

ORDERING INFORMATION

SYS-705	Electro 705 Electrometer
Probe, driven	guard shield and micropipette holder MEH1SF included for
alace microel	ectrodes 0.0 1.0 mm 1.2 mm 1.5 mm or 2.0 mm

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

M3301L	Left-hand Micromanipulator
MM3301R	Right-hand Micromanipulator
M-3	80° Tilting base
RC1T	Reference cell (Ag/AgCl)
2541	Driven guard shield for 705PF Probe
MEH1SF	Microelectrode holder
15804	Replacement Probe Handle & Clip
705PF	Replacement probe (includes calibration)*
*100tru0000	at payed he returned to MDI for free calibration with new probe

*Instrument must be returned to WPI for free calibration with new probe.

See cables and connectors, page 96.

See microelectrode holders,page 98.

See capillary glass, page 104.

^{*} Dependent on measuring technique.

Single Channel Pulse Generator

The accuracy of digital electronics and convenience of analog controls



Features

- · Single channel pulse generator with train capability
- TTL and variable voltage output

Benefits

- Variety of pulses: continuous run, single-shot, train/burst
- Multiple outputs available: monitor, isolator, sync and variable

Applications

Electrophysiology

The **A310** pulse generator/stimulator combines the reproducibility and accuracy of digital electronics with the fine resolution and continuous adjustment possible with analog circuitry. All timing parameters are entered with high resolution, ten-turn potentiometers and six-position range switches. Timing is accurate to within 1% of the set value.

Variety of Pulses

Pulses can be created in continuous run, single-shot or train/burst modes. Duration of the train/burst is controlled using the onboard envelope generator or by using either of two external gating inputs. Used in conjunction with the A360, A365, A385 or A395, constant current pulses and trains can be created easily. A foot switch allows hands-free, manual triggering.

Multiple Outputs Available

Five separate standard BNC outputs are available on the front panel. The isolator output sends full pulse width control signals to any TTL triggered stimulus isolator, such as WPl's A360, A365 or A385 and others. The monitor output sends synchronized large scale full pulse width signals to recording or monitoring instrumentation such as a data acquisition system or oscilloscope. The sync output provides an additional synchronized 5 μs TTL pulse for triggering external instrumentation. A variable voltage output provides two separate full pulse width signals in both positive and negative polarities in two ranges for applications that require a specific output voltage other than TTL.

References

Cha, R., Marescaux, J., & Diana, M. (2014). Updates on gastric electrical stimulation to treat obesity: Systematic review and future perspectives. *World Journal of Gastrointestinal Endoscopy*, 6(9), 419–31. http://doi.org/10.4253/wjge.v6.i9.419

SPECIFICATIONS

TIMING PARAMETERS

 EVENT INTERVAL
 100 μs to 1000 s*

 EVENT DELAY
 10 μs to 100 s *

 PULSE WIDTH
 10 μs to 100 s *

 TRAIN DURATION (ENVELOPE)
 100 μs to 1000 s*

 PULSE INTERVAL
 20 μs to 100 s*

OUTPUTS

SYNC 5 µs, TTL, and 5 V CMOS compatible,

20 mA max.

MONITOR 10-15 V, 50 mA max.

ISOLATOR TTL & 5 V CMOS compatible, 20 mA max.

VARIABLE (Pos or Neg)

 PULSED/DC
 LOW RANGE
 HIGH RANGE

 Range
 0 to ±1 V
 0 to ±10 V

 Resolution
 1 mV
 10 mV

NOISE

Pulsed at 100 kHz bandwidth $<500~\mu V$ $>500~\mu V$ $<500~\mu V$ OUTPUT IMPEDANCE $<1~\Omega$

INPUTS

EXTERNAL SYNC Accepts 1-µs minimum pulses TTL, CMOS compatible

EXTERNAL GATE Accepts 1-µs pulse to continuous

TTL, CMOS compatible

POWER 95-130 V or 190-260 V, switch selectable

single phase, 50/60 Hz

DIMENSIONS 17 x 5.25 x 10 in. (43 x 13 x 25 cm)

SHIPPING WEIGHT 14 lb (6.4 kg)

*Continuously variable in six ranges. All accuracies better than 1% of set value. 50kHz maximum pulse frequency.

ORDERING INFORMATION

SYS-A310 Accupulser™ Signal Generator

Specify line voltage

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

3259	Foot Switch for A310	
2933	Rack Mount Kit, 5¼ in. high	

Constant Current Stimulus Isolator

Automated bipolar pulsing for zero net charge on biological preparations

Features

- · Constant current
- Unipolar and bipolar stimulation modes
- Built-in non-compliance alarm
- · Input is optically isolated
- · Standard TTL triggering
- DC test mode
- Powered by 9V alkaline or rechargeable batteries

Benefits

- Compliance voltage is 100V or better
- Bipolar mode auto generates alternating positive and negative pulses from TTL input
- Test mode simplifies performance verification
- Optical isolation enhances safety of the preparation and reduces noise susceptibility



Applications

- Electrophysiology
- Brain slice stimulation
- · In vivo brain and CNS stimulation

Activated by conventional logic-level commands, Model **A365** can be gated by any pulse generator, stimulator, or computer output; automated bipolar pulsing for zero net charge on biological preparations.

Dual Tone Audible Alarms

A tone sounds when an open electrode circuit is detected or when system compliance is reached. A second optional tone sounds when a signal is applied to the input. A test switch is also provided to check battery charge.

Current Delivery up to 10 mA at More than 100V

Stimulus currents are set using a three-digit control knob and a three-position range switch. Output current tracks control settings to better than 1%. Output current is load independent; voltage sufficient to push the desired current through the load is automatically developed, subject only to compliance limits. Model **A365** produces up to 10 milliampere current, in three ranges, at more than 100 volts compliance.

Bipolar Output Polarity

Output polarity is determined by a push switch on the front panel. Bipolar current is toggled by the command waveform, setting alternating pulses as positive or negative.

References

Lee, E., Hong, J., Park, Y.-G., Chae, S., Kim, Y., & Kim, D. (2015). Left brain cortical activity modulates stress effects on social behavior. **Scientific Reports**, 5, 13342. http://doi.org/10.1038/srep13342

Gindrat, A.-D., Quairiaux, C., Britz, J., Brunet, D., Lanz, F., Michel, C. M., & Rouiller, E. M. (2015). Whole-scalp EEG mapping of somatosensory evoked potentials in macaque monkeys. *Brain Structure & Function*, 220(4), 2121–42. http://doi.org/10.1007/s00429-014-0776-y

Younce, J. R., Albaugh, D. L., & Shih, Y.-Y. I. (2014). Deep Brain Stimulation with Simultaneous fMRI in Rodents. *Journal of Visualized Experiments*, (84), e51271–e51271. http://doi.org/10.3791/51271

Avila, I., & Lin, S.-C. (2014). Motivational Salience Signal in the Basal Forebrain Is Coupled with Faster and More Precise Decision Speed. *PLoS Biology*, 12(3), e1001811. http://doi.org/10.1371/journal.pbio.1001811

Herrera, C., Directores, R., Panetsos, F., Carlos, P., & Trueba, A. (2014). TESIS DOCTORAL Efectos de la estimulación artificial de un nervio periférico seccionado sobre la vía somatosensorial desaferentizada de la rata.

A365 SPECIFICATIONS

OUTPUT WAVEFORM
OUTPUT CURRENT RANGES
CURRENT AMPLITUDE ERROR
CURRENT RESOLUTION
OUTPUT LOAD VOLTAGE
EXCURSION (COMPLIANCE)
EXTERNAL COMMAND VOLTAGE
TRIGGER THRESHOLD
OUTPUT POLARITY
CURRENT RISE TIME & DELAY
CURRENT FALL TIME & DELAY

POWER

Model A365D (Dry Cell)

Model A365R (RECHARGEABLE)

DIMENSIONS

SHIPPING WEIGHT

OPTOCOUPLER

DC or current pulse 0.1, 1.0, and 10 mA 0.5% of full scale, max. 0.1% of full scale, typical

100 V (typically 144 V)

 $\begin{array}{lll} \text{EXTERNAL COMMAND VOLTAGE} & 5.0 \text{ V at } 3.0 \text{ mA (TTL level), } 10 \text{ V max.} \\ \text{TRIGGER THRESHOLD} & 2.0 \text{ V at } 0.5 \text{ mA} \\ \text{OUTPUT POLARITY} & \text{Reversible, manual switch or automatic} \\ \text{CURRENT RISE TIME \& DELAY} & 6 \ \mu\text{s, typical (1 K}\Omega \ \text{load)} \\ \text{CURRENT FALL TIME \& DELAY} & 10 \ \mu\text{s, typical (1 K}\Omega \ \text{load)} \\ \text{OUTPUT TO GROUND RESISTANCE} & 10^{12}\Omega \\ \end{array}$

2500 V, rated min. breakdown voltage

16 alkaline 9 V batteries, included 16 rechargeable NiMH 9 V batteries incl. 8.5 x 3.5 x 5 in. (22 x 9 x 12 cm) 4 lb. (1.8 kg)

	ORDERING INFORMATION
SYS-A365D	High Voltage Isolator, Bipolar, alkaline batteries
A365RC	A365R with charger (A362)
SYS-A365R	High Voltage Isolator, Bipolar, rechargeable
SYS-A362	Battery Charger for A320R, A365R, A395R
	Specify line voltage

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

DRL	Dummy Load Resistor Kit (set of 3)
13347	BNC-to-Double Banana Adapter

A362 Battery Charger

Required for A320R, A365R, A395R

Recharges the high-voltage nickel-cadmium or NiMH battery stack in the A320R, A365R or A395R. LED lamp indicates charging status. Full charge overnight. Dimensions: $2.8 \times 4.1 \times 5$ in. ($7 \times 10 \times 13$ cm). Shipping weight: 4 lb (1.8 kg).



High Current Stimulus Isolator

Constant current stimulus isolator with 100 mA current range

Features

- · Constant current to 100 mA
- Unipolar or bipolar stimulation modes
- Built-in non-compliance alarm
- Input is optically isolated
- Standard TTL triggering
- · DC test mode
- Powered by six rechargeable lead acid
- 36V compliance
- Output polarity and output "on/off" switches

Benefits

- 100 mA current capability
- Bipolar mode automatically generates alternating positive and negative pulses from TTL input
- Test mode simplifies performance verification
- Optical isolation enhances safety of the preparation and reduces noise susceptibility
- Battery charge status LEDs keep the experimenter informed of battery status
- Charger included at discounted price when system is purchases as

Applications

- Muscle electrophysiology
- In vivo/in vitro muscle stimulation

The A385 is an optically isolated current source, which can generate up to 100 mA of unipolar or biphasic constant current pulses or DC. Pulse duration is controlled manually or by an external 5V command. Output current amplitude is determined by a 3-digit 10-turn potentiometer. Maximum output voltage between the stimulating electrodes is +36V.

Delivers Positive, Negative or Bipolar Currents

For bipolar delivery, the polarity of the output can be toggled to the opposite polarity state with each successive pulse presented to the input. Pulse duration is controlled by an externally applied voltage. The input connector is a standard BNC, allowing TTL signals from a data acquisition system to be used.

Excellent Accuracy and Repeatability

The output amplitude is controlled by a 3-digit, ten-turn dial as a percentage of the range selected: for example, a setting of 45.6 in the 0-10 mA range translates to 4.56 mA at the output. Accuracy and repeatability are excellent. Designed for subcutaneous stimulation, maximum output voltage at the stimulating electrodes is 36 volts, reducing the possibility of accidental transcutaneous shocks. A compliance/output alarm sounds when the 36V limit is reached. Internal circuitry ensures electrodes are short-circuited during inactive periods ("electrode exhauster" feature). The ${\bf A385}$ is not appropriate for transcutaneous stimulation.

Rechargeable Battery

The 1.2 amp-hour rating of the six heavy-duty lead-acid rechargeable batteries ensures that all day experiments will not be interrupted by dead batteries when charged daily. Indicator lights and audible alarms keep the user constantly apprised of the battery charge status. The batteries are recharged by the A382 System Charger, which is designed especially for the A385, and included with the A385RC.

References

Lin, C., Disterhoft, J., & Weiss, C. (2016). Whisker-signaled Eyeblink Classical Conditioning in Head-fixed Mice. Journal of Visualized Experiments, (109), e53310-e53310. http://doi.org/10.3791/53310

Li, T., Finch, E. A., Graham, V., Zhang, Z.-S., Ding, J.-D., Burch, J., ... Rosenberg, P. (2012). STIM1-Ca(2+) signaling is required for the hypertrophic growth of skeletal muscle in mice. Molecular and Cellular Biology, 32(15), 3009-17. http://doi.org/10.1128



A385 SPECIFICATIONS

OUTPUT WAVEFORM OUTPUT CURRENT RANGES CURRENT AMPLITUDE ERROR CURRENT RESOLUTION REPEATABILITY

OUTPUT LOAD VOLTAGE EXCURSION (COMPLIANCE) EXTERNAL COMMAND VOLTAGE

EXTERNAL COMMAND VOLTAGE: **OUTPUT POLARITY**

EXTERNAL COMMAND VOLTAGE: CURRENT RISE TIME AND DELAY EXTERNAL COMMAND VOLTAGE: CURRENT FALL TIME AND DELAY EXTERNAL COMMAND VOLTAGE: **OUTPUT TO GROUND RESISTANCE** EXTERNAL COMMAND VOLTAGE: **OPTOCOUPLER**

POWER

A385RC

SYS-A385R

SYS-A382

DIMENSIONS SHIPPING WEIGHT DC or current pulse 1, 10, and 100 mA 0.5% of full scale, max.

0.1% of full scale, typical

36 V (typically 40 V)

5 V at 3 mA minimum, 8.5 V max. Reversible, manual switch, monophasic or electronically switched bipolar delivery

6 μs, typical (1 KΩ load)

10 μs, typical (1 KΩ load)

 $10^{12}\,\Omega$

2500 V, rated minimum breakdown

Six rechargeable lead-acid batteries (Requires companion charger A382) 8.5 x 3.5 x 5 in. (22 x 9 x 12 cm)

ORDERING INFORMATION A385R with A382 Charger High Current Isolator, rechargeable Battery Charger for A385 (see below)

Smart Battery Charger

Specify line voltage

Required for A385

An innovative three-step charger, the A382 employs fast, medium, and trickle charges at a safe, low current, greatly extending battery life. After a fast initial phase, the charger automatically switches to a constant

voltage mode. When charging

is complete, the charger switches to the trickle-charge mode. LED lamps indicate charging status. (For use only in charging the lead acid batteries installed in the A385.)

Linear Stimulus Isolator

Replicates a programmed waveform of any shape or polarity

Features

- Creates a constant current replica of analog waveforms
- Amplitude of the output current is voltage controlled
- Input voltage from –10V to +10V
- 3 current ranges from $100 \mu A$ to 10 mA
- Built-in test resistors
- Digital display shows current being delivered for non-varying currents of adequate duration
- Output offset adjustment
- ±70V compliance range



Benefits

- Amplitude of current is voltage controlled
- Built-in test resistances
- Error LEDs illuminate when current is less than commanded by control voltage

Applications

- Neuroscience
- Muscle Physiology

All WPI stimulus isolators are designed to supply constant current, because current threshold (not voltage) is the most quantitatively reproducible parameter for stimulation of nerve and muscle. Model A395 dispenses current reproducibly from its Output terminals; the amplitude being determined by the selected current RANGE and the input voltage. Current amplitude is "constant", that is, load resistance independent, provided that the I x R (load) product does not exceed the available battery supply voltage. A visual indicator (the compliance LEDs) displays if Ix R reaches this limit. When the unit is out of compliance, one of the two LEDs (labeled - and +) illuminate, depending in which direction the current is flowing. Model **A395D** can generate a voltage of 70Vor more across its OUTPUT terminals. You can be sure that the amplitude of the current is as dialed as long as the voltage drop across the load (stimulus electrode path) does not reach the magnitude of the supply voltage. The compliance LEDs will then be visible. Then, you would know that:

- · Too much current was dialed for a given load or
- · Inter-electrode resistance was too high or the electrode circuit path was

User Defined Output Current of Various Forms

Model A395 generates a user-defined output current of wave shape; DC, AC, pulse and combinations. Battery operated, photoelectrically-isolated from the input voltage drive, the instrument regenerates output currents which are linearly proportional to the analog voltage waveforms provided by your D/A converter or signal generator (see diagram below).

The A395 is ideally suited for data acquisition and stimulator generators.

Current Delivery for Selected Ranges

A 10 V input produces the maximum output current for the current range selected. (For example, 100 µA, 1 mA, or 10 mA) Front panel controls allow DC current to be generated. Externally applied signals can be



superimposed simultaneously (DC offset). Warning lamps indicate open circuit or excessive current conditions.

Digital Meter Shows DC or Average Output

The digital display meter shows the measures DC current or the average output current. Overload lamps indicate when output voltage has reached positive or negative compliance voltage limit.

References

E.D. Zonnevijille, N.N. Somia, g.Perez Abadia, R.W. Stremel, C.J. Maldonado, P.M.N. Werker, M. Kon, J.H. Barker "Sequential Segmental Neuromuscular Stimulation Reduces Fatigue and Imrproves Perfusion in dynamic Gracilolasty" Ann Plast Surg 45. 2000: 292-297

A395 SPECIFICATIONS

OUTPUT CURRENT, Imax 3 ranges: 100 μA, 1 mA, and 10 mA

OUTPUT VOLTAGE RANGE ±70 V

OUTPUT BANDWIDTH 10 kHz (measured across 1KΩ load R)

INPUT RESISTANCE $>20 \text{ m}\Omega$ INPUT VOLTAGE @ Imax ±10 V INPUT/OUTPUT LINEARITY ERROR < 0.5%

RISE, FALL TIME 26 μs @ 10 ΚΩ POWER: Model A395D 17 alkaline 9 V batteries

17 rechargeable NiMH 9 V batteries POWER: Model A395R **DIMENSIONS** 6.5 x 4 x 3.5 in. (16 x 10 x 9 cm)

SHIPPING WEIGHT 4 lb. (1.8 kg)

	ORDERING INFORMATION
A395RC	A395R with Charger (A362)
SYS-A395D	Linear Stimulus Isolator
SYS-A395R	Linear Stimulus Isolator, Rechargeable
SYS-A362	Battery Charger

Specify line voltage

A362 Battery Charger

Required for A320R, A365R, A395R

Recharges the high-voltage nickel-cadmium or NiMH battery stack in the A320R, A365R or A395R.

LED lamp indicates charging status. Full charge overnight.

Dimensions: 2.8 x 4.1 x 5 in. (7 x 10.5 x 12.7 cm). Shipping weight: 4 lb. (1.8 kg).



4-Channel Transducer Amplifier

Amplify output voltage signals

Features

- Use with many different types of resistive based transducers
- WPI resistive force transducers plug in directly
- Supports full resistive bridge or single ended operation
- Output offset control
- Four gain ranges from 1–1000 ×
- Provides ±5V "excitation" voltage for resistive bridge transducers

Benefits


- Provided with blank connectors to rewire to any resistive bridge transducer
- Bridge balance LEDs provide visual cue that unloaded transducers are at zero output state
- · Low noise amplifier

Applications

 Amplify signals from resistive strain gages and other resistive bridge configured transducers

Transbridge (**TBM4M**) is a four-channel analog transducer manifold, specifically designed to amplify output voltage signals from pressure, force, displacement, and temperature transducers as well as a wide variety of other signal sources. Analog output signals are available from each channel for input to a data acquisition system for digital signal processing in a computer. Each channel contains a regulated 10-volt power supply (+5 and -5 volts with respect to signal ground) to provide DC power to transducers, and a precision differential amplifier with selectable voltage amplification and variable position adjustment control.

Transducers can be connected to Transbridge via any of the 8-pin connectors on the front panel. Four spare 8-pin DIN plugs are provided with each instrument to allow you to rewire cables of other manufacturers' transducers and connect them to Transbridge. Each Transbridge channel may be used in either Full Bridge or the Half Bridge mode independently. For transducer types other than resistive bridges, such as active transistor circuits, magnetic, photocell or piezoelectric



devices, the instrument's differential amplifiers may still be used effectively for signal amplification in differential (full bridge) and single-ended (half bridge) modes.

	ORDERING INFORMATION
SYS-TBM4M	Transbridge Transducer Amplifier
	Specify line voltage

OFIIONA	ACCESSORIES/REPEACEMENT PARTS
13024	Single Rack Mount Kit
13025	Dual Rack Mount Kit
500184	BNC-to-BNC cable, 10 ft
3161	8-pin DIN plug
3718	Package of 4, 8-pin DIN (startup kit)

Battery Operated Impedance Measurement

Measure mV and $M\Omega$ impedance of metal or glass microelectrodes

Features

• Battery operated

Benefits

 Determine impedance of electrode during micropipette beveling process for pipette reproducibility

Applications

• Measure impedance of metal or glass capillary microelectrodes

ORDERING INFORMATION

SYS-OMEGAZ Omega-Tip-Z™ with Probe & Holder

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

711P Replacement Probe
5468 Adapter to connect metal microelectrodes to probe,
2 mm socket to .031 in. receptacle

Omega-Tip-Z™ was created especially for measuring impedance in etched tungsten, platinum-iridium* and steel microelectrodes, as well as electrolytefilled micropipettes. The meter's AC impedancemeasuring circuit is unaffected by electrode offset or tip junction potentials. The goldplated miniature probe lets you conveniently monitor microelectrode impedance in electrolytes, and an electrode tip cleaning feature lets you



remove buildup quickly. **Omega-Tip-Z** can also measure DC electrode tip potentials up to 2000 millivolts. The instrument operates for hundreds of hours without battery failure.

*See Metal Microelectrodes, page 94.

Force Transducers

These rigid-lever force transducers transform applied force into proportional voltage. Using balanced strain gages, **FORT** transducers produce linear output voltage vs. applied force input with very little deflection.

To use, clamp the handle of the **FORT** transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.

FORT 5000


FORT 1000

FORT 250

FORT 100

	ORDERING INFORMATION
FORT100	Force Transducer (100 g)
FORT250	Force Transducer (250 g)
FORT1000	Force Transducer (1000 g)
FORT5000	Force Transducer (5000 g)

FORT SPECIFICATIONS							
	FORT100	FORT250	FORT1000	FORT5000			
FORCE RANGES, FULL SCALE	100 g	250 g	1000 g	5000 g			
OUTPUT SENSITIVITY (± 10%)	7 μV/V/g	3 μV/V/g	0.84 μV/V/g	0.38 µV/V/g			
INPUT & OUTPUT RESISTANCE	350 Ω	350 Ω	350 Ω	350 Ω			
RESOLUTION	0.01% of full scale force	0.01% of full scale force	0.01% of full scale force	0.1% of full scale force			
RESONANT FREQUENCY	300 Hz	300 Hz	300 Hz	60 Hz			
LINEARITY ERROR	Less than 0.1% of full scale						
MAX. OPERATING VOLTAGE	10 V AC or DC						
MAXIMUM APPLIED FORCE	3× rated full scale force	3× rated full scale force	3× rated full scale force	3× rated full scale force			
DRIFT	thermally compensated	thermally compensated	thermally compensated	thermally compensated			
DIMENICIONIC	0.3 inch diam × 4 in.						
DIMENSIONS	(7.6 mm diam × 10.2 mm)						
WEIGHT (excluding cable)	0.3 oz (8 g)						



10g & 25 Force Transducers

FORT SPECIFICATIONS						
	FORT10g	FORT25				
FORCE RANGE, FULL SCALE	0-10 g	0-25 g				
OUTPUT SENSITIVITY	2.75 mV/g, nominal	3 mV/g, nominal				
INPUT & OUTPUT RESISTANCE	4000 Ω	1500 Ω				
RESOLUTION	< 1 mg	< 2 mg				
RESONANT FREQUENCY	480 Hz	450 Hz				
LINEARITY ERROR	<0.2% of full scale	<0.2% of full scale				
MAXIMUM OPERATING VOLTAGE	10 V DC (-5V ~ +5V or 0 ~ 10V)	10 V DC (-5V ~ +5V or 0 ~ 10V)				
MAXIMUM APPLIED FORCE	2× rated full scale force	3× rated full scale force				
DRIFT	<30 mg/hr	<50 mg/hr				
DIMENSIONS	46 x 22 x 15 mm Handle 100 mm x 9.8 mm diam	40 x 22 x 17 mm Handle 109 mm x 9.8 mm diam				
WEIGHT	135 g	200 g				

These 10 g and 25 g force transducers are reliable tools for high precision force measurement. Using balanced semiconductor strain gages, both produce linear output voltage vs. applied force input with very little deflection. The rigid lever force transducer transforms the applied force into a proportional voltage. Featuring a temperature-compensated, full-bridge configuration with four high sensitivity semiconductor strain gages, these transducers have broad dynamic measuring range and very high sensitivity.

To use, clamp the handle of the **FORT10** or **FORT25** transducer in a horizontal position and apply the forces to be measured to a rivet or hook mounted in the hole at the end of the flat sensing leaf.

Metal Microelectrodes

Superior microelectrodes for outstanding extracellular recording— tungsten, iridium, platinum-iridium and Elgiloy®

Features

- Available in Tungsten, Platinum/Iridium, Elgiloy and Pure Iridium metal, and are insulated with a thin film of vapor-deposited Parylene-C
- Four different tip profiles also available (Standard, Heat-Treated, Blunt, and Fine tips)
- High corrosion resistance offers consistent long-term performance.

Benefits

- Many standard types available (web) and custom
- Connecting pin fits Amphenol series 220-223 connectors.

Applications

- Type C: Excellent for bipolar stimulation.
- For acute and chronic recording.

EXPOSED TIP DIMENSIONS (nominal)								
Nominal Impedance	Tungsten	Elgiloy	Platinum Iridium	Pure Iridium				
10 kΩ	250 µm			_				
50 kΩ	200 µm			_				
0.1 ΜΩ	100 µm	120 µm	60 µm	45 µm				
0.5 ΜΩ	55 µm	66 µm	18 µm	14 µm				
1.0 MΩ	30 µm	36 µm	10 μm	10 µm				
2.0 ΜΩ	12 µm	15 µm	6 µm	5 µm				
5.0 MΩ	5 μm	6 µm	3 µm	2.5 µm				

Heat Treated Tip is ideal for penetrating tough membranes (not recommended for chronic




implantation). This process is performed using a microforge in which the heating element is positioned in close proximity to the tip in order to melt the Parylene-C distal to the exposed metal. It provides a smooth transition and produces better adherence of the Parylene-C to the metal.

Kapton* tubing, indicated by "KT" in the part number, extends from the connector to within 5 mm of the tip, providing stiffness and additional insulation to the electrode shaft. Kapton-clad electrodes are recommended when the electrode is to be inserted through a cannula for extra deep penetration.

* Parylene is a trade mark of Union Carbide. Kapton is a trade mark of DuPont. Elgiloy is a trade mark of Elgiloy Ltd.

Monopolar



Bipolar

2" for PT/Ir - 3" or 5" for Tungsten or Stainless Steel

| S mm Typical Gold-Plated Connector Pin Taper Length

Shrink Tubing


Polyimide Tubing

NOTE: Electrode diagrams are not shown to scale.

Concentric Bipolar Electrodes

Excellent for shielded macro recording as well as evoked potentials — especially well suited for bipolar stimulation


The tungsten, platinum-iridium and elgiloy electrode is sharpened to a point and is 75 µm in diameter. The outer stainless steel conductor is insulated with Polyimide tubing to within 0.2 mm of the end of the stainless steel tube. Also available without the outer Polyimide insulation.



Insulated metal conductor with exposed concentric surface

Concentric Bipolar, Rhodes Style Tips


This electrode is designed to replicate the Peter Rhodes-produced semi-micro concentric, and features the same length, thickness and tip dimensions. These probes feature a stepped tip with well-defined element lengths, allowing surface area to be easily determined for charge density calculations. Featuring a shaft diameter considerably thinner than the standard biplar electrodes, this electrode is ideal for localized stimulation or population recording in fixed preparations and acute animal experiments where implant size is a critical concern.



Concentric Bipolar with Extension, Rhodes Style Tips

These Concentric Bipolar Electrodes are designed for Neurological recording and stimulation and are ideal for acute applications (part numbers without the X) or chronic applications (part number ends with X).

The Concentric Bipolar Electrodes, designed for neurological recording and stimulation and are ideal for acute applications. The Stainless Steel external tubing (E) provides support for accurate placement and manipulation during acute research. The extra extension may be cut off after insertion leaving only a small protrusion.



CONCENTRIC ELECTRODES* ORDERING INFORMATION									
Metal Probe Outer Tip X dim. w/ (pkg. of									
ltem	Core	Length	lmp	Diameter (total)	Diam.	Core diam.	Y dim.	polyimide	5)
TM33CCNON	Tungsten	3" (76 mm)	10-15 ΚΩ	0.013" uninsulated (325 µm)	3-4 μm	.003" (76 μm)	0.4 mm	.005" (127 μm)	
TM33CCINS	Tungsten	3" (76 mm)	10-15 ΚΩ	0.016" insulated (400 µm)	3-4 μm	003" (76 μm)	0.4 mm	.005" (127 μm)	
TM53CCINS	Tungsten	5" (127 mm)	10-15 ΚΩ	0.018" insulated (450 µm)	3-4 μm	005" (127 μm)	0.4 mm	.008" (203 μm)	
PTM3CC02INS	Pt/Ir NS fine	3" (76 mm)	200 ΚΩ	0.013" insulated (325 µm)	2-4 μm	0.002" (50.8 μm)	.25 mm	.004" (114 µm)	

^{*}All have a stainless steel outer shaft

METAL ELECTRODES ORDERING INFORMATION								
				Nominal				
		Insul.		Impedance				
Item	Length	Thick	Shaft Diam	. (± 20%)	Diam.	Typical Use		
Tungsten — Prof						Package of 10		
TM31A10	76 mm	1 µm	0.127 mm	1.0 ΜΩ	1 µm	Multi unit and single unit recording		
TM31A20	76 mm	1 µm	0.127 mm	2.0 ΜΩ	1 µm	Multi unit and single unit recording		
TM33A05	76 mm	3 µm	0.127 mm	0.5 ΜΩ	1 µm	Multi unit and single unit recording		
TM33A10	76 mm	3 µm	0.127 mm	1.0 ΜΩ	1 µm	Multi unit and single unit recording		
TM33A20	76 mm	3 µm	0.127 mm	2.0 ΜΩ	1 µm	Multi unit and single unit recording		
TM33B01	76 mm	3 µm	0.254 mm	0.1 ΜΩ	1-2 µm	Single and multi unit recording		
TM33B05	76 mm	3 µm	0.254 mm	0.5 ΜΩ	1-2 µm	Single and multi unit recording		
TM33B10	76 mm	3 µm	0.254 mm	1.0 MΩ	1-2 µm	Single and multi unit recording		
TM33B20	76 mm	3 µm	0.254 mm	2.0 ΜΩ	1-2 µm	Single and multi unit recording		
TM33C10	76 mm	1 µm	0.085 mm	1.0 MΩ	1 µm	Single and multi unit recording		
Tungsten — Prof	ile C					Package of 10		
TM33A10KT	76 mm	3 µm	0.216 mm	1.0 MΩ	1 µm	Multi unit and single unit recording		
TM33B01KT	76 mm	3 µm	0.356 mm	0.1 MΩ	1-2 µm	Single and multi unit recording		
Elgiloy®/Stainless	: — Profile A	\				Package of 10		
SSM33A70	76 mm	3 µm	0.229 mm	7.0 MΩ	1-2 µm	Recording and Stimulating (Prussian blue staining)		
SSM33A20KT	76 mm	3 µm	0.356 mm	2.0 ΜΩ	1-2 µm	Recording and Stimulating (Prussian blue staining)		
Tungsten — Prof	ile B					Package of 10		
TST33A001KT	76 mm	3 µm	0.356 mm	10 kΩ	1 µm	Tissue slice stimulation		
TST33A05KT	76 mm	3 µm	0.356 mm	0.5 ΜΩ	1 µm	Stereotrode / Bipolar, differential measurements		
TST33A10KT	76 mm	3 µm	0.356 mm	1.0 ΜΩ	1 µm	Stereotrode / Bipolar, differential measurements		
TST33A20KT	76 mm	3 µm	0.356 mm	2.0 ΜΩ	1 µm	Stereotrode / Bipolar, differential measurements		
TST33C05KT	76 mm	3 µm	0.216 mm	0.5 ΜΩ	1 µm	Stereotrode / Bipolar, diff. meas. — extra fine (75 µm separation)		
TST53A10KT	127 mm	3 µm	0.356 mm	1.0 MΩ	1-2 µm	Stereotrode / Bipolar, differential measurements		
Pure Iridium — P	rofile A					Package of 10		
IRM23E10	50 mm	3 µm	0.106 mm	1.0 ΜΩ	1-2 µm	Single and multiunit recording and stimulation		
IRM23E15	50 mm	3 µm	0.106 mm	1.5 ΜΩ	1-2 µm	Single and multiunit recording and stimulation		
IRM23E25	50 mm	3 µm	0.106 mm	2.5 ΜΩ		Greater selectivity - small cells		
IRM23E30	50 mm	3 µm	0.106 mm	3.0 MΩ	1-2 µm	Greater selectivity - small cells		
Pure Iridium — P	rofile B					Package of 10		
IRM23E20KT	50 mm	3 µm	0.180 mm	2.0 ΜΩ	1-2 µm	Greater selectivity & microstimulation		
IRM23E30KT	50 mm	3 µm	0.180 mm	3.0 MΩ	1-2 µm	Greater selectivity - small cells		

Elgiloy Steel *Cobalt/chromium/nickel alloy. The KT suffix refers to Kapton™ cladding.

All Metal Microelectrodes are available in custom lengths, blunt or heat treated (extra charge).

Assorted kits may be customized.

Ordering:

Add the **B** suffix where blunt electrodes are desired. (For example, an IRM123A10KT ordered as a blunt will be IRM123A10KTB.)

Add the ${\bf H}$ suffix where heat treated electrodes are desired. (For example, an IRM123A10KT ordered with heat treatment will be IRM123A10KTH.)

ORDERING INFORMATION

OPTIONAL ACCESSORIES/REPLACEMENT PARTS


300102	Micromanipulator holder, 4 in., 2mm to 0.031 socket
5468	2 mm receptacle to 0.031-inch jack (for Omega-TipZ)
5469	Adapts mini banana plug (DAM80) to 0.031-inch receptacle (metal microelectrode)
5470	0.031-inch jack, 28 ga. wire, 12 inch (pkg. of 4)
5482*	Pins, 0.031-inch, gold-plated (pkg. of 50)
5483*	Sockets, 0.031-inch gold-plated (pkg. of 50)

*Gold-plated pins (#5482) and sockets (#5483) may be attached to 24-, 26-, or 28-gauge wire.



Gold-plated pins (5482) and sockets (5483) may be attached to 24-, 26-, or 28-gauge wire.





Metal Electrodes, Rhodes Style Tips with Extension

Concentric bipolar electrodes for neurological recording and stimulation

Features

- Rhodes-style tips are great for those that require a precisely controlled electrode surface area
- The stainless steel external tubing provides support for accurate placement and manipulation during acute research
- These are an excellent replacement for long-term Rhodes electrodes (originally made by Peter Rhodes)

Benefits

- Concentric bipolar electrodes are ideal for bipolar stimulation paradigms
- Tools are available for shielded macro recordings, as well as evoked potential
- Electrodes may be acutely reused many times (> 20-30 insertions), if they are properly cared for:
 - Avoid bending the shaft and blunting/hooking the tip
 - Clean electrode carefully between session

Core Electrode Polyimide Tubing Stainless Steel Tube Stainless Steel Tube Pin Connector E - Extension S = 10 mm Mount Insulation

Applications

 Designed for neurological recording and stimulation and are ideal for chronic applications

Connection to the

Core Electrode

D>0.79mm

The Concentric Bipolar Electrodes, designed for Neurological recording and stimulation, are ideal for acute applications (without the X) or chronic applications (ends with X). The stainless steel external tubing (E) provides support for accurate placement and manipulation during acute research. The extra extension may be cut off after insertion leaving only a small protrusion.

Concentric Bipolar Electrodes with Extension



Semi-Micro Concentric Bipolar Electrodes with Extension




ORDERING INFORMATION									
Item	Туре	Metal Core	Length (s)	Probe Total OD (z)	OD Inner sleeve (x)	Extension (E)			
SSM04RC100	Semi micro concentric with extension	Stainless Steel	0.4" (10 mm)	Insulated (411 µm)	140 μm	50 mm			
PTM04RC100	Semi micro concentric with extension	Platinum Iridium	0.4" (10 mm)	Insulated (411 µm)	140 µm	50 mm			
SSM041RC225	Concentric bipolar with extension	Stainless steel	0.4" (10 mm)	insulated (625 µm)	305 μm	50 mm			
PTM041RC225	Concentric bipolar with extension	Platinum iridium	0.4" (10 mm)	insulated (625 µm)	305 μm	50 mm			

^{*}All have a stainless steel outer shaft. Package of 5.

Metal Electrodes, Rhodes Style Tips, No Extension

Concentric bipolar electrodes for neurological recording and stimulation



Concentric Bipolar Electrodes

Choose stainless steel or platinum iridium.



ORDERING INFORMATION								
Item	Metal Core	Length (s)	Probe Total OD (z)	OD of Inner Sleeve (x)	Outer Polyimide	Extension (E)		
SSM321RC225	Stainless steel	3.125" (80 mm)	insulated (625 µm)	225 μm	Yes	30 mm		
PTM321RC225	Platinum iridium	3.125" (80 mm)	insulated (625 µm)	225 μm	Yes	30 mm		

Semi-Micro Concentric Bipolar Electrodes

Choose stainless steel or platinum iridium.



ORDERING INFORMATION								
Item	Metal Core	Length (s)	Probe Total OD (z)	OD of inner sleeve (x)	Outer Polyimide	Extension (E)		
SSM321RC100	Semi micro concentric	Stainless steel	3.125" (80 mm)	insulated (411 μm)	140 µm	30 mm		
PTM321RC100	Semi micro concentric	Platinum iridium	3.125" (80 mm)	insulated (411 µm)	140 µm	30 mm		

^{*}All have a stainless steel outer shaft.

Nerve Cuff Electrodes

For acute/chronic experiments

Features

- For both acute and chronic experiments
- 1-24 electrode sites
- Broad range of available inner diameters, from 5 mm down to as small as 56 µm
- Platinum, stainless steel and platinum/iridium metal electrodes

Benefits

- Customizable electrode arrangements, including concentric and tri-bipolar
- Adaptable to interface with the acquisition system of your choice

Custom length 1 to 24 contacts Custom diameter

Applications

• Suitable for rodent, feline, bird and primate research

Nerve cuff electrodes are designed for reliable recording and/or stimulation of any peripheral nerve. They can be used acutely or for chronic implantation, with a wide selection of inner diameters available to be selected based on the nerve diameter. The nerve cuffs are designed to provide flexibility in electrode contact location for recording and stimulation protocols. Three different electrode variations are available:

- Standard, Micro Cuff and Nanocuff
- Concentric Nerve Cuff
- X-wide Nerve Cuffs


1. Standard, Micro Cuff and Nano Cuff

Many custom combinations of inner diameters and contact arrangements are available. Standard nerve cuffs have a variety of configurable parameters.



AVAILABLE OPTIONS					
Configuration	Metal Type	ID			
Standard	125µm stranded Stainless Steel	1.0 – 5.0 mm			
	100µm Platinum	1.0 - 2.0 mm			
	250µm Platinum*	2.0 - 5.0 mm			
Micro Cuff	100µm Platinum	0.5 - 0.75 mm			
	50µm Platinum/Iridium	0.3 mm			
Nano Cuff	25µm Platinum/Iridium	160 - 250 µm			
	12 5um Platinum/Iridium	56 140 um			

* Recommended for stimulation




2. Concentric Nerve Cuff

Concentric electrodes have multiple contacts around a single point of the nerve, allowing recording or stimulation at different locations around the same point.




Configuration Options for Standard, Micro Cuff, Nano Cuff & Concentric Nerve Cuffs

Many custom combinations of inner diameters and contact arrangements are available for standard, micro and nano nerve cuffs.



- Optional Sutures
- 2 1 to 24 Contacts (Platinum, Platinum/Iridium or Stainless Steel)
- 3 Custom diameter
- 4 Custom length
- 5 Stranded stainless steel leads. The angle is selectable.

3. X-Wide Contact Nerve Cuffs



- Inner diameter
- Spacing between contacts
- Number of contacts n
- Distance from the last contacts to the edge of the cuff У
- Stranded stainless steel leads with Teflon insulation 300 mm long (typical)
- Width of the contact, Pt/Ir 90/10
- Extra large platinum ribbon contacts maximize electrode surface area while maintaining the reliability and customizability
- Ideal for any electrical stimulation studies, especially those using KHz stimulation to block activity of peripheral nerves
- Lower thresholds and increased safety margins provide you with new freedom to develop even more effective stimulation paradigms
- The same reliability enhancement features as our standard cuffs, including reinforced bonds and specially-designed helical leads, providing exceptional implanted longevity
- Fully-customizable specifications, including contact spacing, cuff length and lead features, with a wide range of available inner diameters from over 6 mm down to 300 µm.
- · Configurable for both acute and chronic use
- Connector options for any available stimulation or recording system

The key utility for these cuffs is provided by the very large contact surface area, which is ideal for studies with high current levels and when using kilohertz electrical stimulation to achieve block of activity in peripheral nerves. The reduced impedances ensure that the output of stimulation sources is not attenuated, that no destructive charge builds up, and that true block of nerve activity is achieved.

Types of Nerve Cuff Electrodes

Application	Electrode Type	d (mm)	x (mm)	y (mm)	W (mm)
Sciatic	Rat, 1.4 mm avg. nerve OD, 2.5 mm cuff ID	2.5	1	1.5	0.5
Tibial	Rat, 1.2 mm avg. nerve OD, 1.5 mm cuff ID	1.5	1	1.5	0.5
Cervical Vagus	Rat, 0.5 mm avg. nerve OD, 0.75 mm cuff ID	0.75	1	1.5	0.5
	Rat, 0.3 mm avg. nerve OD, 0.5 mm cuff ID	0.5	1	1.5	0.5
Custom Standard		1.0-	5.0		
Custom Microcuff		0.5 or	0.75		

Sterilization

Nerve cuffs are made entirely of autoclavable materials — silicone rubber, Teflon and stainless steel. They can be steam autoclaved without special precautions. If gas sterilization (EtO) is preferred, pack the nerve cuffs in a gas-permeable bag and allow adequate outgassing time (at least 48 hrs.) to be sure all toxic gases have been desorbed from the silicone rubber.

> Order online at www.wpiinc.com/nce

Nerve Cuff Connectors

Omnetic connectors may be ordered separately. Other connectors must be ordered with nerve cuff electrodes. Connector options include:







- · No connector
- Invivo1 (not sold separately)
- Banana Terminal (not sold separately)
- Touch Proof Connector (not sold separately)
- Omnetic Connectors

Male Implant Connector	Description
\$F96	MP-A11365-001
Ti di	Male Omnetics connector - 8 channel, 10-pin 0.050" pitch
ib.	MP-A12623-001
	Male Omnetics connector- 8 channel, 10-pin, 0.050" pitch with latching mechanism
	MP-48393-001



MP-A8393-001

Male Omnetics connector- 8 channel, 10-pin, 0.025" pitch, 2 guide post



MP-A79000-001 (NPS-09-DD-GS)

Male Omnetics connector- 7 channel, 9-pin, 0.050"



MP-A79014-001 (NPD-18-DD-GS)

Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 6 guides posts



MP-A70242-001

Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 6 guides posts; Nickel free, MRI compatible



MP-A79038-001 (NPD-18-DDGS)

Male Omnetics connector- 16 channel, 18-pin, 0.025" pitch, 2 guide posts



MP-A79022-001 (NPD-36-DD-GS)

Male Omnetics connector -32 channel, 36-pin, 0.025" pitch, 4 guide posts



MP-A72312-001

Male Omnetics connector -32 channel, 36-pin, 0.025" pitch, 4 guide posts; Nickel free, MRI compatible

References

Guillem, K., & Ahmed, S. H. (2017). Preference for Cocaine is Represented in the Orbitofrontal Cortex by an Increased Proportion of Cocaine Use-Coding Neurons. Cerebral Cortex, 28(3), 819-832. http://doi.org/10.1093/cercor/

Uy, J., Yu, M., Jiang, X., Jones, C., Shen, B., Wang, J., ... Tai, C. (2017). Glutamatergic Mechanisms Involved in Bladder Overactivity and Pudendal Neuromodulation in Cats. The Journal of Pharmacology and Experimental Therapeutics, 362(1), 53-58. http://doi.org/10.1124/jpet.117.240895

Khurram, A., Ross, S. E., Sperry, Z. J., Ouyang, A., Stephan, C., Jiman, A. A., & Bruns, T. M. (2017). Chronic monitoring of lower urinary tract activity via a sacral dorsal root ganglia interface. Journal of Neural Engineering, 14(3), 036027. http://doi.org/10.1088/1741-2552/aa6801

Restaino, S. M., Abliz, E., Wachrathit, K., Krauthamer, V., & Shah, S. B. (2014). Biomechanical and functional variation in rat sciatic nerve following cuff electrode implantation. Journal of NeuroEngineering and Rehabilitation, 11(1), 73. http://doi.org/10.1186/1743-0003-11-73

Lissandrello, C. A., Gillis, W. F., Shen, J., Pearre, B. W., Vitale, F., Pasquali, M., ... Gardner, T. J. (2017). A micro-scale printable nanoclip for electrical stimulation and recording in small nerves. Journal of Neural Engineering, 14(3), 036006. http://doi.org/10.1088/1741-2552/aa5a5b

Cables & Connectors

For wiring any laboratory setup



	ORDERING INFORMATION							
Part #	Application/Description	Connector A	Connector B	Cable Length				
2026-10	2 mm Socket, unwired (pkg. of 10) (Not Shown)	2 mm socket	unwired	none				
2851	Standard BNC Cable	BNC (male)	BNC (male)	6 ft. (1.8 m)				
3142	Mini-Banana Adapter	Screw Terminals	Dual Mini-Banana	none				
3161	Connector for input to TBM4M and BP-1	DIN (male)	unwired	none				
3294	Ground Wire for DAM80 probe	Clip	none	3 ft. (0.9 m)				
3417-10	2 mm Plug, unwired (pkg. of 10)	2 mm pin	unwired	none				
3491	Extension for any 8-pin DIN	DIN (male)	DIN (female)	5 ft. (1.5 m)				
3492	Connector, adapts WPI transducers to non-WPI equipment	DIN (female)	unwired	none				
3517	DAM50, DAM60, DAM70, shielded (two cables/pkg)	Modular phone plug, 4 wire	none	3 ft. (0.9 m)				
3578	Adapter cable for Ag/AgCl pellets	2 mm pin	none	5 ft. (1.5 m)				
3670	Double banana plug with solder turret terminals	Dual Banana (male)	Dual Banana (female)	none				
5371	Low-noise cable for microelectrode holders	2 mm gold pin/jack	2 mm gold pin/jack	2 ft. (0.6 m)				
5373	Low-noise cable for microelectrode holders	2 mm gold pin/jack	2 mm gold pin/jack	2 ft. (0.6 m)				
5374	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold pin	4 ft. (1.2 m)				
5375	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold jack	4 ft. (1.2 m)				
5385	Cable, shielded transducer stock	none	none	25 ft. (7.6 m)				
13324	Adapter	Double-banana (female)	BNC (male)	none				
13347	ISO2 (chart recorder adapter)	Double-banana (male)	BNC (female)	none				
13388	Electrode adapter for DAM probes	Miniature banana (male)	2 mm jack	none				
13451	Adapter: Iso-DAM, Iso-DAM8	BNC (female)	two 2 mm pins	6 in. (15 cm)				
13620	Low-noise cable for microelectrode holders	2 mm gold pin	2 mm gold jack	2 ft. (0.6 m)				
13685	SP Series pump-to-pump linking cable	Modular phone plug	Modular phone plug	7 ft. (2.1 m)				
13776	Adapts reference electrode to VF4 ground jack	Banana (male)	2 mm jack	none				
13854	BNC T-connector, male to:	BNC (female)	BNC (female)	none				
14254	BNC Straight Adapter	BNC (female)	BNC (female)	none				
15623	Serial cable and adapter, SP Series pump	SP Pump	IBM 9-pin "D" connector	5 ft. (1.5 m)				
15975	Adapter	2 mm socket	1 mm pin	none				
15976	Adapter	1 mm socket	2 mm pin	none				
300040	Adapter Extension	2 mm socket	2 mm socket	4 in. (10 cm)				
500184	Standard BNC Cable	BNC (male)	BNC (male)	10 ft. (3 m)				
500256	BNC Right Angle Adapter	BNC (male)	BNC (female)	none				
500257	Standard BNC Cable	BNC (male)	BNC (male)	6 in. (15 cm)				
	Standard BNC Cable	BNC (male)	BNC (male)	12 in. (30 cm)				
	Standard BNC Cable	BNC (male)	BNC (male)	18 in. (46 cm)				
	Cable, Extension	8-pin miniDIN (male)	8-pin miniDIN (female)	10 ft. (3 m)				
	Cable, USB	USB (male)	USB (female)	10 ft. (3 m)				
	Cable (red and black pair)	Banana (male)	Banana (male)	36 in. (91 cm)				
	Cable (red and black pair)	Banana (male)	Mini-Gator	36 in. (91 cm)				
504715		Banana (male)	Mini-Clip	36 in. (91 cm)				
	Cable (red and black pair)	Banana (male)	Micro-Clip	36 in. (91 cm)				
	MiniPhone Patch Cable	3.5 mm MiniPhone plug	3.5 mm MiniPhone plug	6 ft. (1.8 m)				
CBL102	DAM Series, PM Series	3.5 mm MiniPhone plug	BNC (male)	6 ft. (1.8 m)				



Micropipette Holders & Half-Cells

Couple fluid filled glass micropipettes to high input impedance amplifiers

Features


- Connector pin, jack, luer ports, Ag/AgCl half-cells, wire selections available
- Optional handles (for some varieties) allow for easy manipulation. Handles are sold separately.
- Screw cap for tight fitting of glass electrode to avoid drifting of electrode

Benefits

- Available for 1.0, 1.2, 1.5 and 2.0 mm OD glass electrodes
- Spare gaskets available
- Optional handles available in two different sizes for some holders

Applications

- Microinjection
- · Electrophysiology recording
- Fluid handling



The most popular micropipette holders are stock items. Custom holders (designated by ◆) can be manufactured on demand but require an additional setup fee. Call for a quote. See all the options at www. wpiinc.com/MEH.


WPI's microelectrode holder-half-cells couple fluid-filled glass micropipettes to high input impedance amplifiers. A Ag/AgCl pellet (or a silver wire) molded into the holder body provides stable potential. Electrical connection is made via male 2 mm pins or female 2 mm sockets. The pipette may be mounted axially or at right angles to the holder. Pipettes are held with screwcaps or rubber gaskets (without caps). Filling WPI microelectrode holders with electrolytes containing chloride results in stable electrode potential. Suitable electrolytes include KCl, NaCl and CaCl₂. Holders are supplied for standard WPI single capillary tubing of 1.0, 1.2, 1.5 and 2.0 mm outside diameters. (Call WPI regarding custom designs for other glass diameters.) The holder style you select will depend on your experimental application, space, and instrumentation.

Hints for selecting and ordering micropipette holders


- Determine the required electrical connection on the holder: for example, if you wish to connect the holder to a 2 mm pin you should select a holder equipped with a 2 mm jack. Most WPI probes require a holder equipped with a 2 mm jack.
- 2. Decide on the required alignment of the electrical connection: either in-line with the glass pipette, or at a right angle to it. Space considerations in your experimental setup and requirements imposed by other pieces of equipment typically determine which alignment is appropriate.
- 3. Determine if you want to hold the glass pipette by a rubber gasket (e.g., MEH1S) or a screwcap (e.g., MEH3S). Rubber gaskets offer easier insertion and removal of glass pipettes whereas screw-caps provide more secure mounts for micropipetters.


- 4. Choose a holder with either a silver wire or a silver/silver chloride pellet for the metal/liquid coupling. Silver/silver chloride pellets provide a more stable low-noise baseline which is important for low-noise DC recording. Pellets require the glass pipette and holder to be free of air bubbles to achieve a good connection. Silver wire holders are durable and are easier to use when the holder is equipped with a pressure port because the fluid in the pipette does not have to be filled to the top of the pipette to achieve a good electrical connection
- 5. Choose a holder equipped with a pressure port only when you want to pressure inject liquid from the pipette. Two types of ports are available: 2.0 mm O.D. and standard "syringe-style" Luer. The Luer port is often recommended because it makes assembly and disassembly much easier. Quick-connect Luer fittings for four common sizes of tubing (1/16", 3/12", 1/8", 5/12" I.D.) are included with each Luerequipped holder.
- 6. Some non-WPI preamplifiers or headstages cannot be mounted on micromanipulators. In such cases, a holder equipped with a rod (e.g., MEH8) permits the holder to be conveniently mounted on a micromanipulator.
- Finally, remember to specify the O.D. of the glass you will be using when you place your order.


MEH6RF/SF is designed primarily for use with the Model 900A Micropressure System; **EHB1** for use in electrode beveling; and **MEH3SW** for microtitration of chloride with a silver wire as the electrode and a solution of silver nitrate filling the holder. **MPH** models do not contain Ag/AgCl half-cells and are used for pressure injection of substances through microelectrodes. **PicoNozzle**, used for pressure injection with PV800 Series PicoPumps, includes an **MPH6S** holder — which may also be used to couple a micropipette to a syringe.















♦ denotes holder manufactured for you as custom orders. Call for price.













































Micropipette Holders & Half-Cells

Couple fluid filled glass micropipettes to high input impedance amplifiers










































• denotes holder sizes manufactured for you as custom orders. Call for price.

	MI	CROELECTRODE	HOLDERS	ORDERI	NG INFOR	MATION	I
♦ = custom order (call for price)							
Order Number Replace XX with glass OD * EHBF	Order Order	Electric Connection Angle Straight	Connector Female	Half-Cell Wire	Pressure Port No Port	Screw Cap N/A	Designed for WPI Products MBS, 48000
FOIMPH		Straight	Fiber Optic	None	No Port	w/Cap	MBS, 48000
MEH1Sxx		Straight	Male	Pellet	No Port	No Cap	
MEH1SFxx		Straight	Female	Pellet	No Port	No Cap	705, 773, 767, 721, FD223
MEH2Rxx	•	Right	Male	Pellet	Male Luer	w/Cap	
ЛЕH2RFxx	•	Right	Female	Pellet	Male Luer	w/Cap	705, 773, 767, 721, FD223
MEH2RFWxx	•	Right	Female	Wire	Male Luer	w/Cap	705, 773, 767, 721, FD223
MEH2RWxx	•	Right	Male	Wire	Male Luer	w/Cap	
ЛЕН2Sxx		Straight	Male	Pellet	Male Luer	w/Cap	
MEH2SFxx		Straight	Female	Pellet	Male Luer	w/Cap	705, 773, 767, 721, FD223
/IEH2SFWxx	•	Straight	Female	Wire	Male Luer	w/Cap	705, 773, 767, 721, FD223
MEH2SWxx		Straight	Male	Wire	Male Luer	w/Cap	
ЛЕН345хх		45°	Male	Pellet	No Port	w/Cap	
/IEH3F45xx	•	45°	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH3FW45xx	•	45°	Female	Wire	Port	w/Cap	
/IEH3Rxx		Right	Male	Pellet	No Port	w/Cap	
/IEH3RFxx		Right	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH3RFWxx	•	Right	Female	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH3RWxx	•	Right	Male	Wire	No Port	w/Cap	
ЛЕН3Sxx		Straight	Male	Pellet	No Port	w/Cap	
//EH3SBWxx	•	Straight	Banana	Wire	No Port	w/Cap	ISO-80, ISO-DAM8A
/IEH3SFxx		Straight	Female	Pellet	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH3SFWxx		Straight	Female	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH3SWxx		Straight	Male	Wire	No Port	w/Cap	
/IEH3W45xx	•	45°	Male	Wire	No Port	w/Cap	705, 773, 767, 721, FD223
/IEH6RFxx		Right	Female	Pellet	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
/IEH6RFWxx	•	Right	Female	Wire	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
/IEH6SFxx		Straight	Female	Pellet	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
/IEH6SFWxx	•	Straight	Female	Wire	2.0-mm Port	w/Cap	705, 773, 767, 721, FD223
ЛЕН7хх		Right	Male	Pellet	2.0-mm Port	w/Cap	
/IEH7Wxx	•	Right	Male	Wire	2.0-mm Port	w/Cap	
ЛЕН8хх		Right	Male	Pellet	No Port	w/Cap	
/IEH900Rxx		Right	Male	Pellet	2.0-mm Port	w/Cap	900A
ЛЕН900Sxx		Straight	Male	Pellet	2.0-mm Port	w/Cap	900A
ЛРН1хх		_	None	None	Female Luer	w/Cap	
ЛРН3хх		_	None	None	Male Luer	w/Cap	
MPH4xx		_	None	None	2.0-mm Port	w/Cap	
MPH6Pxx	•	Right	Male	Pellet	Female Luer	w/Cap	Piconozzle Kit (5430-XX)
MPH6Rxx	•	Right	Male	Wire	Female Luer	w/Cap	Piconozzle Kit (5430-XX)
MPH6Sxx		_	None	None	Female Luer	w/Cap	Piconozzle Kit (5430-XX)

 $[\]star$ Specify O.D. of glass (1.0, 1.2, 1.5 or 2.0 mm) by replacing XX in the Order Number with 10, 12, 15 or 20.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

	Handle 2505 is for use with WPI manipulators. The smaller diameter handle 5444 is required for use with Narishige and Zeiss manipulators.
2505	¹ /4-in (6.3 mm) diameter handle
5444	³/16-in (4.8 mm) diameter handle
GO1-10	Replacement gasket 1.0 mm, Package of 100
GO2-10	Replacement gasket 1.2 mm, Package of 100
GO3-10	Replacement gasket 1.5 mm, Package of 100
GO4-10	Replacement gasket 2.0 mm, Package of 100
1571	Clear Silicone Rubber Sealant (-4.7 oz-)

Precious Metal and Specialty Wire

Bare and coated metal wire for most laboratory applications



Micro coaxial cables (MAXxxxx) are ideal for microelectrode fabrication and construction of similar research tools. The dual shielding eliminates electrical interference caused by radio frequencies (RF), electrostatic and microphonics (e.g., bending and vibration. Available with single or dual (twin) conductors

Teflon-coated stainless steel (type 304) wire (SSTxxxx) is available in 25-ft and 50-ft lengths. The Teflon coating is 150 micro-in. thick (4 µm). The Teflon coating is designed to reduce surface friction, only. It is not insulation.

Carbon wire (C3005) is a single 30-micron fiber of electrochemically activated carbon. This fiber is especially useful in micro-electrochemical experiments.

Platinum/iridium wire — uncoated (PTxxxx) and Teflon-coated (PTTxxxx) is an alloy of 90% platinum and 10% iridium, giving excellent tensile strength and corrosion resistance. Uncoated pure platinum wire (PTPxxx) is 99.95% pure. Indium wire (IN1003) is 99.99% pure, with

Annealed silver wire (AGWxxxx), 99.99% pure, is available in five diameters; three of those sizes are also available with a Teflon coating (AGTxxxx).

a melting point of 156.4°C.

Tungsten wire (TGWxxxx), available in three diameters, is 99.95% pure.

Gold wire (AUWxxxx) is 99.99% pure. Stainless steel wire (SSxxxxx) is type 316.

Wire Cutters and Scissors



ORDERING INFORMATION							
Catalog No.	Metal	Coating	AWG*	Diameter	Precut Length		
AGT0510	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	10 ft (3 m)		
AGT0525	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	25 ft (7.6 m)		
AGT05100	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	100 ft (30 m)		
AGT1010	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	10 ft (3 m)		
AGT1025	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	25 ft (7.6 m)		
AGT10100	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	100 ft (30 m)		
AGT1510	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	10 ft (3 m)		
AGT1530	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	30 ft (9.1 m)		
AGW0510	Silver	_	36	0.005 in. (0.125 mm)	10 ft (3 m)		
AGW0530	Silver	_	36	0.005 in. (0.125 mm)	30 ft (9.1 m)		
AGW1010	Silver	_	30	0.010 in. (0.25 mm)	10 ft (3 m)		
AGW1030	Silver	_	30	0.010 in. (0.25 mm)	30 ft (9.1 m)		
AGW1510	Silver	_	26-27	0.015 in. (0.38 mm)	10 ft (3 m)		
AGW1530	Silver	_	26-27	0.015 in. (0.38 mm)	30 ft (9.1 m)		
AGW2010	Silver	_	24	0.020 in. (0.5 mm)	10 ft (3 m)		
AGW2030	Silver	_	24	0.020 in. (0.5 mm)	30 ft (9.1 m)		
AGW4010	Silver	_	18	0.040 in. (1.0 mm)	10 ft (3 m)		
AUW0170	Gold	_	50	0.001 in. (0.025 mm)	70 ft (21 m)		
AUW201	Gold	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)		
C3005	Carbon	_	49	0.0012 in. (30 µm)	5 ft (1.5 m)		
PT1002	Platinum / Iridium	_	30	0.010 in. (0.25 mm)	2 ft (61 cm)		
PT0402	Platinum / Iridium	_	38	0.004 in. (0.102 mm)	2 ft (61 cm)		
PT0203	Platinum / Iridium	_	44	0.002 in. (0.051 mm)	3 ft (91 cm)		
PT0110	Platinum / Iridium	_	50	0.001 in. (0.025 mm)	10 ft (3 m)		
PTP101	Platinum	_	30	0.010 in. (0.25 mm)	1 ft (30 cm)		
PTP201	Platinum	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)		
PTP401	Platinum	_	18	0.039 in. (1.0 mm)	1 ft (30 cm)		
PTP406	Platinum	_	18	0.039 in. (1.0 mm)	0.5 ft (15.2 cm)		
PTT0502	Platinum / Iridium	Teflon	36	0.005 in. (0.125 mm) ¹	2 ft (61 cm)		
PTT0203	Platinum / Iridium	Teflon	44	0.002 in. (0.051 mm) ¹	3 ft (91 cm)		
PTT0110	Platinum / Iridium	Teflon	50	0.001 in. (0.025 mm) ¹	10 ft (3 m)		
SS31605	Stainless Steel	_	36	0.005 in. (0.125 mm)	50 ft (15.2 m)		
SS31614	Stainless Steel	_	27	0.014 in. (0.36 mm)	30 ft (9.1 m)		
SST30407-25	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	25 ft (7.6 m)		
SST30407-50	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	50 ft (15.2 m)		
TGW0325	Tungsten	_	40	0.003 in. (0.075 mm)	25 ft (7.6 m)		
TGW0515	Tungsten	_	36	0.005 in. (0.125 mm)	15 ft (4.6 m)		
TGW1510	Tungsten	_	26-27	0.015 in. (0.38 mm)	10 ft (3 m)		
		/ICROCO	AXIAL				
MAX3820	Tinned Cu Alloy	Coaxial	0.0450	0.0173 in. (0.44 mm)	20 ft (6 m) ⁴		
MAX4020	Tinned Cu Alloy	Twin Coaxial	0.0158x0	0.024 in. (0.4x0.61 mm)	20 ft (6 m) ⁵		

*Brown & Sharpe

¹ Plus 0.002 in. for Teflon coating

ORDERING INFORMATION 504749 Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) Micro-shear flush cutters for delicate wires 504750 Ergonomic Mini-Scissors, 12.7 cm (5 in.) for cutting fine or delicate items with a clean, square edge. Handle design is advantageous for users with arthritic hands. 504751 Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) for delicate wires. ESD safe. Extra tough high carbon steel blades. Sized for smaller hands and maximum maneuverability


³ Teflon adds 0.00015 in. (4 µm) to diameter

⁴ Impedance: 50 ohm; capacitance: 95 pF/m; resistance: 5 ohm/m

⁵ Impedance: 100 ohm; capacitance: 54 pF/m; resistance: 1.9 ohm/m


Ag/AgCl Half-Cells


Sintering pellets with low resistance and high strength





New, improved sintered pellets with lower resistance and high strength. Stable and well balanced in the presence of current, these small and inexpensive half-cells are easy to work with as bath electrodes.


	ORDERING INFORMATION
RC1	Reference Cell with 1.5 m lead
RC1T	Reference Cell, 1.5 m lead, 2 mm pin
RC2	Reference Cell with 2.0 mm pin
RC2F	Reference Cell with female connector
RC3	Ref. Cell with epoxy body, 4.5 mm diam x 50 mm
RC6	Reference Cell with glass body, 1.5 mm diam x 50 mm
EP05	Ag/AgCl Electrode 0.5 mm diam x 20 mm
EP08	Ag/AgCl Electrode 0.8 mm diam x 20 mm
EP1	Ag/AgCl Electrode 1.0 mm diam x 3 mm
EP2	Ag/AgCl Electrode 2.0 mm diam x 4 mm
EP4	Ag/AgCl Electrode 4.0 mm diam x 1 mm
EP8	Ag/AgCl Electrode 8.0 mm diam x 1 mm
EP12	Ag/AgCl Electrode 12.0 mm diam x 1 mm
3578	Adapter Cable for Ag/AgCl Pellets











Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

Features

- · Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

Benefits

- Superior pricing
- Most glass orders ship within 48 hours

Applications

- Microinjection
- Electrophysiology
- Patch clamp
- Fluid Handling

Fire Polishing

Fire-Polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



Making Uniform, Reproducable Microelectrodes


Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

Filament Glass Capillaries


Single Barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.

Thin Wall Glass Capillaries

Thin Wall single barrel capillaries are offered both with or without inner filaments.



Specialty glass is also available. See page 105.



ORDERING INFORMATION							
	Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	Item
	3 in. (76 mm)	1.0	0.58	V		500	1B100F-3
	3 in. (76 mm)	1.0	0.58			500	1B100-3
	3 in. (76 mm)	1.2	0.68	V		350	1B120F-3
SS	3 in. (76 mm)	1.2	0.68			350	1B120-3
a	3 in. (76 mm)	1.5	0.84	~		225	1B150F-3
Single Barrel Standard Borosilicate Glass	3 in. (76 mm)	1.5	0.84		~	300	1B150-3
cat	4 in. (100 mm)	1.0	0.58	V	V	500	1B100F-4
=	4 in. (100 mm)	1.0	0.58		~	500	1B100-4
ĕ	4 in. (100 mm)	1.2	0.68	~	V	400	1B120F-4
Bo	4 in. (100 mm)	1.2	0.68			350	1B120-4
5	4 in. (100 mm)	1.5	0.84	V	V	300	1B150F-4
qa	4 in. (100 mm)	1.5	0.84		~	300	1B150-4
ian Ean	4 in. (100 mm)	2.0	1.12	V		125	1B200F-4
S	4 in. (100 mm)	2.0	1.12		~	200	1B200-4
re	6 in. (152 mm)	1.0	0.58	~		500	1B100F-6
Bal	6 in. (152 mm)	1.0	0.58			500	1B100-6
<u>o</u>	6 in. (152 mm)	1.2	0.68	~		350	1B120F-6
50	6 in. (152 mm)	1.2	0.68			350	1B120-6
S	6 in. (152 mm)	1.5	0.84	V		225	1B150F-6
	6 in. (152 mm)	1.5	0.84			225	1B150-6
	6 in. (152 mm)	2.0	1.12	V		125	1B200F-6
	6 in. (152 mm)	2.0	1.12			125	1B200-6
	3 in. (76 mm)	1.0	0.75	✓		500	TW100F-3
	3 in. (76 mm)	1.0	0.75			500	TW100-3
_	3 in. (76 mm)	1.2	0.90	~	V	400	TW120F-3
arc	3 in. (76 mm)	1.2	0.90			350	TW120-3
b	3 in. (76 mm)	1.5	1.12	~		225	TW150F-3
ita	3 in. (76 mm)	1.5	1.12		<u> </u>	300	TW150-3
<u> </u>	4 in. (100 mm)	1.0	0.75	/		500	TW100F-4
Ę	4 in. (100 mm)	1.0	0.75		~	500	TW100-4
Ä	4 in. (100 mm)	1.2	0.90	✓		350	TW120F-4
p	4 in. (100 mm)	1.2	0.90			350	TW120-4
ï.	4 in. (100 mm)	1.5	1.12	v		225	TW150F-4
Thin-Wall Single-Barrel Standard	4 in. (100 mm)	1.5	1.12		V	300	TW150-4
Ma	6 in. (152 mm)	1.0	0.75	v		500	TW100F-6
Ë	6 in. (152 mm)	1.0	0.75		V	500	TW100-6
Ę	6 in. (152 mm)	1.2	0.90	~	V	400	TW120F-6
	6 in. (152 mm)	1.2	0.90			350	TW120-6
	6 in. (152 mm)	1.5	1.12	~		225	TW150F-6
	6 in. (152 mm)	1.5	1.12		V	300	TW150-6


Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330.

Tolerance ±10%

Multi-Barrel Glass Capillaries

Multi-barrel configurations are designed especially for microiontophoresis. Because the capillaries are fused together during manufacturing, you will not need to twist them while pulling to seal the tips together. An inner filament in each barrel makes filling easy and fast.

Also see PolyFil for a novel way to connect multi-barrel pipettes




ORDERING INFORMATION

Multi-Barrel Borosilicate Glass Tubing with Filaments


Length	Description	OD/ID (per barrel)	Filament	Qty	Item
4" (102mm)	2-Barrel	1.5/0.84 mm	~	100	2B150F-4
4" (102mm)	3-Barrel	1.2/0.68 mm	✓	100	3B120F-4
4" (102mm)	5-Barrel	1.2/0.68 mm	V	65	5B120F-4
4" (102mm)	7-Barrel	1.0/0.58 mm	~	60	7B100F-4
4" (102mm)	7-Barrel	1.2/0.68 mm	V	75	7B120F-4
6" (152mm)	2-Barrel	1.5/0.84 mm	✓	100	2B150F-6
6" (152mm)	3-Barrel	1.2/0.68 mm	~	100	3B120F-6
6" (152mm)	5-Barrel	1.2/0.68 mm	✓	65	5B120F-6
6" (152mm)	7-Barrel	1.0/0.58 mm	V	60	7B100F-6

Special Configuration Glass Capillaries Septum Theta

Septum Theta offers superior cell impalement. The natural bevel resulting from the prominent spear-like projection of the septum gives microelectrodes a sharp, spear-point tip. This style has low resistance for use as a single microelectrode, and it can be used to make superior double-tipped microelectrodes with low transtip coupling. The natural bevel of Septum Theta also significantly increases the effective tip cross-section. As supplied, the width of the septum is approximately 0.2 mm; wall thickness is approximately 0.2 mm.



Piggyback



Piggyback glass consists of a pair of borosilicate capillaries fused together during manufacturing. One barrel is larger than the other, and both have inner filaments for quick filling. Piggyback glass makes it simple to fabricate two-barrel electrodes with a significant tip diameter differential.

ORDERING INFORMATION

Special Configuration Borosilicate Glass Tubing

Description	OD/ID (mm)	Length	Qty	Item
Septum Theta	1.5/1.02	6 in. (152 mm)	100	TST150-6
Piggyback	1.51/0.84 0.75/0.35	4 in. (102 mm)	50	PB150F-4
Piggyback	1.51/0.84 0.75/0.35	6 in. (152 mm)	50	PB150F-6

Borosilicate Glass Rod

1.0 mm diameter — for making tools, probes, tips

	ORDERING INFORMATION						
Borosilicate Glass Rod							
	Description	OD (mm)	Length	Qty	Item		
	Glass Rod	1.0	4 in. (102mm)	500	GR100-4		
	Glass Rod	1.0	6 in. (152mm)	500	GR100-6		

Multi-barrel pipette coupling kit

Securely couple multi-barrel pipettes to a pressure source



Features

- Complete kit (all in one) for making multi-barrel pipettes
- Secure coupling of multi-barrel micropipettes to a pressure source

Benefits

- Allows for microinjection to multi-barrel micropipette independently for one pressure source
- Pressure safe and convenient luer lock connections

Applications

• Multi-port microinjection

PolyFil™ allows easy and secure coupling of a multi-barrel micropipette to a pressure source. Coupling is achieved by bonding temperature-resistant and flexible MicroFil to the capillary tube with hot melt adhesive. The luer end of each MicroFil is connected to PVC tubing (200 PSI rated). Kits also include a five-port manifold that allows use of a single PV800 Series PicoPump to drive up to six micropipette barrels independently by switching on only the barrels to be injected. All connections are locking luers — pressure safe and convenient.

ORDERING INFORMATION

PolyFil™ Multi-Barrel Micropipette Coupling Kit
 Includes: 1 pipette holder/handle, plastic; 7 pieces MF28G MicroFil;
 7-pieces tubing with male luer lock fittings; 1 flow-thru manifold with five luer lock ports; 1 hot melt glue gun (110V only); 3 glue sticks.
 13316 Mini Glue Gun and (3) glue sticks

Patch Clamp Capillary Glass


Evaluated for best softening temperature, electrical properties, sealability, leachable components

PG52151-4 and **PG52165-4** are prepared from Schott #8250 glass (equivalent to Corning #7052), one of the most widely used patch clamping glasses. This is a specially formulated borosilicate glass with a softening temperature that is 110°C lower than regular borosilicate glass (Corning 7740, or Pyrex). It has excellent sealing properties for most cells. Electrical properties are also very good.

PG10150-4 and **PG10165-4** are composed of Corning #0010 glass, a high lead content (22% PbO) glass. Its thermal and electrical performance is between the Schott #8250 and Corning #8161 glasses described above. It is much more economical than Corning #8161 glass. It has been found that this glass causes much less alteration in channel behavior than Corning #8161 and Schott #8250 glass (Furman and Tanaka, *Biophys. J.* 53, p287, 1988).

Patch clamp capillaries do not have microfilaments.

ORDERING INFORMATION					
Patch Clamp Capillary Glass					
Catalog#	Glass Type	OD/ID (mm)	Dielectric Constant	Softening Point °C	QTY
PG52151-4	#8250	1.5/1.0	4.9	720°	100
PG52165-4	#8250	1.65/1.1	4.9	720°	100
PG10150-4	#0010	1.5/0.75	6.7	625°	100
PG10165-4	#0010	1.65/1.1	6.7	625°	100



Glass Handling Forceps



Borosilicate Glass Micropipettes

Eliminate the cost and trouble of making your own micropipettes U ip™



Features

- Schott Duran borosilicate glass
- 0.5 µm and smaller ID micropipettes include internal glass fiber for easy filling
- Tip inner diameter tolerance ±20%
- Short taper yields high strength
- Nominal length ≈ 50 mm
- OD:ID = 1.33:1
- Standard capillary outer diameters are 1.0 mm (thin-wall) or 1.14 mm
- Every pipette individually tested and inspected
- Vacuum sealed packs of 10

Benefits

• Plain Shank or Luer Fittings

Applications

• Injection of dyes or proteins into cells, oocytes or othe biomedical laboratory applications

WPI can quickly supply your need for consistently sized pre-pulled glass micropipettes.

Tip diameters (ID) range from 0.1 to 15 μm. Silanized Tips (Luer Shank) are available. Silanization waterproofs the glass to retard water when inserting into cell. This will not let the outside fluid run down the pipette and get inside so easily.

			ORDERING I	NFORM <i>A</i>	ATION		
Shank	Tip I.D.	Silanize Coating Length	Glass O.D.	Filament	Fire Polished	Catalog #	
PLAIN	0.1 µm	_	1.0 mm Thin-Wall	Yes	No	TIP01TW1F	
	0.2 µm	_	1.0 mm Thin-Wall	Yes	No	TIP02TW1F	
	0.3 µm	_	1.0 mm Thin-Wall	Yes	No	TIP03TW1F	
	0.4 µm	_	1.0 mm Thin-Wall	Yes	No	TIP04TW1F	
	0.5 µm	_	1.0 mm Thin-Wall	Yes	No	TIP05TW1F	
	1 µm	_	1.0 mm Thin-Wall	No	Yes	TIP1TW1	
	2 µm	_	1.0 mm Thin-Wall	No	Yes	TIP2TW1	
	5 µm	_	1.0 mm Thin-Wall	No	Yes	TIP5TW1	
	10 µm	_	1.0 mm Thin-Wall	No	Yes	TIP10TW1	
	10 μm	_	1.14 mm A203XV glass *	No	Yes	TIP10XV119	
	_	_	1.0 mm Thin-Wall	No	Yes	TIP10LT	Long Taper †
	_	_	1.0 mm Thin-Wall	Yes	Yes	TIP10FLT	Long Taper †
	_	_	1.5 mm Thin-Wall	No	Yes	TIP15LT	Long Taper †
	_	_	1.5 mm Thin-Wall	Yes	Yes	TIP15FLT	Long Taper †
	30 µm	_	1.0 mm Thin-Wall	No	Yes	TIP30TW1	
LUER	0.1 µm	_	1.0 mm Thin-Wall	Yes	_	TIP01TW1F-L	
	0.2 µm	_	1.0 mm Thin-Wall	Yes	_	TIP02TW1F-L	
	0.3 µm	_	1.0 mm Thin-Wall	Yes	_	TIP03TW1F-L	
	0.5 µm	_	1.0 mm Thin-Wall	Yes	_	TIP05TW1F-L	
	1 µm	_	1.0 mm Thin-Wall	No	_	TIP1TW1-L	
	2 µm	_	1.0 mm Thin-Wall	No	_	TIP2TW1-L	
	5 µm	_	1.0 mm Thin-Wall	No	_	TIP5TW1-L	
	10 µm	_	1.0 mm Thin-Wall	No	_	TIP10TW1-L	
	30 µm	_	1.0 mm Thin-Wall	No	_	TIP30TW1-L	
LUER/	5 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS01	
SILANIZED	5 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS02	
	10 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP10TW1LS01	
	10 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP10TW1LS02	
	30 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS01	
	30 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS02	

^{* 10} µm (ID), 1.14 mm capillary pipettes are for use in WPI's Nanoliter 2010.

µTIP SAMPLER ASSORTMENTS

TIPMIX01-05	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, plain shank
TIPMIX05-10	Two each, 0.5, 1, 2, 5, 10 µm ID, plain shank
TIPMIX01-05-L	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Luer
TIPMIX05-10-L	Two each, 0.5, 1, 2, 5, 10 μm ID, Luer



Micropipette Storage Jar

1965 1966 1967

Protect your micropipettes until ready to use them

Stores up to 30 micropipettes, filled or unfilled, up to three inches in length. A gentle sliding action inserts or removes pipettes without damage to the delicate tips.

	ORDERING INFORMATION	
E210	Storage Jar for 1.0 mm OD Micropipettes	
E212	Storage Jar for 1.2 mm OD Micropipettes	
E215	Storage Jar for 1.5 mm OD Micropipettes	
E220	Storage Jar for 2.0 mm OD Micropipettes	
REPLACEM	ENT PARTS	
1965	Foam Ring for 0.75 - 1.0 mm glass	
1966	Foam Ring for 1.2 - 1.5 mm glass	

Foam Ring for 2.0 mm glass

WORLD PRECISION INSTRUMENTS

[†] Long Taper micropipettes are pulled with a 12-15 mm taper which the customer cuts back to obtain the desired tip diameter.

Microelectrode Puller

A compact, versatile and reliable workhorse



Features

- Program sequences up to four steps
- Store up to 95 programs in memory for easy recall
- Supplied with a vial of capillary glass

Benefits

- Tempered glass cover to reduce the effects of humidity on puller reproducibility
- Switchable power supply ensures that line voltage fluctuations don't affect reproducibility

Applications

• Pull your own microelectrodes and micropipettes

PUL-1000 is a microprocessor controlled horizontal puller for making glass micropipettes or microelectrodes used in intracellular recording, patch clamp studies, microperfusion or microinjection. The puller was designed with tight mechanical specifications and precision electronics for complete control of the pulling process and accurate reproducibility. It offers programmable sequences of up to four steps with heating, force, movement and cooling time. This allows graduated cycles for applications like patch clamp recording.

This puller is a reasonably priced, compact, versatile and reliable workhorse. The microprocessor, combined with the LCD display, makes the **PUL-1000** easy to use.

Tempered Glass Cover

The cover of pulling chamber is made with tempered glass to minimize the humidity effect on the reproducibility of pulled pipettes.

Switchable Power Supply

PUL-1000 has a high quality switching power supply for use anywhere in the world without worry about the line voltage differences. Pulling

reproducibility is unaffected by line voltage fluctuation. Heating voltage can be controlled to within 0.1% accuracy even when line voltage fluctuates from 90 to 240 VAC.

References

Plautz, C. Z., Williams, H. C., & Grainger, R. M. (2016). Functional Cloning Using a Xenopus Oocyte Expression System. *Journal of Visualized Experiments*, (107), e53518–e53518. http://doi.org/10.3791/53518

Komarova, Y., Peloquin, J., & Borisy, G. (2011). Components of a microinjection system. *Cold Spring Harbor Protocols*, 2011(8), 935–9. http://doi.org/10.1101/pdb.ip27

PUL-1000 SPECIFICATIONS

HEATER ELEMENT Platinum/Iridium
PULLING FORCE Solenoid, adjustable
CAPILLARY OD RANGE 1.0–2.0 mm
MAXIMUM CAPILLARY LENGTH 170 mm
MINIMUM CAPILLARY LENGTH 55 mm
MEMORY SETS 95

POWER 90-240 VAC, 50/60 Hz DIMENSIONS 34 x 24 x 12 cm

SHIPPING WEIGHT 15 lb.

REPLACEMENT FILAMENTS 2.5 mm Square Box Filament, 2.5 mm wide (13834)

ORDERING INFORMATION

PUL-1000 Micopipette Puller

13834 Filaments, Platinum/Iridium, 2.5mm wide 2.5mm Square Box

504951 Filaments, Platinum/Iridium, 4.5mm wide Trough

Micropipette Puller for Quartz and Glass

CO₂ laser-based heat source for pulling quartz pipettes

Features

- Capable of pulling quartz, borosilicate and aluminosilicate glass
- Fully programmable, including heating filament characteristics
- Pulls electrodes with tip diameters smaller than 0.03 µm
- Preprogrammed sample programs for intracellular and patch pipettes

Benefits

 The laser has no melting point limit as with conventional metal filaments and cannot be burned out

- Optimized velocity sensing circuit for maximized sensitivity and reproducibility
- Operating life of the CO₂ laser is expected to be in excess of ten years with normal use, after which the laser can be refurbished
- Individual programs can be writeprotected in order to secure them from inadvertent changes
- The total time that the heat is on during the pull is displayed for improved program development and troubleshooting
- A date and time stamp is displayed to show the last time that a program has been changed



- Patch clamp, single isolated and whole cell
- · Intracellular recording
- Microinjection
- Nano probe research
- SECN
- P-2000/F is ideal for applications such as nanospray and NSOM

The **SU-P2000** micropipette puller represents a significant advance in the technology of fabrication of micropipettes, optical fiber probes, and nanospray tips. The **SU-P2000** integrates a CO_2 laser-based heat source with the technology derived from our extensive experience with conventional pullers. This system offers capabilities unmatched by other pullers.

While the **SU-P2000** is suitable for working with most conventional glasses, its primary advantage is the ability to work with quartz glass (fused silica). Quartz offers superior material properties for a variety of research applications. Quartz is stronger than other glasses and can facilitate penetration through tough tissues which would normally break conventional pipettes. For applications requiring a low noise glass, users will find that quartz is the lowest noise glass available. Quartz contains none of the metals used in conventional glasses. Optically, quartz is virtually free from fluorescence when illuminated.

A CO_2 laser was selected as the heat source for the SU-P2000 for several reasons:

- The nominal emission wavelength of the laser approximates the resonant frequency of the SiO₂ lattice in glass. Thus, quartz and other conventional glasses can be melted when the appropriate laser power is supplied.
- Laser heat is clean and leaves no metal residue on the pipette as do conventional heating filaments.
- Laser heat can be turned off instantly, leaving no residual filament heat.
- The user can program the amount and distribution of heat supplied to the glass.
- Laser heat source means there are no filaments to burn out or replace.

The **SU-P2000** can store up to 100 separate programs, with each program consisting of up to 8 command lines. Programmable parameters include: laser power level, scan width, trip velocity, delay/laser on time, and hard pull strength.

One important consideration for the use of the **SU-P2000** is the diameter of the glass used. The **SU-P2000** is designed to produce even heating on glass up to 1.2 mm in outside diameter. Larger diameter glass can be used with the **SU-P2000** (up to 1.5 mm quartz and 1.8 mm conventional glass), but the performance is best with glass that is 1.2 mm diameter or less.

SPECIFICATIONS

DIMENSIONS 30 x 14.25 x 13.25 in. (76 x 36 x 33.5 cm)

WEIGHT 90 lbs. (41 kg) ELECTRICAL 115/230 V, 50/60 Hz

*Patent No.4,600,424 CLASS 1 LASER PRODUCT

MODEL P-2000

ORDERING INFORMATION

SU-P2000 Micopipette Puller

References

Munoz, J.L. and Coles, J. Quartz micropipettes for intracellular voltage microelectrodes and ion selective microelectrodes. *Journal of Neuroscience Methods*: 22:57-64. 1987.

Rae, J.L. and Levis, R. A. A method for exceptionally low noise single channel recordings. *European Journal of Physiology* - Pflügers Archiv: 420:618-620, 1992.

Zuazaga, C. and Steinacker, A. Patch-clamp recording of ion channels: Interfering effects of patch pipette glass. *News in Physiological Science*: 5:155-159. 1990.

Levis, R.A. and Rae, J. L. The use of quartz patch pipettes for low noise single channel recording. *Biophysical Journal*: 65:1666-1677, 1993.

Next Generation Micropipette Puller

Flaming/Brown style puller with color touch screen control

Features

- Color touch-screen display
- Help topics and error detection
- Copy and paste functions for writing new programs
- Record of last two pull results
- Two symmetrical pipettes with each pull
- Self-contained air supply with filtration system and humidity control chamber
- Memory storage for up to 100 programs
- Two cooling modes: time and delay

Benefits

- Safe heat mode to protect and extend filament life
- Pre-heat mode improves stability
- Line repeat mode simplifies multi-line programming
- · Pipette Cookbook program directory
- Glossary with micropipette and puller terminology
- Jaw temperature sensor helps define ideal pulling conditions
- Ramp test more easily accessed and can now be stored and referenced within each program. Helps to establish program heat settings and protect filament.

Applications

- Patch pipettes
- Sharp electrodes
- · Pronuclear injection
- · Zebrafish injection
- · Insect egg injection
- · Aspiration pipettes

The **SU-P1000** micropipette puller was developed through years of experience with the Flaming/Brown style micropipette pullers and infused with leading edge technology. The most obvious new feature is the color touchscreen display that provides an intuitive and full-featured interface.

The extensive library of programs found in the popular Sutter Cookbook has been incorporated into the puller and is available to the user via the touchscreen display. You need only specify the glass, filament, and type of pipette you require and a suitable program will be identified and available for installation. This takes the guesswork out of pipette pulling and simplifies programming.

The Pre-heat mode actively heats and controls the jaw temperature and assures that the jaws have reached a specific temperature before the glass is pulled. This can increase the stability of the program from pull to pull. Copy and paste functions assure that programs can be easily written and the line repeat mode simplifies multi-line programming. A safe heat mode is an additional feature that helps the user avoid using heat settings that might damage or burn out the filament. When the safe heat mode is turned on, the puller will "check" the installed heat and alert the user when a given heat value is too low or too high in relation to the ramp value. User notes can be added to each program for annotating important information.

New features

- Diagnostic testing of all puller components
- Built-in error detection of air pressure loss or filament burnout
- · Easy access to ramp test



- Measurement of jaw temperature
- Access to previous pull results with the heat on times for each cycle of the program
- Help topics are pre-loaded to assist with on-site troubleshooting and the built-in glossary includes text, pictures, and diagrams explaining the terms used in micropipette fabrication
- Rotary dial is offered as an alternative to the keypad for numerical entry

SU-P1000 SPECIFICATIONS

 DIMENSIONS
 21 x 14 x 12" (53 x 36 x 30 cm)

 WEIGHT
 41 lb. (18.59 kg)

 POWER
 115/230 V, 50/60 Hz

	ORDERING INFORMATION
SU-P1000	Pipette Puller
OPTIONAL	ACCESSORIES/REPLACEMENT PARTS
505094	2 mm Box Filament, 1.5 mm wide
505095	2 mm Box Filament, 2 mm wide
505096	2 mm Box Filament, 3 mm wide
505097	2.5 mm Box Filament, 4.5 mm wide
13834	2 mm Box Filament, 2.5 mm wide
505098	3 mm Box Filament, 1.5 mm wide
14074	3 mm Box Filament, 2 mm wide
506122	3 mm Box Filament, 3 mm wide
505099	Custom Filament
505101	3 mm Loop Filament, 1.5 mm wide
505103	3 mm Loop Filament, 2 mm wide
505105	3 mm Loop Filament, 2.5 mm wide
505106	Fire Polishing Spacer
505108	Trough Filament, 1.5 mm wide
13835	Trough Filament, 2 mm wide
13836	Trough Filament, 3 mm wide
505109	Trough Filament, 4.5 mm wide

Flaming/Brown Pipette Puller

Create custom micropipettes, patch pipettes and microinjection needles

Features

- · Environmental chamber for humidity control
- Programmable air pressure
- Memory storage for up to 100 programs
- Write protection and date stamp for each program
- Two symmetrical pipettes with each pull
- Two cooling modes: time and delay
- · Vacuum fluorescent display

Benefits

- Preprogrammed sample programs for intracellular and patch pipettes. Special programming on request
- Ramp test self test for establishing program heat settings when a new filament or glass is introduced
- Constant current power supply for filament and pull solenoid
- Looping pull cycle for fabrication of patch type micropipettes
- Self-contained air supply with filtration system and humidity control
- \bullet Consistent and reliable electrodes with tip diameters less than 0.1 μm
- Control over the time and pressure at which the air is delivered
- Optimized velocity sensing circuit for maximized sensitivity and reproducibility
- Quality control, SEM photograph of a tip pulled with each puller.
 Criterion is tip measurement less than 0.1 µm and typically is ~0.06 µm.

Applications

- Patch pipettes
- Sharp electrodes
- Pronuclear injection
- · Zebrafish injection
- Insect egg injection
- · Aspiration pipettes

The **SU-P97** Flaming/Brown type micropipette puller is ideal for fabricating micropipettes, patch pipettes and microinjection needles. While retaining many of the features of earlier models, the **SU-P97** offers improvements in mechanical, electronic and software design. The result is better control of the pulling process and a higher degree of reproducibility. The **SU-P97** combines a proven mechanical system with a sophisticated, programmable microprocessor controller. This programmable control of the pulling parameters allows the investigator to design application specific pipettes from a wide range of glass compositions and sizes.

A number of other features have been incorporated in the design of the SU-P97. Most apparent is the environmental chamber which surrounds the heating filament. This environmental chamber is designed to minimize the effect of changing humidity on the reproducibility of pulled pipettes. A 25% increase in power over the P-87 allows for the use of larger heating filaments, larger diameter glass and multi-barreled glass. The metal jaws that clamp the heating filament have also been redesigned to minimize heat retention. A gas delivery mode switch provides for extended cooling for large diameter and multi-barreled glass. A spring-loaded clamping mechanism has been added for easier loading of glass. A vacuum fluorescent display has been added that allows easy viewing.

Software improvements on the **SU-P97** include a display of the total heaton time to assist in program development and troubleshooting. Up to 100 programs can now be written and stored in memory, which makes the **SU-P97** suitable for multiple users. These programs can now be write-protected, adding security to prevent programs from being changed or altered inadvertently. The display shows the last date and time the



program was written or edited. In addition, the air pressure is included as a programmable parameter.

SU-P97 SPECIFICATIONS

DIMENSIONS 21 x 16 x 12 (53 x 40.6 x 30 cm)
WEIGHT 50 lb. (23 kg)
ELECTRICAL 115/230 V, 50/60 Hz

	ORDERING INFORMATION
SU-P97	Pipette Puller
OPTIONAL	ACCESSORIES/REPLACEMENT PARTS
505094	2 mm Box Filament, 1.5 mm wide
505095	2 mm Box Filament, 2 mm wide
505096	2 mm Box Filament, 3 mm wide
505097	2.5 mm Box Filament, 4.5 mm wide
13834	2 mm Box Filament, 2.5 mm wide
505098	3 mm Box Filament, 1.5 mm wide
14074	3 mm Box Filament, 2 mm wide
506122	3 mm Box Filament, 3 mm wide
505099	Custom Filament
505101	3 mm Loop Filament, 1.5 mm wide
505103	3 mm Loop Filament, 2 mm wide
505105	3 mm Loop Filament, 2.5 mm wide
505106	Fire Polishing Spacer
505108	Trough Filament, 1.5 mm wide
13835	Trough Filament, 2 mm wide
13836	Trough Filament, 3 mm wide
505109	Trough Filament, 4.5 mm wide

Digital Microforge

Microforging, Micropipette Calibration and Microinjection — in a single device!

Features

• Microprocessor-controlled microforge

Digital signal processor technology precisely controls the polish

Unique digital pneumatic pressure feature polishes the tip without changing the size

• W30S-LED Microscope (Optional)

Benefits

- Digital Signal Processor (DSP) technology
- · Complete system package available
- Kohler illuminator and Abbe condenser for less glare and sharper images.
- Pneumatic pressure polishing that allows the preparation of blunt tips without change of tip ID
- Heating filament is attached to the microscope objective so they move together
- Pipette holder sits on the microscope stage to simplify the locating and polishing of the pipette



- · Polishing patch pipettes
- Microforging holding pipettes
- Microforging beveled injection pipettes
- Pipette tip calibration and microinjection

The **DMF1000** is a microprocessor-controlled microforge offering unmatched performance. Designed for fabrication of both small patch clamp glass pipettes and larger injection pipettes, the DMF1000 should find many uses in the laboratory. The DMF1000 is based on a design similar to that first used in WPI's extremely popular microforge model, the MF200. The extensive improvements incorporated into the DMF1000 greatly increase its versatility and performance, making it one of the most powerful microforges on the market.

Digital Signal Processor (DSP) Technology


The **DMF1000** is powered by the latest digital signal processor (DSP) technology. A digital timer is used to precisely control the polish heating time. Ten memories can be used to store settings of the heating power and heating duration. All of the settings are controlled and displayed digitally for better accuracy and reproducibility. Two different operating modes are provided: Manual and Auto. In the Manual mode, the DSP will memorize the duration of the time that is used to achieve a desired polishing. In Auto mode, the heat will be applied for the duration of the timer setting.

Complete System Available


The DMF1000 system includes a specially configured WPI model W30S-LED research grade compound microscope (optional) equipped with a high quality metallurgic 40x long-working distance objective and a pair of 10x eyepieces. The long working


UK: +44 (0)1462 424700 • wpiuk@wpi-europe.com • www.wpi-europe.com


distance objective reduces the danger of damage to the objective lens during the heating process.




stretching)











Filament Holder mounts directly to objective to provide precise control of heating element position.

Kohler Illuminator and Abbe Condenser

DMF1000-M1

Other benefits of the DMF1000 design include the use of a Kohler illuminator and Abbe condenser, which provide the reduced glare and sharper image contrast necessary when polishing pipettes as small as half a micron (0.5 µm) in diameter.

Pressure Polishing

The **DMF1000** incorporates a unique digital pneumatic pressure feature that enables pressurized air to be delivered through the pipette during fire polishing. In the fabrication of patch pipettes, the pressurized air can be used to blunt the taper at the pipette tip without changing the size of the tip opening. This reduces electrical resistance of the tip, leading to lower noise during patch-clamp recordings (Goodman & Lockery, 2000).

The Heating Filament

With a conventional microforge often the most difficult and timeconsuming part of using a high magnification objective is being able to move both the heating filament and the pipette into the same viewing area. Finding and moving both the heating filament and the pipette without collision can be a challenge. However, this difficulty is eliminated with the DMF1000 because the heating filament is directly attached to the microscope's objective. Hence it can be easily adjusted to any position within the viewing area.

The low heat capacity and low thermal coefficient of linear expansion of the filaments are key design features. The low heat capacity of the filament allows it to reach fire-polishing temperatures without excessive heat. This permits you to bring the pipette tip close to the filament during polishing without fear of collapsing the pipette tip. Low heat capacity eliminates the need for an auxiliary air-cooling system. The low coefficient of expansion characteristic of the filament ensures minimal displacement of the filament during heating. This feature eliminates much of the guesswork out of tip placement in relation to the filament.

Two different heating filaments are provided to accommodate various applications. The H5 filament is large gauge and can be reformed into a "U" for fabrication of pipettes, air forming of patch pipettes and other applications. The **H4** is a smaller gauge filament and is ideal for polishing patch clamp pipettes.

Pipette Holder Sits on the Microscope Stage

The pipette rests on a specially designed holder that sits on top of the microscope stage. The position of the pipette, relative to the heating filament, is controlled by the (X, Y, Z) adjustment of the stage. This unique design makes locating and polishing the pipette extremely easy. The stage of the microscope has a high quality rail that gives precise, smooth and stable control of the pipette's movement. This configuration also eliminates the need and expense of an additional micromanipulator to control pipette movement.

Typical applications of the DMF1000

Polishing the Patch Pipettes

Proper fire polishing of patch pipettes is the single most important factor for forming a stable giga-seal in patch clamp recording. This is even more important than the type of glass capillary used. Difficulties often arise in forming giga-seals because the polishing of patch pipettes using a conventional low magnification microforge is inadequate. However, the DMF1000 uses a 40X long-working distance objective. Pipette polishing is much more accurately controlled. Both whole cell patch pipettes and single channel patch pipettes can be conveniently polished to the highest quality and reproducibility achievable with any microforge.

Microforging Holding Pipettes

A holding pipette with a large blunt tip and a small opening is used to hold a floating cell in place prior to microinjection by applying suction to the rear of the pipette. The procedure for making holding pipettes involves three steps: squaring off, large bore flame polishing and tip reducing. These steps are accomplished with a larger heating filament.

Microforging Beveled Injection Pipettes

Occasionally, a beveled large bore pipette is not sharp enough to penetrate a cell without damaging the area around the pipette. With the **DMF1000** and the large heating filament, a sharp point can be formed on the beveled tip to assist the penetration of the cell. This process is referred to as contact stretching.

Pipette Tip Calibration & Microinjection

The integrated digital pneumatic pressure system can be used to calibrate the precise diameter (I.D.) of a micropipette tip, based on a technique described previously (Hagag & Randolph 1990, Bowman & Ruknudin 1999). The pressure system can also be used separately as a simple but highly accurate controller for microinjection applications.

References

Wu, Z.-Z., Chen, S.-R., & Pan, H.-L. (n.d.). Differential Sensitivity of N-and P/Q-Type Ca²⁺ Channel Currents to a Opioid in Isolectin B 4 -Positive and -Negative Dorsal Root Ganglion Neurons. http://doi.org/10.1124/jpet.104.073429



SPECIFICATIONS

AC POWER MODULE	100-240 VAC 50/60 Hz
TIMER RANGE (for heater & timer)	0.01 to 360 sec
NUMBER OF MEMORYS	10

PRESSURE ADJUSTING RANGE 0.5 - 60 PSI (3.5 - 414 kPa) PRESSURE RESOLUTION 0.1 PSI (0.7 kPa)

FILAMENTS: H4 Small filament for working with 40x long working distance objective.

FILAMENTS: H5 Large filament for working with 10x objective. Filament adjustment assembly provided for both objectives.

HEATER AND TIMER CONTROL Auto or Manual via Pushbutton, TTL,

or Optional Foot switch.

4 x 7 x 17 in. (10.2 x 17.8 x 4.8 cm) DIMENSIONS: Control Unit SHIPPING WEIGHT 4 lb. (1.8 kg)

MICROSCOPE See W30S, page 187 MICROSCOPE: SHIPPING WEIGHT 16 lb. (7.3 kg)

ORDERING INFORMATION

DMF1000-1	Microforge, w W30S-LED Microscope (110 V)
DMF1000-2	Microforge, w W30S-LED Microscope (220 V)
DMF1000-M1	Microforge without microscope (110 V)
DMF1000-M2	Microforge without microscope (220 V)

*Above DMF1000 microforges include 40X long working distance objective

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

500329	25x Long Working Distance Objective, 5 mm 0.5 NA
13142	Optional foot switch
800292	40x Long Working Distance Objective, 3mm 0.25 NA
503513	21 mm 10X Eyepiece with 100/10 reticle
DMF1000-H5	Replacement heating filament (large gauge)
MF200-H4	Replacement Heating filament (small gauge)
75050	Replacement Micropipette Slide
75040	Replacement Filament Cable

Analog Microforge

Sometimes the simplest designs work best

Features

- Simple, reliable and economically priced
- Analog temperature controller
- W30S-LED microscope (optional)

Benefits

- Includes 40x long-working distance objective and 10x eyepiece
- Kohler illuminator and Abbe condenser for less glare and sharper images.

Applications

- Patch pipette tip polishing
- Micropipette tip size reduction
- Contact stretching in in vitro fertilization pipette production

The **MF200** Microforge is a versatile instrument designed specifically for the fabrication of glass micropipettes and other related tools. The system was developed in collaboration with Dr. Ming Li of the

Department of Pharmacology, University of South Alabama. The **MF200** is simple, reliable and priced economically.



The **MF200** system includes: An easy to use analog temperature controller, a specially configured WPI model **W30S-LED** research

grade compound microscope (optional),
40x long-working distance objective
and 10x eyepiece. 40x magnification
is essential when polishing pipettes as
small as half a micron (0.5 µm) in diameter.
Compared to a conventional 40x objective,
the long working distance objective reduces
the danger of damage to the pipette and/or
objective lens during the polishing process.

It is the only commercial microforge using the Kohler illuminator and Abbe condenser for illumination.

This provides less glare and sharper image of the pipette than the frosted glass illuminator, which is used on other commercial microforges.

the danger of damage to the piper objective lens during the polishing Kohler Illuminator

References

Guillou, L., Babataheri, A., Puech, P.-H., Barakat, A. I., & Husson, J. (2016). Dynamic monitoring of cell mechanical properties using profile microindentation. Scientific Reports, 6, 21529. http://doi.org/10.1038/srep21529

Vasauskas, A. A., Chen, H., Wu, S., & Cioffi, D. L. (2014). The serine-threonine phosphatase calcineurin is a regulator of endothelial store-operated calcium entry. Pulmonary Circulation, 4(1), 116–27. http://doi.org/10.1086/675641

FEATURE COMPA	RISON	
	MF200	DMF1000
W30S-LED Microscope	V	V
40x Long Working Distance Objective	V	✓
Analog Controller	/	/
Digital Controller		/
Pressurized Air Control		V
Microinjection Capability		V
Optional Foot Switch	V	V
Memory		V
Auto-sense of Filament Type		V
Digital Temperature Control		V



SPECIFICATIONS

AC POWER MODULE 100-240 VAC 50/60 Hz

FILAMENTS (3) H2, H3, H4

FILAMENT ON Pushbutton Controlled or Optional Foot

Switch Controlled

FILAMENT ADJUSTMENT ASSEMBLY For 40x and 25x Long-Working,

Distance Objectives: mounts on objec-

FILAMENT ADJUSTMENT ASSEMBLY: 40x Long-Working Distance (3 mm)

OBJECTIVE

FILAMENT ADJUSTMENT ASSEMBLY: 25x Long-Working Distance (5 mm)

OPTIONAL

EYEPIECE 10x (pair)

EYEPIECE: RETICLE (10x eyepiece 1.25 µm/division (at 40x), 0-90° Angle at

only) (OPTIONAL) 5°/division

EYEPIECE: OPTIONAL EYEPIECE 15x (pair)

GLASS HOLDER Mounts on Microscope Stage

DIMENSIONS: Control Unit $4 \times 7 \times 1.875$ in. (10.2 x 17.8 x 4.8 cm)

SHIPPING WEIGHT 3 lb. (1.4 kg)
MICROSCOPE: Note See W30S-LED
MICROSCOPE: SHIPPING WEIGHT 16 lb. (7.3 kg)

	ORDERING INFORMATION
MF200-1	Microforge System w. W30S-LED Microscope (110 V)
MF200-2	Microforge System w W30S-LED Microscope (220 V)
MF200-M1	MF 200 without microscope (110 V)
MF200-M2	MF 200 without microscope (220 V)
*Above MF200 microforges include 40X long working distance objective	

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

500329	25× Long-Working Distance Objective (fits most
	microscopes with a 160 mm Focal Length)
13142	Optional foot switch
MF200-H3	Replacement heating filament (medium gauge)
MF200-H4	Replacement heating filament (small gauge)
75090	Filament Adjustment Assembly,22mm OD Objectives
75050	Replacement Micropipette Slide
75040	Replacement Filament Cable

Micropressure System

Measure hydrostatic pressure in small vessels and oocytes

Features

- Utilizes a liquid filled micropipette (2–5 µm tip opening) for sensing pressure
- Pressure range from 1–100 mmHg (pressure range to 350 mmHg is available)
- Lower limit 1 mmHg (133 Pa)
- Includes calibration/test chamber
- Tubing and fittings for interconnecting system sub-components are provided
- Pressure in the pipette can be manually set to positive or negative relative to the outside
- Probe holder for mounting on micromanipulator is included.
- 10 pre-pulled pipettes are included

Benefits

- Measures biological pressures in very small liquid (aqueous) filled spaces
- · Pre-configured pipettes are available for convenience

Applications

- · Pressure in kidney tubules
- · Embryonic blood pressure
- Intracellular pressure
- Mouse intraocular pressure

The **900A** system is designed to measure liquid pressures dynamically in aqueous biological micro-environments, such as in kidney tubules or intracellular pressures. A liquid filled micropipette is used as a pressure probe, and pressure external to the pipette is measured at the tip. The outside diameter of the micropipette tip typically measures between 2–7 µm. Pressure measurement is achieved by monitoring the pipette's electrical resistance. The resistance changes according to changes in the pressure outside the pipette tip via displacement of an electrolyte concentration gradient. As the position of the concentration gradient changes, the resistance of the pipette changes. The resistance signal from the pipette is used as feedback to control a pressure source that applies pressure to the inside of the pipette to counterbalance pressure from the outside. The feedback loop forces the gradient to a neutral balance point, which is user-defined at atmospheric pressure beforehand. The internal pressure required to equally balance the external pressure to the neutral point is readily measured, and it is converted into an analog voltage available at the pressure output BNC and displayed numerically on the LED meter.

System requirements

The **900A** requires stable sources of both pressure and vacuum, which are essential for the system to rapidly counteract changing pressures encountered at the pipette tip. Pressure and vacuum sources are not provided with the **900A** instrument because some labs are already equipped with suitable sources of pressure and vacuum. For researchers who do not possess pressure or vacuum sources, a cylinder of compressed air or inert gas with a dual stage regulator serves very well as a pressure source. Vacuum must be very stable. It is often best supplied by a quality vacuum pump. WPl offers a very quiet continuous duty vacuum pump well suited for use with the system. In addition, a vacuum regulation kit is recommended to fine tune the vacuum source to the ideal level (**900A-VAC**).

A manometer or meter for independent pressure measurement is necessary to calibrate the pressure and vacuum sources, as well as for validation of the performance of the **900A** system prior to experimental use. A pressure measurement device capable of measuring within a range of +300 mmHg and -150 mmHg is recommended (**PM015D** or **PM015R**). For system performance validation at pressures well below 100 mmHg, the **PM01D** or **PM01R** is recommended, because it provides higher resolution at low pressure.

For transient response performance evaluation, a rapid burst of air or water pressure is required. WPl's **PV830** or **PV820** series PicoPumps provide this capability. Rapidly occurring transient pressure measurements are typically captured on a data acquisition system. For details about testing and measurement of rapidly occurring pressure phenomena, contact a WPI sales representative for additional information.



References

Inamoto, R., et. al. (2017). The difference in endolymphatic hydrostatic pressure elevation induced by isoproterenol between the ampulla and the cochlea. Auris Nasus Larynx, 44(3), 282–287. http://doi.org/10.1016/J.ANL.2016.07.018

Wei, J., et. al. (2017). Role of intratubular pressure during the ischemic phase in acute kidney injury. American Journal of Physiology - Renal Physiology, 312(6), F1158–F1165. http://doi.org/10.1152/ajprenal.00527.2016

Inamoto, R., et. al. (2017). The difference in endolymphatic hydrostatic pressure elevation induced by isoproterenol between the ampulla and the cochlea. Auris Nasus Larynx, 44(3), 282–287. http://doi.org/10.1016/j.anl.2016.07.018

Wei, J., et. al. (2017). Role of intratubular pressure during the ischemic phase in acute kidney injury. American Journal of Physiology-Renal Physiology, 312(6), F1158–F1165. http://doi.org/10.1152/ajprenal.00527.2016

Petrie, R. J., et. al. (2017). Activating the nuclear piston mechanism of 3D migration in tumor cells. The Journal of Cell Biology, 216(1), 93–100. http://doi.org/10.1083/icb.201605097

Warmerdam, T., et. al. (n.d.). Perilymphatic and endolymphatic pressures during endolymphatic hydrops. European Archives of Oto-Rhino-Laryngology, 260(1), 9–11. http://doi.org/10.1007/s00405-002-0518-2

Petrie, R. J., & Koo, H. (2014). Direct measurement of intracellular pressure. Current Protocols in Cell Biology / Editorial Board, Juan S. Bonifacino ... [et al.], 63, 12.9.1–9. http://doi.org/10.1002/0471143030.cb1209s63

SPECIFICATIONS

PRESSURE RANGE ±0-100 mmHg (900APP100) LINEARITY < ±0.5% from a straight line **STABILITY** ±0.1 mmHg up to 1 hour or more **ACCURACY** ±0.5% of full scale RISE TIME >10 ms (10-90%), based on residual volume OUTPUT ("Pressure Signal") 10 mV/mmHg AMPLIFIER PROBE Input Resistance >1010Ω, Voltage Gain 1.0 **DIMENSIONS (MAIN FRAME)** 43.2 x 13.3 x 25.4 cm (17 x 5.25 x 10 in.) DIMENSIONS (PRESSURE POD) 9.4 x 2.5 x 5.7 cm (3.7 x 1 x 2.25 in.) 110 VAC/220 VAC

ORDERING INFORMATION

SYS-900A	Micropressure System
Specify line vo	oltage. See training videos at www.wpiinc.com/sys-900a.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

900AP	Replacement Probe
CAL900A	Pressure Calibration Chamber
3491	Probe Extension Cable
2933	Rack Mount Kit
5332	Replacement Liquid Trap
MEH6RF	Micropipette Holder (1.0, 1.2, 1.5 or 2.0 mm - Specify OD)
TIPTW900A	Pre-pulled Micropipette for 900A (1 mm thin-wall, 2 μTip)
	(pkg. of 10)
900APP100	Replacement Pressure Pod
SYS-PM015D	Pressure Manometer (15 psi)
SYS-PM01D	Pressure Manometer (1 psi)
801566	Vacuum Pump, 110V
801963	Vacuum Pump, 220V
900A-VAC	Vacuum Pump Regulation Klt

Pressure Manometer

For measuring hydrostatic pressures

Features

- Measure vacuum and pressure in gases
- Ranges available: ±1 PSI, ±15 PSI, ±100 PSI
- Measure in PSI or kPa on the 100 PSI unit or PSI and mmHg on the 15PSI unit
- · Battery powered
- Includes tubing and mini-phone to BNC cable

Benefits

• Easy and accurate measurements

Applications

• Measure pressure of non-corrosive gases

Hand-held and battery operated, PM Series pressure manometers monitor vacuum and pressure in non-corrosive gases. An integral transducer and digital display allow easy and accurate pressure readings. Three versions measure pressures in the range of ± 1 PSI, ± 15 PSI or ± 100 PSI. A range switch allows measurement in units of PSI or kPa for the 100 PSI version, and PSI or mmHg for the 15 psi version. Pressure can be read on the built-in LCD display or relayed to a chart recorder, oscilloscope, or computer.

PM Series pressure manometers come with 4 feet of $\frac{1}{8}$ -in. ID soft vinyl tubing. A mini-phone-to-BNC cable for the recorder output is also available (**CBL102**). Standard versions are equipped with a nine-volt alkaline battery.




	PRESSURE MA	ANOMETER SPECIFICATION	IS
	PM01	PM015	PM100
PRESSURE RANGES	±1 psi (±52 mm Hg)	±15 psi (±775 mm Hg)	±100 psi (±690 kPa)
MAX. PRESSURE	20 psi (1035 mm Hg)	30 psi (1550 mm Hg)	150 psi (1035 kPa)
RESOLUTION	0.001 psi (0.1 mm Hg)	0.01 psi (1 mm Hg)	0.1 psi (1 kPa)
OUTPUT	1 V/psi	100 mV/psi	10 mV/psi
OUTPUT RANGE	±1.0 V	±1.5 V	±1.0 V
LINEARITY		0.5% full-scale	
TEMPERATURE EFFECT		1.0% full-scale (0-70°C)	
ZERO		Screwdriver-adjust	
RESPONSE TIME		30 ms	
POWER	Nine-volt battery		
BATTERY LIFE	Alkaline, 200 hours; rechargeable, 25 hours		
RECORDER OUTPUT	Mini-phone jack, 0.141 inch (3.5 mm)		
OUTPUT IMPEDANCE	1 kΩ		
PNEUMATIC CONNECTORS	Barbed, for ⅙-inch or ⅙-inch ID soft tubing		
DIMENSIONS	3 x 6 x 1 inches (8 x 15 x 4 cm)		
SHIPPING WEIGHT		3 lb (1.4 kg)	

	ORDERING INFORMATION
SYS-PM01D	Pressure Manometer (1 psi)
SYS-PM01R	Pressure Manometer (1 psi), Rechargeable*
SYS-PM015D	Pressure Manometer (15 psi)
SYS-PM015R	Pressure Manometer (15 psi), Rechargeable*
SYS-PM100D	Pressure Manometer (100 psi)
SYS-PM100R	Pressure Manometer (100 psi), Rechargeable*
CBL102	Mini-Phone-to-BNC Cable
	Specify line voltage

^{*}Rechargeable versions come with nickel/cadmium battery and charger

Animal Physiology



Specialized Equipment for Small Animal Research

WPI offers a broad range of instruments and laboratory supplies that are designed for small animal physiology. Let us help you get the laboratory equipment that's right for your research application.

Motorized Stereotaxic Frame

Automatically calculates coordinates and accurately places the probe

Features

- Accurate microstepping motor drive for high resolution placement
- Digital controller
- No computer required
- Touch screen controller for easy, intuitive control with better than 10 μm precision
- Graphic controller display for instant operational feedback
- Hand controller for complete manual control

Benefits

- Reliable positioning reduces errors
- Increased precision and repeatabilty of motion
- Feature rich for your convenience of use
- Multiple modes of operation for precise probe positioning

Applications

Stereotaxic surgery for small animal research

When precision and repeatability of motion are critical, the **MTM-3** Motorized Stereotaxic Frame and digital display outperforms manual

models, and it greatly reduces human error when performing routine stereotaxic surgery. The motorized axis of the **MTM-3** provides precise, controlled, 3-dimensional placement of any probe or accessory within the working space of a stereotaxic frame.

No computer is required. The **MTM-3** 3-axis manipulator supplied with WPI stereotaxic frames is compatible with standard stereotaxic frames and may be adapted to existing frames of other manufacturers.

Single and dual manipulator arm motorized systems are available. This allows you to mount a stereotaxic drill and a probe simultaneously. Both 3-Axis Wheel Control and Touchscreen Controller are included.

Reliable positioning for reduced errors

The controller automatically calculates the brain map coordinates for precise probe placement during a stereotaxic surgery. This eliminates the errors resulting from reading Vernier scales.

Increased precision and repeatability of motion

The MTM-3 offers better than 10 μ m accuracy with the microstepping motor drive for high resolution placement of your probe, drill or electrode. Coordinate distances are automatically calculated.

Convenient to use

Brain atlas coordinates may be entered into the controller, and no computer required. The simple to use touch screen controller offers ease of control. The graphic controller display provides instant operational feedback during a stereotaxic surgery. For example, the "final approach" speed may be set between 2 mm/sec and 0.02 mm/sec.



MTM-6, which includes both manipulator arms and dual controllers.

Multiple modes for positioning


The arms of the motorized stereotaxic frame may be positioned manually using the digital interface or hand controller. Use the Coordinates mode to control the arm positions based on specific mapped coordinates. Set your retracted position, and store up to three reference (origin) positions. Then, display the probe position with respect to any of your references or the absolute coordinates. If you prefer to write short scripts, computer control is also available which uses a simple, text based command set. For example, if your brain research requires that you routinely drill in the same location of the skull and implant a probe to the same prescribed depth, a script could dramatically improve your process.

	SPECIFICATIONS
TRAVEL	80 mm on each axis
PRECISION	10 μm

	ORDERING INFORMATION
MTM-3	3-Axis Stereotaxtic Frame (Left or Right)
MTM-6	6-Axis Stereotaxtic Frame (Left and Right)
MTM	3-Axis Stereotaxtic Arm/Controller only (Left), Upgrade
MTM-R	3-Axis Stereotaxtic Arm/Controller only (Right), Upgrade

Digital Stereotaxic Frame with LCD Display

Battery powered, low noise for electrically quiet operation



Features

- Large LCD display from sealed sensors on each axis
- Adaptors available for use with rats, mice, birds, cats, guinea pigs and other species
- 80 mm of travel in three planes
- 180° rotation and lock at any vertical angle
- 360° rotation and lock at any level angle
- Syringe pump and drill can be attached directly
- Stable and accurate movement
- Targets specific coordinates by zeroing
- 10 µm precision in all directions
- Extended base plate (400 mm x 255 mm) suitable for various animals
- Vertical lock and fixing knob are separated to ensure accurate position at any angle

Benefits

- Quiet, battery-powered operation
- Battery-powered sensors, without electronic noise, are suitable for electrophysiology experiments

Applications

- Stereotaxic surgery for small animal research
- Electrophysiology applications

This new Digital Stereotaxic Frame features sealed electronic sensors and an easy-to-read LCD display with 10 µm resolution in all three axes. A zeroing function aids in targeting specific coordinates.

- The precisely designed rotary knob and U frame allow sufficient space for the anterior-posterior operation
- Laser engraved scales and a darkened rod for easy to read scales
- Curved design of nose clamp fixes the head of the animal more securely

	ORDERING INFORMATION
502300	Digital Stereotaxic Frame with 18°Ear Bars
502350	Digital Stereotaxic Frame with 45°Ear Bars
502303	Dual Manipulator Digital Stereotaxic Frame, 18°Ear Bars
502353	Dual Manipulator Digital Stereotaxic Frame, 45°Ear Bars
TAXIC-300	Digital Stereotaxic Frame with 18° Ear Bars & UMP3T-1
TAXIC-303	Dual Manipulator Digital Stereotaxic Frame with 18° Ear
	Bars with UMP3T-1 syringe pump
TAXIC-350	Digital Stereotaxic Frame with 45° Ear Bars & UMP3T-1
TAXIC-353	Dual Manipulator Digital Stereotaxic Frame with 45° Ear
	Bars with UMP3T-1 syringe pump
UPGRADE MANUAL STEREOTAXIC TO DIGITAL	

502360 Manipulator, 3 Axes, Right Hand, Digital Display502361 Manipulator, 3 Axes, Left Hand, Digital Display

Precision Stereotaxic Instruments

For small animal research

Features

- 5 mm linear movement per revolution on each axis
- Absolute lock at 90° (vertical)
- Entire frame is electrically continuous, ideal for electrophysiology
- Accessories available for wide variety of small animals

Benefits

- · Versatile positioning
- · Easy to read scales
- Convenient for electrophysiology because the frame may be grounded

Applications


- Stereotaxic surgery for small animal research
- Electrophysiology research

WPI's Precision Stereotaxic Instrument is built around the time-proven U-frame design concept, providing stability and adaptability to most species. Precision alignment ensures accurate placement of electrodes, micropipettes and other devices. It is ideal for researchers in need of a versatile, reliable instrument for stereotaxic procedures with small animals

The manipulator arm controls mediolateral and vertical positioning via lead screws with 80 mm of travel. This allows the fastest positioning possible, consistent with lining up the scales easily at a given coordinate. The anteroposterior movement is controlled via a dovetail slide movement, with 80 mm of travel possible in each direction. A universal joint allows the investigator to change the angle of the probe up to 90° in either the anteroposterior or mediolateral planes. The locking mechanism will hold any angle position without drift or creep. It also provides an absolute lock at 90° vertical.

All scales are oriented to be read easily from the open end of the "U." This is the position from which most scientists prefer to work. The numerals on the scales are clear and easy to read. Precise alignment with facing Vernier scales gives accurate resolution to 0.1 mm.

The entire stereotaxic frame, including the dovetails, manipulator arms and base are electrically continuous. Grounding of the entire frame including the base plate can be accomplished by connecting the provided grounding stud to earth. This is ideal for electrophysiological studies where the animal and surrounding structures need to be grounded to reduce electrical noise.






	ORDERING INFORMATION
502600	Precision Stereotaxic Frame with 18°Ear Bars
502650	Precision Stereotaxic Frame with 45°Ear Bars
502603	Dual Manipulator Stereotaxic Frame with 18°Ear Bars
502653	Dual Manipulator Stereotaxic Frame with 45°Ear Bars
502612	Stereotaxic Frame, Single Manipulator, for Rat, 18° Ear Bars
502614	Stereotaxic Frame, Single Manipulator, Rail Mounted, for
	Rat, 18° Ear Bars
TAXIC-600	Stereotaxic Frame with 18°Ear Bars plus UMP3T-1
TAXIC-650	Stereotaxic Frame with 45°Ear Bars plus UMP3T-1
TAXIC-603	Dual Manipulator Stereotaxic Frame with 18°Ear Bars
	plus UMP3T-1 System
TAXIC-653	Dual Manipulator Stereotaxic Frame with 45°Ear Bars
	nlus LIMP3T-1 System

Mice and Neonatal Rats

Specialized equipment for the smallest subjects



Mouse/Neonatal frame with the optional digital controller

Parallel Rail Stereotaxic Instruments

For large animal research

Features

- 5 mm linear movement per revolution on each axis
- Entire frame is electrically continuous, ideal for electrophysiology
- Includes the U-frame base plate, manipulator(s), cat/monkey adaptor, 18° ear bars and swivel mount

Benefits

- · Versatile positioning
- Easy to read scales
- Convenient for electrophysiology because the frame may be grounded
- Ability to expand to accommodate a variety of large animals

Applications

- Stereotaxic surgery for large animals
- Electrophysiology applications



WPI's Parallel Rail Stereotaxic Frame systems are heavy-duty research instruments for large laboratory animals such as cats, monkeys and dogs. The solid, large frame and superior rigidity ensure the precise alignment of animals for stereotaxic surgery, injection and recording. The system can accommodate up to four manipulators with 100 μm resolution on each axis. Each manipulator can smoothly move to and lock at any location on both parallel rails in a range of 20 cm.

Parallel Rail Stereotaxic Frame system for large animals includes the Parallel Rail Frame, base plate, manipulator(s), cat/monkey or dog adaptor and swivel mount.

	ORDERING INFORMATION
502227	Stereotaxic Frame, 1 Manipulator for Cat/Monkey
502228	Stereotaxic Frame, 2 Manipulators for Cat/Monkey
502229	Stereotaxic Frame, 3 Manipulators for Cat/Monkey
502230	Stereotaxic Frame, 4 Manipulators for Cat/Monkey
502231	Stereotaxic Frame, 1 Manipulator for Dog
502232	Stereotaxic Frame, 2 Manipulators for Dog
502233	Stereotaxic Frame, 3 Manipulators for Dog
502234	Stereotaxic Frame, 4 Manipulators for Dog



502616MRI Stereotaxic Instrument for Dogs and Monkeys

ORDERING INFORMATION	
MRI Stereotaxic Frame, for Dogs and Monkey	'S


502616







0	RDERING INFORMATION (ADAPTORS)
502063	Mouse and Neonatal Rat Adaptor
502213	Platform, Gas Anesthesia, with Mouse Mask (use with 502063)
502265	Mouse Adaptor with Nose Clamp
502266	Rat Adaptor
502226	Cat/Monkey Adaptor for 502600 series
502259	Spinal Adaptor for Rat
502060	Guinea Pig Adaptor for 502600 series
502241	Dog/Monkey Adaptor for Parallel Rail Instruments, with Ear Bars, 18°
502259	Rat Spinal Adaptor, for Microscope Stage







	ORDERING INFORMATION (EAR BARS)
502055	Ear Bars, Rat, 18° (pair)
502056	Ear Bars, Rat, 45° (pair)
502224	Ear Bars, Cat, 18º (pair)
502225	Ear Bars, Cat, 45° (pair)
502235	Ear Bars, Mouse, 60°, Non-rupture (pair)
502242	Ear Bars, Rat, Hollow, 1.5 mm Hole for Auditory Stimulation
502249	Non-Rupture, Rubber Ear Bars for Mice
502250	Non-Rupture, Cuff Ear Bars for Mice

ORDERING INFORMATION (OTHER ACCESSORIES)

502257	MRI Rat Adaptor,36 x 12.5 x 9.6 cm
502201	V-Clamp, with 10/32 Screw for Mounting UMP3
504608	V-Clamp, with 10/32 Thumb Screw
502213	Platform, Gas Anesthesia with Mouse Mask (use with 502063)
502243	Adjustable Stage Platform for 502600 series (5 cm high)
503598	Micro-Drill, 35K rpm, with a set of bits, 110/220VAC
503599	Micro-Drill, 35K rpm, with a set of bits, 240VAC, UK plug
61840	Heating Plate for 502063, 4×15 cm (use with ATC2000)

Stereotaxic Anesthesia Facemask

Ensures a snug fit for mice and rats with zero leakage and zero dead space

Features

- Facemask custom fitted to a toothbar that fits a Cunningham Stereotaxic Frame
- Ports swivel 360°
- Optically clear
- Includes OC-tubing and OC-adapter (requires **EZ-1130**)

Benefits

- Anatomically correct design for a tight fit for mice and rats
- · Zero leakage
- Zero dead space
- Provides access to eyes and head of the rodent
- Mask provides stability of the animal for examinations and surgery
- · Reduces animal stress by limiting handling
- Reduces technician's exposure to gas

Applications

• Small animal stereotaxic surgery

Our custom toothbar for the Cunningham adaptor is designed to accommodate the use of a soft gas mask. The included mask securely fits the toothbar to eliminate gas leakage.

Choose either small or medium size. The **502063-SGM** and **502063-MGM** are customized versions of our standard rodent gas masks.

The OC gas masks allow you access to the eyes and head of the rodent while sealing the nose and mouth. The gas masks provide stability of the animal for examinations and surgery, and they can be easily incorporated in existing examination procedures. The wasteline connects to gravity or vacuum filters, and swiveling ports allow for greater flexibility in positions. They are easily adapted for stereotaxic applications.

Thanks to the anatomical face mask design, you have head-to-tail access. The face mask reduces the technician's exposure to gas by sequestering and filtering waste, minimizing dead space and leakage.

You can reduce rodent stress by limiting handling required when using injectable anesthesia. By choosing these gas masks with isoflurane and thereby reducing the use of Schedule II drugs, you can simplify your lab management. Setting up and changing masks is simple, and the swivel adapters allow you to make adjustments quickly and even modify the mask for stereotaxic uses. Gas anesthesia reduces recovery time, making for a faster and more efficient lab.



	ORDERING INFORMATION
502063-SGM Small Mask for Cunningham Stereotaxic Fram	
	Includes OC-TUBING and OC-ADAPTER
502063-MGM	Medium Mask for Cunningham Stereotaxic Frame
	Includes OC-TUBING and OC-ADAPTER

PLEASE NOTE: This part is made to fit the WPI neonatal frame only. If you have a neonatal frame of another make, please contact WPI.



	RODENT FACE	MASKS
*Requires the i	use of OC-TUBING, OC-ADAPT	ER and EZ-1130.
OC-SFM	Small Mask	Small Mouse 0.5 – 2 months
		8 mm Aperture x 25 mm overall length
OC-MFM	Medium Mask	Mouse 2+ months and Small Rat
		10 mm Aperture x 30 mm overall length
OC-LFM	Large Mask	Adult Rat up to 250 g
		13 mm Aperture x 35 mm overall length
OC-XLFM	Extra Large Mask	Adult Rat from 250 – 500 g
		18 mm Aperture x 42 mm overall length
502063-SGM	Small Mask for Cunningham	Small Mouse 0.5 – 2 months
	Stereotaxic Frame	8 mm Aperture x 25 mm overall
	*Includes OC-TUBING and OC-ADAPTER	length
502063-MGM	Medium Mask for Cunningham Stereotaxic	Mouse 2+ months and Small Rat
	Frame	10 mm Aperture x 30 mm

*Includes OC-TUBING and overall length



Stereotaxic Frame Anesthesia Masks

Active and passive anesthesia masks for use with stereotaxic equipment

Two kinds of stereotaxic frame nosecone masks are available:

- · Passive masks have collinear inlets and outlets.
- Active masks have separate (non-collinear) inlets and outlets on both

The passive mask is used with a charcoal canister only, and the active mask can be used with an evacuation apparatus.

Active Anesthesia



502046 (above) — 20 mm ID 502047 (not shown) — 22 mm ID



502043 (not shown) — 11.5 mm ID



ORDERI	NG INFORMATION (ACTIVE SCAVENGING)
502042	Stereotaxic Anesthesia/Evacuation Mask, Mice or
	Neonatal Rats, < 30 g
502043	Stereotaxic Anesthesia/Evacuation Mask, Mice or
	Neonatal Rats, 30-70 g
502046	Stereotaxic Anesthesia/Evacuation Mask, Rats ,< 300 g
502047	Stereotaxic Anesthesia/Evacuation Mask, Rats, > 300 g
502262	Mouse Gas Anesthesia Head Holder with Rubber Ear Bars

Passive Anesthesia





ORDERING INFORMATION (PASSIVE SCAVENGING)

502213	Platform, Gas Anesthesia, with Mouse Mask (use with 502063)
502053	Stereotaxic Mask, Gas Anesthesia, Mouse or Neonatal Rat (< 70g)
502054	Stereotaxic Mask, Gas Anesthesia, Rat (< 300 g)

Animal Temperature Controller

Maintain optimal animal temperature during research procedures

Features

- Animal Temperature Controller with Adaptive mode–auto adjust PID regardless of animal size
- Plate's internal temperature sensor prevents localized hot spots under animal with maximum temperature stability
- Compatible with RTD (resistive temperature device) and thermocouple probes
- Extremely quiet DC heater to facilitate electrophysiological recordings
- Three temperature sensor inputs
- Auto shutdown if the plate reaches 45°C

Benefits

- Low noise temperature control system
- · Adapts to different size animals automatically
- · Accepts RTD and thermocouple probes

Applications

- · Small animal surgery
- Stereotaxic surgery

The ATC2000 Animal Temperature Controller is a low-noise heating system for maintaining animal body temperature during experimental procedures. The DC heater is extremely quiet in terms of electromagnetic radiation. This is essential in electrophysiological recordings which are very sensitive to electromagnetic interference.

The controller uses proportional, integral, derivative (PID) technology to provide precise and stable control of a subject's temperature. Compared with switched on/off type controllers, PID controllers provide a much more precise and stable control of temperature. The PID approach is also more immune to the variation of the experimental conditions such as change in animal size and unexpected disturbances. Our unique adaptive mode technology automatically senses and adjusts the PID values based on the size of the animal.

Multiple temperature sensing inputs

The ATC2000 has three temperature sensing inputs.

- The resistive temperature device (RTD) probe input can be used to monitor an RTD rectal probe to control the animal temperature or to monitor ambient temperature, induction chamber temperature or any other temperature.
- When using a thermocouple probe, the thermocouple (TC) probe input can be used just like the RTD input. (A T type thermocouple must be used)
- The heater plate also has an internal temperature sensor. The ATC2000 monitors this sensor to prevent the localized hot spots under animal.

Operational modes

The controller has three operational modes:

- Normal mode uses the configured sensor (RTD or TC) or the plate sensor to drive the control loop.
- Adaptive mode uses the temperature of the heated plate and the temperature of the subject to control. This approach is less prone to overshoot, but somewhat slower than the normal mode, depending on the sampling rate used.
- Shutdown is a fail safe mode used if the plate temperature ever exceeds 45°C.

Required Accessories

The **ATC2000** is tuned at the factory. However, the PID parameters may also be set manually. The temperature resolution of the controller is 0.1°C. The rectal temperature probe has a 5-ft. ultra-flexible, shielded cable and an RTD sensor.



The metal heating plates (available separately) have built-in temperature sensors. The smaller heating plate (WPI# 61840) is compatible with stereotaxic systems; the rigid, flat surface of the warming pad fits under the U-frame. Our homeothermic warming pads are washable with water and detergent.

A heating plate and a probe are required for use with this unit.

References

Nguyen Chi, V., Müller, C., Wolfenstetter, T., Yanovsky, Y., Draguhn, A., Tort, A. B. L., & Brankačk, J. (2016). Hippocampal Respiration-Driven Rhythm Distinct from Theta Oscillations in Awake Mice. *Journal of Neuroscience*, 36(1).

Okun, M., Lak, A., Carandini, M., Harris, K. D., Buzsaki, G., Stevenson, I., ... Kaufman, M. (2016). Long Term Recordings with Immobile Silicon Probes in the Mouse Cortex. *PLOS ONE*, 11(3), e0151180. http://doi.org/10.1371/journal.pone.0151180

Gaylo, A., Overstreet, M. G., & Fowell, D. J. (2016). Imaging CD4 T Cell Interstitial Migration in the Inflamed Dermis. *Journal of Visualized Experiments*, (109), e53585–e53585. http://doi.org/10.3791/53585



	ORDERING INFORMATION
ATC2000	Animal Temperature Controller
RECOMME	NDED ACCESSORIES (select one sensor and one plate)
61800	Heating Plate with Built-in RTD Sensor, 15 x 25 cm
61830	Heating Plate with Built-in RTD Sensor, 15 x 10 cm
61840	Heating Plate with Built-in RTD Sensor*, 15 x 4 cm
	*for stereotaxic frame

OPTIONAL ACCESSORIES/REPLACEMENTS

RTD Rectal Temp Probe, 1.25 mm Shaft 2.5 mm Ball Tip
TC Rat Rectal Temp Probe, 1 mm Shaft, 3.2 mm Ball Tip
TC Mouse Rectal Temp Probe, 1 mm shaft, 1.6 mm Ball Tip
Silicone Pad for ATC2000 (10 x 15 cm)
29 ga 1 cm Needle Microprobe, 5-ft Cable
23 ga 3 cm Needle Microprobe, 5-ft Cable
Needle Microprobe, 5-ft Cable
Flexible Implantable Probe, 0.025-inch Diam., 3-ft Cable
Flexible Implantable Probe, 0.009-inch Diam., 3-ft Cable
Flexible Implantable Probe, 0.025-inch Diam., 3-ft Cable

Small Animal Anesthesia

Complete systems and all the accessories



Features

- Safe for surgical personnel, 90% below OSHA isoflurane limit
- Time efficient and cost-effective
- Virtually stress free for the animals
- Speedy recovery time

Benefits

- Compact and portable
- Easy to set up and use, simplifying the training of new staff and reducing the risk of human error
- Turnkey (all-in-one) system or plug-and-play system

Applications

· Small animal surgery

EZ-Anesthesia is the system of choice for anesthetizing small animals, and it comes with a variety of choices. Animals to be anesthetized are placed in the acrylic induction chamber, and the system delivers a precisely blended mixture of oxygen and isoflurane. An activated charcoal air filter canister at the top of the chamber releases safe, filtered air back into the room. A water-heated cage warmer or warming plate (ATC2000) is used to retain the animal body temperature while in the induction chamber. After the initial anesthetizing, the animal may be moved to the heated surgical water bed and positioned properly in the snugly fitted nose cone. A highly sensitive valve regulated by the animal's breathing works with the nose cone to ensure non-rebreathing efficiency. It allows safe anesthesia for up to several hours.

The breathing device also includes an air filter that releases safe, filtered air back into the room.

Oxygen and liquid isoflurane are not supplied. These are required to operate the system. Each system comes with all necessary components and connections for immediate use, including:

- · Oxygen regulator
- Vaporizer unit
- Induction chamber (standard **EZ-178**)
- Breathing device (standard EZ-103)
- Case of charcoal filters
- · Connecting tubing

Other available options include:

- Heated water beds
- Additional breathing devices
- Chambers in a variety of sizes
- Key fill isoflurane vaporizer is available separately or may be substituted in a complete system (call for details)

EZ-FF9000 Fixed Flow System

The Fixed Flow System is our most advanced system and provides preset, fixed flow rates with no need for adjustment. The system offers five gas outlets, each with an individual ON/OFF switch. The system ensures consistent, precise gas flow when connected to any pressurized gas source higher than 7 psi.

- Induction chamber output is fixed at 1 LPM
- Four breathing circuit outputs are fixed at 0.5 LPM

This unit features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

EZ-7000 Classic System

The Classic System has a user-friendly flow meter. The system offers five gas outlets that can supply a single induction chamber and four breather circuits simultaneously.

The system also features an oxygen flush button that purges the induction chamber with pure oxygen, thereby protecting personnel from exposure to anesthesia gas when opening the chamber. A handle on the vaporizer makes the unit easy to safely transport.

EZ-B800 Basic System

The Basic System is designed for use with a single animal. It utilizes one output directly from the vaporizer into a Y-splitter which creates a dual feed that can be directed to the induction chamber or the breather circuit.

This unit incorporates an oxygen flush system that purges the induction chamber with pure oxygen, thereby eliminating personnel exposure to anesthesia gas when opening the chamber.

	ORDERING INFORMATION
EZ-FF9000	Fixed Flow System
EZ-7000	Classic System
EZ-B800	Basic System

Anesthesia Accessories

Induction Chambers

These chambers incorporate a positive seal O-ring gasket for containment during use.

	ORDERING INFORMATION
EZ-177	Sure-Seal Mouse Chamber (5"L × 4.75"W × 4.38"H)
EZ-178	Sure Seal Mouse/Rat Chamber (9.75"L × 4.75"W × 4.38"H)
EZ-1785	Large Mouse/Rat Chamber (7"W × 11"D × 6"H inside)
EZ-179	Rabbit/Guinea Pig Chamber (18.75" × 12.75" × 12.75")

Ventilator Connection Kit

This kit is used to connect the **SAR-830** ventilators with EZ-Anesthesia systems. It includes all required components, pre-assembled for simple connection between the ventilator and the anesthesia system.

ORDERING INFORMATION

EZ-830 Ventilator Connection Kit

Connection Kit

This kit is used for connecting additional components to the EZ-Anesthesia Systems. It includes 6 ft of clear PVC tubing with a quick-disconnect fitting.

ORDERING INFORMATION

EZ-1130 Connection Kit

EZ Oxygen Regulators

Regulators are preset to 50 psi.

- **EZ-320** utilizes a CGA-540 connection for large "H" tanks.
- **EZ-330** utilizes a CGA-870 connection for small "E" tanks.

	ORDERING INFORMATION
EZ-320	Oxygen Regulator for large tanks
EZ-330	Oxygen Regulator for small tanks

Custom Hose Assemblies Available

Custom built hose assemblies for wall or ceiling outputs specific to facility needs: Chemetron, Ohmeda, Schrader or DISS.

Other Accessories



ORDERING INFORMATION
Versaflex Non-Rebreathing Unit
Microflex Non-Rebreathing Unit
Rat Stereotaxic Non-Rebreathing Unit
Multi-Animal Non-Rebreathing Unit
Mouse/Rat Thin-Line Heated Waterbed

Mobile Workstations


Two mobile workstations, constructed of heavy-duty stainless steel with locking casters, integrate all your EZ-Anesthesia components into one portable unit. Open side shelves accommodate 20 lb. cylinders, and convenient 2" port holes allow for easy rigging of gas and electrical lines. Below the work surface of each mobile workstation is an open shelf and a locking cabinet. **EZ-E25000** provides a 42"x24" work surface and holds up to four cylinders. **EZ-E27000** has a 22"x21" work surface and holds up to two cylinders. These systems are easy to set up and provide maximum flexibility and mobility.

	ORDERING INFORMATION
EZ-E25000	Mobile Workstation, 42" x 24" Top
EZ-E27000	Mobile Workstation, 22" x 21" Top



gas source. A foam lid gasket ensures a good seal on the cage. Multiple lids may be used to treat several cages at once.

	ORDERING INFORMATION
EZ-20027	Small Lid (13" x 9") Fits old-style mouse cages
EZ-20028	Small Lid (16" x 10") Fits new-style mouse cages
EZ-20030	Square Lid (13" x 13") Fits Thorn cages
EZ-20032	Medium Lid (20.5" x 11") Fits rat cages
EZ-20034	Large Lid (23" x 16.5") Fits guinea pig cages
EZ-20029	Lid Storage Bracket (wall-mounted, holds up to four lids)
EZ-20027G	Small Lid Gasket (13" x 9")
EZ-20028G	Small Lid Gasket (16" x 9")
EZ-20030G	Square Lid Gasket (13" x 13")
EZ-20032G	Medium Lid Gasket (20.5" x 11")
EZ-20034G	Large Lid Gasket (23" x 16.5")



Blood Pressure Monitor and Transducer

Audible monitor with variable pitch provides feedback

Features

- Variable amplification
- Displays systolic, diastolic or average blood pressure
- Audio monitor
- · BNC output for data logging

Benefits

• BLPR2 can be used for the direct arterial and venous pressure measurement in animal blood vessels

Applications

· Invasive monitoring of animal arterial or venous blood

BP1 accepts WPI's BLPR2 blood pressure transducer (below) as well as other blood pressure transducers. An audio monitor provides a signal with variable pitch and amplitude, allowing you to hear changes in blood pressure. Digital LCD display provides average or peak signal values from 0 to 1999 mV. With an optional pressure gauge (not provided — see **PM015D**, page 116), you may calibrate the display to read mmHg. Recorder output connector allows direct connection to a pen recorder, oscilloscope or computer via a data acquisition system.

Supplied sterile, **BLPR2** is accurate, linear and stable with temperature. May be sterilized cold with Rapicide OPA or a similar bactericide.

BLPR2 is equipped with a 12-ft cable and connector compatible with WPI's four-channel signal conditioning unit, TBM4M Transbridge, and the single-channel BP1 blood pressure monitor. The cable has a moisture-resistant locking connector. A continuous, uniform lumen eliminates places for bubbles to form and lodge. The clear fluid path is easy to inspect. Easy to mount — slotted transducer body accepts Velcro strap.

To facilitate setup and operation, a four-way stopcock that allows easy filling, flushing and zeroing of the transducer is included. Typically, the stopcock is located between the transducer and the animal catheter where it can be used to quickly zero, flush or de-bubble the transducer.



Stopcock 14036 included with BLPR2

Micro Cannula KZ1101 0.4 mm O.D., 0.2 mm I.D. tubing

- Autoclavable
- Biocompatible perfluorocarbon tubing material

See description on page 221.

CE Pressure Monitor BP-1 SYS-BP1

BP1 SPECIFICATIONS

AMPLIFICATION x1, x10, x100, variable (x5 to x1000)

OUTPUT VOLTAGE SWING ± 5 V MAXIMUM OUTPUT CURRENT 2 mA

INPUT IMPEDANCE, EACH INPUT 100 kΩ || 0.01 μF

TRANSDUCER APPLIED VOLTS 10 V nominal, varies with load.

25 mA, maximum

POWER 95-135 V or 220-240 V, 50/60 Hz 8.5 x 5.12 x 10 in. (21.6 x 13 x 25.44 cm) DIMENSIONS

SHIPPING WEIGHT 11 lb. (5 kg)

BLPR2 SPECIFICATIONS

WORKING PRESSURE -50 to +300 mm Hg -400 to +4000 mm Hg **OVERPRESSURE EXCITATION VOLTAGE** 1-10 VDC or RMS to 5 kHz

SENSITIVITY 5 uV/V/mm Hg DYNAMIC RESPONSE 100 Hz

EIGHT HOUR DRIFT 1 mm Hg after 5 minute warm-up

Total combined effects of Sensitivity, Linearity, MAXIMUM FRROR

Hysteresis (at 25°C and 5 μ V/V/mm Hg) do not exceed ±2% or 1 mm Hg, whichever is greater.

SHIPPING WEIGHT

ORDERING INFORMATION

Pressure Monitor (transducer & cable not included) SYS-BP1 Specify line voltage

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

BLPR2	Blood Pressure Transducer & Cable
BPCABLE2	Cable (12 ft) with DIN connector for BLPR2
503067	BLPR2 Transducer without cable
13024	Single Rack Mount Kit
13025	Dual Rack Mount Kit
3491	Extension Cable, 5 ft
500184	BNC-to-BNC Cable, 10 ft
14036-15	4-Way Luer Stopcock, Blue Tint (package of 15)
KZ1101	Micro Cannula, 3 inch
NOTE DIDDO	

NOTE: BLPR2 is intended for animal research only and may not be used for human blood pressure measurement.

Blood Pressure Measurement for Rodents

Unique environmental control system and single tail cuff allows tail movement, minimizing stress in the animal subjects



Features

- Test up to 200 animals at a time
- Sensor is MRI compatible
- Quick and accurate blood pressure measurement at temperatures as low as 32°C
- Highly sensitive photoelectric sensor for blood pressure detection
- Monitor, record, store or export real time systolic, diastolic, mean and heart rate

Benefits

- Automatic evaluation of data
- Touch screen or software control
- Data can be collected and stored on and displayed from a USB flash drive, which is included

Applications

• Non-invasive blood pressure monitoring of rodents

This revolutionary design brings non-invasive blood pressure testing to a new level — a true turn-key system for accurate, consistent blood pressure measurement on mice, rats or any other laboratory animal test subject

It is a compact, simple yet versatile system that can test from one to 200 animals at a time with independent control of each channel. Simple daisy-chaining allows expansion of up to 200 independently controlled systems.

All components are built into one small unit — controls inflation of tail cuff, warming environment with whisper-quiet fans — providing an ideal system for teaching facilities and for the pharmaceutical industry when high throughput is a must.

Single animal systems are controlled from the touch screen, which allows keying in all necessary test setups. Touch screen control allows ease of operation, supplying automatic evaluation of test results — systolic, diastolic, mean and heart rate.

Data is collected, stored, displayed and can be transferred to the supplied memory stick. The USB interface allows for software control of multichannel systems. For single animal systems, built-in software lets you view and export data. Reports are in an Excel-style format and may be easily exported.

No computer is required. However, the analog output may be interfaced with your own data acquisition software.

System is easily cleaned. Removable trays are included with each system.

Tail Cuff Sensors have the ability to detect mouse blood pressure or rat blood pressure via the tail photoelectrically. The photoelectric sensor, light source and sensor are built into a unit which attaches to a restrainer with nylon hand screws. Sensors are nonmagnetic. The ambient temperature that a reading can be taken starts at 30°C. The maximum temperature needed for a test is usually 32°C for rats and 34°C for mice. The systolic readings are within five mmHg of Cannulation and Telemetry.

In addition to the standard one-year warranty on the system, tail cuff sensors have a lifetime warranty.



ORDERING INFORMATION

II-MRBP-M	Mouse Blood Pressure System
II-MRBP-R	Rat Blood Pressure System
Includes	one tail cuff cancer and one restrainer (select when or

Includes one tail cuff sensor and one restrainer (select when ordering).

Call for pricing on multiple channels.

OPTIONAL ACCESSORIES

II-B64-1	Tail Cuff Sensor, 650-800 grams, ID 1 in.
II-B63-075	Tail Cuff Sensor, 500-650 grams, ID ¾ in.
II-B60-044	Tail Cuff Sensor, 150-380 grams, ID ¾6 in.
II-B60-038	Tail Cuff Sensor, 70-160 grams, ID ¾ in.
II-B60-025	Tail Cuff Sensor, 20-80 grams, ID ¼ in.
II-78	Rodent Restrainer, 750 to 1000 grams, ID 4" (105 mm), 13" Long
II-79	Rodent Restrainer, 600 to 800 grams, ID 3½" (88 mm), 13" Long
II-80	Rodent Restrainer, 400 to 600 grams, ID 3" (75 mm), 12" Long
II-81	Rodent Restrainer, 250 to 400 grams, ID 2½" (63 mm), 9" Long
II-82	Rodent Restrainer, 180 to 270 grams, ID 2" (50 mm), 8" Long
II-83	Rodent Restrainer, 70 to 170 grams, ID 1½" (38 mm), 6" Long
II-84	Rodent Restrainer, adult mice and neonate rats, 20 to 45
	grams, ID 1" (25 mm), 4" Long
II-84XL	Rodent Restrainer, 45 to 70 grams, ID 11/4" (30 mm), 5" Long

Microprobe Thermometers

Instrument of choice for biological and laboratory use





Features

- Super accuracy
- Fast response
- · Analog output signal
- Multiple inputs
- Differential temperature measurement

Benefits

- Portable
- A variety of probes available

Applications

- · Temperature monitoring
- · Cryogenic measure for blood banking
- Skin temperature measurement in exercise experiments

A Microprobe Thermometer is the instrument of choice for biological and laboratory temperature measurements. These thermometers are very versatile, providing fast response, high accuracy and stability with digital display and analog signal for connection to a computer or recorder. With the wide selection of probes, the instruments can be used in almost any application.

BAT-12: This thermometer has a sealed construction making it water, dust and fume resistant. The **BAT-12** has a single microprobe input and a single range with the same high accuracy as the **BAT-10**. Comes complete with carrying case.

The thermometers can be used with any "Type T" thermocouple. Select a temperature microprobe on the following page for your specific application.

BAT-10: This is the most versatile thermometer available. The instrument has a wide temperature range and fast response with most microprobes. The **BAT-10** accuracy is NIST traceable and in each of the two temperature ranges, the accuracy is the same as the resolution. There are three microprobe inputs, 1 and 2 can be selected as separate inputs while 2 and 3 will read the differential temperature measurement between the two. The instrument has automatic warnings for low battery or faulty probes on the digital display. The linearized analog output (LOP) signal allows ease of connection to a data acquisition system or recorder.

	BAT-10	BAT-12
TEMPERATURE RANGE & RESOLUTION	-200°C to +400°C, 1°C resolution -100°C to +199.9°C, 0.1°C resolution	-100°C to +199.9°C, 0.1°C resolution
DIFFERENTIAL TEMP. RANGE	-19.99°C to +19.99°C Linearization centered at 40°C 0.01°C resolution	N/A
ACCLIRACY		

MICROPROBE THERMOMETER SPECIFICATIONS

1° Range 1°C ±1 least significant digit 0.1°C ± 1 digit between 0-50°C 0.1° Range 0.1°C ± 1 least significant digit 0.1% ± 1 digit over full range Diff. Range 0.01°C ± 1 least significant digit REPEATABILITY ± 1 least significant digit Follows NIST thermocouple tables CALIBRATION CONFORMITY Conforms to NIST tables within 1 digit DISPLAY 3½ Digit LCD 3½ Digit LCD Miniature, quick disconnect, copper-Miniature, quick disconnect, copper-INPUT SOCKET constantan constantan Non-linearized set at 1.6 V, ANALOG OUTPUT corresponding to temperature of ≈ 10 mV per degree C 401°C BAT-10: 4 alkaline "C" cells (life: 1000 hr) BAT-12: 9V cell POWER SUPPLY/BATTERIES BAT-10R: 4 Ni-Cad "C" cells (recharge- BAT-12R: 9V Ni-Cad with charger able unit) SENSORS Three Type T thermocouple inputs One Type T thermocouple input AMBIENT OPERATING Auto-compensated to 0.1°C from 15-45°C RANGE 0°C to 50°C

ORDERING INFORMATION

21.6 x 22.9 x 8.9 cm (8.5 x 9 x 3.5 in.) 12.7 x 6.4 x 15.2 cm (5 x 2.5 x 6 in.)

1.6 kg (3.5 lb), including carrying case 1 kg (2 lb), including carrying case

BAT-10R/LOP Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)

BAT-10R/LOP220 Multiple Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)

BAT-12R Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 110 VAC adapter (microprobes ordered separately)

BAT-12R-220 Single Input Type T Thermocouple Thermometer, rechargeable NiCad batteries and 220 VAC adapter (microprobes ordered separately)

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

EXT-6	Probe Extension Lead, 180 cm long
501608	Tripod Stand for BAT-12

DIMENSIONS

WEIGHT

Temperature Probes

A variety of probes for any application

Features

- Flexible Teflon microprobes are used for implantation in tissue, in spectrophotometer cuvettes, rectally in neonatal mice, in water baths, PCR thermal cyclers, etc.
- Monitors animal rectal temperatures during surgical procedures and pyrogen testing


Benefits

 Thermocouple wire used meets NIST standards (certification is available for an additional cost)

Applications

- · Rectal temperature monitoring
- Instant read in tissue, semi-solids, small specimen and other materials

When precise temperature measurements are required, WPI can provide you with a very accurate monitor and thermocouple microprobes. WPI monitors have both resolution and accuracy of 0.1°C in the 0-50°C range and are traceable to NIST standards, whereas, other competitive electronic thermometers have an accuracy that is usually to 0.5°C or worse. Furthermore, all our type T clinical probes are guaranteed accurate to 0.1°C, due to our stringent wire standards. These are five times more accurate than competitive probes made with regular "Special Limits"



NEEDLE MICROPROBES ORDERING INFORMATION

Fast-response needle probes for instant readings in tissue, semi-solids, liquids, very small specimens, powders and materials. Needle tip is sealed to ensure only stainless steel contacts specimen. Maximum autoclavable temperature is 135°C.

			Time		Max.	Lead	
Probe Type	Size	Style	Constant	Isolated	Temp.	Length	Description
MT-29/1HT	29 ga / 1 cm	Α	0.125 sec	No	200°c	5 ft	29 gauge approximately 0.013 in.
MT-29/2HT	29 ga / 2 cm	Α	0.125 sec	No	200°c	5 ft	-
MT-29/3HT	29 ga / 3 cm	Α	0.125 sec	No	200°c	5 ft	
MT-29/5HT	29 ga / 5 cm	Α	0.125 sec	No	200°c	5 ft	-
MT-26/2HT	26 ga / 2 cm	Α	0.1 sec	No	200°c	5 ft	26 gauge approximately 0.018 in.
MT-26/4HT	26 ga / 4 cm	Α	0.1 sec	No	200°c	5 ft	-
MT-26/6HT	26 ga / 6 cm	Α	0.1 sec	No	200°c	5 ft	
MT-23/3HT	23 ga / 3 cm	Α	0.15 sec	No	200°c	5 ft	23 gauge approximately 0.125 in.
MT-23/5HT	23 ga / 5 cm	Α	0.15 sec	No	200°c	5 ft	
MT-23/8HT	23 ga / 8 cm	Α	0.15 sec	No	200°c	5 ft	-
MT-4	29 ga / 1 cm	А	0.025 sec	No	200°c	5 ft	Similar to MT-29/1 but has a blunt tip. Good for instant skin and surface temperatures, liquids
MT-D	_	C	0.025 sec	No	200°c	5 ft	Fast response surface probe (stainless steel for locating inflammation, arteries, etc. Also for dental use.

FLEXIBLE IMPLANTABLE PROBES ORDERING INFORMATION


Designed for high accuracy on extremely small specimens such as insects, seeds, etc. Maximum insertion depth 1/8". Totally sheathed in chemical resistant Teflon.


Probe Type	Sensor Lead Diameter	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
IT-14	0.050" dia	D	0.3 sec	Yes	150°c	3 ft	-
IT-18	0.025" dia	D	0.1 sec	Yes	150°c	3 ft	-
IT-18EXLONG	0.025" dia.	D	-	Yes	150°c	5 ft	-
IT-21	0.016" dia	D	0.08 sec	Yes	150°c	1 ft	-
IT-23	0.009" dia	E	0.005 sec	Yes	150°c	3 ft	For ultra fast measurements and for use on micro- size specimens. Tissue implantable with 23ga. Needle (supplied). Rather fragile. Teflon coated.
IT-1E	0.025" dia	F	0.005 sec	Yes	150°c	3 ft	As IT-18 sensor except bead exposed. Combines ultra-fast reponse of IT-23 with sheath strength of IT-18.


RECTAL PROBES ORDERING INFORMATION							
Probe Type	Sensor Lead Diameter	Style	Time Constant	Isolated	Max. Temp.	Lead Length	Description
RET-2	-	G	0.8 sec	No	125°c	5 ft	Rectal probe for rats typically for fast intermittent mea- surements. Smooth ball tip (0.125 in. dia.) with 1" long (0.59 in. dia) stainless steel shaft.
RET-3	-	G	0.5 sec	No	125°c	5 ft	Rectal probe for mice similar to RET-2 except tip diam. 1.6mm (0.063") and shaft $\frac{3}{4}$ " long (0.028" diam.)

	GENERAL PURPOSE ORDERING INFORMATION								
Probe Type	Probe Type Sensor Lead Diameter Style Constant Isolated Temp. Length Description								
HT-1		Н	0.5 sec	No	400°c	5 ft	"Workhorse" probe for liquids, gases, semi-solids. Plastic handle with straight stainless steel shaft. Not good for surface temperatures.		
HT-2		Н	0.5 sec	No	400°c	5 ft	Like HT-1 except shaft length is 9 in.		

Due to the fragility of these probes, warranty is limited: Defective probes may be returned within 14 days of receipt.









These probes may also be used with **ATC2000** Animal Temperature Contoller.

See page 127.

Neuroscience Cannulas

For in vivo investigation of rodents

Features

- · Exceptional quality
- · Best prices
- Rapid order response

Benefits



Cannula assembly: fixing screw, internal cannula with attached tubing, and guide cannula.

• Beveled inside and out and polished to remove burrs and ensure that and ensures smooth operation.

· Exceptional quality at the best prices

Applications

- · Drug administration
- Optogenetics

the inside diameter is perfectly cylindrical. This limits trauma to tissue

Internal cannula secured with fixing screw.

These cannulas for neuroscience study and pre-clinical research include an entire range of cannula options, including the Guide Cannula, the Internal Cannula and the Dummy Cannula (cap).

Guide Cannulas



The Guide Cannula is a surgical grade, stainless steel tube that is implanted into a rodent's skull and cemented into place using dental cement and screws. It guides the Internal Cannula to the specific injection site.


Internal Cannula

↓ 0.7mm	Gauge 22	OD	ID
V 0.7111111	22	0.41	0.26
	24	0.36	0.21
	26	0.30	0.16

Insert the Internal Cannula into the Guide Cannula to sample or inject

Dummy Cannula

The Dummy Cannula has a stainless steel wire core, and it is placed in the Guide Cannula when the Internal Cannula is removed. It seals the



it while grooming. Note that the DUM26 cap is used for the 24 ga cannulas also.

Ordering

Order the Guide Cannulas based on the gauge and length from the base and the Internal or Dummy Cannulas based on the length of the Guide Cannula and the projection from the Guide Cannula tip.

Understanding Part Numbers

onderstanding rait is	GC22-37
Cannula Type	
GC -Guide Cannula	Cannula Type 🥕 🁌 🧳
INC-Internal Cannula	Gauge of Guide Cannula
DUMC -Dummy Cannula	Length of Cannula or Projection

Gauge-Choose the gauge of the Guide Cannula to be used, even if you are ordering an Internal or Dummy Cannula. Choices include 22, 24 and 26 gauge.


Length–Guide Cannulas can be ordered in lengths from 1.0–9.9 mm. Length can be specified to 0.1 mm, with a tolerance of ±0.07mm. For example, a GC24-60 is a 24 gauge cannula that is 6.0 mm long.

For other cannulas, the length is determined by the desired projection length and the guide cannula length. The projection can extend beyond the tip of the Guide Cannula up to 1.0 mm.

DUMC Length = GC Length + Projection

INC Length = GC Length + Projection

For example, for a Dummy Cannula flush with the end of a 26 ga Guide Cannula that is 6.0 mm long, order a **DUMC26-60**. For an internal cannula that projects 0.5 mm beyond the 6.0 mm Guide Cannula (6.5 mm), order an INC26-65.



* DUMC Length=GC Lenth + Projection. If the cap is screwed on too tightly, the projection will be longer than expected.
** INC Length=GC Lenth + Projection. Internal cannula mounts flush and does not screw into place.

	ORDERING INFORMATION
Cannulas	
GC22-X*	Guide Cannula, 22 ga, X.0 mm
INC22-X*	Internal Cannula, 22 ga, X.0 mm
DUMC22-X*	Dummy Cannula, 22 ga, X.0 mm
GC24-X*	Guide Cannula, 24 ga, X.0 mm
INC24-X*	Internal Cannula, 24 ga, X.0 mm
GC26-X*	Guide Cannula, 26 ga, X.0 mm
INC26-X*	Internal Cannula, 26 ga, X.0 mm
DUMC26-X*	Dummy Cannula, 26 ga, X.0 mm
* V indicator th	as length of the cannula or projection. Can Understanding Part

 st X indicates the length of the cannula or projection. See Understanding Part

Flexible PE Tubing

504278	0.25 mm ID, 0.5 mm OD, 1 m long		
504279	0.42 mm ID, 0.85 mm OD, 1 m long		
504280	0.6 mm ID, 1.1 mm OD, 1 m long		
Miscellaneous Accessories			
504281	Fixing Scraw-connect INC & GC		

Flexible PE tubing is recommended to connect to Internal Cannula (INC) only, NOT Guide Cannula (GC) or Dummy Cap (DUMC).

504278 matches INC26-XX 504279 matches INC24-XX 504280 matches INC22-XX

Electronic von Frey Anesthesiometer

Assessing mechanical allodynia with 15 SuperTips

Features

- LCD readout (floating or last maximum/minimum)
- Rigid tips up to 800 g
- Supertips™ up to 65 g
- 1,000 g probe available
- · Optional analog output cable for chart recorder
- Pipette tips can be customized to any specification
- Microprocessor electronics 0.1 g plug-in probes

Benefits

- Plug up to three probes into a single unit
- Independent from temperature

Applications

• Small animal analgesia testing

To assess mechanical allodynia, which is a painful response to a light touch or pressure from a stimulus that is not normally painful, the Electronic von Frey Anesthesiometer was developed. The Electronic von Frey meter uses one of 15 different flexible von Frey hairs called "SuperTips™" (or rigid tips up to 800 grams). Each hair, regardless of model chosen, is exactly 0.8 mm in diameter. This uniformity of design eliminates false readings and allows for comparison of test results. The Electronic von Frey can be used with chart recorders and analog/digital converters, and it never needs calibration. This system includes either a 90, 800 or 1,000 gram probe. Mesh stands are available in a variety of sizes for large group studies.

1 13 Superrips	
383	
180	II-2391

	ORDERING INFORMATION
II-2390	Electronic von Frey Anesthesiometer, Rigid Tips, 90 g Range
II-2391	Electronic von Frey Anesthesiometer, Rigid Tips, 800 g Range
II-2392	Electronic von Frey Anesthesiometer, Rigid & 15 Supertips, 90 g Range
II-2393	Electronic von Frey Anesthesiometer, Rigid & 15 Supertips, 800 g Range
II-23931	Electronic von Frey Anesthesiometer, Custom Rigid Tips, 1000 g Range
II-2394	von Frey Probe, 90 g Range
II-2395	von Frey Probe, 800 g Range
II-2396	von Frey Probe, 1000 g Range
II-2397	MRI Probe Option (add to price of probe above)
II-2400	Analog Output Cable

Quattro

Four test systems



This special package offers four tests, including Electronic von Frey, Plethysmometer, Randall Selitto and the Grip Strength Meter. You get all four test modules and the electronic controller that is interchangeable with all four systems. The electronic controller has up to three inputs, so you can perform up to four unique tests with only one electronic system. If you prefer, you may build your system as you grow. Because of the modular design of these four systems, you need to order only one complete system. Then, the modules for the other three tests, which integrate into the system, can be purchased separately, as needed. The stand and sling for the Randall Selitto test are sold separately.

ORDERING INFORMATION

II-2889 Quattro 4-in-1 System

Trio

Three test systems



Get three test systems in one package with the Trio, featuring the Electronic von Frey, Plethysmometer and Randall Selitto Meters. Just like the Quattro package, the modular design allows these three test systems to communicate with the same electronic controller. The stand and sling for the Randall Selitto test are sold separately.

ORDERING INFORMATION

II-2888 Trio 3-in-1 System

WORLD PRECISION INSTRUMENTS

Plethysmometer (Paw Volume Meter)

Test effectiveness of anti-inflammatory agents

Features

- · No wetting solution needed
- One calibration/year
- · Battery-operated or line powered controller
- One-year warranty

Benefits

- A sensor in the water notes a pressure change when the paw is
- Pressure is calibrated in 0.1 milliliters and displays on the batterypowered monitor

Applications

• Measure the effectiveness of anti-inflammatory agents and agents to reduce edemic conditions

effectiveness of anti-inflammatory agents and agents to reduce edemic conditions with the Plethysmometer. Simply insert the rat or mouse paw into water. The acrylic stand offers



full visibility of the subject throughout the testing.

ORDERING INFORMATION

II-520MR

Paw Volume Meter for Mouse & Rat

Randall Selitto Paw Pressure Meter

Get instantaneous, live readings **Features**

- Visible force limit indicator
- Portable electronic display
- Battery-operated or line powered controller
- · No calibration required

Benefits

- Hands-free operation with foot switch
- 3-oz. hand held probe with accuracy of 0.5%

Applications

• Analgesic drug studies

The Randall Selitto Paw Pressure Meter for analgesia testing is digitally controlled. Use the handheld instrument to apply force to an animal's extremity and get instantaneous, live readings. You can even view the last maximum force



800 g pressure applicator

applied during the test. The new limit indicator lets you select the maximum force limit, and then indicates with a warning light when the system reaches that limit. This unit comes with an acrylic stand to allow for easy viewing of the display. Stand and sling are sold separately.

ORDERING INFORMATION

II-2500

Randall Selitto Paw Pressure Meter

Grip Strength Meter

Measure muscle hyperalgesia in rodents **Features**


- · Mouse and rat wire mesh grid plates included
- Maximum force range is 2,000 g (10% over range allowed). Higher ranges are available by special order.
- 1 gram increment readings
- Suction feet on the heavy, anodized base plate resist even large pulling forces

Benefits

- Battery-operated or line powered controller
- · One-year warranty

Applications

· Measure muscle hyperalgesia in rats and mice



Measure muscle hyperalgesia in rats and mice with the Grip Strength Meter, which gauges the forelimb grip force using a digital force transducer. Simply hold the animal by the tail and gently dangle it over the wire mesh plate until the animal grasps the plate with its forepaws. The force transducer, connected with the wire mesh plate, measures the strength of the animal at the time of the test. The battery-operated, electronic control device calculates the average of three measurements, and it holds the last maximum force in a "peak and hold" type readout until you reset it.

ORDERING INFORMATION

Grip Strength Meter for Mouse & Rat

Incremental Hot Cold Plate Analgesia Meter

Latency and threshold-based nocicetption

Features

- Heat or cool, from 0-70°C
- Ramping temperatures for threshold & latency results
- Rapid increase or decrease in temperature
- Printout of data
- Includes clear animal enclosure
- Plate size 4" x 8"
- Two-year warranty

Benefits

- Precise programmable digital control
- Temperature stability is 0.1 °C

Applications

- · Analgesia research
- Latency and threshold-based nocicetption

This safe, humane device for rats and mice is used for latency and threshold-based nocicetption, ramping temperatures for 0-70°C. Because this hot cold plate is incremental, it measures latencies of much more than just the strong narcotic agents, broadening dramatically the range of analgesia research with devices of this type. Microprocessor-controlled, the Incremental Hot Cold Plate can heat or cool in increments of 0.1°C, at a rate of 1-10°C per minute. With



uniform heating and cooling and upper/lower cut-off limits, this device is predictable and safe. It can also function as a constant temperature plate with great stability (0.1°C). As soon as a reaction is observed from the chosen paw, the unit reverses to the standby temperature.

ORDERING INFORMATION

II-PE34 Incremental Hot Cold Plate Analgesia Meter

Incapacitance Meter for Mouse & Rat

Test pain and inflammation in small animal hind limbs

Features

- Precise programmable digital control
- Start, stop and reset test from controller's front panel
- 180-270 gram holder included (other sizes available)
- · Reaction detected automatically
- Manual override of all timer functions
- · Alphanumeric readout
- 5 to 999 seconds test period
- Alphanumeric readout
- All functions and parameters entered via keypad
- Two-year warranty

Benefits

 Overcomes the drawback of placing unnecessary stress on the animals

Applications

• Test pain and inflammation in the hind limbs of mice, rats or birds

The Incapacitance Meter uses a technique called dual channel weight averaging, which tests both hind limbs. This gives you a clean, stress-free correlation of the paw pressure test. Conduct control and testing of the animal at the same time. Place the animal in



the holder with its hind limbs resting on the two weight-averaging platform pads. The controller records the average weight (grams) over the test period as the animal shifts its weight from each pad.

ORDERING INFORMATION

II-600MR

Incapacitance Meter for Mouse & Rat

Hot Plate Analgesia Meter

Latency testing in rodents

Features

- · Includes both mouse and rat enclosure
- Temperature indicated in 0.1 °C increments
- Holding accuracy is ±0.1°C
- · Digitally controlled
- Two-year warranty

Benefits

 Plate has consistent temperature throughout a procedure, ensuring accurate testing

Applications

· Latency testing in rats and mice

To use the Hot Plate Analgesia meters, simply place the animals on the black anodized aluminum plate (11" x 10.5" x 0.75", 275 x 263 x 15 mm) and set the plate's surface temperature



to the desired set point (up to 75°C). The plate maintains a consistent temperature throughout the test.

ORDERING INFORMATION

Hot Plate Analgesia Meter for Mouse & Rat

Plantar Test Apparatus/Tail Flick Test Analgesia

II-39

Meter to test narcotics and strong non-narcotic drugs

Features

- Includes three acrylic animal enclosures that each hold two rats or four mice
- Precise programmable digital control
- User-defined humane cutoff feature
- Adjustable beam intensity in 1% increments up to 250°C
- Reaction is detected automatically
- Alphanumeric readout
- Manual override of all timer functions
- All functions and parameters entered via keypad
- · Heated glass option
- Tail temperature monitor option (for use with the Tail Flick Meter)

Benefits

- Plantar test and tail flick test applications are combined in one
- Stimulate other body parts by adjusting the height of the glass

Applications

• Test the properties of narcotics on unrestrained mice/rats in Plantar

This unit, which is designed for testing narcotics and strong non-narcotic drugs, offers both Plantar (Hargreaves Method) and Tail Flick testing with a single unit. Either testing system is also available individually. In plantar mode, the visible light/heat source is directed at the paw or other desired body part, and in tail flick mode it is directed at the subjects' tails. Test up to 12 mice or 6 rats simultaneously. If desired, other animals like cats and rabbits may also be used. Tests are simple to setup. The focused, radiant heat/light source creates a 4 x 6 mm intense



spot. Because the light is visible, you know when the test starts and ends. The equipment is silent (no whining or clicking sounds) to avoid triggering an automatic response in conditioned animals. You can set a humane cutoff timer that automatically shuts off the heat if no response is observed during the designated time frame.

When an animal is placed on cold glass, its reaction time may be slower. This unique system offers a heated glass option that prevents the glass enclosure from acting as a heat sink, giving a more accurate reading. An optional tail temperature monitor can also be selected for use with the Tail Flick Meter. This option actually preheats the tail before experimentation. Once the preset tail temperature is reached, the test and timer automatically begin. A glass stand is also available in two sizes for large group studies.

	ORDERING INFORMATION
II-336T	Combination Plantar/Tail Flick Meter (non-heated glass and tail temperature for mouse and rat)
II-336TG	Combination Plantar/Tail Flick Meter (tail temperature and heated glass for mouse and rat)
II-390	Plantar Test Analgesia Meter (non-heated glass for mouse and rat)
II-390G	Plantar Test Analgesia Meter (heated glass for mouse and rat)

Adhesives

	Α	DHESIVES APPLICA	TION GUIDE	
WPI Part #	Description	Curing Time	Useful Applications and Characteristics	
Epoxies	Form strong bonding. Used in wire bonding applications.			
4898	Silver Filled Conductive Epoxy	12 hr @ 50° C – 5 min @ 150° C	Connecting conductors that can't be soldered. Constructing or connecting Ag/AgCl pellets.	
7335	Carbon Filled Conductive Epoxy	48 hr @ 25° C – 5 min @ 150° C	Constructing carbon electrodes	
4886	High Performance Structural Epoxy	12 hr @ 25° C	Forms a strong and slightly flexible bond on plastic, metal, & glass. Bonds some low surface.	
Hot melt (EVA)	Easy to use for bonding, need	ls large gap filling		
13316	Mini Glue Gun with glue sticks	As soon as it cools down	Bonds wood, glass, metals, and many plastics	
Silicone Adhesiv	es/Sealants/Primers	Good moisture resistar	nt and elastic. Low toxic.	
1571	Room Temperature Vulcanizing (RTV) adhesive. Acyloxy/moisture cure system. Acetic acid is cure by-product.	24 hr @ 25° C	Has the best adhesion property in this silicone family. Will bond to many materials.	
7128	RTV sealant. Alkoxy/moisture cure system. Methanol is cure by-product.	72 hr @ 25° C	Good for bonding or sealing electronics circuits (metal)	
SYLG184	Sylgard. Two-part vinyl/platinum cure sealant. Hydrogen is cure by-product. Very low toxicity.	24 hr @ 25° C – 15 min @ 150° C	Coating patch clamp electrodes, cell culture dish, making dissection pads	
KWIK-SIL	Two-part adhesive. Vinyl/platinum system. Hydrogen is cure by-product. Very low toxicity.	< 5 min @ 25° C	Live tissue and nerve studies. Medium strength adhesion.	
KWIK-CAST	Two-part sealant. Vinyl/platinum cure system. Hydrogen is cure by-product. Very low toxicity.	< 5 min @ 25° C	Sealant for live tissues. Embedding peripheral nerves with electrodes.	
6820	Primer for silicone	N/A	Enhances adhesion of silicone adhesives for difficult to bond plastic surfaces	
Cyanoacrylate	Forms an instantaneous bond	ding.		
7341	Ethyl cyanoacrylate, low viscosity 90-120 cps	<10 seconds	Mounting rat/mouse brain slices. Ideal for relatively small gaps and smooth surfaces. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.	
7342	Ethyl cyanoacrylate, high viscosity 1100-1600 cps	<30 seconds	Use in brain slice experiments. Ideal for larger gaps, allows slightly longer bonding time. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.	
VETBOND	Butyl cyanoacrylate, low toxicity	<10 seconds	Bonds tissues, alternative to suture, helps small wound healing. Antimicrobial effect. Used in forensic science.	
503763	Octyl cyanoacrylate, low toxicity	<15 seconds	Suitable for surface wound bonding, protection, holding a sensor or other device on the tissue.	

Kwik-Gard[™]

Kwik-Gard™ is specially packaged Sylgard 184 silicone for quicker and easier application, eliminating the messy procedure of preparing the mixture before application. Its special cartridge controls the precise mixing ratio to ensure proper curing. The disposable tip mixes resin



and hardener as they are dispensed. Since no air is introduced during mixing, the resin does not need degassing for most applications. The mixed silicone is applied directly to the site, reducing preparation time and material waste.

Each Kwik-Gard cartridge contains 37 mL of resin and hardener. The dispensing tip has a dead volume of 0.75 mL.

ORDERING INFORMATION			
KWIKGARD	Kwik-Gard™ Start-up Kit (incl. dispenser, 1 cartridge, 5 tips)		
KWIKGLUE	Kwik-Gard™ Refill (2 cartridges, 10 dispensing tips)		
KWIKMIX	Dispensing Tips (pack of 10)		
KWIKGUN	Kwik-Gard™ Dispenser		

Scotch-Weld 2216 Structural Epoxy

Scotch-Weld 2216
remains the best epoxy
for bonding plastic, often
used as the benchmark
for testing the binding
strength of other
adhesives. The slightly
rubbery texture also makes it

less easy to break off. It is the only epoxy known that can bond PEEK.

Color: gray Cures at room temperature. Shipping weight: 1 lb. (0.5 kg)

	ORDERING INFORMATION
4886	Scotch-Weld 2216 (2 oz.)

2216 B/A

4886

Low Toxicity Adhesive

Ideal for neuroscience applications, nerve studies and more

600022

Features

- Bio-compatible adhesive for live tissue and nerve studies
- · Pre-mixing tips simplify use
- Medium strength adhesion
- Low toxicity
- · Rapid curing silicone adhesive, cure on contact
- Cures without producing heat
- Includes 10 mixing tips
- Volume discounts save up to 15%!



KWIK-SIL

Benefits

- Low toxicity
- Rapid cure time

Applications

- Neuroscience and nerve studies
- Biomedical applications

Kwik-Sil and Kwik-Cast silicones have very low toxicity before, during and after curing. The by-product of curing is a small amount of hydrogen gas, which is much less toxic to cells than acetic acid or alcohol from traditional RTV silicone sytems.

Kwik-Sil and **Kwik-Cast** curing speed is hundreds of times faster than traditional RTV silicones. A curing time of a few minutes at room temperature is especially useful for encapsulation of live tissue or implanting into a live animal.

Unlike many vinyl-based silicones in which the platinum complex catalysts are easily poisoned by contamination from amines and animal tissue, **Kwik-Sil** and **Kwik-Cast** are not sensitive to contamination from animal tissue

Kwik-Sil™ is a translucent, medium-viscosity silicone adhesive, developed for chronic peripheral nerve studies such as anterograde tracing with fluorescent indicators or electrode recording. Good adhesion and mechanical properties (tear strength and elongation) allow days of study without breaking of the bonding. Curing speed is very reproducible.

Kwik-Cast™ is a very low viscosity silicone sealant developed to embed peripheral nerves with electrodes for acute multi-fiber recordings. It flows easily, filling the small spaces around the nerve and leaving no channels through which peritoneal fluid can travel and thus short the nerve/electrode contact. Equally important is the ability of the material to flow into itself and create one continuous mass from underneath the nerve all the way to the top of the nerve/electrode contact to ensure long-term recording stability. Kwik-Cast is color-coded to make the mixing foolproof. The catalyst is yellow and the base is blue. When uniformly mixed, it is green. Kwik-Cast can be applied and cured underneath mineral oil. After recording, electrodes are easily recovered due to the low tear strength.



KWIK-CAST

References

Maura Casciola, Shu Xiao & Andrei G. Pakhomov. Damage-free peripheral nerve stimulation by 12-ns pulsed electric field. *Scientific Reports 7*, Article number: 10453. https://www.nature.com/articles/s41598-017-10282-5

Tiran, E., Ferrier, J., Deffieux, T., Gennisson, J.-L., Pezet, S., Lenkei, Z., & Tanter, M. (2017). Transcranial Functional Ultrasound Imaging in Freely Moving Awake Mice and Anesthetized Young Rats without Contrast Agent. *Ultrasound in Medicine & Biology*, 43(8), 1679–1689. http://doi.org/10.1016/j.ultrasmedbio.2017.03.011

Lee, D., & Lee, A. K. (2017). Efficient Method for Whole-Cell Recording in Freely Moving Rodents Using Ultraviolet-Cured Collar-Based Pipette Stabilization. *Cold Spring Harbor Protocols*, 2017(4), pdb.prot095810. http://doi.org/10.1101/pdb.prot095810

Welle, T., Hoekstra, A. T., Daemen, I. A. J. J. M., Berkers, C. R., & Costa, M. O. (2017). Metabolic response of porcine colon explants to in vitro infection by Brachyspira hyodysenteriae: a leap into disease pathophysiology. *Metabolomics*, 13(7), 83. http://doi.org/10.1007/s11306-017-1219-6

Lee, D., & Lee, A. K. (2017). In Vivo Patch-Clamp Recording in Awake Head-Fixed Rodents. *Cold Spring Harbor Protocols*, 2017(4), pdb.prot095802. http://doi.org/10.1101/pdb.prot095802

Suwabe, T., Fukami, H., & Bradley, R. M. (n.d.). Synaptic Responses of Neurons Controlling the Parotid and von Ebner Salivary Glands in Rats to Stimulation of the Solitary Nucleus and Tract. http://doi.org/10.1152/jn.01115.2007

KWIK-CAST & KWIK-SIL SPECIFICATIONS			
Kwik-Sil	Kwik-Cast		
1 to 1	1 to 1		
< 5 minutes*	4 minutes		
5-10 minutes**	<10 minutes		
~15 minutes			
15,000	10,000		
6 months	6 months		
5 mL	5 mL		
10	10		
<0.12 mL	<0.12 mL		
90	44		
650	60		
30	36		
translucent	green		
1x10 ¹⁵	1x10 ¹⁵		
	Kwik-Sil 1 to 1 < 5 minutes* 5-10 minutes** ~15 minutes 15,000 6 months 5 mL 10 <0.12 mL 90 650 30 translucent		

^{* 3} minutes average with about 90 seconds of liquidity

^{**} no longer mixable at this point

	ORDERING INFORMATION
KWIK-SIL	Silicone Adhesive Compound (two 5-mL syringes)
KWIK-CAST	Silicone Casting Compound (two 5-mL syringes)
600022	Replacement KWIK Mixing Tips (pack of 10)
	Ouantity discounts available

"Super" Adhesives for Life Science Research

Four times stronger than butyl cyanoacrylate and less toxic

Cyanoacrylate adhesives have been on the market since 1958. Most industrial or household grade cyanoacrylate is made of shorter alkyl chain derivatives such as methyl or ethyl cyanoacrylate (WPI's 7341 and 7342). They are very useful for temporarily holding tissues such as mounting specimens for microtome sectioning. However, they are not suitable for bonding wounds on live animals. The difficulties of using cyanoacrylate for bonding live animals are: (1) a strong, irritating odor; (2) quick loss of bonding strength due to breakdown of the bonding by hydration; (3) the breakdown products, cyanoacetate

and formaldehyde, are toxic and



VETROND

can cause inflammatory reactions; and (4) they have low flexibility and tend to be brittle.

To overcome these problems, several longer alkyl chain cyanoacrylates have been developed especially for veterinary and human use. The first longer alkyl chain product is butyl cyanoacrylate. This product has been used for animal and human applications outside the USA since 1970. It is much less toxic and has a lower odor than the methyl or ethyl cyanoacrylate. The butyl cyanoacrylate offered by WPI is **Vetbond**[™].

A family of adhesives containing octyl cyanoacrylate, a plasticizer and stabilizer, was developed In the 1990's (one of them approved by FDA for human use). When bonding to tissue, these new adhesives are four times stronger and less toxic than butyl cyanoacrylate. Compared with the traditional suture, the new super adhesive has several advantages. On average, it takes only one-tenth of the time to close an incision. The bonding strength is equal to 5-0 monofilament suture. It also has a

mysterious antimicrobial effect that can decrease infection rates in contaminated wounds. Bonding will slough off naturally in 5 to 7 days. Cosmetic appearance of the healed incision is also better

Gluture Topical Tissue Adhesive **503763** forms a strong and flexible film and is thus more suitable for surface wound bonding, protection, and holding a sensor or other device on the tissue.



Setting time is about 10 seconds, which gives ample time for application. It can also be used for temporarily holding a live tissue. For example, there is a report of using it to hold nematodes on a glass slide for patchclamp neurons recording.

All of the products offered by WPI are veterinary grade (not suitable for human application). Though very similar to the grade for human use, they are not sterile and do not have FDA approval.

	ORDERING INFORMATION
503763	Gluture Topical Tissue Adhesive (10 tips), 1.5 mL
7341	Cyanoacrylate Adhesive, low viscosity—90-120 cps (pack of 10 vials, each approx. 1.5 mL)
7342	Cyanoacrylate Adhesive, high viscosity—1100-1600 cps (pack of 10 vials, each approx. 1.5 mL)
VETBOND	3M Vetbond™ Adhesive (3 mL)

Sylgard



A two-part silicone elastomer, ideal for potting and encapsulating applications. Very low dielectric constant sealing compound used in patch clamping and many other lab applications. After cure, will withstand -55° to 200 °C.

ORDERIN	IG INF	ORMAT	ION.
OILDEILII			

SYLG184 Sylgard, 0.5 kg (1.1 lb)

Silicone Dissecting Pad Kit



Make your own silicone dissecting pads easily and quickly. Mix the 2-part silicone right in the plastic Petri dishes and allow to cure 24 hours at room temperature. Kit includes enough two-part Sylgard silicone elastomer to prepare 20 dishes; pins; and 20 plastic Petri dishes with lids, 65 mm.

ORDERING INFORMATION

501986 Silicone Dissecting Pad Kit

WORLD PRECISION INSTRUMENTS

www.wpiinc.com

Electrically Conductive Silver Epoxy



Two-component silver-filled epoxy for electrical connections which cannot be soldered, such as Ag/AgCl pellets. This widely used silver-filled epoxy features low viscosity and smooth flowing character. Pure silver is dispersed in both resin and hardener. Cures in 15 minutes at 120 °C. Mix ratio 1:1. May be premixed and frozen for later use.

ORDERING INFORMATION

4898

Silver Epoxy, 28 g (1 oz)

Electrically Conductive Carbon Epoxy



Two-component carbon-epoxy, curable at room and elevated temperatures. Ideal for electrostatic discharge protection and EMI/RFI shielding. 1:1 mix ratio. May be premixed and frozen for later use.

ORDERING INFORMATION

Carbon Epoxy, 56 g (2 oz)

Silicone RTV Adhesive (non-acidic)



Because it is non-corrosive, this material is ideal for use on metal, for encapsulating small circuits on connectors. After cure, will withstand -55° to 200°C. No mixing required.

ORDERING INFORMATION

7128 RTV Coating, 90 mL (3 fl oz)

Silicone RTV Adhesive



Clear silicone sealant provides good bonding to plastic. After cure, will withstand -55 to 200 °C. No mixing required. A handy, general purpose laboratory sealant. (Releases acetic acid during curing.)

ORDERING INFORMATION

1571 RTV Sealant, 139 mL (4.7



RTV Prime Coat

Enhances adhesion of silicone adhesives to many difficult-to-bond plastic surfaces.

ORDERING INFORMATION

RTV Prime Coat, 400 mL (13.5 fl oz)




ORDERING INFORMATION

13316 Mini Glue Gun

WORLD PRECISION INSTRUMENTS

Biosensing



Selection includes high selectivity and low detection limit sensors

WPI's biosensors are unique, because they offer a high selectivity and low detection limit (down to nM concentration) with a broad dynamic range; covering physiological concentration of species with different sizes from nM to mM. The majority of our sensors are the only commercially available sensors in the world. Scientists across a variety of disciplines have relied on our sensors for over 25 years. These scientists use WPI's sensors for research performed in universities, hospitals, biomedical research labs, pharmaceutical companies, food/ nutrition research labs, environmental monitoring centers and military labs. Our popular biosensors are listed in thousands of publications.

Leading Scientist Heads WPI Biosensing Division

About Dr. Xueji Zhang

The notable expert in the area of electrochemical sensors, Xueji Zhang, PhD, Sr. Vice President and Chief Scientist at WPI, is leading the development for microsensors (including the ISO-H2S-100-Cxx) for many applications. His research includes bio-analysis, electrochemical sensors and biosensors, microelectrodes, ultramicroelectrodes and nanoelectrodes, nanosensors, free radical sensors, nitric oxide sensors, cancer diagnostics, design and application of biomedical instrumentation.


Dr. Zhang has written hundreds of published research papers. Listed below is a select list of Dr. Zhang's papers covering the detection of Hydrogen Sulfide. We are very proud to share the research findings with you. Please email wpi@wpiinc.com for a copy of the paper(s).


- Zhang XJ, et al., Electrochemical Sensors Biosensors and their Applications, Elsevier 2008.
- Zhang XJ, et al., Electrochemical hydrogen sulfide biosensors. Analyst. 2016;141:1185-95.
- Zhang XJ, et al., A novel enzymatic method for determination of homocysteine using electrochemical hydrogen sulfide sensor, Frontiers in Bioscience, 2007, 12, 3774


Dr. Zhang is a national chair professor and director at the University of Science & Technology, Beijing, China. He is also an adjunct professor of University of Tokyo, University of South Florida, Chinese Academy of Sciences, PLA 301 hospital. He is president of the Chinese Alliance of Biodetection & Biomonitoring. He is the honorary president of MITBJ. He has published over 400 papers in peer-reviewed journals, authored 8 books and 80 patents. He has developed numerous biomedical sensors, instruments and devices for commercialization. Dr. Zhang also serves as the editor-in-chief of Frontiers in Bioscience (www.bioscience.org) and The American Journal of Biomedical Sciences (www.nwpii.com/ajbms) and has been a member of the advisory editor board of 19 other international journals. He has received numerous awards and honors, including academician of Russian Academy of Engineering (2013), Fellow of American Institute for Medical & Bioengineering (2016), Fellow of Royal Chemistry Society (2014), Fellow of Chinese Chemical Society (2014), Scientist of the Year in China (2015), Engineer of the Year in China (2016), National Innovation Award (2017) and W. Simon Fellow of ICSC-World lab, United Nations (1996).











Biosensors Designed for Selective Detection

Selecting microsensors for your application

Choosing the right sensor for your application is critical for successful research. The manufacturer's specifications provide valuable information for selecting your sensor. Consider the following five performance factors in reference to your application:

Response–Electrochemical electrodes produce changes in current in response to changes in concentration. "Response" is most often specified in terms of the amount of current per concentration unit: nA/ μ M or pA/nM, etc. The larger the current per concentration the higher the sensitivity of the sensor.

Detection Limit–Detection limit is the minimum change in concentration that can be reliably seen. This specification is directly related to the noise of the sensor. A sensor with a 100nA/ μ M response but a 3 μ M detection limit is not as good as a 10nA/ μ M response sensor with a 1 μ M detection limit.

Free Radical Detection–The best sensors have low detection limits and high sensitivity. A sensor can have a low detection limit and a good response, however, to be useful in long term studies it must be stable

when temperature and concentration are constant. A drifting baseline, if monotonic, can be corrected, but wandering baselines limit the utility of sensors to short experiments.

Selectivity–It is a rare instance that the ion species of interest is the only ion in the medium to be measured. In a perfect world your sensor would respond ONLY to the ion of interest. In reality there is always some contribution from competing species. The lower the contribution the better.

Linearity–For an electrode to be useful and easy to calibrate the response must be "linear" with changes in concentration over the range of interest. Non-linear behavior requires special curve fit software to calibrate the sensors. This approach is more time consuming and can be unreliable. "Good" linearity is expressed by a R² of 0.980 or higher. (1.00 is perfect.) All of the electrochemical sensors made by WPI are 0.98 or better.

The table below presents the specifications of WPI's macro sensors (2 mm sensor) and microsensors.

MACRO SENSORS					
SPECIES	Carbon Monoxide	Nitric Oxide	Hydrogen Peroxide	Oxygen	Hydrogen Sulfide
Order Number	ISO-COP-2	ISO-NOP	ISO-HPO-2	ISO-OXY-2	ISO-H2S-2
Available Diameters	2 mm	2 mm	2 mm	2 mm	2 mm
Response Time	< 10 sec	< 5 sec	< 5 sec	< 10 sec	< 5 sec
Detection Limit/Range	10 nM to 10 μM	1 nM to 40 μM*	< 100 nM to 100 μM	0.1%-100%	< 5 nM to 100 μM
Sensitivity	~0.5 pA/nM	≤ 2 pA/nM	8 pA/μM	0.3-0.6nA/%	2 pA/nM
Drift	< 1 pA/min	< 1 pA/min	< 1 pA/min	≥ 1%/min	< 1 pA/min
Temperature Dependent	Yes	Yes	Yes	Yes	Yes
Physiological Interference	nitric oxide	NaNO ₂ (10 ⁻⁶ or better)	None	None	None
Replacement Sleeves (pkg. of 4)	95620	5436	600012	5378	600016
Filling Solution	95611	7325	100042	7326	100084
Start-up Kit	95699	5435	600011	5377	600015

^{*}Higher detection limit available on request — call for custom pricing.

MICRO SENSORS												
	ISO-NOPF200	ISO-NOPF200-Lxx³	ISO-NOPF100	ISO-NOPF100-Lxx³	ISO-NOPF500-Cxx	ISO-NOP3005	ISO-NOP3020	ISO-NOP007	ISO-NOPNM	ISO-HPO-100	ISO-HPO-100-L	ISO-H2S-100-Cxx
Species	Species Nitric Oxide H ₂ O ₂ H ₂ S							H ₂ S				
Package Quantity	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 2)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 3)	(pkg. of 2)	(pkg. of 2)
Fiber Diameter (µm)	200	200	100	100	500	30	30	7	7 Conical tip: 100nm	100	100	100
Tip Length ² (mm)	1-5 ¹	1-10 ¹	1-5 ¹	1-5 ¹	5–10	0.5	2	2	2	1-5 ¹	1-5 ¹	2-5 ¹
Response Time (sec.)	< 5	< 5	< 5	< 5	< 10	< 3	< 3	< 3	< 3	< 5	< 5	< 5
Lowest Detection Limit/Range (nM)	0.2	0.2	0.2	0.2	0.2	1	1	0.5	0.5	1	1	< 5
Nominal Sensitivity-New Sensor ² (pA/nM)	≥ 10	≥ 20	≥ 6	≥ 20	≥ 20	≥ 0.2	≥ 1	≥ 0.5	≥ 0.5	≥ 1	≥1	1-4
Baseline Drift (pA/min.)	none	none	none	none	none	none	none	none	none	< 2.0	< 2.0	< 2

¹Sensor available in 1 mm length increments (for example, 1 mm, 2mm, 3mm...)

Any 100 µm sensor can be purchased with a hypodermic sheath. Add -H to the end of the part number (for example, ISO-HPO-100-H).

Some nitric oxide sensors are available in custom lengths. When ordering custom lengths, use the part numbers **ISO-NOPF100-Cxx** or **ISO-NOPF200-Cxx** and replace the **xx** with the desired length. For example, for a 1 mm flexible sensor tip, the part number should be **ISO-NOPF200-C01**. Sensors can be ordered in the following custom lengths: 1 mm, 2mm, 3mm, 4mm or 5 mm (**ISO-H2S-100**: 2-5mm only).

²Sensor sensitivity varies with length and diameter.

³L-shaped sensor for use with a tissue bath.

Nitric Oxide Sensors

For routine detection of NO at ultra low concentrations

Features

- Excellent selectivity to NO
- Rapid response time
- Highly sensitive

Benefits

 NO sensors in different sizes from 100 nm to a few mm can be used for many NO detection applications, such as single cell measurement, in vivo measurement in tissues (even in animals), and NO release from drugs

Applications

- Cell culture, cell suspensions
- Measurements at the cellular level
- Arteries, microvessels, in vivo applications
- Tissue bath applications

WPI offers the most extensive range of nitric oxide (NO) sensors available. Developed over a decade of extensive research in the field of nitric oxide, the result is a superior range of NO sensors that enable routine detection of nitric oxide at ultra low concentrations.

The ideal NO sensor should be insensitive to other reactive species likely to be present within the measurement environment. The conventional Nafion coated carbon fiber NO micro sensor exhibits a large response to such species. WPI's unique NO sensor technology utilizes a novel surface membrane which amplifies the response to NO while eliminating responses to a vast range of reactive species, including nitrite, ascorbic acid, hydrogen peroxide, catecolamines, and much more.

ISO-NOP

The original nitric oxide probe - ideal for cell cultures, cell




suspensions and many other applications

The **ISO-NOP** is a popular, robust and high performance sensor encased within a 2 mm diameter disposable stainless steel protective sleeve. The tip of the sleeve is covered with a NO-selective membrane. Replacement membrane sleeves can be purchased separately (5436) and require an internal electrolyte (7325).

A simple change in experimental protocol enables the ISO-NOP to be

conveniently used for indirect rapid accurate determination of nitrite (NO₂) and nitrate (NO₃) concentration in samples. Using this method a detection limit for NO₂ or NO₃ as low as 1 nM is routinely possible.



Abdominal X-ray showing the appratus consisting of two customized ISO-NOP nitric oxide probes, 4-channel pH catheter, and Teflon nasogastric tube. (Courtesy Prof. K.E.L. McColl, University Department of Medicine and Therapeutics, Western Infirmary, Glasgow, Scotland.) lijima, K., et al. Gastroenterology 2002: **122**: 1248-1257.

Brazil: (013) 406-29703 · info@brazil.wpiinc.com · www.wpiinc.com

References

Zhang X.J. Real time and *in vivo* monitoring of nitric oxide by electrochemical sensors—From dream to reality, *Frontiers in Bioscience*, 9,3434-3446(2004)

Zhang X.J., Ju H.X., Wang J. Electrochemical Sensors, *Biosensors and Their Biomedical Applications*, Elsevier, 2008

Zhang X.J., et al. Nanometer size Electro for nitric oxide and s-nitrosothiols measurement, *Electrochem. Commun.* 4, 11-16 (2002)

Zhang X.J. Construction and Characterization of a new *in vivo* nitric oxide microsensor, *Electroanal.* 9, 640-643, (2004)


Zhang X.J. An integrated nitric oxide sensor based on carbon fiber electrode coated with selective membranes, 12. 1113-1117(2000)

ISO-NOPNM

The world's smallest nitric oxide NanoSensor, designed for measurement of NO at the cellular level.


The ${\bf ISO-NOPNM}$ NanoSensor has a tip diameter of just 100 nm (0.1 μ m) and a detection limit for NO of less than 0.5 nM — making it indisputably the smallest and most sensitive NO sensor in the world!

The ISO-NOPNM is based on a novel design in which an electrochemically



"activated" composite graphite nanofiber is used as the NO-sensing element. The surface of the NanoSensor is then modified using a

unique multi-layered NO-selective membrane. Figure at right illustrates the response of the **ISO-NOPNM** following successive additions of nanomolar concentrations of NO. The ultra-low noise of the **ISO-NOPNM** (0.5 pA) enables a detection limit of just 0.5 nM NO. The response time of **ISO-NOPNM** is less than 3 s.




Amperometric response of the NO nanosensor (ISO-NOPNM) to the successive additions of 2 nM, 4 nM, 8 nM NO into 0.1 M CuCl₂.


ISO-NOPF

Unique flexible NO sensor! Designed for arteries, microvessels, in vivo applications and similar applications.

ISO-NOPF electrodes are available in 100 µm and 200 µm diameters. Utilizing the latest advances in nano-technology and material science, scientists at WPI's Sensor Laboratory have created these completely flexible and virtually unbreakable NO sensors. The new sensors are based on a composite graphite NO-sensing element combined with a



reference electrode. The surface of the sensor is then coated with a unique multi-layered NO-selective membrane.





NITIRC OXIDE SENSOR SPECIFICATIONS								
ISO-NOPF ISO-NOP ISO-NOPF100-L10 ISO-NOPF200-L10 ISO-NOP30 ISO-NOP007 ISO-NOPNM								
APPLICATION	In Vivo	Cell Cultures, NOz, NOz	Tissue Bath	Tissue Bath	Microvessels	Microvessels	Single Cell	
SENSOR DIAMETER	100, 200, or 500 μm	2 mm	30 µm	70 μm	30 µm	7 μm	100 nm	
RESPONSE TIME	< 5 sec	< 5 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec	< 3 sec	
LOWEST DETECTION LIMIT	0.2 nM	1 nM	1 nM	1 nM	1 nM	0.5 nM	0.5 nM	
TEMPERATURE SENSITIVITY	yes	yes	yes	yes	yes	yes	yes	
DRIFT	none	none	none	none	none	none	none	
SENSITIVITY	8 pA/nM	1-2 pA/nM	≥1 pA/nM	≥1 pA/nM	1.4 pA/nM	0.5 pA/nM	0.5 pA/nM	
PHYSIOLOGICAL INTERFERENCE	none	none	none	none	none	none	none	

All WPI NO sensors are 100% compatible with the TBR4100 and TBR1025 Free Radical Analyzers.
They are also compatible with the discontinued ISO-NO Mark II (NOMK2), APOLLO 4000 and APOLLO 1000.

ISO-NOP007, ISO-NOP30xx

7 and 30 micron sensors with exceptional performance — ideal for tissues and microvessels






The response of a 7 µm NO sensor (ISO-NOP007) to successive additions of NO (100 nM). Inset shows the linearity of the resulting calibration plot.


The **ISO-NOP007** has a tip diameter of just 7 microns and a length of 2 mm. The ISO-NOP30 has a tip diameter of 30 microns and is available in two different tip lengths (i.e., ISO-NOP3020 has tip length of 2 mm, ISO-NOP3005 has tip length of 0.5 mm). The response of the ISO-NOP007 and ISO-NOP30 is linear over a wide dynamic concentration range of NO. The design of both electrodes is based on a single carbon fiber coated with WPI's NO-selective membrane. A detection limit of approximately 1 nM NO makes these electrodes ideal for use in tissues and microvessels.

ISO-NOPF100-L10, ISO-NOPF200-L10, ISO-HPO-100-L



L-shaped sensors for tissue bath & cell culture studies

These unique L-shaped nitric oxide sensors are designed specifically for use in tissue bath studies and similar applications. The shape of the sensor has been engineered to facilitate placement of the electrode within the lumen of the tissue vessel under study. The ISO-NOPF200-L10 is designed specifically for cell culture studies.



Simultaneous measurement of force (top trace) and changes of NO concentration (lower trace) in (A) the rat superior mesenteric artery relaxed with ACh and (B) a small human artery relaxed with ACh and SNAP using and ISO-NOP30. In this artery oxyhaemoglobin(oxyHb) partly reversed the increase in NO concentration, with only a small change in force. [U. Simonsen, et al., J. Physiol., 1999, **516**: 271-282.]

RDERING INFORMATION
Replacement 2 mm shielded sensor and cable
100 nm NanoSensor, pkg. of 3 (requires cable 91580)
100 μm Flexible NO Sensor, pkg. of 2
200 μm Flexible NO Sensor, pkg. of 2
500 μm Flexible NO Sensor, pkg. of 2
7 μm Nitric Oxide Sensor (pkg. of 3)
30 μm Sensor Tips (2 mm length), pkg. of 3 (requires 91580)
30 μm Sensor Tips (0.5 mm fiber), pkg. of 3 (requires 91580)
NO Sensor, 100 μm Flexible L-shaped (pkg. of 2)
NO Sensor, 200 µm Flexible L-shaped (pkg. of 2)
HPO Sensor, L-Shaped 100 μm (pkg. of 2)
ISONOP Startup Kit (recommended with first purchase)
Replacement Sleeve Kit for 2mm sensor, pkg. of 4
ISO-NO Electrolyte (10 mL)
ISO-NO Electrolyte, CO2-insensitive (10 mL)
T-Adapter Kit (pkg. of 3) for ISO-NOP
Nitrite Standard Solution, 1 gram/liter (100 mL)
Microsensor Adapter Cable
SNAP, 50 mg vial

*Tip diameters given do not include the coatings. Tips are bare wire.

www.wpiinc.com

Hydrogen Sulfide Sensors

Low detection limit sensors to record in vitro

Features

- High selectivity and sensitivity to H₂S with fast response time
- Broad linearity


Benefits

- A dry sensor that works like a 2 mm wet sensor
- Use like an ISO-H2S-2 without sleeves and filling solution


Applications

- Measure biological H₂S in vivo or in vitro in tissues or biological media
- Measure H₂S released from drugs
- Cell culture, cell suspensions, arteries, in vivo applications

Although hydrogen sulfide (H₂S) is generally thought of as a poisonous gas, it is endogenously produced in many mammalian tissues. It has



Stepped response to increasing concentrations of H₂S



The sensor is insensitive to competing species such as ascorbic acid, dopamine, nitrate, nitrite and

been detected in micromolar amounts in blood and brain tissue. Hydrogen sulfide is reported as having a broad range of biological functions and although its potential to participate in cell signaling is clear, this biological role is not well understood. H₂S is strongly analogous to nitric oxide (NO) because they share several physical and metabolic properties.

Like NO, H₂S is a potent vascular signal that can mediate vasoconstriction or vasorelaxation depending on the O₂ level and tissue. In the rat aorta, H₂S concentrations that mediate rapid constriction at one O₂ level will cause rapid relaxation at lower O_2 levels.

The ISO-H2S sensor is a low detection limit sensor to record H₂S in vitro. This is the only sensor available that measures H₂S. The **ISO-H2S-100** is a hydrogen sulfide sensor with a 100 µm diameter tip. It is designed like the dry microsensors, however, it works like a traditional 2 mm sensor. The sensor can be ordered in a variety of lengths from 2-5 mm. It incorporates WPI's proprietary combination electrode technology in which the hydrogen sulfide-sensing element and separate reference electrode are encased within a single shielded sensor design. The ISO-H2S-100 offers several advantages:

- Requires no sleeves or filling solutions like 2mm macrosensors
- Durable for long-term use, because of its platinum wire construction
- · Rapid response time
- Broad linear range (Range is based on the length of the sensor tip)

ORDERING INFORMATION

ISO-H2S-100-Cxx Hydrogen Sulfide Micro Sensor (pkg. of 2) 2 mm Shielded Hydrogen Sulfide Macro Sensor

Hydrogen Peroxide Sensors

Direct quantitative measurements in biological samples



Benefits

• Non breakable integrated hydrogen peroxide sensor with tip dimension of 100 µm, detection limit down to nM range. Tip size can be custom made as small as a few micrometers.

Applications

- Cell culture, cell suspensions
- Cell tissue measurements

Despite the recognized importance of this oxidant in biology, real-time measurements at low concentration have been difficult. The hydrogen peroxide sensors developed by WPI are designed to compliment existing high sensitivity fluorescent approaches with direct quantitative measurement in biological samples in the low nM range.

The ISO-HPO-2 is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (600012) and an internal refillable electrolyte (100042). It is designed for use in cell cultures and similar applications.

The ISO-HPO-100 is a 100 µm tip diameter hydrogen peroxide micro

sensor designed for use in tissues and similar applications. The design is based on a platinum wire sensing electrode coated with a proprietary membrane to enhance H₂O₂ detection.

These sensors incorporate WPI's proprietary combination electrode technology whereby the hydrogen peroxide sensing element and separate reference electrode are encased within a single Faradayshielded probe design. This design has been shown to enhance performance during measurements and minimizes overall sensor size.

Our hydrogen peroxide (H2O2) sensors work with the TBR4100 and TBR1025 free radical analyzers.

H ₂ O ₂ SENSOR SPECIFICATIONS							
	ISO-HPO-2	ISO-HPO-100	ISO-HPO-100 H	ISO-HPO-100-L			
APPLICATION	Cell Cultures, Cell Suspensions	Tissue/Microvessels	Hypodermic Sheath	Tissue Bath			
SENSOR DIAMETER	2.0 mm	100 μm	100 μm	100 μm			
RESPONSE TIME	< 5 SEC (90%)	< 5 SEC (90%)	< 5 SEC (90%)	< 5 SEC (90%)			
DETECTION LIMIT	$<$ 100 nM to 100 μM	1 nM to 1 mM	< 10 nM to 1 mM	1 nM to 100 μM			
DRIFT	< 0.2 pA/min	< 2.0 pA/min	< 2.0 pA/min	< 2.0 pA/min			
SENSITIVITY	8 pA/μM	1 pA/nM	1 pA/nM	1 pA/nM			
PHYSIOLOGICAL INTERFERENCE	none	Contact WPI	Contact WPI	Contact WPI			

	ORDERING INFORMATION
600011	ISO-HPO Startup Kit (recommended with first purchase)
ISO-HPO-2	2 mm Shielded HPO Sensor & Cable
ISO-HPO-100	100 μm HPO Sensor*, pkg. of 3
ISO-HPO-100-L	100 μm HPO Sensor, L-shaped*, pkg. of 2
ISO-HPO-100H	100 μm HPO Sensor in hypodermic sheath*, pkg. of 3
600012	Replacement Sleeve Kit for ISO-HPO-2, pkg. of 4
100042	ISO-HPO-2 Electrolyte (10 mL)
91580	Microsensor Adapter Cable
-1 1	Paguiras 01590 Microsopsor Adapter Cable

* Requires 91580 Microsensor Adapter Cable

Oxygen Sensors

Make direct quantitative measurements in biological samples



Features


- Gas permeable polymer membrane sleeve blocks liquids, ions and particulate matter
- Incorporated reference electrode

Benefits

• Clark type of oxygen sensor with tip size of 2 mm, can be used for very small volume O2 measurement

Applications

- Cell culture, cell suspensions
- Cell tissue measurements



This sensor incorporates WPI's proprietary combination electrode technology whereby the oxygen-sensing element and separate reference electrode are encased within a single shielded sensor design. A gaspermeable polymer membrane is fitted over the end of the sleeve, which allows oxygen to pass while blocking liquids, ions and particulate matter.

Oxygen diffuses through the membrane. The voltage applied to the sensor is held at -0.7V when the monitoring device is on and the sensor is properly connected. The magnitude of the generated electrical current is determined by the rate of diffusion through the membrane. The rate is proportional to the partial pressure of oxygen outside the membrane. The current serves as a measure of the partial pressure of O_2 .

The ISO-OXY-2 is a 2.0 mm stainless steel sensor, with replaceable membrane sleeves (5378) and an internal refillable electrolyte (7326).

O₂ SENSOR SPECIFICATIONS

APPLICATION	Cell cultures, cell suspensions.
SENSOR DIAMETER	2.0 mm
RESPONSE TIME	<10 SEC (90%)
DETECTION LIMIT	0.1 % to 100%
DRIFT	<1%/min
SENSITIVITY	N/A
PHYSIOLOGICAL INTERFERENCE	None

ISO-OXY Startup Kit (recommended with first purchase)Replacement Electrode Sleeve Kit, pkg. of 4		ORDERING INFORMATION
5378 Replacement Electrode Sleeve Kit, pkg. of 4	ISO-OXY-2	2 mm Shielded Oxygen Sensor & Cable
	5377	ISO-OXY Startup Kit (recommended with first purchase)
T005	5378	Replacement Electrode Sleeve Kit, pkg. of 4
7326 ISO2 Filling Solution (electrolyte)	7326	ISO2 Filling Solution (electrolyte)

Implantable Glucose Sensor

Measure glucose in vivo or in vitro over long term

IGS100

Features

- Implantable sensor for long term studies
- Incorporated reference electrode

Benefits

 Implantable microsensor,



· Long term monitoring

Applications

• In vivo long term measurement of glucose in animals or tissues

Measuring glucose in vivo over the long term is challenging and difficult. Previous measurement systems were limited to acute studies or a few days at best. WPI introduces a new kind of implantable glucose sensor based on a patented technology. This sensor provides a tool for researchers to directly detect glucose in chronic studies in vitro or in vivo. The sensor is fully compatible with WPI's TBR systems.

GLUCOSE SENSOR SPECIFICATIONS

IN VITRO PRECISION Coefficient of Variation (CV) ≤5% **GLUCOSE RANGE** 36-450 mg/dl (or 2-25 mM/L) RESPONSE TIME (SEC) 100-300 s

IN VIVO CALIBRATION In vivo calibration

INTERFERENCE SPECIES Acetaminophen, ascorbic acid, uric acid

LENGTH 5 cm SENSOR SIZE $0.6 \times 0.7 \text{ mm}$ REFERENCE ELECTRODE Ag/AgCl

POLARIZATION VOLTAGE (V) 0.65 -0.7 V vs. Ag/AgCl

3-4 months in solutions at room temperature SENSOR LIFE under continuous polarization; 15-30 days in vivo

SHELF LIFE

OPERATION CONDITIONS 20° to 40° C (68° to 104° F) STORAGE CONDITIONS 10° to 25° C (50° to 77° F)

	ORDERING INFORMATION
100	Implantable Glucose Sensor (pkg. of 2)

IGS1 91580 Microsensor Adapter Cable

Four-Channel Free Radical Analyzer

Fast, reliable, real-time detection— Measure redox-reactive species



Features

- Real-time detection using electrochemical sensors
- Integrated system includes one temperature sensor, your choice of two additional sensors and a start-up kit
- Current measurement range from 300 fA to 10 μA (four ranges) permits wide dynamic range for detection
- Wide bandwidth allows recording of fast events
- Measure carbon monoxide from 10 nM to 10 μM
- Measure nitric oxide from < 0.3 nM to 100 μM
- Measure hydrogen peroxide < 10 nM to 100 mM
- Measure hydrogen sulfide
- Measure glucose
- Measure oxygen from 0.1% to 100%
- Isolated architecture allows Lab-Trax interface to simultaneously measure free radical and independent analog data (for example, ECG, BP, etc.) on any channel

Benefits

- Measure up to four different species and temperature in the same preparation or simultaneous measurement in four different preparations
- Lab-Trax data acquisition system is flexible

Applications

• Free radical detection (NO, H₂O₂, H₂S, CO, O₂ and glucose)

Real-time detection and measurement of a variety of redox-reactive species is fast and easy using the electrochemical (amperometric) detection principle employed in the **TBR4100**. This optically isolated four-channel free radical analyzer has ultra low noise and independently operated channels.

Measure four species simultaneously

For use with WPI's wide range of nitric oxide, hydrogen peroxide, hydrogen sulfide and oxygen sensors, the **TBR4100** can measure four different species simultaneously in the same preparation. Simply plug a sensor into any one of the input channels on the front panel and select the current range. Poise voltage can be selected from a range of values tuned for optimal response from WPI sensors. An independent output for real-time monitoring of temperature is also included.

Lab-Trax data acquisition system is flexible

The **TBR4100** analyzer utilizes PC-based data acquisition via our **Lab-Trax** interface. Data traces are displayed and recorded in realtime. The LabScribe software comes pre-configured for single or multiple electrode recording; filters, gains and smoothing are all set for optimal results. Data can be viewed making adjustments to smoothing and filter settings without affecting the original stored raw data. Electrode calibration from multiple concentration readings can be input into the software's Multipoint Calibration utility quickly, providing a plot and slope calculation for electrode sensitivity determination. Alternately, the **Lab-Trax** data interface can be used for providing simultaneous acquisition of Free Radical data along with other physiological data (ECG, HR, BP, *etc.*) as each of the four input channels has its own independent input, filters and 24-bit converter. *See www.wpiinc.com/TBR4100 for more information on Lab-Trax data acquisition*.

References

Zhang X.J. Real time and *in vivo* monitoring of nitric oxide by electrochemical sensors—From dream to reality, *Frontiers in Bioscience*, 9,3434-3446(2004) **Zhang XJ. Ju HX. Wang J.** Electrochemical Sensors, *Biosensors and Their*

Biomedical Applications, Elsevier, 2008 Murine strain differences in inflammatory angiogenesis of internal wound in diabetes. (2017). Biomedicine & Pharmacotherapy, 86, 715–724. https://doi.

org/10.1016/J.BIOPHA.2016.11.146

Xie, L., Feng, H., Li, S., Meng, G., Liu, S., Tang, X., ... Ji, Y. (2016). SIRT3

Mediates the Antioxidant Effect of Hydrogen Sulfide in Endothelial Cells


Mediates the Antioxidant Effect of Hydrogen Sulfide in Endothelial Cells. Antioxidants & Redox Signaling, 24(6), 329–343. http://doi.org/10.1089/ars.2015.6331 Nguyen, T.-K., Selvanayagam, R., Ho, K. K. K., Chen, R., Kutty, S. K., Rice, S.

A., ... Boyer, C. (2016). Co-delivery of nitric oxide and antibiotic using polymeric nanoparticles. *Chem. Sci.*, 7(2), 1016–1027. http://doi.org/10.1039/C5SC02769A Chen, G., Yang, L., Zhong, L., Kutty, S., Wang, Y., Cui, K., ... Bin, J. (2016). Delivery of Hydrogen Sulfide by Ultrasound Targeted Microbubble Destruction Attenuates Myocardial Ischemia-reperfusion Injury. *Scientific Reports*, 6, 30643. http://doi.org/10.1038/srep30643

Xu, T., Scafa, N., Xu, L.-P., Zhou, S., Abdullah Al-Ghanem, K., Mahboob, S., ... Zhang, X. (2016). Electrochemical hydrogen sulfide biosensors. *The Analyst*, 141(4), 1185–1195. http://doi.org/10.1039/C5AN02208H

Zong, Y., Huang, Y., Chen, S., Zhu, M., Chen, Q., Feng, S., ... Jin, H. (2015). Downregulation of Endogenous Hydrogen Sulfide Pathway Is Involved in Mitochondrion-Related Endothelial Cell Apoptosis Induced by High Salt. Oxidative Medicine and Cellular Longevity, 2015, 1–11. http://doi.org/10.1155/2015/754670

Cho, Y., Park, Y. M., Barate, A. K., Park, S.-Y., Park, H. J., Lee, M. R., ... Holden, D. (2015). The role of rpoS , hmp , and ssrAB in Salmonella enterica Gallinarum and evaluation of a triple-deletion mutant as a live vaccine candidate in Lohmann layer chickens. *Journal of Veterinary Science*, 16(2), 187. http://doi.org/10.4142/jvs.2015.16.2.187



TBR4100 SPECIFICATIONS

POWER 100 ~ 240 VAC, 50-60 Hz, <15 W OPERATING TEMPERATURE (ambient) 0-50°C (32-122°F) OPERATING HUMIDITY (ambient) 15 - 70% RH non-condensing WARM UP TIME <5 minutes

DIMENSIONS

135 X 419 X 217 mm (5.25" X 16.5" X 8.16")

WFIGHT 1.35 kg (3 lb.)

DISPLAY FUNCTIONS 18 mm (0.7") LCD readout, 4.5 digit

Polarization Voltage (mV) Current input (nA, µA) Power (on/off) Current Input Range

Polarization Voltage ANALOG OUTPUT RANGE ±10 V (continuous)

ANALOG OUTPUT IMPEDANCE 10 KΩ CHANNEL TO CHANNEL ISOLATION >10 GΩ CHANNEL TO OUTPUT ISOLATION >10 GΩ >100 MΩ POWER SUPPLY TO AC LINE ISOLATION ANALOG OUTPUT DRIFT <10 pA/h

TEMPERATURE INPUT

CONTROLS

NUMBER OF CHANNELS

Platinum RTD, 1000 Ω SENSING ELEMENT

RANGE 0-100°C ACCURACY +1°(0.1°C RESOLUTION

31.25 mV/°C (continuous) ANALOG OUTPUT

AMPEROMETRIC INPUT

NUMBER OF AMPEROMETRIC CHANNELS	4
SIGNAL	0-3 Hz
BANDWIDTH	U-3 ПZ
POLARIZATION VOLTAGE (SELECTABLE VIA ROTARY SWITCH)	
Nitric Oxide	865 m\

150 mV Hydrogen Sulfide Hydrogen Peroxide 450 mV Glucose 600 mV 700 mV Oxygen ADJ (user adjustable) ± 2500 mV POLARIZATION VOLTAGE ACCURACY $\pm 5 \, \text{mV}$ POLARIZATION VOLTAGE DISPLAY RESOLUTION ±1 mV

CURRENT MEASUREMENT PERFORMANCE

Range	Analog Output	Noise @ 3Hz	Noise @ 0.3 Hz
± 10 nA	1 mV / 1 pA	< 1 pA	< 0.3 pA
± 100 nA	1 mV / 10pA	< 7 pA	< 3 pA
± 1 μA	1 mV / 100pA	< 70 pA	< 30 pA
+ 10 µA	1 mV / 1 µA	< 700 nA	< 300 nA

^{*}Instrument performance is measured as the (max-min) over 20-second period with open input. Typical values are given at 3 Hz and 0.3 Hz bandwidth.

Typical sensor performance with TBR4100: ISO-NOPF100 NOISE 0.2 nM NO (<2 pA)** **Sensor noise is measured as the (max-min) over a 20-second period with the sensor immersed in 0.1 M CuCl₂ solution.

Temperature Sensor



The temperature sensor (ISO-TEMP-2) is based on a 2.0 mm tip diameter high quality miniature platinum RTD (Resistance Temperature Detector) electrode. This design has been shown to provide greater accuracy, stability and interchangeability during temperature measurements than traditional thermistor and thermocouple sensors. The ISO-TEMP-2 is included with the purchase of a system.



TBR1025

Don't need four channels? The single-channel TBR1025 packs the power of the full-size 4-channel unit in a small, economical package.

ORDERING INFORMATION

TBR4100-416 Four-Channel Free Radical Analyzer with Lab-Trax 4/16 Data Acquisition System

Includes TBR4100 analyzer & power cord, Lab-Trax-4/16 data acquisition system & USB cable, 4 BNC cables, 1 electrode adapter cable, 1 temperature probe, 2 sensors of your choice, and sensor start-up kit(s), if applicable.

Single-Channel Free Radical Analyzer **TBR1025** Includes 1 sensor of your choice & 1 start-up kit

RECOMMENDED ACCESSORIES / REPLACEMENT PARTS

SNAP50 SNAP S-Nitroso-N-acetyl-D-penicillamine, 50 mg vial ISO-TEMP-2 2 mm Platinum RTD Temperature Sensor (requires 91580) 91580 Microsensor Adapter Cable

Optical Oxygen Sensor

Oxygen detection in small samples

Features

- Oxygen detection using the phase shift between a reference signal and measured signal
- Common sterilization procedure via gamma radiation, ethanol, paraacetic acid, ethylene oxide
- Easy cleaning procedure with 3% H₂O₂, ethanol or common soap solution
- Compatibility with ethanol, methanol and isopropanol
- No cross-sensitivity in static aqueous solutions with pH 1-14, CO₂, H₂S, SO₂ or any other ionic species
- Cross-sensitivity for organic solvents such as acetone, toluene, chloroform or methylene chloride and any chlorine gas
- Available with SMA or ST connector
- Small sensor tip diameter ≤ 50 µm

Benefits

- Designed for use in gaseous and aqueous solutions
- Easy calibration via a two-point calibration with an oxygen-free environment (nitrogen, sodium sulfite) and an air-saturated environment
- Fast response time t₉₀ ≤ 3 seconds (in gas phase)
- Retractable up to 10 mm
- No oxygen consumption during measurements

Applications

- · Oxygen measurement in small samples
- Bioprocess control

BioOxy is a new and innovative technology for measuring oxygen in gaseous and aqueous phase. **BioOxy** is an optical oxygen sensor with important advantages over using common Clark type electrodes.

- BioOxy does not consume oxygen while making measurements.
- Similar performance overall at a lower cost. This helps lower research or operational expenses compared to competitive consumable sensors.
- **BioOxy** may be gamma sterilized. Sensors can be irradiated, as well as ETOH and PAA, which may be ideal for biological applications.
- Very good resistance to long term photobleaching allows for extended use, a periodic re-calibration is recommended for optimal results.
- Sensor-to-sensor repeatability because of the very small unit-to-unit variation allows much more confident data from test to test using different sensors.

The BioOxy is ideally adapted for examination of small sample volumes thanks to its miniature sensor tip diameter of less than 50 μm . The BioOxy is well suited for long-term measurements and biotechnological applications. The signal does not depend on the flow rate of the sample and the measurement principle requires no contact.

The physical principle is based on quenching of luminescence, which is caused by the collision between molecular oxygen and a luminescent dye molecule in the excited state. There is a direct relationship between the measured phase angle and oxygen concentration. The phase angle decreases when oxygen is present in your sample. The transfer of the emitted light signal is done optically using high quality fused silica fibers, ensuring a constant performance from lot-to-lot.



BIOOXY SPECIFICATIONS

OXYGEN MEASUREMENT Gaseous & Dissolved O_2 PROBE SIZE 240 μ m with active area ~50 μ m Retractable 10 mm

MEASUREMENT TEMPERATURE

RANGE

RESPONSE TIME (T₉₀) < 3 sec. (in gas phase)

ORDERING INFORMATION

0 to + 80 °C

02-T0050L250-12012 BioOxy Sensor, 50 µm Tip, ST Connector, pack of 2 **02-T0050L250-12022** BioOxy Sensor, 50 µm Tip, SMA Connector, pack of 2

COMPATIBLE WITH:

 OXY-MICRO-AOT
 Fiber Optic Oxygen Meter with ST Connection

 FLOX-PROBE
 In Situ Oxygen Monitoring Kit with SMA Connection

Oxygen Monitoring Kit

Detect O₂ for real-time analysis

Features

- Fluorescence-sensing detector for optical sensors, a viable alternative to traditional chemical sensing devices
- Measure fluorescence lifetime, phase and intensity
- LED based photometry makes this unit affordable-Half the cost of comparable phase measurement systems
- · Excellent stability, extremely low drift and phase noise
- · Simple calibration, setup and control

Benefits

- Self-contained, benchtop system-invariant to fiber bending or stray
- Easy setup and control, combined with its stability and sensitivity make this unit perfect for long term studies
- Three coatings for probes and patches give you many options.

Applications

- Oxygen sensing applications requiring stability and sensitivity to drift which must be undisturbed for lengths of time
- Monitor O₂ partial pressure in gas and aqueous solutions and in nonaqueous vapors and solutions
- Monitor traces of oxygen in gas and liquids

FLOX is a device for measuring fluorescence lifetime, phase and intensity. It uses LED excitation and photodiode detection with filterbased wavelength selection for easy experimental set-up and control.

Self contained unit

The compact, self-contained unit makes it invariant to fiber bending and stray light. It also has a wide dynamic range of optical intensity, as well as low optical and electronic crosstalk, and low drift and phase noise.

FLOX is especially useful for oxygen sensing applications where stability and sensitivity to drift is important and where sample set-ups must be left undisturbed for long periods of time.

Easy setup, great for long term studies

The new oxygen sensing system measures fluorescence lifetime, phase and intensity, using LED excitation and photodiode detection with filterbased wavelength selection. The system is simple to set-up and control. When stability and sensitivity to drift are important in your oxygen sensing experiment, this unit is ideal. It is perfect for applications where sample set-ups must be left undisturbed for long periods of time.

Three coatings available for probes and patches

The optical sensors consist of transducer materials, applied to the tips of optical fibers or to substrates such as patches or cuvettes, which change optical properties in response to specific analytes in their immediate environment.

OXY-The standard oxygen sensor designed for monitoring oxygen partial pressure in gas and aqueous solutions is a fiber optic fluorescence probe with a proprietary oxygen sensing coated tip.

HIOXY–Designed for monitoring oxygen partial pressure in non-aqueous vapors and solutions. The sensor coating chemistry is compatible with oils, alcohols, and hydrocarbon-based vapors and liquids.

FOSPOR-A new generation of highly sensitive sensor coating can be used for monitoring traces of oxygen in gas and liquids.

The oxygen sensor probes are low-power and offer high sensitivity, reversibility and stability, ideal for remote monitoring. The thin coating on the probe tips consumes no oxygen, allowing for continuous contact with the sample. They are ideal for viscous samples and are immune to interference caused by pH, ionic strength or salinity fluctuations or



fiber, stainless steel ferrule; 3.175 mm OD, 63.5 mm length

R Sensor Probe

1000 µm fiber in a stainless steel 1/16" ferrule; 1.587 mm OD, 152.4 mm length

AL300

A 500 µm OD (300 µm core diameter) aluminum-jacketed optical fiber probe for applications that require fine spatial resolution. • 300 μm alμminμm-jacketed fiber assembly; 500 µm OD, 1 m length

The OR125 is a 1/8" OD optical fiber probe used as a direct replacement for 1/8" OD O2 electrodes. • 100 μm optical fiber, stainless steel ferrule; 3.175 mm OD, 63.5 mm length

P1600

The PI600 is a silicone-jacketed, polyimidecoated optical fiber probe used in environments where a non-metallic probe is required. • 200 or 600 µm optical fiber with silicone jacketing; 710 µm OD, 2 m length



ORDERING INFORMATION


FLOX-PATCH Non-Invasive Oxygen Monitoring Kit, including phase measurement system, temperature probe; select sensor patches when ordering In Situ Oxygen Monitoring Kit, including phase measurement

FLOX-PROBE

system, temperature probe; select sensor when ordering

Carbon Fiber Microelectrodes

Electrochemical detection of oxidizable compounds



Features

- Sensitive, renewable/durable and economical carbon fiber electrodes for electrochemical detection of oxidizable compounds
- Excellent linearity to the oxidizable compounds

Benefits

- Precision tip size and length of the CF electrodes
- Renewable

Applications

 Detection of neurotransmitters and oxidizable compounds in vivo or in vitro

Carbon fiber microelectrodes have been used in both the detection of oxidizable compounds (Gonon, et al., 1978; Cahill and Wightman, 1995) and extracellular single-unit recording (Armstrong-James and Millar, 1979). WPI's ultra-sensitive and low-noise carbon fiber (CF) electrodes can be applied, with our Micro-C Potentiostat or similar instruments, in the electrochemical detection of catecholamines (epinephrine, norepinephrine and dopamine), indolamines (serotonin, 5-HT or melatonin), ascorbic acid, Fe (II) and other oxidizable compounds. CF electrodes (diameter of 10 or 30 µm) respond with an excellent linearity to the oxidizable compounds (see figure below) and can detect the compounds as low as 0.2 nM. While the shorter (25-100 μ m) CF electrodes are suitable for in vivo amperometric and voltammetric measurements, the longer CF electrodes provide higher sensitivity and are especially useful for the in vitro studies (amperometric or differential pulse voltammetry). When used with the Micro-C Potentiostat, these CF electrodes can be activated and renewed

Dopamine Concentration/Response Curve




Fig. 1 — Excellent linearity in the response of carbon fiber electrode (CF30-500) to dopamine recorded on Micro-C. Courtesy: Drs. D. Yeomans and X.-T. Wang, University of Illinois at Chicago.

in sensitivity for multiple use. The selective detection of catecholamines can be achieved with our Nafion-coated CF electrodes. For selective detection of 5-HT and ascorbic acid, please contact WPI for more information.

References

P. S. Cahill, R. M. Whightman, Anal. Chem., 67, 2599-2605 (1995).
F. Gonon, et al., Hebd Seances Acad. Sci. Ser. 286, 1203 (1978).
M. Armstrong-James, J. Millar, J. Neurosci. Methods, 1, 279 (1979).




Fig. 2 — Extracellular recording using a carbon electrode in CA1 region of the hippocampus in an anesthetized rat shows ultra-low noise (<5 μV). Courtesy: Dr. Carolyn Harley of Memorial University, Newfoundland, Canada.

ORDERING INFORMATION

CARBON FIBER MICROELECTRODES, UNCOATED

	Diameter	Length	(pack of 5)
CF10-100	10 μm	100 µm	
CF10-250	10 μm	250 µm	
CF10-500 *	10 μm	500 µm	
CF30-50 *	30 µm	50 µm	
CF30-100	30 µm	100 μm	
CF30-500 *	30 µm	500 µm	
CF30-1000 *	30 µm	1000 µm	

CARBON FIBER MICROELECTRODES, NAFION-COATED

	Diameter	Length	(pack of 5)
CFN10-50 *	10 μm	50 μm	•
CFN10-100 *	10 μm	100 µm	
CFN10-250 *	10 μm	250 µm	
CFN30-50 *	30 µm	50 μm	
CFN30-100 *	30 µm	100 µm	
CFN30-250 *	30 µm	250 µm	
CFN30-1000 *	30 µm	1000 µm	

^{*} Built to order — allow up to 4 weeks manufacturing time.

Fiber Optic Oxygen Meter and Sensors

A generation of sensors based on luminescence lifetime

DXY MICRO

Features

- Oxygen is not consumed during the experiment
- Immune to electrical and magnetic interference
- Excellent long-term stability
- No lengthy polarization necessary (like Clark-type O₂ electrodes)
- Fast response time < 0.5 s for MicroTip sensors
- Probe size of MicroTip sensors as small as 50 μm
- Measurement is feasible in dry gas
- Optical isolation of sensor tip available for fluorescent or photosynthetically active samples

Benefits

- Unaffected by light source stability and intensity fluctuations, because it is based on luminescence lifetime detections.
- Compact, portable meters may be used inside or outside

Applications

- Process control like bottling plant in breweries and quality control of packages (OxyMini)
- Biotechnology like control of cell culture media and non-invasive control of bioreactors (OxyMini)
- Implantation of oxygen sensors into soil and trees (OxyMini)
- Oxygen profiles of marine sediment, soils, or tissue (OxyMicro)
- Implantation in living tissue like heart or muscle tissue (OxyMicro)
- Control of cell culture media in Biotechnology (OxyMicro)

The measurement principle of the sensor system is based on the detection of oxygen concentration as a function of luminescence lifetime, either in dissolved or gaseous phase environments.

Unaffected by light source stability/intensity

Conventional fiberoptic oxygen sensor systems based on intensity measurements are limited in their accuracy by light source stability and ambient light fluctuations. Using a luminescence lifetime detection, measurements are not affected by light source stability, intensity fluctuations caused by fiber bending or changes of the optical properties of the sample (turbidity, refractive index, coloration, etc.).

Calibration: The sensors can be calibrated by a simple two point calibration, 100% air-saturation and 0% air saturation.

Compact, portable meters, inside and outside

The **OxyMini** and **OxyMicro** fiber optic oxygen meters are compact and easy to transport. Designed for in/outdoor use, they can be connected to a PC via a RS232 interface. Data can be visualized, analyzed and stored with supplied software. A full range of sensors are available.

OxyMini systems

The OxyMini is a one channel fiber optic O_2 meter for fiber optic O_2 minisensors. These sensors are based on 2 mm polymer optical fibers and have a length of 2.5 m.

MiniTip – This dipping probe (501641) has a tip diameter of 4 mm and consists of a polymer optical fiber with an O_2 sensitive coating. Its range is 0–100%. It has a response time (t_{90})

MiniFlow – The MiniFlow probe (501642) is a miniaturized fiber optic chemical sensor integrated in a standard T-shape flow through cell which can be easily connected via Luer-Lock adapters to external tubings. Liquids can be pumped through the cell. It has a response time (t₉₀) of approximately 40 s and an excellent long-term stability.

MiniFoil – Sensor material on a 1 cm² support disk made of polyester can be glued inside glass vials. Measure oxygen concentration non-invasively and non-destructively from



outside through the wall of the vial. Illuminate the sensor foil with a plastic fiber optic cable (501644, 01645). The wall of the flask must be transparent/non-fluorescent (response time (t_{90}) of ~50s). The material can be implanted into animal tissues or custom-made housings.

OxyMicro systems

OxyMicro is a one channel O_2 meter for fiber optic O_2 microsensors.

MicroTip – The MicroTip (501656) is a needle-type (27 ga.) oxygen micro sensor designed for applications where a small tip



size $(50 \ \mu m)$ and fast response time (t_{90}) of 1 s are required. The oxygen sensitive sensor tip consists of 140 μm fiber tapered to a 50 μm tip. Housed inside a stainless steel needle 22 mm long and 0.4 mm diameter, it can penetrate a septum rubber or similar material. They are ideal for oxygen profiling in sediments and biofilms.

MicroFlow – The MicroFlow oxygen sensor (501657) is a miniaturized fiber optic chemical sensor optimized for



fast response time (t_{90} <1 sec in gases, <5 sec in liquids). The tiny probe has a tip size of 50 μ m and is integrated in a T-shape flow cell for easy connection via Luer-Lock adapters to external tubings. Liquids (like water, blood, etc.) can be pumped through the cell.

MicroImplant – The MicroImplant oxygen sensor (501658) is an implantable probe (IMP) with a probe



tip size 50 µm, an exposed fiber length of 5 mm and a jacket diameter of 900 µm. It has been successfully implanted in crabs, fishes and soil.

SPECIFICATIONS									
MiniTip MiniFlow MiniSpot MicroTip MicroFlow MicroImplant									
MEASURE RANGE DISSOLVED/ GASEOUS	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg	0-45 ppm, 0-100% 0-760 mmHg			
RESPONSE TIME [T ₉₀] DISSOLVED/ GASEOUS	40 s 10 s								
STERILIZATION ETOH, H ₂ O ₂ AUTOCLAVABLE*	Y N	Y Y	Y Y	Y N	Y Y	Y			
DRIFT**	< 0.1%	< 0.1%	< 0.1%	< 0.3%	< 0.3%	< 0.3%			
ACCURACY***	ACCURACY*** 0.2%								
RESOLUTION***	2.75 ±0.01 ppm, 9.00 ±0.05 ppm, 220 ±0.15 ppm, 45.0 ±0.25 mmHg, 150 ±0.75 mmHg, 375 ±2.6 mmHg								
TEMP RANGE	-10°C to 50°C								
PROBE ASSY LENGTH	2.5 m								

*130°C, 1.5 atm **100,000 data points, 20°C ***20°C

ORDERING INFORMAT	ΓΙΟΝ
MINICENCOD CYCTEM	

MINAISEIA	30K 3131 EWI
OXY-MIN	II-AOT Fiber-optic Oxygen Meter for Minisensors *
MINISEN	SORS (not interchangeable with Microsensors)
501641	MiniTip, fiber-optic oxygen sensor
501642	MiniFlow, flow-through cell with integrated planar oxygen sensor
503090	MiniSpot, planar oxygen-sensitive spot, 5 mm diam. (includes 10)
	Requires 501644
501644	Polymer optical fiber with 1 SMA connector
MICROSE	INSOR SYSTEM
OXY-MIC	RO-AOT Fiber-optic Oxygen Meter for Microsensors *

OXY-MICR	O-AOT Fiber-optic Oxygen Meter for Microsensors *
MICROSE	NSORS (not interchangeable with Microsensors)
501656	MicroTip, needle-type housing fiber-optic oxygen sensor, 50 μm tip
501656-C	MicroTip, needle-type housing, 50 μm tip, optical isolation
501656-F	MicroTip, needle-type housing, 140 μm flat tip
501657	MicroFlow, flow-through housed oxygen microsensor
501658	MicroImplant, implantable oxygen microsensor, 50 μm tip
501658-F	MicroImplant, 140 μm flat tip
*/	Meter contains two analog outputs and one trigger input

4-Channel Data Acquisition System with Software

Low noise, high resolution system with 8 analog input and 3 analog output channels

Features

- Powerful low-noise (<1 mV RMS) and high-resolution (16 bits) data acquisition system for sampling up-to 8 analog input channels and 3 analog output channels simultaneously, using standard **BNC** connections
- MDAC software provides easy to use interface controlling, with extensible standard and customized Data Processing and Analysis Tools



Benefits

- · Online Channel Math operations, general purpose Fast Fourier analysis (FFT) and digital filtering of Analog In channels
- Numerous basic signal forms can be combined to design experimental protocols, for most physiological applications
 - · Factory designed standard or customized protocols
 - Semi-automated data analysis toolbox
- Protocol repeat function to avoid time consuming protocol programming of extended experiments

Applications

- Muscle physiology (Can be used with SI-MTM Muscle Testing Platform, SI-CTS200 Cell Tester System, SI-HTB2 Horizontal Tissue Bath and SI-BF-100 Biofluorometer)
- Stand alone general data recorder for Spectroscopy, Neuroscience and Electrophysiology (Can be used with TBR4100 Free Radical Analyzer, Extracellular Bioamplifiers like SYS-DAM50, SYS-DAM80, SYS-900A, ISO-80, EVOM2™ Volt Ohm Meter, ATC2000 Animal Temperature Controller, BP-1 Blood Pressure Monitor or the BAT-12 Microprobe Thermometer)
- Instrument control for software triggered devices like A365/A385/A395 Constant Current Stimulators, MPS-2 Perfusion System, SYS-PV820/ SYS-PV830 Pneumatic PicoPumps, Duo 773 Intracellular amplifiers and the SYS-TBM4M Transbridge Transducer Amplifier (e.g. for FORT force

Knowledge of the physiological characteristics of muscle tissue can be useful to quantify beneficial or adverse effects of drug supply on muscle function in pre-clinical and toxicological studies, evaluating muscle dystrophies, training effects in sports and rehabilitation (disuse vs. overuse) and advanced physiology and biomedical research.

This is usually achieved by quantifying the contractile and/or the elastic properties of muscle tissue. This needs the programming of different and specific experimental protocols (isometric, concentric and eccentric, isokinetic or isotonic), so that the physiological structure of interest can be quantified. LabTrax-MDAC data acquisition software was designed for use with WPI's Muscle Physiology line to test physiological characteristics of muscle tissues in various conditions, using factory designed standard or customized protocols. The semi-automated Data Analysis Toolbox of standard protocols gives quick access to user-friendly, readable and interpretable results of the experiments.

Variety of muscle physiology applications

The physiological response of muscle tissue to training, disuse, nutrition, drug supply and others factors may be studied by adding accessories to the system, like:

• Study of the muscle's force production capacities in combination with the Ca²⁺ release from the sarcoplasmic reticulum (SR) and ATPase consumption. The perfect instrument for this is WPI's Biofluorometer (SI-BF-100) in combination with any system of WPI's Muscle Physiology line, controlled via LABTRAX-MDAC.

• Study of the muscle's force production capacities from direct muscle or peripheral nerve stimulation. For this experiment, use WPI's programmable isolated current stimulators (A365, A380 or A395), controlled via LABTRAX-MDAC.

LABTRAX-MDAC provides easily used continuous stimulation protocols, so that especially cardiac cells/tissue remain intact during experimental resting periods.

LABTRAX-MDAC is also well suited for other software triggered instruments or as a stand-alone general data recorder for selected WPI



The back panel of the Lab-Trax-8/16 has four analog outputs, digital inputs or outputs, a USB port, power socket and power switch.

LAB-TRAX-8/16 SPECIFICATIONS

ANALOG INPUTS 8 BNC connections **INPUT RANGE** ± 10V SYSTEM NOISE < 1 mV RMS ISOLATION 1,500V OPERATING CURRENT 800 mA maximum ANALOG OUTPUTS 4 BNC connections OUTPUT RANGE

cut-off frequency OUTPUT IMPEDANCE 1000

OUTPUT CURRENT

IMPLEMENTED FILTER

16/16 TTL (BNC or DB-9 Connector) DIGITAL I/O LOGIC HIGH VOLTAGE 3.3 V minimum

5th order low-pass Bessel filter with 3dB

LOGIC LOW VOLTAGE 1.0 V maximum ANALOG & DIGITAL INPUTS Operating voltage protected to ±30V

PC INTERFACE **USB 2.0**

RESOLUTION 16 bits **POWER SOURCE**

ORDERING INFORMATION

LABTRAX-MDAC Lab-Trax-8/16 with MDAC software LAB-TRAX-4 4-Channel General Data Acquisition System BNC to BNC Cable

Novel Fiber Optic pH System

Referenced measurements with single excitation

Features

pHOptica Meter

- Single-channel, compact, easy to transport fiber-optic meter for pH measurements with miniature sensors
- Two 12-bit, programmable analog outputs, with electrical isolation
- One external trigger input, with electrical isolation
- Computer with RS232 interface required for operation

pH Optical Sensors

- Immune to electrical interferences and magnetic fields
- Low drift
- High spatial resolution due to small tip size
- Measurement in very small sample volumes
- Additional optical isolation of the sensor tip is available for measurements in colored or photosynthetically active samples

Benefits

- User-friendly software saves and visualizes measured values
- Several pHOptica meters can be connected to one computer
- Temperature variation is recorded using a temperature sensor
- No reference electrode is needed

Applications

pHOptica Micro System

The pH Optica micro system is a single channel pH system for use with fiber optic micro sensors. The applications include:

- Penetration or implantation into living tissue (heart, muscle or animal blood vessels).
- Soil implantation for pH measurement.
- Implantation into customer-made housing.

TIP: To protect the small glass fiber tip against breaking, suitable housings and tubings around it, depending on the respective application, were designed.

Optical fiber

god to was

with transparent window

pHOptica Mini System

The pHOptica mini system is a one channel pH system for use with fiber optic mini sensors, foil and spot surface sensors for applications like:

- Non-invasive and non-distractive pH measurements from outside through flask walls (cell culture)
- Online pH monitoring by flow through cells
- Dipping probe pH measurements

pH Micro Sensors

- Tip size 140 micrometer
- Drift of 0.1 pH units for 2000 measurements (16 hours measurement in the 30 sec data update mode)



Implantable sensor—without any housings implantation into animal blood circuits; soil implantation; implantation in custom-made housings



Needle-Type Housing Sensor—the glass-fiber with its pH-sensitive tip is protected inside a stainless steel needle (18 ga.); fiber has to be extended during measurement; penetration through septum.

pH Mini Sensors

• OD of the dipping sensor is 4 mm

PH-OPTICA-MINI

pH**Optica**

- Sterilization of the pH sensor spots via gamma radiation
- pH mini sensor meter is based on 2 mm PMMA waveguides
- Drift of 0.1 pH units for 10,000 measurements (4 days measurement in the 30 sec data update mode).

pHOptica™ is a pH measuring system which uses fiber optic sensors and patented DLR technology. This method allows referenced measurements with single excitation to be implemented.

Two different housings and sensor spots (sensorfoils) are available.



POF Coated with a pH-Sensitive Foil— Small and robust pH dipping sensor; no reference electrode needed.

POWER SUPPLY



Flow-Through Cell with Integrated pH Sensor—On-line monitoring; can be easily connected via Luer-Lock adapters.

SPECIFICATIONS

DATA INTERFACE RS232

SAMPLE RATE 1 sample per sec

MEASURING PH RANGE 5 - 9

RESOLUTION (20°C) ± 0.03 (microsensors) ± 0.01 (minisensors)

RESPONSE TIME 1 in in.

DIMENSIONS 185 x110 x 45 mm

WEIGHT 630 g

ORDERING INFORMATION

100-220 VAC

MINISENSOR SYSTEM (cannot be used with microsensors)

PH-OPTICA-MINI Fiber Optic pH Meter for Minisensors, Foils and Spots
 503538 pH MiniTip, fiber optic pH sensor dipping probe, disposable (4 mm OD), pkg. of 3. Requires cable 503110.
 501644 Polymer Optical Fiber with 1 SMA connector
 503110 Fiber Optic Cable with 1 SMA connector

MICROSENSOR SYSTEM (cannot be used with minisensors)

PH-OPTICA-MICRO Fiber Optic pH Meter for Microsensors
 502124 pH MicroImplant, fiber optic pH implantable sensor (140 μm OD), pkg. of 3

Reference Electrodes

Low electrolyte leakage, stable potential with low resistance

Features

- Extremely low electrolyte leakage
- Stable, reproducible potential with low resistance
- · Chemically resistant to strong acids and bases

Benefits

• May be used with ion selective electrodes without contamination from the reference electrode

Applications

• Small volume, low salt concentration measurement (SDR)

Dri-Ref™ reference electrodes were developed by WPI to have extremely low electrolyte leakage properties, hence the name "Dri-Ref." In addition to this key feature, these electrodes exhibit stable and reproducible potential and low resistance. Stored in KCl when not in use, they have a long life expectancy.

May be used with ion selective electrodes


Although the internal filling solution contains KCI, the low fluid leakage means Dri-Ref may be used in combination with ion selective electrodes, including those for K^{\star} and Cl^{\star} , without significant contamination from the reference electrode.

The Dri-Ref electrodes are chemically resistant to strong acids and alkalines. Dri-Ref electrodes are not suitable for use in organic solvents. In addition, the long, thin FLEXREF may be easily manipulated to accommodate a difficult experimental setup.

SUPER-Dri-Ref – With a diameter of 2mm, SUPER-Dri-Ref *does not leak electrolyte at all.* Exhibiting the electrical stability of a classic flowing junction reference cell, this electrode exhibits low resistance and a stable half-cell potential essentially independent of sample electrolyte concentration. SUPER-Dri-Ref is ideal for small volume and low salt concentration measurements.

Micro-Reference Electrode – Only 450 µm in diameter and one inch long, the new **DRIREF-450** reference electrode can be used along with other sensors in space-restricted areas and very small sample volumes.

Luer-Tip Reference – The male luer fitting at the front of the **DRIREF-L** allows it to be easily connected to a female luer port (see WPI's luer fittings kit, page 219 to form a tight seal — a very convenient installation for a flow-through system.



L	KI-KEI	. ELEC	. I KUL	E 251	CIF	ICATIC	2אוי	
							DRIREF-	DRIREF-
	DRIREF-450	DRIREF-5	DRIREF-2	FLEXREF	SDR	DRIREF-L	5SH	2SH
LENGTH	2.54 cm	9 cm	13 cm	13 cm	9 cm	7.5 cm	3.5 cm	2.54 cm
DIAMETER	450 µm	4.7 mm	2 mm	1.5 mm	2 mm	Standard Luer	4.7 mm	2 mm
CONSTRUCTION	Coated Glass	Ероху	Isoplast™	Teflon™	PVC	Polypropylene	Ероху	Isoplast™
RESISTANCE (TYPICAL)	<5 K Ω	~500 Ω	~2.7 KΩ	~2.7 KΩ	<5 KΩ	~500 Ω	~500 Ω	~2.7 KΩ
ELECTROLYTE LEAKAGE (ML/HR)	_	~7.4×10 ⁻⁷	~5.7×10 ⁻⁸	~5.7×10 ⁻⁸	-	7.4×10 ⁻⁷	~7.4×10 ⁻⁷	~5.7×10 ⁻⁸
LEAD LENGTH				30 in. (76	cm)			
CONNECTOR	2 mm pin							
Filling Solution				KCl				

	ORDERING INFORMATION
FLEXREF	Flexible Dri-Ref, 1.5 mm diam.
DRIREF-2	Dri-Ref, 2 mm diam.
DRIREF-2SH	Dri-Ref, 2 mm diam. (Short)
DRIREF-5	Dri-Ref, 4.7 mm diam.
DRIREF-5SH	Dri-Ref, 4.7 mm diam. (Short)
SDR2	SUPER-Dri-Ref, 2 mm diam.
DRIREF-450	Micro-Dri-Ref, 450 μm diam.
DRIREF-L	Reference Electrode with Luer Tip
	-

Calcium Calibration Solutions

CALBUF-1

For use with calcium electrodes

A set of eight calcium buffers covering the range of concentration from 10^{-1} to 10^{-8} M Ca $^{++}$. Each buffer contains 20 mL of solution and enough potassium chloride to set the ionic strength to 0.1 M. Limited shelf life; use within 30 days.

Concentration: $1x10^{-1}$, $1x10^{-2}$, $1x10^{-3}$, $1x10^{-4}$, $1x10^{-5}$, $1x10^{-6}$, $1x10^{-7}$, $1x10^{-8}$ M at 20° C. Limited shelf life; use within 30 days.

ORDERING INFORMATION

CALBUF-1 Kit of 8 Calcium Buffer Solutions

CAI BUF-2

Use with calcium fluorescent indicators

CALBUF-2 is especially suitable for calibrating fluorescent Ca $^{++}$ indicators. It provides eleven buffer standards in the 10^{-4} to 10^{-8} M Ca $^{++}$ range, whereas other commonly used fluorescent Ca $^{++}$ indicators have the apparent $\rm K_d$ in the range of 100 to 300 nM. As with any ionic sensitive indicator, the sensitivity range of these indicators is about 1.0 log unit above and below the K_d. CALBUF-2 provides seven calibration points in this sensitivity range. It has an osmolarity of 0.305, which is isotonic with most mammalian cells.

Concentration: $1x10^{-8}$, $4x10^{-8}$, $1x10^{-7}$, $2.5x10^{-7}$, $5x10^{-7}$, $7.5x10^{-7}$, $1x10^{-6}$, $4x10^{-6}$, $1x10^{-5}$, $4x10^{-5}$, and $1x10^{-4}$ M at 20° C. Ionic strength: 0.150 M. 11 bottles, 20 mL each. Limited shelf life; use within 30 days.

ORDERING INFORMATION

CALBUF-2 Kit of 11 Calcium Buffer Solutions

Ion Selective Electrodes

Accurately measure calcium, potassium, hydrogen or TPP

Features

- Fast, accurate, economical
- Superior, stable PVC membrane
- Fast response
- 2mm diameter tips
- Interchangeable tip holder
- Each kit includes 3 electrode tips and MicroFil filling syringe

Benefits

- Inexpensive
- Use to measure various ions in biological media

Applications

• Detection of ions in vivo or in vitro for biological applications

These highly stable electrodes accurately measure calcium, potassium, hydrogen and TPP (Tetraphenylphosphonium) ion activity. Tips consist of 2 mm diameter plastic tubes sealed at one end with an ion-sensitive membrane. After filling with electrolyte solution, you can insert the tube into the holder and connect it to a pH meter. Tips and holders are interchangeable, so one tip may be replaced with another sensitive to a different ion. Replacing a tip takes less than a minute. Electrode tips normally last several months, when stored properly in saline solution. When replacement is necessary, only the tip needs be replaced.

Kwik-Tip electrodes are available separately and as kits. Each "KWIK" Electrode Holder kit includes a reusable holder and three removable tips. In addition to a 4-foot BNC cable and an electrolyte filling syringe; "TIP" Electrode Kits contain three electrode tips for a specific ion. A separate reference electrode, such as WPI's Dri-Ref™, is also required.



KWIKH-2 Holder & 3 Hydrogen Electrodes KWIKPOT-2 Holder & 3 Potassium Electrodes KWIKTPP-2 Holder & 3 TPP (Tetraphenylphosphonium) Electrodes HOLDERS AND REPLACEMENT TIPS KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)		ORDERING INFORMATION
KWIKH-2 Holder & 3 Hydrogen Electrodes KWIKPOT-2 Holder & 3 Potassium Electrodes KWIKTPP-2 Holder & 3 TPP (Tetraphenylphosphonium) Electrodes HOLDERS AND REPLACEMENT TIPS KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	COMPLETE	KITS
KWIKPOT-2 Holder & 3 Potassium Electrodes KWIKTPP-2 Holder & 3 TPP (Tetraphenylphosphonium) Electrodes HOLDERS AND REPLACEMENT TIPS KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	KWIKCAL-2	Holder & 3 Calcium Electrodes
KWIKTPP-2 Holder & 3 TPP (Tetraphenylphosphonium) Electrodes HOLDERS AND REPLACEMENT TIPS KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	KWIKH-2	Holder & 3 Hydrogen Electrodes
HOLDERS AND REPLACEMENT TIPS KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	KWIKPOT-2	Holder & 3 Potassium Electrodes
KWIK-2 Electrode Holder with BNC cable TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	KWIKTPP-2	Holder & 3 TPP (Tetraphenylphosphonium) Electrodes
TIPCA Calcium Electrode Tips (3) TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	HOLDERS A	ND REPLACEMENT TIPS
TIPH Hydrogen Electrode Tips (3) TIPK Potassium Electrode Tips (3)	KWIK-2	Electrode Holder with BNC cable
TIPK Potassium Electrode Tips (3)	TIPCA	Calcium Electrode Tips (3)
	TIPH	Hydrogen Electrode Tips (3)
TIPTPP TPP+ (Tetraphenylphosphonium) Electrode Tips (3)	TIPK	Potassium Electrode Tips (3)
	TIPTPP	TPP+ (Tetraphenylphosphonium) Electrode Tips (3)

KWIK-TIP ELECTRODE SPECIFICATIONS									
Part #	Recommended Min. Slope/ Concentration Selectivity Coefficients rt # Electrode Color Code Filling Solution Decade Range (-log)								
TIPCA	Calcium	Green	0.1 M CaCl ₂	28 mV	0.1 M - 10 ^{-6.75} M	Na+ 5.5, K+ 5.4, Mg++ 4.9			
TIPH	Hydrogen	Red	1 M Citric Acid, 0.01 M NaCl, pH 5.6	54 mV	pH 5.0 - 12	Na ⁺ 10.4, K ⁺ 9.8, Ca ⁺⁺ 11.1			
TIPK	Potassium	Yellow	0.1 M KCI	54 mV	0.1 M - 10 ^{-4.5} M	Na+ 4.0, Ca++ 3.9, Mg++ 3.0			
TIPTPP	TPP+	Purple	10 mM TPP+	54 mV	0.001 M - 10 ⁻⁴ M	K+ 6.0			

Liquid Ion Exchangers

Make micropipettes to record cellular concentrations

WPI's Liquid Ion Exchangers (LIX), for use with the **FD223A** Electrometer, allow intracellular measurements to be made for cations (hydrogen, potassium and calcium).

p - 1 - 1 - 1 - 1 - 1 - 1 - 1 - 1 - 1 -							
SPECIFICATIONS							
ION	H+	K+	Ca++				
CATALOG NO.	IE 010	IE 190	IE 200				
SELECTIVITY COEFFICE	ENTS*						
Na ⁺ Mg ⁺⁺	12.7 —	1.97 2.95	5.5 4.9				
K ⁺ Ca ⁺⁺	_	<u> </u>	5.4				
USEFUL pH RANGE SLOPE	2-10 56 mV	4-10 58 mV	4-10 28 mV				
LINEAR RANGE APPROX. EQUIV.	pH 4-12 —	pK 0-3 Corning 477317	pCa 1-7 ETH1001				

^{*}Selectivity Coefficients are expressed here as -log Kij or pKij.



When used in micropipettes to record cellular ion concentrations, consider using WPI's Duo 773 electrometer (channel A).

	ORDERING INFORMATION
IE010	Hydrogen Ion Exchanger (0.1 mL)
IE190	Potassium Ion Exchanger (1.0 mL)
IE200	Calcium Neutral Ion Exchanger (0.1 mL)

Multi-Port Measurement Chamber

4-port closed chamber for measurements of NO, O_2 , H_2O_2 & other species in cell culture, temperature stabilized

Features

- Four port (NOCHM-4) chamber accommodates WPI's 2 mm sensors for nitric oxide (ISO-NOP), oxygen (ISO-OXY-2), hydrogen peroxide (ISO-HPO-2) and WPI's KWIK-TIP ion selective electrodes in combination with WPI's 2 mm Dri-Ref™ reference electrodes
- Two additional top ports for injection of reagents using WPI's MicroFil™ syringe needles
- Temperature control through an external circulating bath
- \bullet The chamber can be used for nitric oxide and other species calibration at temperatures from 4-40 $^{\circ}\text{C}$

Benefits

- Closed chamber design greatly reduces the surface area of the solution exposed to air
- One top port and up to three side ports configuration provides adequate space for convenient sample and electrode manipulation

Applications

 Simultaneously measurement of free radicals such as NO, H₂O₂, H₂S, O₂ and other ions at controlled conditions for cultured cell, cell suspensions or biological media

References


Olson, K. R., et. al. (2017). Catalase as a sulfide-sulfur oxido-reductase: An ancient (and modern?) regulator of reactive sulfur species (RSS). Redox Biology, 12, 325–339. http://doi.org/10.1016/j.redox.2017.02.021

Zhou, D., et. al. (2017). Oxygen binding and nitric oxide dioxygenase activity of cytoglobin are altered to different extents by cysteine modification. FEBS Open Bio, 7(6), 845–853. http://doi.org/10.1002/2211-5463.12230

Liu, X., et. al. (2017). Cytoglobin regulates blood pressure and vascular tone through nitric oxide metabolism in the vascular wall. Nature Communications, 8, 14807. http://doi.org/10.1038/ncomms14807

NSA Pre-polarizer Keep extra NO sensors ready to use

Achieve a stable background current quickly. This small battery-powered device applies a potential to the NO electrode equivalent to the potential applied by the ISO-NO meter. A sensor, which has been connected to the activator, may be transferred to the meter for immediate use. **For use with all WPI NO electrodes.**



ORE	DERIN	IG IN	IFOR	MAT	ION

NSA-3 ISO-NO Activator



Santos, S. S., et. al. (2017). NO production and potassium channels activation induced by Crotalus durissus cascavella underlie mesenteric artery relaxation. Toxicon, 133, 10–17. http://doi.org/10.1016/j.toxicon.2017.04.010

DeLeon, E. R., et. al. (2016). A case of mistaken identity: are reactive oxygen species actually reactive sulfide species? American Journal of Physiology-Regulatory, Integrative and Comparative Physiology, 310(7), R549–R560. http://doi.org/10.1152/aipregu.00455.2015

	ORDERING INFORMATION
NOCHM-4	Four-Port Closed Chamber, for use with WPI's 2.0 mm elec-
	trodes (e.g., ISO-NOP and ISO-OXY-2, etc.)
NOCHM-P	Spare Plug-adapter for ISO-NOP nitric oxide electrode
800100-5	Spare Center Chamber Gasket (package of 5)

SNAP

S-Nitroso-N-acetyl-D-penicillamine

SNAP is a stable green crystalline S-nitrosothiol compound that mimics the action of nitric oxide in vivo. It has vasodilatory properties and has been shown to relax isolated bovine coronary artery rings by activating soluble granulate cyclase. This reagent also actuates apoptosis in mouse thymocytes and has been accounted for reversible inactivation of protein Kinase C. SNAP can be used for calibration of all WPI NO sensors.

M.W. 220.2 • Purity > 98% by NMR or TLC

	ORDERING INFORMATION
SNAP25	SNAP, 25 mg vial
SNAP50	SNAP, 50 mg vial
SNAP100	SNAP, 100 mg vial

GSNO


S-nitrosoglutathione

GSNO has been identified *in vivo* as a potential storage and transport vehicle for NO in the body. GSNO has been used in clinical trials to treat a form of preeclampsia and to prevent platelet aggregation. It also has considerable potential as an NO donor in medicine.

M.W.336.3 • C_{10} H₁₆N₄O₇S • Purity > 98% • Soluble in water or DMSO • Storage: -20°C

	ORDERING INFORMATION
GSNO-50	GSNO 50 mg vial
GSNO-100	GSNO 100 mg vial

Optical Measurement Systems



All you need for optical measurements

WPI now provides a full range of components and solutions, from proprietary qualified fibers to optical assemblies and fully integrated spectroscopy measurement systems. WPI knows the importance of making measurements where the science happens, and WPI has a full range of reflectance and transmission dipping probes, flowcells and long pathlength liquid waveguide capillary cells (LWCCs) with up to 500 cm pathlength that eliminate the need to take a sample to an instrument. Our Photonics Center of Excellence in Friedberg, Germany also allows us to build proven custom configurations quickly. Whether that is a new optical design for a probe, an optical assembly or a customized measurement system.

Solutions for All Your Optical Needs

Optical measurement is important to researchers. Whether it is fluorescence studies or measuring absorbance of trace compounds, WPI is your solution partner who has the tools you need, and the technical understanding to help you solve your application challenges.

Starting in 1993 with our first patented liquid waveguide capillary cell (**LWCC**) for precise absorbance measurement of very low concentrations, WPI now provides a full range of components and solutions, including proprietary qualified fibers, optical assemblies and fully integrated spectroscopy measurement systems.

WPI's Photonics Center of Excellence in Friedberg, Germany, allows WPI to build proven custom configurations quickly. The close collaboration with key universities such as the THM University of Applied Sciences gives WPI a leading edge.

Features

- High quality silica/silica fibers for DUV/UV/VIS/NIR
- High quality silica/polymer clad fiber for high NA and broad wavelength range
- High flexible PMMA fibers for VIS wavelength range, large core fiber and high NA
- Industry leading solarization resistance at DUV-UV wavelengths: 180-1200 nm
- UV-enhanced broad range: 260-2200 nm
- High laser damage resistance
- · Broad temperature range
- Easy customization by: core size, connector, jacketing, and configuration
- Assemblies manufactured with biocompatible materials (connectors, epoxies and furcation tubing)
- · Radiation resistance dose dependent
- Sterilizable by common methods

Benefits


- Robust fiber and assemblies for UV use. WPI's fiber resists degradation due to UV exposure better than virtually any fiber on the market. This means less changeouts and better data comparison run to run.
- Unique extended range fiber that covers the high UV to long NIR wavelength range
- Qualified fiber allows true "plug and play" replacement
- Qualified QVLUX fiber is ideal for OEM and applications requiring field service or where multiple instruments must be used for a battery of tests
- Wide selection of fiber materials and sizes to manufacture the optimum fiber assembly for your application
- WPI's fiber optic assemblies are uniquely designed for biology and physiology applications with small size and excellent photon efficiency
- Custom configurations for your specific application are easy to specify and quick to deliver
- Customer support with in-house engineering and manufacturing for optical assemblies
- For critical biological work, many of our assemblies can be sterilized with gamma radiation, H₂O₂, ETOH or autoclave. Ask for details.

The QVLUX proprietary fiber series is WPI's most advanced fiber product line. QVLUX-DUV silica fibers are very resistant to UV light down to 185 nm but also have the uniformity and performance that OEM applications require. QVLUX-NIR fibers perform exceptionally well in the 700 to 2200 nm wavelength region. WPI is the only company offering DIN 58145:2018-01 certified low solarization UV fibers and assemblies for 185 m to 340 nm. For less demanding applications, we offer our standard WLUX fiber in both silica and polymer cores. While WPI offers QVLUX and WVLUX assemblies in many standard configurations, custom configurations can be quickly manufactured in house.

Basics of using optical fibers in science

Different than the fiber used for telecommunications, scientific optical fibers efficiently transmit light at many wavelengths. The amount of light that can enter a fiber and be transmitted is determined mostly by three factors:

- · Core diameter size
- Material of the core
- Material of the cladding.



The difference in the refractive index between core and cladding sets the numerical aperture, i.e. the maximum angle that light can enter the fiber and transmits through the fiber.

In addition, fibers are specified by wavelength region: typically, UV-VIS (190-800 nm) or VIS-NIR (400-2200 nm). For sensing and light delivery applications we offer fused silica, polymer cladded fused silica and solid polymer fibers, whichever are most suitable for your application. Silica fiber has a lower numerical aperture, but more efficiently transports light, especially at UV wavelengths. Polymer core fiber has a higher numerical aperture allowing more light to enter the fiber and may be a better choice for normal temperature conditions in a laboratory, especially in fluorescence detection.

Qualified fibers for UV applications

Our **QVLUX** fiber has been engineered to not only provide superior transmission in the UV but have excellent consistency over length. This makes **QVLUX** a perfect choice for critical measurements or where "plug-and-play" capability is required without having to recalibrate after a fiber assembly change. Standard silica fiber exhibits a "solarization impact" when transmitting light in the UV below ~280 nm. Over time, depending on the optically transmitted power, the fiber's ability to transmit UV light drops off. Our QVLUX fiber offers superior resistance to solarization, lowest basic attenuation at 200 nm (< 1 dB/m) and the lowest solarization at 214 nm on the market. WPI is the only company offering DIN 58145:2018-01 certified low solarization UV fibers and assemblies for 185 nm to 340 nm. Our standard WVLUX fused silica optical fibers are designed for value, but keep in mind that they lose their transmission over time and at high power due to the solarization effects of UV light below 280 nm.

Qualified fibers for NIR applications

We also offer **QVLUX** fiber for the 700-2200 nm range that works well in demanding petro-chemical, food and pharmaceutical applications. This region is very difficult for standard fibers from traditional suppliers since they cannot control their processes or glass chemistry sufficiently to have the purity required for high-standard of uniformity and superior optical throughput that these applications demand.

Match your Fiber to your Needs

With WPI you can match the fiber, the assembly and even the instrument exactly to your application needs:

Fiber	QVLUXDUV	QVLUXNIR	WVLUXDUV
FIBER MATERIAL	Qualified Silica	Qualified Silica	Specialty Silica
WAVELENGTH RANGE (nm)	180-1200	260-2200	200-1200
AVAILABLE CORE DIAMETERS (µm)	400, 600	200, 400, 600	100, 200, 300, 400, 600
KEY SPECIFICATIONS	 Lowest attenuation @ 180 nm Lowest solarization @ 214 nm NA = 0.22 @ 633 nm 	• Low solarization @ 248 nm • NA = 0.22 @ 633 nm	Low solarization @ 214 nmNA = 0.22 @ 633 nm



Fiber	WVLUXPV	WVLUXUVIS	WVLUXNXR	WVLUXHNA
FIBER MATERIAL	Silica/PolyClad	Silica	Silica	PMMA
WAVELENGTH RANGE (nm)	300-1200	260-1200	320-2200	365-720
AVAILABLE CORE DIAMETERS (µm)	200, 300, 400, 600, 800, 1000	100, 200, 300, 400, 600, 800, 1000	200, 400, 600	250, 500, 750, 1000, 1500, 2000, 3000
KEY SPECIFICATIONS	• 365 nm for low power UV	• Great general purpose fiber	Extended NIR range	 Solid polymer fiber (PMMA)
	applications	for UV-VIS spectroscopy	• NA=0.22 @ 633 nm	 Highly flexible
	• NA = 0.37 @ 633 nm	• NA = 0.22 @ 633 nm	 Custom core diameters available 	 Largest range of core diameters

Fiber Selection Guide

WPI offers standard, bifurcated and cross fiber assemblies. For Y assemblies and X assemblies, the split point is approximately at 50 cm. The standard jacketing is silicone monocoil.



1x1 Standard – S Optical fibers with 2 connectors and furcation tubing



1x2 Bifurcated – Y Split or combines similar intensities by mixing different or same fiber type




2x2 Cross – X Cross fiber assembly with mixed fiber types (call for details)

Ordering assemblies with our fiber is simple! Standard lengths are 100 cm and terminated with SMA connectors.

- 1 Choose fiber type.
- 2 Choose style.
- 3 Choose diameter (in microns).
 The available diameters are limited by the fiber type chosen in step 1.
 Refer to the tables above for the options.
- **4** Select length (in centimeters)

 If your assembly uses all the same fibers, start with the standard part number defined above. The standard fiber is 100 cm.
- **5** Choose connector type.

For example, the part number for a standard QVLUX-UV fiber with a Y assembly and 600 µm core diameter would be **QVLUXDUV-Y-0600**. For The part number for a custom cable of 250 cm with FC connectors using 750 µm high NA polymer fiber would be **WVLUXHNA-S-0750-250-F125**. For a mixed fiber, mixed connectors, X-cross assemblies or a special configuration, contact WPI (wpi@wpiinc.com).



Standard Option

Core diameters available depend on fiber chosen.

High Performance Spectrophotometer

A unique multiple long pathlength sample cell for absorbance spectroscopy



Features

- Process Control & Oceanography
- · Rugged system for laboratory and onboard measuring
- Portable & easy to use
- User-selected optical path lengths: 2, 10, 50 & 200 cm
- Highly sensitive and stable

Benefits

- Designed with NASA for colored dissolved organic matter in seawater and fresh water
- Highly sensitive with extended dynamic range for UV and VIS absorbance measurements, with fewer complications associated with standard long pathlength systems
- Portable system for field operations and mobile lab environments

Application

- · CDOM Colored dissolved organic matter
- QFT Quantitative Filter Technique

UltraPath™ is a unique high-performance spectrophotometer system offering selectable optical path lengths of 2, 10, 50 and 200 cm. The instrument operates in the wavelength range of 250 to 730 (**UPUV**) or 380 to 730 nm (**UPVIS**) and has an exceptional dynamic range. Designed for the detection of low absorbing species in aqueous solutions, **UltraPath** is an ideal tool for any study requiring precise and highly sensitive spectroscopic determination of analytes, either in the lab or in the field.

Designed with NASA for CDOM

UltraPath was developed by WPI under a collaborative agreement with NASA (Stennis Space Center) for the spectroscopic determination of colored dissolved organic matter (CDOM) in seawater and fresh water environments. It can be used in the laboratory and in the field (i.e., at sea). CDOM concentrations vary significantly between open ocean samples with low CDOM (e.g., 0.007 m-1 at 380 nm), and high CDOM freshwater environments (e.g., 10-20 m-1 at 380 nm). To address these problems, the design requirements of **UltraPath** mandated the development of a rugged portable system capable of high sensitivity measurements across a wide dynamic range. The **UltraPath** system meets these stringent design criteria and enables reliable measurement of CDOM in the range of 0.002 m-1 to 200 m-1 (250 to 730 nm).

Highly Sensitive, Extended Dynamic Range

UltraPath has four optical pathlengths contained within a single sample cell (i.e., 2 cm, 10 cm, 50 cm and 200 cm). The pathlengths are selectable, offering a very high sensitivity and an extended dynamic range for UV and VIS absorbance measurements. The fluid path of the sample cell is optimized to produce a laminar flow that is virtually free of interference from trapped air bubbles and adherence of dissolved substances to the cell wall. In particular, the design greatly minimizes the problems commonly found with flow cells of long optical pathlengths.

- · Reduces the risk of trapping dust particles.
- Reduces contamination of fibers or particulate matter inside the cell.

The UV/VIS UltraPath system (UPUV) includes a low noise photodiode array-based **TIDAS S-300 UV/VIS** spectrophotometer with a deuterium/ halogen light source. The VIS UltraPath system includes a low noise photodiode array-based TIDAS E-Base UV/VIS spectrometer module and a FO-6000 visible light source. Light is coupled from the light source to the sample cell and from the sample cell to the detector using two fused silica fibers. A peristaltic pump (PeriPro-4LS) draws the sample into the UltraPath sample cell.

A standard PC or laptop (not included) is connected to the TIDAS units using an RJ-45 Ethernet interface.

Portable system for field environments

The system is designed for mobility. The components of the **UltraPath** system are designed to function over a broad range of laboratory and field environments.



The Mini Star is a compact, lightweight peristaltic pump that fits just about anywhere. It can be mounted directly on the bench, in a regular rack or to a post. See page 53.

References

M. G. Novak, A. Mannino, N. Nelson, E. D'Sa, R. Miller, J. Werdell, R. Del Vecchio, C. Del Castillo, J. Chaves, J. F. Berthon, E. Boss, M. Tzortziou, A. Neeley, S. Freeman, A. Bricaud, R. Röttgers, A. Matsuoka, M. Belz, N. Blough, Measurements of CDOM absorption spectra using different Instruments and techniques: A round robin exercise and extensive field data set., Ocean Optics XXII, (2014)

N. B. Nelson, D. A. Siegel, C. A. Carlson, C. Swan, W. M. Smethie Jr. and S. Khatiwala. 2007. Hydrography of chromophoric dissolved organic matter in the North Atlantic. *Deep-Sea Res.* I. 54: 710 – 731.

V. Kitidis, A. P. Stubbins, G. Uher, R. C. Upstill Goddard, C. S. Law, E. M. S. Woodward, "Variability of chromophoric organic matter in surface waters of the Atlantic Ocean", *Deep Sea Research Part II: Topical Studies*, Vol. 53, Issue 14-16, 2006, p. 1666-1684.

R. L. Miller, M. Belz, C. Del Castillo, R. Trzaska, "Determining CDOM Absorption Spectra in Diverse Coastal Environments Using a Multiple Pathlength, Liquid Core Waveguide System", *Continental Shelf Research*, July 2002, 22:9, p 1301-1310.

"System Analyzes Water Samples at Sea", NASA Aerospace Technology Innovation, 2001, 9 (5). http://nctn.hq.nasa.gov/innovation/innovation95/3-techtrans2 html

R. L. Miller and E. D'Sa. "Evaluating the influence of CDOM on the remote sensing signal in the Mississippi River Bight". In *Eos Transactions AGU Ocean Sciences*, 2002. Honolulu, HI, p. 171.

E. D'Sa, R.L. Miller and R. Trzaska. "Aparent Optical Properties in Waters Influenced by the Mississippi River", *Proceedings of the Seventh Thematic Conference, Remote Sensing for Marine and Coastal Environments*, 2002, 6 pg, Miami, FL.

R. L. Miller, C. Hall, C. Del Castillo, B. McKee and M. Dagg. "Bio-optical Properties of the Mississippi River Plume and Adjacent Shelf." *ASLO Aquatic Sciences*, Albuquerque, NM, 2001.

R. L. Miller, M. Belz and S. Y. Liu, "Measuring the absorption of CDOM in the field using a multiple pathlength liquid waveguide system", *Ocean Optics XV*, paper 1308, Monaco, October 2000.

ULTRAPATH SPECIFICATIONS

WAVELENGTH RESOLUTION (FWHM) 5 nm

NOISE (PEAK TO PEAK) < 0.2 mAU

DRIFT < 1 mAU/h

OPTICAL PATHLENGTH 2, 10, 50 & 200 cm (user selectable)

SAMPLE CELL INNER DIAMETER 2 mm

CELL VOLUME 10 mL (at 200 cm pathlength)

SAMPLE INLET / OUTLET 1/8"
FIBER INPUT/OUTPUT 600 µm

SOLVENT RESISTANCE Most organic and inorganic solvents

SHIPPING WEIGHT UPUV: 44 lb. (20 kg)
UPVIS: 33 lb. (15 kg)

ORDERING INFORMATION

UPVIS Ultrapath System, Visible Light

UPUV Ultrapath System, Ultraviolet & Visible Light

The UltraPath system includes Multiple pathlength cell, Tidas S300 (UPUV) with integrated light source or Tidas E Base (UPVIS) with TidasDAQ3 software, two WVLUXUVIS-S-600 optical fibers, PeriStar Pro peristaltic pump, silicone tubing, sample injector and Waveguide Cleaniing Kit.

501609 Waveguide Cleaning Kit

89575 QFT1, Fiber Optic Holder for Glass Fiber Filters

WLUXUVIS-S-1000 Fiber Optic Cable 1 m, SMA, 1000 µm core, UV enhanced

Applications

Simple measurements for particulate absorption

WPl's filter holder for particulate absorption measurements is specially designed for field use. It is rugged and portable. It performs as well with a laboratory-based spectrophotometer. It can be directly connected to WPl's line of fiber optic spectrometers and light sources. Improve your particulate absorption measurements on site, instead of loss of spectral information by transporting your sample to a laboratory (Sosik, Limnol Oceanogr. Vol 44, 1139-1141, 1999).

How does it work?

Particulate absorption of fresh and seawater can be determined by filtering a known amount of sample through a Glass Fiber Filter (GF/F) and measuring the particulate absorption coefficient concentrated on the filter. This technique is called quantitative filter technique (QFT) and corrects for pathlength amplification, an effect of scattering (Mitchell, SPIE Vol 1302, 137-148, 1990).

Advantages of QFT

WPI's filter holder for glass fiber filters is designed for particulate absorption or cuvette measurements in a single instrument. A further advantage of the filter holder is its large beam diameter of 5 mm, resulting in "averaging out" of larger non-organic particles frequently found on the filter pad when using natural samples. The removable filter fixture allows simple filter alternation and cleaning.




QFT for particulate absorption

Particulate absorption can be measured by the well established Quantitative Filter Technique (QFT). WPI now offers a fiber optic filter holder for Glass Fiber Filters (QFT1, page 168) which can be used with the spectrometer (Tidas E Base)



and light source (**D4H** or **FO6000**) supplied with the **UltraPath**. With this accessory, particulate absorption can be measured on site, avoiding loss of spectral information due to freezing and shipping particulate samples to a laboratory.



Ponceau S absorption measured with UltraPath (200 cm cell). Ponceau S was dissolved in Millipore water.

Special attention should be drawn to the exceptional sensitivity of UltraPath enabling detection of CDOM absorption below 0.03 m $^{\rm 1}$. To exemplify the performance of the UltraPath in laboratory chemistry and process control, Ponceau S absorbance was measured with the 200 cm pathlength of an UltraPath. Normalizing the Ponceau absorbance graph to AU/cm, the range of this measurement is 150 μ AU with a noise level below 2 μ AU peak to peak. Sub-nanomolar concentration of this dye can clearly and reliably be detected, which is a novelty in absorbance based spectroscopy.

Liquid Waveguide Capillary Cell

Long pathlengths for small sample volumes

Features

- 10-500 cm pathlength
- 10–500 fold sensitivity improvement in comparison to 1cm cuvette
- 0.55 mm ID for low sample volume sampling
- 2 mm ID for unfiltered liquid samples
- SMA 905 fiber optic connections
- 250 nm 720 nm wavelength range with MilliPore water (see datasheet for details)

Benefits

- Adapts to most fiber optic detection systems
- 20 years of manufacturing experience
- Low UV drift

Applications

- Trace detection of nutrients (nitrite, nitrate, phosphate, iron) in seawater
- Environmental and oceanographic monitoring
- Drinking water analysis
- Colored dissolved organic matter (CDOM)
- Process control

UV/VIS/NIR absorbance spectroscopy is governed by Beer's Law, where the absorbance signal is proportional to chemical concentration, light path length and the compound's specific molar absorption coefficient. Typical optical pathlengths of cuvettes and flow cells are between 0.1 cm and 10 cm. Longer pathlengths are difficult to achieve due to mechanical constraints. Liquid Waveguide Capillary Cells (LWCCs) fill this gap. LWCCs are fiber optic flow cells that combine an increased optical pathlength (10–500 cm) with small sample volumes ranging from 2.4 μL to about 3 mL. Compared with a standard 1 cm cell, a 1 mAU signal is enhanced one hundred fold with a 100 cm flowcell to 100 mAU, using WPl's patented aqueous waveguide technology.* They can be connected via optical fibers to a spectrophotometer with fiber optic capabilities. Ultra-sensitive absorbance measurements can be performed in the ultraviolet (UV), visible (VIS) and near-infrared (NIR) to detect low sample concentrations in a laboratory or process control environment.

Your sample is the core of a light guide

WPl's Liquid Waveguide Capillary Cells are made of fused silica tubing with an outer coating of a low refractive index polymer. Your liquid sample is guided through the capillary and represents the core of the waveguide. The hydrophilic character of the fused silica capillary inner wall results in high signal stability and easy removal of air bubbles trapped in the flow cell. However, the transmission of the **LWCC** is mainly dependent on the intrinsic attenuation of the sample liquid. In case of water, a usable wavelength range from 250 nm to 720 nm wavelength can be observed in a 100 cm pathlength **LWCC**. Using a 500 cm pathlength **LWCC** will reduce that transmission range from 300 nm to about 700 nm. However, when switching from water to methanol as a solvent, transmission into the NIR are possible with suitable light sources and detectors.

Connections

The **LWCC-3xxx** series of flow cells uses traditional HPLC type 10-32 coned port fittings with 1/32 inch tubing for liquid connection and 500 μm SMA fiber optic adapters for light input and output. The **LWCC-4xxx** series of flow cells uses 1/4-28 flangless flat bottom fittings with 0.125" tubing and 600 μm SMA fiber optic adapters. Liquid can be pumped into the flow cells using (in the simplest case) a sample injector (**58006**) and a ministar peristaltic pump (**MINISTAR**). The **LWCC** may be connected directly to a fluid injection analysis (FIA) system or to a gas segmented fluid injection analysis (GFIA) system via a debubbler. Finally, for routing discrete measurements, WPl's LWCC Injection system (**89372**) may be used where the sample is injected into a constant flow via an injection loop of 3–4 times the internal flow cell volume to ensure a stable baseline and avoid the introduction of micro air bubbles into the flow cell.

Applications

When the entire spectral shape of an absorbance curve is required for analysis, WPI's TIDAS E Base spectrometer with a **D4H** or a **FO-6000**, or the TIDAS S300 spectrophotometer can be used. **LWCC**s have



been used in a variety of applications, such as liquid chromatography, stopped-flow and colormetric detection, drinking water analysis, as well as environmental and oceanographic monitoring systems.

Accessory: LWCC Injection System

For flow analysis, including simple fluid injection analysis (FIA) setups, add WPI's **LWCC** injection system (**89372**). A selection valve provides baseline or cleaning solutions to the sample stream. The injection valve injects a sample into the stream, avoiding the introduction of air bubbles or changes of flow rate.



89372 LWCC Injector System (pump and LWCC not included)

Related Patents

Micro Chemical Analysis Employing Flow Through Detectors, 1995, U.S. Patent No. 5.444.807.

Aqueous Fluid Core Waveguide, 1996, U.S. Patent No. 5,507,447. Long Capillary Waveguide Raman Cell, 1997, U.S. Patent No. 5,604,587. Chemical Sensing Techniques Employing Liquid-Core Optical Fibers, U.S. Patent No. 6,016,372

References

Yuhan Liu, Keding Lu, Huabin Dong, Xin Li, Peng Cheng, Qi Zou, Yusheng Wu, Xingang Liu, Yuanhang Zhang, (2016) In situ monitoring of atmospheric nitrous acid based on multi-pumping flow system and liquid waveguide capillary cell. *Journal of Environmental Sciences Volume 43*, 2016, Pages 273–284 http://dx.doi.org/10.1016/j.ips.2015.11.034

Ye, C., Zhou, X., Pu, D., Stutz, J., Festa, J., Spolaor, M., ... Knote, C. (2016). Rapid cycling of reactive nitrogen in the marine boundary layer. *Nature*, 532(7600), 489–491. http://doi.org/10.1038/nature17195

Catelani, T. A., Tóth, I. V., Lima, J. L. F. C., Pezza, L., & Pezza, H. R. (2014). A simple and rapid screening method for sulfonamides in honey using a flow injection system coupled to a liquid waveguide capillary cell. *Talanta*, 121, 281–7. http://doi.org/10.1016/j.talanta.2013.12.034

LWCC SPECIFICATIONS							
	LWCC-3050	LWCC-3100	LWCC-3250	LWCC-3500	LWCC-4010	LWCC-4050	LWCC-4100
OPTICAL PATHLENGTH	50 cm	100 cm	250 cm	500 cm	10 cm	50 cm	100 cm
INTERNAL VOLUME	125 μL	250 μL	625 µL	1250 μL	0.31 mL	1.57 mL	3.1 mL
FIBER CONNECTION	500 µm SMA	500 µm SMA	500 µm SMA	500 µm SMA	600 µm SMA	600 µm SMA	600 µm SMA
TRANSMISSION @254nm*	20	10	5	-	5	3	2
TRANSMISSION @540 nm*	35	30	25	20	5	4	4
NOISE [mAU]**	<0.1	<0.2	<0.5	<1.0	<0.1	<0.2	<0.2
MAXIMUM PRESSURE	100 PSI	100 PSI	100 PSI	100 PSI			
WETTED MATERIAL PEEK, Fused Silica, PTFE							

LIQUID INPUT LWCC-3xxx series: 10-32 coned port fittings for 1/16" tubing

LWCC-4xxx series: 1/4-28 flangeless flat bottom fitting for 1/8" tubing

* Referenced using coupled 500 µm fibers or LWCC-3xxx series and 600 µm fibers for LWCC-4xxx series

^{***} A one-meter waveguide of 550 µm internal diameter requires approximately 1.5PSI for water flow of 1.0 mL/min.



Waveguide Cleaning Kit (#501609), above, includes the most commonly needed cleaning solutions for the LWCC waveguides. The LWCC Start-up Kit (#KITLWCC), at right, includes two fiber optic cables (#WVLUXUVIS-S-600), Sample Injector Assembly (#58006), MiniStar™ Peristaltic Pump, and Waveguide Cleaning Kit (#501609).

	ORDERING INFORMATION
LWCC-3050	Liquid Waveguide Capillary Cell, 50 cm pathlength, 0.55 mm ID
LWCC-3100	Liquid Waveguide Capillary Cell, 100 cm pathlength, 0.55 mm ID
LWCC-3250	Liquid Waveguide Capillary Cell, 250 cm pathlength, 0.55 mm ID
LWCC-3500	Liquid Waveguide Capillary Cell, 500 cm pathlength, 0.55 mm ID
LWCC-4010	Liquid Waveguide Capillary Cell, 10 cm pathlength, 2 mm ID
LWCC-4050	Liquid Waveguide Capillary Cell, 50 cm pathlength, 2 mm ID
LWCC-4100	Liquid Waveguide Capillary Cell, 100 cm pathlength, 2 mm ID

Accessories

A sample injector assembly can be used to conveniently fill an LWCC with sample solution using a peristaltic pump. The LWCC requires two optical fibers to connect to spectrophotometer system. Choose between anti-solarized 400 µm core or UV-enhanced cables (ordered in 1 or 3 meter lengths).

89372	LWCC Injection System
58006	Sample Injector Attachment
PERIPRO-4LS	Peri-Star™ Pro Peristaltic Pump (See page 52.)
MINISTAR	Miniature Peristaltic Pump, 1-channel (See page 53.)
WVLUXUVIS-S-60	NO Fiber Optic Cable, 1 m, SMA, 600 μm core,
	UV-enhanced
501609	Waveguide Cleaning Kit (available only in USA)
KITLWCC	LWCC Start-up Kit*
58450	Kit, Adapter Syringe, LWCC

*includes WVLUXUVIS-S-600 (two), 58006, MINISTAR, 501609

Low Volume Flow Cell

For FIA, HPLC and process analysis

Features

- UV/VIS flow cell for absorbance
- Low internal volume
- Fits 500 and 600 µm fibers
- High UV transmission
- 0–10 mL/min flow rate

Benefits

- · High efficiency coupling
- · Low refractive index offset
- Fits WPI **TIDAS** systems

Applications

- FIA, GFIA, HPLC, Optofluidics
- · Process control

MicroLWCC is a new fiber optic low volume flow cell for UV/VIS/ NIR absorbance analysis. Based on WPI's established liquid core waveguide technology, the analyte solution functions as the core of a fluid filled light waveguide. Wetted parts in the sample cell light path are PEEK, fused silica and PTFE. Optical fibers are used to transport light to and from the sample cell. The cell can be used in biochemistry for DNA, RNA & protein quantification, colorimetric nutrient and trace metal analysis, drug discovery and dissolution testing, process control, and HPLC analysis.

LWCC-M SPECIFICATIONS						
LWCC-M-10	LWCC-M-50	LWCC-M-100				
10 mm	50 mm	100 mm				
2.4 μL	12 µL	24 µL				
< 15 mAU	< 20 mAU	< 30 mAU				
25%	20%	15%				
40%	35%	30%				
WAVELENGTH RANGE 200 – 1000 nm						
BER CONNECTION [µm] 500 (SMA)						
50 Bar						
PEEK, Fused Silica, PTFE						
	10 mm 2.4 μL < 15 mAU 25% 40%	LWCC-M-10 LWCC-M-50 10 mm 50 mm 2.4 μL 12 μL < 15 mAU				

^{*} Reference: 2 * 600 µm Fiber, butt-coupled

WPI U.S. Patents: 5,444,807; 5,570,447; 5,604,587; 6,603,556; 6,385,380.


ORDERING INFORMATION				
LWCC-M-10	Low Volume Flow Cell, 10 mm pathlength			
LWCC-M-50	Low Volume Flow Cell, 50 mm pathlength			
LWCC-M-100	Low Volume Flow Cell, 100 mm pathlength			

^{**} Measured using ASTM E685-93

^{**} Measured using ASTM E 685 - 93

Measure Colored Dissolved Organic Matter

Detect lower concentrations of solutes & conduct broader range of measurement



Colored Dissolved Organic Matter (CDOM) is organic matter whose optical properties are measurable using WPl's Liquid Waveguide Capillary Cell (LWCC). CDOM occurs naturally in water systems and is derived from organic tannins. CDOM concentration depends on the location where samples are taken, with coastal waters showing higher CDOM concentrations compared to open-ocean waters. In addition, CDOM absorption depends on open-ocean water depth.

LWCCs are fiber optic flow cells that combine an increased optical pathlength (1–500 cm) with small sample volumes ranging from 2.4 μL to about 3 mL.

Features

Select models:

- 10–250 cm pathlength (depending on LWCC model)
- 0.31–3.1 mL internal volume (depending on LWCC model)
- Up to 230–730 nm wavelength range (depending on LWCC model)

All models:

- 10 to 500 fold sensitivity improvement in comparison to 1 cm cuvettes
- 0.55 mm ID for low sample volume sampling
- 2 mm ID for filtered and unfiltered liquid samples
- SMA 905 fiber optic connections

Benefits

- Improved dynamic range for broader range of absorbance measurements
- Improved sensitivity of measurement
- · Detect low concentrations of solutes
- · Measurements can be made using small sample volumes
- Compact, portable system for real time measurement aboard ship

Determining LWCC pathlength

Selecting the proper equipment is imperative when setting up your system. Here are a few considerations:

- Select the LWCC pathlength based on the desired absorption range.
 Here are some reference values for different water types:
 - Fresh water with absorption range > 4.0 m⁻¹
 - Coastal-ocean waters with absorption range 1.0–4.0 m⁻¹
 - Open ocean water with absorption range < 1.0 m⁻¹
- · Select the usable internal volume.

Next, you can select components to complete your CDOM analysis system, depending on the selected LWCC. The effective pathlength of

WPI's LWCC is defined as the equivalent pathlength of the cell, if it is assumed that the LWCC strictly follows the Beer-Lambert law.

Typically, the longer LWCC pathlength is used to increase the sensitivity when the maximum absorbance values are supposed to be <0.1 AU (Absorbance Unit). Inversely, when absorbance measurements are above 1.4 AU, the LWCC pathlength should be decreased to ensure that measurements still remain within the linear range of the LWCC detection system.

LWCC Type	Pathlength (cm)	Noise (mAU)	Absorbance Range (mAU)	Absorption Range (m ⁻¹)	Internal Volume (mL)
LWCC-4010	10	< 0.1	0.5 - 1400	0.012 - 32.0	0.31
LWCC-4050	50	< 0.2	1.0 - 1400	0.005 - 6.4	1.57
LWCC-4100	100	< 0.5	2.5 - 1400	0.006 - 3.2	3.1
LWCC-3050	50	< 0.1	0.5 - 1400	0.002 - 6.4	0.125
LWCC-3100	100	< 0.2	1.0 - 1400	0.002 - 3.2	0.250
LWCC-3250	250	< 0.5	2.5 - 1400	0.002 - 1.2	0.625


The useful absorption range calculation is based on the absorbance detection limits of the LWCC, considering a wavelength range of 300—700 nm.


Absorbance measurements obtained with WPI's LWCC and Tidas S300 UV/VIS spectrophotometer are linear up to 1.4 AU. The measured absorbance can be converted to the spectral absorption coefficient $\alpha(\lambda)$, commonly used in oceanography for CDOM measurements. Absorbance and spectral absorption are related by the formula: $\alpha(\lambda)=2.303$ A (λ) / L, where 2.303 is the conversion factor from decimal to natural logarithmic, A (λ) is the absorbance at wavelength λ and L the LWCC pathlength.


Common applications of CDOM detection

- Determine the biogeochemical cycles, e.g., the organic carbon-based cycle in the ocean.
- Monitor and map surface-water masses.
- Measure the UV light penetration into the ocean to determine:
 - Photosynthesis reaction with effects on phytoplankton population
 - Effect on oceanic food chains
 - Atmospheric oxygen concentration
- Monitor the light absorption of CDOM as it relates to heat storage and the decline of sea ice.

Two typical absorption spectra recorded with an **UltraPath** (**UPUV**) of a seawater and a fresh water sample collected in November 2007 are shown in the figures. Due to their high absorbance, both samples were analyzed in the 10 cm pathlength. The CDOM sample labeled Mayagüez Bay is from oligotrophic, low productive waters with high salinity collected off the west coast of Puerto Rico in the Mayagüez Bay and was measured with a 200 cm pathlength.







CDOM Sample "Mayagüez Bay" was collected from the high salinity oligotrophic waters of Mayagüez Bay on the west coast of Puerto Rico (2001). Data courtesy of NASA Stennis Space Center.

System Configurations

Refer to the table on the previous page for selecting the LWCC pathlength.

The CDOM-FRESH System > 4 m ⁻¹						
LOW VOLUME HIGH VOLUME						
Product Description Item		Product Description	ltem #			
CDOM-FRESH-LV System Includes:	CDOM-FRESH-LV	CDOM-FRESH-HV System Includes:	CDOM-FRESH-HV			
LWCC, 50 cm pathlength	LWCC-3050	LWCC, 10 cm pathlength	LWCC-4010			
Photo Diode Array (PDA) Spectrophotometer System, UV/VIS (190-720 nm), with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System, UV/VIS (190-720 nm), with integrated deuterium/halogen light source	505067			
(2) UV-Enhanced Fiber Optic Cables, 1 m, 600 µm Core	WVLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables, 1 m, 600 μm Core.	WVLUXUVIS-S-600			
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS			
LWCC Injection System	89372	Injector Kit	72100			
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120			
UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200			
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395			

	The CDOM-CO	DAST System 1-4 m ⁻¹			
LOW VOLUME	HIGH VOLUME				
Product Description	Item #	Product Description	Item #		
CDOM-COAST-LV System Includes:	CDOM-COAST-LV	CDOM-COAST-HV System Includes:	CDOM-COAST-HV		
LWCC, 100 cm pathlength	LWCC-3100	LWCC, 50 cm pathlength	LWCC-4050		
PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System, UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067		
(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 μm Core.	WVLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 μm Core.	WVLUXUVIS-S-600		
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS		
LWCC Injection System	89372	Injector Kit	72100		
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120		
UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200		
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395		

The CDOM-Ocean System < 1 m ⁻¹							
LOW VOLUME HIGH VOLUME							
Product Description	Item #	Product Description	Item #				
CDOM-OCEAN-LV System Includes:	CDOM-OCEAN-LV	CDOM-OCEAN-HV System Includes:	CDOM-OCEAN-HV				
LWCC, 250 cm pathlength	LWCC-3250	LWCC, 100 cm pathlength	LWCC-4100				
PDA Spectrophotometer System. UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067	PDA Spectrophotometer System. UV/VIS (190-720 nm) with integrated deuterium/halogen light source	505067				
(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 μm Core	WVLUXUVIS-S-600	(2) UV-Enhanced Fiber Optic Cables. 1 m, 600 μm Core.	WVLUXUVIS-S-600				
Ministar Peristaltic Pump	504011	PeriStar Pro Pump	PeriPro-4LS				
LWCC Injection System	89372	Injector Kit	72100				
TTL Control Module for Ministar and/or Peristar	503120	TTL Control Module for Ministar and/or Peristar	503120				
UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200	UV-Enhanced Fiber Optic Cable	WVLUXUVIS-S-200				
(2) SMA Bulkhead feed through Connector/Coupler	13395	(2) SMA Bulkhead feed through Connector/Coupler	13395				

Filter Holder for Glass Fiber Filters

Designed for field use

Features

- Simple measurements for particulate absorption
- · Rugged and portable
- · Performs as well as a laboratory based spectrophotometer

Benefits

- · QFT corrects for the effect of scattering
- · Averages out larger, non-organic particles

Applications

Measurement of UV/VIS absorption of particulate matter in seawater

WPI's filter holder for particulate absorption measurements is specially designed for field use. It is rugged and portable. It performs as well as a laboratory based spectrophotometer. It can be directly connected to WPI's line of fiber optic spectrometers and light sources. Instead of collecting your samples, transporting them to a laboratory, and accepting the loss of spectral information associated with it (Sosik, 1999), particulate absorption can now be measured on site.



Particulate absorption of fresh water and seawater can be determined by filtering a known amount of sample through a Glass Fiber Filter (GF/F) and measuring the particulate absorption coefficient $a_{\text{P}}(\lambda)$ concentrated on the filter. This technique is called quantitative filter technique (QFT) and corrects for the pathlength amplification, an effect of scattering. The correction of the pathlength amplification and the correction of the non-linear relationship between the optical density of samples on a Whatman GF/F filter and in suspension are discussed in Mitchell (1990).

Averages out larger, non-organic particles

A significant advantage of the filter holder is its large beam diameter of 5 mm, resulting in "averaging out" of larger non-organic particles frequently found on the filter pad when using natural samples. The removable filter fixture allows simple filter alternation and cleaning.

System Requirements

The optical throughput of QFT1 equipped with a classical GF/F filter is very low and requires a matched light source/spectrometer system. TIDAS E Base UV/VIS (504718) in combination with a D4H UV/VIS light source







89575

or a **FO-6000** VIS light source, as it is offered in the Ultrapath system, are ideally suited, as well as TIDAS S300 UV/VIS (**505067**). Light should be coupled into the **QFT1** with a 1000 µm fiber (**WVLUXUVIS-S-1000**) and from the **QFT1** to a spectrometer with a 600 µm fiber

(WVLUXUVIS-S-600).

SPECIFICATIONS

GF/F FILTER DIAMETER	25 mm
WAVELENGTH RANGE	280-730 nm *
FIBER OPTIC CONNECTION	Input: 1000 μm Output: 600 μm
MATERIAL IN CONTACT WITH FILTER PAD	Delrin
WEIGHT	0.5 kg (1 lb)

^{*} Using a Tidas E Base spectrometer and D4H UV/VIS light source.

References

M. Belz, K. Larsen, K.-F. Klein, "Fiber optic sample cells for polychromatic detection of dissolved and particulate matter in natural waters", *Proc. SPIE*, Vol. 6377, Oct 2006, 63770X

ORDERING INFORMATION

89575 QFT1, Fiber Optic Holder for Glass Fiber Filters

In-Line Fiber Optic Filter Holder

Insert optical filters into a fiber optic pathway

Features

- In-line filter for SMA terminated fibers
- Collimators for UV/VIS/NIR
- Filter diameter 8–25.4 mm
- Filter thickness 2-10 mm

Benefits

- · No stray light
- Solid design

Applications

- · Stray light filtering
- Removal of excitation light in fluorescence detection



This In-Line Fiber Optic Filter Holder allows the insertion of optical filters within a fiber optic pathway. The connectors of the filter holder assembly are compatible with WPI's range of fiber optic jumper cables and can be coupled using SMA connectors.

Filters with outer diameters from 8 to 25.4 mm and thicknesses from 2 to 10 mm can be accomodated. The design limits lateral and axial movement of the filter when secured in the holder.

Two fiber optic collimators are internally mounted in the holder to pass collimated light through the filter and then refocus the filtered light into the aperture of the output fiber. Spectral range will be

largely limited by the bandpass of the optical fibers (from UV to near IR using WPI UV-enhanced cables).

ORDERING INFORMATION

56200 In-Line Fiber Optic Filter Holder (SMA)

Smallest Fiber Optic Dipping Probe

Perfect for UV/Vis spectroscopy



Features

- · Compact and efficient dipping probe
- Tip diameter 1.8 mm

Benefits

- · Perfect for mobile applications

Applications

- · Protein and DNA sample measurements
- · Dissolution system

Mini DipTip™ is a miniature transmission probe for microliter spectroscopic sampling. Mini DipTip's tip diameter is only 1.5 mm—the size of a 17-gauge needle. It will fit into all micro centrifuge tubes on the market. Microliter samples can be analyzed cost effectively when you combine the Mini DipTip™ with one of the following:

- TIDAS E Base with FO-6000 or D4H
- TIDAS S300 series
- · Compatible with most fiber coupled spectrometers
- Ideal for multi channel applications with LEDSpec



DIPTIP SPECIFICATIONS

TIP DIAMETER	1.8 mm
LIGHT PATHLENGTH	2, 5, 10 mm
WAVELENGTH RANGE (nm)	200-1000
SAMPLE VOLUME REQUIRED	20-50 μL
DISTANCE FROM TIP TO UPPER EDGE OF SAMPLE WINDOW	7 mm
FIBER LENGTH	1.0 m
FIBER OPTIC CONNECTION	SMA 905
LAUNCH FIBER BUNDLE (7 x 200 μ m)	680 µm*
RETURN FIBER BUNDLE (7 x 200 μm)	680 µm*

*Circular packaging of the fiber bundle results in an active area equivalent to a fiber with a core diameter of 680 µm. Using a 600 µm connection is recommended and will result in negligible light loss.

ORDERING INFORMATION DIP-UV-MINI-2 Mini DipTip™ for UV/VIS/NIR (2 mm path) DIP-UV-MINI-5 Mini DipTip™ for UV/VIS/NIR (5 mm path) DIP-UV-MINI-10 Mini DipTip™ for UV/VIS/NIR (10 mm path)

Fiber Optic Collimator

Collimate an emitted light beam or couple light into an optical fiber



Features

- Maximum coupling efficiency with fused silica fibers
- Easily adjust the distance between the lens and the optical fiber for focusing
- · Not suitable for single mode laser applications

Benefits

• UV/VIS/NIR light collimation

Applications

 Generates plan-parallel beam with a 5mm beam diameter from UV/ VIS light guided in an fiber

WPI's Fiber Optic Collimator can be used for both collimating a light beam emitted by an optical fiber or coupling light from a collimated light beam into an optical fiber. The numerical aperture of the collimator is optimized for maximum coupling efficiency into typical fused silica fibers. The collimator can, for example, be used to guide a parallel light beam through a sample cuvette or an optical filter with virtually no optical losses. In this application, one collimator collimates the light into a parallel beam 5 mm in diameter, enabling it to pass a long distance without losing the energy. After the light passes the sample media, a second collimator can be used to collect the beam into the receiving fiber. A unique design feature of this collimator is that the distance between the lens and the optical fiber can be easily adjusted. This permits it to be used as a focusing device or for fine-tuning the color balance when coupling light from a light source into multimode fibers.

COLLIMATOR SPECIFICATIONS

LENS DIAMETER 5 mm

LENS FOCAL DISTANCE 10 mm

LENS MATERIAL Ultraviolet grade synthetic fused silica (KU-1)

WAVELENGTH RANGE 220 nm-2 µm

MOUNTING THREADS 3/8-24 UNF

DIVERGENCE < 0.1 rad for 1 mm core fiber

FIBER CONNECTOR INTERFACE SMA or ST

	ORDERING INFORMATION
300051	Fiber Optic Collimator (SMA)
300052	Fiber Optic Collimator (ST)
OPTIONAL	ACCESSORIES/REPLACEMENT PARTS
13395	SMA Bulkhead feed through connector/coupler, D-hole
13370	SMA half-length bulkhead coupler/connector
CC-3-UV	Cosine Corrector



Z-Dimensions Are Not Created Equal

Cuvettes come in a variety of shapes and sizes, but one of the most important specifications of a cuvette is its Z-dimension. The Z-dimension of an instrument (cuvette holder or spectrometer) is the distance from the bottom of the cuvette chamber floor to the center

of its light beam (see image). A cuvette's Z-dimension must match the Z-dimension of the instrument with which it will be used.



Each manufacturer designs its instruments

with a specific Z-dimension. Common Z-dimensions include 8.5 and 15 mm, and sometimes 20 mm. When purchasing small volume cuvettes, the correct Z-dimension becomes critical. Matching the Z-dimension of the cuvette to the Z-dimension of the instrument ensures that the light beam passes through the center of small samples. The table below shows the standard Z-dimension of the spectrometer sample compartments for many manufacturers.

Manufacturer	Z-Dimension
Agilent®	15 mm
Avantes®	15 mm
Beckman®	8.5 mm
Bio-Rad®	8.5 mm
Cecil®	15 mm
Eppendorf®	8.5 mm
Hewlett – Packard®	15 mm
Hitachi®	8.5 mm
Jasco [®]	11 mm
J & M®	8.5 mm
Ocean Optics®	15 mm
Perkin – Elmer®	15 mm
Pharmacia®	15 mm
Shimadzu®	15 mm
Spectronics®	8.5 mm
Stellarnet®	15 mm
Turner®	8.5 mm
Varian®	20 mm
WPI	15 mm

To determine the Z-dimension of a cuvette holder:

- Use strips of heavy paper that will fit neatly into a cuvette (for example, 12 mm x 50 mm) and not allow light to pass through the cuvette.
- Poke a tiny hole in each paper "sample." For example, one paper sample could have a hole at 8.5 mm, one at 15 mm, one at 20 mm.
- One at a time, insert the paper samples into the cuvette and place the cuvette into the cuvette holder. The paper sample with the pin hole at the instrument's Z-dimension will allow light to pass. The other paper samples will not allow light to pass.

If you have an instrument that is not on the list and need to know its Z-dimension, please contact WPI at (941) 371-1003 or wpi@wpiinc.com.

Optical Glass and Quartz Cuvettes

For spectrophotometry and fluorometry

Features

- Absorbance and fluorescence cuvettes
- 2-100 mm path length
- Fused silica and glass cuvettes
- · Flow cuvettes
- Microliter cuvettes

Benefits


- · Eight different styles
- Inexpensive single use glass cuvettes

Applications


- Absorbance spectrophotometers
- Fluorescence fluorometers
- · Chemistry, Biochemistry
- · FIA, GFIA, Quality control

WPI's glass and synthetic quartz cuvettes are ideal for UV/VIS/NIR absorbance or fluorescence experiments.

Synthetic quartz can be used in deep UV applications and is recommended for fluorescent applications, as it does not exhibit background fluorescence. Quartz cuvettes (absorbance, fluorescence and flow) are shipped individually packaged, glass cuvettes are shipped in packages of 10 cuvettes. These economic quartz and glass cuvettes are ideal for precision measurements because of their high quality materials used and their low manufacturing tolerances. Typical transmission curves of glass and synthetic quartz cuvettes are shown below.



Transmission curves of Glass and Synthetic Quartz Cuvettes. The cuvettes were empty, thickness 1.25 mm x 2, including surface reflections, measured with a TIDAS II against air as reference.



A complete transmission spectrum of glass and synthetic quartz cuvettes from 190 nm to 4 mm is shown. Cuvettes were empty, thickness 1.25 mm x 2, including surface reflections, measured with a TIDAS II.

5	Р	E	CI	FΙ	C	Αī	Ц	O	N	S	

Cuvette Material	Spectral Range (>80%)	Transmission Difference Between Different Cuvettes
OPTICAL GLASS	350 - 2500 nm	Less than 1%
SYNTHETIC QUARTZ	200 – 2500 nm	Less than 1%

Style G



ORDERING INFORMATION								
WPI PN	Style	material	Polished windows	path [mm]	Dimensions [mm]	volume [mL]	Beam width [mm]	
			STAND	ARD RECTAI	NGULAR CUV	ETTES		
CUV2101-1*	В	Quartz	2	1	3.5x12.5x45	0.35	10	
CUV2102-1*	В	Quartz	2	2	4.5x12.5x45	0.7	10	
CUV2011-1*	В	Quartz	2	5	7.5x12.5x45	1.7	10	
CUV1022-10	C	Optical Glass	2	10	12.5x12.5x45	3.5	10	pack of 10
CUV2012-1	C	Quartz	2	10	12.5x12.5x45	3.5	10	
CUV2105-1	C	Quartz	2	20	22.5x12.5x45	7	10	
CUV2106-1	C	Quartz	2	30	32.5x12.5x45	10.5	10	
CUV2107-1	C	Quartz	2	40	42.5x12.5x45	14	10	
CUV2108-1	С	Quartz	2	50	52.5x12.5x45	17.5	10	
4								

^{*89341} Cuvette spacer for 1-mm cuvettes (part **CUV2101-1**)

^{*89337} Cuvette spacer for 5-mm cuvettes (part CUV2011-1, CUV2023-1, CUV2063-1)

05007		er for 5 fillin cav	ettes (part	.012011 1, 0012023	1, 2012005 17				
	SELF MASKING SEMI MICRO CELL CUVETTE								
CUV2023-1	l * D	Quartz	2	5	7.5x12.5x45	0.7	4		
CUV2031-1	I D	Quartz	2	10	12.5x12.5x45	1.4	4		
CUV2025-1	I D	Quartz	2	20	22.5x12.5x45	2.8	4		
CUV2032-1	I D	Quartz	2	10	12.5x12.5x45	1	3		
CUV2033-1	D	Quartz	2	10	12.5x12.5x45	0.7	2		
CUV2034-1	l D	Quartz	2	10	12.5x12.5x45	0.35	1		
		SE	LF MASH	KING CONTINU	OUS FLOWTHR	ROUGH	CELL		
CUV2063-1	* E	Quartz	2	5	7.4x12.5x45	0.035	Ø3		
CUV2061-1	I E	Quartz	2	10	12.5x12.5x45	0.07	Ø 3		
CUV2065-1	I E	Quartz	2	20	22.6x12.5x45	0.14	Ø3		
CUV2066-1	I E	Quartz	2	30	32.6x12.4x45	0.21	Ø 3		
CUV2062-1	F	Quartz	2	10	12.5x12.5x45	0.48	4x12		
SEL	F MASKII	NG CONTIN	UOUS FI	OW THROUGH	CELL, SMALL	INPUT,	LARGE OUTPUT Z=8.5 MM		
CUV2614-1	I Н	Quartz	2	10	12.4x12.4x35.6	0.03	Ø 2		
			N	IICRO CELL WI	TH BLACK WAL	LS.			
CUV2674-1	l J	Quartz	2	10	12.5x12.5x45	0.05	2		
				FLUOR	ESCENCE				
CUV2051-1	I A	Quartz	4	10	12.5x12.5x45	3.5	10		
CUV2052-1	I A	Quartz	4	10	12.5x12.5x45	1.4	4		
				LONG PAT	TH CUVETTE				
CUV2071-1	l G	Quartz	2	100	102.5 x 22 Ø	28	19		



^{*89342} Cuvette spacer for 2-mm cuvettes (part CUV2102-1)

Photodiode Array Based Spectrophotometers

High performance fiber optic spectrometer systems





Features

Plug and Play – with help of standard accessories the TIDAS® E sits ready to be immediately operational

Increased Flexibility – the patented cuvette holder allows a variety of measurement configurations that are seamlessly integrated within one instrument

Diode Array Technology – enables quick and precise measuring with a low level of scattered light

Software – modular and user friendly, permits the data acquisition and evaluation of individual tasks

- · Fiber-optic diode-array spectrometer
- Measures in transmission and absorption with external light source
- Supports external immersion probes, temperature-controlled cuvette holders and flow cells
- Supports external reflection measurement heads
- Perfoms film thickness measurement using white light interference
- Suitable for color measurement and analysis
- Chemometric modelling
- Certified instrument quality (wavelength and photometric accuracy)
- Software for instrument control, data acquisition and data evaluation included

Benefits

- Compact, rugged and transportable instruments
- · Multimodal and combined methods
- · High dynamic range
- Fast and reliable results
- Powerful, yet easy to use software

Applications

- · Standard chemistry and biochemistry
- · Trace level nutrient analysis
- · Detection in flow injection analysis (FIA)
- HPLC analysis with WPI's LWCC-M flow cells
- CDOM detection with UltraPath

The **TIDAS** spectrometer and spectrophotometer series is based on a high performance monolithic spectrometer module manufactured by Zeiss for UV, UV/VIS and VIS/NIR applications. The **Tidas E Base** spectrometers are ideal for modular spectrometer systems with separate detectors, light sources and sample cells. The **TIDAS S300** spectrophotometers already include a light source.


The **TIDAS** series instruments are specialized for research and instructional concepts. This Diode Array based technology allows for fast and precise measurements, and the instruments are universally applicable. The **TIDAS E Base** and the **TIDAS S300** are suitable for routine measurements in the laboratory and outperform conventional benchbased spectrophotometers and CCD-based spectrometer modules, when it comes to high precision fiber optic sampling. They rely on a monolithic optical bench made by Zeiss, which is optimized for fiber optic applications. Most cuvette-based standard spectrometers lose more than 90% of light through expensive prism decoupling. The Tidas systems are designed for fiber optic sampling cells. Using suitable light sources and sample cells, spectral detection in the wavelength range of 190 – 720 nm can be performed.

Applications


The Tidas systems are ideally suited for WPI's fiber optic sampling equipment. High sensitivity detection systems for flow analysis can be assembled using WPI's Liquid Waveguide Capillary Cells (LWCC) with effective pathlengths ranging from 50 to 500 cm. These setups are frequently used in fluid injection analysis systems for nutrient analysis (nitrite, nitrate, phosphate, iron) in oceanographic applications. Microliter sampling systems for UV/VIS applications can be assembled using WPI's DipTip™ dipping probes.

Software

TidasDAQ 3 software is included with each instrument for data collection and data analysis. TidasDAQ is used to run the spectrometer module, collect spectra in either single or continuous mode, control the digital I/ Os, save the experimental data to disk, and analyze the data. Further, TidasDAQ can export data directly into GRAMS/AI, a feature very useful for advanced data analysis for pharmaceutical applications and requirements.



TIDASDAQ acquisition window, showing an absorbance baseline.



Spectra may be displayed and analyzed in 2D and 3D format. This allows the user to conveniently interpret "time acquisition" data typically done with a TIDAS-E-BASE-LWCC flow system.

TidasDAQ: Data Collection & Instrument Control


With TidasDAQ, high precision intensity, absorbance, transmittance or normalized spectra can be obtained in less than a second. Only a few parameters need to be adjusted to obtain spectral data. Sampling of single scans, continuous full spectra scans or triggered scans is possible. Chromatograms can be displayed and logged to disk at up to four wavelengths. Data Export of 2D and 3D Spectrograms, as well as Chromatograms is supported in ASCII, Spectralys/SpectraView, Excel and Grams/AI formats. Light sources and other sampling instrumentation can be controlled via the TTL level digital outputs, and data collection can be triggered by TTL leveled external inputs of the Tidas E Base.

Spectra can be recorded in 2D and 3D view. Mathematical computation, Derivation, Smoothing, Quantification and other functions are available to work with your data. The Quantification module allows single point and multiple point analysis, multiple linear regression, partial least square and principle component analysis. Data can be exported out of a 3D analysis file into separate scans. Furthermore, chromatograms as well as spectrograms can be copied directly into Excel for further data analysis.

Typical Usage


These high performance spectrometer modules and spectrophotometers combine modularity, flexibility, ease of use and high quality in one unit. They are ideally suited for WPI's fiber optic sampling equipment. High sensitivity detection system for flow injection analysis (FIA) and gas segmented flow injection analysis (GFIA) systems can be assembled with WPI's proprietary Liquid Waveguide Capillary (LWCC) flow cells with effective pathlengths ranging from 50 to 500 cm. Typical applications include colorimetric trace detection of nitrite, nitrate, phosphate and iron in oceanographic research.

Small flow volume samples, like those that may be found in Optofluidics and High Performance Liquid Chromatography (HPLC), can be measured conveniently and accurately with WPl's **LWCC-M** series of flowcells. LWCC-M flow cells feature pathlengths of 10, 50 and 100 mm with corresponding sample cell volumes of 2, 4, 12 and 24 μ L. Finally, small discrete samples of 2-40 μ L volume often found in biochemistry (like, Protein or DNA) may be analyzed using WPl's **DIP-UV-MINI** miniature dipping probes with 2, 5 and 10 mm pathlength. Other applications for



Mean measured absorption spectra for marine samples, measured with TIDAS1 and a 100-cm pathlength LWCC (adapted from Beine et al., J Geophys Res, Vol. 117, D00R15).

the TIDAS series spectrometers and spectrophotometers include material science (forensics, semiconductor technology), quality control (automotive and food), as well as new research in the field of nanotechnology and rapid kinetics.



Absorbance spectrum displayed using TIDASDAQ software

	SPECIFIC	CATIONS				
	TIDAS S300/E BASE UV	TIDAS S300/E BASE UV/VIS	TIDAS S300/E BASE VIS/NIR			
WAVELENGTH RANGE	190-390 nm	190-720 nm	300-1100 nm			
SPECTRAL RESOLUTION	<3 nm	<7 nm	103 nm			
WAVELENGTH ACCURACY	±1 nm	±1 nm	<3 nm			
PHOTOMETRIC ACCURACY UV	±10 mAU	±10 mAU	-			
PHOTOMETRIC ACCURACY VIS	-	±10 mAU	-			
WAVELENGTH REPRODUCIBILITY	<0.1 nm	<0.1 nm	<0.1 nm			
BASELINE DRIFT @250NM	5x10 ⁻⁴ AU/h according to ASTM E685	5x10 ⁻⁴ AU/h according to ASTM E685	5x10 ⁻⁴ AU/h according to ASTM E685			
SIGNAL-TO-NOISE RATIO	<3x10 ⁻⁵ AU according to ASTM E685	<3x10 ⁻⁵ AU according to ASTM E685	<3x10 ⁻⁵ AU according to ASTM E685			
NUMBER OF DIODES	256	256	256			
LIGHT SOURCE Not included for Tidas E Base	Deuterium	Deuterium & Halogen	Halogen			
OPTICAL BENCH Monolithic spectrometer module with concave aberration corrected holographic grating;						
DETECTOR ARRAY	Hamamat	su photodiode arra	ay, 256 pixel			
A/D RESOLUTION		16 Bit				
FIBER OPTIC CONNECTION SMA905, 600 μm						
SOFTWARE		Tidas DAQ3 (include	ed)			
SCRIPT LANGUAGE METHOD DEVELOP		Yes				
SYSTEM REQUIREM	ENT	Windows 7, 8, 10)			

	ORDERING INFORMATION
504717	TIDAS E Base, UV 190-390 nm
504718	TIDAS E Base, UV/VIS 190-720 nm
504719	TIDAS E Base, VIS/NIR 300-1100 nm
505066	TIDAS S300, UV 190-390 nm with deuterium lamp
505067	TIDAS S300, UV/VIS 190-720 nm with deuterium/halogen lamp
505068	TIDAS S300, VIS/NIR 300-1100 nm with halogen lamp
505069	Tidas S300 VIS/NIR 300-1100 nm (without halogen lamp)
Systems	includes power supply, TIDAS DAO software, RI-45 cable and manual.

Deuterium Halogen Fiber Light Source

For a continuous spectrum in the UV, VIS and NIR range

Features

- Continuous spectrum from 200 nm-1700 nm
- Integrated shutter with switch and TTL control
- SMA fiber optic connection
- Separate UV and VIS bulb control

Benefits

- Matched deuterium and halogen bulbs
- Matched optical output for WPI's flow cells and fiber optic probes
- Low drift < 1 mAU/h @ 254 nm
- Simple exchange of light bulbs

Applications

- UV/VIS/NIR absorbance spectroscopy
- · Excitation light source for fluorescence applications
- Colored Dissolved Organic Matter (CDOM) detection in seawater
- Nutrient analysis in fresh water, drinking water and seawater using WPI's LWCCs
- Protein detection using WPI's DIP UV MINI fiber optic probes

The **D4H** is a combined deuterium and halogen light source for UV/ VIS and NIR applications. This light source is ideally suited to work with WPI's spectrometer modules and sample cells. It supplies a continuous spectrum in the UV, VIS and NIR range from 200 nm to 1100 nm. The **D4H** is equipped with an integrated electrical shutter, which can be controlled by a switch or a TTL signal.



D4F



Replacement Deuterium Lamp **503847**



Replacement Halogen Lamp **503848**

LIGHT SOURCE SPECIFICATIONS					
	D4H	FO-6000			
APPLICATION	UV/VIS/NIR	VIS/NIR			
SPECTRAL RANGE	200—1100 nm	380—1700 nm			
DEUTERIUM LAMP LIFE	2000 hr	NA			
TUNGSTEN/HALOGEN LAMP LIFE	2000 hr	3000* hr			
STABILITY	1-2 mAU/h	<0.5 mAU/h			
POWER CONSUMPTION	140 W	6 W			
POWER REQUIREMENTS	110/240V, 50-60 Hz, 1A	12VDC/1A			
SHUTTER/TTL TRIGGER	Yes	Yes			
MAX. FIBER OUTPUT	1000 μm	1000 μm			
CONNECTIONS	SMA	SMA			
SHIPPING WEIGHT	13.2 lb (6 kg)	1.3 lb (0.6 kg)			
DIMENSIONS (W/H/L)	7 x 6.2 x 9.8 in. (17.8 x 15.7 x 25 cm)	4.8 x 2.8 x 7.5 in. (12 x 7 x 19 cm)			

*Lamp life is dependent upon internal power settings.

		D4	IH Spec	trum		
	60000					
	50000 -	_	Deuterium Halogen	1	^	
Intensity [counts]	40000 -	_		ı + Halogen		٨
0)/	30000 -	"\ \\				
nsity	20000 -	ľ		W/		
Inte	10000 -			w	\mathcal{N}	
		200	400	F00	600	700
	200	300	400	500	600	700
		W	<i>l</i> avelength	[nm]		

ORDERING INFORMATION				
D4H	Deuterium Halogen Light Source (200 nm-1100 nm)			
503848	Halogen Replacement Lamp for D4H			
503847	Deuterium Replacement Lamp for D4H (> 215 nm)			
OPTIONAL ACCESSORIES/REPLACEMENT PARTS				
TIDAS-D2	Replacement Deuterium Lamp, for Tidas II			
TIDACII	Deple coment Helegen Lemp (Type 1) for Tides II			

TIDAS-D2 Replacement Deuterium Lamp, for Tidas II
TIDAS-H Replacement Halogen Lamp (Type 1), for Tidas II
TIDAS-H2 Replacement Halogen Lamp (Type 2), for Tidas II
D2H-DB Replacement Deuterium Lamp, for D2H
D2H-HB Replacement Halogen Lamp, for D2H
D2H-HBER Replacement Deuterium Lamp, Extended Range, for D2H

Tungsten Fiber Light Source

High color temperature



Features

- · Visible light source
- 380 nm 1700 nm
- Low drift < 0.5 mAU/h
- SMA fiber connection
- Electrical shutter control (switch & TTL)


Benefits

- Bulb with 10,000-hour lifetime
- Fits all of WPI's fiber optic probes and flow cells
- Design matched to **LWCC** applications
- Temperature controlled optical bench

Applications

- Low noise VIS absorbance measurements
- Trace analysis of nutrients in seawater and freshwater
- · CDOM analysis in seawater
- · Analytical chemistry, environmental science and life science

The **FO-6000** is a continuous fiber optic light source featuring an extended visible part of the light source (380 nm - 1700 nm). It has a SMA fiber optic connector. The shutter and lamp can both be controlled via a switch or external TTL triggering. This light source offers a wide assortment of applications. A special feature of the **FO-6000** is its color balancing optics, which shifts the usable range of the light source from traditional 420 nm down to 380 nm wavelength. Due to its thermally controlled optical bench, it is particularly suitable for low noise and low drift applications.



FO-6000 Fiber Optic Light Source 800120 Replacement Lamp for FO-6000

Absorbance Detection

Detection of organic compounds in water analysis

Absorption of light correlates to the energy of a photon that is taken-up by electrons of the substance atom. The electromagnetic energy is transformed into internal energy of the absorbent substance. The absorbance of a substance quantifies how much of the incident light is absorbed by it (instead of being reflected or refracted). Precise measurements of the absorbance at many wavelengths allow the identification of a substance via absorption spectroscopy, where a sample is illuminated from one side, and the intensity of the light that exits from the sample in every direction is measured (see Fig. 1). A few examples of absorption are ultraviolet–visible (UV-Vis) spectroscopy or infrared (IR) spectroscopy.




Fig 1. Concept of absorbance spectroscopy using white light and optical components to filter out light of a specific wavelength that interacts with molecules in the solution. Absorbance at this specific wavelength by the molecules in the solutions is detected as a decrease in light intensity (Spectrophotometer-Source: http://chemwiki.ucdavis.edu/).

Absorption is the amount of light that a substance takes in and does not allow to pass through it. Spectrophotometers actually measure transmission, the amount of light that passes through a sample, but this is converted into absorption by comparing the bulb output to the light that has passed through the sample. Light sources that can be used for absorbance spectroscopy depend strongly on the used substance to label a specific molecule and can span the entire electromagnetic spectrum of light.

Protein detection uses the UV- spectrum (typically, 260 nm and 280 nm, while further information is obtained at 230 nm and 320 nm, but compensated by a selection of ratios and background corrections) and is most commonly used to estimate DNA or RNA concentration and to analyze the purity of the preparation. Further application spans the measurement of the light scatter at 600 nm to monitor the growth rate of a cultured bacterial population and to identify the peak concentration.

Industrial applications that also cover the UV spectrum for protein detection is linked to food analysis (for example, characterizing the grading of olive oil as extra virgin, virgin oil, etc., as set out by European regulations) or quality control in the pharmaceutical industry. In addition, industrial applications of absorbance spectrophotometry cover the characterization of water purity or waste water analysis, in addition to the determination of specific organic molecules like nitrate, nitrite or phosphor. This last application is interesting when using the visible light (Vis) spectrum.

Biofluorometer

Reliable, 2-channel, LED-based fluorometer

Features

- Two photomultiplier inputs
- · Light excitation with high power LEDs
- Modes: single excitation & single emission, single excitation & dual emission, dual excitation & single emission
- 2 channel mode (two single excitation & single emission)
- Optical connections: Liquid Light Guides (LLGs) and SMA terminated fibers
- Sampling rates up to 1 kHz (1000 ratios/second)
- Automatic LED light drift correction for long term measurements
- Automatic room light correction
- Optional fiber optic probes for horizontal tissue bath applications
- Optional imaging probes for Langendorff systems
- Optional attachments for direct connection to fluorescence microscopes via epifluorescence port (excitation) and C-Mount (Emission) via liquid light guides

Benefits

- Versatile instrument for horizontal tissue bath, Langendorff and microscope applications
- Warmup time < 1 minute
- Low bleaching mode (5% LED On time)
- Can be combined with imaging based Sarcomere Detection System (OptiSarc)
- SMA and Liquid Light Guide (LLG) may be exchanged (optional)
- Customized analysis techniques in WPI's MDAC data acquisition software


Applications

• Auto-fluorescence imaging

The **SI-BF-100** is an LED-based fluorometer for life science applications. It is ideally suited for ratiometric calcium detection (FURA-8) and ATPase detection (via NADH fluorescence). With up to three LED modules (wavelengths), the **SI-BF-100** covers many fluorometric applications in neuroscience and cell biology.

The **SI-BF-100** enables the detection and analysis of fluorescence signals in four different modes:

- Single excitation/single emission—In this classical mode, a fluorophore
 is excited at one wavelength and the fluorescence signal is detected
 at a single higher wavelength using one photomultiplier. The
 concentration of the analyte is directly proportional to the intensity of
 the detected signal.
- Dual excitation/single emission—A fluorophore is excited at two wavelengths and the fluorescence signal is detected at one wavelength using one photomultiplier. The concentration of the analyte is proportional to the ratio of the two detected fluorescence signals. This ratiometric concept minimizes the effect of indicator dye bleaching and motion artifact in experiments. A typical example is the detection of free calcium in muscle tissue using the indicator dye
- Single excitation/dual emission-A fluorophore is excited at one wavelength and the fluorescence signal is detected at two wavelength using two photomultipliers.
- Dual excitation/dual emission–Two separate fluorophores are excited at different wavelengths and the fluorescence signal of each fluorophore is detected at two separate wavelengths using two photomultipliers.



The Biofluorometer can be connected to the epifluorescence port of an inverted microscope and its high intensity LED light source is used for excitation of the corresponding dye.

Fluorescence Probe

The principle of a fluorescence probe is to bring light to a sample and to detect reflected or fluorescent light from the sample surface. These probes may be used for the detection of the transient response of free ion concentrations, like calcium, potassium, sodium, magnesium. Fluorescence probes can also be designed to detect pH and membrane potential. Auto-fluorescence like the detection of ATPase activity via NADH or FAD is another application. If a fluorescent dye is used to label a molecule, a custom fluorescence probe can be designed to detect the emission.

To select the appropriate fluorescent probe, consider the detection principle and the dye to be used. All probes use fibers with 300 μm core diameter. Excitation fibers have 1000 μm SMA connectors for excitation and double emission probes, and 1500 μm SMA connectors for excitation and single emission probes. The standard probes come with a rectangular head (3 mm \times 0.7 mm) or round head (2 mm diameter). The fibers that are used for such probes are highly flexible plastic fibers or fused silica fibers. Custom fibers for your application may be designed with different core fiber diameters.



(Left) Customized single emission round probe with a 5 mm diameter for one excitation and one emission wavelength.


(Right) Customized rectangular probe.

Application	Probe	Excitation	Emission
Free Calcium via Indo-1	Double Emission	365nm	405nm/486nm
Free Calcium via Fura-8	Single Emission	365nm/410nm	525nm
ATPase via NADH/TAMRA	Double Emission	365nm	472/572nm
FAD	Single Emission	470nm	525nm


Auto-fluorescence imaging

Genetically modified frogs

Fluorescence signaling with single emission probe on living frogs of wildtype tissue (wt) and genetically modified tissue from frog expressing a protein of "auto-fluorescence" called EGFP. Data below showing an increase in autofluorescence signal for EGFP frogs.



- **A** Wildtype (wt, left) and EGFP (right) frog strain, showing the difference in color between both strains. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.
- **B** Experimental setup for frog dorsum in vivo fluorescence measurement with single emission probe. Observe the LED beam/spot as measurement area of frog dorsum. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.




Fluorescence response of frog dorsum in vivo with single emission probe from whole frog dorsum measurements for EGFP (red) vs. wildtype (wt, green). Note the distinguishable increase in auto-fluorescence signal for EGFP frog dorsum but confounded with tissue color backreflection

Genetically modified mice isolated hearts

Fluorescence signaling with LLG Imaging probe of wildtype tissue (wt) and genetically modified tissue from mice isolated heart expressing a protein of "auto-fluorescence" called EGFP. Data below showing an increase in autofluorescence signal for EGFP heart at two LED current intensity.



Experimental setup for isolated mouse heart fluorescence measurement with LLG Imaging Probe.



Fluorescence signal from isolated mouse heart with LLG Imaging Probe for EGFP (red) vs. wildtype (wt, green). Note that auto-fluorescence signals can be distinguished between EGFP vs. wildtype, although the signals and difference are weak. Note that back reflection can be considered as constant, as heart and dish support colors are the same. Photo courtesy of Dr. Xu-Feng, Qi Jinan University.

Dye loading concept for measuring fluorescence signal on a leaf


Often calcium fluorescence is measured in leaves via specific dye loading techniques, like microinjection. A simplified concept for dye loading can also be the use of the root fluid transport system. In this example Fluoresceine dye was used to exemplify this technique. Measurements

of Fluoresceine in the leaf was done with the BF-100 and an SMA single emission probe for direct sensing under a dark condition. Backreflection is considered as constant as booth leaves had the same coloring.

For the setup, the plant with root was incubated in Fluoresceine solution (using physiological transportation) for ~1hr or 2~hrs in Fluoresceine concentration of 16 mg/l and 100 mg/l. Measurements were taken on the side without vein. Clear distinguishable effects of the different Fluoresceine concentration can be observed, notably regarding the control measurement with Milli-Q water.

Example of plant dye loading via the plant root for two different

Fluoresceine concentration. Photo courtesy of Professor Yao Youli of Yangzhou University, China.



Plant dye loading of fluoresceine at two different concentrations (green, red) via the plant root. For comparison see only Milli-Q water (blue). Note that dye loading time was 1 hr (green) and 2 hrs (red). Fluorescence measurement was conducted by shining directly on the leaf under a dark condition.

SI-BF-100 SPECIFICATIONS

MEASUREMENT PRINCIPLE Fiber optic fluorometer with 2 inputs and 1 output DETECTOR INPUTS 2 PMTs

EXCITATION High Power LED Modules: 365 nm, 420 nm, 470 nm,

530 nm (select any 3 modules, when ordering)
ANALOG OUTPUT RANGE 0–10V

OPTICAL CONNECTIONS Choice of Liquid Light Guide (LLG) or SMA connections
POWER 12 V/2 A (includes external 100 – 240 V / 50 – 60 Hz

power supply)

DIMENSIONS (h x w x d) 3.5 x 17 x 13 in. (88 x 431 x 330 mm)

	ORDERING INFORMATION					
SI-BF-100L	SI-BF-100LLG Biofluorometer with LLG Optical Connections					
SI-BF-100S	SI-BF-100SMA Biofluorometer with SMA Fiber Optic Connections					
99261	C-Mount Microscope Attachment for 2x PMTs					
99259	C-Mount microscope attachment for 1x camera & 2x PMTs (includes 1x camera C-Mount adapter with adjustable aperture,					

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

SI-BF-SMA-U	SI-BF-SMA-UPGRADEBiofluorometer Upgrade Kit for SMA optical probes				
SI-BF-LLG-U	SI-BF-LLG-UPGRADE Biofluorometer Upgrade Kit for LLG connections				
802407	Liquid Light Guide (LLG), 3mm diameter, 6' long				
M3301	Manual Manipulator for securing the probe				
M10	Magnetic Base				
94650	Single Emission, Small Tissue Probe				
94689	Dual Emission, Small Tissue Probe				

Check our website for new LED modules, emission filters and dichroic mirrors for specific applications.

Precision Stereo Zoom Microscope

Modular, two parallel beam path design with high quality optics

Features

- Planachromatic objectives, no optical distortion
- Large zoom ratio: 8:1
- Large zoom range: 0.62x-5.0x
- Step and continuous zoom

Benefits

- True trinocular version available
- High-contrast imaging, ideal for observing transparent, low-contrast objects
- · Long working distance option available

Applications

• Integrated optical and biological research

The fourth generation of WPI's precision stereo zoom microscopes uses modular, two parallel beam path design and high quality optical system. The advanced optical design with planachromatic objectives provides a sharp and distortion-free contrast image throughout the entire zoom range and comes with an impressive 5-year warranty.

The **PZMIV** is available in a binocular or trinocular version. In addition, an extensive list of optional accessories is available that makes the PZMIV suitable for integrated optical and biological research.

The microscope comes with a track stand, standard 10x eyepieces (wide-field, distortion-free and high eye point) and a 1x planachromatic distortion-free objective.


All **PZMIV** and **PZMTIV** microscopes come with 10x eyepieces, built-in 1x auxiliary lens and light ring adapter.

Long working distance option

A common application for the PZMIV and the **PZMIII** is to add a 0.5X objective to allow for more working distance under the lens. This increases the PZMIV working distance from 80 mm to 189 mm. The PZMIII working distance increases from 100 mm to 177mm.

The use of this auxiliary lens drops the magnification range to half, so we also recommend acquiring the 20X eyepieces to return the microscope to its standard magnifications.

See the specification tables marked areas for more details. On trinocular versions of the microscope with the LWD option, the camera view is nearly the same as the eyepiece visual field of view.



The trinocular version is a true trinocular, with continuous operation of both eyepieces and photo tube simultaneously. There is no need to block the right eyepiece to use the photoport.



	PZMIV & PZMTIV Eyepieces and Objectives								
	10x Eyepiece 16x Eyepiece 20x Eyepiece 25x Eyepiece								
Objective	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Mag	Field (mm) (Video Field)	Working Distance
0.32x	1.9x - 16x	106 -13.1 <i>(49.8 - 6.1)</i>	3.2x 25.6x	70.6 - 8.8 (49.8 - 6.1)	3.9x - 32x	55.4 - 6.9 <i>(49.8 - 6.1)</i>	5x - 40x	45.4 - 5.6 (49.8 - 6.1)	296 mm
0.5x	3.1x - 25x	67.7 - 8.4 (31.8 - 3.95)	5x - 40x	45.2 - 5.6 (31.8 - 3.95)	6.2x - 50x*	35.5 - 4.4 (31.8 - 3.95)	7.8x - 62.5x	29 - 3.6 <i>(31.8 - 3.95)</i>	189 mm*
0.63x	3.9x - 31.5x	53.8 - 6.7 (25.3 - 3.15)	6.2x - 50.4x	35.8 - 4.4 <i>(25.3 - 3.15)</i>	7.8x -63x	28.2 - 3.5 (25.3 - 3.15)	9.8x - 78.8x	23 -2.9 <i>(25.3 - 3.15)</i>	149 mm
1.0x (inc)	6.2x - 50x	33.9 - 4.2 <i>(15.9 - 1.97)</i>	9.9x - 90x	22.6 - 2.8 (15.9 - 1.97)	12.4x - 100x	17.7 - 2.2 (15.9 - 1.97)	15.5x - 125x	14.5 - 1.8 <i>(15.9 - 1.97)</i>	80 mm

The Video Field is based on a 1/2-inch CCD (8 mm diagonal) and a 0.5x camera adapter. *Long working distance (LWD) configuration



PZMIV SPI	ECIFICATIONS	
EYEPIECES	WFH 10×	
AUXILIARY LENSES	1×	
ZOOM RANGE	0.62× - 5×	
TOTAL MAGNIFICATION	6.2× - 50×	
ZOOM RATIO	8:1	
FIELD OF VIEW	Ø33.9- Ø4.2 mm	
WORKING DISTANCE	80 mm	
BINOCULAR TUBE	Inclined 45°	
INTERPUPILARY DISTANCE	50 – 75 mm	
DIOPTER ADJUSTMENT	± 5 Diopter	
MICROSCOPE BODY	Rotatable 360°	
OPTIONAL ACCESSORIES		
Eyepieces	16×, 20×, 25×	
Auxiliary lenses	0.32×, 0.5×, 0.63×	
Total Magnification	1.9× - 125×	
Field of view	Ø106 - Ø1.8 mm	
Working Distance	80 – 296 mm	
SHIPPING WEIGHT	23 lb.	

	on Track Stand
	BS PZMTIV Microscope on Boom Stand
	DIG50 PZMTIV Microscope System, including PZMTIV,
US	SBCAM50 USB Computer Camera, 0.5× CCD Camera Coupler,
Z-	LITE Optical Illuminator, Bifurcated Optical Fiber Light Guide
► PZMTIV-	BS-DIG50 PZMTIV Microscope System (80 mm WD), including
PZ	ZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD Camera
	oupler, Z-LITE Optical Illuminator, Bifurcated Optical Fiber Light
	uide, Boom Stand
	BS-LWD-DIG50 PZMTIV Microscope System (189 mm WD),
	cluding PZMTIV, USBCAM50 USB Computer Camera, 0.5× CCD
	amera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical
	ber Light Guide, Boom Stand, 0.5× Objective, 20× Eyepieces for
	ong Working Distance Viewing
	See website for complete configurations.
	* Z-LITE-Z may be substituted in EU countries.
	* Z-LITE-Z May be substituted in EU countries.
OPTION	AL ACCESSORIES/REPLACEMENT PARTS
04:	I Dadia
Optiona	
502000	PZMIV Binocular Body, 10× Eyepieces, 1× Objective, Eye
	Guards
502001	PZMTIV Trinocular Body, 10× Eyepieces, 1× Objective, Eye
	Guards
Optiona	l Stands, Mounts, Bases
502002	76mm Rectanglular Track Stand
502002	Boom Stand (Heavy) without Focus Mount
302004	(requires 502009 Focus Mount for PZMIV)
ENZONE	
502005	Ball Bearing Boom Stand (Heavy) without Focus Mount
F0200C	(requires 502009 Focus Mount for PZMIV)
502006	Boom Clamp Stand (Heavy)
F00000	(requires 502009 Focus Mount for PZMIV)
502009	Universal Focus Mount for 76 mm PZMIV
F04400	(Required for BS, AAC, BBS, and BCS) (5/8" pin)
504123	Extension for Heavy Clamp Stand
504929	LED Microscope Stand, 12.5-in. Post
Optiona	l Eyepieces
502010	10× Wide Field Eyepiece for PZMIV (pair)
502011	16× Wide Field Eyepiece for PZMIV (pair)
502012*	
502013	25× Wide Field Eyepiece for PZMIV (pair)
500264	10× Eyepiece with Reticle (matches 10× eyepiece #502010)
500266	20× Eyepiece with Reticle (matches 20× eyepiece #502012)
	l Objectives
502016	0.32×, Planachromatic Objective (Distortion-free) (278 mm WD
502017*	
502018	0.63×, Planachromatic Objective (Distortion-free) (138 mm WD
502019	1.0×, Planachromatic Objective (Distortion-free) (73mm WD)
Other A	ccessories
500028	1× CCD Camera Coupler
500261	0.35× CCD Camera Coupler, C-Mount (Use with USBCAM33
500262	0.5× CCD Camera Coupler, C-Mount (Use with COLCAM,
300202	USBCAM50)
Z-LITE	
	Z-Lite Fiber Optic Illuminator 7 Lite Fiber Optic Illuminator with Rifurcated Light Guide
Z-LITE-18	1
F00406	and Lenses
500186	Bifurcated Light Guide with Lenses
R-8-8-WF	
502015	Ring Light Adapter for PZMIV
	(For R-8-8-WPI01 Ring Light Guide)
505131	Replacement Trinocular Tube
505132	Replacement Eye Guards for 30 mm Eyepieces (pair)
.1.1	1: 1: 1: 1: 1: 1: 1: 1: 1: 1: 1: 1: 1: 1

*Long working distance (LWD) option with the 502012 and 502017 on a trinocular microscope shows nearly the same scene on the screen as the

ORDERING INFORMATION

PZMIV Microscope on Boom Stand

on Track Stand

on Track Stand

Precision Stereo Zoom Binocular Microscope (Model IV),

Precision Stereo Zoom Trinocular Microscope (Model IV),

PZMIV

PZMIV-BS

PZMTIV

viewer sees in the eyepieces.

Precision Stereo Zoom Microscope

Quality and precision to improve your vision

Features

- Stereo viewing with ample working distance when used with the long working distance (LWD) option
- · Affordably priced

Benefits

- Advanced optics
- True trinocular version available
- Perfect for a boom stand or articulating arm

Applications

- · Tissue Dissection
- Cell injection
- Specimen manipulation
- · Electrode inspection

WPI's third-generation stereo microscope, PZMIII, is an ideal tool for tissue dissection, cell injection, specimen manipulation, electrode inspection, and many other applications that require a magnified, stereo viewing and ample working distance. It offers the leading brand's quality and performance at an affordable price. Advanced optics provide the sharpest image that can only be found among the best of this class. It is superior to many stereo microscopes costing almost twice as much. Zooming is achieved by a spring-loaded knob that is smooth and effortless. The compact size and light weight make it more stable and easily manipulated on the boom stand. A specially designed photo/video module is used in the trinocular version of the microscope (PZMTIII) for photo, video, or digital imaging. In addition, an extensive list of optional accessories is available that can make the PZMIII suitable for almost any bio-research applications requiring a stereo microscope. See next page for options.

The PZMTIII trinocular version of this microscope offers a true trinocular view. Both eyepieces and the trinocular ports are active all the time. There is no need for the right eye piece to switch off. The left eye piece dims slightly, and the right eye piece is 50/50.



PZMTIII-BS-DIG50 includes Boom Stand and USB Digital Camera



503051 Manual Stage — Mounts in the circular opening in the PZMIII base. XY travel distance: 75 x 56 mm. Glass size: 116 x 96 mm. Active diameter: 37.6 mm. Dimensions: 180 x 155 x 27 mm. **Fits 503102 base only.**



PZMIII Eyepieces and Objectives									
	10× Eyepiece 15× Eyepiece 20× Eyepiece 25× Eyepiece								
Objective	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	Mag.	Field (mm) (Video Field)†	Working Distance (mm)
0.3×	2× - 13.5×	114 - 17 (53.6 - 8)	3× - 20.3×	84 -13 (53.6 - 8)	4× -27×	69 - 10.3 (53.6 - 8)	5× - 33.8×	44.8 - 6.7 (53.6 - 8)	287 mm
0.5×	3.4× - 22.5×	69 - 10 (32.4 - 4.7)	5× - 33.8×	51 - 7 (32.4 - 4.7)	6.7× - 45×*	42 - 6.2 (32.4 - 4.7)	8.4× - 56.3×	26.9 - 4.0 (32.4 - 4.7)	177 mm*
0.75×	5× - 33.8×	45 - 7 (21.1 - 3.3)	7.5× - 50.6×	34 - 5 (21.1 - 3.3)	10× - 67.5×	28 - 4.2 (21.1 - 3.3)	12.6× - 84.4×	17.9 - 2.7 (21.1 - 3.3)	117 mm
1.0×	6.7x - 45x	34 - 5 (16 - 4.7)	10× - 67.5×	25 - 3.7 (16 - 4.7)	13.4× - 90×	21 - 3.1 (16 - 4.7)	16.8×-112.5×	13.4 - 2.0 (16 - 4.7)	100 mm
1.5×	10× - 67.5×	23 - 3.4 (10.8 - 1.6)	15× - 101.3×	17 - 2.5 (10.8 - 1.6)	20.1×-135×	14 - 2.1 (10.8 - 1.6)	25.1×-168.8×	9.0 - 1.3 (10.8 - 1.6)	47 mm
2.0×	13.4× - 90×	12 - 2.5 <i>(5.6 - 1.17)</i>	20.1× - 135×	13 - 1.8 <i>(5.6 - 1.17)</i>	26.8×-180×	10 - 1.5 <i>(5.6 - 1.17)</i>	33.5× - 225×	6.7 - 1.0 (5.6 - 1.17)	26 mm

† The video field of view is based on a 1/2-inch (8 mm diagonal) CCD camera and a 0.5× camera adapter. *Long working distance (LWD) configuration



PZMIII SPECIFICATIONS

EYEPIECES	WFH 10×			
ZOOM RANGE	0.67× - 4.5×			
TOTAL MAGNIFICATION	6.77× - 45×			
FIELD OF VIEW	Ø 34 mm – Ø 5 mm			
WORKING DISTANCE	100 mm			
BINOCULAR TUBE	Inclined 45°			
INTERPUPILLARY DISTANCE	Adjustable 47-70 mm			
DIOPTER ADJUSTMENT	±5 Diopter (both eyepieces)			
MICROSCOPE BODY	Rotatable 360°			
ALIVILIA DV LENICEC	Total magnification 2x - 21			

AUXILIARY LENSES Total magnification 2x - 225xBiggest Field of View \emptyset 110 mm

Working Distance

SHIPPING WEIGHT 23 lb.

ODD	EDINIC II	LEODRAGTION	
OKD	EKING II	NFORMATION	

PZMIII Precision Stereo Zoom Microscope (Model III), on Post Stand
PZMIII-BS PZMIII Microscope on Boom Stand
PZMIII-AAC PZMIII Microscope on Articulated Arm with Table Clamp

PZMTIII Precision Stereo Zoom Trinocular Microscope (Model III)

*PZMTIII-DIG50 PZMTIII Microscope System (100 mm WD)

Including PZMTIII, USBCAM50 Computer Camera, 0.5× CCD
Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical
Fiber Light Guide

*PZMTIII-LWD-DIG50 PZMTIII Microscope System (177 mm WD)
Including PZMTIII, Pair of 10× Eyepieces and Eye Guards, Ring
Light Adapter, USBCAM50 Camera, 0.5× CCD Camera Coupler,
Z-Lite Optical Illuminator With Bifurcated Optical Fiber Light
Guide, 0.5× Objective, 20× Eyepieces

*PZMTIII-BS-DIG50 PZMTIII Microscope System (100 mm WD)
Including PZMTIII, USBCAM50 Computer Camera, 0.5× CCD
Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical
Fiber Light Guide, Boom Stand

*PZMTIII-BS-LWD-DIG50 PZMTIII Microscope System (177 mm WD)
Including PZMTIII, USBCAM50 Computer Camera, 0.5× CCD
Camera Coupler, Z-LITE Optical Illuminator, Bifurcated Optical
Fiber Light Guide, 0.5× Objective, 20× Eyepieces for Long
Working Distance Viewing

PZMTIII-BS PZMTIII Microscope on Boom Stand

PZMTIII-AAC PZMTIII Microscope on Articulated Arm with Table Clamp
All PZMIII and PZMTIII microscopes come with 10x eyepieces and

built-in 1x auxiliary lens.

PZMTIII-AAC-DIG50 PZMTIII Microscope System

Including PZMTIII, 10× Eyepieces and Eyeguards, USBCAM50 Camera, 0.5× CCD Camera Coupler, Ring Light Adapter, and Articulated Arm With Table Clamp

* **Z-LITE-Z** may be substituted in EU countries.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

Optional	Bodies
501352	PZMIII Binocular Body, 10x eyepieces (pair), eye guards 13338 Ring Light Adapter NOT included
501379	PZMTIII Trinocular Body, 10× eyepieces & eye guards True trinocular view — 13338 Ring Light Adapter NOT included
Optional	Stands, Mounts, Bases
501353	Fan Post Stand with 76mm Focus Mount
502004	Boom Stand (Heavy) without Focus Mount
502005	Ball Bearing Boom Stand (Heavy) without Focus Mount
502006	Boom Clamp Stand (Heavy)
	(requires 502009 Focus Mount)

(requires 502009 Focus Mount)
 502007 Articulated Arm, 40 cm Table Clamp, without Focus Mount
 502009 Universal Focus Mount, 76 mm ID for PZMIII, PZMIV body
 502163 Wall-Mount Plate, 6" x 6" (or 15.24 cm x 15.24 cm)

503051 Manual Stage for PZMIII503102 76 mm Rectangular Base Post Stand for PZMIII

504123 Post Extension for Heavy Boom Stand 504929 LED Microscope Stand, 12.5-in. Post

Optional Eyepieces

501509	wide Field TOX Eyepieces (pair)
501370	Wide Field 15× Eyepieces (pair)
501371*	Wide Field 20× Eyepieces (pair)
501372	Wide Field 25× Eyepieces (pair)
504128	10× Eyepiece with Reticle (matches 10× eyepiece on PZMIII)
504129	20× Eyepiece with Reticle (matches 20× eyepiece 501371)

Optional Objectives

501373	0.3× Long Working Distance Objective Lens
501375*	0.5× Long Working Distance Objective Lens
501376	0.75× Long Working Distance Objective Lens
501377	1.5× Long Working Distance Objective Lens
501378	2.0× Long Working Distance Objective Lens
	•

Other Accessories

13338	Ring Light Adapter for PZMIII Series			
	(included with all microscope configurations on previous page)			
501381	0.5× C-Mount CCD Camera Coupler			
504596	76 mm Halogen-Halogen Dual Illuminated Track Stand			
504597	Replacement Lamp for 504596			
	the transfer of the transfer teachers			

*Long working distance (LWD) option with the 501371 and 501375 on a trinocular microscope shows nearly the same scene on the screen as the viewer sees in the eyepieces.

26-287 mm

Stereo Microscope with LED Illuminated Base

Articulating mirror ideal for Brightfield and Darkfield applications

Features

- Illuminated base
- Articulating mirror
- Includes base, microscope head and focus mount
- Trinocular version available
- Includes opaque, black/white contrast stage plate

Benefits

- Oblique, transmitted light for low magnification view of tissues, cells or embryo transfer
- Adjustable brightfield and darkfield illumination
- Pseudo DIC illumination

Applications

- Microinjection
- · Nematode research

PZMIII-MI Microscope with Illuminated Base and Articulating Mirror is perfect for microinjection and transfection. It includes our standard stereo microscope head mounted on a research grade Brightfield/Darkfield pole type stand. It has a large stable work surface and a rotatable lens/mirror system which provides transmitted LED intensity illumination.

The sliding mirror is gimbaled, allowing for a full range of movement front to back, as well as rotation. The mirror rotates 360° on one axis and can slide for further lighting effect directionally, front to back. A knob on the right of the base adjusts the mirror and a locking ring holds the desired mirror position.

The articulating mirror lets you vary the microscope illumination from Brightfield LED to Darkfield LED at an appropriate angle. It is an effective tool for viewing live bacteria. At low magnifications, view tissues, cells or embryo transfer where oblique, transmitted illumination is critical.

Microscope options

This unit is sold with our standard PZMIII binocular stereo microscope head. Other options are available. Speak with a specialist today to configure a microscope for your needs





Left: If you want to add a camera or connect with a video monitor, select a PZMTIII trinocular head for your lighted microscope base.

Right: This trinocular microscope head has a camera attached. It also shows an M3301 mounted to the base using a M4C coupler for microinjection studies. A setup like this facilitates classroom or collaborative environments so everyone can see on the remote screen what the researcher is viewing through the microscope..



Systems include

The scope mounting pole diameter is 32 mm.

This platform includes the PZMIII microscope, a focus mount, two stage clips, 94.5 mm glass stage plate, opaque black/white contrast stage plate and an automatic voltage sensing power supply and a 25.6 cm vertical post. The base may be fitted to other microscope heads.

PZMTIII-MI SPECIFICATIONS

DIMENSIONS 13 x 11 x 3.37 in. (33 x 28 x 8.5 cm)

10 in. (25.6 cm) PILLAR HEIGHT GLASS STAGE PLATE 3.75 in. (9.5 cm) diameter

STAGE CLIPS 75 mm clip length, 4 mm pin, paired ILLUMINATION PORT

Built in incident and transmitted variable intensity

ROTATABLE MIRROR Dual-reflection lens/mirror system provides trans-

mitted brightfield/ pseudo-darkfield illumination STAND MANUFACTURER Made in Japan

WARRANTY Limited lifetime warranty

SHIPPING WEIGHT 10 lb. (4.5 kg)

ORDERING INFORMATION

PZMIII-MI Stereo Microscope with LED Illuminated Base Stand with

Articulating Mirror

Includes PZMIII microscope, a focus mount, 2 stage clips, 94.5 mm glass stage plate, black/white stage plate, automatic

voltage sensing power supply

PZMTIII-MI Trinocular Microscope with LED Illuminated Base Stand

with Articulating Mirror

Includes PZMTIII microscope, a focus mount, 2 stage clips, 94.5 mm glass stage plate, black/white stage plate, automatic **Precision Surgical Microscope**

Ideal for small animal surgery

Features

- · Motorized focusing system, allows hands-free operation
- Optional video adapter
- Improved optics 119 lp/mm
- Convenient handles
- · New head tilting mechanism

Benefits

- Light weight, compact and easy to maneuver—weighs only 70 lb.
- Dual bulbs prevent illumination failure during surgery
- · Five magnification steps

Applications

- Small animal surgery
- · Veterinary surgery

WPI's popular precision SurgioScope with five magnification steps is a portable, high quality surgical microscope offering outstanding image quality and value. Incorporating an agile extension arm and excellent working distance objectives, the SurgioScope provides convenient movement and maneuverability necessary for accurate positioning. These important features, together with a high quality optical system, provide sharp image contrast and enhanced large field of vision. The SurgioScope comes fully equipped with a foot-controlled motorized focusing system, normally only found in more expensive surgical microscopes. A unique dual lamp housing enables safe and rapid changing of the lamp during an operation, without the need to power down. The optional video port on the trinocular version permits operational procedures to be monitored or recorded simultaneously using a video recorder and a COLCAM-HD video camera or digital stills with USBCAM50.

SURGIOSCOPE SPECIFIC	CATIONS

TOTAL MAGNIFICATION (F200) 3.2x - 25x± 6 Diopter ADJUSTABLE DIOPTER

ADJUSTABLE INTERPUPILLARY DISTANCE $\,$ min. 50 $\,$ mm - max. 70 $\,$ mm

FYFPIFCF 12 5x FINE FOCUS ADJUSTMENT RANGE 30 mm WORKING HEIGHT (Arm Movement Range Above Floor)

> 89 cm Post.......Focus from 34.5" (88 cm) to 51" (130 cm) above floor * 103 cm Post......Focus from 40.5" (103 cm) to 57" (146 cm) above floor *

* Subtract Working Distance for height above specimen, 103 cm post recommended for F350 objective.

RANGE OF MOTION

Maximum Stretch Radius of Arm .. 870 mm

Vertical Movement Range700-1100 mm

ILLUMINATION

Spot = 42 mm

Dual lamp housing with quick-change spare and internal coaxial fiber

optic cable.

HALOGEN-TUNGSTEN LAMP 12 V, 100 W, with cold reflection OPTIONAL CAMERA COLCAM-HD, COLCAM-HD1080P,

> USBCAM50 (1/2" CCD) USBCAM33 (1/3" CCD)

POWER 110 V, 50-60 Hz, or 220 V, 50-60 Hz

SHIPPING WEIGHT 94 lb. (43 kg)



	ОВ	ECTIVES		
Working Distance	Magnification step	Visual Field of view (mm)	Camera field 1/2" CCD (mm)	Camera field 1/3" CCD (mm)
90 mm	6.4, 10, 16, 26, 40x	25, 15.5, 10, 6, 4	25, 15.5, 10, 6, 4.5	17.5, 11.5, 7, 4.6, 2.8
190 mm	3.2, 5, 8, 13, 20x	50, 31, 20, 12, 8	50, 31, 20, 12, 8	35, 23, 14, 9, 5.5
240 mm	2.6, 4, 6.4, 10.4, 16x	65, 40, 25, 16, 10	63, 40, 25, 16, 10	45, 28, 18, 11, 7
290 mm	2.1, 3.3, 5.3, 8.7, 13x	75, 46.5, 30, 18, 12	75, 46.5, 30, 18, 12	52.5, 34.5, 21, 13.5, 8.3
340 mm	1.8, 2.9, 4.6, 7.4, 11x	91, 57, 36, 22, 14	88, 55, 35, 21, 13	60, 38, 24, 15, 9.5
	90 mm 190 mm 240 mm 290 mm	Working Distance Magnification step 90 mm 6.4, 10, 16, 26, 40x 190 mm 3.2, 5, 8, 13, 20x 240 mm 2.6, 4, 6.4, 10.4, 16x 290 mm 2.1, 3.3, 5.3, 8.7, 13x 340 mm 1.8, 2.9, 4.6,	Distance Magnification step view (mm) 90 mm 6.4, 10, 16, 26, 40x 25, 15.5, 10, 6, 4 190 mm 3.2, 5, 8, 13, 20x 50, 31, 20, 12, 8 240 mm 2.6, 4, 6.4, 10.4, 65, 40, 25, 16x 16, 10 290 mm 2.1, 3.3, 5.3, 75, 46.5, 30, 8.7, 13x 18, 12 340 mm 1.8, 2.9, 4.6, 91, 57, 36,	Working Distance Magnification step Visual Field of View (mm) Camera field 1/2" CCD (mm) 90 mm 6.4, 10, 16, 26, 40x 25, 15.5, 10, 6, 4 25, 15.5, 10, 6, 4.5 190 mm 3.2, 5, 8, 13, 20x 50, 31, 20, 12, 8 50, 31, 20, 12, 8 240 mm 2.6, 4, 6.4, 10.4, 16x 65, 40, 25, 16, 10 63, 40, 25, 16, 10 290 mm 2.1, 3.3, 5.3, 87, 13x 75, 46.5, 30, 75, 46.5, 30, 18, 12 75, 46.5, 30, 18, 12 340 mm 1.8, 2.9, 4.6, 91, 57, 36, 88, 55, 35, 21, 13 88, 55, 35, 21, 13

For additional objectives and specifications, see www.wpiinc.com

	ORDERING INFORMATION
PSMB5N	Binocular SurgioScope, F200 objective (Specify post height)
PSMT5N	Trinocular SurgioScope, beam splitter, standard video adapter, F200 objective (Specify post height)
	Spacify 90 cm or 102 cm pact

Specify 89 cm or 103 cm post Specify line voltage

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

501636	1/2" CS-mount Adapter (requires Beam Splitter 501637)
501637	Beam splitter
504284	F100 Objective
504285	F250 Objective
504286	F300 Objective
504287	F350 Objective
500162	Replacement lamp, 12V, 100 W

Inverted Trinocular Microscope

Excellent for video recording and photography



- Phase contrast at 10X and 20X
- 160 x 260 mm stage
- Infinity Optical System
- Trinocular head for attaching camera
- Abbe condenser
- · Coarse and fine focusing
- · Multiple objectives included

Benefits

- Stage includes inserts for holding standard culture dishes and slides
- · Halogen illumination, offering an image with true to life coloring
- Fixed Stage—optics move during focusing—excellent for patch clamp and brain slice recording

Applications

- Video recording/photography
- Patch clamping
- · Brain slice recording

The INV-101 is an affordable inverted microscope for pathologists, biologists and medical researchers. Able to perform intricate and varied applications for vital tissue cultures, it is perfect for observation and education or for professional research in a clinical lab.

Accessories included: Green, blue and neutral filters, dust cover, immersion oil, manual and warranty card.

TRI	NOC	:UL	.AR	HE/	AD
	=\A/ 10×	/22	ovtra	wido	fiolo

EYEPIECES

CONDENSER

Abbe condenser with 10× / 20× phase

FOCUS

Coarse adjustment range +8 to -3 mm

Fine adjustment: 0.002 mm

STAGE 160×250 mm

INSERTS 35 mm round, 50 mm round, 87×46 rectangular

STAGE DRIVE Coaxial drive controls

MECHANICAL STAGE X-Y coaxial control; 120×78mm range of tra-

verseNosepiece

NOSEPIECE Quintuple nosepiece
OBJECTIVES Infinity optical system
Plan 4× and 40×,
Plan Phase 10× and 20×

6V/30 W halogen bulb

FIXED STAGE Optics move during focusing—excellent for patch

clamp and brain slice recording

Centering eyepiece for phase objective

	ORDERING INFORMATION
INV-101	Trinocular Inverted Microscope
503510	30 mm 10× Eyepiece with 100/10 reticle
503520	Replacement lamp

INV-101

Video camera available

separately

WORLD PRECISION INSTRUMENTS

ILLUMINATION

LED Illuminated Microscopes

Superb optics, durable and high performance microscopes



performance, value and resolution to meet the exacting standards of life science professionals and students. With a newly designed infinity plan optical system, a best-in-class 20 mm field of view and a super-bright 3-watt LED illuminator, the 3000-LED Series provides high contrast images with outstanding resolution, precision design and enhanced illumination.

3000-LED SERIES SPECIFICATIONS OPTICAL SYSTEM Infinity Optical System, f=180 mm, Anti-Mold Siedentopf type, inclined 30°; interpupillary distance

VIFWING HFAD adjustment 48-75 mm

HWF Plan 10× eyepiece, 20 mm field of view with built-in **EYEPIECES** diopter adjustment; a pointer is standard in one eyepiece

NOSEPIECE Rear facing quadruple

Infinity Plan achromat 4× (N.A. 0.10), 10× (N.A. 0.25), **OBJECTIVES** 40×R* (N.A. 0.65), 100×R oil* (N.A. 1.25) are standard

*spring-loaded

Stage size 140 mm × 132mm with X-Y movement range of STAGE

76mm × 50 mm

PHASE CONTRAST Phase sliders for 10×/BF/40× (OPTION) CONDENSER N.A. 1.25 Abbe condenser

ILLUMINATION 3 watt LED with variable intensity control

Dust cover, immersion oil and instruction manual, universal **ACCESSORIES**

power supply 110v-240v

STAND Cast alloy aluminum; coaxial coarse and fine focusing controls.

DIMENSIONS 15.25 × 7.75 × 15.4 in. (387 × 196 × 391 mm) WEIGHT

16 lb. (7.26 kg)

	ORDERING INFORMATION
504221	Binocular, Infinity Plan Achromat 4×, 10×, 40×R and 100×R oil objectives
504443	Trinocular, Infinity Plan Achromat 4×, 10×, 40×R and 100×R oil objectives
504445	Binocular, Plan Phase Contrast 10× and 40×R objectives
504444	Trinocular, Plan Phase Contrast 10× and 40×R objectives
504416	0.50× C-mount Adapter for 1/2" Sensors, adjustable focus
504417	0.35× C-mount Adapter for 1/3" Sensors, adjustable focus

W30S SPECIFICATIONS

Professional-Grade Microscope

Best-seller in universities, medical schools, research laboratories

Features

- · Titanium finished DIN or semi-plan optic
- 30-year antifungal coating

Benefits

· Affordable research grade microscope

Applications

Universities, medical schools, and research laboratories

The W30S professional-grade microscope is a best-seller in universities, medical schools, and research laboratories. Equipped for performance, its features include titanium-finished DIN or Semi-Plan optics and a 30-year

antifungal coating. The W30S is the choice for superior performance at a great price.



HEAD	Binocular (Seidentopf), True Trinocular Inclined 30°, rotates 360° Dual diopter adjustment, Interpupillary distance range 55-75 mm 10×/18 wide field eyepieces
NOSEPIECE	Quadruple forward-facing nosepiece
OBJECTIVES	DIN Plan, antifungal 4×, 10×, 40×, 100×R (oil) Parfocal, parcentric, color-coded
STAGE	Mechanical stage (140 mm × 140 mm) Coaxial drive controls XY Movement: 73 mm × 43 mm
FOCUS	Coarse adjustment: range of 30 mm Fine adjustment: graduation of 2 µm Tension control knob
ILLUMINATION	Moveable Abbe condenser, NA 1.25, Iris diaphragm Variable LED light source (3W bulb) 110 V/220 V switchable electronics
ACCESSORIES INCLUDED	Replacement 0.5 amp fuses, mirror attachment (for field use), blue and green filters, dust cover, immersion oil
DIMENSIONS AND WEIGHT	15" (38 cm) × 9" (23 cm) × 7" (17.8 cm) 14 lb. (6.4 kg)

	1
	ORDERING INFORMATION
W30S-LED	Binocular Microscope
W30ST-LED	Trinocular Microscope
503513	21 mm 10X Eyepiece with 100/10 reticle
500828	Stage Micrometer, 1 mm scale, 200 div. at 10 µm
504606	Stage Micrometer, 50-0.5 mm scale and 10 µm scale

High Definition Camera & Monitoring System

HD imaging for scientific and industrial applications

Features

- High definition camera and monitoring system
- 30 fps (HDMI) and 15 fps (USB 2.0), Aptina sensor and an 4 GB SC card
- Includes HD camera, 1/2.5" (7.182 diagonal), 11.6" HD display screen, 4 GB SD card, USB 2.0 mouse, CaptaVision PC Imaging Software, HDMI Cable, USB 2.0 Cable, 12 V Power Adapter, HDMI Adapter, Y Power Splitter, Mounting Brackets and Hardware
- · Specify line voltage

Benefits

- · Connect your way: use directly, USB to PC or HDMI to a projector
- · Superior performance and color
- Built-in mouse control software for HDMI viewing and recording without a computer
- · Fluorescence imaging with 3D noise reduction
- · Many built-in camera functions

Applications

Scientific and industrial imaging

The new **PRO-300HDS** sets a new standard for excellence in high definition imaging for scientific and industrial applications. This full featured HD camera offers super fast frame rates in video preview, with unrivaled color fidelity and on-board image capturing capability.

PRO-300HDS lets you view and capture images and video directly to the supplied SD card without the need for a computer or separate monitor. The 11.6-in. HD display offers beautiful, crystal-clear image quality, and vibrant, true-to-life color with exceptional viewing from all angles.

Connect your way

PRO-300HDS can be used as a stand-alone system or connected to a PC via USB cable so images can be displayed simultaneously on a PC. For even more flexibility—and ideal for teaching environments—the HDMI architecture allows the device to be connected to an HDMI-enabled projector as well. Each camera is supplied with the PC/Mac-compatible on-board image capture software, as well as full-featured CaptaVision PC image analysis software.

Superior performance and color

PRO-300HDS incorporates dual FPGA processors and unique algorithms that produce perfect color reproduction. It offers full 1080p video preview for HD output with no lag time or compression.

Built-in mouse control software for HDMI viewing and recording without a computer

An external HDMI port allows you to connect directly to an HD monitor for live/real time viewing, capturing and saving of images to an SD card without being connected to a computer. On-board software lets you control the camera with a click of the mouse rather than searching for buttons on the camera, making **PRO-300HDS** the ideal choice for teaching, group presentations or when a computer set-up is not an option.

Fluorescence imaging

The ultra-high signal-to-noise ratio sensor lets you set exposure time from 1 ms to up to 10 seconds and adjust 20 scales of gain value. The incredibly efficient 3D noise reduction performance delivers detailed low light images, making **PRO-300HDS** a great choice for fluorescence applications.





Built-in camera functions

Auto or Manual Exposure Time Capability–You are in full control of exposure and gain. Use the Auto Exposure function or set the exposure time from 1 ms up to 10s and adjust 20 scales of gain value.

 ${\bf 3D}$ ${\bf Noise}$ ${\bf Reduction}\text{-}{\bf Longer}$ exposure times increase image noise. The integrated 3D noise reduction keeps the image clean and sharp.

1080P Video Recording–Just click on the video record icon to start recording 1080P videos at 30fps. The recorded video files are saved directly to the high speed SD card. You can playback videos directly from the SD card.

Magnify (ROI), Rotate and Flip Images-Image operation buttons on the right side of the screen allow you to select a ROI (Region of Interest), as well as flip or rotate an image.

Side by Side Image Comparison–The image comparison function allows you to choose one image, move the image position or select the ROI area to compare with the live image.

Browse Captured Images and Video–Easily browse images on the SD card, zoom in on images or delete them. You can even playback video files saved to the SD card from the Browse feature.

On Board Imaging Tools–No driver is needed when you connect to a PC or Mac via the USB 2.0 camera port. You'll be up and running in no time with the on-board software that features basic imaging tools:

- · Flip Horizontally
- Flip Vertically
- · Zoom in
- 200111111

· Crop

- Zoom Out
- Divison
- CancelInsert a Crossline
- · Compare Images
- Browse Images

ORDERING INFORMATION

PRO-300HDS High Definition Camera with Screen Specify line voltage

Includes HD camera, 11.6" HD display screen, 4GB SD card, USB 2.0 mouse, CaptaVision PC Imaging Software, HDMI Cable, USB 2.0 Cable, 12V Power Adapter, HDMI Adapter, Y Power Splitter, Mounting Brackets and Hardware

Digital Microscope Cameras

Photograph or share your research images

USBCAM133/USBCAM152/USBCAM202

Features

 Ultra-Compact USB cameras with color CCD

Benefits

- Rear mini-USB connector
- Hardware & Software Trigger

Applications

 Image capturing from a microscope



Options

- USBCAM133 with a 1/3-inch color CCD, 6mm diagonal
- USBCAM152 with a 1/2-inch color CCD, 8mm diagonal
- USBCAM202 with a 1/1.8-inch color CCD, 8.9 mm diagonal

All three cameras are cased models with a rear mini-USB connector, hardware and software triggering image capture, digital zoom and a feature-rich, user based menu setup and control. Sentech USB cameras include a SDK, DirectX, Twain and Linux driver, as well as the Sentech Viewing Software.

SPECIFICATIONS USBCAM133 USBCAM152 USBCAM202 IMAGE SENSOR 1/3" Interline SXVGA color 1/2" Interline SXGA color 1/1.8" Interline UXGA color progressive CCD progressive CCD progressive CCD CELL SIZE 3.75 (H) x 3.75 (V) µm, 4.65 (H) x 4.65 (V) µm, 4.40 (H) x 4.40 (V) µm, 8.9 6mm diagonal mm diagonal 8mm diagonal SCANNING SYSTEM Progressive Progressive Progressive 1360 (H) x 1024 (V) 1600 (H) x 1200 (V) 1280 (H) x 960 (V) RESOLUTION MIN. SCENE 11 Lux at F1.2 18 Lux at F1.2 7.7 Lux at F1.2 ILLUMINATION 22.4 Frames per Second 19.26 Frames per Second 15.3 Frames per Second SPEED ELECTRONIC Auto/Manual (software Auto/Manual (software Auto/Manual (software SHUTTER selectable) selectable) selectable) GAIN Auto/Manual (software Auto/Manual (software Auto/Manual (software selectable) selectable) selectable) GAMMA Manual (software Manual (software Manual (software selectable) selectable) selectable) WHITE BALANCE Auto/Manual/One shot Auto/Manual/One shot Auto/Manual/One shot (software selectable) (software selectable) (software selectable) INPUT/OUTPUT USB 2.0 High Speed USB 2.0 High Speed USB 2.0 High Speed POWER +5 VDC through USB +5 VDC through USB +5 VDC through USB connector, < 300 mA connector, < 420 mA connector, < 450 mA DIMENSIONS 28 (W) x 28 (H) x 37 (D) mm 28 (W) x 28 (H) x 42 (D) mm 28 (W) x 28 (H) x 42 (D) mm (excluding connector) (excluding connector) (excluding connector) LENS MOUNT CS mount C mount WFIGHT Approximately 45 g Approximately 45 g Approximately 45 g INTERFACE CON-USB: mini-B USB connector USB: mini-B USB connector USB: mini-B USB connector **NECTOR** IO signal: 6 pin connector 10 signal: 6 pin connector 10 signal: 6 pin connector (HR10A-7R-6PB or (HR10A-7R-6PB or (HR10A-7R-6PB or equivalent) equivalent) equivalent) RoHS RoHS Compliant RoHS Compliant RoHS Compliant

	ORDERING INFORMATION
USBCAM133	Digital Microscope Camera
USBCAM152	Digital Microscope Camera
USBCAM202	Digital Microscope Camera
504570	Replacement USB Cable, 3 m (10 ft.)

USBCAM33/USBCAM50

Features

· Automatic and manual controls

C/CS mount

Software included

Benefits

Connect via USB

Applications

Image capturing from a microscope



Record images directly to your computer. These digital microscope cameras offer fle

microscope cameras offer flexibility, with a range of configurations for image capture, a choice of mount option (C or CS) and file output alternatives. Since both cameras connect via the USB port, installing the image capture software is simple. Either camera can be used on WPI's stereo microscopes PZMTIV, PZMTIII, compound microscopes W30ST and GPL-T and also the PSMT5 Surgical Microscope. Choose from the one third-inch CCD with 1024×768 resolution and 30 frames per second (USBCAM33) or one half-inch CCD with 1280×960 resolution and 15 frames per second (USBCAM50).

These cameras include IC Imaging Control software that has:

- · Real-time video preview
- Text and graphics can be drawn on a live video stream
- · Scroll and Zoom
- Acquisition of single frames
- Capture pause, for intermittent image capture
- Timestamps

SPECIFICATIONS		
	USBCAM33	USBCAM50
IMAGE SENSOR	1/3" Sony CCD, progressive scan	1/2" Sony CCD, progressive scan
MAX RESOLUTION	1024 x 768	1280 x 960
SIZE	4.65 μm x 4.65 μm, 6mm diameter	4.65 μm x 4.65 μm, 8mm diameter
SPEED (PC DEPENDENT)	30 fps, 15 fps, 7.5 fps or 3.75 fps	15 fps, 7.5 fps or 3.75 fps
SENSITIVITY	0.5 lux @ 1/15 s	0.5 lux @ 1/7.5 s
EXPOSURE, SHUTTER CONTROL, WHITE BALANCE	Automatic/Manual	Automatic/Manual
INTERFACE	USB 2.0 cable	USB 2.0 cable
SYSTEM REQUIREMENT	Windows Vista (32 & 64 bit), Windows 7 (32 & 64 bit), Windows 8 (32 & 64 bit) or Windows 10 (32 & 64 bit)	Windows Vista (32 & 64 bit), Windows 7 (32 & 64 bit), Windows 8 (32 & 64 bit) or Windows 10 (32 & 64 bit)
SOFTWARE	IC Imaging Control Software	IC Imaging Control Software
LENS MOUNT	C/CS-Mount	C/CS-Mount
CAMERA BODY	50.6 x 50.6 x 50 mm	50.6 x 50.6 x 50 mm
WEIGHT	265 g (9.5 oz.)	265 g (9.5 oz.)

	ORDERING INFORMATION
USBCAM33	Digital Microscope Camera, 1/3-in. CCD
USBCAM50	Digital Microscope Camera, 1/2-in. CCD
503536	Cable, USB Extension (male-female)

Color CCD Video Camera for Microscopy

Watch live video images

COLCAM-HD1080P

Features

- 16:9 Aspect Ratio, 1:1, no scaling
- DVI signal output via HDMI cable
- Improved design of COLCAM-HD with better low light sensitivity
- Low cost alternative to existing HD and 3 CCD cameras, with outstanding image quality



COLCAM-HD1080P

Benefits

• Directly connect to HDTV, no PC required and no software to load

Applications

Microscopy video output to HDTV

The **COLCAM-HD1080P** is a CMOS-based camera that outputs a true HD 1080P or 720P image at 60fps in the 16x9 format, with the capability to program individual DSP profiles accessed via remote hand held controller.

If still image or HDTV recording is required, contact WPI.

COLCAM-HD1080P SPECIFICATIONS 1/2.33" 14 MP CMOS **IMAGER** Progressive HD ACTIVE PICTURE ELEMENT 1920 (H) x 1080 (V) CHIP SIZE 7.8 mm diagonal MINIMUM SCENE ILLUMINATION 650 mV/Lux-sec SYNC SYSTEM Internal DVI 1.0 conformity 1080P VIDEO OUTPUT RGB, 1920 H x 1080 V 60/50/30 Hz GAIN AGC * SHUTTER SPEED Auto * GAIN AGC or Fixed gain * **GAMMA** presets or manual gamma * WHITE BALANCE Auto/Manual/Push-to-set * **POWER** 12V power jack; 5.5x2.1 mm **DIMENSIONS** 40 (W) x 45 (H) x 41.1 (D) mm OPTICAL FILTER IR cut filter included LENS MOUNT C/CS mount VIDEO OUTPUT HDMI connector WEIGHT **ROHS** RoHS compliant

ORDERING INFORMATION

COLCAM-HD1080P Color Video Camera, 1920 x 1080 HD resolution

Includes 3m DV Cable, Power Supply, C/CS Mount

504136 3 meter HDMI Cable

C-Mount Eyepiece Adapters



For 1/3-inch (6 mm diagonal) and 1/2-inch (8 mm diagonal) video cameras and eyepiece camera conversion, this lens and its accessories make it possible to connect a typical video or C-mount camera to almost any microscope on the market. The lens fits right into the ocular socket of standard 23.2mm microscopes and the 30 mm adapter allows for use on the typical stereo zoom microscope. If you already have a trinocular microscope you can add the included C-Adapter to the top of an existing 1X C-mount (no lens) adapter.

ORDERING INFORMATION

503097

C-Mount to Eyepiece Adapter Kit, 0.45X for 1/3-in. and 1/2-in. Video Cameras, 30 mm Stereo Adapter, 1X C-Mount Adapter

SLR Digital Camera-to-Microscope Eyepiece Adapter



This adapter connects T-mount SLR digital cameras to almost any microscope on the market. The adapter is built to 23.2 mm ocular tubes that are found on most high magnification (upright, inverted, standard) microscopes. The 30 mm adapter allows mounting on most stereo zoom microscopes that use 30 mm oculars. If you already have a trincocular microscope, you can add this adapter to the top of an existing 1X C-mount adapter.

The 2X magnification of this microscope adapter yields an approximate 65% field of view from the visual field as measured on a Canon 80D Digital camera. (CCD Sensor size = $22.7 \times 15.1 \text{ mm}$). 35 mm film reference size is $24 \times 36 \text{mm}$.

Please contact your camera dealer for a suitable T-mount to bayonet adapter for your camera.

ORDERING INFORMATION

503099

C-Mount to Eyepiece Adapter Kit, 2.0x for SLR digital Cameras Includes 30 mm Stereo adapter, 1X C-mount adapter

^{*} Selectable via the UART communication

High Intensity Fiber Optic Illumination Source

Uninterrupted, directed light for microscopes and other applications

Features

- · Reliable, uninterrupted high-intensity light
- Use with microscopes
- · Intensity controlled by rotary knob
- · Use with ring light, or single/bifurcated guides, sold separately

Benefits

· Contains full spectrum of visible light

Applications

- · General laboratory use
- Microscopy illumination

The **Z-LITE** Fiber Optic Illuminator provides reliable, uninterrupted high-intensity light for microscopes. **Z-LITE** allows a continuous range of subdued or concentrated lighting controlled by a rotary dimmer on the front panel. **Z-LITE** may be used with a ring light and single or bifurcated flexible fiber bundles, enabling the light beam to be placed exactly where it is needed. Forced air cooling prolongs the lamp life. Lamp color temperature is 3350°K. An interlock switch automatically cuts off power when the front panel is opened to replace the bulb.

Light guides available separately



Z-LITE SPECIFICATIONS

LAMP 150 W quartz halogen (EKE lamp)
SIZE 30.5×25×25 cm (12×10×10 in.)
POWER 115 VAC, 50/60 Hz, 3 A
WEIGHT 5.9 kg (13 lb.)

	ORDERING INFORMATION
Z-LITE-186	Z-Lite & Bifurcated Light Guide (115 V, 60 Hz, beige case)
Z-LITE-Z186	Z-Lite & Bifurcated Light Guide (230 V, 50 Hz, black case) *
Z-LITE	Z-Lite Fiber Optic Illuminator (115 V, 60 Hz, beige case)
Z-LITE-Z	Z-Lite Fiber Optic Illuminator (230 V, 50 Hz, black case) *

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

EKE

500186	Bifurcated Light Guide (with lenses)
504930	Flexible Light Guide with Focusing Lens, 18" (46cm)
504931	Flexible Light Guide with Focusing Lens, 24" (61cm)
	Ring Light Guide for PZM and PZMIII Series
Ring LigI	ht Guide requires adapter #13338 for use with PZM, PZMII and
PZMI	II, included with each PZMIII and PZMIV microscope system.
13338	Ring Light Adapter (48 mm Ø) for PZM, PZMII, PZMIII
502015	Ring Light Adapter for PZMIV
EJA	Replacement Halogen Lamp, 150W, 3350°K, 40-hour

Replacement Lamp, 150W, 3250°K, 200-hour

* Not available in EU countries







Features

- "White" light illumination 72 LED bulbs
- Maximum opening 61 mm
- Ring light is divided into four areas and each area is turned on and off separately
- · Brightness adjustable
- ESD safe
- Power supply AC 90-264V, 50/60 Hz, US plug only

Benefits

- Direct top illumination plus four zones for shadows
- Low cost

504134

- Less bulky than Halogen
- · Cold light with no heat

Applications

Stereo microscope use





Microdissection Instruments



Quality you can count on

WPI surgical instruments are manufactured with the highest quality materials and craftsmanship to provide you with reliable instruments at a cost-effective price. We offer a wide variety of instruments, including surgical kits. Whether you are looking for a pair of quality European surgical instruments, marine grade stainless steel forceps or precision scissors, the quality of our surgical instruments is backed by our 100% satisfaction guarantee.



Swiss Tweezers

Our Swiss surgical tweezers are manufactured from fine Swiss steel. We offer a variety of tweezers and forceps.

WPI SWISS TWEEZERS

- · Used in microsurgical procedures to grasp fine tissue
- Sharp tips may puncture tissue
- 10.5 cm (4.1 in.) Long
- Extra fine 90° angled tips

504508	Stainless Steel
504509	Acid-Resistant/Anti-magnetic



WPI SWISS TWEEZERS

- 11 cm (4.3 in.) Long
- · Handcrafted to perfection in Switzerland
- For use in general laboratory and medical applications, microscopy, precision manufacturing, electronic industry, watch making
- 0.1 x 0.06 mm Straight tips

504506	Stainless Steel
504507	Acid-Resistant/Anti-magnetic





504506 504507

WPI SWISS TWEEZERS

- 11 cm (4.3 in.) Long
- Diamond Like Carbon coated tips to improve hardness, toughness and wear characteristics

504524	Ni-Cr-Mo Superalloy
504516	Ni-Cr-Mo Superalloy, Diamond Coated Tips
504515	Anti-magnetic, Diamond Coated Tips





WPI SWISS TWEEZERS

- 11.5 cm (4.5 in.) Long
- Stainless Steel

504505 Acid-Resistant/Anti-magnetic, Curved Tips





WPI SWISS CUTTING TWEEZERS

- 11.5 cm (4.5 in.) Long
- Angled blades
- Stainless Steel

• Stalliless	s Steel
504744	Stainless Steel
504745	High Precision, Tungsten Carbide Tips







504745

WORLD PRECISION INSTRUMENTS

Dumont Forceps

High quality Dumont tweezers are available in a variety of metals.

Dumostar – Anti-magnetic. Dumostar is more elastic and more corrosion-resistant than the best stainless steel.

Dumoxel – Extremely flexible. Dumoxel is 95% anti-magnetic and stain resistant.

Stainless Steel - This is the standard metal for most surgical instruments.

MINI DUMONT #M5S

- 8.2 cm (3.25 in.) Long
- 0.09 x 0.05 mm Tips
- Stainless Steel

501764 Mini Dumont #M5S



DUMONT #3C

- 11 cm (4.3 in.) Long
- 0.08 x 0.04 mm Tips
- · Non-magnetic
- · Non-corrosive
- Dumostar Steel

500064 Dumont #3C



DUMONT #4

- 11 cm (4.3 in.) Long
- 0.06 x 0.02 mm Tips
- · Non-magnetic
- · Non-corrosive
- · Dumostar Steel

500339 Dumont #4



DUMONT #4

- 11 cm (4.3 in.) Long
- · Stainless Steel

500340	0.13 x 0.08 mm Tips
500231	0.02 x 0.06 mm Tips



DUMONT #5

- 11 cm (4.3 in.) Long
- · Non-magnetic
- Non-corrosive
- Dumostar Steel

500085	0.025 x 0.015 mm Tips
500233	0.1 x 0.06 mm Tips





500340

500231

Dressing Forceps


Dressing forceps are used when dressing wounds. They hold gauze and other dressings. They may also be used during wound debridement to remove infected or necrotic tissue or debris from the wound. They may also be used for suturing.

ADSON FORCEPS

- Used to grasp tissue during stapling of the skin
- 12 cm (4.75 in.) Long
- Serrated
- · Stainless Steel

14226	1.5 mm Tips
14226-G	1.5 mm Tips, German
500222	Tungsten Carbide Jaws
500222-G	Tungsten Carbide Jaws, German

Standard and German versions may differ slightly in appearance.



IRIS FORCEPS

- Used for grasping tissue and dressing applications in microscopy
- 10 cm (4 in.) Long
- 0.8 mm Tips
- 1x2 Teeth
- · Stainless Steel

15916	Straight
15916-G	Straight, German
15917	Curved
15917-G	Curved, German

Standard and German versions may differ slightly in appearance.









15916 15917

Hemostatic Forceps

Ring forceps (also called hemostatic forceps) are hinged and look like ring scissors. Frequently, hemostatic forceps have a locking mechanism, called a ratchet, which is used for clamping. The jaws of the locking forceps gradually come together as each increment of the ratchet is employed.


- Used to compress blood vessels or other tubular structures to obstruct the flow of blood or fluids
- Used for small, delicate vessels
- Full serrated jaws
- Used in ophthalmology, biopsy, skin grafting or general surgery and lab procedures

HALSTED MOSQUITO HEMOSTATIC FORCEPS

- Serrated
- Stainless Steel

15920	12.5 cm (5 in.), Straight
15920-G	12.5 cm (5 in.), Straight, German
15921	12.5 cm (5 in.), Curved
15921-G	12.5 cm (5 in.), Curved, German
501705	9 cm (3.5 in.), Straight
501705-G	9 cm (3.5 in.), Straight, German
501291	9 cm (3.5 in.), Curved
501291-G	9 cm (3.5 in.), Curved, German

Standard and German versions may differ slightly in appearance.



Titanium Forceps & Needle Holders

100% anti-magnetic, corrosion-resistant, lightweight and strong, titanium alloy is ideal for biological and medical applications. Titanium has the tensile strength of carbon steel and is completely resistant to corrosion from nitric acid, chloride, salt water, and industrial and organic chemicals. Titanium is more flexible and 40% lighter than lnox. When heated or cooled, the dimensions of titanium alloy change less than half of what stainless steel alloys will, making titanium surgical instruments much more durable. Titanium is stain-free and temperature resistant up to 430°C. Titanium tools are the premium choice for corrosive environments or MRI applications.

- 100% non-corrosive (great for sea water procedures)
- · 40% lighter than stainless steel (reduces hand fatigue)
- 100% non-magnetic (MRI compatible)
- · Anodized, non-glare blue finish

FORCEPS

- 8.5 cm (3.3 in.) Long
- 12 mm Tips
- 4.5 mm Tying platform
- · Straight
- · Titanium

WP5000 Forceps



FORCEPS WITH TYING PLATFORM

- 8.5 cm (3.3 in.) Long
- 12 mm Tips
- 5 mm Tying platform
- Straight with 0.12 mm teeth (1x2)
- Titanium

WP2000 Forceps with Tying Platform



TROUTMAN NEEDLE HOLDER

- 11.5 cm (4.5 in.) Long
- Curved
- 10 mm Jaws
- Blunt
- Titanium


WP3620 Troutman Needle Holder



NEEDLE HOLDER, SERRATED JAWS

- 16 cm (6.3 in.) Long
- Straight
- Serrated Jaws
- Titanium

503391 Needle Holder



WORLD PRECISION INSTRUMENTS

Spring Scissors

Vannas scissors – small, sharp blades designed for researchers working in tight spaces.

- · May be used with either left or right hand
- · Designed for neurosurgical, vascular and ophthalmic procedures
- Cuts by shearing, reduces crushing of tissue
- Fine tips are ideal for use in very restricted spaces
- Curved tips are useful to avoid cutting of underlying tissue
- · Heavier blades are helpful when a thick tissue or vessel is to be cut
- Length of the scissor tips can help with the depth of the incision: Shorter tips for near-to-surface incision Longer tips for deep incision

VANNAS SCISSORS

- Small, long blades for uses in ophthalmology, neurosurgery, microsurgery
- Super fine
- 8 cm Long
- · 3 mm Blades
- 0.015 mm Tips
- · Stainless Steel

501778	Straight	
501839	Curved	



VANNAS SCISSORS

- Diamond shaped tips
- 9 cm (3.5 in.) Long
- 5 mm Blades
- 0.1 mm Tips
- Stainless Steel

14003	Straight
14003-G	Straight, German
Standard and	German versions may differ slightly in annearance

Standard and German versions may differ slightly in appearance.



NOYES SCISSORS

- 12 cm Long
- 15 mm Tips
- Originally designed to remove the iris (iridectomy)
- · Fine for multi-purpose use
- · Best for trimming tissue from nerves

500228	Sharp/Sharp Tips, Straight
500228-G	Sharp/Sharp Tips, Straight, German
501236	Sharp/Sharp Tips, Curved
501236-G	Sharp/Sharp Tips, Curved, German
503305	Sharp/Blunt Tips, Straight
503306	Blunt/Blunt Tips, Straight
501237	14 cm, 20 mm Tips, Straight
501238	14 cm, 20 mm Tips, Curved

Standard and German versions may differ slightly in appearance.



501238



503306

501237

KATENA-VANNAS SCISSORS

- 11 cm (4.3 in.) Long
- · Stainless Steel

504024	7 mm Thin Blades, Straight
504025	7 mm Thin Blades, Curved
504026	10 mm Thin Blades, Straight
504027	10 mm Thin Blades, Curved



MICRO SCISSORS – DIAMOND LIKE CARBON COATED BLADES

Diamond like carbon coated blades increase the longevity of the cutting edge and improve the quality of the cut. Used to reduce wear on razor blades and industrial metal cutting tools, this coating is now being used in biomedical applications. Some studies have shown that even after double the normal usage of the cutting implement, no blade wear was detected.

- · Hardness of coated surface resists wear and increases life of cutting edge
- Anti-reflective coating resists glare and assists in better view of surgical field
- Coating is corrosion resistant
- Reduced friction at cutting site decreases necrosis of surrounding tissue and bone
- 12 cm (4.75 in.)
- Sharp coated tips
- 8 mm Blades
- · Stainless Steel,

503365	Straight	
503364	Curved	




SPRING SCISSORS

- 12 cm (4.75 in.) Long
- 12 mm Extra fine, long blades
- Stainless Steel

14125	Straight
14125-G	Straight, German
14126	Curved

Standard and German versions may differ slightly in appearance.



Titanium Scissors

100% anti-magnetic, corrosion-resistant, lightweight and strong, titanium alloy is ideal for biological and medical applications. Titanium has the tensile strength of carbon steel and is completely resistant to corrosion from nitric acid, chloride, salt water, and industrial and organic chemicals. Titanium is more flexible and 40% lighter than lnox. This reduces hand fatigue during long procedures. When heated or cooled, the dimensions of titanium alloy change less than half of what stainless steel alloys will, making titanium surgical instruments much more durable. Titanium is stain-free and temperature resistant up to 430°C. Titanium instruments are the premium choice for corrosive environments or MRI applications.

- 100% non-corrosive (great for seawater procedures)
- 40% lighter than stainless steel (reduces hand fatigue)
- 100% non-magnetic (MRI compatible)
- · Anodized, non-glare blue finish

VANNAS SCISSORS

- 8.2 cm (3.2 in.) Long
- Sharp 9.5 mm tips
- Titanium

WP5050	Straight	
WP5070	Angled Up	



IRIS SCISSORS

- 9 cm (3.5 in.) Long
- Curved
- Sharp 11 mm tips
- Titanium

WP4410 Iris Scissors



VANNAS CAPSULOTOMY SCISSORS

- 9.5 cm (3.7 in.) Long
- Angled
- Sharp
- 10.5 mm Tips
- Titanium

WP5070L Vannas Capsulotomy Scissors



VANNAS CAPSULOTOMY SCISSORS

- 9.5 cm (3.7 in.) Long
- Curved
- Sharp
- 12.5 mm Tips
- Titanium

WP5060 Vannas Capsulotomy Scissors



Ring Scissors

Our standard scissors all have ring handles. You can typically recognize the construction style of scissors based on the following color code:

- Black handled rings indicate SuperCut scissors (with one serrated blade).
- Gold handles indicate that the scissors have tungsten carbide inserts.
- One gold handle and one black handle indicates that the scissors have serrations in one blade (black handle) and tungsten carbide on the other (gold handle).
- · Blue instruments are typically titanium.

IRIS DISSECTING SCISSORS WITH PROBE TIPS

- · Multi-purpose
- Ideal for detailed dissection of delicate tissues and can be used in suture removal
- 10 cm (4 in.) Long
- Blunt probe tips
- Stainless Steel

500366	Straight
500366-G	Straight, German
500367	Curved
500367-G	Curved, German

Standard and German versions may differ slightly in appearance.




OPERATING SCISSORS

- · Recommended for cutting and dissecting tissue
- Also used for cutting surgical drapes
- 16 cm (6.25 in.) Long
- Stainless Steel

501223	Sharp/Blunt, Straight
501223-G	Sharp/Blunt, Straight, German
501224	Blunt/Blunt, Straight
501224-G	Blunt/Blunt, Straight, German
501225	Sharp/Sharp, Straight
501225-G	Sharp/Sharp, Straight, German
501226	Sharp/Blunt, Curved
501226-G	Sharp/Blunt, Curved, German
501227	Blunt/Blunt, Curved
501227-G	Blunt/Blunt, Curved, German
501228	Sharp/Sharp, Curved
501228-G	Sharp/Sharp, Curved, German

Standard and German versions may differ slightly in appearance.



SuperCut Scissors

Our black handled surgical scissors designate our SuperCut scissors. These scissors have one on razor sharp blade and one micro-serrate blade. This construction offers a couple of advantages for surgeons.

- · Sharp edge gives a clean cut with minimal tissue damage
- · Serrated edge actually holds the tissue to prevent it from slipping while you are making an incision


IRIS SUPERCUT SCISSORS

 One edge micro serrated, one edge honed to the sharpness of a knife edge—prevents tissue slippage

اد	Stai	229	Stain	•

14218	10 cm (4 in.), Straight
14218-G	10 cm (4 in.), Straight, German
14219	10 cm (4 in.), Curved
14219-G	10 cm (4 in.), Curved, German
503259	11.5 cm (4.5 in.), Straight
503260	11.5 cm (4.5 in.), Curved
14225	12.5 cm (5 in.), Straight
14225-G	12.5 cm (5 in.), Straight, German
14224	12.5 cm (5 in.), Curved
14224-G	12.5 cm (5 in.), Curved, German

Standard and German versions may differ slightly in appearance.



Tungsten Carbide SuperCut Scissors

Our black handled surgical scissors designate our SuperCut scissors. These scissors have one on razor sharp blade and one micro-serrate blade. Surgical scissors with one black handle and one gold handle are SuperCut scissors with tungsten carbide inserts. Tungsten carbide is harder than steel, which means the instruments last longer. The SuperCut edges offer a couple of advantages for surgeons.

- Sharp edge gives a clean cut with minimal tissue damage
- · Serrated edge actually holds the tissue to prevent it from slipping while you are making an incision

IRIS SUPERCUT SCISSORS, TUNGSTEN CARBIDE

- Gold blade has tungsten carbide insert black blade is serrated
- · Produces clean cut without tearing the tissue
- Used for cutting delicate or heavy tissue
- 11 cm (4.3 in.) Long
- Sharp tips
- Tungsten Carbide blades
- Stainless Steel

501263	Straight
501263-G	Straight, German
501264	Curved
501264-G	Curved, German

Standard and German versions may differ slightly in appearance.



Ring-Handle Needle Holders

OLSEN-HEGAR NEEDLE HOLDER WITH SUTURE SCISSORS

- Straight
- Serrated
- Tungsten Carbide jaws
- · Stainless Steel

501312	12 cm (4.75 in.)
500227	14 cm (5.5 in.)
500227-G	14 cm (5.5 in.), German
501725	16.5 cm (6.5 in.)
501725-G	16.5 cm (6.5 in.), German

Standard and German versions may differ slightly in appearance.



OLSEN-HEGAR NEEDLE HOLDER

- 14 cm (5.5 in.) Long
- · Straight
- Serrated jaws with suture scissors
- Stainless Steel

500023	Olsen-Hegar Needle Holder
500023-G	Olsen-Hegar Needle Holder, German

Standard and German versions may differ slightly in appearance.



CRILE-WOOD NEEDLE HOLDER

- · Straight
- Serrated
- · Tungsten Carbide jaws
- · Stainless Steel

500224	14.5 cm (5.75 in.)
500224-G	14.5 cm (5.75 in.), German
503741	15 cm (5.9 in.), Left-Handed

Standard and German versions may differ slightly in appearance.



MAYO-HEGAR NEEDLE HOLDER

- 15 cm (6 in.) Long
- Straight
- Stainless Steel

V503382	Mavo-Hegar Needle Holder



Black Instruments

JAFFE FORCEPS

- 10 cm (3.9 in.)
- · Curved 0.52 mm Tip
- · Stainless Steel with black coating

WPB505210 Jaffe Forceps, Black



IRIS SUPERCUT SCISSORS, TITANIUM POWDER COATED

- More resistance to rust and body fluids
- · Anti slippage serrated blade
- Used in microsurgery
- 11 cm (4.3 in.) Long
- · Anti-reflective
- · One serrated blade
- Stainless Steel

504487	Straight	
504500	Curved	

BLACK EPOXY COATED TWEEZERS WITH PATENTED EZ LEVER

- 11 cm (4.25 in.) Long
- Strong blades
- · Stainless Steel with black coating

WPB315AB 4.5 mm from tip to bend WPB315AS 6.5 mm from tip to bend







WPB315AB

WPB315AS

Skin Stapler

The **WPSS01** skin stapler is an affordable, reliable solution for ideal veterinary wound closures.

- Made in the U.S.A.
- · Designed specifically for veterinarians
- · Clear handle makes it easy to spot contaminants
- Uses reloadable cartridges containing 35 stainless steel staples
- · Autoclavable handle
- Improves patient safety by considerably reducing surgery and anesthesia time
- Angled tip design helps you maintain a clear view of the incision line at all times
- · Lightweight but durable with a comfortable ergonomic grip design
- Simple to use simple to reload
- Designed to reduce problems with rotated staples

WPSS01	Skin Stapler (one staple cartridge included)
WPSS02	Staple Cartridge (35 staples)





Stainless steel skin staples are 3.6 mm high by 6.9 mm wide after closure.

Reusable Biopsy Punches with Plunger

Features

- · Individual Sterile Packs, Limited Reuse
- · Tip material: 304 stainless steel
- · Sterilization: ETO (ethyline oxide)

Applications

- Electrophysiology—Patch recording (various regions of the specimen can be punched prior to separation).
- Biomedical Pharmacology—Collect samples from the brain neurotransmitters for analyzing the metabolic changes in response to active substances, or skin samples to analyze the response of anti aging medications.
- Anatomy—Sample brain regions for analysis of neurotransmitter concentrations of mRNA levels.
- Forensic Sampling—WPI punches can easily cut, retrieve, store and eject material like skin, gel, films, paint chips and paper. (Better



504529

samples are collected using one of our self-healing cutting mats $\bf 504620$ or $\bf 504621.)$

 Microfluidics—Punch flow inlets in microfluidics chambers molded of PDMS (polydimethylsiloxane).

Surgeons and laboratory technicians use biopsy punches and curettes to remove tissue samples. These punches are for use in research laboratories.

504638	ID 0.35 mm, OD 0.63 mm, Wall 0.15 mm
504528	ID 0.50 mm, OD 0.80 mm, Wall 0.15 mm
504529	ID 0.75 mm, OD 1.07 mm, Wall 0.15 mm
504646	ID 1.0 mm, OD 1.26 mm, Wall 0.15 mm
504530	ID 1.2 mm, OD 1.5 mm, Wall 0.15 mm
504647	ID 1.5 mm, OD 1.9 mm, Wall 0.20 mm
504531	ID 2.0 mm, OD 2.4 mm, Wall 0.20 mm
504648	ID 2.5 mm, OD 2.9 mm, Wall 0.20 mm

504649	ID 3.0 mm, OD 3.4 mm, Wall 0.20 mm
504650	ID 3.5 mm, OD 3.9 mm, Wall 0.20 mm
504651	ID 4.0 mm, OD 4.4 mm, Wall 0.20 mm
504652	ID 5.0 mm, OD 5.5 mm, Wall 0.25 mm
504653	ID 6.0 mm, OD 6.5 mm, Wall 0.25 mm
504654	ID 7.0 mm, OD 7.5 mm, Wall 0.25 mm
504655	ID 8.0 mm, OD 8.5 mm, Wall 0.25 mm

The Rapid-Punch consists of a razor sharp cutting tip designed to cut, retrieve and store cored samples from source materials such as tissue, gels, paper, cloth, leaves, paint chips, films or other thin substrates. The Rapid-Punch is recommended as the special DNA sampling tool. It is ideal for forensic and most laboratory sampling purposes. Designed for heavy, numerous and precise sampling operations. It is durable, lasting up to 2000 punches.

Reusable Rapid-Punch Biopsy Kit

504639	Reusable Rapid Punch Kit with 0.5 mm Tip
WP1010	Reusable Rapid Punch Kit with 1.0 mm Tip
WP1212	Reusable Rapid Punch Kit with 1.2 mm Tip
WP2020F	Reusable Rapid Punch Kit with 2.0 mm Tip
WP3030	Reusable Rapid Punch Kit with 3.0 mm Tip
504640	Replacement Tip for Rapid Punch size 0.5 mm
504641	Replacement Tip for Rapid Punch size 1.0 mm
504642	Replacement Tip for Rapid Punch size 1.2 mm
504643	Replacement Tip for Rapid Punch size 2.0 mm
504644	Replacement Tip for Rapid Punch size 3.0 mm
504735	Replacement Plunger for Rapid Punch, 0.5 mm
504736	Replacement Plunger for Rapid Punch, 1.0 mm
504737	Replacement Plunger for Rapid Punch, 1.2 mm
504738	Replacement Plunger for Rapid Punch, 2.0 mm
504739	Replacement Plunger for Rapid Punch, 3.0 mm



Replative to the second

Rapid-Punch with 0.5 mm tip

Clips and Clamps

Save time in surgery using our reflex clips. They are a fast, effective suture alternative for wound closure. Choose the 9 mm clips for use in rats or 7 mm clips for mice and young rats. The spring-loaded clip applier is easy to use. You get minimal tissue trauma and maximum holding power. Use the clip remover to extract the clips when the wound heals. The clips, applier and remover are all autoclavable.

Vessel clips are used to clamp arteries or veins during surgical procedures. These vascular clamps come in a variety of shapes and sizes.

REFLEX CLIP APPLIER

11 cm (4.3 in.) Long

500343	For 7 mm Clips	
500345	For 9 mm Clips	



REFLEX CLIP

- 100 per box
- Non-sterile
- Stainless Steel

500344	7 mm	
500346	9 mm	



500344

VESSEL CLIPS

- 5 per pack
- · Stainless Steel


15911	0.8 x 5 mm Jaws, 10 g Pressure
15913	1 x 6 mm Jaws, 15 g Pressure
14120	1.5 x 5 mm Jaws, 30 g Pressure
14121	1.5 x 8 mm Jaws, 60 g Pressure



Disposable Scalpels

- · Carbon steel blades
- · Non-slip handles
- · Used for incisions and underlying tissue
- Packed 10 per box
- Sterile
- Individually wrapped

500348	#10	
500349	#11	
500350	#12	
500351	#15	
500352	#20	
500353	#21	
500354	#22	
500355	#23	



Rodent Accessory

Features

- · Large, stable base
- · Hardened blades for long service
- · Ambidextrous configuration
- 3 sizes
- · Cuts through bone


Benefits

- Blades are drawn together by magnetic force to ensure a clean and precise cut through very strong bones and skin
- · Spring action, locking device prevents blades from accidentally falling
- Fluoropolymer coated base for easy cleanup

Applications

• Decapitator

DCAP	For Rodents/Other Small Animals, 1.5 x 1.5 in. Opening
DCAP-M	For Large Rodents/Other Medium Animals, 2.5 x 2.5 in. Opening
DCAP-L	For Larger Animals, 4.0 x 4.0 in. Opening




The small animal guillotine is easy to use with extra added safety features. There is a large base for stability, long handle for extra leverage, spring action so the blades can not fall down unexpectedly, hardened stainless blades for endurance, simplified construction for easy maintenance. The fluoropolymer coated surface on the base makes cleaning easy.

The guillotine is considered one of the most humane methods to sacrifice a subject.

Rodent Brain Matrices

WPI offers a large selection of acrylic and stainless steel brain matrices. These matrices are sturdy and can be heated, chilled, autoclaved, scrubbed and stand up to rigorous daily use. Coronal matrices have the

additional feature of a mid-line sagittal cut to facilitate splitting of the left and right hemispheres. Sections can be as fine as 1 mm. The olfactory/spinal/notch is cut into each matrix.





		SI	PECIFICATI	ONS						
Order #	Subject	Material	Section	Α	В	С	D	E	Cavity Depth	Weight
RBMA-200C	Adult Mouse, 40-75 g	Acrylic	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb.
RBMA-200S	Adult Mouse, 40-75 g	Acrylic	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	0.5 lb.
RBMA-300C	Rat, 175-300g	Acrylic	Coronal	4.76	15.9	12.7	36.6	23.8	7.61	0.5 lb.
RBMA-300S	Rat, 175-300g	Acrylic	Sagittal	4.76	15.9	12.7	36.6	23.8	7.61	0.5 lb.
RBMA-600C	Rat, 300g-600g	Acrylic	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb.
RBMA-600S	Rat, 300g-600g	Acrylic	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	0.5 lb.
RBMS-200C	Adult Mouse	Stainless Steel	Coronal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb.
RBMS-200S	Adult Mouse	Stainless Steel	Sagittal	3.18	11.1	8.73	19.1	12.2	7.4	1.0 lb.
RBMS-300C	Rat, 175-300g	Stainless Steel	Coronal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb.
RBMS-300S	Rat, 175-300g	Stainless Steel	Sagittal	4.76	15.9	12.7	36.6	23.8	7.61	1.0 lb.
RBMS-600C	Rat, 300g-600g	Stainless Steel	Coronal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb.
RBMS-600S	Rat, 300g-600g	Stainless Steel	Sagittal	4.76	19.8	14.7	36.6	24.7	10.91	1.0 lb.

Surgical Kits

PHYSIOLOGY KIT I

Kit Includes:

- SuperCut Tenotomy Scissors, Curved (14396)
- Dumont Forceps (500342)
- Rongeur, 3 mm jaw (14091)
- Filter Forceps, 11 cm, Straight, Flat Jaw (500456)
- Utility Scissors (501322)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Operating Scissors, Straight, Sharp/Blunt (14192)

KIT-PHYSIO-I

Physiology Kit I (7 instruments)



PHYSIOLOGY KIT II

Kit Includes:

- Vannas Scissors (14003)
- Stevenson Retractor (14131)
- SuperCut Iris Scissors, Straight (14218)
- Iris Forceps, Curved, Serrated (15915)
- Rongeur, 1.3 mm Jaw (14292)
- Adson Forceps, 1x2 Teeth (500092)
- Utility Scissors, Straight, Sharp/Blunt (501322)
- Olsen-Hegar Needle Holder (500227)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Dumont Forceps (500342)
- Operating Scissors, Straight, Sharp/Blunt (14192)
- Portfolio (501319)

PHYSIO-II Physiology Kit II (11 instruments plus zipper case)



MOUSE DISSECTING KIT


Kit Includes:

- Dumont Forceps (500342)
- Wire Retractor (14130)
- Vannas Scissors (14003)
- Needle Holder (14109)
- · Iris Forceps, Curved, Serrated (15915)
- Probe, 1.0 mm Diameter, Blunt (501313)
- Dissecting Scissors, 10 cm, Straight (14393)
- Portfolio (503294)

MOUSEKIT Mouse Kit (7 instruments plus zipper case) 504741 Mouse Kit for Telemetry Training (includes 500453)



Vessel Cannulation Forceps (500453) Included only with Kit 504741



500342

15915


14240

500342

RAT DISSECTING KIT

Kit Includes:

- Dumont Forceps 500342
- Alm Retractor 14240
- · Vannas Scissors 14124
- Needle Holder 14110
- Iris Forceps, Curved, Serrated 15915
- Probe, 1.0 mm Diameter, Blunt 501313
- SuperCut Iris Scissors, Straight 14218



15915

14110

14109

14393

Included only with Kit 504740

STORAGE PORTFOLIO

JIONAC	E I OKII OLIO	W I	
501319	Large Surgical Instrument Storage Portfolio, Holds up to 40 Instruments	*****	
503294	Small Surgical Instrument Storage Portfolio, Holds up to 10 Instruments		1
		503294	501319

OmniDrill35 Micro Drill

Drill, grind or finish bone or teeth Features

- Excellent tool for grinding, finishing, and drilling bone and other material
- High-torque motor (35,000 rpm) is quiet
- · Includes foot switch

Benefits

- · Removable nose cone that can be cleaned and sterilized
- · Consistant power for the duration of its use
- Variable speeds so you control the amount of heat generated

Applications

• Grinding, finishing, cutting and drilling bone, teeth and other material This line-powered micro drill will make easy work of grinding, finishing, cutting and drilling bone, teeth and other material. The high-torque 35,000 rpm (maximum) motor is quiet and has minimal vibration which reduces wear on the motor and provides greater comfort for you. It also features a forward and reverse switch, "E Type" handpiece and handpiece holder. The handpiece has a removable nose cone that can be cleaned and sterilized. It accepts 3/32" and 2.33mm bur shanks. Unlike battery-powered drills, this unit maintains consistent power for the duration of use. The wide range of speeds allows you to control the amount of heat generation.

α	$\mathbf{M} \mathbf{N} \mathbf{H}$	וומח	I DE CI	SECIFI	CATI	ONIC
w	VI NI	UKIL	L35 SI	449121	LAII	UIV5

110 V, 50/60 Hz (**503598**) 240 V, 50/60 Hz (**503599**)

OUTPUT 0-32 VDC FUSE 1 A

OPERATING SPEED 0-35,000 RPM

DIMENSIONS 178 x 114 x 89 mm (7 x 4.5 x 3.5in.)

WEIGHT 1.7 kg (3.75 lbs.)



	ORDERING INFORMATION
503598	OmniDrill35 Micro Drill System, 110 V
503599	OmniDrill35 Micro Drill System, 220 V
Stand	tips, cutoff disk, mandrels and foot switch are all included.
	Replacements available below.

OPTIONAL ACCESSORIES/REPLACEMENT PARTS

501851	Abrading Tip, Stone, pkg. of 5
501852	Accessory Stand
501853	Ball Mill, Carbide, #1, .031" Diameter, pkg. of 5
501854	Ball Mill, Carbide, #2, .039" Diameter, pkg. of 5
501855	Ball Mill, Carbide, #3, .047" Diameter, pkg. of 5
501856	Ball Mill, Carbide, #4, .055" Diameter, pkg. of 5
501857	Ball Mill, Carbide, #5, .063" Diameter, pkg. of 5
501858	Ball Mill, Carbide, #6, .071" Diameter, pkg. of 5
501842	Ball Mill, Carbide, #7, .083" Diameter, pkg. of 5
501860	Ball Mill, Carbide, #1/4, .019" Diameter, pkg. of 5
501861	Ball Mill, Carbide, #1/2, .027" Diameter, pkg. of 5
501862	Cutoff Disk, pkg. of 20
501863	Mandrel, Screw, pkg. of 5
501864	Mandrel, Threaded, pkg. of 5
504459	Foot switch
502237	Stereotaxic Holder for OmniDrill35 Microdrill

Economy Electrosurgical Unit

Cut and coagulate with a press of the foot switch

Features

- 10 levels of output intensity
- Three operational modes (cut, coagulate, cut/ coagulate)
- · Choice of electrodes

Benefits

• Ready to use out of the box.

• Includes a variety of electrodes, which are also sold individually

Applications

• Cutting and coagulation of tissue

Electrosurgery utilizes alternating current at radio frequencies to cut and coagulate. Using this method, the current enters the subject's body and the subject becomes part of the electrical circuit. This requires the use of a return, or ground, plate. This economically priced electrosurgical

unit has 10 levels of output intensity, three operational modes (cut, coagulate, and cut/coagulate) and various choices of electrodes. The unit comes complete and ready-to-use with a handpiece, ground plate, foot switch and one of each electrodes. All accessories can also be ordered separately.

ELECTROSURGICAL UNIT SPECIFICATIONS

OPERATION FREQUENCY 1.5 MHz
STABLE & FINE POWER SETTING 10 Steps

POWER SUPPLY $\begin{array}{c} 115 \text{ V} \pm 10\% \text{ - } 50/60 \text{ Hz } 1.8 \text{ A, } 210 \text{ VA} \\ 230 \text{ V} \pm 10\% \text{ - } 50/60 \text{ Hz } 0.9 \text{ A, } 210 \text{ VA} \\ \end{array}$

 $\begin{array}{ll} \text{OUTPUT POWER} & 70 \text{ W} \pm 5\% \\ \text{WORKING FRQUENCY} & 1.5\text{-}1.7 \text{ MHZ} \pm 5\% \\ \end{array}$

DIMENSIONS 24 cm x 22 cm x 8.5 cm (lxwxh)

SHIPPING WEIGHT 10 lb. (4.5 kg)

www.wpiinc.com

	ORDERING INFORMATION
501274	Electrosurgical Unit, 110 V
501284	Ball Electrode, Ø 1.6 mm shaft

501274

Dry Sterilizer

Heat sterilize instruments in seconds

Features

- · Large LED control panel
- · Adjustable temperature control up to 300°C
- · Sterilizes forceps, needles and other instruments
- Choose Bead capacity: 300 g (ST5193), 150 g (ST5191)

Benefits

- No chemicals
- · No flames
- No risk of burns
- · No disinfectant fluids

Applications

 Sterilize your microdissecting and tissue culture instruments, thoroughly and conveniently, in seconds

WPI's **Steri-Lite Micro Bead Sterilizer** is ideal for the sterilization of small research instruments (forceps, scissors, etc.). When the chamber is filled with the included glass beads, high temperatures (up to 300°C) can be used to eliminate bacteria, spores and other microorganisms.

Chamber temperature is displayed on the large LED control panel and may be adjusted with the control knob. The outer casing of the unit stays cool to the touch, even when it is used continually. The supplied safety cover can be placed over the chamber to conserve energy and reduce the heating/ramping rate. The Steri-Lite Sterilizer does not use gas, flame or chemicals and is safe for use in hoods.

With adjustable temperature, the sterilizer can be used for incubations and the decontamination of plastic objects/tools.

This product is designed for RESEARCH and DEVELOPMENT USE ONLY.



STER	I I IT	E CD	ECIEI	CATI	ONIC
SIEK		ESF	ECIFI	CAII	CINO

TEMPERATURE RANGE: 100 to 300°C

PREHEAT TIME: Approximately 25 min. to 300°C

CHAMBER DIAMETER: 40 mm

CHAMBER HEIGHT: ST5191: 80 mm

ST5193: 140 mm

DIMENSIONS: ST5191: 14x14.5x15.5 cm (5.5x5.7x6.1 in.)

ST5193: 14X14.5X21.5 cm (5.5X5.7X8.7 in.)

BEAD CAPACITY: ST5191: 150 g

ST5193: 300 g

WEIGHT: ST5191: 2 Kg


ST5193: 2.5 Kg

ELECTRICAL: ST5191: 115 V or 230 V, 50-60 Hz, 120 W

ST5193: 115 V or 230 V, 50-60 Hz, 250 W

WARRANTY: 2 Years





Request the surgical instrument catalog to see the full line of instruments:

www.wpiinc.com/catalog

Ultrasonic Cleaning Systems

Quantrex® bath with a timer

Features

- · 60 Minute Timer
- · Vinyl-clad steel and stainless steel
- 14 quality inspection steps for strength and durability

Benefits

- · Provides super-strength cleaning every time
- When used with L&R's specialty formulated solutions, the selfcontained, compact unit offers efficient trouble-free cleaning
- Each Quantrex machine comes standard with increased power strength you can see as soon as you turn the unit on
- Stainless steel drain with multi-positional outlet for easy removal of solution

Applications

· Versatile enough for a variety of cleaning applications

	SPECIFICATIONS
DIMENSIONS	26.0 x 16.5 x 21.0 cm (10.25 x 6.5 x 8.25 in.)
SHIPPING WEIGHT	4.5 kg (10 lb.)
TANK CAPACITY	3.2 L (0.85 gal.)
INTERNAL DIMENSIONS	23.8 x 13.7 x 10.2 cm (9.38 x 5.38 x 4.0 in.)



504216

	ORDERING INFORMATION
504216	Quantrex Ultrasonic Cleaning System
	Includes Timer, Drain and Cover. Heater Optional.

Economy ultrasonic cleaner



Features

- · Half liter stainless steel tank
- · Durable and compact

Benefits

· Robust all-metal construction allows for continuous duty


Applications

Microdissection and veterinary instrument cleaning

	·
U	BATH SPECIFICATIONS
INPUT POWER	22W
PEAK OUTPUT	70W, 55 kHz
TANK CAPACITY	0.53L (18 oz.)
TANK I.D.	12.1 x 8.6 x 6.6 cm, (4.75 x 3.375 x 2.625 in.)
TANK O.D.	13.7 x 10.5 x 12.1 cm, (5.375 x 4.125 x 4.75 in.)
SHIPPING WEIGHT	1.8 kg (4 lb.)



	ORDERING INFORMATION
UBATH-Y	Ultrasonic Cleaner (110 V)
13740	Ultrasonic Detergent (4 lb)
504766	Sterilization mesh casette 40x40x20 mm
504767	Sterilization mesh casette 80x80x34mm
504768	Sterilization mesh casette 105x70x25 mm



Request a complete catalog of sterilization baskets and trays:

www.wpiinc.com/catalog

Lab Supplies



Affordable laboratory equipment

Looking for a quality source for affordable laboratory equipment? From capillary glass and adhesives up to micromanipulators, pumps and microscopes, we offer the full range for the well stocked scientific laboratory.

Which Adhesive is Right for Me?

	А	DHESIVES APPLICA	ATION GUIDE
WPI Part #	Description	Curing time	Useful Applications and Characteristics
Epoxies	Form strong bonding. Used in	wire bonding application	s.
4898	Silver filled conductive Epoxy	12 hr @ 50°C; 5 min @ 150°C	Connecting conductors that can't be soldered. Constructing or connecting Ag/AgCl pellets.
7335	Carbon filled conductive Epoxy	48 hr @ 25°C; 5 min @ 150°C	Constructing carbon electrode.
4886	High performance Structural Epoxy	12 hr @ 25°C.	Forms a strong and slightly flexible bond on plastic, metal, & glass. Bonds some low surface.
Hot Melt (EVA)	Easy to use for bonding. For la	rge gap filling.	
13316	Mini Glue Gun with glue sticks	As soon as it cools down	Bonds wood, glass, metals, and many plastics.
Silicone Adhesive	es/Sealants/Primers	Good moisture resistan	nce and elasticity. Low toxicity.
1571	Room temperature vulcanizing (RTV) adhesive. Acyloxy/moisture cure system. Acetic acid is cure by-product.	24 hr @ 25°C	Has the best adhesion property in this silicone family. Will bond to many materials.
7128	RTV sealant. Alkoxy/Moisture cure system. Methanol as cure by-product.	72 hr @ 25°C	Good for bonding or sealing electronics circuits (metal).
SYLG184	Sylgard, two parts, vinyl/platinum cure sealant. Hydrogen as cure by-products. Very low toxic	24 hr @ 25°C, 15 min. @ 150°C	Coating Patch Clamp electrodes, Cell culture dish, making dissection pads.
KWIK-SIL	Two part, adhesive. Vinyl/platinum system, Hydrogen as cure by-products. Very low toxic.	< 5 min. @ 25°C	Live tissue and nerve studies. Medium strength adhesion.
KWIK-CAST	Two part sealant. Vinyl/platinum cure system. Hydrogen as cure by-products. Very low toxic.	< 5 min. @ 25°C	Sealant for live tissues. Embedding peripheral nerves with electrodes.
6820	Primer for silicone	N/A	Enhances adhesion of silicone adhesives for difficult to bond plastic surfaces
Cyanoacrylate	Forms an instantaneous bond.		
7341	Ethyl Cyanoacrylate, low viscosity 90-120 cps	<10 seconds	Mounting rat/mouse brain slices. Ideal for relatively small gaps and smooth surfaces.Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
7342	Ethyl Cyanoacrylate, high viscosity 1100-1600 cps	<30 seconds	Use on brain slice exp. Ideal for larger gaps, allows slightly longer bonding time. Bonds plastic, metals and rubber. Package of 10 vials, each approximately 1.5 mL.
VETBOND	Butyl Cyanoacrylate, low toxic	<10 seconds	Bonds tissues, alternative to suture, helps small wound healing. Antimicrobial effect. Used in forensic science.
503763	Octyl Cyanoacrylate, low toxic	<15 seconds	Suitable for surface wound bonding, protection, holding a sensor or other device on the tissue.

Easy Application Silicone Structural Epoxy

Kwik-Gard™ is specially packaged Sylgard 184 silicone for quicker and easier application, eliminating the messy procedure of preparing the mixture before application. Its special cartridge controls the precise mixing ratio to ensure proper curing. The disposable tip mixes resin



and hardener as they are dispensed. Since no air is introduced during mixing, the resin does not need degassing for most applications. The mixed silicone is applied directly to the site, reducing preparation time and material waste.

Each Kwik-Gard cartridge contains 37 mL of resin and hardener. The dispensing tip has a dead volume of 0.75 mL. Perfect for PDMS. Biopsy punches (0.5–8mm) are available online at www.wpiinc.com/punch.

	ORDERING INFORMATION
KWIKGARD	Kwik-Gard Start-up Kit (incl. dispenser, 1 cartridge, 5 tips)
KWIKGLUE	Kwik-Gard Refill (2 cartridges, 10 dispensing tips)
KWIKMIX	Dispensing Tips (pkg. of 10)
KWIKGUN	Kwik-Gard Dispenser

Scotch-Weld 2216
remains the best epoxy
for bonding plastic,
often used as the
benchmark for testing
the binding strength of
other adhesives. The slightly
rubbery texture also makes it less easy to break
off. It is the only epoxy known that can bond PEEK.

Color: gray
Cures at room temperature.
Shipping weight: 1 lb. (0.5 kg)

	ORDERING INFORMATION	
4886	Scotch-Weld 2216 (2 oz.)	

2216 B/A

Low Toxicity Adhesive

Ideal for neuroscience applications, nerve studies and more

Kwik-Sil"

Kwik-Sil"

KWIK-SIL

Features

- Bio-compatible adhesive for live tissue and nerve studies
- · Pre-mixing tips simplify use
- · Medium strength adhesion
- · Low toxicity
- Rapid curing silicone adhesive, cure on contact
- Cures without producing heat
- Includes 10 Mixing Tips
- · Volume discounts Save up to 15%!

Benefits

- Low toxicity
- Rapid cure time

Applications

- Neuroscience and nerve studies
- Biomedical applications



Kwik-Sil and **Kwik-Cast** silicones have very low toxicity before, during and after curing. The by-product of curing is a small amount of hydrogen gas, which is much less toxic to cells than acetic acid or alcohol from traditional RTV silicone sytems.

Kwik-Sil and **Kwik-Cast** curing speed is hundreds of times faster than traditional RTV silicones. A curing time of a few minutes at room temperature is especially useful for encapsulation of live tissue or implanting into a live animal.

Unlike many vinyl-based silicones in which the platinum complex catalysts are easily poisoned by contamination from amines and animal tissue, **Kwik-Sil** and **Kwik-Cast** are not sensitive to contamination from animal tissue

Kwik-Sil™ is a translucent, medium-viscosity silicone adhesive, developed for chronic peripheral nerve studies such as anterograde tracing with fluorescent indicators or electrode recording. Good adhesion and mechanical properties (tear strength and elongation) allow days of study without breaking of the bonding. Curing speed is very reproducible.

 $Kwik-Cast^m$ is a very low viscosity silicone sealant developed to embed peripheral nerves with electrodes for acute multi-fiber recordings. It

flows easily, filling the small spaces around the nerve and leaving no channels through which peritoneal fluid can travel and thus short the nerve/electrode contact. Equally important is the ability of the material to flow into itself and create one continuous mass from underneath the nerve all the way to the top of the nerve/electrode contact to ensure long-term recording stability. **Kwik-Cast** is color-coded to make the mixing foolproof. The catalyst is yellow and the base is blue. When uniformly mixed, it is green. **Kwik-Cast** can be applied and cured underneath mineral oil. After recording, electrodes are easily recovered due to the low tear strength.

KWIK-CAST

KWIK-CAST & KWIK-SIL	. SPECIFICAT	IONS
	Kwik-Sil	Kwik-Cast
MIX RATIO	1 to 1	1 to 1
WORKING TIME	< 5 minutes*	4 minutes
SETTING TIME (ROOM TEMP., 1:1 RATIO)	5–10 min- utes**	<10 minutes
CURE TIME	~15 minutes	
VISCOSITY, CPS	15,000	10,000
GUARANTEED SHELF LIFE AT 23 °C	6 months	6 months
VOLUME	5 mL	5 mL
NUMBER OF MIXING TIP	10	10
DEAD VOLUME OF THE MIXING TIP	<0.12 mL	<0.12 mL
AFTER CURING 24 HOURS:		
TEAR STRENGTH, PPI	90	44
ELONGATION %	650	60
DUROMETER (SHORE A-2)	30	36
COLOR	translucent	green
VOLUME RESISTIVITY, W/CM	1x10 ¹⁵	1x10 ¹⁵
+ 2 :	I CI: : I:	

^{* 3} minutes average with about 90 seconds of liquidity

^{**} no longer mixable at this point

	ORDERING INFORMATION
KWIK-SIL	Silicone Adhesive Compound (two 5 mL syringes)
KWIK-CAST	Silicone Casting Compound (two 5 mL syringes)
600022	Replacement KWIK Mixing Tips (pkg. of 10)
	Quantity discounts available

MicroFil™ Non-metallic Syringe Needle

Perfect for filling micropipettes

WPI's MicroFil™ fills micropipettes easily and reliably. Its long and fine tip allows you to start the filling very close to the pipette tip, eliminating both air bubble formation and clogging due to the washing down of dust particles. The transparent amber MicroFil™ needle is constructed from a combination of plastic and fused silica — no metal components are used. The MicroFil™ needle can be stored for days with the filling solution inside without clogging.

The MicroFil's tip elasticity is sturdy and very flexible though not unbreakable. Since it is more flexible than stainless steel needles, moderate bending will not block or damage the MicroFil™ needle. The combination of plastic and fused silica in the MicroFil™ tip is sturdier than plastic tips, allowing easy and repeated insertions into micropipettes. MicroFil's luer fitting allows easy coupling to syringes and syringe filters.



	ORDERING INFORMATION
MF34G-5	MicroFil™, 34 ga., 67 mm long (pkg. of 5)
MF28G-5	MicroFil™, 28 ga., 97 mm long (pkg. of 5)
MF28G67-5	MicroFil™, 28 ga., 67 mm long (pkg. of 5)

Custom MicroFil™

All MicroFil™ products, including custom orders, can be shipped immediately. Custom orders for special needs can be made using nine sizes of MicroFil™ tubing in lengths up to 50 cm — except for CMF90U which has a maximum length of 10 cm because of its high resistance to flow. Quantity discounts available. Specify length when ordering.

CME20G MicroFil™ 20 ga 700 μm ID, 850 μm OD (pkg. of 4)

CIVIF20G	Microfil''', 20 ga, 700 µm iD, 850 µm OD (pkg. 01 4)
CMF22G	MicroFil™, 22 ga, 530 μm ID, 700 μm OD (pkg. of 4)
CMF23G	MicroFil™, 23 ga, 530 μm ID, 665 μm OD (pkg. of 4)
CMF26G	MicroFil™, 26 ga, 320t μm ID, 430 μm OD (pkg. of 4)
CMF28G	MicroFil™, 28 ga, 250 μm ID, 350 μm OD (pkg. of 4)
CMF31G	MicroFil™, 31 ga, 100 μm ID, 238 μm OD (pkg. of 4)
CMF34G	MicroFil™, 34 ga, 100 μm ID, 164 μm OD (pkg. of 4)
CMF35G	MicroFil™, 35 ga, 75 μm ID, 144 μm OD (pkg. of 4)
CMF90U	MicroFil™, ~36 ga, 20 μm ID, 90 μm OD (pkg. of 4)

"Super" Adhesives for Life Science Research

Four times stronger than butyl cyanoacrylate and less toxic

Cyanoacrylate adhesives have been on the market since 1958. Most industrial or household grade cyanoacrylate is made of shorter alkyl chain derivatives such as methyl or ethyl cyanoacrylate (7341 and 7342). They are very useful for temporarily holding tissues such as mounting specimens for microtome sectioning. However, they are not suitable for bonding wounds on live animals. The difficulties of using cyanoacrylate for bonding live animals include:

- · A strong, irritating odor
- Quick loss of bonding strength due to breakdown of the bonding by hydration
- Breakdown products (cyanoacetate and formaldehyde) are toxic and can cause inflammatory reactions
- · Low flexibility and tend to be brittle.

To overcome these problems, several longer alkyl chain cyanoacrylates have been developed especially for veterinary and human use. The first longer alkyl chain product is butyl cyanoacrylate. This product has been used for animal and human applications outside the USA since 1970. It is much less toxic and has a lower odor than the methyl or ethyl cyanoacrylate. The butyl cyanoacrylate offered by WPI is **Vetbond™**.

A family of adhesives containing octyl cyanoacrylate, a plasticizer and stabilizer, was developed In the 1990's (one of them approved by FDA for human use). When bonding to tissue, these new adhesives are four times stronger and less toxic than butyl cyanoacrylate. Compared with the traditional suture, the new super adhesive has several advantages. On average, it takes only one-tenth of the time to close an incision.



VETBOND

The bonding strength is equal to 5-0 monofilament suture. It also has a mysterious antimicrobial effect that can decrease infection rates in contaminated wounds. Bonding will slough off naturally in 5 to 7 days. Cosmetic appearance of the healed incision is also better.

Gluture Topical Tissue Adhesive **503763** forms a strong and flexible film and is thus more suitable for surface wound bonding,



protection, and holding a sensor or other device on the tissue. Setting time is about 10 seconds, which gives ample time for application. It can also be used for temporarily holding a live tissue. For example, there is a report of using it to hold nematodes on a glass slide for patch-clamp neurons recording.

All of the products offered by WPI are veterinary grade (not suitable for human application). Though very similar to the grade for human use, they are not sterile and do not have FDA approval.

	ORDERING INFORMATION
503763	Gluture Topical Tissue Adhesive (10 tips), 1.5 mL
7341	Cyanoacrylate Adhesive, Low Viscosity—90-120 cps
	(package of 10 vials, each approximately 1.5 mL)
7342	Cyanoacrylate Adhesive, High Viscosity—1100-1600 cps
	(pkg. of 10 vials, each approximately 1.5 mL)
VETBOND	3M Vetbond™ Adhesive (3 mL)

Sylgard



A two-part silicone elastomer, ideal for potting and encapsulating applications. Very low dielectric constant sealing compound used in patch clamping and many other lab applications. After cure, will withstand -55° to 200 °C.

	ORDERIN	NG INFOR	MATION
--	---------	----------	--------

SYLG184 Sylgard, 0.5 kg (1.1 lb)

Silicone Dissecting Pad Kit



Make your own silicone dissecting pads easily and quickly. Mix the 2-part silicone right in the plastic petri dishes and allow to cure 24 hours at room temperature. Kit includes enough 2-Part Sylgard silicone elastomer to prepare 20 dishes, pins and 20 plastic petri dishes with lids, 65 mm.

ORDERING INFORMATION 501986 Silicone Dissecting Pad Kit

WORLD PRECISION INSTRUMENTS

www.wpiinc.com

Electrically Conductive Silver Epoxy



Two-component silver-filled epoxy for electrical connections which cannot be soldered, such as Ag/AgCl pellets. This widely used silver-filled epoxy features low viscosity and smooth flowing character. Pure silver is dispersed in both resin and hardener. Cures in 15 minutes at 120 °C. Mix ratio 1:1. May be premixed and frozen for later use.

ORDERING INFORMATION

4898 Silver Epoxy, 28 g (1 oz)

Electrically Conductive Carbon Epoxy



Two-component carbon-epoxy, curable at room and elevated temperatures. Ideal for electrostatic discharge protection and EMI/RFI shielding, 1:1 mix ratio. May be premixed and frozen for later use.

ORDERING INFORMATION

7335 Carbon Epoxy, 56 g (2 oz)

Silicone RTV Adhesive (non-acidic)



Because it is non-corrosive, this material is ideal for use on metal, for encapsulating small circuits on connectors. After cure, will withstand -55° to 200 °C. No mixing required.

ORDERING INFORMATION

7128 RTV Coating, 90 mL (3 fl oz)

Silicone RTV Adhesive



Clear silicone sealant provides good bonding to plastic. After cure, will withstand -55 to 200 °C. No mixing required. A handy, general purpose laboratory sealant. (Releases acetic acid during curing.)

ORDERING INFORMATION

1571 RTV Sealant, 139 mL (4.7 fl oz)



RTV Prime Coat

Enhances adhesion of silicone adhesives to many difficult-to-bond plastic surfaces.

ORDERING INFORMATION

6820 RTV Prime Coat, 400 mL (13.5 fl oz)



WORLD PRECISION INSTRUMENTS

Glass Capillaries

Quality glass, superior prices for microinjection/microelectrodes

Features

- Quality borosilicate glass capillaries
- Large variety available, including fire polished, filaments, thin wall, specialty glass and multi-barrel

Benefits


- · Superior pricing
- Most glass orders ship within 48 hours

Applications

- Microinjection
- · Electrophysiology
- · Patch clamp
- Fluid Handling

Fire Polishing

Fire-Polished glass capillaries are easier to insert into microelectrode holders without damaging the gasket. More importantly, fire-polished glass won't scratch the chloridized wire used in a recording electrode. Fire-polishing does not affect the glass's mechanical or electrical properties.



Making Uniform, Reproducable Microelectrodes


Borosilicate glass capillaries: Close dimensional tolerances assure microelectrode uniformity and reproducibility. Capillaries are available in 1, 2, 3, 5 and 7-barrel configurations, complete range of single barrel thin-wall sizes and a variety of special configurations. Capillaries with filaments contain a solid filament fused to the inner wall, which speeds filling of electrodes. Capillaries with or without inner filaments are available for making microelectrodes in a wide range of diameters.

Filament Glass Capillaries


Single Barrel standard wall thickness capillaries are offered either with or without inner filaments for quick filling in a variety of lengths and diameters.

Thin Wall Glass Capillaries

Thin Wall single barrel capillaries are offered both with or without inner filaments.



Specialty glass is also available. See page 105.



ORDERING INFORMATION							
	Length	OD(mm)	ID(mm)	Filament	Fire-Polished	Quantity	Item
	3 in. (76 mm)	1.0	0.58	~		500	1B100F-3
	3 in. (76 mm)	1.0	0.58			500	1B100-3
	3 in. (76 mm)	1.2	0.68	V		350	1B120F-3
SS	3 in. (76 mm)	1.2	0.68			350	1B120-3
<u> </u>	3 in. (76 mm)	1.5	0.84	V		225	1B150F-3
ë	3 in. (76 mm)	1.5	0.84		V	300	1B150-3
cat	4 in. (100 mm)	1.0	0.58	V	V	500	1B100F-4
=	4 in. (100 mm)	1.0	0.58		✓	500	1B100-4
õ	4 in. (100 mm)	1.2	0.68	v	V	400	1B120F-4
80	4 in. (100 mm)	1.2	0.68			350	1B120-4
5	4 in. (100 mm)	1.5	0.84	V	V	300	1B150F-4
g	4 in. (100 mm)	1.5	0.84		~	300	1B150-4
an	4 in. (100 mm)	2.0	1.12	V		125	1B200F-4
S	4 in. (100 mm)	2.0	1.12		~	200	1B200-4
<u>e</u>	6 in. (152 mm)	1.0	0.58	V		500	1B100F-6
3ar	6 in. (152 mm)	1.0	0.58			500	1B100-6
e E	6 in. (152 mm)	1.2	0.68	V		350	1B120F-6
뗼	6 in. (152 mm)	1.2	0.68			350	1B120-6
Single Barrel Standard Borosilicate Glass	6 in. (152 mm)	1.5	0.84	v		225	1B150F-6
	6 in. (152 mm)	1.5	0.84			225	1B150-6
	6 in. (152 mm)	2.0	1.12	V		125	1B200F-6
	6 in. (152 mm)	2.0	1.12			125	1B200-6
	3 in. (76 mm)	1.0	0.75	V		500	TW100F-3
	3 in. (76 mm)	1.0	0.75			500	TW100-3
_	3 in. (76 mm)	1.2	0.90	V	V	400	TW120F-3
arc	3 in. (76 mm)	1.2	0.90			350	TW120-3
ğ	3 in. (76 mm)	1.5	1.12	V		225	TW150F-3
ţ	3 in. (76 mm)	1.5	1.12		V	300	TW150-3
S	4 in. (100 mm)	1.0	0.75	V		500	TW100F-4
ırr	4 in. (100 mm)	1.0	0.75		V	500	TW100-4
-Ba	4 in. (100 mm)	1.2	0.90	~		350	TW120F-4
ė O	4 in. (100 mm)	1.2	0.90			350	TW120-4
ing	4 in. (100 mm)	1.5	1.12	~		225	TW150F-4
S	4 in. (100 mm)	1.5	1.12		V	300	TW150-4
Thin-Wall Single-Barrel Standard	6 in. (152 mm)	1.0	0.75	V		500	TW100F-6
<u>-</u>	6 in. (152 mm)	1.0	0.75		V	500	TW100-6
2	6 in. (152 mm)	1.2	0.90	~	V	400	TW120F-6
Τ.	6 in. (152 mm)	1.2	0.90			350	TW120-6
	6 in. (152 mm)	1.5	1.12	~		225	TW150F-6
	6 in. (152 mm)	1.5	1.12		V	300	TW150-6

Single barrel glass is Kimble N51A. All thin wall glass is Schott Duran 8330. Tolerance ±10%

Borosilicate Glass Micropipettes

Eliminate the cost and trouble of making your own micropipettes



Features

- · Schott Duran borosilicate glass
- 0.5 µm and smaller ID micropipettes include internal glass fiber for easy filling
- Tip inner diameter tolerance ±20%
- · Short taper yields high strength
- Nominal length ≈ 50 mm
- OD:ID = 1.33:1
- Standard capillary outer diameters are 1.0 mm (thin-wall) or 1.14 mm
- Every pipette individually tested and inspected
- Vacuum sealed packs of 10

Benefits

• Plain Shank or Luer Fittings

Applications

 Injection of dyes or proteins into cells, oocytes or other biomedical laboratory applications

WPI can quickly supply your need for consistently sized pre-pulled glass micropipettes.

Tip diameters (ID) range from 0.1 to 10 µm.

Silanized Tips (Luer Shank) are available. Silanization waterproofs the glass to retard water when inserting into cell. This will not let the outside fluid run down the pipette and get inside so easily.

ORDERING INFORMATION							
Shank	Tip I.D.	Silanize Coating Length	Glass O.D.	Filament	Fire Polished	Catalog #	
PLAIN	0.1 µm	_	1.0 mm Thin-Wall	Yes	No	TIP01TW1F	
	0.2 µm	_	1.0 mm Thin-Wall	Yes	No	TIP02TW1F	
	0.3 µm	_	1.0 mm Thin-Wall	Yes	No	TIP03TW1F	
	0.4 µm	_	1.0 mm Thin-Wall	Yes	No	TIP04TW1F	
	0.5 µm	_	1.0 mm Thin-Wall	Yes	No	TIP05TW1F	
	1 µm	_	1.0 mm Thin-Wall	No	Yes	TIP1TW1	
	2 µm	_	1.0 mm Thin-Wall	No	Yes	TIP2TW1	
	5 µm	_	1.0 mm Thin-Wall	No	Yes	TIP5TW1	
	10 µm	_	1.0 mm Thin-Wall	No	Yes	TIP10TW1	
	10 µm	_	1.14 mm A203XV glass *	No	Yes	TIP10XV119	
	_	_	1.0 mm Thin-Wall	No	Yes	TIP10LT	Long Taper †
	_	_	1.0 mm Thin-Wall	Yes	Yes	TIP10FLT	Long Taper †
	_	_	1.5 mm Thin-Wall	No	Yes	TIP15LT	Long Taper †
	_	_	1.5 mm Thin-Wall	Yes	Yes	TIP15FLT	Long Taper †
	30 µm		1.0 mm Thin-Wall	No	Yes	TIP30TW1	
LUER	0.1 µm	_	1.0 mm Thin-Wall	Yes	_	TIP01TW1F-L	
	0.2 µm	_	1.0 mm Thin-Wall	Yes	_	TIP02TW1F-L	
	0.3 µm	_	1.0 mm Thin-Wall	Yes	_	TIP03TW1F-L	
	0.5 µm	_	1.0 mm Thin-Wall	Yes	_	TIP05TW1F-L	
	1 µm	_	1.0 mm Thin-Wall	No	_	TIP1TW1-L	
	2 µm	_	1.0 mm Thin-Wall	No	_	TIP2TW1-L	
	5 µm	_	1.0 mm Thin-Wall	No	_	TIP5TW1-L	
	10 µm	_	1.0 mm Thin-Wall	No	_	TIP10TW1-L	
	30 µm		1.0 mm Thin-Wall	No		TIP30TW1-L	
LUER/	5 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS01	
SILANIZED	5 μm	2 inch	1.0 mm Thin-Wall	No	_	TIP5TW1LS02	
	10 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP10TW1LS01	
	10 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP10TW1LS02	
	30 µm	1 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS01	
	30 µm	2 inch	1.0 mm Thin-Wall	No	_	TIP30TW1LS02	

^{* 10} μ m (ID), 1.14 mm capillary pipettes are for use in WPI's **Nanoliter 2010**.

μTip Sampler Assortments

TIPMIX01-05	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Plain Shank
TIPMIX05-10	Two each, 0.5, 1, 2, 5, 10 µm ID, Plain Shank
TIPMIX01-05-L	Two each, 0.1, 0.2, 0.3, 0.4, 0.5 µm ID, Luer
TIPMIX05-10-L	Two each, 0.5, 1, 2, 5, 10 µm ID, Luer

Slides

These clean glass microscope slides are 25 x 75 mm, 1.0~1.2 mm thick with 90° ground edges. They are available as frosted and red ended. The frosted end slides feature a fine 20 mm frosted area on both sides of one end for easy



marking. The red frosted slides feature a 20 mm colored end useful for identifying hazardous materials.

	ORDERING INFORMATION
503506	Frosted Glass Microscope Slides, Box of 144
503507	Red Frosted Glass Microscope Slides, Box of 144

Cover Slips

These cover slips (made of German glass) can be used for growing and culturing cells that normally have poor adhesion to plastic surfaces. They are small enough to be placed in the micro plate or other cell culture devices. The 5 mm size will fit inside the 96-well culture plate and leave enough room to pick it up from the bottom of the well with forceps. The 8mm size fits inside the 24-well plates.



ORDERING INFORMATION									
Order#	Diameter	Thickness	Quantity						
502040	5 mm	#1.5 (0.16 - 0.19 mm)	100						
502041	8 mm	#1.5 (0.16 - 0.19 mm)	100						
503508	25 mm	#1.5 (0.16 - 0.19 mm)	100						

[†] Long Taper micropipettes are pulled with a 12-15 mm taper which the customer cuts back to obtain the desired tip diameter.

Luer Valve Assortment Kit

Build your own liquid flow experiment

Features

- Over 300 assorted parts
- Luer fittings for quick and easy connect and disconnect

Benefits

• Sold individually or in kits

Applications

Liquid flow experimental setups

A useful kit (right) for building your own liquid flow experiment. It provides the means to start, stop, add, divide and control a flow of liquid or gas. Included in the kit are *over 200 assorted parts* such as one-way and threeway stopcocks, manifolds, Y-connectors, injection sites, male and female luer caps,

check valves, syringe-activated check valves, slide clamps, roller clamps, and pinch clamps. All (except clamps) have a luer fitting for quick and easy connecting and disconnecting. Includes assorted luer fittings for use with flexible tubing.



ORDERING INFORMATION

14011 Luer Valve Assortment Kit

Luer-to-Tubing Coupler Assortment Kit

Quick connects in nylon and polypropylene

Features

- Over 250 assorted parts in each kit
- Valves are polycarbonate, and the valve handles are polyethylene. Do not autoclave those parts.

Benefits

- Polypropylene parts (504954) can be autoclaved repeatedly at 121°C/15PSI, 15 min. cycle
- Polypropylene fittings are chemically inert and resistant to most organic and inorganic solvents
- Nylon fittings are strong and can be bonded with adhesive.

Applications

• Liquid flow experimental setups



504954



504955

Assemble quick-disconnect luer fittings for use with flexible tubing with internal diameters of 1/16", 3/32" and 1/8". A variety of quick-disconnect connectors can be quickly made for connecting small diameter flexible tubing; 3-way connections can be made with the use of the 3-way luer tee; luer plugs, tees, connectors, bulk-head mounts, color coding rings, locking nuts, male and female luers—are all included to enhance the versatility of this kit. The kit has over 250 assorted parts and is offered in two different types of materials. Nylon parts are not autoclavable.

ORDERING INFORMATION

504954 Luer-to-Tubing Coupler Assortment Kit (Polypropylene)504955 Luer-to-Tubing Coupler Assortment Kit (Nylon)

Luer Valve Assortment Kit 14011 Parts

Kit parts are also available individually



14034-40 Injection Site Male luer lock Pack of 40



14039-10 Check Valve Pack of 10



14044-5 Syringe Activated Dual Check Valve Pack of 5



14045-20 Syringe Slip Luer Valve Activated Check Pack of 20



13822-10 0.135"/3.4 mm **OD Tubing** Pack of 10



14041-60 Roller Clamp 3/16" Tubing Pack of 60



7465-20 Pinch Clamp Large Bore Pack of 20



Pinch Clamp for 7mm Tubing Pack of 50



Female T Luer Pack of 20



14047-10 4-Port Infusion Y Swivel Thread Pack of 10



14048-20 3-Port Infusion Y Swivel Thread Pack of 20



4-Way Stopcock, Luer Lock Pack of 10



14036-15 4-Way Luer Stopcock Pack of 15



14058-10 4-Way Stopcock, Luer Lock Pack of 10



14035-10 4-Way Stopcock, Luer lock, Pack of 10



14051-100 Pinch Clamp for 5 mm Tubing Pack of 100



14038-10 1-Way Stopcock Luer Lock, Pack of 10



14054-10 1-Way Stopcock, Luer Slip Pack of 10



14055-2 4-Port Manifold (6 Female Ports) Pack of 2



3-Port Manifold (5 Female Ports) Pack of 2



13156-100 Female Luer Fitting for 1/16" ID Tubing Pack of 100



13157-100 Female Luer Fitting for 3/32" ID Tubing Pack of 100



13158-100 Female Luer Fitting Female Luer 1/8" ID Tubing Pack of 100



13159-100 Fitting for 5/₃₂" ID Tubing Pack of 100



13160-100 Male Luer Fitting for 1/16" ID Tubing Pack of 100



13161-100 Male Luer Fitting for 3/₃₂" ID Tubing Pack of 100



13162-100 Male Luer Fitting for 1/8" ID Tubing Pack of 100



13163-100 Male Luer Fitting for 5/₃₂" ID Tubing Pack of 100



14061-60 Male/Female Luer Plug Pack of 60



14042-100 Slide Clamp for 2.5 mm O.D. Tubing Pack of 100

Parts in kit may differ slightly in appearance from those pictured.

BNC Cables & Connectors

For wiring any laboratory setup



	ORDERING INFORMATION							
PART #	APPLICATION/DESCRIPTION	CONNECTOR A	CONNECTOR B	CABLE LENGTH				
2851	Standard BNC cable	BNC (male)	BNC (male)	6 ft (1.8 m)				
5374	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold pin	4 ft (1.2 m)				
5375	Low-noise cable for microelectrode holders	BNC (male)	2 mm gold jack	4 ft (1.2 m)				
13324	Adapter	Double-banana (female)	BNC (male)	none				
13347	ISO2 (chart recorder adapter)	Double-banana (male)	BNC (female)	none				
13451	Adapter: Iso-DAM, Iso-DAM8	BNC (female)	two 2 mm pins	6 in. (15 cm)				
13854	BNC T-connector, male to:	BNC (female)	BNC (female)	none				
14254	BNC Straight Adapter	BNC (female)	BNC (female)	none				
500184	Standard BNC Cable	BNC (male)	BNC (male)	10 ft (3 m)				
500256	BNC Right Angle Adapter	BNC (male)	BNC (female)	none				
500257	Standard BNC Cable	BNC (male)	BNC (male)	6 in. (15 cm)				
500258	Standard BNC Cable	BNC (male)	BNC (male)	12 in. (30 cm)				
500259	Standard BNC Cable	BNC (male)	BNC (male)	18 in. (46 cm)				
CBL102	DAM Series, PM Series	3.5 mmMiniPhone plug	BNC (male)	6 ft (1.8 m)				

Micro Cannula

Quickly connect to existing experimental plumbing

Features

- 0.4 mm O.D., 0.2 mm I.D. tubing
- Autoclavable
- Biocompatible Perfluorocarbon tubing material

Benefits

 May be used with a pressure tranducer (BLPR2) or a micro syringe injection system (UMP3 or MMP)

Applications

 Cannula for carotid or femoral arteries of rodents and small animal blood vessels

This micro cannula is ideal for placement in the carotid or femoral artery of mice, rats, and other small animal blood vessels. It can be used with a pressure transducer (WPI'S **BLPR2**) for blood pressure measurement, or in conjunction with a micro-syringe injection system (like WPI's **UMPIII** or

MMP pumps). The incorporated standard female luer fitting makes connecting to existing experimental plumbing quick and easy. The cannula is provided with a contoured-tip stainless steel stylet (trocar) to facilitate placement using established techniques. A movable "shoulder" ring provides a tie-in point to prevent accidental removal. The cannula may be left in place for 2 hours or more, and with proper care and cleaning, may be re-used multiple times. Instructions for use included.

ORDERING INFORMATION

KZ1101 Micro Cannula, 3-inch

www.wpiinc.com

Precious Metal and Specialty Wire

Bare and coated metal wire for most laboratory applications



Micro coaxial cables (MAXxxxx) are ideal for microelectrode fabrication and construction of similar research tools. The dual shielding eliminates electrical interference caused by radio frequencies (RF), electrostatic and microphonics (e.g., bending and vibration. Available with single or dual (twin) conductors.

Teflon-coated stainless steel (type 304) wire (SSTxxxx) is available in 25-ft and 50-ft lengths. The Teflon coating is 150 micro-in. thick (4 µm). The Teflon coating is designed to reduce surface friction, only. It is not insulation.

Carbon wire (C3005) is a single 30-micron fiber of electrochemically activated carbon. This fiber is especially useful in micro-electrochemical experiments.

Platinum/iridium wire — uncoated (PTxxxx) and Teflon-coated (PTTxxxx) - is an alloy of 90% platinum and 10% iridium, giving excellent tensile strength and corrosion resistance. Uncoated pure platinum wire (PTPxxx) is 99.95% pure. Indium wire (IN1003) is 99.99% pure, with a melting point of 156.4°C.

Annealed silver wire (AGWxxxx), 99.99% pure, is available in five diameters; three of those sizes are also available with a Teflon coating (AGTxxxx).

Tungsten wire (TGWxxxx), available in three diameters, is 99.95% pure. Gold wire (AUWxxxx) is 99.99% pure. Stainless steel wire (SSxxxxx) is type

Wire Cutters and Scissors



ORDERING INFORMATION							
Catalog No.	Metal	Coating	AWG*	Diameter	Precut Length		
AGT0510	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	10 ft (3 m)		
AGT0525	Silver	Teflon	36	0.005 in. (0.125 mm) ¹	25 ft (7.6 m)		
AGT05100	Silver	Teflon	36		100 ft (30 m)		
AGT1010	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	10 ft (3 m)		
AGT1025	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	25 ft (7.6 m)		
AGT10100	Silver	Teflon	30	0.010 in. (0.25 mm) ¹	100 ft (30 m)		
AGT1510	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	10 ft (3 m)		
AGT1530	Silver	Teflon	26-27	0.015 in. (0.38 mm) ¹	30 ft (9.1 m)		
AGW0510	Silver	_	36	0.005 in. (0.125 mm)	10 ft (3 m)		
AGW0530	Silver	_	36	0.005 in. (0.125 mm)	30 ft (9.1 m)		
AGW1010	Silver	_	30	0.010 in. (0.25 mm)	10 ft (3 m)		
AGW1030	Silver	_	30	0.010 in. (0.25 mm)	30 ft (9.1 m)		
AGW1510	Silver	_	26-27	0.015 in. (0.38 mm)	10 ft (3 m)		
AGW1530	Silver	_	26-27	0.015 in. (0.38 mm)	30 ft (9.1 m)		
AGW2010	Silver	_	24	0.020 in. (0.5 mm)	10 ft (3 m)		
AGW2030	Silver	_	24	0.020 in. (0.5 mm)	30 ft (9.1 m)		
AGW4010	Silver	_	18	0.040 in. (1.0 mm)	10 ft (3 m)		
AUW0170	Gold	_	50	0.001 in. (0.025 mm)	70 ft (21 m)		
AUW201	Gold	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)		
C3005	Carbon	_	49	0.0012 in. (30 μm)	5 ft (1.5 m)		
PT1002	Platinum / Iridium	_	30	0.010 in. (0.25 mm)	2 ft (61 cm)		
PT0402	Platinum / Iridium	_	38	0.004 in. (0.102 mm)	2 ft (61 cm)		
PT0203	Platinum / Iridium	_	44	0.002 in. (0.051 mm)	3 ft (91 cm)		
PT0110	Platinum / Iridium	_	50	0.001 in. (0.025 mm)	10 ft (3 m)		
PTP101	Platinum	_	30	0.010 in. (0.25 mm)	1 ft (30 cm)		
PTP201	Platinum	_	24	0.020 in. (0.5 mm)	1 ft (30 cm)		
PTP401	Platinum	_	18	0.039 in. (1.0 mm)	1 ft (30 cm)		
PTP406	Platinum	_	18	0.039 in. (1.0 mm)	0.5 ft (15.2 cm)		
PTT0502	Platinum / Iridium		36	0.005 in. (0.125 mm) ¹	2 ft (61 cm)		
PTT0203	Platinum / Iridium		44	0.002 in. (0.051 mm) ¹	3 ft (91 cm)		
PTT0110	Platinum / Iridium	letion	50	0.001 in. (0.025 mm) ¹	10 ft (3 m)		
SS31605	Stainless Steel	_	36	0.005 in. (0.125 mm)	50 ft (15.2 m)		
SS31614	Stainless Steel	— T-fl	27	0.014 in. (0.36 mm)	30 ft (9.1 m)		
SST30407-25	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	25 ft (7.6 m)		
SST30407-50	Stainless Steel	Teflon	33	0.007 in. (0.18 mm) ³	50 ft (15.2 m)		
TGW0325	Tungsten	_	40	0.003 in. (0.075 mm)	25 ft (7.6 m)		
TGW0515 TGW1510	Tungsten	_	36	0.005 in. (0.125 mm)	15 ft (4.6 m)		
19441210	Tungsten	_ 1ICROCO <i>l</i>	26-27	0.015 in. (0.38 mm)	10 ft (3 m)		
MAX3820	Tinned Cu Alloy	Coaxial	TAIAL C	0.0173 in. (0.44 mm)	20 ft (6 m) ⁴		
MAX4020	Tinned Cu Alloy	Twin Coaxial	0.0158×0	0.024 in. (0.4x0.61 mm)			
011020	*Pro	wn 9 Charno	0.013000	1 111. (0. 170.01 111111)	20 10 (0 111)		

*Brown & Sharpe

¹ Plus 0.002 in. for Teflon coating

³ Teflon adds 0.00015 in. (4 µm) to diameter

4 Impedance: 50 ohm; capacitance: 95 pF/m; resistance: 5 ohm/m 5 Impedance: 100 ohm; capacitance: 54 pF/m; resistance: 1.9 ohm/m

ORDERING INFORMATION

504749 Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) Micro-shear flush cutters for delicate wires 504750 Ergonomic Mini-Scissors, 12.7 cm (5 in.) For cutting fine or delicate items with a clean, square edge. Handle design is advantageous for users with arthritic hands.

Ergonomic Micro-Shear Flush Cutters, 12.7 cm (5 in.) For delicate wires. ESD safe. Extra tough high carbon steel blades. Sized for smaller 504751 hands and maximum maneuverablity


FrameWorks

Non-Magnetic Bases, Stainless Steel Rods, & Clamps

These high quality components are made of stainless steel and polymer that resist organic solvents and corrosion. They can be easily assembled to make a stand-alone setup for student labs or to make a complicated frame for research labs.











503086

14073-4

503082-4





503078-4





503079-4

503080-4

	ORDERING INFORMATION			
503041	Large Clamp with Rod (157 mm), opens up to 85 mm			
503042	Medium Clamp with Rod (157 mm), opens up to 45 mm			
503086	Small Clamp with Rod (157 mm), opens up to 16 mm			
14073-4	Open-sided Frame Clamp (pkg. of 4)			
503082-4	Board Frame Clamp, opens to 8.5 mm (pkg. of 4)			
503080-4	Frame Clamp with Parallel Surface Mount (includes mounting screws) (pkg. of 4)			
502193-4	Parallel Frame Clamp (pkg. of 4)			
503078-4	T-joint Frame Clamp (pkg. of 4)			
503079-4	In-line Frame Clamp (pkg. of 4)			
502190	Heavy Rectangular Base (with M8 Thread mount and thumbscrew mount), 23×15.6 cm, 4 lb.			
503083	Light Rectangular Base (with M8 Thread mount and thumbscrew mount), 23×15.6 cm, 0.5 lb.			
503085	Large 10-in. V-base with M8 Thread Mount			
503084	Small 6-in. V-base with M8 Thread Mount			
503081-4	Vertical Surface Mount, M8 Threaded (pkg. of 4)			
503070	Polished Stainless Steel Post, 12 mm OD, 25 cm long, No Thread			
503071	Polished Stainless Steel Post, 12 mm OD, 50 cm long, No Thread			
503072	Polished Stainless Steel Post, 12 mm OD, 75 cm long, No Thread			
503073	Polished Stainless Steel Post, 12 mm OD, 25 cm long, M8 Thread			
502191	Polished Stainless Steel Post, 12 mm OD, 50 cm long, M8 Thread			
503075	Polished Stainless Steel Post, 12 mm OD, 60 cm long, M8 Thread			
503076	Polished Stainless Steel Post, 12 mm OD, 75 cm long, M8 Thread			
503077	Polished Stainless Steel Post, 12 mm OD, 80 cm long, M8 Thread			



The high quality electronic digital caliper is a useful tool — no laboratory should be without one because it is more accurate and easier to use than the traditional analog devices. Measure in either inches or millimeters at the touch of a button. The floating zero feature allows

	ORDERING INFORMATION	
501601	Digital Caliper	
502157	Replacement Battery (package of 10)	

Non-rotating Spindle Digital Micrometer Head

Build your own precision micro-positioning device



The new non-rotating spindle digital micrometer head allows you to create your own micro-positioning instrument. With micron-level accuracy, it gives higher precision than a normal micromanipulator. Since the spindle does not rotate as it advances, instruments can be directly attached without the need for a complicated decoupling device. The digital display eliminates the need to squint at the notational scale. Readings can be clearly seen in either inches or millimeters. You can read both absolute position and the increment relative to a previously chosen point.

502102 SPECIFICATIONS

 TOTAL TRAVEL DISTANCE
 25 mm

 RESOLUTION
 0.001 mm

 ACCURACY
 ± 0.003 mm

 SPINDLE
 Ø 8 mm

 MOUNTING
 Ø 12 mm x 10 mm

TOTAL LENGTH 166 mm

MEASUREMENT MODE Absolute and incremental DIGITAL READOUT mm or inch

ANALOG READOUT mm
DATA OUTPUT RS232
ENVIRONMENTAL PROTECTION IP54

SHIPPING WEIGHT 0.51 kg (1.12 lb.)

ORDERING INFORMATION

502102 Non-Rotating Spindle Micrometer Head

Powerful Ball Joint Rare Earth Magnet

Construct holding devices for small parts/equipment

- Small but very powerful: holds 2 kilograms (~5 pounds)!
- Steel ball rotates freely 360° on a 180° axis
- M3 mounting screw on ball for attachment to equipment
- Magnet base threaded (M3) for mounting onto a base or equipment

This novel magnetic ball joint has phenomenal holding power for up to 2kg of attached weight while permitting the ball a full 360° rotation on a 180° axis. You can freely orient your equipment to an infinite number of positions within this rotation. This is made possible by the combination of a steel ball (10 mm diameter) and a powerful rare earth magnet contained in the magnet cylinder (ϕ 10 x

20 mm). Convenient M3 attachment sites are provided on both the ball (male) and the magnet base (female). For use with micromanipulators for the positioning and holding of optical instruments including various lighting sources and lasers, pipettes and any small parts that would benefit from the flexibility offered by this new magnetic ball joint.



ORDERING INFORMATION

500871 Magnetic Ball Joint

Back by popular demand: Black Wall

FluoroDish

Culture Dishes with Optical Glass Bottom,

Cover-glass bottom for observing & growing cells for imaging

Features

- Optical quality glass bottom for better imaging quality (RI=1.525)
- · Low sample volume for expensive chemicals
- Lowest access angle for micropipette
- Black wall available for low background fluorescent measurement

Benefits

- Multiple sizes and designs to suit your application
- · Optional Poly-D-Lysine coating for neurons
- · Dishes designed for low volumes or large growth areas

Applications

- · High resolution image analysis
- Microinjection
- Electrophysical recording of fluorescent-tagged cells
- Black wall available for use with confocal microscopes

WPI's FluoroDish™ tissue culture dishes provide exceptional imaging quality for many applications requiring the use of inverted microscopes such as high resolution image analysis, microinjection and electrophysical recording of fluorescent-tagged cells. Taking advantage of WPI's extensive experience with low toxicity adhesives, FluoroDish™ uses a specially formulated adhesive that is optically clear, durable and with extremely low toxicity. Tests by an independent laboratory have shown that the 96-hour surviving rate of embryos is 100% when kept in FluoroDish™, substantially better than other brands. The bottom glass has superior UV transmission (30% transmission at 300 nm, compared to less than 7% for the most popular German glass). Stringent quality control ensures that glass thickness is 0.17 ±0.01 mm.

Conventional plastic dishes and chambers limit the use of the inverted scope, because the thick plastic bottom requires a long working distance objective available only in lower magnifications. Each WPI dish has a flat (0.17mm thick), optical quality glass bottom, allowing the use of a much shorter working distance, larger numerical aperture (NA) and higher magnification (up to 100X). The larger NA and higher magnification provide superior quality imaging for both classical and fluorescence microscopy. Higher effective NA yields brighter images for fluorescence and higher resolution in image analysis. The glass bottom permits the use of immersion objectives with media such as water, glycerin or oil for the highest magnification possible. To optimize heat-exchange, WPI's glass-bottom dish is designed to be flush (flat) with the microscope stage or heating unit, eliminating the air gap that exists with modified plastic dishes where a cover slip was inserted.

Multiple sizes

We have a 50 mm diameter dish and two types of 35 mm diameter dishes. An inner well is created within the dish by the glass bottom and the tissue culture grade polystyrene which forms the sides of the dish. They are individually packed and gamma sterilized.

The 35 mm dish has outside dimensions similar to that of a Corning 35 mm dish and has $\emptyset 23.5$ mm glass window (**FD35**) or $\emptyset 10$ mm glass window (**FD3510**). Most heaters and perfusion adapters designed for the Corning 35 mm dish also fit this dish.

Poly-D-Lysine coating or uncoated

The 23.5 mm glass dish is available uncoated or poly-D-lysine-coated. Some cell lines (e.g., PC3 and HEK) adhere well to the uncoated glass bottom dish. The poly-D-lysine coating reportedly improves the adhesion of neuron cells. You may also apply any special coating that is best for your cell line to uncoated dishes.

Low volume dishes or large growth areas

The 10 mm glass dish (**FD3510**) has low sidewalls for easy microelectrode access and low solution volume. The low microelectrode access angle is the lowest among all of 35 mm dishes on the market (very close to that of a 50 mm dish). The dish needs only about 115 μL to cover the bottom well, an important feature when using expensive drugs and chemicals.

The 50 mm dish (**FD50**) has a large growth area (35 mm well diameter), a low access angle for microelectrodes, and grips for easy handling.




Petros Batsios, Xiang Ren, Otto Baumann, Denis A. Larochelle, Ralph Gräf, Cells. 2016 Mar; 5(1): 13. Published online 2016 Mar 18. doi: 10.3390/cells5010013 PMCID: PMC4810098

P. Calcraft, A. Arredouani, M. Ruas, Z. Pan, X. Cheng, X. Hao, J. Tang, K. Rietdorf, L. Teboul, K. Chuang, P. Lin, R. Xiao, C. Wang, Y. Zhu, Y. Lin, C. Wyatt, J. Parrington, J. Ma, A. M. Evans, A. Galione, M. Zhu. Nature. 2009 May 28; 459(7246): 596–600. Published online 2009 April 22. doi: 10.1038/nature08030 PMCID: PMC2761823

T Mitra-Ganguli, I. Vitko, E. Perez-Reyes, A. Rittenhouse. J Gen Physiol. 2009 November; 134(5): 385–396. doi: 10.1085/jgp.200910204 PMCID: PMC2768804 ArticlePubReaderPDF–3.0 m

S. Khuon, L. Liang, R. Dettman, P. Sporn, R. Wysolmerski, T. Chew. *J Cell Sci.* 2010 February 1; 123(3): 431–440. Published online 2010 January 12. doi: 10.1242/jcs.053793 PMCID: PMC2816188

	SPECIFICATIONS									
Style	Style ID (mm) OD (mm) Glass Ø height Height Access (mm) (inside) (outside) Angle									
FD35	33	35.5	23.5	7.8	9	29°				
FD5040	47.5	49.82	35	7.25	7.4	17°				



	ORDERING INFORMATION
FD35-100	FluoroDish™ Sterile Culture Dish, clear wall, 35 mm,
	23 mm well, box of 100
FD35B-100	FluoroDish™ Sterile Culture Dish, black wall , 35 mm,
	23 mm well, box of 100
FD35PDL-100	FluoroDish™ Sterile Culture Dish, Poly-D-Lysine Coated,
	clear wall, 35 mm, 23 mm well, box of 100
FD3510-100	FluoroDish™ Sterile Culture Dish, clear wall, 35 mm,
	10 mm well, low sidewall, box of 100
FD3510B-100	FluoroDish™ Sterile Culture Dish, black wall , 35 mm,
	10 mm well, low sidewall, box of 100
FD5040-100	FluoroDish™ Sterile Culture Dish, clear wall, 50 mm,
	35 mm well, box of 100
FD5040B-100	FluoroDish™ Sterile Culture Dish, black wall , 50 mm,
	35 mm well, box of 100

Product Index

tem P	age Ite	em Pa	age	Item	Page	Item	Page	Item	Pag
B100-3 42, 47, 104	. 217 35	78 97, 100,	103	13620	97	15920-	G 195	500199	
B100-4 47, 104	. 217 35	79-20		13685	63, 97		195	500200	
B100-6 47, 104		69	9	13740	211	15921-	G 195	500222	
B100F-3 47, 104	217 36	70	97	13776	97	15975.	97	500222-G	
B100F-4 47, 104	. 217 37	18 79	, 88	13822-10	220	15976.	97		
B100F-6 47, 104	. 217 37	42-20	220	13834	48, 108, 110, 111		59		
B120 47, 104	. 217 39	55	9	13835	110, 111	53330-)1 5	500227	202, 2
B120-4 47, 104	. 217 39	60	9		110, 111		170	500227-G	
B120-6 47, 104	. 217 39	93		13854	97, 221	58006.	167	500228	
B120F-3	.217 48	86			63		167		
B120F-4	217 48	98 139, 142, 213,	216	14003	197, 207, 208				
B120F-6	217 50	52	37	14003-G	197		127		
B150-3 47, 104	217 50	53	0		72, 219		127		97, 2
B150-4 47, 104	217 51	10			220		124, 127		97, 2
B150-6 47, 104	217 52	33			220				97, 2
B150F-3	. Z17 JZ	.32			130, 220				,
D150F-5 47, 104	217 53								97, 2
B150F-4	. 217 53	61			220				
B150F-6	. 217 53	62			220		79		
B200-4 47, 104	21/ 53	71 97,			220		79		
B200-6 47, 104		73			220		113, 114		
B200F-4	. 21/ 53	74 97,			220		113, 114		
B200F-647, 104	. 217 53	75 97,	221		220		210 45		
B150F-4	105 53	77 145,	149		220		310 45		
B150F-6	105 53	78 145,	149	14047-10	220	75122-	410 45		
mm Plug	. 97 53	85		14048-20	220	77020	106		207,
-Axis Micropositioner	. 26 53	89		14051-100			173		
-Axis Stereotaxtic Frame		99			220		173		
B120F-4		30-10		14055-2			173		
B120F-6		30-12		14057-10			170, 173		
-Axis Motorized Micromanip		30-15		14058-10			147, 148, 149, 151		
							, , ,		
B120F-4		30-20		14059-2			4		
3120F-6		30-ALL 44			220		4		
lb Weight		35 145,		14073-4			4		
-Axis Stereotaxtic Frame		36 145,		14074	110, 111		179, 21		
B100F-4		40	105	14091	207	94689.	179, 21	500353	
B100F-6		44	101	14109	208	95611.	145	500354	
B120F-4	105 54	47	75	14110	208	95620	145	500355	
-pin DIN Plug 7	9, 88 54	.64	25	14120	205	95699.	145	500366	
Ó5		.68	. 91	14121	205	96824.	37	500366-G	
05PF		69 74, 75, 91			208		11, 17		
11P		70 74, 91			198				
12P		.79			198		37		
15P		82			198		42		
73		83			208				25, 29
00A		 89			207		179, 21		25, 29
							145, 148		23, 29
00AP				14192					
00APP100	115 /1	28 139, 142, 213,			201, 207, 208			500828	
00A-VAC	115 /3	145,			201		59		30,
571 101, 139, 142, 213		26 145,			201				
965		35 139, 142, 213,		14219-G	201		171		
966	107 73	41 139, 141, 213,	215	14224	201	300052	171		
167	107 73	42 139, 141, 213,	215	14224-G	201		75, 91	501224-G	
126-10		57		14225	201	300411	37	501225	
33 74, 7	5, 78 74	65-20	220	14225-G	201	300521	59	501225-G	
34 74, 7		21			195		4		
35 7	,	024 88,			195		35, 37		
05	,	025 88,			208		75		
41		142 55, 59, 113,			97, 221				
47 80, 8		156-100			207		-G		
51 15, 17, 156, 97, 15		157-100			208		181		
		158-100			208				
32									
33 45, 81, 84,		159-100			37				
72		160-100			63, 97		207		
61 79, 8		161-100			81				
59 8		162-100			83		79, 88, 97, 130, 221		
160		163-100			25, 29	500186	181, 191	501263-G	
94 74, 97	. 100 13	316 105, 139, 142, 213,	216	15911	205	500191	71		
16		324 97,		15913			71		
117-10		338 183,			207, 208				
184		347 97, 100,			195		71		
		370			195				
491 97, 100, 115		388 74, 75, 97,	100		195				
192192									
	. 7/ 13	395	17.1	ט-/וענו	195	20019/	71	DU1512	

ltem	Page	Item	Page	Item	Page	Item	Page	Item	Page
501319	207, 208	502043	. 126	503070	223	504444	187	505062	9
	207	502046			223		187		9
	183	502047			223	504459	209	505066	175
	183	502053			223		203		175
	183	502054			223		203		175
	183	502055			223		193		27, 28, 175
	183	502056			223		193 193		28 27, 28
	183 183	502060 502063			223 223		193		27, 28
	183	502063-MGM			223		193		27
	183	502063-SGM			223		193		27, 28
	183	502067			223		193	505094	110, 111
501378	183	502068		503083	223	504524	193		110, 111
	183	502070			223		204		110, 111
501381	183	502102			223		204		110, 111 110, 111
	224	502105			223				110, 111
501608	167, 171	502109 502110			155 190				110, 111
501622	,	502124			190		61, 62		110, 111
501623		502125			183		61, 62	505105	110, 111
501624		502157			157		61, 62		110, 111
	185	502163		503120	52, 53		70		110, 111
	185	502190	. 223		53		183		110, 111
	155	502191			53		183		181 181
501642		502193-4			55		187		110, 111
	155, 157 32	50220155 50221055			201 201				145, 148
501652		502213 123, 124			201		204	600012	145, 148
501653		502224			55		204		145
	155	502225			197		204		145
501656-C	155	502226	. 123	503306	197	504642	204		140, 214
	155	502227			198		204		160
	155	502228			198		204		177
	155	502229			196		204		11, 41, 71, 115
	155	502230			218		204		11, 41, 71, 115
	195 195	502231 502232			218 218		204		179, 21
	202	502233					204	μManager	18
501725-G		502234			113, 187		204	μTip	107, 218
	194	502235			186				A
501778	197	502236	. 124	503536	97, 189	504653	204	A310 Accumulse	r 84
	197	502237 124			157		204		87, 88
	209	502241			32		204		85
	200	E000 10	. 174	503569	32	50/1/13	97		
	209	502242						A365R	85
501053	209	502243	. 124	503570	32	504714	97	A365R A365RC	85
		502243 502244	. 124	503570 503571		504714 504715		A365R A365RC A385R	85 86
501854		502243 502244 502245	. 124 . 124 . 124	503570 503571 503573		504714 504715 504716		A365R A365RC A385R A385RC	
501854		502243 502244	. 124 . 124 . 124 . 124	503570503571503573503598	32	504714 504715 504716 504717	97 97 97 97 97 175	A365R A365RCA385R A385RC	
501854 501855 501856		502244 502245 502249	. 124 . 124 . 124 . 124 . 124	503570503573503598503599503741		504714 504715 504716 504717 504718		A365R	
501854 501855 501856 501857 501858	209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124	503570 503571 503573 503598 503599 503741 503763		504714 504715 504716 504717 504718 504719 504735	97 97 97 175 175 175 204	A365R	
501854	209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123	503570 503571 503573 503598 503599 503741 503763		504714 504715 504716 504717 504718 504719 504735	97 97 97 175 175 175 204	A365R	
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126	503570		504714 504715 504716 504717 504718 504719 504735 504737	97 97 97 175 175 175 204 204	A365R	
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123	503570 503571 503573 503598 503599 503741 503763 503847 503848 504011		504714 504715 504716 504717 504718 504719 504736 504737 504738	97 97 97 175 175 204 204 204	A365R	
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123	503570		504714 504715 504716 504717 504718 504735 504736 504737 504738 504739	97 97 97 175 175 204 204 204 204	A365R	
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119	503570		504714 504715 504716 504717 504718 504735 504736 504737 504738 504739 504740	97 97 97 175 175 175 204 204 204 204 204 204 204	A365R	
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119	503570		504714 504715 504716 504717 504718 504735 504735 504736 504737 504738 504740 504741	97 97 97 175 175 175 204 204 204 204 204 204 208	A365R	85 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 126 . 123 . 119 . 119	503570		504714 504715 504716 504717 504718 504736 504736 504738 504738 504739 504740 504744	97 97 97 175 175 175 204 204 204 204 204 204 204	A365R	85 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119 . 119 . 119	503570		504714 504715 504716 504717 504718 504735 504736 504737 504739 504740 504741 504745	97 97 97 97 175 175 204 204 204 204 204 208 208 193	A365R	85 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119 . 119 . 119 . 119	503570		504714 504715 504716 504717 504718 504735 504736 504737 504739 504740 504741 504745	97 97 97 97 175 175 204 204 204 204 204 208 208 193 193 193 102, 222	A365R	85 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 123 . 119 . 119 . 119 . 119 . 119	503570		504714 504715 504716 504717 504719 504736 504736 504737 504739 504741 504744 504744 504745 504750	97 97 97 97 175 175 204 204 204 204 208 208 3193 193 102, 222 102, 222	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 123 . 119 . 119 . 119 . 119 . 119 . 120 . 120	503570		504714 504715 504716 504717 504719 504736 504736 504737 504739 504741 504744 504745 504750 504756	97 97 97 97 175 175 204 204 204 204 208 208 3193 193 102, 222 102, 222	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 123 . 123 . 123 . 119 . 119 . 119 . 119 . 119 . 119 . 119 . 120 . 120	503570		504714 504715 504716 504717 504718 504735 504735 504737 504740 504741 504744 504745 504751 504766 504766	97 97 97 97 175 175 175 204 204 204 204 208 208 193 193 192 2102, 222 211 211	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 124 123 123 123 119 . 119 . 119 . 119 . 119 . 120 . 120 . 120	503570	32 32 32 127 124, 209 124, 209 124, 209 202 139, 141, 213, 215 176 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504718 504735 504735 504737 504740 504741 504744 504745 504751 504766 504766	97 97 97 97 175 175 175 204 204 204 204 208 308 3193 3193 3102, 222 3102, 222 211 211	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 124 123 123 123 119 . 119 . 119 . 119 . 120 . 120 . 120 . 120	503570		504714 504715 504716 504717 504719 504736 504736 504738 504740 504741 504744 504745 504750 504767 504765 504767	97 97 97 97 175 175 175 204 204 204 204 204 208 208 193 193 102, 222 102, 222 211 211 211	A365R	85 86 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 124 123 123 123 119 . 119 . 119 . 119 . 120 . 120 . 120 . 120 . 120	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504718 504735 504736 504736 504740 504741 504745 504756 504766 504766 504768 504768 504926	97 97 97 97 175 175 175 204 204 204 204 208 308 3193 3193 3102, 222 3102, 222 211 211	A365R	85 86 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 123 126 123 123 129 . 119 . 119 . 120 . 12	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 199 198 199 198 181, 183 59 183 183 191 190 211 187 134	504714 504715 504716 504717 504718 504735 504736 504736 504740 504741 504745 504745 504766 504766 504766 504767	97 97 97 97 97 97 97 175 175 175 204 204 204 204 204 208 193 193 192 208 208 208 208 208 2102, 222 222 2102, 222 211 211 211 211	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 123 126 123 123 129 . 119 . 119 . 120 . 12	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 198 199 198 191 190 1181, 183 183 191 190 191 191 191 191 191 191 191 191	504714 504715 504716 504717 504718 504719 504736 504736 504739 504741 504744 504745 504766 504766 504767 504926 504929 504929	97 97 97 97 97 175 175 175 204 204 204 204 208 208 193 193 102, 222 102, 222 211 211 211 121 121 181, 183	A365R	85 86 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 124 124 124 124 124 124 122 126 123 123 123 119 . 119 . 119 . 120 . 120 . 120 . 121 . 120 . 120 . 122 . 120 . 120	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 191 190 211 187 134 134 134	504714 504715 504716 504717 504718 504739 504736 504739 504741 504741 504745 504750 504766 504766 504767 504768 504926 504927 504930 504930	97 97 97 97 97 175 175 175 204 204 204 204 208 208 208 193 102, 222 102, 222 211 211 211 211 121 181, 183 191 191	A365R	85 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 123 . 123 . 119 . 119 . 119 . 120 . 120	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 198 198 191 197 198 181, 183 183 183 181 181 181 181 183 183 183	504714 504715 504716 504717 504719 504735 504735 504739 504740 504741 504745 504765 504766 504766 504768 504926 504927 504930 504930	97 97 97 97 97 175 175 175 204 204 204 204 208 208 208 193 193 102, 222 102, 222 211 211 211 211 121 121 181, 183 191 191	A365R	85 86 86 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 123 . 123 . 119 . 119 . 119 . 119 . 120 . 223 . 223 . 223	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504719 504735 504735 504737 504740 504741 504745 504751 504766 504766 504927 504929 504930 504930 504930	97 97 97 97 97 175 175 175 204 204 204 204 208 208 208 193 193 102, 222 102, 222 102, 222 102, 222 11 211 211 211 121 121 121 121 121	A365R	85 86 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119 . 119 . 119 . 119 . 119 . 120 . 223 . 223 	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504718 504735 504736 504736 504740 504741 504745 504756 504766 504766 504767 504763 504929 504930 504931 504931 504931	97 97 97 97 97 97 97 97 175 175 175 204 204 204 204 204 208 193 193 102, 222 102, 222 102, 222 102, 222 11 211 211 211 121 181, 183 191 191 191 191 191 59 59 59	A365R	85 86 86 87 87 87 87 209 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119 . 119 . 119 . 120 . 223 . 223 253 253	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504718 504735 504736 504736 504740 504741 504745 504766 504766 504766 504763 504926 504927 504930 504931 504949 504951	97 97 97 97 97 97 97 97 175 175 175 204 204 204 204 204 208 193 193 193 102, 222 102, 222 102, 222 102, 222 112, 211 211 211 211 121 181, 183 191 191 191 191 191 191 191 59 59 59 48, 108 72, 219	A365R	85 86 87 87 87 87 299 ection
501854	209 209 209 209 209 209 209 209 209 209	502243	. 124 . 124 . 124 . 124 . 124 . 124 . 124 . 123 . 126 . 123 . 123 . 119 . 119 . 119 . 119 . 120 . 223 . 223 223 52 . 188	503570	32 32 32 127 124, 209 124, 209 202 139, 141, 213, 215 176 176 53 198 198 198 198 198 198 198 198 198 198	504714 504715 504716 504717 504718 504735 504736 504736 504741 504744 504745 504766 504765 504765 504767 504929 504930 504931 504951 504955	97 97 97 97 97 97 97 97 175 175 175 204 204 204 204 204 208 193 193 102, 222 102, 222 102, 222 102, 222 11 211 211 211 121 181, 183 191 191 191 191 191 59 59 59	A365R	85 86 86 87 87 87 87 209 ection

tem	Page	Item	Page	Item	Page	Item	Page	Item	Pa
.GW10301	102, 222	Bipolar Stimulation	90	CHM1	9	CUV2108-1	173	CX-KIMAKIT-PC	
GW15101	102, 222	Black Epoxy Coated Tweezers .		CHM2	9	CUV2614-1	173	CX-KIMA-P	
GW15301		Black Instruments		CHM3		CUV2674-1		CX-MIRUS	
GW20101		Blood Pressure Monitor		CHM4				CX-MIRUS-MULTI8	
						Cuvettes			
GW2030 1		Blood Pressure, Rodent		CHM5		CUVxxx-1		CX-MIRUS-PRO	
GW40101		Blood Pressure System		CHM6		CX-001		CX-UNIGO-80	
GWxxxx 1	102, 222	Blood Pressure Transducer	130	CHM7	9	CX-002	69	CX-UNIGO-EC-80	
rTherm-SAT		BLPR2	130	CHM8	9	CX-003	69	CX-UNIGO-ECI-80	
RTHERM-SAT-1W		BLPR2 Transducer		Chopstick Electrode Set .		CX-004		CX-UNIGO-ECP-80	
RTHERM-SAT-2W		BNC Cable 15, 17, 9						CX-UNIGO-I-80	
				Circulating Bath		CX-4UECIPAD-4X80FS			
Therm Satellite		BNC Right Angle Adapter 9	37, 221	Circulation Reservoirs		CX-4UECPC-4X80FS	68	CX-UNIGO-P-80	
rTherm Satellite Heater	36	BNC Straight Adapter	221	Clamps	33, 223	CX-4UNCIPAD-4X80FS	68	Cyanoacrylate Adhesive .	141,
Therm-SMT	. 35, 37	BNC T-Connector9	97, 221	Clamp with Rod		CX-4UNCPC-4X80FS	68	_	
RTHERM-SMT-2W		BNC to BNC Cable 23, 75,		Classic System, Anesthes		CX-005		D	
-300		130, 156	73,00,					D2H-DB	
			05	Clip Applier		CX-007		D2H-HB	
1000		BNC-to-Double Banana Adapt		Clips and Clamps	205	CX-008		D2H-HBER	
1000HP		Board Frame Clamp	223	Closed Chamber	160	CX-009	69		
-1000HP2	65	Boom Clamp Stand	181	Closed Loop	37	CX-010	69	D4H	
-2000	65	Boom Stand 18	31 183	CMF20G		CX-018		DAM50	7
4000		Borosilicate 43, 47, 104, 10						DAM80	74
			J, 107,	CMF22G		CX-023		DAM80P	
-6000		217, 218		CMF23G	214	CX-024	68	DAM Series Bioamplifier	rc
6000H	65	Box Filament 11	,	CMF26G	214	CX-026	69		
8000	65	BP1	130	CMF28G		CX-027		Data Acquisition	
ddin Syringe Pump		BPCABLE2		CMF31G		CX-028		DC Amplifier System	
odynia		Brain Matrices						DCAP	
				CMF34G		CX-029		DCAP-L	
M Retractor		BRIDGE8		CMF35G	214	CX-030	69		
plifier	76	Buffer Solutions		CMF90U	214	CX-031	69	DCAP-M	
plifier System	13	Butyl Cyanoacrylate	139	CMOS Camera with Opti		CX-032		Deuterium Halogen Ligh	
algesia Meter 1						CX-032		Deuterium Lamp	
esthesia		C		C-Mount Adapter				Deuterium Replacement	
		C3005 10	12 222	C-Mount CCD Camera Co	oupler. 183	CX-034	69	Diamond Coated Micro	
esthesia Accessories	129			C-Mount Eyepiece Adap	ters 190	CX-035	69		
sthesia Facemask	125	Cables 74,		C-Mount Microscope Att		CX-036	69	Digital Caliper	
esthesia Masks	126	CAL900A	115	21, 179		CX-037		Digital Filtering	
esthesiometer		CALBUF-1	158	,	100			Digital Microforge	
gled Electrode Holder		CALBUF-2		C-Mount to Eyepiece Ad		CX-038		Digital Micrometer Head	
		Calcium Calibration Solutions.		Coaxial Cables, Micro	102, 222	CX-039	69	Digital Microscope Came	
imal Physiology				COLCAM-HD1080P	190	CX-040	69		
imal Temperature Controll	ler 127	Calcium Detection		Collimator		CX-041		Digital Mouse Stereotax	
ti-Oscillation Module	. 13, 17	Calcium Electrode	159			CX-042		Digital Stereotaxic Frame	5
iculated Arm 1		Calcium Neutral Ion Exchange	r. 159	Colored Dissolved Organ	iic ivialler			Dipping Probe	
C2000		Caliper, Digital		168		CX-043		DipTip	
		Camera, Digital		Color Video Camera	190	CX-044	69	DIP-UV-MINI-2	
to-Fluorescence Imaging		•		Compact Magnetic Stand	d 32	CX-045	68		
tomated Tissue Resistance)	Camera, High Definition		Concentric Bipolar Electr		CX-046		DIP-UV-MINI-5	
Measuring System	7	Cannula, Micro	130		0000			DIP-UV-MINI-10	
xiliary Knob Controller		Cannulas	134	90, 92, 93		CX-047		Dispensing Tips	
W201 1		Carbon Epoxy 14	12, 216	Concentric Nerve Cuff		CX-048	69	Disposable Scalpels	
		Carbon Fiber Microelectrodes		Condenser	9	CX-049	69	Dissecting Pad Kit	
Wxxxx1		Carbon Filled Conductive Epox		Connectors	96	CX-050	69	Dissecting Scissors	
s Adjustment Tool	25	the state of the s	,	Constant Load		CX-051			
_		Carbon Monoxide	145					DMF1000	
В		Carbon Wire 10)2, 222	Controller		CX-052		DMF1000-1	
ll Bearing Boom Stand 1	181. 183	Cat/Monkey Adaptor	123	Cosine Corrector	1/1	CX-053	69	DMF1000-2	
Il Electrode		CBL100		Cover Slips, Glass	218	CX-054	69	DMF1000-H5	
				Crile-Wood Needle Hold	ler 202	CX-055	69	DMF1000-M1	
l Joint 25, 29		CBL102 74, 75, 78, 97, 11		Cross Fiber Assemblies		CX-056			
l-Joint Holder Attachment	33	CC-3-UV						DMF1000-M2	
l Joint, Magnetic	224	CCD Camera Coupler	181	CS-mount Adapter		CX-057		DMP	
l Mill		CDOM		Culture Dishes	225	CX-058		Dog/Monkey Adaptor	
e		CDOM-COAST-HV		Current Sources	85	CX-059	69	Double Banana Plug	
ic System, Anesthesia		CDOM-COAST-LV		Cutoff Disk	209	CX-060		Double-Holder Microma	
				CUV1022-10		CX-061			
-10R/LOP		CDOM-FRESH-HV						Dovetail Extension	
-10R/LOP220		CDOM-OCEAN-HV	169	CUV2011-1		CX-062		Dressing Forceps	
-12R	132	CDOM-OCEAN-LV	169	CUV2012-1	173	CX-063		Drill, Micro	
-12R-220		Cell Tester		CUV2023-1	173	CX-067	68	Dri-Ref	
		Cell & Tissue		CUV2025-1		CX-073		DRIREF-2	
tery				CUV2031-1					
tery Charger 85		Center Chamber Gasket				CX-084		DRIREF-2SH	
m Splitter	185	CF10-100	154	CUV2032-1		CX-087		DRIREF-5	
ırcated Fibers		CF10-250	154	CUV2033-1	173	CX-EC	68	DRIREF-5SH	
rcated Light Guide 1		CF10-500		CUV2034-1	173	CX-EXIGO	66	DRIREF-450	
		CF30-50		CUV2051-1		CX-EXIGO-I		DRIREF-L	
ocular Body				CUV2052-1					
ocular Microscope		CF30-100				CX-EXIGO-I-7		Driven Guard Shield 8	
ocular SurgioScope	185	CF30-500	154	CUV2061-1		CX-EXIGO-IM		DRL	
-amplifier		CF30-1000	154	CUV2062-1	173	CX-EXIGO-IM-7	66	Dropchip	
chip		CFN10-50		CUV2063-1		CX-EXIGO-M		Dry Sterilizer	
. '				CUV2065-1		CA EAICO M			
fluorescence		CFN10-100				CX-EXIGO-P	00	Dual Channel Differentia	
fluorometer 14, 18-20, 1	1/8, 179	CFN10-250	154	CUV2066-1		CX-EXIGO-P-7		Electrometer	
Oxy Sensor		CFN30-50	154	CUV2071-1	173	CX-EXIGO-PM		Dual Emission	
psy Punch		CFN30-100		CUV2102-1		CX-EXIGO-PM-7		Dual Manipulator Digita	
				CUV2105-1					
sensing		CFN30-250				CX-KIMA		Frame	
	1/15	CFN30-1000	154	CUV2106-1	1/3	CX-KIMA-I	hh	Dual Manipulator Stered	maxic
osensors oolar Electrodes		Chamber, Measurement		CUV2107-1		CX-KIMAKIT-IPOD		Duai Manipulator Steret	otanic

em	Page	Item	Page	ltem	Page	Item F	age	Item	Pa
ual Microprobe System	80	EZ-830	129	FORT25	79. 89	Hydrogen Sulfide	145	Iris Forceps	195
ual Rate Syringe Pump		EZ-1130		FORT100		Hydrogen Sulfide Macro Senso		Iris Scissors	
ual Tool-Holder Micromar		EZ-1785		FORT250		Hydrogen Sulfide Micro Sensor		Iris SuperCut Scissors	
JMC22		EZ-7000		FORT1000		Hydrogen Sulfide Sensor		IRM23E10	
JMC26		EZ-20027		FORT5000		Hydrostatic Pressure	116	IRM23E15	
mmy Load Resistor Kit	85	EZ-20027G	129	FORT Transducers	89			IRM23E20KT	
mont	194	EZ-20028	129	Frame Clamp	223	l l		IRM23E25	
mont #3C		EZ-20028G		Frame, Stereotaxic		IE010	159	IRM23E30	
mont #4		EZ-20029				IE190			
				FrameWorks				IRM23E30KT	
mont #5		EZ-20030	129	Free Radical Analyzer.	150, 151	IE200		ISDB	
mont Forceps	194, 208	EZ-20030G	129	_		IGS100		ISDB Blank Panels	
nostar		EZ-20032	129	G		II-39	138	ISO2 Adapter	
nostar Forceps		EZ-20032G		Gas Anesthesia	123 124 126	II-78	131	ISO2 Filling Solution	
						II-79			
0773		EZ-20034		Gasket				ISO-80	
773 Electrometer		EZ-20034G		Gaskets for PicoNozzle		II-80		ISO-80P	
Loading	179	EZ Anesthesia	128	GC22	134	II-81	131	ISO-COP-2	
		EZ-B800		GC24	134	II-82	131	Iso-DAM8A	
E		EZ-E25000		GC26	134	II-83	131	ISO-H2S-2	
0	107					II-84			
		EZ-E27000		Glass Capillaries 4				ISO-H2S-100-Cxx	
2	107	EZ-E28000	129	Glass Fiber Filters	170	II-84XL	131	ISO-HPO-2	145
5	107	EZ-FF900		Glass Handling Forcer	os 106	II-336T	138	ISO-HPO-2 Electrolyte	
0		EZ Lever		Glass Micropipettes		II-336TG		ICO UDO 100	1 //
Bars		LT FEACI	∠UJ			II-390		ISO-HPO-100	
		F		Glass Rod				ISO-HPO-100H	
nomy Manual Micromar				Glucose Sensor		II-390G		ISO-HPO-100-L	14
-H5	35	Facemask	125	Glue Gun	105, 142, 216	II-520MR	136	ISO-HPO Startup Kit	
-HC		FC-H5		GN-NET7		II-600MR			
-HOC		FC-HC				II-2200		Isolated Bioamplifier	
				GN-NET25				Isolated Differential Am	
1		FC-HOC		GN-PC7		II-2390		ISO-NO Activator	
F	99, 101	FC-SMT	35	GN-PC25	65	II-2391	135	ISO-NOP1	
	,	FD35-100		GN-TTL		II-2392			
		FD35B-100				II-2393		ISO-NOP007	
				G01-100				ISO-NOP3005	145
		FD35PDL-100	225	G02-100	101	II-2394	135	ISO-NOP3020	14
Cartridges	9	FD223A	82	G03-100	101	II-2395	135	ISO-NOPF	
		FD223AP	87	GO4-100		II-2396	135		
								ISO-NOPF100	
		FD3510-100		Gold Plated Pins, Sock		II-2397		ISO-NOPF100-L10	
Cartridges	9	FD3510B-100	225	Gold Wire	102, 222	II-2400	135	ISO-NOPF200	140
tro 705	76. 83	FD5040-100	225	GR100-4	105	II-2500	136	ISO-NOPF200-L10	
trode Adapter		FD5040B-100		GR100-6		II-2888			
								ISO-NOPF500	
trode Holder 33, 74		Fiber Optic Assemblies		Grip Strength Meter		II-23931		ISO-NOPF500-Cxx	
trolyte, ISO-NO	147	Fiber Optic Cables 15	7, 165	Ground Wire	97	II-B60-025	131	ISO-NOPNM1	145 146
trometer	83	Fiber Optic Collimator	171	GSNO	160	II-B60-038	131	ISONOP Startup Kit	
trophysiology		Fiber Optic Filter Holder		GSNO-50		II-B60-044			
trosurgical Unit		Fiber Optic Holder 16				II-B63-075		ISO-0XY-2	143
				GSNO-100				ISO-TEMP-2	
oy Microelectrodes		Fiber Optic Illuminator 17	9, 181	Guide Cannula	134	II-B64-1		IT-1E	12
ohm	5	Fiber Optic Light Source	177	Guillotine	206	II-MRBP-M	131	IT-14	
0HM-6G	579	Fiber Optic Oxygen Meter. 15	2 155	Guinea Pig Adaptor		II-MRBP-R		IT-18	
0HM-12G		Fiber Optic pH Meter	157	Gairiea i igridaptor		II-PE34.			
OTTAL 24CCNAD				H				IT-18EXLONG	
OHM-24GSNAP		Fiber Optic pH System				ILS005LT		IT-21	
ronmental Control	34	Filament 11	0, 111	Half-Cells	98, 100	ILS010LT	55	IT-23	12
74	1 97 103	Filament Adjustment Assembly	/ 114	Halogen-Halogen Dua		ILS025LT	55	11-ZJ	12
		Filament Cable	,	Track Stand		ILS Gas-Tight Syringes		1	
								J 5	
		Filament for Puller		Halogen Lamp		Impedance Measurement		Jaffe Forceps	
5	103	Filaments, Platinum/Iridium	108	Halogen Replacement		Implantable Probes 127, 133	, 155	Joystick Control	
3	103	Filling Solution		Halsted Mosquito For		INC22		Julabo Circulating Bath .	
2		Filter Holder		HD imaging		INC24		jaiano circulatilig Datili.	
								K	
helial Tissue Voltohmme		Fire-Polished 47, 10		HDMI Cable		INC26			
ky 139, 142,	213, 216	Fixed Flow System		Heated Waterbed		Incapacitance Meter		Katena-Vannas Scissors	
l Cyanoacrylate		Flaming/Brown Pipette Puller		Heater Controller	37	Incremental Hot Cold Plate Ana	lgesia	KG Transducers	
anex Lids		Flexible Dri-Ref	150	Heating Filament		Meter		Kima Pump	
				0				,	
M2		Flexible Implantable Probe		Heating Plate		Induction Chambers		KITE-L	
0	66	Flexible Light Guide	191	Hemostatic Forceps		Infinity Plan Achromat	187	KITE-M3-L	
б	132	Flexible NO Sensor		High Current Isolator.	86	Infusion Pump		KITE-M3-R	
nsion		FLEXREF		High Performance Str		Injection Assembly Parts Kit		Kite Manual Manipulato	
					1 1			KITE-R	,
nsion Cable		Flow Cell, Low Volume		11:-b \/-la ll-k		injector			
Guards		FLOX-PATCH	153	High Voltage Isolator .		Inlet and Outlet Hoses		KITLWCC	
iece Adapter, SLR	190	FLOX-PROBE 15	2, 153	Horizontal Tissue Bath	11	In-Line Fiber Optic Filter Holder	170	KIT-PHYSIO-I	
ieces 181, 182,		Fluorescence Spectroscopy		Hose Assembly		In-line Frame Clamp		KWIK-2	
iece with Reticle		FluoroDish4		Hot Plate Analgesia M		Input Kit		KWIKCAL	
03A	129	Flush Cutters	102	HPO Sensor	147, 148	In Situ Oxygen Monitoring 152	., 153	KWIK-CAST 139, 1	
04A	129	FO-6000	177	HT-1	133	Internal Cannula	134	KWIKGARD	139
07A		Foam Ring		HT-2		Intracellular Amplifier 8		Kwik-Gard Dispenser	
09		FOIMPH		HTS Electrodes		INV-10118		Kwik-Gard Refill	
77	129	Foot Switch45,	55, 81,	Humidifier	35, 37	Inverted Microscope	18	Kwik-Gard Start-up Kit	
78		84, 113, 114, 209	. ,	Humidifier Kit		Inverted Trinocular Microscope		KWIKGLUE	
			100						
79		Forceps		Humidity Probe		IO-KIT		KWIKGUN	
.11	129	Forceps with Tying Platform	196	Hydrogen Electrodes .	159	Ion Selective Electrodes	159	KWIKH-2	
		Force Transducers 15, 22,		Hydrogen Ion Exchang		Iridium Microelectrodes	90	KWIKMIX	
320			,	, u. u gui i i oi i Lauidii)					1
20 30		FORT10g		Hydrogen Peroxide Se		Iris Dissecting Scissors	200	KWIK Mixing Tips	

item P	age ite	em	Page	item	Page	item	Page	item	Page
KWIKPOT-2		agnetic Clamping System	32	Microdialysis	63	MK-BM-8175	38	NanoSensors	. 146, 14
NIK-SIL 139, 140, 213,	214 Ma	gnetic Holding Devices	30, 32	Micro Drill	124, 209	MK-BM-8200	38	Needle Holder	. 196, 20
wik-Tip		gnetic Stand		Micro-Dri-Ref		MK-BM-8225	38	Needle Microprobes	127 13
VIKTPP-2		andrel		Microelectrode Hold		MK-BM-8250		Nerve Cuff Connectors	
						MM1			
		anipulator Clamps		Microelectrode Pulle				Nerve Cuff Electrodes	
L		anipulator, Motorized		Microelectrodes		MM1-3		Neuroscience Cannulas	
	Ma	anipulator Mounting Clam	р 31	MICRO-ePORE		MM1-A		NF26BV-2	
ab Supplies		anometer	116	MicroFil	41, 214	MM1-ALL	26	NF33-36BL	
AB-TRAX-423,	156 Ma	anual Manipulator 21,	25. 179	Microflex Non-Rebre		MM1-C	26	NF33-36BV	
ab-Trax 4/16		anual Micromanipulator		MicroFlow	0	MM3	24, 26	NF33BL-2	
abTrax 8/16 10, 11, 15, 23,		anual Microsyringe Pump .		Microfluidic Biochips		MM3-3		NF33BV-2	
ABTRAX-MDAC 11, 15, 23,	156	, , , ,				MM3-A		NF34BL-2	
ED Illuminated Base	104	anual Stage		Microfluidic Pumps		MM3-ALL		NF34BV-2	
	107	ask for Cunningham Stereo		Microfluidic Syringe I					
ED Illuminated Microscopes		Frame	125	Microforge	. 112, 114, 115	MM3-C		NF35BL-2	
ED Microscope Stand 181,		ask, Gas Anesthesia	126	Microforge System	114	MMJL	24, 29	NF35BV-2	
ED Ringlight 189,	191 MA	X38201	02, 222	MicroImplant	155, 157	MMJR	29	NF36BL-2	
egato	C 1	AX40201		Microinjection		MMP	46	NF36BV-2	
id	170	ayo-Hegar Needle Holder .		Microinjection Syring		MMP-KIT	46	NFGSK-5	
ight Guides				, , ,		Mobile Workstations		NFINHLD	
ight Source	176	32		Microinjection Syster		Mobity Magnetic Clam			
inear Motor 13	0 1/1	04		Microinjection & Tran	•			NFINHLD-G10	
	Y IVIL	04-M3-L		Molds	46	Modular Amplifier Sys		NFQ34-5	
inear Motor Controller 13		04-M3-R	29	MicroLWCC	167	Modular Cables		Nitric Oxide Sensors 145	
inear Stimulus Isolator 87		04R		Micromanipulator	24, 25, 82, 83	Modular DC Amplifier.	79	Nitrite Standard Solution	14
iquid Ion Exchangers		DAC Software 11, 15,		Micrometer		Monopolar Electrodes	90	NL2010MC2T	40, 5
quid Light Guide 21,		echanical Allodynia		Micrometer Head		Mosquito Hemostatic		NOCHM-4	
iquid Trap	115					Motorized Micromania		NOCHM-P	
quid Waveguide Capillary Cell	100	egohm Meter		Micrometer Slide Mi		Motorized Stereotaxic			
	27 IVIL	H1S		Micropipette Calibra				Non-Rotating Micrometer.	
ve Cell Imaging	2/I	H1SF 83,		Micropipette Holder.	45	Mounting Adapter		Non-Rupture	
ve Cell Microscopy		H2R	98, 101	Micropipette Holders	5. 98, 100, 115	Mounting Adapter Plat		Nose Clamp	12
ong Taper Glass Pipettes		H2RF	99. 101	Micropipette Puller		Mouse Adaptor	123	Noyes Scissors	19
ong Working Distance Objective		H2RFW		Micropipettes		Mouse & Neonatal Rat	: Adaptor 123	NSÁ-3	16
113, 114, 183		H2RW				Mouse Dissecting Kit	208		
oop Filament 110,				Micropipette Slide		Mouse Gas Anesthesia		0	
ow Volume Flow Cell	1C7 IVIL	H2S		Micropipette Storage		Holder		02-T0050L250-12012	10
-Shaped Sensors	1.17 IVIE	H2SF		Micropressure Syster	n 115			02-T0050L25012022	
'	O IVIE	H2SFW		Microprobe Thermor	neters 132, 133	MOUSEKIT			
uer Fittings	IVIE	H2SW	99, 101	Micro-Reference Elec	trode 158	Mouse/Rat Chamber		Objective	
uer Stopcock	130 MF	H3F45	99 101	Micro Scissors, Diam		MPH1	100, 101	OC	
uer-to-Tubing Coupler 72,		H3FW45		Microscope		MPH3	100, 101	OC-LFM	12
uer Valve Kit	. 72	H3R	00, 101		100, 102, 104,	MPH4	100, 101	OC-SFM	
uer Valves 219,	220		99, 101	185, 186, 187	Ch	MPH6P	100, 101	Octyl Cyanoacrylate 139	9, 141, 21
.WCC	162 IVIE	H3RF		Microscope Environn		MPH6R		OC-XLFM	
WCC-3050 167,	171 IVIE	H3RFW				MPH6S		Olsen-Hegar Needle Hold	
.WCC-3100 167,	171 IVIE	H3RW		Microscope, Inverted	l 186			Omoga Tip 7	.1 20
WCC-3100 107,	171 ME	H3S	98, 101	Microscope on Boom	Stand 181	MPS-2		Omega-Tip-Z	
WCC-3250 167,	17 I ME	H3SBW	99 101	Microscope Slides		MRI Probe		OmniDrill35	
WCC-3500	167 ME	H3SF	00 101	Microscope Stage		MRI Rat Adaptor	124	Oocyte injection	
WCC-4010 167,	1/1					MRI Stereotaxic Frame	122	Open-Sided Frame Clamp	22
WCC-4050 167,	171	H3SFW	99, 101	Microscope Stage Ad		MT-4	133	Operating Scissors	20
WCC-4100 167,	171 IVIE	H3SW		Microsensor Adapter	Cable	MT-23/3HT		Optical Glass Bottom Dish	
WCC Injection System		H3W45 1	00, 101	147, 148, 151		MT-23/5		Optical measurement	
		H6RF 100, 1	01, 115	Micro Sensors	145				
WCC-M-10		H6RFW 1	00. 101	Microsyringe Pump	46	MT-23/8		Optical Sarcomere Spacing	0
WCC-M-50	107	H6SF 1		MicroTip		MT-26/2		Optical Transducer Amp	
WCC-M-100	107					MT-26/4	133	Optics	16
WCC Start-Up Kit	107	H6SFW 1		MicroTip Oxygen Mic		MT-26/6	133	Optisarc	<i>'</i>
		H7 1		Microvolume Syringe		MT-29/1HT	127 133	OptiSarc Software	
M	ME	H7W 100, 1	01, 104	Millivolt & Megohm N	Meter 88	MT-29/2		Optogenetics	
И1	. 31 ME	H8 1	00, 101	Miniature Microposit	ioners 26	MT-29/3		O-Rings	
11L		H345	99, 101	Miniature Peristaltic	Pump 11, 53,			Oxygen Monitoring Kit	
12		H900R 1		167		MT-29/5		, , ,	
		H900S 1		Mini-Banana Adapte	r 97	MT-D		Oxygen Regulators	
A-3 25, 41, 82	-, 03	etal Electrode Extension				MTM	118	Oxygen Sensing System	1
14C 25, 33	,, _,			Mini-Banana Plug		MTM-3	118	Oxygen Sensors 145	
45 3°	.,	etal Electrodes		Mini DipTip		MTM-6	118	Oxygen Sensors, Optical	1
Лб	. 33 Me	etal Microelectrode Adapte	er 75	Mini Dumont #M5S		MTM-R	118	OXY-MICRO-AOT	
л7	. 32 Me	etal Microelectrodes	90	MiniFlow	155	Multi-Animal Non-Reb		OXY-MINI-AOT	
18		-28G-5	214	MiniFoil	155		0	OX1 WIIN 7.01	
19		-28G67-5		Mini Glue Gun		Unit		Р	
						Multi-Barrel Micropipe			2
И10 30, 179		-34G-5		Mini Micropositioner		Coupling	105	Parallel Frame Clamp	
/10L		-200-1		MiniPhone Patch Cab		Multichannel Perfusio	n System 51	Parallel Rail instruments	
111		200-2		Mini-Phone-to-BNC (Muscle Hyperalgesia		Parallel Rail Stereotaxic In:	st 1.
1325 24	1, 29 MF	E200-H3	114	Mini-Scissors, Ergono	omic 102	Muscle Research Syste		Patch Clamp	1
/13301 21, 24, 41,	, -	- - - - - - - - - - - - - - - - - - -		MINISTAR		mascie nescaren syste	10	Paw Pressure Meter	
лээот 21, 24, 41, Л3301ЕН 25, 74		-200-M1		MiniTip		N		Paw Volume Meter	
				Mini Vacuum Pump			A1 E7		
/3301L 25, 82		-200-M2				NANOFIL		PB150F-4	
/13301-M3-L		200 Microforge		Mixing Tips		NANOFIL-100		PB150F-6	
/13301-M3-R	. 25 Mi	copipette Puller	48	MK-BM-825	38	NanoFil Application Kit	ts 58	Perfusion	39,
42201D 2F 01	2, 83 MI	CRO2T	55	MK-BM-850	38	NanoFil Needle		Perfusion Chambers	
/I3301K 25, 82		cro Bead Sterilizer		MK-BM-8100		Nanoliter 2010		Perfusion System	
	145 1911								
Macro Sensors			30, 221	MK-BM-8125	38	NANOI ITFR2010	59	PFRIPRO-2HS	11 ¹
W3301R	224 Mi	cro Cannula 1 cro Cuff		MK-BM-8125 MK-BM-8150		NANOLITER2010 Nanoliter Injector		PERIPRO-2HS PERIPRO-4HS	

Page

Item

Page

Item

Page

Page

Item

Item

Page

Item

Item	Page	Item	Page	Item	Page	Item	Page	Item	Page
PERIPRO-4LS	11 52 167	Pre-Pulled Pipettes	<i>I</i> 11	Rat Stereotaxic Non-R	ehreathing	RTV Coating	1/12 216	SI-TCM2	16
PERIPRO-8LS		Pressure Calibration Chan		Unit	U	RTV Prime Coat		SI-TCM2B	
Peristaltic Pump		Pressure Manometer		RBMA-200C		RTV Sealant		SI-TM1	,
Peri-Star Pro		Pressure Monitor		RBMA-200S		Rubber Ear Bars		SI-TM2	
		Pressure Pod		RBMA-300C				SI-TM3	
PE Tubing, Flexible		PRO-300 HDS				S		SI-TM4	
PG10150-4				RBMA-300S		Sample Cell			
PG10165-4		PRO-300HDS		RBMA-600C		Sample Injector Atta		SI-TM5	
PG52151-4		Probe 7!		RBMA-600S		Scalpels, Sterile		SI-TM6	
PG52165-4		Probe Extension Cable		RBMS-200C		Scissors		SI-TM7	
pH Implantable Senso		Probe Extension Lead		RBMS-200S		Scotch-Weld		SI-TM8	
pH Meter		Probe Handle		RBMS-300C		SDR2		SI-TM9	
pH Meter, Fiber Optic.		Programmable Filter Mod		RBMS-300S				SI-TM10	
pH Mini Sensors		Programmable Syringe Pu		RBMS-600C		Semi-Micro Concen		SI-TM11	
Phone Plug-to-BNC Ca	able 75	PSMB5N		RBMS-600S	206	Electrodes		Skin Staplers	203
PH-OPTICA-MICRO	157	PSMT5N		RC1		Septum Theta		Slides, Glass	218
PH-OPTICA-MINI	157	PT0110		RC1T	83, 97, 103	Serial Cable		SMA Bulkhead	
Photodiode Array Base	ed	PT0203	. 102, 222	RC2	103	SGE001RN		Small Tissue Probe	179
Spectrophotometers	s 174	PT0402	. 102, 222	RC2F	103	SGE0005RN		SMARTouch	
Photonics Center of Ex		PTM3CC02INS	91	RC3	103	SGE010RNS		SNAP	151
pH System		PTM04RC100	92	RC6		SGE025RN		SNAP25	
PHYSIO-II		PTM041RC225		Rectal Probes		SGE050RN	55	SNAP50	
Physiology Kit I		PTM321RC100		Rectangular Base		SGE050TLL	55	SNAP100	
Physiology Kit II		PTM321RC225		Rectangular Base Post		SGE100RN	55	S-nitrosoglutathione	
		PTP201		Reference Cell		SGE100TLL	55	S-Nitroso-N-acetyl-D,I	
PicoNozzle Kit	41, 45	PTP401				SGE250TLL	55		periiciliamine
PicoPump		PTP406	,	Reference Electrodes.		SGE Syringes		160	07
Piggyback Glass 47,		PTPxxx	,	Reflex Clip		Shielded Electrode		Socket	
Pinger			,	Reflex Clip Applier		Shielded Modular (SO-OXY Startup Kit	
Pinpoint Injection		PTT0110	,	REG2		SI-AOSUB		SP100i	
Pinpoint Injector		PTT0203		REG5K		SI-BAM21-LC		SP100iZ	
Pipetters	70	PTTxxxx		REG10	70	SI-BAM21-LCB		SP101i	63
Pipetter Stand	70	PTxxxx		REG10K	70			SP101iZ	63
Pipette Tips	71	PUL-1000		REG20	70	SI-BF-100		SP120p	63
Planachromatic Object		Puller 43, 48, 108, 109	9, 110, 111	REG50	70	SI-BF-100LLG	21, 179	SP120pZ	63
Plan Phase Contrast		PV820	40, 44, 45	REG100	70	SI-BF-100SMA	21, 179	SP200i	
Plantar/Tail Flick Mete		PV830	40, 43, 44	REG200		SI-BF-LLG-UPGRADI		SP200iZ	
Plantar Test		PZMIII	41, 183	REG1000		SI-BF-SMA-UPGRAD		SP210c	
Plantar Test Analgesia		PZMIII-AAC	183	Regal Pipetters		SI-BMFA		SP210cZ	
Platinum/Iridium Micr		PZMIII-BS	183	REGS5		SI-COLUB	, ,	SP210iw	
Platinum/Iridium Teflo		PZMIII-MI	41, 184	REGS6		SI-CTS200		SP210iwZ	
222	011 WIIC 102,	PZMIV	181	REGS7		SI-CTS200A	12	SP220i	
Platinum/Iridium Wire	102 222	PZMIV-BS	181	REMS-24		SI-CTS200B	12	SP220iZ	
Platinum Wire		PZMTIII	41, 183	REMS-24M		Signal Conditioning	Amplifiers 13	SP230iw	
		PZMTIII-AAC	183	REMS-96		Signal Conditioning	Amplifier	SP230iwZ	
Plethysmometer		PZMTIII-AAC-DIG50	183	REMS-96C		System	13, 14, 16	SP250i	
PM01		PZMTIII-BS		REMS AutoSampler		Signal Generator			
PM01D		PZMTIII-BS-DIG50				SI-HTB2		SP250iZ	
PM01R		PZMTIII-BS-LWD-DIG50		Replacement Electrod		SI-HTB4	11	SP260p	
PM015		PZMTIII-DIG50			25, 29	SI-KG2	15, 22	SP260pZ	
PM015D		PZMTIII-MI		Replacement Electrod	e Sleeve 149	SI-KG2B	15, 22	Spare Plug-Adapter	
PM015R		PZMTIV		Replacement Hoses	3/	SI-KG4	15, 22	Specialty Micromanip	
PM100		PZMTIV-BS		Replacement Lamp	177	SI-KG7		Spectrophotometers	174
PM100D		PZMTIV-BS-DIG50		Replacement Needles		SI-KG7a	15, 22	SPLG100	
PM100R	116			Replacement Sleeve K	lit 147, 148	SI-KG7b		SPLG101	61
PM Series	116	PZMTIV-BS-LWD-DIG50		Replacement Sleeves.	145	SI-KGxx		SPLG110	61
Pneumatic PicoPump	42, 43, 44	PZMTIV-DIG50	101	Replacement Tubing	52	SILFLEX-2		SPLG200	61, 62
Polishing Patch Pipette	es 113	Q		Restrainer	131	SilFlex Tubing		SPLG210	61
PolyFil	105		165 170	RET-2	127, 133	Silicone Adhesive		SPLG212	61
Polymer Core Fiber		QFT1 QUAD		RET-3		Silicone Casting Cor		SPLG270	
Polymer Optical Fiber .		Qualified Optical Fiber		Rhodes Style Tips	92, 93	Silicone Dissecting I		SPLG272	
Portfolio				Ring-Handle Needle H		Silicone Gasket		SPLG-G88PLUS	64
Post		Quantrex		Ring Light		Silicone Pad for ATO		SPL Syringe Pump	
Post Extension		Quartz Tubing		Ring Light Adapter		Silicone RTV		Spring Scissors	
Post Mounting Kit		Quattro		Ring Light Guide				SS316xx	
Post Stand		Quattro 4-in-1 System		Ring Scissors		Silicone Rubber Sea		SS31614	
Potassium Electrodes .		QVLUX -DUV		Ringstand Mounting K		Silicone Tubing		SSM04RC100	
Potassium Ion Exchang		QVLUX-NIR	163	RN001		Silver Epoxy		SSM33A20KT	
Power Frame Enclosur		R		RN0005		Silver Filled Conduc		SSM33A70	
			101 101			Silver Wire			
Precious Metal		R-8-8-WPI01		RN010		SI-MOTDB		SSM041RC225	
Precision Current Sour		Rabbit/Guinea Pig Chamb		RN025		SI-MOT-MT		SSM321RC100	
Precision Stereotaxic F		Rack Mount Kit 45, 75,	/9, 81, 84,	Rodent Brain Matrices		SI-MT-L		SSM321RC225	
Precision Stereo Zoom		88, 115, 130		Rodent Restrainer		SI-MTM-L		SST30407-50	
Microscope		Randall Selitto Paw Pressu		Rod Holder		Single-Barrel Stand		SSTxxxxx	
Precision Stereo Zoom		Meter		Room Temperature Vi		Glass Tubing	42, 47, 104, 217	ST5191	
		Rapid-Punch Biopsy Kit	204	Rotating Base	27, 28	Single Emission		ST5193	
Precision SurgioScope		Rat Adaptor	123	RPE-KIT		Sintering Pellets		Stage Micrometer	
Pre-Polarizer		Rat Dissecting Kit	208	RTD Rectal Temp Prob	e 127	SI-OSARC	18	Stagetop Environment	
Pre-Pulled Glass Pipet	ttes 42	RATKIT	208	RTD Temperature Sen	sor 151	SI-PF100	11, 17	Stainless Steel Rods	
Pre-Pulled Micropipet		Rat Spinal Adaptor		RTV	139	SI-PF100B		Stainless Steel Wire	

tem	Page	Item	Page	Item	Page	Item	Page	Item	Р
tand Alone Optical Force Tra	nsducer	SYS-A365R	85	TIP5TW1LS02	107, 218	TTL Control Module	52, 53		V
Amplifier		SYS-A382		TIP10FLT		Tubing Cartridge		VENDOD	•
		SYS-A385R			,			V503382	
tandard Cuff				TIP10LT		Tubing Kit		Vacuum Pump.	
teel Base Plate		SYS-A385RC		TIP10TW1	,	Tungsten Carbide SuperC		Vacuum Pump I	tegulation Klt
tereo Microscope	184	SYS-BP1	130	TIP10TW1-L	107, 218	Scissors		Vannas Capsulo	tomy Scissors,
tereo Microscope with LED		SYS-DAM50	75	TIP10TW1LS01	107, 218	Tungsten Light Source	177		
Illuminated Base	184	SYS-DAM80	75	TIP10TW1LS02	107. 218	Tungsten Microelectrodes	90	Vannas Scissors	
tereotaxic		SYS-OMEGAZ		TIP10XV119		Tungsten Wire			
		SYS-PM01D				•		Vaporizer Pole I	
ereotaxic Accessories				TIP12FLT		TW100-3 4		Vascular Hooks	
tereotaxic Anesthesia/Evacu	ation	SYS-PM01R		TIP15FLT		TW100-4 4		V-Base	
Mask, Mice or Neonatal Ra	ats 126	SYS-PM015D	115, 116	TIP15LT	107, 218	TW100-6 4	7, 104, 217	V-Clamp	
ereotaxic Anesthesia/Evacu		SYS-PM015R	116	TIP30TW1		TW100F-3 4	7. 104. 217	Ventilator Conn	
		SYS-PM100D		TIP30TW1-L		TW100F-4 47, 10			
Mask, Rats								Versaflex Non-F	
ereotaxic Anesthesia Facem		SYS-PM100R		TIP30TW1LS01		TW100F-6 4		Vertical Surface	Mount
ereotaxic Frame 1	20, 121	SYS-PV820		TIP30TW1LS02	107, 218	TW120-3 4		Vessel Cannulat	ion Forceps
ereotaxic Holder	209	SYS-PV830	42, 45	TIPCA	159	TW120-4 4	7, 104, 217	Vessel Clips	
ereotaxic, Mice and Neona		SYS-REMS	7	TIPH		TW120-6 4	7 104 217		
								VETBOND	
Rats		SYS-TBM4M	00	TIPK	159	TW120F-3 4		Vibration-Free 1	ables
ereo Zoom Microscope 1	80, 182	т		TIPMIX01-05	107, 218	TW120F-4 4	7, 104, 217	Vibration-Free \	Vorkstation
eri-Lite Micro Bead Sterilize				TIPMIX01-05-L	107 218	TW120F-6 4	7, 104, 217	Vibration Isolati	on Platform
erilization Mesh Casette		T-Adapter Kit	147	TIPMIX05-10		TW150-3 4	7, 104, 217		
		Tail Cuff Sensor	131			TW150-4 4		Vibration Platfo	
erilizers		Tail Flick Meter		TIPMIX05-10-L	,		, ,	Video Camera	
-H5	36			TIPTPP	159	TW150-6 4		VIS-NIR	
-HC		Tail Flick Test		TIPTW900A		TW150F-3 4	7, 104, 217	Von Frey Anesth	
-H0C		TAXIC-300		Tissue Bath, Horizonta		TW150F-4 47, 10	4, 105. 217		
		TAXIC-303	119			TW150F-6 4		Von Frey Probe	•••••
mulator		TAXIC-350		Tissue Mounting Hoo					14/
mulus Isolator 86	, 87, 88			Tissue Probe	21	Tygon Tubing	68		W
opcock	220	TAXIC-353		Titanium Forceps				W30S-LED	
		TAXIC-600	120			U		W30ST-LED	
orage Jar		TAXIC-603	120	Titanium Scissors		UBATH-Y	211		
orage Portfolio	208	TAXIC-650		T-Joint Frame Clamp	223	UltraMicroPump		Wall-Mount Pla	
uctural Epoxy	213			TM31A10	91			Waterbed	
X2		TAXIC-653		TM31A20		UltraPath		Waveguide Clea	ning Kit 167
		TBM4M	88, 89			Ultrapath System	165	Wide Field	
X3		TBR1025	151	TM33A05		Ultrasonic Cleaner	211	Wire	
X100		TBR4100		TM33A10	91	Ultrasonic Cleaning System			
X100C	6			TM33A10KT	91			Wire Cutters, Er	gonomic
X100C96		TBS		TM33A20		UMP3		Wire Cutters	102
		TC Mouse Rectal Temp	Probe 127			UMP3-3		Wire Retractor.	
X100F		TC Rat Rectal Temp Pro	obe 127	TM33B01		UMP3-4	55	WP1010	
X100M	6	TEER Measurement		TM33B01KT	91	UMP3T-1			
b-Microliter Injection	56			TM33B05	91	UMP3T-2		WP1212	
J-P9741,		Telemetry Training		TM33B10	91			WP2000	
		Temperature Control N				UMS		WP2020F	
J-P1000 41,		Temperature Probes	37, 133	TM33B20		UniGo	67	WP3030	
J-P2000	109	Temperature Sensor	151	TM33C10	91	Universal Focus Mount	179 181		
per Adhesive 1	41, 215	Test Resistor		TM33CCINS	91	Universal Manipulator Sta		WP3620	
perCut Iris Scissors				TM33CCNON				WP4410	
		TGW0515	102, 222			Universal Pipette Tips	/1	WP5000	
perCut Scissors		TGW1510	102, 222	TM53CCINS		Univette	10	WP5050	
JPER-Dri-Ref	158	TGWxxxx	102, 222	TPP+ (Tetraphenylpho	osphonium)	UPUV	165		
pport Stand	9	Thermocouple Thermo		Electrodes	159	UPVIS		WP5060	
				TPP Electrode				WP5070	
I-QUAD		Thin-Line Heated Water				USBCAM33		WP5070L	
-QUAD-L		Thin Wall	. 47, 104, 217	Transbridge		USBCAM50	189	WPB315AB	
-QUAD-R	28	TIDAS-D2		Transbridge Transduc	er Amplifier 88	USBCAM133	189		
re-Seal Chambers	129	TidasDAQ		Transducer Amplifier .	79, 88. 89	USBCAM152		WPB315AS	
rgical Instruments				Transducer Amplifier				WPB505210	
•	174	TIDAS E Base				USBCAM202		WPSS01	
rgical Instrument Storage		TIDAS-H		Transmission Probes .		USS1L		WPSS02	
Portfolio	208	TIDAS-H2	176	Trinocular	187	USS1S	9		
rgical Kits	207	TIDAS S300		Trinocular Body	183	USS2L		WVLUX-DUV	
rgioScope				Trinocular Microscope				WVLUXHNA	
• .		Tilting Base		,		USS2S		WVLUXNXR	
-TRIO		Tinned Cu Alloy	102	Trinocular Microscope		USS3L		WVLUXPV	
-TRI0235		TIP1TW1	. 107, 108, 42	Illuminated Base	184	USS3S	9	//////////////////////////////////////	
-TRI0245	24	TIP01TW1F		Trinocular SurgioScop		USS4L		WVLUXUVIS	
-TRIO245-L				Trinocular Tube		USS4S			V
		TIP01TW1F-L							X
-TRIO245-R		TIP1TW1-L		Trio Analgesia Meters		USS5L		X-Wide Contact	Nerve Cuffs
riss Tweezers		TIP2TW1	107, 218	TRIO 3-Axis Manipula	tor 27	USS5S			
.G184 139, 141, 2	13, 215	TIP02TW1F	107. 218	Trio 3-in-1 System	135	USS6L			Z
gard 139, 1				Tripod Stand		USS6S		Z-Axis Vertical E	
		TIPO2TW1F-L		_ '					
ringe Multi Rack		TIP2TW1-L		Trough Filament		USS7L		Z-Dimensions	
ringe Pumps 60, 62	, 63, 65	TIP03TW1F	107, 218	Troutman Needle Hol	der 196	USS7S	9	Zebrafish	
ringes with Luers		TIP03TW1F-L		TST33A001KT		USS8L		Zhang, Xueji	
5-705		TIP04TW1F		TST33A05KT		USS8S		Z-LITE	
S-773	81	TIP5TW1	107, 218	TST33A10KT	91	Ussing Current Electrode	9	Z-LITE-186	18
S-900A	115	TIP05TW1F	107. 218	TST33A20KT	91	Ussing Electrode Kit		Z-LITE-Z	
S-A310		TIP05TW1F-L		TST33C05KT				Z-LITE-Z186	
						Ussing System			
	۲۷, ۵۵	TIP5TW1-L	107, 218	TST53A10KT	91	Ussing Voltage Electrode.	9	Z-MOLDS	4
S-A362 S-A365D		TIP5TW1LS01		TST150-6		UV-VIS		Zoom Microsco	

TERMS & CONDITIONS

WPI values its customers and will make every effort to maintain the highest professional and ethical standards. However, WPI regards as invalid any special terms and conditions set forth by any other party, including purchasing, quotations, warranties or other statements. WPI's response to a quotation, purchase order or inquiry shall not be construed as acceptance of any conditions stated therein.

Credit Terms

Net 30 days from date of invoice to customers with satisfactory credit references.

Standard Equipment & Accessories

To insure against misunderstandings or delays, always include model number, name, and significant specifications of systems and accessories when ordering. Also specify any options desired. WPI retains the right to change specifications or discontinue products without prior notice and without responsibility for incorporating such changes into equipment previously sold. Please state the desired operating line voltage and frequency.

Warranty

WPI strives to maintain the highest quality standards in all products. We warranty these products against defects in workmanship or materials, and all WPI products retain a warranty that is stated in the product's manual. (Typically, one year for parts and labor except for consumable items such as electrodes, glass, etc., for which warranty varies from 30 days to six months, or as otherwise specified.) However, WPI hereby specifies that no use or application of any product may be suitable or successful for a specific application and therefore waives all liability in that regard. It is the sole responsibility of the buyer/user to apply the product in a manner consistent with its intended use.

WPI products are not approved for human use unless specifically stated in writing with accompanying FDA (or applicable) documents. WPI assumes no liability or legal, moral, ethical, or fiduciary responsibility for a use of any product in the WPI catalog in treating or studying humans.

Prices

WPI maintains competitive world wide pricing. Prices listed for standard items as described are net and subject to change without notice. Prices (in U.S. dollars) do not include shipping charges. Formal quotations are valid for a period of 30 days unless otherwise specified. Purchases may be charged to your VISA, MasterCard, or American Express account.

Taxes

We reserve the right to add any taxes required by local law or ordinance. Institutions operating under tax-free conditions must state applicable license or contract number and be prepared to furnish necessary exemption certificate.

Packaging and Shipping

Equipment is shipped in adequate commercial packing. All packages should be checked thoroughly before being discarded. Shipping and handling expenses, including adequate insurance coverage, are paid by the customer.

Claims and Returns

Inspect all shipments upon receipt. Missing cartons or obvious damage to cartons should be noted on the delivery receipt before signing. Concealed loss or damage should be reported at once to the carrier and an inspection requested. All claims for shortage or damage must be made within ten (10) days after receipt of shipment. Claims for lost shipments must be made within thirty (30) days of receipt of invoice or other notification of shipment. Pleass save damaged or pilfered cartons until claim is settled. In some instances, photographic documentation may be required. Some items are time-sensitive; WPI assumes no extended warranty or any liability for use beyond the date specified on the container

Do not return any goods to us without obtaining prior approval and instructions from our Returns Department. Goods returned unauthorized or by collect freight may be refused. Goods accepted for restocking will be exchanged or credited to your WPI account and are subject to a restock fee of up to 25%. Goods returned which were ordered by customers in error are subject to a 25% restock fee. Equipment which was built as a special order cannot be returned.

Repair

Contact our Returns Department for assistance in the repair of apparatus. Do not return goods until instructions have been received. Returned items must be securely packed to prevent further damage in transit. The Customer is responsible for paying shipping expenses, including adequate insurance on all items returned for repairs. Identification of the item(s) by model number, name, serial number and proof of purchase (packing slip number or invoice number) as well as complete description of the difficulties experienced should be written on the RMA request form and a copy of the form must be included with the item.

ORDERING INFORMATION

Orders from North America should be directed to WPI's main office:

World Precision Instruments, Inc. 175 Sarasota Center Boulevard Sarasota, Florida 34240-9258

Telephone: (941) 371-1003 (collect calls not accepted). Facsimile copies of purchase orders, inquiries or other correspondence may be transmitted to WPI's Fax number: (941) 377-5428. When sending written purchase orders to confirm a telephone order, please clearly mark "CONFIRMING" on the order to avoid it being shipped a second time. Our Sales Department's e-mail address is: sales@wpiinc.com

North American customers may also order many products online through WPI's Web site — www.wpiinc.com. All customers may use the online Quote Request Form and online Purchase Order Form.

Orders from Belgium, Denmark, Eire, Finland, Luxembourg, Norway, Portugal, Spain, Sweden, Switzerland (FR) and the United Kingdom should be directed to

WPI UK: 1 Hunting Gate
Hitchin, Hertfordshire SG4 0TJ England
Tel: +44 (0)1462 424700 Fax: +44 (0)1462 424701
E-mail: wpiuk@wpi-europe.com

Customers in **France** may call +33 (0)970 44 9000 or E-mail: wpifr@wpi-europe.com

Orders from Austria, Germany, Greece, Italy, Malta, Netherlands, Russian States, Switzerland (DE) and Turkey should be directed to

WPI Germany: Saarstrasse 23 D-61169 Friedberg, Germany

Tel: +49 (0)6031 1602171 Fax: +49 (0)6031 1602180 E-mail: wpide@wpi-europe.com

Orders from China, Hong Kong and Taiwan should be directed to

WPI China: WPI Shanghai Trading Co., Ltd. Rm 25e, No8 Dong Fang Rd., Lu Jia Zui Financial District Shanghai PRC

Tel: + 86 688 85517

Email: chinasales@china.wpiinc.com

Orders from **Brazil** should be directed to

WPI Brazil: Av: Conselheiro Nebias, 756 sala 2611

Santos - Sao Paulo - Brasil CEP: 11045-002

Email: Info@brazil.wpiinc.com

Orders from the countries listed below should be directed to the appropriate WPI distributors. All other international orders should be directed to WPI's main office and factory in Sarasota, Florida.

To expedite shipment of orders outside North America (other than to areas served by the international offices above), payment must be made in advance via wire transfer, check payable in U.S. dollars, or charged to a credit card. Proforma Invoices will gladly be furnished upon request. Please specify proper electrical current for your order.

