



WORLD  
PRECISION  
INSTRUMENTS

# EVOM™ AUTO


Accelerate Your Drug Discovery with Our **New** EVOM™ System  
With Both 24 and 96 HTS Multiwell Plate Capability



# Automated TEER Measurement System

EVOM™ Auto with both 24 and 96 HTS multiwell plate capability is a high throughput screening (HTS) platform offering fast, non-invasive sample scanning by comparing electrical resistance measurements (transepithelial/transendothelial electrical resistance: TEER). TEER measurement experiments are simple to set up and less time consuming than more complex molecular studies. TEER measurement can be used as the primary scanning method to identify any physiological changes that can then be further evaluated by other methods. EVOM™ Auto can capture TEER measurements in 24 and 96 HTS multiwell plates from Corning, Millipore, or MatTek. The sample preparation time in these HTS plates is efficient, allowing for fast, multi-channel pipetting options.

Additionally, the EVOM™ Auto electrode disinfection capability during measurement minimizes sample cross-contamination. Its wireless device control offers the convenience of operating the instrument from a distance, and the small footprint of EVOM™ Auto enables you to use it inside a sterile cell culture hood or an incubator. The EVOM™ Auto provides a fast and efficient platform for early drug discovery, by narrowing down drug targets and drug concentrations through automated, non-invasive sample scanning.



## EVOM™ AUTO EXPEDITES DRUG DEVELOPMENT AND LIFE SCIENCE RESEARCH BY ACCELERATING:

- Hit Discovery Process
- Hit To Lead
- Target Screening, Identification and Validation
- Formulation Optimization & Improving Bioavailability
- Assay Development
- Safety Assessment: Toxicity
- Quality Control of 2-D and 3-D
- In Vitro Tissue Models in Drug Discovery

## APPLICATIONS

- Confluence of Monolayer
- Toxicity
- Blood Brain Barrier (BBB)
- Epithelial and Endothelial Barrier Studies
- Intestinal Drug Absorption: ex. Caco-2, 3-D Tissue Function, Primary Cell
- Lung In Vitro Models for COVID Study, Lung Viral Infection
- Cancer Tissue Studies
- Antibody-Antigen Binding

## SPECIFICATIONS

Autosampler Dimensions (WxDxH)	16x10x8.4"	Electrode Array for 96 HTS Plate	Array of 8 pair of (1.2mm Φ) electrodes
Autosampler Weight	15 lbs.	Electrode Array for 24 HTS Plate	Array of 4 pair of (1.2mm Φ) electrodes
Compatibility	Corning, Millipore, and MatTek 96-well HTS Plates with 96-well electrode array (EVA-EL-03-01) Corning 24-well HTS plates with 24C electrode array (EVA-EL-03-02) Millipore 24-well HTS plates with 24M electrode array (EVA-EL-03-03)	Minimum Sample Reading Time	1 Second
Rinse Stations	3	Control Device to Run Software	Tablet, Laptop, Desktop with Wi-Fi adapter
		Output Data	CSV/Microsoft® Excel
		Resistance Range	10KΩ, 50KΩ, 100KΩ
		CE Certified	Yes

## SYSTEM (CONFIGURATION) OPTIONS

- EVA-MT-03-01 for 96-well HTS (Corning, Millipore, MatTek)
- EVA-MT-03-02 for Corning 24-well HTS
- EVA-MT-03-03 for Millipore 24-well HTS

**NOTE:** The system configurations can be used for any of the listed plate types (24 or 96) by easily switching the electrode arrays and plate adapters. The electrode array heads and plate adapter switch mechanism enables you to analyze samples in different plate types.

**SAVE TIME BY AUTOMATING YOUR PROCESS AND MOVE THROUGH A PLATE QUICKLY**

**FLEXIBILITY TO MANAGE YOUR DATA**

**MINIMIZE HUMAN ERRORS**

**HARDWARE SETUP IS EASY AND REQUIRES NO CONFIGURATION**

**MINIMIZE PROBE DAMAGE AND AVOIDS COSTLY REPAIRS**

**WITH COMPLETE CONTROL OF THE SYSTEM, YOU CAN FINE TUNE THE PROGRAMMING AS DESIRED**

**EASY-TO-NAVIGATE SYSTEM SAVES TIME WHEN CONFIGURING SEQUENCES**

# ELECTRODE ARRAY


- Specially designed electrode array fits precisely in the HTS plates from Corning, Millipore, and MatTek, ensuring consistent placement
- Perform resistance measurements directly in the HTS plates, common or divided, reducing the possibility of contamination and mechanical damage to your cultured cells
- Array of multiple electrode pairs enhances the throughput by saving time to read the samples due to shorter sampling time and shorter time requirement for disinfecting electrodes in between measurements
- Measures a whole column of wells each time the electrode array moves once

## 96 Well Array


- Measures a Column of 8 Wells Sequentially
- Array of 8 pair of (1.2mm  $\Phi$ ) electrodes

## 24 Well Array

- Measures a Column of 4 Wells Sequentially
- Array of 4 pair of (1.2mm  $\Phi$ ) electrodes



EVA-EL-03-01



EVA-EL-03-03

# PART NUMBERS

## EVOM™ Auto System

EVA-MT-03-01	EVOM™ Auto for TEER Measurement in 96 HTS Plate
EVA-MT-03-02	EVOM™ Auto for TEER Measurement in 24C HTS Plate
EVA-MT-03-03	EVOM™ Auto for TEER Measurement in 24M HTS Plate
EVA-EL-03-01	EVOM™ Auto 96 HTS Electrode Array for TEER
EVA-EL-03-02	EVOM™ Auto 24C HTS Electrode Array for TEER
EVA-EL-03-03	EVOM™ Auto 24M HTS Electrode Array for TEER

## Warranty

EVA-MT-03-EX1	EVOM™ Auto 1 Year Premium Warranty
EVA-MT-03-EX2	EVOM™ Auto 2 Year Premium Warranty
EVA-MT-INST	EVOM™ Auto On-Site Premium Installation

## System Includes

EVOM™ Auto TEER Measurement System (EVA-MT-03-xx) includes the autosampler, an electrode array with accessories (for 96 or 24 HTS multiwell plates), an interface unit and its cable, an iPad control tablet with software, and a power cord.



WORLD PRECISION INSTRUMENTS