




Microscope Incubator

Live Cell Microscopy



Incubator on a Nikon Eclipse TE-2000 Microscope. Incubators for Zeiss and Olympus microscopes are also available, as well as versions that accommodate confocal modules. All incubators are compatible with all commercially available cameras, light sources, filter wheels, motorized stages, and motorized nosepieces.

This Live Cell Microscope Incubator was extensively tested in laboratories. When compared with other systems, it offers dramatic advantages. For example, other incubators for live cell microscopy rely on passive, random diffusion of heated air from a single source to maintain the desired temperature setpoint.

With no hot air return vent, the heated air escapes from the system through cracks at the microscope/incubator junction in an uncontrolled, random fashion. These systems offer no temperature uniformity, suffer from focus drift and often experience electrical and vibrational interference from the heater. You will also notice dramatic temperature drifts when the imaging environment is disturbed.

Features

- Unique, diffusion grid, combined with air input and return vents provide an air flow pattern for consistent, even heating, with no hot or cold spots in the chamber
- External heater that can be placed far enough from the system to eliminate electrical and vibrational interference from the heater
- High degree of temperature precision and stability
- Minimal focal drift after equilibrium is achieved— Accuracy $\pm 0.1^{\circ}\text{C}$ at the sample itself, and 0.2°C across the microscope stage (allowing for uniform heating of multiwell dishes)
- Airflow pattern and temperature uniformity eliminate dramatic changes in environmental temperature when the incubator door opens
- Ergonomic design for ease of use— The focus and x/y stage controls are outside of the incubator itself. Large doors allow easy access to the specimen and small ones for cords, tubing, etc.
- Precision, shielded temperature probe
- Simple, one person setup of the system



Microscope Incubator

Live Cell Microscopy

Consistent Air Flow

Air flow affects the temperature uniformity of incubators. The red arrows on Fig. 1 and Fig. 3 indicate air flow. The Live Cell Microscope Incubator uses a diffuser grid and proper venting to insure consistent air flow. Traditional incubators with poor air flow suffer with hot and cold spots in the incubator, as seen in thermal images (Fig. 2 and Fig. 4). Warmer temperatures are indicated by red and cooler temperatures by blue.

Temperature Stability

Thermal stability of the Live Cell Microscope Incubator is clearly shown in Fig. 5 (below). Temperatures were recorded with a digital temperature probe placed at the center of the stage. High precision temperature stability is maintained over both short (left) and long (right) durations.

ORDERING

When ordering a system, please have the following information ready.

Scope _____
 Stage _____
 Stage-Up _____
 Perfect Focus _____
 Camera _____
 Left Port _____
 Right Port _____
 Analyzer _____
 Fluor Attachment _____
 Tif _____
 White Light Tif _____
 Binocular D Head _____
 Tilting Head _____
 Filter Wheels: _____
 Excitation _____
 Emission _____
 Dual Lamphouse _____
 Transmitted Light Shutter _____
 Where: _____
 Front Of Lamphouse/
 Attached To Field Light _____
 Cells: _____
 35mm/60mm _____
 Wells _____
 Coverslips _____





Fig. 5—The temperature remains stable both short- and long-term.




Fig. 1—Single air input and no venting causes random air flow in a traditional incubator.




Fig. 2—Hot and cold spots result from inconsistent flow.




Fig. 3—A diffusion grid with air input and exhaust vents yields consistent air flow.




Fig. 4—Consistent air flow means uniform heating.

World Precision Instruments, Inc.

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA
 Tel: 941-371-1003 • Fax: 941-377-5428 • E-mail: info@wpiinc.com • Internet: http://www.wpiinc.com
Germany: Zossener Str. 55, 10961 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wptide@wpi-europe.com
China & Hong Kong: WPI Shanghai Trading Co., Ltd. • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com
UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com

