



# OxyMini Oxygen Meter

*A new generation of oxygen meters for fiber-optic minisensors*



**Small in size:** The OxyMini is a compact, easy to transport and PC-controlled oxygen meter for small fiber-optic sensors.

**Novel technology:** This technology is superior to conventional intensity based sensors in creating very stable, internal referenced values.

**Stable signals:** This enables a more flexible use of optical oxygen sensors in many different fields of interest.

**Temperature compensation:** The OxyMini compensates for variations in the oxygen content of the sample caused by temperature variations.

**Control:** TTL trigger input and analog output.


## Measuring Principle

The OxyMini measures the **luminescence lifetime** of the immobilized luminophore as the oxygen dependent parameter to avoid problems that are inherent with intensity based measurements. Lifetime based measurements are *not* affected by bending of the fiber or the optical properties of the sample (turbidity, refractive index, coloration).

## Possible Applications

### Biotechnology

- control of cell culture media
- online and non-invasive control of bioreactors



### Implantation

- implantation of oxygen sensors into soil, trees, etc.

### In-Process Control

- bottling plant in breweries
- sewage plants
- water-recycling plants
- quality control of packages



# Fiber-Optic Oxygen Minisensors

*Sensors and housings*

WPI fiber-optic oxygen minisensors are based on 2 mm polymer optical fibers (POF). Depending on the respective application, WPI offers a set of different standard designs. Of course, it is possible to build customer-specific designs.



## Flow-Through Cell with Integrated Oxygen Sensor

- very robust
- excellent long-term stability
- online monitoring
- stands CIP (Cleaning In Place) conditions
- **sterilizable** (autoclave (130°C, 1.5 atm), EtOH, H<sub>2</sub>O<sub>2</sub>)
- measurement is possible in both **liquids and dry gases**




## POF Coated with an Oxygen-Sensitive Foil

- POF with a polished distal tip coated with a planar oxygen-sensitive foil
- usable for process application
- very robust sensor
- sterilizable (H<sub>2</sub>O<sub>2</sub>, EtOH)
- not autoclavable (POF does not stand autoclaving conditions)
- excellent long-term stability

## Planar Oxygen Sensitive Foils

- non-invasive and non-destructive measurement from outside through the wall of a flask or a package
- excellent long-term stability
- online monitoring
- stands CIP (Cleaning In Place) conditions
- sterilizable (autoclave (130°C, 1.5 atm), EtOH, H<sub>2</sub>O<sub>2</sub>)



	Dissolved Oxygen	Gaseous Oxygen
<b>Measuring Range</b>	0 - 45 mg/L 0 - 1013 hPa	0 - 500 % air-saturation 0 - 1013 hPa
<b>Response time (t<sub>90</sub>)</b>	< 40 s (< 60 s with optical isolation)	< 10 s (< 15 s with optical isolation)
<b>Resolution</b> (at 20° C)	60 ± 0.3 hPa; 200 ± 0.9 hPa; 500 ± 3.6 hPa; 2.75 ± 0.01 ppm; 9.00 ± 0.04 ppm; 22.0 ± 0.17 ppm 45.0 ± 0.25 mmHg; 150 ± 0.75 mmHg; 375 ± 2.6 mmHg	
<b>Accuracy</b> (at 20° C)	± 1% air-saturation	
<b>Temperature Range</b>	-10 to 50° C	
<b>Chemical Resistance</b>	Sensors can be used in methanol, ethanol and alcohol-water mixtures. Not useful in organic solvents, such as acetone or chloroform.	
<b>Cross-Sensitivity</b>	<b>No interference</b> to carbon dioxide (CO <sub>2</sub> ), hydrogen sulfide (H <sub>2</sub> S), ammonia (NH <sub>3</sub> ), pH, and any ionic species like sulfide, sulfate or chloride. Only affected by gaseous sulfur dioxide (SO <sub>2</sub> ) and gaseous chlorine (Cl <sub>2</sub> ).	
<b>Calibration</b>	Two-point calibration <b>100% air-saturation</b> (air-saturated water, or water-saturated air) <b>0 % air saturation</b> (deaerated water)	
<b>Drift</b>	< 0.5 % air-saturation (continuous sensor illumination; <b>100000 data points</b> ; 20° C, 100 % air-saturation)	

## World Precision Instruments, Inc.

**USA:** International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA  
**Tel:** 941-371-1003 • **Fax:** 941-377-5428 • **E-mail:** info@wpiinc.com • **Internet:** http://www.wpiinc.com

**UK:** Astonbury Farm Business Centre • Aston, Stevenage, Hertfordshire SG2 7EG England • Tel: 01438-880025 • Fax: 01438-880026 • E-mail: wpiuk@wpi-europe.com  
**Germany:** Liegnitzer Str. 15, D-10999 Berlin, Germany • Tel: 030-6188845 • Fax: 030-6188670 • E-mail: wpide@wpi-europe.com