




## Protocol

---

### VenaT4™ Biochip

Distributed in Europe by World Precision Instruments [www.wpi-europe.com](http://www.wpi-europe.com)  
UK Office +44 (0)1462 424700 [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com)  
German Office +49 (0)6031 67708-0 [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com)


## VenaT4 Biochip, Protocol #1: coating VenaT4 biochips

**Step 1**

The VenaT4 biochip microwells are sealed with a thin film strip. The microchannels of the VenaT4 biochip are coated using a standard yellow tip pipette by dispensing approximately 50  $\mu\text{L}$  of protein (e.g. rhICAM) into each microchannel. Note the excess of liquid on the entrance and exit ports.

**Step 2**

The VenaT4 biochip is then placed in a humidified box and incubated at 4°C for overnight coating

**Step 3**

After the incubation period, turn the biochip upside-down and remove the thin-film strips. Again using a standard yellow tip pipette, add approximately 30  $\mu\text{L}$  of Type I Bovine collagen gel solution with chemoattractant into the wells.

Place the biochip into a humidified box kept in the CO<sub>2</sub> incubator for 15–20 minutes at 37°C. Once gel solidifies, re-seal the microwells with thin-film strips. The biochip is now ready to run the assay.

## VenaT4 Biochip Protocol #2: trans endothelial migration assays under shear flow with VenaT4 biochips (single channel version)



### Step 1:

Suspension cells (e.g. T cells) are re-suspended in culture medium at an appropriate concentration (typically  $5 \times 10^6/\text{mL}$ ) in an Eppendorf tube. Cells are stained with a suitable dye.



### Step 2:

Using the Cellix Mirus Evo nanopump or the ExiGo pump, 30  $\mu\text{L}$  of media is dispensed from pump output cable. Following this, the output cable is inserted into a specified channel on the VenaT4 biochip.



### Step 3:

Then using the Cellix Mirus Evo nanopump, or the ExiGo pump, 40  $\mu\text{L}$  of the media is injected through the channel at a shear stress of 40 dynes/cm<sup>2</sup>. This is done to wash the channel. The waste is aspirated from the microwell of VenaT4 biochip with a pipette.




**Step 4:**

Cell sample is placed into the microwell of this channel on the VenaT4 biochip.



**Step 5:**

Cells are introduced into the channel, by specifying the desired shear stress using VenaFlux Assay software or SmartFlo application. The flow rate will be automatically calculated.



**Step 6:**

Time-lapse fluorescent images are recorded as the microscope objective is positioned over the microwell. The rate of image capture is 6 frames per minute for 30 minutes.

Distributed in Europe by World Precision Instruments [www.wpi-europe.com](http://www.wpi-europe.com)

UK Office +44 (0)1462 424700 [wpiuk@wpi-europe.com](mailto:wpiuk@wpi-europe.com)

German Office +49 (0)6031 67708-0 [wpide@wpi-europe.com](mailto:wpide@wpi-europe.com)