APPLICATION NOTE



What is Microfluidics? Behavior of Fluids Through Microchannels

Microfluidics is both the science which studies the behavior of fluids through micro-channels, and the technology of manufacturing microminiaturized devices containing chambers and tunnels through which fluids flow or are confined. Microfluidics deal with very small volumes of fluids, down to femtoliters (fL) which is a quadrillionth of a liter. Fluids behave very differently on the micrometric scale than they do in everyday life: these unique features are the key for new scientific experiments and innovations.

The key concept related to microfluidics is to integrate in a simple micro-sized system operations that commonly require a whole laboratory.

MICROFLUIDIC CHIPS

A microfluidic chip is a pattern of microchannels, molded or engraved. This network of microchannels incorporated into the microfluidic chip is linked to the macro-environment by several holes of different dimensions hollowed out through the chip. It is through these pathways



that fluids are injected into and evacuated from the microfluidic chip. Fluids are directed, mixed, separated or manipulated to attain multiplexing, automation, and high-throughput systems. The microchannels network design must be precisely elaborated to achieve the desired features (lab-on-a-chip, detection of pathogens, electrophoresis, DNA analysis, etc.). To accurately manage fluids inside the microchannels, specific systems are required. These elements can either be found embedded inside the microfluidic chip, such as Quake valves, or outside of it, like in the case of pressure controllers.

Microfluidic devices exploit the physical and chemical properties of liquids and gases at a microscale. Microfluidic devices offer several benefits over conventionally sized systems. Microfluidics allow the analysis and use of less volume of samples, chemicals, and reagents, reducing the global fees of applications. Many operations



can be executed at the same time thanks to their compact size, shortening the time of experiment. They also offer an excellent data quality and substantial parameter control which allows process automation while preserving the performances. They have the capacity to both process and analyze samples with minor sample handling. The microfluidic chip is elaborated so that the incorporated automation allows the user to generate multi-step reactions requiring a low level of expertise and a lot of

functionalities. The microsystems execute functions that extend from detecting toxins to analyzing DNA sequences or creating inkjet printing devices.

Today, microfluidics provide efficient tools for multiple research areas, and more specifically for biological analysis:

- Whole biological process integrated and simplified for the end-users
- High-throughput, multiplexed and highly parallelled assays
- Faster analyses due to the shorter reactions and/or separation times
- Portable devices for point-of-care applications
- Low reagent consumptions
- Global cost reduction per analysis
- Accurate measurement, microfluidics allowing to increase the measurement resolution in given applications

LAB-ON-A-CHIP

A lab-on-a-chip (LOC) is a device performing on a miniaturized scale one or several analyses commonly carried out in a laboratory. It integrates and automates multiple high-resolution laboratory techniques such as synthesis and analysis of chemicals or fluid testing into a system that fits on a chip. There are many advantages to



operating at this scale. Samples analysis can occur on location, where the samples are generated, rather than being carried to an extensive laboratory facility. Fluids' behavior at this scale makes it easier to control the movement and interaction of samples, causing reactions to be much more potent, and minimizing chemical waste. It also allows reduces exposure to dangerous chemicals.

Lab-on-a-chip devices are a subdivision of micro-electromechanical systems (MEMS) devices and often indicated by

> www.wpiinc.com 04/03/25

"Micro Total Analysis Systems" (µTAS). MEMS are information sensors that relay information to a microcontroller that performs the analysis. These mechanical miniaturized systems are composed of some substrates (silicon, glass etc.). They make use of various types of technology. The one encompassing the whole concept is nanofluidics. Nanofluidics examines the behavior, manipulation, and control of fluids that are limited to structures of nanometric (10⁻⁹ m) dimensions. Nanosensors are a primary component of many lab-on-a-chip systems. Sensors have been elaborated using nanomaterials such as carbon nanotubes, suited to apprehend very low volumes like a single cell sample or even smaller samples. These are highly advantageous since they allow a high degree of analytical flexibility in a lab-on-a-chip system while keeping the small size of the device.

ORGAN-ON-CHIP

Organ on chips are 3D cell culture microdevices aiming to reproduce the key functions of living organs on a computer chip. These microfluidic devices are more efficient than conventional cell culture techniques, since they are able to mimic microenvironments,



as well as their influence on organ function. This allows you to research the human physiology for a specific organ and initiate advancements in artificial disease models.

Organs on chips use microfluidics and microfabrication technologies to better replicate the functionalities of living organs. Among them we can find models like gut on a chip, heart on a chip, liver on a chip, lung on a chip, tumor on a chip, muscle on a chip, multiple organs on a chip, etc.

HISTORY OF MICROFLUIDICS

Microfluidics is an integrative domain originating from a mix of microanalysis, biodefense and microelectronics down to the micrometric and sub-micrometric scale. The history of microfluidics started in the 1950's, with the invention of the first transistor and the apparition of microtechnology. The inkjet printhead technology in



which tiny tubes transport the ink for printing was also invented during this period.

In the 60's computers were miniaturized in order to allow space exploration. The first integrated circuits and microprocessors were created. Technologies like photolithography were developed and allowed to miniaturize and integrate thousands of

transistors on semiconductor wafers. These technologies were then used in pressure sensor production.

In 1979 the first miniaturized gas chromatograph containing mechanical micro-elements



integrated on a silicon wafer was created, thanks to the use of silicon etching procedures developed for the microelectronics industry.

Over the 80s, a new sort of devices called Micro Electro Mechanical Systems (MEMS) emerged, allowing industrial applications such as pressure sensors and printheads. By the end of the 1980s the first microvalves and micro-pumps based on silicon micro-machining were manufactured. During the following



years, several silicon-based analysis systems were developed. At this time microfluidic devices were still made of silicon or glass substrates and required heavy industry facilities and hefty skills.

In the 90s researchers spent a lot of time investigating the applications of MEMS in biology, chemistry, and biomedical fields in order to allow the control of liquids in microchannels. They developed microfluidic fluid handling devices for integrated analysis systems such as micropumps, thermic flow sensors and massic flow sensors and mechanisms for fluid transport, metering, mixing and others.

In the middle 1990s, MEMS technologies implemented tools for genomics research. The military, especially the Defence Advanced Research Projects Agency (DARPA) supported research due to their interest in portable biological and chemical warfare agent detection systems. A leading research area was launched around the elaboration of concepts that could transform the sensor function in a full laboratory analysis on a single microfluidic chip. In the late 90s, the use of soft-lithography allowed the production of cheap microfluidic devices by using polymer molds.

In the early 2000s, technologies based on molding microchannels in polymers like PDMS encountered a great expansion. A large number of laboratories were able to conduct microfluidic research thanks to the reduction of costs and production time of these devices. Nowadays, a lot of micro-pumps, mixers, microvalves and other devices are available and employed in microfluidic applications.

REFERENCES

- 1. Nguyen, N. T., Wereley, S. T., & Shaegh, S.A.M. (2006). Fundamentals and Applications of Microfluidics. *Artech House*.
- Sia, S. K., & Whitesides, G. M. (2003). Microfluidic devices fabricated in Poly(dimethylsiloxane) for biological studies. *Electrophoresis*, 24(21), 3563-3576. https://doi.org/10.1002/elps.200305584
- Weibel, D. B., & Whitesides, G. M. (2006). Applications of microfluidics in chemical biology. *Current Opinion in Chemical Biology*, *10*(6), 584-591. <u>https:// doi.org/10.1016/j.cbpa.2006.10.016</u>
- Beebe, D. J., Mensing, G. A., & Walker, G. M. (2002). Physics and applications of microfluidics in biology. *Annual Review of Biomedical Engineering*, *4*, 261-286. <u>https://doi.org/10.1146/annurev.bioeng.4.112601.125916</u>
- Thorsen, T., Maerkl, S. J., & Quake, S. R. (2002). Microfluidic large-scale integration. *Science*, *298*(5593), 580-584. <u>https://doi.org/10.1126/ science.1076996</u>
- Unger, M. A., Chou, H. P., Thorsen, T., Scherer, A., & Quake, S. R. (2000). Monolithic microfabricated valves and pumps by multilayer soft lithography. *Science*, 288(5463), 113-116. <u>https://doi.org/10.1126/science.288.5463.113</u>
- Duffy, D. C., McDonald, J. C., Schueller, O. J. A., & Whitesides, G. M. (1998). Rapid prototyping of microfluidic systems in poly(dimethylsiloxane). *Analytical Chemistry*, 70(23), 4974-4984. <u>https://doi.org/10.1021/ac980656z</u>

USA: 941-371-1003 UK: +44 1462 424700 Germany: +49 6031 67708-0 China & Hong Kong: +86 688 85517