# CAPSTAR-100 CO<sub>2</sub> Analyzer Fast response end-tidal Carbon Dioxide monitor for small animals



End-tidal CO<sub>2</sub> has long been known to be a reliable measure of respiratory status. The CAPSTAR-100 Carbon Dioxide Analyzer provides continuous measurement of expired CO<sub>2</sub> over the range 0-10 percent. Its rapid response time, small sample flow requirements, and long-term stability make it ideal for respiratory gas measurement of rodents and larger animals.

The heart of the instrument is a miniature infrared analyzer with exceptional stability, low maintenance, and a linear output voltage for recording. The heated measurement cell prevents water condensation. The SENSORPOD, which contains the infrared measurement cell, attaches to the front panel by a 1.5m cable. This cell can be removed for cleaning, inspection, or replacement in minutes. The digital display shows CO2 concentration in percent or mmHg, and allows breath-bybreath peak end-tidal reading or continuous rapid

response. Calibration is easily performed with a single calibrated gas and room air.

The reliability and ease of operation of the CAPSTAR-100 make it the perfect instrument for routine respiratory CO<sub>2</sub> monitoring. It is the natural companion to our SAR-830 series Small Animal Ventilators for monitoring respiratory status. It can also be used to monitor spontaneously breathing animals. Gas sampling rates of less than 30ml/min can be used, making it possible to monitor animals as small as rats. The fast response time of the CAPSTAR-100 allows end-tidal response plateaus with respiratory rates of up to 100 breaths/min.

A range of accessories is available to ensure easy setup and convenient operation. An accessory pack containing sample tubing, moisture traps, and Nafion® moisture-absorbing sample tubes is included with the instrument.

#### Standard Features:

- End-tidal peak or continuous readings
- Low sample flow requirements
- Accurate and stable monitoring
- Linear output signal for recording

#### Applications:

- Respiratory monitoring
- Use with ventilated or unassisted animals
- Suitable for rats and larger animals
- Control of breathing & sleep studies



ACQUIRE • VENTILATE • ANESTHESIA RESPIRATION

## CAPSTAR-100 CO<sub>2</sub> Analyzer

### **Specifications**

Analyzer typesingle beam infrared Sample cell materialssapphire and stainless steel Sample inlet connector		
Ordering Information		
•		
11-10000 CapStar-100 CO <sub>2</sub> Analyzer, complete with accessory kit , 120/240V 50/60Hz 11-01100 Replacement CO <sub>2</sub> sample cell		

11-10000 CapStar-100 CO <sub>2</sub> Analyzer, comple	ete with accessory kit	. 120/240V 50/60Hz
---	------------------------	--------------------

11-01101 Accessory kit; 2 moisture traps, 2 13mm hydrophobic filters,

1 sample line, 2 Nafion®. sample lines (6" and 2')

11-01102 Sample line, plain, 5', with male Luer-loc connectors

11-01103 Sample line, 2' Nafion® tubing (absorbs and evaporates conder

moisture), with Luer connectors

11-01111 Sample line, 6" Nafion® tubing, with Luer connectors

11-01104 Sample tubing, .040" ID Tygon™, 10' coil

11-01105 Luer-loc adapter, .040" - .062" ID flex tubing to male Luer, pack

11-01106 Luer-loc adapter, .125" ID flex tubing to male Luer, pack of 5

In-line miniature (13mm) hydrophobic filter, Luer fittings, pack c. -11-01108

Calibration gas, 4 liter aerosol can, 5% CO<sub>2</sub>, 20% O<sub>2</sub>, balance N<sub>2</sub> 11-01110

11-10001 Serial Port, RS232 Serial data port, with cable

Note 1: This instrument is not to be used for Human Life Support applications or Critical Care monitoring. No claims are made regarding the diagnostic properties of any measurements made using the diagnostic properties of any measurements made using this device.

ACQUIRE • VENTILATE • ANESTHESIA RESPIRATION



SensorPod CUTAWAY VIEW