

BONE WAX

Product

SMI Bone Wax is indicated to control bleeding on bone surfaces. It acts as a mechanical barrier to stop bleeding, exerting a haemostatic function.

Composition

A mixture of whitened bee wax (80%) and isopropyl palmitate (20%). It is slightly yellowish and has a characteristic bee wax odour.

Features

- ✓ Easy to shape and apply "hand-warm"
- ✓ Remains in place and does not swell
- ✓ Ready to use out of package
- ✓ Offers controlled, precise application

Advantages

- ✓ Stops bleeding upon application
- ✓ Not soluble in corporal fluids
- ✓ Biochemically inert and biocompatible



Packaging

Supplied in sterile, individually wrapped sticks of 2.5 g (0.09 oz)
Packs of 12 pcs

Contraindication

Do not use where a rapid regeneration and bone melt is required.

Possible adverse effects

May include foreign body reaction with inflammatory tissue reaction or increased susceptibility to infection.



Product	Unit	Order no
Bone Wax (12 foil packs of 2.5 g)	1	501771-EU

*For quality
&
safety*

