



WORLD
PRECISION
INSTRUMENTS

TABLETOP UV STERILIZER

CLEANRAY FOR QUICK, EASY STERILIZATION USING UV LIGHT

WPI's portable ultraviolet (UV) sterilization cabinet is perfect for quick disinfection of small instruments, masks, cell phones, keys, money and other small items. The UV light automatically turns off when the cabinet door is opened, ensuring safe operation.

BENEFITS

- Rapid sterilization with high power UV light and ozone
- No added reagents required for sterilization
- Large sterilization capacity
- Easy access to the sterilization chamber
- Safety lock turns off the UV light before the cabinet opens
- Safe for users because dark blue lid blocks ultraviolet light
- CE certified
- 1 year warranty

FEATURES

- Professional counter top UV sterilizer
- Perfect for masks, towels and tools
- Automatically sterilizes when the lid closes
- Automatically shuts off when the lid opens



SPECIFICATIONS

Type	UV sterilizer
UVC Wavelength	254 nm
Safety	Interlock
Interior Volume	~ 6.5 L
Product Size	40 × 30 × 14 cm
Light Source	9 W , 3000 hour ave. bulb life
Voltage	110-120 V/220-240 V, 50/60 Hz
Plugs	EU, AU, US, UK
Material	ABS, AS cover, UV bulb
Weight	2 kg (4.4 lbs.)
Certifications	CE



WORLD PRECISION INSTRUMENTS

USA: 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA

Tel: (941) 371-1003 • Fax: (941) 377-5428 • E-mail: wpi@wpiinc.com • Internet: www.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com

Germany: Pfingstweide 16, D-61169 Friedberg (Hessen) • Tel: +49 (0)6031 67708-0 • Fax: +49 (0)6031 67708-80 • E-mail: wpide@wpi-europe.com

China & Hong Kong: Rm 27A, No 8 Donfang Rd., Pudong New District, Shanghai PRC • +86 6888 5517 • 400 688 5517 • ChinaSales@china.wpiinc.com

Brazil: Conselheiro Nabias, 756 sala2611, Santos-São Paulo 11045-002 Brazil • E-mail: info@brazil.wpiinc.com