

MOORFIELD FORCEPS

- Thumb Forceps
- 11 cm (4.33") Long
- 0.9mm Tips
- Black Coated

WPB5056 Moorfield Forceps



WPB5056

FORCEPS

WPB315AS Black Epoxy Coated Cutter, 11 cm (4.3 in.), Patented EZ Lever, strong blades, 6.5 mm from tip to bend

WPB315AB Black Epoxy Coated Cutter, 11 cm (4.3 in.), Patented EZ Lever, 4.5 mm from tip to bend



WPB315AS



WPB315AB



WORLD PRECISION INSTRUMENTS

USA: International Trade Center, 175 Sarasota Center Boulevard, Sarasota FL 34240-9258 USA

Tel: (941) 371-1003 • **Fax:** (941) 377-5428 • **E-mail:** wpi@wpiinc.com • **Internet:** www.wpiinc.com

UK: 1 Hunting Gate, Hitchin, Hertfordshire SG4 0TJ England • Tel: 44 (0)1462 424700 • E-mail: wpiuk@wpi-europe.com

Germany: Saarstraße 23, D-61169 Friedberg (Hessen), Tel: +49 (0)6031 1602171 • Fax: +49 (0)6031 1602180 • E-mail: wpide@wpi-europe.com

China & Hong Kong: Rm 18a, No8 Donfang Rd., Pudong District, Shanghai 200120 PRC • Tel: +86 688 85517 • E-mail: ChinaSales@china.wpiinc.com


Brazil: Conselheiro Nabias, 756 sala2611, Santos-São Paulo 11045-002 Brazil • E-mail: info@brazil.wpiinc.com



WORLD
PRECISION
INSTRUMENTS



NON-REFLECTIVE INSTRUMENTS FOR MICROSCOPY



BLACK COATED

Our black instruments are coated with titanium nitride (TiN), an extremely hard ceramic material. The TiN coating hardens and protects the cutting edge. Ceramic coated anti-reflective instruments are perfect for microscopy and microsurgical applications. Coating surgical instruments with a black ceramic adds a thin layer to the metal instrument, making the instrument harder and giving you greater precision. This anti-glare surface minimizes reflections off the surface of your instruments. The incredibly smooth coating improves the instruments' resistance to corrosion and minimizes friction. The ceramic coating is virtually impenetrable, because the raw material is bonded to the instrument both physically and chemically. These instruments are much more resilient to the pressure of daily use and chemical processing. Coated instruments last considerably longer.

- Corrosion resistant
- Non-reflective
- Low friction
- Wear and tear resistant
- Biocompatible

STARBISMUS SCISSORS

- Round tips
- Ceramic coated

WPB407400 12 cm (4.7 in.), Straight
WPB409000 11 cm (4.5 in.), Curved



WPB409000 **WPB407400**



WPB409000

IRIS SCISSORS

- 11 cm Long
- SuperCut
- Straight
- Ceramic coated

504487 Iris scissors



504487

VANNAS TUBINGEN SCISSORS

- 8 cm Long
- Ceramic coated

504488 Straight
504499 Curved



504488

504499

JAFFE FORCEPS

- 10cm Long
- 0.3mm Tips
- Curved
- Ceramic coated
- Stainless Steel



WPB505210

WPB505210 Jaffe Forceps

MICRO BARRAQUER NEEDLE HOLDER

- 12 cm Long
- 0.45 mm Fine Tips
- Curved
- Ceramic coated
- Stainless Steel


WPB6036 Micro Barraquer Needle Holder



WPB6036

ST MARTIN MICRO CORNEAL FORCEPS

- 11cm Long
- Fine Tips
- Ceramic coated
- Stainless Steel



WPB5035

WPB5035 St Martin Micro Corneal Forceps

WESTCOTT SCISSORS, BLACK COATED

- 11 cm Long, Flat Handle
- Curved
- Ceramic coated
- Stainless Steel

WPB4041 0.14 mm Sharp Tips

WPB4030 0.9 mm Blunt Tips



WPB4041

LIEBERMANN SPECULUM V

- Wire lid tension is adjustable with thumbscrew control
- 12 mm Blades

WPB105675 Liebermann Speculum V



WPB105675