



SUTTER INSTRUMENT MICROMANIPULATORS

Features	MPC-200/MPC-385/MPC-325		MPC-365	MP-285	MP-225
Manipulators per Controller	2		2	1	1
Manipulators per ROE	4		4	1	1
Memory Moves	2 (Home & Work)		2 (Home & Work)	500 & Home	2 (Home & Work)
Travel	25mm on all three axes		25mm on X and Z axes 12.5mm on Y axis	25mm on all three axes	25mm on all three axes
Interface	USB		USB	Serial	None
Expandable	Yes		Yes	No	No
Step Size	62.5nm (w/MP-225 manipulator)	62.5nm (w/MP-285 manipulator)	62.5nm	40nm and 200nm	62.5nm
Selectable Speeds	6		6	14	14
Max. Speed	3mm/sec		5mm/sec	2.9mm/sec	2mm/sec
Type of Bearing	Ball Bearing		Cross Roller	Cross Roller	Ball Bearing
PRICE¹	\$11,900 (MPC-325-2)	\$13,900 (MPC-385-2)	\$14,600 (MPC-365-2)	\$7,900	\$6,500

1. US prices. International prices are 5% higher.

SUTTER INSTRUMENT

One Digital Drive • Novato • CA 94949 • Phone 415.883.0128
Fax 415.883.0572 • Web www.sutter.com • Email info@sutter.com