

THUMB FORCEPS

THUMB FORCEPS

SCALPELS & KNIVES

WOUND CLOSURE

BONE INSTRUMENTS

BONN FORCEPS WITH TYING PLATFORM
504475

- 7.2 cm (2.8 in.) Long
- 0.3 mm Tips
- Stainless Steel



MICRO DISSECTING IRIS FORCEPS
504478

- 10.2 cm (4 in.) Long
- 0.5 mm Tips
- Stainless Steel



MICRO DISSECTING FORCEPS
504479 Angled, Serrated 0.8 mm Tips
504482 Angled, Fine Sharp Tips
504483 Straight Fine Sharp Tips

- 10.2 cm (4 in.) Long
- Stainless Steel



MICRO ADSON TISSUE FORCEPS
501244 1x2 Teeth
501244-G 1x2 Teeth, German
501245 Serrated Tips
501245-G Serrated Tips, German

- 12 cm (4.75 in.) Long
- 0.8 mm Fine point tips
- Stainless Steel



ADSON FORCEPS
500092 1.5 mm Tips
500092-G 1.5 mm Tips, German
500178 Tungsten Carbide Jaws
500178-G Tungsten Carbide Jaws, German

- 12 cm (4.75 in.) Long
- Stainless Steel



TISSUE FORCEPS
503311 4 x 5 Teeth
503312 3 x 4 Teeth
503313 2 x 3 Teeth

- For grasping heavy tissue
- 14.5 cm (5.7 in.) Long
- Stainless Steel



RUBBER-TIPPED FORCEPS
503370 Rubber-Tipped Forceps

- For rodent handling
- 21 cm (8.25 in.)
- Stainless Steel



FOERSTER TISSUE FORCEPS
504476 Straight
504477 Curved

- 10.2 cm (4 in.) Long
- 0.5 mm Tips
- Stainless Steel



IRIS FORCEPS
15916 Straight
15916-G Straight, German
15917 Curved
15917-G Curved, German

- 10 cm (4 in.) Long
- 0.8 mm Tips
- Stainless Steel




GRAEFE FIXATION FORCEPS
501988 Graefe Fixation Forceps

- 11.5 cm (4.5 in.) Long
- 4 mm Tips
- Stainless Steel



ADSON-BROWN FORCEPS
500177 Adson-Brown Forceps
500177-G Adson-Brown Forceps, German

- 12 cm (4.75 in.)
- Tungsten carbide maws
- Stainless Steel



TISSUE FORCEPS
15918 Straight
15918-G Straight, German
501216 Curved
501216-G Curved, German
501742-G Straight, German
501741-G Straight, German

- 12.5 cm (5 in.) Long
- 1 x 2 Teeth
- Stainless Steel



RUSSIAN FORCEPS
501987 Russian Forceps

- 15.25 cm (6 in.) Long
- 7 mm Diameter oval tips
- Serrated
- Stainless Steel




GILLIES DISSECTING FORCEPS
501265 Tungsten Carbide Jaws, Standard

- 15.5 cm (6 in.) Long
- 1 x 2 Teeth
- Stainless Steel



STUDENT TWEEZER, STYLE 1
501974 Style 1
501974-6 Style 1, 6 per pack

- 11 cm (4.3 in.) Long
- 0.45 x 0.55 mm Tips
- Stainless Steel



STUDENT TWEEZER, STYLE 3
501976 Style 3
501976-6 Style 3, 6 per pack

- 11 cm (4.3 in.) Long
- 0.20 x 0.40 mm Tips
- Stainless Steel



STUDENT TWEEZER, STYLE 4
501978 Style 4
501978-6 Style 4, 6 per pack

- 11 cm (4.3 in.) Long
- 0.40 x 0.45 mm Tips
- Stainless Steel



STUDENT TWEEZER, STYLE 2
501975 Style 2
501975-6 Style 2, 6 per pack

- 11 cm (4.3 in.) Long
- 0.40 x 0.55 mm Tips
- Stainless Steel



STUDENT TWEEZER, STYLE 3C
501977 Style 3C
501977-6 Style 3C, 6 per pack

- 11 cm (4.3 in.) Long
- 0.40 x 0.45 mm Tips
- Stainless Steel



STUDENT TWEEZER, STYLE 5
501979 Style 5
501979-6 Style 5, 6 per pack

- 11 cm (4.3 in.) Long
- 0.40 x 0.45 mm Tips
- Stainless Steel



SCALPELS & KNIVES

SCALPEL HANDLE #3
500236 Scalpel Handle #3
500236-G Scalpel Handle #3, German

- 13 cm (5 in.) Long
- Reusable, non-sterile
- Stainless Steel



SCALPEL HANDLE #4
500237 Scalpel Handle #4
500237-G Scalpel Handle #4, German

- 14 cm (5.5 in.) Long
- Reusable, non-sterile
- Stainless Steel



SCALPEL HANDLE #7
500238 Scalpel Handle #7
500238-G Scalpel Handle #7, German

- 16 cm (6.25 in.) Long
- Stainless Steel



RAZOR BLADE HOLDER
14134 Razor Blade Holder
14134-G Razor Blade Holder, German


- 12 cm (4.75 in.) Long
- Round Tips
- Convex/Concave Jaws



FEATHER® BLADES

- Sterile
- 100 per box


504169 #10 Feather Blades for Handles #3 or #7
504170 #11 Feather Blades for Handles #3 or #7
504171 #12 Feather Blades for Handles #3 or #7
504172 #15 Feather Blade for Handles #3 or #7
504173 #20 Feather Blades for Handle #4
504174 #21 Feather Blades for Handle #4
504175 #22 Feather Blades for Handle #4
504176 #23 Feather Blades for Handle #4
504177 #24 Feather Blades for Handle #4



STANDARD BLADES

500239 #10 Standard Blades for Handles #3 or #7
500240 #11 Standard Blades for Handles #3 or #7
500242 #15 Standard Blades for Handles #3 or #7
500243 #20 Standard Blades for Handle #4
500244 #21 Standard Blades for Handle #4
500245 #22 Standard Blades for Handle #4
500246 #23 Standard Blades for Handle #4

- Surgical quality
- Exceptionally strong and uniformly sharp
- A good alternative to Feather
- 100 per Box



SCALPELS & KNIVES

BLADE SAFE SURGICAL BLADE REMOVER
501335 Blade Safe Surgical Blade Remover

- Remove blades easily and safely from any style surgical blade handle.



MICRO DISSECTING KNIFE
14135 Micro Dissecting Knife

- 13 cm (5.1 in.) Long
- 1 mm Blade, 90° Angle



DISPOSABLE SCALPELS

500348 #10
500349 #11
500350 #12
500351 #15
500352 #20
500353 #21
500354 #22
500355 #23

- Carbon steel blades
- Packed 10 per box
- Sterile
- Individually wrapped



DISPOSABLE SAFETY SCALPELS

- Retractable Blade
- Ambidextrous
- 10 per box



501337 #10
501338 #11
501339 #15

WOUND CLOSURE

DERF NEEDLE HOLDER
501246 Derf Needle Holder
501246-G Derf Needle Holder, German

- 12 cm (5 in.) Long
- Straight
- Serrated jaws
- Stainless Steel



BRUCE CLARKE NEEDLE HOLDER
503219 Bruce Clarke Needle Holder

- 12.5 cm (5 in.) Long
- Straight
- Serrated jaws,
- Stainless Steel



HALSEY NEEDLE HOLDER
15926 Halsey Needle Holder
15926-G Halsey Needle Holder, German

- 14 cm (5.5 in.) Long
- Straight
- Serrated
- Tungsten Carbide jaws
- Stainless Steel



MICRO OLSEN-HEGAR NEEDLE HOLDER
501989 Micro Olsen-Hegar Needle Holder

- 12 cm (4.75 in.) Long
- Serrated jaws with suture scissors
- 10 mm Jaws
- 2 mm Tips
- Stainless Steel



RYDER NEEDLE HOLDER
500226 Ryder Needle Holder
500226-G Ryder Needle Holder, German

- 13 cm (5 in.) Long
- Serrated
- Tungsten Carbide Jaws
- Stainless Steel



WEBSTER NEEDLE HOLDER
501843 Smooth Jaw,
14109 Serrated/Extra Delicate
14109-G Serrated/Extra Delicate, German
501929 Straight

- 12.5 cm (5 in.) Long
- Straight
- 1.2 mm Tips with lock
- Stainless Steel



WOUND CLOSURE

HALSEY NEEDLE HOLDER
141110 Halsey Needle Holder
141110-G Halsey Needle Holder, German

- 14 cm (5.5 in.) Long
- Straight
- Serrated jaws
- Stainless Steel



CRILE-WOOD NEEDLE HOLDER
500224 14.5 cm (5.75 in.)
500224-G 14.5 cm (5.75 in.), German
503741 15 cm (5.9 in.), Left-Handed


- Straight
- Serrated
- Tungsten Carbide jaws
- Stainless Steel



OLSEN-HEGAR NEEDLE HOLDER WITH SUTURE SCISSORS

501312 12 cm (4.75 in.)
500227 14 cm (5.5 in.)
500227-G 14 cm (5.5 in.), German
501725 16.5 cm (6.5 in.)
501725-G 16.5 cm (6.5 in.), German

- Straight
- Serrated
- Tungsten Carbide jaws
- Stainless Steel



DIEFFENBACH BULLDOG CLAMPS

14117 Curved
14117-G Curved, German
14116 Straight
14116-G Straight, German

- 6.5 cm (2.5 in.) Long
- Serrated jaws
- Stainless Steel



BONE INSTRUMENTS

BONE CUTTER SCISSORS
503253 Serrated, Straight
503255 Serrated, Curved
503252 Smooth, Straight
503254 Smooth, Curved

- 10.5 cm Long
- Strong Tips
- Stainless Steel



MINI BLUMENTHAL RONGEUR

14090 Mini BLumenthal Rongeur
14090-G Mini BLumenthal Rongeur, German

- 11 cm (4.3 in.)
- Durable, strong handles
- Double spring handles
- Cup shaped tips
- 2 mm Jaw
- Stainless Steel



BONE INSTRUMENTS

LISTON BONE CUTTING FORCEPS
500462 14 cm (5.5 in.), Straight
500358 19 cm (7.5 in.), Straight

- Spring handles
- Stainless Steel



NAIL NIPPERS
500362 Nail Nippers
500362-G Nail Nippers, German

- 14 cm (5.5 in.)
- Concave jaws
- Double spring handles
- Stainless Steel



BLUMENTHAL RONGEUR
14091 Blumenthal Rongeur
14091-G Blumenthal Rongeur, German

- 15 cm (6 in.)
- 3 mm Jaw
- Double spring handles
- Stainless Steel



NAIL NIPPER
500137 Nail Nipper
500137-G Nail Nipper, German

- 11 cm (4.3 in.) Long
- Sharp tips
- Stainless Steel




WIRE CUTTERS, TUNGSTEN CARBIDE
503674 Wire Cutters

- 15 cm (6 in.) Long
- Tungsten carbide tip
- No spring
- Round head
- Stainless Steel



SINGLE BLADE BONE CUTTER
501990 Single Blade Bone Cutter

- 14 cm (5.5 in.) Long
- Double spring handles
- 23 mm (0.9 in.) Jaw
- 1 mm Blade tip
- 2 mm Block tip
- Stainless Steel



MICRO FRIEDMAN RONGEUR
501332 0.8 mm Tip
501333 1.0 mm Tip
14292 1.3 mm Tip
14292-G 1.3 mm Tip, German

- 15 cm (5.9 in.) Long
- Double spring handles
- Stainless Steel



BONE CUTTER
501269 Bone Cutter

- 11 cm (4.3 in.) Long
- Oval Cutting Surface
- 1 cm Cup
- Sharp tips
- Stainless Steel



CORNEOSCLERAL PUNCH
500079 2 mm Wide
500143 1.5 mm Wide

- 8 cm (3 in.)



FRIEDMAN RONGEUR
14089 13 cm (5 in.)
14089-G 13 cm (5 in.), German
503295 15 cm (6 in.)

- Double spring handles
- 2.5 mm Jaw
- Stainless Steel



For more surgical instruments log on and view our full catalog at si.wpiinc.com

SCISSORS

POTTS-SMITH SCISSORS

- 503251** Angled 60°
503250 Angled 45°
 • 18 cm (7 in.) Long
 • Stainless Steel

OPERATING SCISSORS

- 501753** Sharp/Blunt
501753-G Sharp/Blunt, German
501754 Sharp/Sharp
501754-G Sharp/Sharp, German
501743 Blunt/Blunt
501743-G Blunt/Blunt, German
 • Blunt/blunt
 • 11.5 cm (4.5 in.) Long
 • Curved
 • Stainless Steel

OPERATING SCISSORS

- 14192** Sharp/Blunt
14192-G Sharp/Blunt, German
501218 Sharp/Sharp
501218-G Sharp/Sharp, German
501219 Blunt/Blunt
501219-G Blunt/Blunt, German
 • 14 cm (5.5 in.) Long
 • Straight
 • Blunt/blunt
 • Stainless Steel

IRIS SCISSORS

- 501758** Straight
501758-G Straight, German
501759 Curved
501759-G Curved, German
503225 45° Angled, Serrated
503224 35° Angled, Serrated
 • 11.5 cm (4.5 in.) Long
 • Stainless Steel

METZENBAUM SCISSORS

- Rounded tips
 • Stainless Steel

WILMER IRIS SCISSORS

- 503248** Wilmer Iris Scissors
 • 10.5 cm (4 in.) Long
 • 45° Angled blades
 • Stainless Steel

IRIS SCISSORS

- 555646S** Iris Scissors, Titanium
 • 11 cm (4.3 in.) Long
 • Pointed 32 mm tips
 • Straight
 • Titanium

OPERATING SCISSORS

- 501755** Sharp/Sharp
501755-G Sharp/Sharp, German
501756 Sharp/Blunt
501756-G Sharp/Blunt, German
501757 Blunt/Blunt
501757-G Blunt/Blunt, German
 • 11.5 cm (4.5 in.) Long
 • Blunt/blunt
 • Straight
 • Stainless Steel

OPERATING SCISSORS

- 501220** Sharp/Sharp
501220-G Sharp/Sharp, German
501221 Sharp/Blunt
501221-G Sharp/Blunt, German
501222 Blunt/Blunt
501222-G Blunt/Blunt, German
 • Blunt/blunt
 • 14 cm (5.5 in.) Long
 • Curved
 • Stainless Steel

GOLDMAN FOX SCISSORS

- 503256** Goldman Fox Scissors
 • 13 cm (5 in.)
 • Serrated Tips
 • Curved
 • Stainless Steel

- 501747** 11.5 cm (4.5 in.), Straight
501747-G 11.5 cm (4.5 in.), Straight, German
501748 11.5 cm (4.5 in.), Curved
501748-G 11.5 cm (4.5 in.), Curved, German
501252 14.5 cm (5.75 in.), Straight
501252-G 14.5 cm (5.75 in.), Straight, German
501253 14.5 cm (5.75 in.), Curved
501253-G 14.5 cm (5.75 in.), Curved, German
501254 17.5 cm (7 in.), Straight
501254-G 17.5 cm (7 in.), Straight, German
501255 17.5 cm (7 in.), Curved
501255-G 17.5 cm (7 in.), Curved, German

ANGLED SCISSORS

- 503225** Angled Scissors
 • 11 cm (4.3 in.) Long
 • 45° Angled to side
 • Sharp/sharp
 • Serrated
 • Stainless Steel

SCISSORS

MAYO HARRINGTON SCISSORS

- 555659S** Mayo Harrington Scissors
 • 14.5 cm (5.75 in.) Long
 • 32 mm Tips
 • Straight
 • Titanium

IRIS SUPERCUT SCISSORS

- 500046** Iris Supercut Scissors
 • 10 cm (3.9 in.) Long
 • 30° Angled blades
 • Stainless Steel

TENOTOMY SUPERCUT SCISSORS

- 14395** 10 cm (4 in.), Straight
14395-G 11.5 cm (4.5 in.), Straight, German
14396 10 cm (4 in.), Curved
14396-G 10 cm (4 in.), Curved, German
 • Stainless Steel

KILNER SCISSORS

- 503307** Kilner Scissors
 • 14 cm (5.5 in.) Long
 • Flat tips
 • Curved
 • Stainless Steel

IRIS SUPERCUT SCISSORS

- 14218** 10 cm (4 in.), Straight
14218-G 10 cm (4 in.), Straight, German
14219 10 cm (4 in.), Curved
14219-G 10 cm (4 in.), Curved, German
503259 11.5 cm (4.5 in.), Straight
503260 11.5 cm (4.5 in.), Curved
14225 12.5 cm (5 in.), Straight
14225-G 12.5 cm (5 in.), Straight, German
14224 12.5 cm (5 in.), Curved
14224-G 12.5 cm (5 in.), Curved, Germansnip
 • Stainless Steel

METZENBAUM SUPERCUT SCISSORS

- 503269** 11.5 cm (4.5 in.), Straight
503270 11.5 cm (4.5 in.), Curved
14214 12 cm (5 in.), Straight
14214-G 12 cm (5 in.), Straight, German
14215 12 cm (5 in.), Curved
14215-G 12 cm (5 in.), Curved, German
503271 16 cm (6.2 in.), Straight
503272 16 cm (6.2 in.), Curved
503273 20 cm (7.75 in.), Straight
503274 20 cm (7.75 in.), Curved
503275 25 cm (9.75 in.), Straight
 • Not recommended for cutting sutures or drains
 • Stainless Steel

ANGLED SCISSORS

- 500047** Standard
500047-G German
 • 14 cm (5.5 in.) Long
 • Sharp/sharp
 • 30° Angled to side
 • Stainless Steel

MAYO SCISSORS

- 501749** Straight
501749-G Straight, German
501750 Curved
501750-G Curved, German
 • 14 cm (5.5 in.) Long
 • Stainless Steel

MAYO SCISSORS

- 501751** Straight
501751-G Straight, German
501752 Curved
501752-G Curved, German
 • 17 cm (6.75 in.) Long
 • Stainless Steel

IRIS SCISSORS, TUNGSTEN CARBIDE

- 500216** Straight
500216-G Straight, German
500217 Curved
500217-G Curved, German
 • 11 cm (4.3 in.) Long
 • Stainless Steel

SCISSORS

MAYO SUPERCUT SCISSORS

- 14216** 14 cm (5.5 in.), Straight
14216-G 14 cm (5.5 in.), Straight, German
14217 14 cm (5.5 in.), Curved
14217-G 14 cm (5.5 in.), Curved, German
503277 16 cm (6.3 in.), Straight
503278 16 cm (6.3 in.), Curved
503279 20 cm (7.8 in.), Straight
503280 20 cm (7.8 in.), Curved
503281 25 cm (9.8 in.), Straight
503282 25 cm (9.8 in.), Curved
503384 17 cm (6.75 in.), Straight
 • Beveled blade, rounded tips to avoid punctures
 • Smooth and clean cuts with minimal tissue tearing
 • Stainless Steel

METZENBAUM SCISSORS, TUNGSTEN CARBIDE

- 501739** 11.5 cm (4.5 in.), Straight
501739-G 11.5 cm (4.5 in.), Straight, German
501740 11.5 cm (4.5 in.), Curved
501740-G 11.5 cm (4.5 in.), Curved, German
501733 14.5 cm (5.75 in.), Straight
501733-G 14.5 cm (5.75 in.), Straight, German
501734 14.5 cm (5.75 in.), Curved
501734-G 14.5 cm (5.75 in.), Curved, German
501744 17.75 cm (7 in.), Straight
501744-G 17.75 cm (7 in.), Straight, German
501745 17.75 cm (7 in.), Curved
501745-G 17.75 cm (7 in.), Curved, German

- Tungsten carbide inserts help retain blade sharpness longer
 • Commonly used for dissecting or cutting heavy tissue
 • Beveled blade, rounded tips to avoid punctures
 • Stainless Steel

IRIS SUPERCUT SCISSORS, TUNGSTEN CARBIDE

- 501263** Straight
501263-G Straight, German
501264 Curved
501264-G Curved, German
 • Gold blade has tungsten carbide insert — black blade is serrated
 • Produces clean cut without tearing the tissue
 • Used for cutting delicate or heavy tissue
 • 11 cm (4.3 in.) Long
 • Sharp tips
 • Tungsten Carbide blades
 • Stainless Steel

METZENBAUM SCISSORS, SUPERCUT

- 501261** 15 cm (5.9 in.), Curved
501261-G 15 cm (5.9 in.), Curved, German
503359-G 15 cm (5.9 in.), Straight, German
501262 18 cm (7 in.), Curved
501262-G 18 cm (7 in.), Curved, German
 • Beveled blade, rounded tips to avoid punctures
 • Tungsten Carbide blades
 • Stainless Steel

MAYO SCISSORS, TUNGSTEN CARBIDE

- 500218** 14 cm (5.5 in.), Straight
500218-G 14 cm (5.5 in.), Straight, German
500219 14 cm (5.5 in.), Curved
500219-G 14 cm (5.5 in.), Curved, German
500220 17.5 cm (7 in.), Straight
500220-G 17.5 cm (7 in.), Straight, German
500221 17.5 cm (7 in.), Curved
500221-G 17.5 cm (7 in.), Curved, German
 • Tungsten carbide inserts help retain blade sharpness longer
 • Used to cut heavy tissue or suture
 • Beveled blade, rounded tips to avoid punctures
 • Stainless Steel

SCISSORS

MAYO SUPERCUT SCISSORS

- 501259** Straight
501259-G Straight, German
501260 Curved
501260-G Curved, German
 • Gold blade has tungsten carbide insert — black blade is serrated
 • Used for cutting heavy tissue
 • Produces clean cut with minimal tearing
 • 17 cm (6.75 in.) Long
 • 5.5 mm Tips
 • Tungsten Carbide blades
 • Stainless Steel

SPENCER STITCH SCISSORS

- 500380** 9 cm (3.5 in.)
500380-G 9 cm (3.5 in.), German
501703 10.75 cm (4.25 in.)
501703-G 10.75 cm (4.25 in.), German
501704 13.25 cm (5.5 in.)
501704-G 13.25 cm (5.5 in.), German
 • Used for suture removal
 • A hook-shaped blade is used to grab the suture and lift it for cutting
 • Stainless Steel

ROGER WIRE CUTTING SCISSORS

- 501321** Roger Wire Cutting Scissors
501321-G Roger Wire Cutting Scissors, German
 • One serrated blade and two notched blades
 • For cutting fine wires
 • 12 cm (4.75 in.) Long
 • Stainless Steel

METZENBAUM SCISSORS

- WP4099** Metzenbaum Scissors
 • 14 cm (5.5 in.) Long
 • Curved
 • Stainless steel

LISTER BANDAGE SCISSORS

- 500444** 15 cm (6 in.)
500444-G 15 cm (6 in.), German
501701 8.75 cm (3.5 in.)
501702 10.75 cm (4.25 in.)
501702-G 10.75 cm (4.25 in.), German
 • Angled blade with smooth tips
 • Used for cutting dressings and bandages
 • Stainless Steel

CHECK OUT
OUR SURGICAL
INSTRUMENTS
SITE:
SI.WPIINC.COM

SURGICAL INSTRUMENTS

QUICK SELECTION GUIDE

| SCISSORS | THUMB FORCEPS | SCALPELS & KNIVES
 | WOUND CLOSURE | BONE INSTRUMENTS



(866) 606-1974

SI.WPIINC.COM